2026-27 CATALOGUE

- General Laboratory Chemicals & Solvents
- Bioreagents, Biochemicals & Speciality Fine Chemicals
- Animal Cell Culture
- Nanomaterials and Renewable Energy
- Dehydrated Culture Media Products
- BioLit® DNA & Protein Tools and Kits
- Laboratory Accessories
- Bio-Indenta®

endless process of science...


Quality Policy

We, at Sisco Research Laboratories Pvt. Ltd., shall always strive to provide the best quality Laboratory Reagents and Biochemicals, Microbiological Culture Media, Supplements & Kits, and BioLit™ range of protein and DNA tools required for teaching, research, analysis and manufacture of high quality products for better quality of life.

We shall always maintain and upgrade our products and services to meet the highest International Quality Standards.

It is our constant endeavor to achieve maximum customer satisfaction by implementing an efficient customer relationship management system to meet the requirements and expectations of our valued customers.

To achieve our goals, we have instituted a quality management system conforming to ISO 9001:2015, ISO 17025:2017, CE Certification for our Microbiology range of products and an FDA Approval for select range of Pharmacopoeial products; and are firmly committed to continually improving its effectiveness through technological advancement, enhancement of skills, cultivation of a suitable work environment and dedicated team work



SRL

An ISO 9001:2015 Company

From the desk of the

Chairman

Dear Friends and Partners in Science

As we present this new edition of the SRL catalogue, I reflect with immense gratitude on our shared journey over the past five decades. The year 2025 marked a significant milestone—SRL's 50th anniversary—a celebration not merely of our company's legacy, but of the enduring trust and partnership of scientists, researchers, educators, and institutions across the world who have shaped our progress.

Guided by your evolving research needs, we have significantly expanded and strengthened our portfolio this year. I am pleased to introduce our advanced range of Cell Culture Media, thoughtfully developed to support optimal cell growth and reproducibility; our expanded portfolio of Antibiotic Susceptibility Test discs for more accurate and dependable susceptibility testing; a broader selection of specialized biochemicals, organics, and inorganics designed to simplify and enhance experimental workflows; We also present the long awaited range of high-purity LC-MS Grade Solvents, engineered to deliver exceptional sensitivity, precision, and consistency for advanced analytical and research applications.



Our commitment to the diagnostic and life sciences industry remains stronger than ever. Through our PrimeDx range, we continue our endeavour to provide high-quality solutions that empower laboratories, healthcare professionals, and researchers with reliability and confidence.

As SRL continues to expand its global footprint, we remain deeply anchored in the spirit of "Make in India" — proudly building world-class scientific solutions in India for the global scientific community. We believe true innovation begins with strong foundations at home and meaningful impact across the world.

This catalogue represents far more than a compilation of products. It is a reflection of our unwavering commitment to supporting your discoveries, your ambitions, and the countless scientific breakthroughs that shape a better future for humanity.

On behalf of the entire SRL family, I extend my heartfelt appreciation for your continued trust and partnership.

Together, let us continue shaping the future of science — one discovery at a time.

With warm regards and sincere gratitude,

S. K. Agarwal

S. K. Agarwal
Chairman
15th May 2026

Quality Control







Quality control facilities in our manufacturing plants

- Anaerobic Culture Jars
- Analytical Balances
- Antibiotic Zone Reader
- Atomic Absorption Spectrophotometer (AAS)
- Autoclaves

- Bacteriological Incubator
- B.O.D. Incubator

- Calorimeter
- Conductivity Meters
- Colony Counter
- Cyclomixer

- Electrophoresis Equipments

- FTIR Spectrophotometer

- Gas Chromatograph (w/Head space)
- Gel-Strength Indicator

- High Speed Centrifuges
- HPLC Instrument
- Hot Air Ovens

- IR Spectrophotometer

- Karl-Fischer Apparatus
- Kjedhal's Apparatus

- Lab Furnace
- Laminar Air Flow Bench

- Melting Point Apparatus
- Micro-Ovens
- Microscopes

- pH Meters
- Polarograph
- Potentiometric Titrator

- Refractive Index Apparatus
- Rotary Flask Shaker

- TLC Equipments
- Transilluminator
- Turbidometer

- UV Cabinet
- UV Spectrophotometers

- Viscometer

- other support instruments

It is the precise instrumental, chemical, biochemical and microbiological analysis which enables **SRL** to guarantee that its products shall always meet the highest international standards.

Contents

15 **Part A**
General Laboratory Chemicals & Solvents

77 **Part B**
Bioreagents, Biochemicals &
Speciality Fine Chemicals

153 **Part B1**
Animal Cell Culture

159 **Part C** Nanomaterials and Renewable Energy
(Carbon Nanotubes, Nanodispersants and
Nanodispersions, Nanowires, Nanofilms and Nanopowders)

167 **Part D** Dehydrated Culture Media
• Product Listing • Media Supplements
• ChroMed™ • Media Kits


211 **Part E** BioLit®
• DNA & Protein Tools

233 **Part L**
Laboratory Accessories

239 Bio-Indenta™
Global searching station for Unlisted Compounds

International Business

You will find the **SRL** brand in 6 different regions around the world and in more than 100 countries. Our main motive has always been to provide our products to customers in different corners of the world, wherever there is an industry or an institution that is working on improving the quality of our lives through research, innovation, quality control and manufacturing. Having more than 40 authorized agents and distributors and more than 100 trade partners, **SRL** is always on the lookout to make brighter and long lasting partnerships with companies that have a similar goal in mind, i.e. to be the best in their business. Before we appoint trade partners, we make it a point to conduct systematic checks to ensure that our partners are in a position to offer better sales support and post-sales services to their valued customers who choose an **SRL** product. We welcome new partnerships for growth and prosperity.



Our global reach

Europe & Eurasia

- Austria
- Azerbaijan
- Belgium
- Bulgaria
- Cyprus
- Czech Republic
- Denmark
- Finland
- France
- Germany
- Greece
- Hungary
- Ireland
- Italy
- Malta
- Moldova
- Netherlands
- Norway
- Poland
- Portugal
- Russia
- Spain
- Sweden
- Switzerland
- Turkey
- Turkmenistan
- United Kingdom
- Uzbekistan

Africa

- Benin
- Botswana
- Chad
- Cameroon
- Djibouti
- Ethiopia
- Ghana
- Kenya
- Madagascar
- Malawi
- Mali
- Mauritius
- Mozambique
- Namibia
- Nigeria
- Rwanda
- Seychelles
- South Africa
- Tanzania
- Uganda
- Zambia
- Zimbabwe

Middle East

- Bahrain
- Egypt
- Jordan
- Kuwait
- Lebanon
- Morocco
- Oman
- Qatar
- Saudi Arabia
- Tunisia
- United Arab Emirates

Central Asia

- Bangladesh
- Bhutan
- Kazakhstan
- Maldives
- Nepal
- Srilanka
- Tajikstan

North & South America

- Argentina
- Brazil
- Canada
- Chile
- Costa Rica
- Ecuador
- Guatemala
- Jamaica
- Mexico
- Peru
- Trinidad & Tobago
- Venezuela
- United States of America

East Asia & Pacific

- Australia
- Brunei
- Burma
- Cambodia
- China
- Fiji
- Indonesia
- Japan
- Laos
- Malaysia
- Mongolia
- New Zealand
- Philippines
- Singapore
- South Korea
- Thailand
- Vietnam

overseas customers can reach us at
export@srlchem.com

Depots

Hyderabad

Plot No.08, H. No. 5-9-234/1,
Rajiv Gandhinagar,
Prashanti Nagar, Kukatpally,
Hyderabad 500 072.
Telangana, India.
Telefax: +91-40-7981540907
e-mail: hyderabad@srlchem.com

Chennai

Plot No. 8B,
Self Help Industrial Estate,
Keeelkattalai, Kovilambakkam,
Chennai 600 117,
Tamil Nadu, India.
Telephone: +91-9962005587
Email: chennai@srlchem.com

Delhi

69/6A, (Part A)
Rama Road Industrial Area,
Najafgarh Road,
New Delhi 110 015, India.
Telephone: +91-11-4591 4530
Email: delhi@srlchem.com

Mumbai

Raj Rajeshwari Logistic Park,
Gala No. 3-8, Bldg.No. B-4,
Behind Mini Punjab Grill,
Vehale, Bhatala,
Bhiwandi 421 302,
Maharashtra, India.
Telephone: +91-93229 31625 /
9619147231
Email: bhiwandi@srlchem.com

Kolkata

46, N.S.C. Bose Road,
P.O. Malancha Mahinagar,
South 24 Paraganas,
Before Malancha Bazaar,
Kolkata 700 145,
West Bengal, India.
Telephone: +91-8334800445
Email: kolkata@srlchem.com




Administrative Office:
608-B, Satellite Gazebo,
Near Solitaire Corporate Park,
Andheri Ghatkopar Link Road,
Chakala, Andheri (East), Mumbai 400 099. India
Tel.: +91-22-4268 5800
E-mail: marketing@srlchem.com

Website: www.srlchem.com

Label Explanation



for Laboratory & Research Chemicals



12052 – 500 Gms

Potassium Hydroxide Pellets
ExiPlus, Multi-Compendial, 85%
meets compendial specs
of USP ,BP, & Ph.Eur

Specifications
Assay (HCl titration) ..... min. 85%
pH (10% aq. solution)..... min.10.5
Content of Potassium (K) min.59.9%

Maximum Limits of Impurities
Chloride (Cl) 0.005%
Sulphate (SO4) 0.003%
Iron (Fe) 0.001%
Phosphate (Po4)..... 0.0005%
Heavy Metals (Pb)..... 0.001%
Carbonate (K2CO3)..... 2.0%
Aluminium (Al) 0.00002%
Sodium (Na) 1.0%

Sisco Research Laboratories Pvt. Ltd.
608-B, Satellite Gazebo, Chakala, Andheri (E), Mumbai - 400099, Maharashtra, India. +91-22-42865800. Email: customercare@srlchem.com

Mfg. Date:
Exp./Re-test Date:
Batch No.:
Not for Retail Sale.
Industrial/Institutional Pack.


Opened on: _____
Use by:_____

For laboratory, industrial, academic, and research use only. No liability is accepted for accidents, injuries or damages from handling or use. Dispose containers and contents according to local, regional or international regulations.

1310-58-3
KOH
M.W. 56.11

www.srlchem.com

for Dehydrated Culture Media, Ingredients, Supplements & Kits



REF **63971 – 500 Gms**

Nutrient Agar

For the cultivation and maintenance of a wide variety of microorganisms.

Composition	gm/lit.
Agar	15.00
Peptone	5.00
Sodium chloride	5.00
Yeast extract	2.00
Beef extract (ex. Buffalo)	1.00

Directions :
- Add 28.00 gm powder to 1.0 litre purified water and mix thoroughly.
- Gently heat and bring to boiling.
- Autoclave at 15 psi pressure at 121°C for 15 minutes.
- Cool to 45-50°C.
- If desired, aseptically add sterile 5-10% defibrinated blood.

Sisco Research Laboratories Pvt. Ltd.
608-B, Satellite Gazebo, Chakala, Andheri (E), Mumbai - 400099, Maharashtra, India. +91-22-42865800. Email: customercare@srlchem.com

Mfg. Date:
Best Before:
LOT **Batch No.:**
Not for Retail Sale.
Industrial/Institutional Pack.

Opened on: _____
Use by:_____

IVD CE


8°C 25°C

Final pH at 25°C : 7.4 ± 0.2
Formula standardized/ adjusted to meet performance criteria.

For laboratory, industrial, academic, and research use only. No liability is accepted for accidents, injuries or damages from handling or use. Dispose containers and contents according to local, regional or international regulations.

www.srlchem.com

for BioLit™ DNA & Protein Tools



92972 – 50 µg **BioLit™**

ProxiO 100bp DNA Ladder


Mfg. Date:
Exp./Re-test Date:
Batch No.:
Not for Retail Sale.
Industrial/Institutional Pack.

Opened on: _____
Use by:_____

Sisco Research Laboratories Pvt. Ltd.
608-B, Satellite Gazebo, Chakala, Andheri (E), Mumbai - 400099, MH, India. +91-22-42865800. Email: customercare@srlchem.com

For laboratory, industrial, academic, and research use only. No liability is accepted for accidents, injuries or damages from handling or use. Dispose containers and contents according to local, regional or international regulations.


www.srlchem.com



Website

Consistently we follow in all aspects of our products and services continues even more comprehensively on our website. Some of the most revolutionary information we offer,

1. **Interactive Product Pages**
   - Technical product specifications & print-ready Product Data sheets
     - (containing technical & application for each product)
   - LIVE Stock Checks
     - Check the stocks of your desired product packing in real time at each of our
       - Manufacturing Sites
       - Distribution Depots
   - Avail prices in more than 17 Global Currencies
   - Access to H.S. Codes & MSDS of all products
2. **Quick Links for viewing and downloads**
   - Product search by keywords, CAS, application, grade
   - Certificate of Analysis (COA)
   - Material Safety Data Sheet (MSDS)
   - Latest Catalog Download
   - Technical Product Brochures
   - Picture gallery for product packaging
3. **Online Transactions**
   - Get quotations
   - Place Orders
   - Ask us a Technical Query
   - Apply for a position with us
4. **Locate your nearest SRL Stockists & Distributor**
   - Indian Stockist network provided City-wise & State-wise
   - Global Distributors provided country-wise
5. **Multi-Lingual Display**
   - Read content in more than 80 Global languages of your preference
6. **News and Information**
   - Follow our developments and road shows
   - Path-breaking research developments in the Bioscience industry
   - Latest SRL product offerings and highlights
7. **Regulatory Certifications, now available online.**
   Click and Download our
   - ISO 9001-2015 Quality Management Certificate
   - FDA (Food & Drug Association) Certificate
   - CE (European Conformity) Certificate
   - GMP (Good Manufacturing Practices) Certificate



www.srlchem.com

is ideal for researchers, scientists and customers who are always looking for new molecules, product information and availabilities at the click of a button!

Trademark & Grades of SRL

No part of this book may be reproduced by any mechanical, photographic or electronic process or in the form of a phonographic audio/audio visual recording, nor may it be stored in a retrieval system, transmitted or otherwise copied for public or private use, without written permission from the publisher.

Trademarks owned by Sisco Research Laboratories

Registered Trademarks

SRL® Bio-Indenta® BioLit® ChroMed® ExiPlus® FluroBronze®
PrimeDx™ SafeDye™ Serabeads® Seralite® Seralose®

"All trademarks mentioned in the catalogue other than those mentioned here, are not property/owned by Sisco Research Laboratories Pvt. Ltd. (SRL)"

Grades of Purity

SRL products adhere to highest grades of quality. However, we understand that our multiple grades may tend to overlap and disorient users. The explanations below give an approximate representation of the **SRL** grades and most closely suited applications of our products.

AAS	Products for Atomic Absorption Spectrocopy specially suited for AAS instrumentation.
ACS/AR/Reagent Grade	Complying with the specifications of American Chemical Society (ACS) reference manuals/ANALAR specifications, for critical laboratory analysis and quality production applications.
Dried/Specially Dried	Containing highly reduced water content for 'anhydrous' reactions, synthesis and assay studies.
Cell Culture	Low endotoxin products specially suited for cell culture DNA, RNA, Genetics applications.
Electronic Grade	Very low ionic impurities and suited for electronics and high-end applications.
ExiPlus™	Complying with the compendial specifications of USP (United States Pharmacopoeia) or BP (British Pharmacopoeia) or Ph.Eur (European Pharmacopoeia) or all. For production, quality control and laboratory applications.
Extrapure	Purified grade suited for general laboratory or synthesis.
GC-HS	For Gas Chromatography assay. Suited for use in GC instrumentation with Head Space.
High Purity/Ultrapure	Having very high levels of assay/purity (typically 3N (99.9%), 4N (99.99%), 5N (99.999%) or more) as compared to similar products in their class/category. Specially synthesized and purified for critical applications with low elemental impurities.
HPCL/UV	For High Pressure Liquid Chromatography or Ultra Violet Spectroscopy. Suited for use in HPLC or UV instrumentation.
Microbiology/Bacteriology	For microbiological application and bacteriological studies.
Microscopy	Typically dyes, stains and indicators (powders and solutions) used for microscopy and histological applications.
Molecular Biology	For molecular biology, cell culture, tissue culture, genetics, and all such applications, tested 'non-detected' for DNAse, RNAse and Protease.
NMR Spectroscopy	Solvents for NMR studies suited to most advanced NMR instrumentation.
OAS	Organic Analytical Standard for reference studies.
PrimeDx™	Refers to high-quality reagents optimized for precision and reliability in diagnostic applications.
Pure/Technical/Practical/Grade	General, commercial grade for synthesis and routine applications.
OAS/CHR/Reference Standard Grade	Organic Analytical Standard/Chromatographically Harmonized Reagent to be used as in reference studies as an analytically pure standard.
Scintillation	Ideal for scintillation studies.
Tissue Culture	Ideal for plant tissue culture application & studied.

For Pharmaceutical Grade Labelled Products

IP	Indian Pharmacopoeia.
BP	British Pharmacopoeia.
USP-NF	United States Pharmacopoeia & The National Formulary.

For Dehydrated Culture Media

Harmonized	Dehydrated media prepared as mentioned in ICH Guidelines.
Bacteriology	For microbiological application and bacteriological studies.
BactoBio/Bacto	For microbiological preparation and cell culture media.
(I/P)	Dehydrated media prepared as mentioned in IP (Indian Pharmacopoeia).
(U/P)	Dehydrated media prepared as mentioned in USP (United States Pharmacopoeia).
(B/S)	Dehydrated media prepared as mentioned in BIS (Bureau Of Indian Standards).
(I)	Dehydrated media prepared as mentioned in ISO (International Standards Organization) Guidelines.

Production Capabilities







We invite you to tour and audit our plant and discuss any queries with our technical team

Our Bulk Manufacturing facilities follow cutting edge bio-techniques for producing many of our specialized products in bulk volumes. Multi-purpose process vessels and unit operations ensure consistent high quality batch-after-batch in a timely manner.

Our production facilities include –

- SS & Glass-lined Reactors with distillation columns
- All Glass Reactors
- SS Centrifuges
- SS Nutsche Filters
- SS Sparkler Filters
- SS Pulverizers, Sifters and Blenders
- Vacuum Tray Dryers
- SS Rotocone Vacuum Dryers, etc.

To support these unit operations, our utilities include –

- Chilled Brine Plant
- Cooling Water and Chilled Water Plant
- Boiler and Oil Heating System
- High Pressure Vacuum System
- Micron-filter Air Purification System
- DM Water Plant, Nitrogen Plant, Scrubbers, etc.

Contact info@srlchem.com for more info.

Keeping strong the SRL quality promise.

Bulk Manufacturing & Custom Synthesis (BMCS) Business

Developing bulk specialty chemicals and customized new molecules in compliance with your in-house specifications supported by complete tech-packs.

Today, **SRL** is not only known for its high quality laboratory chemicals, but it is also been recognized as an important supplier of various fine chemicals as raw materials. Our company is manufacturing reagents and powders in bulk scale to meet the critical requirements of several industries that we serve. As a supplier of key technologies, we want to grow along with our customers in the pharmaceutical and biopharmaceutical industries. This expansion of capacity and modernization is our response to the growing global demand for and high requirements placed on quality and product safety for final pharmaceutical and biopharmaceutical formulations. SRL has an extensive product range and supports its customers in all phases of the drug manufacturing process, from research, quality control and commercial production, through to providing service in concrete safety and regulatory issues.

SRL's Bulk Manufacturing and Custom Synthesis (BMCS) business has broadened its production at our Taloja MIDC site to include existing and new bulk fine chemicals which are being manufactured at commercial scales. Production of the raw materials have been on since 2011 at this site and complements the existing offering of highly regulated active and non-active pharmaceutical ingredients manufactured here. Most of the products manufactured at the Taloja MIDC bulk plant facility meet and exceed IP, BP, USP & Ph. Eur regulation; follow Current Good Manufacturing Practices (cGMPs) for biopharmaceutical production; and comply with International Consortium for Innovation and Quality in Pharmaceutical Development.



Target Audience:

- Active Pharmaceutical Ingredients (API) (Human & Veterinary)
- Antibodies
- Biosimilars
- Biocides, Insecticides & Crop Sciences
- Bio-Fermentation (Upstream & Downstream)
- Blood Plasma
- Diagnostic Kit
- Excipients & Intermediates (Organic & Inorganic)
- Formulations (New & Generic)
- Home Care & Sanitation
- Hormones
- Injectables
- Neutraceuticals
- Peptides
- Personal Care Products
- Protein Purification
- Vaccines (Human & Veterinary)

Reaction Capabilities:

- Acid Chloride Preparations
- Aldol & Crossed Aldol Reaction
- Assymetric Reduction
- Bromination
- Caroxylation
- Cryogenic Reaction
- Cyclization
- Cyclo Addition (4 + 2, 4 + 3) Reaction
- Dehydration
- Esterification
- Formylation
- Friedel – Craft
- Grignard Reaction
- Heck Reaction
- Hydrochloride preparation
- Iodination
- Ion exchange Reactions
- Jones Oxidation
- Lithiation Chemistry
- Mignonac Reaction
- Mitosonobu Esterification
- Nitration
- Nitrosation
- Oximination
- Pachmann Condensation
- Perkin Reaction
- Sandmeyer Reaction
- Sonogashira Coupling
- Sulfonation
- Suzuki Coupling
- Swern Oxidation
- Williamson's Ether Synthesis
- Writing-Homer Reaction

We offer:

The BMCS division at the MIDC Taloja facility include difficult-to-source materials available from research to commercial-scale volumes. The program takes advantage of **SRL**'s purification, re-crystallization and distillation competencies to manufacture a short list of catalog products in the commercial product segment.

You will find our ExiPlus grade of products that are complying to the compendia specifications of Indian Pharmacopoeia (IP), British Pharmacopoeia (BP), United States Pharmacopoeia (USP), National Formulatory (NF) and European Pharmacopoeia (Ph. Eur). We are well versed in providing these products for production requirements. Alternately, our highly researched ACS and AR grades are anyway much better and stringent in their impurity limits and are also being used as raw materials in various industries. We also are highly skilled in providing additional value through quality documentation packages, readily-available inventory and multiple batches to support customer qualification activities

Tech-packs:

SRL not only provides the raw materials through its BMCS program, but we also cater to the strict documentation requirements which our customers ask. Since the last decade, India has become a hub for the regulated as well as the non-regulated pharmaceutical markets, and to serve them we offer complete vendor questionnaire documents with all the technical, commercial and ethical information that we can provide.

Post Sales Support:

We at **SRL** take special care to value our customers with utmost seriousness post-sales. Customers requiring additional services of more document support, pre-shipment sampling for repeat orders, production planning for scale-up and addressing non-conformity issues are some of the areas which we address. In times of product rejections due to storage or any other issues, we provide a complete CAPA (Corrective And Preventive Action) analysis document and such documents are recorded for future product improvements.

Today we provide a helping hand to the world's top 100 pharmaceutical companies to file successful ANDA (Abbreviated New Drug Application) filings in both regulated and non-regulated markets. We are ready to partner with you, so do get in touch with us on **marketing@srlchem.com**

* Document support:

- Aflatoxin free statement
- Allergen free certificate
- Brief description of manufacturing locations, capabilities
- Certificate of analysis (COA) / Method of analysis (MOA)
- Company Organogram
- Dioxan free certificate
- GC, HPLC, NMR graphs
- Gluten free certificate
- Kosher & Halal certificate
- Latex free statement
- List of equipment and instrumentation
- Material safety data sheet (MSDS)
- Melamine free certificate
- Metal catalyst statement
- Microbial contamination certificate
- Quality management certificate (QMS)
- Residual solvent & impurity profiling
- Route of synthesis (ROS) / Process flow chart
- Stability data
- TSE/BSE Declaration
- Vendor Questionnaire (VQ)



* Document support provided depends on individual product and its chemistry. Not all documents are available for all BMCS products.

ExiPlus Range

(Products confirming to the reagent specifications of USP, BP & Ph. EUR)

Exiplus range of products has been carefully designed to meet the requirements of scientists working in the manufacture of drugs, active pharmaceutical ingredients, bio-pharmaceutical & neutraceutical products in semi-bulk or bulk scale. ExiPlus products are the safest and most reliable products to meet the needs of students, researchers working in quality control, research and development laboratories in premier institutions all over the world.

Offering more than 300 ExiPlus World-class reagents and biochemicals, that are manufactured in compliance with the reagent specifications of various pharmacopeia such as United States Pharmacopeia (USP), British Pharmacopeia (BP), European Pharmacopeia (Ph.Eur).

Take the ExiPlus advantage, why loose a good chance:

- ExiPlus for use as an excipient
- ExiPlus for use as an active pharmaceutical ingredient
- ExiPlus for use in biopharmaceutical production
- ExiPlus for use in neutraceutical production
- ExiPlus for use in critical research and development

The quality of our raw materials and our extensive documentation support enable you to do your work quickly and cost-effectively. With ExiPlus SRL provides its customers with a high level of safety in QC, R&D and production activities. We provide you with all necessary technical documents such as Product Data Sheets, CoA, MoA and other documentary support can be made available upon request.

ExiPlus will simplify your life so let's make life simple.



ACS Grade Reagents

(Products confirming to the reagent specifications set by the American Chemical Society Committee)

Our ACS grade of reagents and chemicals are manufactured and tested to meet the specifications set by the American Chemical Society (ACS) Committee on Analytical Reagents. The specifications listed in this catalog are based on the ACS Reagent Chemicals Tenth Edition.

SRL offers more than 200 reagents and chemicals complying with ACS grade reagent specifications of very high quality for use in laboratory for both quality control in academics and industries. They are superior grade reagents that are contaminant-free that are being preferred world-wide because of the stringent specifications set by the American Chemical Society (ACS) committee. Many customers of SRL want to purchase an ACS grade product instead of an AR grade product, however both grades of products are considered to be superior for routine analytical testing applications.

Some common Application of ACS grade products:

- ACS reagents for industrial cleaning purposes
- ACS reagents for routine quality control
- ACS reagents for chemical process development
- ACS reagents for commercial scale-up of manufacture
- ACS reagents for analytical applications wherein UV absorption is not critical
- ACS reagents for the metal & steel industry

We provide you with all necessary technical documents such as Product Data Sheets, CoA, MoA and other documentary support can be made available upon request. We provide many products in various pack sizes and also at affordable prices so that students and scientists can benefit from the ultrapure purity of the ACS grade reagent manufactured, tested and packed for the most stress-free and reliable experience.

Come and try our ACS grade, because it's worth it!

New Introductions

Part A - General Laboratory Chemicals & Solvents

Alkali Blue 6B
Aspartame ExiPlus
Benzyl Alcohol GC-HS
Biebrich Scarlet (C.I. 26905)
Bromine Number Titration Solvent
Calcium Stearate
Chromium (III) Picolinate
Citric Acid Anhydrous
Crotonaldehyde
Crystal Violet (Gentian Violet) PrimeDx
Cupric Chloride Dihydrate ACS
DL-Cysteine Hydrochloride Monohydrate
Diethanolamine ACS
Diethylene Glycol Dimethyl Ether (DEGDME) Dried
Dimethyl Sulphate
N,N-Dimethylurea (1,3-Dimethylurea) (DMU) PrimeDx
Disodium Tartrate Dihydrate ACS
1-Dodecanol (Dodecyl Alcohol, Lauryl Alcohol) ACS
Ethyl Chloroformate
Ferric Tartrate
Formic Acid for HPLC
D-Fructose ExiPlus, Multi-Compendial
Gallic Acid ACS
Iodine Monochloride ACS
Lead Subacetate ACS
Lead (II) Hydroxide Acetate
Mesityl Oxide (MOX)
Methanol for LC-MS
Methyl Methacrylate (Stabilized)
2-Methyl Tetrahydrofuran Dried
Natamycin (NAT)
Nonylphenol Ethoxylate (NPE)
Orthophosphoric Acid ExiPlus, Multi-Compendial
Paraformaldehyde (PFA)
Petroleum Ether 40-60 ExiPlus, Multi-Compendial
DL-Glutamic Acid
p-Phenylenediamine
2-Phenylethanol (PEA)
Phytic Acid (InsP6)
Platinum (II) Chloride
Platinum (IV) Chloride
Platinum (IV) Oxide Anhydrous
Potassium Bicarbonate ACS
Pyridine ACS
Pyridium p-Toluenesulfonate (PPTS)
Sodium Antimonate
Sodium Diacetate
Sodium Perchlorate Monohydrate
Sodium Pyroantimonate
Sodium Tellurite
Sorbitan Sesquioleate (SPAN 83)
Strontium Hydroxide Octahydrate
Sucralose
Tetrabutylammonium Acetate (TBAA)
Tetrabutylammonium Chloride (TBAC)
Tetraethylammonium Chloride (TEAC)
Tetramethylammonium Hydrogen Sulphate for HPLC
Tetramethylammonium Hydrogen Sulphate Monohydrate
Trimethyl Phosphate (TMeP)
Triton X-100 ExiPlus, Multi-Compendial
Triton X-114
Vitamin A (Retinol)
Water for LC-MS

Part B - Bioreagents, Biochemicals & Specialty Fine Chemicals

2-Acrylamido-2-Methyl-1-Propanesulfonic Acid Sodium Salt Solution
Albumin Bovine (BSA) (pH 7) fraction V PrimeDx
6-Aminoisoquinoline
Bellidifolin (BLF)
Bis(Triphenylphosphine) Palladium(II) Dichloride
5-Bromo-1-Pentene
Carveol
D-Carvone
Chloro(1,5-cyclooctadiene) Rhodium (I) Dimer
2-(2-((6-Chlorohexyl)oxy)Ethoxy) Ethylamine Hydrochloride
Cholesterol Oxidase ex. Microorganism
p-Coumaric Acid
D-Cysteine Hydrochloride Monohydrate
DBCO Acid
Diacetato [(R)-(+)-2,2'-bis(Diphenylphosphino)-1,1'-Binaphthyl]Ruthenium (II))
Diacetone Alcohol
Diatrizoic Acid Dihydrate
Diatrizoic Acid Dihydrate ExiPlus, Multi-Compendial
Dichloro (Cycloocta-1,5-Diene) Ruthenium (II)
3,5-Dichlorophenyldiazonium Tetrafluroborate
Dihydrocarvone
Dihydrocarvyl Acetate (mix of Isomers)
Dysprosium Oxide
Ethyl Bromodifluoroacetate
Fluorescein Diacetate
5-Fluorouracil (5-FU) Exiplus, Multi-Compendial
Gelrite Gellan Gum (High Acyl)
Glycidyl Trimethyl Ammonium Chloride (GTMAC)
Hematoxylin Monohydrate PrimeDx
Hexachloroplatinic Acid
Hexamethylenediamine (HMDA)
Lutein

New Introductions

2-Naphthoyl Chloride
Nickel (II) Acetylacetonate Anhydrous
ß-Nicotinamide Mononucleotide Disodium Salt (Reduced) (ß-NMNH)
Nicotinamide Riboside Chloride
Nicotinamide Riboside Malate
Oregano Oil
Palladium Sponge
Peroxidase ex. Horseradish RZ >3.0, w/o Stab., Salt free (HRP Type 1B)
Phosphoenolpyruvate Carboxylase (PPC) ex. Microorganism
Platinum Sponge
Potassium Amylxanthate (PAX)
Protocatechuate-3,4-Dioxygenase (PCD) ex. Pseudomonas Sp.
D-Pyroglutamic Acid
Pyrroloquinoline Quinone (PQQ)
Pyrroloquinoline Quinone Disodium Salt (PQQ Disodium)
Rhodium (II) Acetate Dimer
Ruthenium AcetylAcetonate
Ruthenium (III) Chloride
Anhydrous
Safranine O PrimeDx
Sebacoyl Chloride
Sodium Isopropyl Xanthate (SIPX)
Sucrose Octasulphate Potassium Salt
1,2,4-Trichlorobenzene (1,2,4-TCB)
2.6M Tris Hydrochloride Buffer pH-8 for molecular biology
Tris(Triphenylphosphine) Rhodium (I) Chloride
Tris(Triphenylphosphine) Ruthenium (II) Dichloride
Tyloxapol
Urolithin A
Xylidine Ponceau 2R
Zirconium Hydride

Part B1 - Animal Cell Culture

Ceftriaxone Disodium Salt Hemiheptahydrate (CFTZ) for cell culture, Endotoxin (BET) 0.2EU/mg
Dulbecco's Modified Eagles Medium 1X w/ 4500mg/lit Glucose, L-Glutamine, Sodium Pyruvate, Sodium Bicarbonate, Phenol Red
Dulbecco's Phosphate Buffered Saline 1X w/o Calcium & Magnesium
Dulbecco's Phosphate Buffered Saline 10X w/o Calcium and Magnesium
RPMI 1640 Medium 1X
RPMI 1640 Medium w/ L-Glutamine w/ 25mM HEPES

Part C - Nanomaterials and Renewable Energy

Boron Micronpowder

Part D - Dehydrated Culture Media

Ampicillin (AMP) AST Disc, 10mcg
Ampicillin (AMP) AST Disc, 25mcg
Azithromycin (AZM) AST Disc, 15mcg
Cefazolin (CZ) AST Disc, 30mcg
Cefixime (CFM) AST Disc, 5mcg
Cefoperazone (CPZ) AST Disc, 30mcg
Cefotaxime (CTX) AST Disc, 30mcg
Ceftazidime (CAZ) AST Disc, 30mcg
Chloramphenicol (CMP) AST Disc, 30mcg
Ciprofloxacin (CIP) AST Disc, 5mcg
Colistin (CL) AST Disc, 10mcg
Erythromycin (E) AST Disc, 15mcg
Gentamicin (GEN) AST Disc, 10mcg
Gentamicin (GEN) AST Disc, 120mcg
Nalidixic Acid (NA) AST Disc, 30mcg
Netilmicin (NET) AST Disc, 30mcg
Polymyxin B (PB) AST Disc, 50Units
Tetracycline (TE) AST Disc, 30mcg
Vancomycin (VA) AST Disc, 30mcg
Bromocresol Purple Agar
Haemoglobin Supplement
Letheen Broth (USP Basis)
Letheen Broth, Modifed (Twin Pack)
Schuberts Arginine Broth

Part A

General Laboratory Chemicals & Solvents

Sodium Lauryl Sulphate

A highly versatile product for Drug Dissolution Testing, Quality Control, Research and Electrophoresis applications.



For Drug Dissolution Testing Studies

Sodium Lauryl Sulphate (Sodium Dodecyl Sulphate, SLS, SDS) is a strong component of the media used for Dissolution Testing in In Vitro methods for water-soluble and slightly-water-soluble drugs. It's importance as a simulated gastric/intestinal media ingredient makes it's quality a critical factor in the accuracy of results.

Principally used as a QC Test, Dissolution Testing, SDS of AR & High Purity grades are the most prevalent reagents in use due to constant production improvements, testing and stability studies. SRL is one of the first developers of the High Purity SLS (99.5%) which is fast-becoming the industry standard across the Pharma QC Testing, with many scientists now adopting the newer improved grade as a Protocol standard.

Our Premier Product, Sodium Lauryl Sulphate High Purity (99.5%) (35825) can,

- Eliminate further purification by the customer by foaming & other processes, demonstrate exceptionally lower lauryl alcohol (LA) impurities
- Improve testing accuracy and confirm lot-to-lot bioequivalence
- Ensure consistent and reliable results in In Vitro Testing and Intrinsic Dissolution Rate (IDR) studies
- Secure accurate Dissolution profiles (suitable to fast-release and slow-release drugs)
- Assist in detection of relevant product changes so as to ensure the quality and consistent performance of final products

For SDS Polyacrylamide Gel Electrophoresis (SDS-PAGE) Studies

Primarily used in the denaturing of proteins, Sodium Lauryl Sulphate of high purity and molecular biology grades play a critical role in PAGE electrophoresis techniques as well as SDS Gradient Gel Electrophoresis to,

- Unfold normally buried polypeptide chains while retaining their normal (native) structure
- Bind to difficult regions (hydrophobic regions) within the protein chains
- Induce a predominant negative charge for the subsequent electrophoresis
- Ensure a electrophoretic separation for a large number (and types) of proteins

Types of Sodium Lauryl Sulphate we offer:

Product code	Product Name	Assay	Pack Size
54468	Sodium Lauryl Sulphate extrapure AR, ACS	(C12) min. 99% (GC)	100g, 250g, 500g, 1kg
99671	Sodium Lauryl Sulphate for HPLC	(C12) min. 99% (GC)	10g, 25g, 100g, 1kg
32096	Sodium Lauryl Sulphate for molecular biology	(C12) min. 99.5% (GC)	25g, 100g, 500g
35825	Sodium Lauryl Sulphate High Purity	(C12) min. 99.5% (GC)	25g, 100g, 500g
85369	Sodium Lauryl Sulphate (10%) solution	10%	100ml
87545	Sodium Lauryl Sulphate (20%) solution	20%	250ml
34378	Sodium Lauryl Sulphate BP	min. 85%	1kg
93117	Sodium Lauryl Sulphate IP	min. 85%	1kg
14374	Sodium Lauryl Sulphate pure	min. 85% (dry basis)	500g, 5kg

Larger pack sizes 1kg, 5kg, 25kg available on request

Code	Name	Packing	Price
95237 [60-35-5]	**Acetamide pure, 98%**	500 Gms	850
81658 [103-84-4]	**Acetanilide extrapure, 99%**	500 Gms	950
68175 [103-84-4]	**Acetanilide ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP	500 Gms	1350
18439 [64-19-7]	**Acetic Acid 1N Aq. Solution**	1000 ml	610
72690 [64-19-7]	**Acetic Acid 2.5M Aq. Solution**	1000 ml	650
90868 [64-19-7]	**Acetic Acid Glacial extrapure, 99.5%**	500 ml 2500 ml 25 Ltr	240 1100 9000
85801 [64-19-7]	**Acetic Acid Glacial ACS, 99.9%**	500 ml 2500 ml	260 1250
93602 [64-19-7]	**Acetic Acid Glacial ACS, ExiPlus, Multi-Compendial, 99.9%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	290 1350
73215 [64-19-7]	**Acetic Acid Dried, 99.9%, water 0.05%**	500 ml 1000 ml	380 740
32532 [64-19-7]	**Acetic Acid for HPLC, 99.9%**	500 ml 1000 ml 2500 ml	520 980 2100
59788 [64-19-7]	**Acetic Acid Glacial for molecular biology, 99.9%**	100 ml 500 ml	520 1600
15868 [1186-52-3]	**Acetic Acid-d4 for NMR spectroscopy, 99.5 Atom %D** (Acetic-d3 Acid-d,Tetradeuteroacetic Acid)	5 ml	14000
31566 [67-64-1]	ⓘ **Acetone pure, 99%**	500 ml 2500 ml 25 Ltr	270 1100 10000
66951 [67-64-1]	ⓘ **Acetone extrapure, 99%**	500 ml 2500 ml 25 Ltr	285 1150 10800
15168 [67-64-1]	ⓘ **Acetone ACS, 99.5%**	500 ml 2500 ml	310 1250
27648 [67-64-1]	ⓘ **Acetone ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP	500 ml 2500 ml	340 1350
35007 [67-64-1]	ⓘ **Acetone Dried, 99%, water 0.005%**	500 ml 2500 ml	380 1500
89140 [67-64-1]	ⓘ **Acetone GC-HS, 99.9%**	250 ml 500 ml 1000 ml	450 700 1050
11340 [67-64-1]	ⓘ **Acetone for HPLC & UV Spectroscopy, 99.9%**	500 ml 1000 ml 2500 ml	500 800 1550
12312 [67-64-1]	ⓘ **Acetone electronic grade, 99.9%**	1000 ml 2500 ml	850 1900
55245 [67-64-1]	ⓘ **Acetone VLSI Semiconductor Grade, 99.5%**	1000 ml 2500 ml	12000 26000
27498 [67-64-1]	ⓘ **Acetone for molecular biology, 99.8%**	500 ml 2500 ml	550 2200
12972 [666-52-4]	🌡 **Acetone-d6 for NMR spectroscopy, 99.5 Atom %D**	10 ml 25 ml	14000 35000
75185 [666-52-4]	🌡 **Acetone-d6 (with 1% TMS) for NMR spectroscopy, 99.5 Atom %D**	10 ml	15500
49967 [75-05-8]	**Acetonitrile (ACN) extrapure, 99%**	500 ml 2500 ml 25 Ltr	500 2350 21000
84578 [75-05-8]	**Acetonitrile (ACN) ACS, 99.5%**	500 ml 2500 ml	600 2500
82745 [75-05-8]	**Acetonitrile (ACN) ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	700 2700
90135 [75-05-8]	**Acetonitrile (ACN) GC-HS, 99.9%**	250 ml 500 ml 1000 ml	650 1000 1700
58209 [75-05-8]	**Acetonitrile (ACN) for HPLC & UV Spectroscopy, 99.9%**	500 ml 2500 ml 2.5 Ltr (Wide-Mouth) 4 Ltr (Wide-Mouth)	800 3500 5000 8000
24899 [75-05-8]	**Acetonitrile (ACN) Gradient grade for HPLC, 99.9%**	1000 ml 2500 ml 2.5 Ltr (Wide-Mouth) 4 Ltr (Wide-Mouth)	1800 4000 6000 9500
34396 [75-05-8]	**Acetonitrile (ACN) for LC-MS & HPLC, 99.9%**	1 Ltr (Wide-Mouth) 2.5 Ltr (Wide-Mouth) 4 Ltr (Wide-Mouth)	12000 24000 37000
81506 [75-05-8]	**Acetonitrile (ACN) Ultra Gradient for HPLC, UHPLC, 99.9%**	1 Ltr (Wide-Mouth) 2.5 Ltr (Wide-Mouth) 4 Ltr (Wide-Mouth)	9000 20000 31000
62006 [75-05-8]	**Acetonitrile (ACN) for molecular biology, 99.9%**	500 ml	1200

Code	Name	Packing	Price
62076 [75-05-8]	**Acetonitrile (ACN) for DNA synthesis, 99.9%**	1000 ml 2500 ml	2400 5600
26483 [2206-26-0]	**Acetonitrile-d3 for NMR spectroscopy, 99.8 Atom %D** (Methyl-d3 Cyanide, Trideuteroacetonitrile)	10 Gms	18000
35845 [2206-26-0]	**Acetonitrile-d3 for NMR spectroscopy, 99.8 Atom %D** (Methyl-d3 Cyanide, Trideuteroacetonitrile)	10 ml	16000
67418	**Aceto Orcein, 1% in Dil Acetic Acid**	100 ml	580
15873 [98-86-2]	**Acetophenone extrapure, 99%**	500 ml 2500 ml	580 2650
22506 [98-86-2]	**Acetophenone Reagent Grade, 99%**	500 ml 2500 ml	750 3570
96169 [123-54-6]	**Acetylacetone pure, 99%**	500 ml 2500 ml	1150 4830
24212 [123-54-6]	**Acetylacetone extrapure, 99.5%**	250 ml 1000 ml	810 2190
39522 [3244-88-0]	**Acid Fuchsin** (C.I. No. 42685, Acid Violet 19, Fuchsin S, Fuchsin Acid, Rubine S)	25 Gms 100 Gms 500 Gms	700 2300 11400
50190 [123334-10-1]	**Acid Fuchsin Calcium Salt (Acid Rubin), 60%** (C.I. No. 42685) (Acid Magenta, Rubin S)	25 Gms 100 Gms	3800 13900
22794 [79-06-1]	**Acrylamide extrapure, 99%**	500 Gms 5 Kg 25 Kg	500 3700 15800
15657 [79-06-1]	**Acrylamide 3x cryst. extrapure AR, 99.9%**	500 Gms 1 Kg 5 Kg	1600 3000 13700
61346 [79-06-1]	**Acrylamide 3x cryst. for molecular biology, 99.9%**	100 Gms 500 Gms 1 Kg	900 2500 4000
29787 [79-06-1]	**Acrylamide 40% Aq. Solution for molecular biology**	250 ml 500 ml	2200 4000
24263	**Acrylamide / Bis-acrylamide Premix Powder, Ratio 19:1**	30 Gms 200 Gms	1200 6300
10762	**Acrylamide / Bis-acrylamide Premix Powder, Ratio 29:1**	30 Gms 200 Gms	1200 6300
74379	**Acrylamide / Bis-acrylamide Premix Powder, Ratio 37.5:1**	30 Gms 200 Gms	1200 6300
67394	**30% Acrylamide / Bis-acrylamide Mix Solution (Ratio 29:1)**	100 ml 500 ml	1200 4500
92035	**40% Acrylamide / Bis-acrylamide Mix Solution (Ratio: 37.5:1)**	100 ml 500 ml	1200 4500
60452 [79-10-7]	**Acrylic Acid (Stabilized) pure, 99%**	500 ml 2500 ml	550 2500
41666 [107-13-1]	**Acrylonitrile (Stabilized) pure, 99%**	500 ml 2500 ml	350 1500
60064 [124-04-9]	**Adipic Acid pure, 99%**	500 Gms 5 Kg	500 4000
38573 [124-04-9]	**Adipic Acid ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	580 4500
19661 [9002-18-0]	**Agar powder regular grade for bacteriology**	100 Gms 250 Gms 500 Gms 1 Kg 5 Kg	1030 2220 4260 8400 40830
85473 [9002-18-0]	**Agar powder for tissue culture**	100 Gms 250 Gms 500 Gms 5 Kg	1400 2450 4670 44330
24970 [9002-18-0]	**Agar powder Bacto grade**	100 Gms 500 Gms	1640 6540
26741 [8001-99-8]	**Ajowan Seed Oil extrapure** (Ajwain Oil)	100 ml 500 ml	1290 5330
48863 [107-95-9]	**b-Alanine extrapure CHR, 99%**	25 Gms 100 Gms 1 Kg	300 1010 8630
71178 [56-41-7]	**L-Alanine extrapure CHR, 99%**	100 Gms 500 Gms 1 Kg 5 Kg	540 1800 3300 14500
71777 [56-41-7]	**L-Alanine for tissue culture, 99%**	100 Gms 500 Gms	650 2200
50683 [56-41-7]	**L-Alanine ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP	100 Gms 500 Gms 5 Kg	740 2400 17500
54479 [56-41-7]	**L-Alanine PrimeDx, 99%**	1 Kg 5 Kg	5000 22000
87279	**Alberts Stain A Solution**	125 ml 500 ml	190 480
35652	**Alberts Stain B Solution**	125 ml 500 ml	190 480
64199 [1324-76-1]	**Alkali Blue 6B** (Acid Blue 110, Acid Blue 119, C.I. 42765)	5 Gms 25 Gms 100 Gms	600 2300 8000
17996 [63393-96-4]	**Aliquat 336, 88%** (Tricaprylylmethylammonium chloride)	250 ml 1000 ml 2500 ml 25 Ltr	1500 6330 14950 138000

Code	Name	Packing	Price
43805 [97-59-6]	**Allantoin ExiPlus, 98.5%** (5-Ureidohydantoin, Glyoxylic Acid Diureide) meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms 5 Kg	580 2300 21280
99564 [2244-11-3]	**Alloxan Monohydrate extrapure, 98%** (2,4,5,6-Tetraoxypyrimidine, 5,6-Dioxyuracil)	10 Gms 25 Gms 100 Gms	430 940 3450
19093 [107-18-6]	**Allyl Alcohol extrapure, 99%**	500 ml 2500 ml	810 3450
69098 [107-11-9]	**Allylamine pure, 99%** (2-Propen-1-Ylamine, 3-Amino-1-Propene)	25 ml 100 ml 250 ml	690 1150 2300
45953 [124-02-7]	**Diallylamine extrapure, 99%**	100 ml 500 ml	1380 5060
73056 [106-95-6]	**Allyl Bromide pure, 98%** (3-Bromo-1-Propene)	250 ml 500 ml	2880 5520
92459 [107-05-1]	**Allyl Chloride pure, 98%** (3-Chloro-1-Propene)	500 ml 2500 ml	460 1730
45262 [8007-69-0]	**Almond Oil pure**	500 ml	7840
54568 [142-03-0]	**Aluminium Acetate Basic Hydrate (Aluminium Subacetate) extrapure, 98%** (Aluminum Hydroxide Acetate Hydrate)	100 Gms 500 Gms	1360 3320
64132 [7784-26-1]	**Aluminium Ammonium Sulphate Dodecahydrate pure, 99%** (Ammonium Alum, Alum Ammonia, Ammonium Aluminium Sulphate)	500 Gms	300
83147 [7784-26-1]	**Aluminium Ammonium Sulphate Dodecahydrate ACS, 98%** (Ammonium Alum, Alum Ammonia, Ammonium Aluminium Sulphate)	500 Gms	430
61315 [7784-26-1]	**Aluminium Ammonium Sulphate Dodecahydrate ExiPlus, Multi-Compendial, 99%** (Ammonium Alum, Alum Ammonia, Ammonium Aluminium Sulphate) meets compendial specs of USP	500 Gms 5 Kg	460 4490
47689 [7803-63-6]	**Ammonium Sulphate Monobasic extrapure, 98%** (Ammonium Bisulphate, Ammonium Hydrogen Sulphate, Monoammonium Sulfate)	500 Gms	1350
93606	**Aluminium AAS Standard Solution in 1N HCl (Traceable to NIST)**	100 ml 250 ml	2070 3340
90015 [7784-13-6]	**Aluminium Chloride Hexahydrate pure, 97%**	500 Gms 5 Kg	670 6330
23700 [7784-13-6]	**Aluminium Chloride Hexahydrate extrapure AR, 99%**	500 Gms 5 Kg	750 7130
35543 [21645-51-2]	**Aluminium Hydroxide Gel ExiPlus, Multi-Compendial** (Dried Aluminium Hydroxide) meets compendial specs of USP, BP	500 Gms 5 Kg	520 4830
12283 [555-31-7]	**Aluminium Isopropoxide (powder) extrapure, 98%** (AIP, Aluminum Isopropylate, Aluminum Triisopropoxide, 2-Propanol Aluminium Salt)	500 Gms 5 Kg	1330 12080
57395 [7429-90-5]	**Aluminium Metal powder, 99%, -325 mesh**	500 Gms	980
92309 [7784-27-2]	**Aluminium Nitrate Nonahydrate extrapure, 98%**	500 Gms 5 Kg	440 3680
19302 [7784-27-2]	**Aluminium Nitrate Nonahydrate ACS, 98.5%**	500 Gms	920
58095 [24304-00-5]	**Aluminium Nitride (AIN) pure, 98%, ~635 mesh**	5 Gms 50 Gms 250 Gms	3450 4700 22080
66289 [1344-28-1]	**Aluminium Oxide Activated (Acidic) (Alumina)**	500 Gms 5 Kg 25 Kg	630 5930 28430
71461 [1344-28-1]	**Aluminium Oxide Activated (Acidic) (Alumina) ExiPlus, Multi-Compendial** meets compendial specs of USP, BP	500 Gms 5 Kg	650 6280
34321 [1344-28-1]	**Aluminium Oxide Activated (Basic) (Alumina)**	500 Gms 5 Kg 25 Kg	500 4620 21760
61492 [1344-28-1]	**Aluminium Oxide Activated (Basic) (Alumina) ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	500 Gms 5 Kg	530 4780
39882 [1344-28-1]	**Aluminium Oxide activated (Neutral) (Alumina)**	500 Gms 5 Kg 25 Kg	520 4870 23200
62690 [1344-28-1]	**Aluminium Oxide G (Neutral) (Alumina G)**	500 Gms 5 Kg	690 6530
34641 [7784-24-9]	**Aluminium Potassium Sulphate Dodecahydrate extrapure, 99%** (Potash Alum, Potassium Aluminium Sulphate)	500 Gms 2.5 Kg	280 1180
42889 [7784-24-9]	**Aluminium Potassium Sulphate Dodecahydrate ACS, 99.5%** (Potash Alum, Potassium Aluminium Sulphate)	500 Gms 2.5 Kg	330 1520
42550 [7784-24-9]	**Aluminium Potassium Sulphate Dodecahydrate ExiPlus, Multi-Compendial, 99%** (Potash Alum, Potassium Aluminium Sulphate) meets compendial specs of USP, BP, Ph.Eur	500 Gms 2.5 Kg	380 1840
53016 [7784-31-8]	**Aluminium Sulphate Octadecahydrate extrapure, 98%**	500 Gms	260
11820 [569-58-4]	**Aluminon extrapure** (Aurintricarboxylic Acid Ammonium Salt, ATA, Ammonium Aurintricarboxylate)	10 Gms 25 Gms	950 2000
55745 [915-67-3]	**Amaranth (Acid Red 27) (C.I. 16185), 85-95%**	25 Gms 100 Gms 1 Kg	210 730 6330

Code	Name	Packing	Price
57432 [83-07-8]	**4-Aminoantipyrine (Ampyrone) ACS, 98%** (4-Amino-2,3-Dimethyl-1-Phenyl-3-Pyrazolin-5-One)	25 Gms 100 Gms 250 Gms	1370 4950 12080
69657 [83-07-8]	**4-Aminoantipyrine (Ampyrone) PrimeDx, 98%** (4-Amino-2,3-Dimethyl-1-Phenyl-3-Pyrazolin-5-One)	1 Kg 5 Kg	55200 258750
82623 [60-32-2]	**6-Aminocaproic Acid extrapure, 99%** (6-Aminohexanoic Acid)	25 Gms 100 Gms 500 Gms	1370 4350 21120
85996 [60-32-2]	**Aminocaproic Acid IP** (6-Aminohexanoic Acid, 6-Aminocaproic Acid IP)	1 Kg 5 Kg	36800 172500
99093 [60-32-2]	**Aminocaproic Acid BP, 98.5-101%** (6-Aminohexanoic Acid, 6-Aminocaproic Acid BP)	1 Kg 5 Kg	36800 172500
50445 [60-32-2]	**Aminocaproic Acid USP, 98.5-101%** (6-Aminohexanoic Acid, 6-Aminocaproic Acid USP)	1 Kg 5 Kg	41400 184000
78719 [1336-21-6]	**Ammonia Solution Reagent Grade, 25%** (Ammonium Hydroxide 25%)	500 ml 2500 ml 25 Ltr	210 640 6000
14388 [1336-21-6]	**Ammonia Solution ACS, 30%** (Ammonium Hydroxide 30%)	500 ml 2500 ml 25 Ltr	250 710 6500
44683 [631-61-8]	**Ammonium Acetate extrapure AR, 98%**	500 Gms 5 Kg	310 2900
62778 [631-61-8]	**Ammonium Acetate ACS, ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	340 3100
78844 [631-61-8]	**Ammonium Acetate for HPLC, 99%**	500 Gms	700
37829 [631-61-8]	**Ammonium Acetate for molecular biology, 98%**	100 Gms 500 Gms	310 870
26757 [1066-33-7]	**Ammonium Bicarbonate (Ammonium Hydrogen Carbonate) extrapure AR, 99%**	500 Gms 5 Kg 25 Kg	250 2200 7000
38586 [7789-09-5]	**Ammonium Dichromate pure, 98%** (Ammonium Bichromate)	500 Gms 5 Kg	990 9520
67077 [7789-09-5]	**Ammonium Dichromate extrapure AR, ACS, 99.5%** (Ammonium Bichromate)	500 Gms 5 Kg	1080 10360
13664 [1341-49-7]	**Ammonium Bifluoride pure, 98%**	500 Gms 5 Kg	900 8070
95407 [31886-41-6]	**Ammonium Bismuth Citrate extrapure AR, 43-46% Bi**	100 Gms 500 Gms	3030 13440
46813 [12124-97-9]	**Ammonium Bromide pure, 99%**	500 Gms 5 Kg	1100 10500
40553 [16774-21-3]	**Ammonium Ceric Nitrate ACS, 99%** (Ammonium Cerium (IV) Nitrate, Ceric Ammonium Nitrate)	100 Gms 250 Gms	900 2000
42077 [16774-21-3]	**Ammonium Ceric Nitrate ExiPlus, Multi-Compendial, 99%** (Ammonium Cerium (IV) Nitrate, Ceric Ammonium Nitrate) meets compendial specs of USP, BP	100 Gms 250 Gms	1200 2500
42725 [16774-21-3]	**Ammonium Ceric Nitrate 0.05N Solution**	500 ml	300
23757 [10378-47-9]	**Ammonium Ceric Sulphate Dihydrate ACS, 94%** (Ammonium Cerium (IV) Sulphate Dihydrate)	100 Gms 500 Gms	1200 3800
25103 [12125-02-9]	**Ammonium Chloride ACS, 99.5%**	500 Gms 5 Kg 25 Kg	290 2700 12980
96452 [12125-02-9]	**Ammonium Chloride ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	340
16992 [12125-02-9]	**Ammonium Chloride for molecular biology, 99.5%**	500 Gms	660
35309 [12125-02-9]	**Ammonium Chloride for tissue culture, 99.5%**	500 Gms 5 Kg	320 3000
11505 [3458-72-8]	**Ammonium Citrate Tribasic (Triammonium Hydrogen Citrate) extrapure AR, 98.5%**	500 Gms 5 Kg	750 7000
35442 [1185-57-5]	**Ammonium Ferric Citrate pure, 16.5-22.5%** (Ammonium Iron (III) Citrate)	500 Gms 5 Kg	600 5750
63571 [7783-83-7]	**Ammonium Ferric Sulphate Dodecahydrate extrapure, 98%** (Ammonium Iron (III) Sulphate Dodecahydrate, Ferric Alum)	500 Gms 5 Kg 25 Kg	310 2900 13000
29168 [7783-83-7]	**Ammonium Ferric Sulphate Dodecahydrate ACS, 98.5%** (Ammonium Iron (III) Sulphate Dodecahydrate, Ferric Alum)	500 Gms 5 Kg 25 Kg	370 3200 14500
89944 [7783-85-9]	**Ammonium Ferrous Sulphate Hexahydrate extrapure, 98%** (Diammonium Iron (II) Sulphate Hexahydrate)	500 Gms 1 Kg 5 Kg 25 Kg	220 390 1700 7000
61732 [7783-85-9]	**Ammonium Ferrous Sulphate Hexahydrate ACS, 99%** (Diammonium Iron (II) Sulphate Hexahydrate)	500 Gms 5 Kg 25 Kg	300 2200 8000
11605 [13826-83-0]	**Ammonium Fluoroborate pure, 98%** (Ammonium fluoroborate, Ammonium tetrafluoroborate)	100 Gms 500 Gms 5 Kg	500 1400 12000
49813 [12125-01-8]	**Ammonium Fluoride pure, 98%**	500 Gms	1070
38431 [12125-01-8]	**Ammonium Fluoride ACS, 98%**	500 Gms	1240
25142 [540-69-2]	**Ammonium Formate pure, 98%** (Formic Acid Ammonium Salt)	500 Gms 5 Kg	390 3360
27858 [540-69-2]	**Ammonium Formate extrapure AR, 98%** (Formic Acid Ammonium Salt)	500 Gms 5 Kg	420 3590

Code	Name	Packing	Price
50504 [540-69-2]	**Ammonium Formate extrapure AR, ExiPlus, Multi-Compendial, 98%** (Formic Acid Ammonium Salt) meets compendial specs of BP, Ph.Eur	500 Gms 5 Kg	460 3920
14969 [16941-11-0]	**Ammonium Hexafluorophosphate extrapure, 95%** (Ammonium Hexafluorophosphate V, Ammonium Phosphorus Hexafluoride)	5 Gms 25 Gms 100 Gms 250 Gms	1900 6300 15600 35000
51775 [12027-06-4]	**Ammonium Iodide extrapure, 99%**	50 Gms 100 Gms 250 Gms	2760 5060 11620
23185 [7803-55-6]	**Ammonium Metavanadate extrapure, 98%**	100 Gms 500 Gms	1040 3800
15068 [7803-55-6]	**Ammonium Metavanadate ACS, 99%**	100 Gms 500 Gms	1210 4030
94314 [7803-55-6]	**Ammonium Metavanadate ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	1340 4260
69429 [12054-85-2]	**Ammonium Molybdate Tetrahydrate pure, 98%** (Ammonium Heptamolybdate)	100 Gms 250 Gms 500 Gms	1840 4550 8050
10299 [12054-85-2]	**Ammonium Molybdate Tetrahydrate ACS, 81-83% MoO3** (Ammonium Heptamolybdate)	100 Gms 250 Gms 500 Gms	2070 4950 8630
82615 [12054-85-2]	**Ammonium Molybdate Tetrahydrate ExiPlus, Multi-Compendial, 99.3%** (Ammonium Heptamolybdate) meets compendial specs of USP, BP, Ph.Eur	100 Gms 250 Gms 500 Gms	2190 5060 9200
92887 [12054-85-2]	**Ammonium Molybdate Tetrahydrate PrimeDx, 99%** (Ammonium Heptamolybdate)	1 Kg 5 Kg	10350 48300
32910 [27546-07-2]	**Ammonium Dimolybdate extrapure, 99%, 56.5% Mo**	100 Gms 250 Gms	2200 5200
85044 [7785-20-8]	**Ammonium Nickel Sulphate Hexahydrate pure, 98%** (Diammonium Nickel (II) Sulphate Hexahydrate)	500 Gms	1840
55517 [7785-20-8]	**Ammonium Nickel Sulphate Hexahydrate extrapure AR, 98%** (di-Ammonium nickel (II) sulphate hexahydrate)	500 Gms	2760
85084 [6009-70-7]	**Ammonium Oxalate Monohydrate extrapure, 98%**	500 Gms	410
72913 [6009-70-7]	**Ammonium Oxalate Monohydrate ACS, 99%**	500 Gms	440
53621 [6009-70-7]	**Ammonium Oxalate Monohydrate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	510
82381 [7727-54-0]	**Ammonium Persulphate (APS) pure, 98%** (Ammonium Peroxodisulphate)	500 Gms 5 Kg	440 4000
84569 [7727-54-0]	**Ammonium Persulphate (APS) ACS, 98%** (Ammonium Peroxodisulphate)	500 Gms 5 Kg	490 4110
29534 [7727-54-0]	**Ammonium Persulphate (APS) ExiPlus, Multi-Compendial, 99%** (Ammonium Peroxodisulphate) meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	520 5010
28575 [7727-54-0]	**Ammonium Persulphate (APS) for electrophoresis, 99%** (Ammonium Peroxodisulphate)	25 Gms 100 Gms	260 790
65553 [7727-54-0]	**Ammonium Persulphate (APS) for molecular biology, 99%** (Ammonium peroxodisulphate)	25 Gms 100 Gms 500 Gms	280 690 2250
30468 [7783-28-0]	**Ammonium Phosphate Dibasic ACS, 98%** (Diammonium Hydrogen Orthophosphate)	500 Gms 5 Kg	600 5750
65984 [7783-28-0]	**Ammonium Phosphate Dibasic ExiPlus, Multi-Compendial, 98%** (Diammonium Hydrogen Orthophosphate) meets compendial specs of USP, Ph.Eur	500 Gms 5 Kg	680 6560
61563 [7722-76-1]	**Ammonium Phosphate Monobasic extrapure AR, 99%** (Ammonium Dihydrogen Orthophosphate)	500 Gms 5 Kg	580 5580
58611 [7722-76-1]	**Ammonium Phosphate Monobasic extrapure AR, ACS, ExiPlus, Multi-Compendial , 99%** (Ammonium Dihydrogen Orthophosphate) meets compendial specs of USP	500 Gms 5 Kg	660 6330
24678 [7722-76-1]	**Ammonium Phosphate Monobasic for HPLC, 99%** (Ammonium Dihydrogen Orthophosphate)	500 Gms	1380
96648 [5108-96-3]	**Ammonium Pyrrolidine Dithiocarbamate (APDC) extrapure, 98.5%** (1-Pyrrolidinecarbodithioic Acid Ammonium Salt, Ammonium pyrrolidinecarbodithioate, PDC)	10 Gms 25 Gms	2990 6560
86689 [7773-06-0]	**Ammonium Sulphamate ACS, 98%**	100 Gms 500 Gms	350 1380
62419 [7773-06-0]	**Ammonium Sulphamate Ultrapure, 99.5%**	100 Gms 500 Gms	2300 7590
50012 [7773-06-0]	**Ammonium Sulphamate Ultrapure, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	2880 8630
88064 [7783-20-2]	**Ammonium Sulphate ACS, 99.5%** (Ammonium Sulfate)	500 Gms 5 Kg 25 Kg	290 2250 10930
88004 [7783-20-2]	**Ammonium Sulphate for tissue culture, 99.5%**	500 Gms 5 Kg	330 2420

Code	Name	Packing	Price
26339 [7783-20-2]	**Ammonium Sulphate ExiPlus, Multi-Compendial, 99%** (Ammonium Sulfate) meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	350 2530
82126 [7783-20-2]	**Ammonium Sulphate for molecular biology, 99.5%** (Ammonium Sulfate)	500 Gms 5 Kg	1380 11500
95876 [3164-29-2]	**Ammonium Tartrate Dibasic pure, 99%**	500 Gms	2760
56357 [3164-29-2]	**Ammonium Tartrate Dibasic extrapure AR, 99%**	250 Gms 500 Gms	1960 3450
81016 [1762-95-4]	**Ammonium Thiocyanate ACS, 97.5%**	500 Gms	690
70158 [1762-95-4]	**Ammonium Thiocyanate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	790
67550 [123-92-2]	**Isoamyl Acetate pure, 98%**	500 ml 2500 ml	640 2530
28605 [71-41-0]	**n-Amyl Alcohol (1-Pentanol) pure, 98%**	500 ml 2500 ml	810 3910
31906 [71-41-0]	**n-Amyl Alcohol (1-Pentanol) ACS, 98%**	500 ml 2500 ml	980 4600
83344 [28983-56-4]	**Aniline Blue (water soluble) (Methyl Blue)** (Acid Blue 93)	25 Gms 100 Gms	810 2800
22285 [62-53-3]	**Aniline ACS, 99.5%**	500 ml 2500 ml	650 2950
46774 [142-04-1]	**Aniline Hydrochloride pure, 99%**	250 Gms 500 Gms	820 1450
69980 ✔ [542-16-5]	**Aniline Sulphate pure, 98%** (Aminobenzene, Phenylamine, Benzenamine)	250 Gms 500 Gms	900 1670
36381 [104-94-9]	**p-Anisidine pure, 98%** (4-Aminoanisole, 4-Methoxyaniline)	100 Gms 500 Gms	290 1040
41010 [120-12-7]	**Anthracene pure, 97%**	25 Gms 100 Gms 500 Gms	400 1200 4000
61273 [120-12-7]	**Anthracene (High Purity) scintillation grade, 99%**	25 Gms 100 Gms 250 Gms	1550 5800 14210
78531 [84-65-1]	**Anthraquinone pure, 98%** (9,10-Anthraquinone, Anthradione)	25 Gms 500 Gms	340 1300
97400	**Antimony AAS Standard Solution in 2N HCl (Traceable to NIST)**	100 ml 250 ml	2720 4420
66105 [7440-36-0]	**Antimony Metal Lumps, 98%**	100 Gms 500 Gms	3000 14000
12042 [7440-36-0]	**Antimony Metal Powder extrapure, 99%, -200 mesh**	100 Gms 250 Gms	3400 8000
67577 [7440-36-0]	**Antimony Metal Ingots ultrapure, 99.9%**	10 Gms 25 Gms	1100 2000
47448 [28300-74-5]	**Antimony Potassium Tartrate Hemihydrate pure, 98%**	100 Gms 250 Gms 500 Gms	1500 3300 6200
83896 [28300-74-5]	**Antimony Potassium Tartrate Hemihydrate extrapure AR, 99.5%**	100 Gms 500 Gms	1650 6800
60274 [10025-91-9]	**Antimony Trichloride pure, 98%**	100 Gms 500 Gms	2300 10000
63365 [10025-91-9]	**Antimony Trichloride ACS, 99%**	100 Gms 500 Gms	2600 11000
97628 [10025-91-9]	**Antimony Trichloride ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	2800 12000
37157 [1309-64-4]	**Antimony Trioxide extrapure, 99%**	25 Gms 100 Gms 500 Gms	600 2300 10200
98063 [1309-64-4]	**Antimony Trioxide extrapure AR, 99%**	25 Gms 100 Gms 500 Gms	700 2500 11000
82829 [1345-04-6]	**Antimony Trisulfide pure, 98%**	100 Gms 500 Gms	1850 8200
24766 [68425-36-5]	**Arachis Oil** (Peanut Oil)	500 ml	600
66637 [74-79-3]	**L-Arginine (free base) extrapure CHR, 99%**	100 Gms 1 Kg 5 Kg	600 5000 23500
16178 [74-79-3]	**L-Arginine for tissue culture, 99%**	100 Gms 1 Kg	720 5330
17863 [74-79-3]	**L-Arginine (free base) ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP	25 Gms 100 Gms 1 Kg 5 Kg	300 880 6000 27500
40495 [74-79-3]	**L-Arginine (free base) PrimeDx, 99%**	1 Kg 5 Kg	6500 30000
55453 [1119-34-2]	**L-Arginine Hydrochloride extrapure CHR, 99%**	100 Gms 1 Kg 5 Kg	600 4800 21000
23006 [50-81-7]	**L-Ascorbic Acid ACS, 99%** (Vitamin C)	100 Gms 500 Gms 5 Kg	300 800 6800
40704 [50-81-7]	**L-Ascorbic Acid ExiPlus, Multi-Compendial, 99%** (Vitamin C) meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms 5 Kg	300 900 8000
14116 [50-81-7]	**L-Ascorbic Acid for molecular biology, 99.7%** (Vitamin C)	25 Gms 100 Gms 500 Gms	300 600 1100

Code	Name	Packing	Price
32488 [50-81-7]	**L-Ascorbic Acid for tissue culture, 99.7%** (Vitamin C)	100 Gms 500 Gms	600 1500
12535 [5743-28-2]	**L-Ascorbic Acid Calcium Salt Dihydrate extrapure, 99%** (Calcium-L-Ascorbate Dihydrate)	100 Gms 500 Gms	900 3500
65265 [134-03-2]	**L-Ascorbic Acid Sodium Salt extrapure, 99%** (Sodium-L-Ascorbate, Vitamin C Sodium Salt)	100 Gms 500 Gms	380 1850
64848 [617-45-8]	**DL-Aspartic Acid extrapure CHR, 99%**	100 Gms 500 Gms 1 Kg	560 2600 4500
84025 [56-84-8]	**L-Aspartic Acid extrapure CHR, 99%**	100 Gms 500 Gms 1 Kg 5 Kg	400 1200 2100 9500
81245 [56-84-8]	**L-Aspartic Acid ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP	100 Gms 1 Kg	460 2500
86229 [56-84-8]	**L-Aspartic Acid PrimeDx, 99%**	1 Kg 5 Kg	2300 11000
61560 ☐ [1115-63-5]	**L-Aspartic Acid Potassium Salt extrapure, 98%** (L-Asparaginic Acid Potassium Salt, Potassium L-Aspartate, Potassium Aspartate)	100 Gms 500 Gms 5 Kg	1800 6000 25000
90156 🌡 [22839-47-0]	**Aspartame ExiPlus, 98%** (Asp-Phe Methyl Ester, N-(L-a-Aspartyl) -L-Phenylalanine Methyl Ester, Asp-Phe-OMe, Aspartame)	5 Gms 25 Gms	1200 3200
78087 [2465-27-2]	**Auramine O** (C.I. No. 41000, Basic Yellow 2)	25 Gms 100 Gms 500 Gms	170 440 1620
41574 [25641-18-3]	**Azocarmine G (C.I. 50085) for histology** (Acid Red 101, Rosinduline)	5 Gms 25 Gms	600 1700
16411 [67-52-7]	**Barbituric Acid extrapure, 99%** (2,4,6-Trihydroxypyrimidine, Malonylurea)	100 Gms 500 Gms	800 3000
21154 [67-52-7]	**Barbituric Acid extrapure AR, 99%** (2,4,6-Trihydroxypyrimidine, Malonylurea)	25 Gms 100 Gms	340 1350
10013 [543-80-6]	**Barium Acetate extrapure, 98%**	500 Gms 5 Kg	690 6330
16985 [543-80-6]	**Barium Acetate ACS, 99%**	500 Gms 5 Kg	750 6620
45569 [543-80-6]	**Barium Acetate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	870
56009	**Barium AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml 250 ml	2070 3340
32005 [7791-28-8]	**Barium Bromide Dihydrate extrapure, 98%**	500 Gms	1840
35145 [513-77-9]	**Barium Carbonate pure, 98.5%**	500 Gms 5 Kg	910 8770
58434 [513-77-9]	**Barium Carbonate ACS, 99%**	500 Gms	1100
66323 [513-77-9]	**Barium Carbonate ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	500 Gms	1280
78118 [10326-27-9]	**Barium Chloride Dihydrate pure, 98%**	500 Gms 5 Kg	410 3680
68554 [10326-27-9]	**Barium Chloride Dihydrate ACS, 99%**	500 Gms 5 Kg	460 4350
92996 [10326-27-9]	**Barium Chloride Dihydrate Exiplus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	520 4720
49669 [10294-40-3]	**Barium Chromate extrapure, 98%**	500 Gms	870
86732 [7787-32-8]	**Barium Fluoride pure, 97%**	500 Gms	1040
85499 [12230-71-6]	**Barium Hydroxide Octahydrate pure, 98%**	500 Gms 5 Kg	350 3220
43469 [12230-71-6]	**Barium Hydroxide Octahydrate ACS, 98%**	500 Gms	500
84482 [12230-71-6]	**Barium Hydroxide Octahydrate ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	670
10290 [10022-31-8]	**Barium Nitrate extrapure, 98%**	500 Gms 5 Kg	690 6440
88030 [10022-31-8]	**Barium Nitrate ACS, 99%**	500 Gms	980
80355 [7727-43-7]	**Barium Sulphate extrapure, 97%**	500 Gms 5 Kg	320 2760
90115 [7727-43-7]	**Barium Sulphate extrapure AR, 98%**	500 Gms 5 Kg	480 4490
40251 [7727-43-7]	**Barium Sulphate extrapure AR, ExiPlus, Multi-Compendial, 98%** meets compendial specs of Ph.Eur	500 Gms 5 Kg	560 5180
24031 [8012-89-3]	**Bees Wax pure**	500 Gms	920
82949	**Benedict`s Reagent Quantitative**	500 ml 2500 ml	440 2020
35003	**Benedict`s Reagent Qualitative**	500 ml 2500 ml	330 870
92118 [1302-78-9]	**Bentonite pure** (Aluminium Silicate Hydrate)	500 Gms 5 Kg	260 2400
37278 [100-52-7]	**Benzaldehyde pure, 98%**	500 ml 2500 ml 25 Ltr	580 2530 20700
34087 [100-52-7]	**Benzaldehyde Reagent Grade, 99%**	500 ml 2500 ml	720 2880

Code	Name	Packing	Price
82149 [100-52-7]	**Benzaldehyde ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	870 3220
43648 [68391-01-5]	**Benzalkonium Chloride (BKC), 50% soln in water**	500 ml 5 Ltr	490 3680
87688 [55-21-0]	**Benzamide pure, 98%**	100 Gms 500 Gms	350 1610
99515 [71-43-2]	**Benzene (Sulphur free) extrapure, 99%**	500 ml 2500 ml 25 Ltr	260 1040 10010
16410 [71-43-2]	**Benzene ACS, 99.5%**	500 ml 2500 ml	310 1250
58857 [71-43-2]	**Benzene ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	350 1440
95582 [71-43-2]	**Benzene GC-HS, 99.9%**	500 ml 1000 ml	610 1160
36823 [71-43-2]	**Benzene for HPLC & UV Spectroscopy, 99.8%**	500 ml 1000 ml	470 860
69981 [71-43-2]	**Benzene Dried, 99.7%, water 0.005%**	1000 ml	940
31332 [1076-43-3]	**Benzene-d6 for NMR spectroscopy, 99.5 Atom %D**	10 ml 25 ml	12000 28000
27094 [121-54-0]	**Benzethonium Chloride extrapure, 98%** (Hyamine 1622, (Diisobutylphenoxyethoxyethyl)dimethyl benzylammonium Chloride, Phemerol Chloride)	25 Gms 100 Gms 250 Gms	800 2400 5000
36805 [65-85-0]	**Benzoic Acid pure, 99%**	500 Gms 5 Kg	330 2750
64239 [65-85-0]	**Benzoic Acid ACS, 99.5%**	500 Gms 5 Kg	720 6380
10105 [65-85-0]	**Benzoic Acid ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	940 8250
56626 ✔ [93-97-0]	**Benzoic Anhydride pure, 98%**	100 Gms 500 Gms	1980 6930
56518 [119-53-9]	**Benzoin extrapure, 99%**	100 Gms 250 Gms	550 850
51448 [119-61-9]	**Benzophenone extrapure, 99%** (Diphenylketone)	500 Gms 5 Kg	660 6380
42810 [119-61-9]	**Benzophenone ExiPlus, Multi-Compendial, 99%** (Diphenylketone) meets compendial specs of USP, BP, Ph.Eur	500 Gms	770
98101 [95-14-7]	**Benzotriazole pure, 97%**	250 Gms 1 Kg	660 2420
92300 [100-51-6]	**Benzyl Alcohol extrapure, 99%**	500 ml 1000 ml 2500 ml 25 Ltr	680 1300 3000 28000
97005 [100-51-6]	**Benzyl Alcohol ACS, 99%**	500 ml 1000 ml 2500 ml	720 1400 3150
84758 △ [100-51-6]	**Benzyl Alcohol ExiPlus, Multi-Compendial, 99%, aldehyde 500ppm** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	1500 3600
18353 △ [100-51-6]	**Benzyl Alcohol GC-HS, 99.9%**	1000 ml 2500 ml	2200 5000
84199 [120-51-4]	**Benzyl Benzoate pure, 99%**	500 ml 2500 ml	810 3680
28482 [100-44-7]	**Benzyl Chloride extrapure, 99%**	500 ml 2500 ml	550 2180
31041 ✔ [4196-99-0]	**Biebrich Scarlet (C.I. 26905) for microscopy** (Acid Red 66, Ponceau BS, Biebrich Scarlet WS, Biebrich Scarlet Sodium Salt)	25 Gms 100 Gms	1150 3200
41322 [5892-10-4]	**Bismuth (III) Carbonate Basic pure, 80%** (Bismuth Subcarbonate)	100 Gms 500 Gms	2420 10120
51451 [7787-60-2]	**Bismuth (III) Chloride extrapure, 98%** (Bismuth Trichloride)	100 Gms 500 Gms	3910 18400
82797 [99-26-3]	**Bismuth (III) Gallate Basic ExiPlus, Multi-Compendial, 52-57%** (Bismuth Subgallate, Gallic Acid Bismuth Basic Salt) meets compendial specs of USP	100 Gms 500 Gms	2650 12080
80053 [7440-69-9]	**Bismuth Metal Lumps, 99.5%**	100 Gms 500 Gms	1730 6900
24008 [7440-69-9]	**Bismuth Metal Granular ultrapure, 99.99%**	10 Gms 50 Gms 250 Gms	1100 3450 9200
72217 [7440-69-9]	**Bismuth Metal Powder extrapure, 99%, -325 mesh**	100 Gms 500 Gms	2990 13800
29444	**Bismuth AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml 250 ml	2070 3340
72918 [10035-06-0]	**Bismuth (III) Nitrate Pentahydrate extrapure, 98%** (Bismuth Trinitrate Pentahydrate)	100 Gms 500 Gms	1610 7480
19737 [10035-06-0]	**Bismuth (III) Nitrate Pentahydrate ACS, 98%** (Bismuth Trinitrate Pentahydrate)	100 Gms 500 Gms	1840 8340
84869 [10035-06-0]	**Bismuth (III) Nitrate Pentahydrate ExiPlus, Multi-Compendial, 98%** (Bismuth Trinitrate Pentahydrate) meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	2070 8970
89049 [1304-76-3]	**Bismuth (III) Oxide extrapure, 99%** (Dibismuth Trioxide)	100 Gms 500 Gms	3300 14000

Code	Name	Packing	Price
32915 [1304-85-4]	**Bismuth (III) Nitrate Basic extrapure AR, 79% Bi2O3** (Bismuth (III) Subnitrate, Bismuth (III) Oxynitrate)	100 Gms 500 Gms	2760 12650
18473 [7787-68-0]	**Bismuth (III) Sulphate pure, 90%** (Bismuth Trisulphate)	100 Gms 500 Gms	2650 11500
17798 [108-19-0]	**Biuret extrapure AR, 97%** (Allophanic Acid Amide, Carbamoylurea)	25 Gms 100 Gms	1330 4720
18952 [10043-35-3]	**Boric Acid extrapure, 99.5%**	500 Gms 5 Kg	580 5520
80266 [10043-35-3]	**Boric Acid ACS, 99.5%**	500 Gms 5 Kg	660 6240
61871 [10043-35-3]	**Boric Acid ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	700 6540
22311 [10043-35-3]	**Boric Acid for molecular biology, 99.5%**	250 Gms 500 Gms 1 Kg	750 1240 2310
64765	**Boron AAS Standard Solution in H2O (Traceable to NIST)**	100 ml 250 ml	2070 3340
52519 [109-63-7]	**Boron Trifluoride Diethyl Etherate pure, 45-50% BF3** (Boron Trifluoride Ethyl Etherate)	250 ml 500 ml 2500 ml	1270 2070 9780
25529 [1303-86-2]	**Boron Trioxide Anhydrous (Boric Anhydride) extrapure AR, 98.5%** (di-Boron Trioxide, Boric Oxide)	100 Gms 250 Gms 1 Kg	2070 4030 13800
24167 [1303-86-2]	**Boron Trioxide Anhydrous Granular (Boric Anhydride) extrapure, 97%, -60 mesh** (di-Boron Trioxide, Boric Oxide)	100 Gms 250 Gms 500 Gms	1790 3220 5980
82076	**Bouins Picro Formal Fixing Solution (Bouins Fluid & Reagent)**	100 ml 500 ml	280 720
98188 [3844-45-9]	**Brilliant Blue FCF (Erioglaucine Disodium Salt)** (C.I. No. 42090, Acid Blue 9, Alphazurine FG)	25 Gms 100 Gms 500 Gms	580 2070 6900
49617 [633-03-4]	**Brilliant Green** (CI No. 42040)	25 Gms 100 Gms 1 Kg	280 810 6890
93883 [9002-92-0]	**Brij-30® (Polyethylene Glycol Dodecyl Ether, Brij-L4)**	500 ml 1000 ml	1440 2650
89429 [9002-92-0]	**Brij-35® (Main Component) (Polyoxyethylene Lauryl Ether, Brij-L23)**	250 Gms 500 Gms 5 Kg	870 1500 13800
97698 [9002-92-0]	**Brij-35® (30% Aq. solution) (Polyoxyethylene Lauryl Ether, Brij-L23)**	100 ml 500 ml	710 1120
43808 [9004-95-9]	**Brij-58® (Polyethylene Glycol Cetyl Ether, Cetomacrogol 1000, C16H20)**	500 Gms 5 Kg	1380 12370
67687	**Bromine Number Titration Solvent (Dichloromethane (MDC) Formulation)**	1 Ltr 2.5 Ltr	1670 3220
64062 [76-60-8]	**Bromocresol Green (BCG) ACS**	5 Gms 25 Gms 100 Gms 500 Gms	880 3450 11500 51750
52512 [76-60-8]	**Bromocresol Green (BCG) ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	5 Gms 25 Gms 100 Gms 500 Gms	950 3800 12650 55200
11456 [62625-32-5]	**Bromocresol Green Sodium Salt (water soluble) ACS** (3',3",5',5-Tetrabromo-m-Cresolsulfonp hthalein Sodium Salt)	5 Gms 25 Gms	810 3390
17375 [115-40-2]	**Bromocresol Purple**	5 Gms 25 Gms	350 1180
24836 [62625-30-3]	**Bromocresol Purple Sodium Salt for molecular biology** (5,5-Dibromo-o-Cresolsulfonphthalein Sodium Salt)	5 Gms 25 Gms	350 870
55551 [115-40-2]	**Bromocresol Purple ACS, ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	5 Gms 25 Gms	450 1330
77644 [75-25-2]	**Bromoform extrapure, 98%** (Tribromomethane)	250 ml	8050
44669 [75-25-2]	**Bromoform Specially Purified extrapure, 99%** (Tribromomethane)	250 ml	8630
11458 [115-39-9]	**Bromophenol Blue ACS**	5 Gms 25 Gms 100 Gms 500 Gms	260 940 2760 12650
93676 [115-39-9]	**Bromophenol Blue ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	5 Gms 25 Gms 100 Gms	330 1130 2990
42888 [16574-43-9]	**Bromopyrogallol Red extrapure AR** (Dibromopyrogallolsulfonphthaleine)	1 Gms 5 Gms 10 Gms	950 3680 7050
48214 [76-59-5]	**Bromothymol Blue**	5 Gms 25 Gms 100 Gms 500 Gms	410 1330 4720 21850
72752 [76-59-5]	**Bromothymol Blue ACS, ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	5 Gms 25 Gms 100 Gms	460 1440 5060
44798 [76-59-5]	**Bromothymol Blue 0.04% Aq. Solution**	100 ml 500 ml	190 550
45899	**Buffer solution pH 4**	500 ml	450

Code	Name	Packing	Price
97409	**Buffer solution pH 7**	500 ml	450
52031	**Buffer solution pH 9.2**	500 ml	450
74887 [2386-54-1]	**Butane Sulphonic Acid Sodium Salt Anhydrous for HPLC, 99%**	25 Gms 100 Gms 250 Gms	1500 5290 13000
99919 [2386-54-1]	**Butane Sulphonic Acid Sodium Salt Monohydrate for HPLC, 99%**	25 Gms 100 Gms 250 Gms	1500 5290 13000
28118 [123-86-4]	**Butyl Acetate pure, 99%**	500 ml 2500 ml 25 Ltr	410 1610 14950
80870 [123-86-4]	**Butyl Acetate ACS, 99.5%**	500 ml 2500 ml	520 1960
30389 [123-86-4]	**Butyl Acetate ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	600 2250
78487 [123-86-4]	**Butyl Acetate GC-HS, 99.9%**	250 ml 500 ml 1000 ml	810 1440 2530
47364 [141-32-2]	**n-Butyl Acrylate pure, 99%**	500 ml 2500 ml	640 2880
90527 [540-88-5]	**tert-Butyl Acetate (TBAc) extrapure, 99%**	500 ml 2500 ml	1380 5700
19147 [71-36-3]	**n-Butyl Alcohol (1-butanol, n-butanol) pure, 99%**	500 ml 2500 ml 25 Ltr	270 1320 10800
72768 [71-36-3]	**n-Butyl Alcohol (1-Butanol, n-Butanol) ACS, 99.5%**	500 ml 2500 ml	340 1560
15612 [71-36-3]	**n-Butyl Alcohol (1-Butanol, n-Butanol) ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	400 1860
83903 [71-36-3]	**n-Butyl Alcohol (1-butanol, n-butanol) GC-HS, 99.9%**	250 ml 1000 ml	480 1320
43706 [71-36-3]	**n-Butyl Alcohol (1-butanol, n-butanol) for HPLC & UV Spectroscopy, 99.9%**	500 ml 1000 ml	600 1140
95725 [78-92-2]	**sec-Butyl Alcohol pure, 98%** (2-Butanol, sec-Butanol)	500 ml 2500 ml 25 Ltr	540 2100 19620
58308 [78-92-2]	**sec-Butyl Alcohol Reagent Grade, 99%** (2-Butanol, sec-Butanol)	500 ml 2500 ml	600 2370
99336 [78-92-2]	**sec-Butyl Alcohol ExiPlus, Multi-Compendial, 99%** (2-Butanol, sec-Butanol) meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	660 2620
21835 [75-65-0]	**tert-Butyl Alcohol extrapure, 99%** (2-Methyl-2-Propanol Trimethyl Carbinol)	500 ml 2500 ml 25 Ltr	350 1460 13200
70816 [75-65-0]	**tert-Butyl Alcohol ACS, 99%** (2-Methyl-2-Propanol Trimethyl Carbinol)	500 ml 2500 ml	400 1590
41609 [75-65-0]	**Tert-Butyl Alcohol ExiPlus, Multi-Compendial, 99%** (2-Methyl-2-Propanol Trimethyl Carbinol) meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	460 1850
90248 [75-65-0]	**tert-Butyl Alcohol for HPLC, 99.7%**	250 ml 1000 ml	740 1460
76035 [109-73-9]	**n-Butylamine pure, 99%**	500 ml 2500 ml	520 2300
54997 [109-73-9]	**n-Butylamine Reagent Grade, 99%**	500 ml 2500 ml	670 2880
12302 [75-64-9]	**tert-Butylamine pure, 99%**	500 ml 2500 ml 25 Ltr	690 3110 28750
55005 [75-64-9]	**tert-Butylamine Reagent Grade, 99.5%**	500 ml 2500 ml	810 3680
60004 ✔ [25013-16-5]	**Butylated Hydroxyanisole (BHA) ExiPlus, Multi-Compendial, 98%** (2-tert-Butyl-4-Methoxyphenol) meets compendial specs of USP, BP	100 Gms 500 Gms	980 3020
32227 [109-69-3]	**n-Butyl Chloride (1-Chlorobutane) extrapure, 99%**	500 ml 2500 ml	920 4490
95594 [109-69-3]	**n-Butyl Chloride (1-Chlorobutane) for HPLC & UV Spectroscopy, 99.5%**	500 ml 2500 ml	3220 5980
38033 [128-37-0]	**Butylated Hydroxytoluene (BHT) pure, 99%** (2,6-di-tert-Butyl-p-Cresol)	500 Gms 5 Kg	1070 10290
82010 [128-37-0]	**Butylated Hydroxytoluene (BHT) ExiPlus, Multi-Compendial, 99%** (2,6-di-tert-Butyl-p-Cresol) meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	1170 11290
32701 [1948-33-0]	**tert-Butyl Hydroquinone (TBHQ) extrapure AR, 98%**	250 Gms 2.5 Kg	1580 15050
79927 [1634-04-4]	**tert-Butyl Methyl Ether pure, 98%**	500 ml 2500 ml 25 Ltr	400 1610 13800
19465 [1634-04-4]	**tert-Butyl Methyl Ether ACS, 99%**	500 ml 2500 ml	460 1960
65701 [1634-04-4]	**tert-Butyl Methyl Ether ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	560 2360
91087 [1634-04-4]	**tert-Butyl Methyl Ether for HPLC, 99.5%**	1000 ml 2500 ml	1210 2760
63952 [590-01-2]	**Butyl Propionate extrapure, 99%** (Propionic Acid Butyl Ester)	500 ml 2500 ml	1730 8280

Code	Name	Packing	Price
81503 [107-92-6]	**n-Butyric Acid extrapure, C4-99% (GC)**	500 ml	610
		2500 ml	2480
55964 [107-92-6]	**n-Butyric Acid Reagent Grade, C4-99% (GC)**	250 ml	370
		1000 ml	1140
82749 [5743-04-4]	**Cadmium Acetate Dihydrate pure, 98%**	100 Gms	520
		500 Gms	1960
84627 [5743-04-4]	**Cadmium Acetate Dihydrate extrapure AR, 99%**	100 Gms	610
		500 Gms	2280
98469	**Cadmium AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml	2070
		250 ml	3340
25442 [513-78-0]	**Cadmium Carbonate pure, 97%**	100 Gms	690
		500 Gms	2760
15411 [10108-64-2]	**Cadmium Chloride Dried Anhydrous pure, 97%**	100 Gms	1800
		500 Gms	7500
99643 [10108-64-2]	**Cadmium Chloride Dried Anhydrous ACS, 99%**	100 Gms	2000
		500 Gms	8000
48076 [35658-65-2]	**Cadmium Chloride Monohydrate extrapure, 99%**	100 Gms	510
		500 Gms	1800
87355 [7790-80-9]	**Cadmium Iodide pure, 99%**	25 Gms	1730
		100 Gms	5980
		500 Gms	28520
48392 [7440-43-9]	**Cadmium Metal Rods, 99%**	100 Gms	640
85545 [7440-43-9]	**Cadmium Metal Powder extrapure, 99.5%, -200 mesh**	25 Gms	1380
		100 Gms	4600
70434 [10022-68-1]	**Cadmium Nitrate Tetrahydrate ACS, 99%**	500 Gms	1500
29996 [1306-19-0]	**Cadmium Oxide pure, 99%**	100 Gms	500
		250 Gms	1040
		500 Gms	1610
39818 [1306-19-0]	**Cadmium Oxide extrapure AR, 99.7%**	100 Gms	690
		250 Gms	1270
15806 [7790-84-3]	**Cadmium Sulphate 8/3Hydrate pure, 98%**	100 Gms	460
		500 Gms	1730
43019 [7790-84-3]	**Cadmium Sulphate 8/3Hydrate ACS, 98%**	100 Gms	580
		500 Gms	2070
98002 [1306-23-6]	**Cadmium Sulphide Orange extrapure, 99%**	100 Gms	580
		500 Gms	1840
22905 [1306-23-6]	**Cadmium Sulphide Orange ultrapure, 99.995%**	5 Gms	4600
		25 Gms	10350
19267 [58-08-2]	**Caffeine anhydrous pure, 98%**	100 Gms	900
		250 Gms	1730
		500 Gms	2880
		5 Kg	26450
47633 [58-08-2]	**Caffeine Anhydrous ExiPlus, Multi-Compendial, 98.5%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	1040
		250 Gms	1840
		2.5 Kg	14380
44509 [8011-96-9]	**Calamine pure, 98%**	500 Gms	750
		1 Kg	1420
28162 [8011-96-9]	**Calamine ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP	500 Gms	1050
70214 [5743-26-0]	**Calcium Acetate Monohydrate ACS, 99%**	500 Gms	1380
21606 [5743-26-0]	**Calcium Acetate Monohydrate ACS, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP	500 Gms	1500
58056 [114460-21-8]	**Calcium Acetate Hydrate extrapure AR, 99%**	500 Gms	490
		5 Kg	4600
60934 [114460-21-8]	**Calcium Acetate Hydrate ultrapure, 99.9965%**	5 Gms	5750
96493 [9005-35-0]	**Calcium Alginate ex. Brown Algae pure, 15% moisture** (Alginic Acid Calcium Salt)	100 Gms	410
		500 Gms	1500
		5 Kg	12650
36920	**Calcium AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml	2070
		250 ml	3340
67162 [471-34-1]	**Calcium Carbonate extrapure, 98%**	500 Gms	230
		5 Kg	1680
		25 Kg	8050
87293 [471-34-1]	**Calcium Carbonate ACS, 99%**	500 Gms	290
		5 Kg	1950
		25 Kg	8080
37683 [471-34-1]	**Calcium Carbonate ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	580
		5 Kg	4520
70650 [10035-04-8]	**Calcium Chloride Dihydrate ACS, 99%**	500 Gms	420
		5 Kg	3220
		25 Kg	14380
31944 [10035-04-8]	**Calcium Chloride Dihydrate ACS, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	650
		5 Kg	5060
97080 [10035-04-8]	**Calcium Chloride Dihydrate for molecular biology, 99.5%**	100 Gms	360
		500 Gms	1640
56949 [10035-04-8]	**Calcium Chloride Dihydrate for tissue culture, 99.5%**	500 Gms	530
		5 Kg	4830
84336 [10043-52-4]	**Calcium Chloride Fused pure, 90-95%**	500 Gms	250
		5 Kg	1730
77838 [5785-44-4]	**Calcium Citrate Tribasic Tetrahydrate pure, 98%**	500 Gms	460
97388 [5785-44-4]	**Calcium Citrate Tribasic Tetrahydrate ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP	500 Gms	580
57055 [7789-75-5]	**Calcium Fluoride extrapure, 97%** (Fluorite, Fluorospar)	100 Gms	360
		500 Gms	1320
52434 [544-17-2]	**Calcium Formate pure, 98%** (Formic Acid Calcium Salt)	500 Gms	660

Code	Name	Packing	Price
85610 [1305-62-0]	**Calcium Hydroxide ACS, 95%**	500 Gms 5 Kg	610 4780
74471 [1305-62-0]	**Calcium Hydroxide ExiPlus, Multi-Compendial, 95%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	870 6400
31549 [7789-80-2]	**Calcium Iodate Anhydrous (Lautarite) pure, 98%**	100 Gms 250 Gms	3220 6560
96342 [7440-70-2]	**Calcium Metal Granular extrapure, 98%, 10-20 mesh**	100 Gms 500 Gms 5 Kg	410 1380 9200
94212 [10102-68-8]	**Calcium Iodide Tetrahydrate pure, 98%**	10 Gms 25 Gms 100 Gms	1730 3680 10350
56571 [71626-98-7]	**Calcium Iodide Hydrate pure, 98%**	25 Gms 100 Gms	4370 11270
79225 [13477-34-4]	**Calcium Nitrate Tetrahydrate ACS, 99%**	500 Gms 5 Kg	380 2760
16807 [1305-78-8]	**Calcium Oxide Powder pure, 95%**	500 Gms 5 Kg	300 2050
51007 [10031-30-8]	**Calcium Phosphate Monobasic Monohydrate extrapure, 95%** (Calcium Tetrahydrogen Diorthophosphate Monohydrate)	500 Gms 5 Kg	610 5440
93239 [7757-93-9]	**Calcium Phosphate Dibasic Anhydrous extrapure, 98%** (Calcium Hydrogen Phosphate Anhydrous, di-Calcium Phosphate Anhydrous)	500 Gms 5 Kg	430 3030
77793 [7789-77-7]	**Calcium Phosphate Dibasic Dihydrate extrapure, 98%** (Calcium Hydrogen Phosphate Dihydrate, di-Calcium Phosphate Dihydrate)	500 Gms 5 Kg	520 4830
57744 [7789-77-7]	**Dibasic Calcium Phosphate Dihydrate, ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP	500 Gms 5 Kg	800 7250
34878 [7758-87-4]	**Calcium Phosphate Tribasic (TCP) pure, 90%** (Tricalcium Phosphate)	500 Gms 5 Kg	650 5980
59781 [4075-81-4]	**Calcium Propionate pure, 97%**	500 Gms	600
46133 [1592-23-0]	**Calcium Stearate, 6.6-7.4% Ca** (Stearic Acid Calcium Salt)	500 Gms 5 Kg	690 5870
94044 [10101-41-4]	**Calcium Sulphate Dihydrate pure, 98%**	500 Gms 5 Kg	380 3390
69095 [10101-41-4]	**Calcium Sulphate Dihydrate ACS, 98%**	500 Gms 5 Kg	550 4960
64813 [76-22-2]	**Camphor pure, 95%**	100 Gms 500 Gms	650 3140
23764 [334-48-5]	**Capric Acid (n-Decanoic Acid) pure, C10-99%** (n-Capric Acid, n-Decoic Acid)	250 ml 500 ml 2500 ml	810 1500 6900
51840 [110-42-9]	**Capric Acid Methyl Ester (Methyl Decanoate) extrapure, 99%** (Methyl Caprate)	250 ml 1000 ml	2100 6710
89028 [142-62-1]	**Caproic Acid extrapure, C6-99%** (n-Hexanoic Acid)	100 ml 500 ml 2500 ml	350 1100 4600
88889 [124-07-2]	**Caprylic Acid pure, C8-99% (GC)**	500 ml 2500 ml	1500 6330
15343 [67-56-1]	**Carbinol pure, 99%**	500 ml 2500 ml 25 Ltr	290 780 6620
82911 [67-56-1]	**Carbinol extrapure AR, 99.8%**	500 ml 1000 ml 2500 ml	320 600 840
39058 [4197-24-4]	**Carbol Fuchsin**	25 Gms 100 Gms 500 Gms	440 1380 6210
93049 [9003-01-4]	**Carbopol 934 (Carboxy Vinyl Polymer 934) extrapure**	500 Gms 1 Kg 5 Kg	1440 2760 12650
51159 [9003-01-4]	**Carbopol 940 (Carboxy Vinyl Polymer 940) extrapure**	500 Gms 1 Kg 5 Kg	1440 2760 12650
93797 [9004-32-4]	**Carboxymethyl Cellulose Sodium Salt (CMC Sodium LV), Low Viscosity 250-350CPS** (Carmellose Sodium)	500 Gms 5 Kg	750 6500
59938 [9004-32-4]	**Carboxymethyl Cellulose Sodium Salt (CMC Sodium MV), Medium Viscosity 400-800 CPS** (Carmellose Sodium)	500 Gms 5 Kg	800 7500
35298 [9004-32-4]	**Carboxymethyl Cellulose Sodium Salt (CMC Sodium HV), High Viscosity 1100-1900 CPS** (Carmellose Sodium)	500 Gms 5 Kg	980 9200
97048 [1390-65-4]	**Carmine** (C.I. No. 75470)	5 Gms 25 Gms 100 Gms	1500 5500 20000
66993	**Carnoy`s Solution (Fixative)** (Carnoy`s Fluid (Fixative), Carnoy Fluid (Fixative))	100 ml 500 ml	370 1700
32385 [11114-20-8]	**Carrageenan (Kappa) Refined ex. Eucheuma Cottonii** (k-Carrageenan, Irish Moss)	25 Gms 100 Gms 500 Gms	1300 2600 8100
16145 [8001-79-4]	**Castor Oil**	500 ml 2500 ml	680 3200
26080 [120-80-9]	**Catechol (Pyrocatechol) extrapure, 99%**	100 Gms 500 Gms	450 900
31018 [8000-27-9]	**Cedarwood Oil Natural for microscopy** (Immersion Oil)	30 ml 125 ml 500 ml	350 900 1900
60705 [1562-90-9]	**Celestine Blue B, 98%-102%** (Mordant Blue 14, C.I. 51050)	5 Gms 25 Gms	4300 20500

Code	Name	Packing	Price
82523 [91053-39-3]	**Celite 545, size 20-40 micrometer**	500 Gms 1 Kg 5 Kg	450 800 3700
61767 ■ [477243-06-4]	**Ceramide ex. Chicken Egg 99%** (C16 Ceramide, N-(Hexadecanoyl) -Sphing-4-Enine), n-Hexadecanoyl-d-Erythro-Sphingosin e	5 Mg 25 Mg	5600 28000
23505 [1306-38-3]	**Ceric Oxide extrapure AR, 99.9%** (Cerium (IV) Oxide)	100 Gms 500 Gms	500 2000
92425 [10294-42-5]	**Ceric Sulphate Tetrahydrate extrapure AR, 80%** (Ceric (IV) Sulfate Tetrahydrate, Cerium Disulfate Tetrahydrate)	25 Gms 100 Gms 500 Gms	550 1300 4200
10731 [10294-42-5]	**Ceric Sulphate Tetrahydrate specially purified, 98%** (Ceric (IV) Sulfate Tetrahydrate, Cerium Disulfate Tetrahydrate)	25 Gms 100 Gms	7800 24000
61414 [13590-82-4]	**Ceric Sulphate 0.1N Solution**	500 ml	620
66676 [18618-55-8]	**Cerous Chloride Heptahydrate extrapure AR, 99%** (Cerium (III) Chloride Heptahydrate)	25 Gms 100 Gms 500 Gms	690 1700 6000
57320 [10294-41-4]	**Cerous Nitrate Hexahydrate extrapure AR, 99.9%** (Cerium (III) Nitrate Hexahydrate, Cerium Trinitrate Hexahydrate)	100 Gms 500 Gms	1020 4000
81586 [1119-97-7]	**Cetrimide (Tetradecyltrimethyl Ammonium Bromide) extrapure AR, 99%** (Tetradecyltrimethyl Ammonium Bromide)	100 Gms 500 Gms 5 Kg	740 2070 19550
33001 [8005-44-5]	**Cetostearyl Alcohol pure**	500 Gms 5 Kg	950 9200
40629 ✔ [36653-82-4]	**Cetyl Alcohol pure, 98%** (1-Hexadecanol)	500 Gms 5 Kg	750 6670
83438 [7440-44-0]	**Charcoal Activated (280) ExiPlus, Multi-Compendial** (Activated Carbon) meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms 5 Kg	280 590 5020
48088 [7440-44-0]	**Charcoal Activated (950) Granular extrapure** (Activated Carbon)	500 Gms 5 Kg	640 5650
64905 [7080-50-4]	**Chloramine-T Trihydrate extrapure AR, 98%** (N-Chloro-p-Toluenesulphonamide Sodium Salt)	100 Gms 250 Gms 500 Gms 2.5 Kg	520 1040 2020 9200
29611 [95-51-2]	**o-Chloroaniline extrapure, 99%** (2-Chloroaniline)	500 ml 2500 ml	820 3910
89370 [108-90-7]	**Chlorobenzene (MCB) extrapure, 99%** (Monochlorobenzene)	500 ml 2500 ml 25 Ltr	350 1500 13800
46364 [108-90-7]	**Chlorobenzene (MCB) ACS, 99.5%** (Monochlorobenzene)	500 ml 2500 ml	460 1730
10124 [108-90-7]	**Chlorobenzene (MCB) ExiPlus, Multi-Compendial, 99.5%** (Monochlorobenzene) meets compendial specs of BP	500 ml 2500 ml	580 2070
35680 [3114-55-4]	**Chlorobenzene-d5 for NMR spectroscopy, 99 Atom %D**	1 ml 5 ml	8000 24000
33023 [59-50-7]	**p-Chloro-m-Cresol (PCMC) pure, 98%**	250 Gms 500 Gms	1300 2100
80908 [88-04-0]	**4-Chloro-3,5-Dimethylphenol (PCMX) extrapure, 99%** (4-Chloro-3,5-Xylenol, Chloroxylenol, p-Chloro-m-Xylenol)	100 Gms 500 Gms	880 3060
96712 [67-66-3]	**Chloroform extrapure, 99%**	500 ml 2500 ml 25 Ltr	300 980 7500
84155 [67-66-3]	**Chloroform ACS, 99.8%**	500 ml 2500 ml	340 1100
32022 [67-66-3]	**Chloroform ExiPlus, Multi-Compendial, 99.8%** meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	400 1200
96764 [67-66-3]	**Chloroform for molecular biology, 99.8%**	100 ml 500 ml	450 1800
30792 [67-66-3]	**Chloroform for HPLC & UV Spectroscopy, 99.8%**	500 ml 1000 ml 2500 ml	450 650 1400
14423 [67-66-3]	**Chloroform GC-HS, 99.9%**	250 ml 500 ml 1000 ml	550 900 1500
51126 [67-66-3]	**Chloroform Dried, 99.5%, water 0.005%**	500 ml 1000 ml 2500 ml	500 970 2180
90603 🌡 [865-49-6]	**Chloroform-d (CDCl3) (w/o TMS, Stab w/ Ag) for NMR spectroscopy, 99.8 Atom %D**	100 ml	10200
20996 🌡 [865-49-6]	**Chloroform-d (CDCl3) (with 1% TMS, Stab w/ Ag) for NMR spectroscopy, 99.8 Atom %D**	100 ml	10900
57034 🌡 [865-49-6]	**Chloroform-d (CDCl3) (with 0.03% TMS, Stab w/ Ag) for NMR spectroscopy,99.8 Atom %D**	100 ml	10500
22288 🌡 [865-49-6]	**Chloroform-d (CDCl3) (w/o TMS, Stab w/ Ag) for NMR spectroscopy, 99.8 Atom %D**	100 Gms	8700
98718 🌡 [865-49-6]	**Chloroform-d (CDCl3) (with 0.03% TMS, Stab w/ Ag) for NMR spectroscopy,99.8 Atom %D**	100 Gms	9200
72017 🌡 [865-49-6]	**Chloroform-d (CDCl3) (with 1% TMS, Stab w/ Ag) for NMR spectroscopy, 99.8 Atom %D**	100 Gms	9700
39657 [937-14-4]	**3-Chloroperoxybenzoic Acid (mCPBA) pure, 55-75%**	25 Gms 100 Gms 500 Gms	425 1300 6250
59271 [128-09-6]	**N-Chlorosuccinimide extrapure, 98%**	250 Gms 1 Kg	750 2880

Code	Name	Packing	Price
47048 [1667-99-8]	**Chromeazurol S (Chromazurol S)** (Mordant Blue 29)	10 Gms 25 Gms	1100 2000
86183	**Chromium AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml 250 ml	1800 2900
72232 [10060-12-5]	**Chromium (III) Chloride Hexahydrate extrapure AR, 98%** (Chromic Chloride, Chromium Trichloride)	500 Gms 5 Kg	700 6500
74767 [7789-02-8]	**Chromium (III) Nitrate Nonahydrate pure, 96%** (Chromic Nitrate Nonahydrate)	500 Gms 5 Kg	1200 11000
68501 △ [14639-25-9]	**Chromium (III) Picolinate pure, 98%** (Chromium Picolinate)	100 Gms	2250
21295 [39380-78-4]	**Chromium (III) Sulphate Basic extrapure, 22.5% Cr** (Chromosal B, Basic Chromic Sulphate)	500 Gms 2.5 Kg	380 1750
57555 [7440-47-3]	**Chromium Powder, 99% -100 mesh**	100 Gms 500 Gms	900 3500
21817 [1333-82-0]	**Chromium (VI) Oxide (Chromium Trioxide) pure, 97%** (Chromic Acid, Chromic Anhydride)	500 Gms 1 Kg	850 1600
61365 [1333-82-0]	**Chromium (VI) Oxide (Chromium Trioxide) ExiPlus, Multi-Compendial , 98%** (Chromic Acid, Chromic Anhydride) meets compendial specs of USP, BP, Ph.Eur	500 Gms 1 Kg	900 1750
36505 [4438-16-8]	**Chrysoidine R** (Basic Orange 1, Methylchrysoidine) (C.I. No. 11320)	25 Gms 100 Gms	210 590
31126 [532-82-1]	**Chrysoidine Y** (4-Phenylazo-m-Phenylenediamine Monohydrochloride, Basic Orange 2, Chrysoidine G)(C.I. No. 11270)	25 Gms 100 Gms	210 590
55990 ⓘ [104-55-2]	**Cinnamaldehyde pure, 98%** (Cinnamic Aldehyde, Cinnamal)	500 ml 2500 ml	1150 6330
29955 [140-10-3]	**Cinnamic Acid extrapure, 99%**	250 Gms 1 Kg	1010 3630
23277 [140-10-3]	**Cinnamic Acid extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	250 Gms 1 Kg	1610 4030
78213 [77-92-9]	**Citric Acid Anhydrous extrapure, 99%**	500 Gms 5 Kg 25 Kg	300 2400 11000
24670 [77-92-9]	**Citric Acid Anhydrous ACS, 99.5%**	5 Kg 500 Gms	2550 320
36890 [77-92-9]	**Citric Acid Anhydrous for tissue culture, 99%**	500 Gms 5 Kg	350 2800
16842 [77-92-9]	**Citric Acid Anhydrous ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	340 2700
49897 [5949-29-1]	**Citric Acid Monohydrate extrapure, 99%**	500 Gms 5 Kg 25 Kg	280 2550 12000
43667 [5949-29-1]	**Citric Acid Monohydrate ACS, 99%**	500 Gms 5 Kg 25 Kg	300 2800 13000
53884 [5949-29-1]	**Citric Acid Monohydrate ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	350 3000
66112	**Cobalt AAS Standard Solution in 0.5N HCl (Traceable to NIST)**	100 ml 250 ml	2400 4100
49393 [6147-53-1]	**Cobalt (II) Acetate Tetrahydrate ACS, 98%** (Cobaltous Acetate Tetrahydrate)	100 Gms 250 Gms 500 Gms	1100 2400 4000
75141 [6147-53-1]	**Cobalt (II) Acetate Tetrahydrate ExiPlus, Multi-Compendial, 98%** (Cobaltous Acetate Tetrahydrate) meets compendial specs of USP, BP	100 Gms 250 Gms 500 Gms	1200 2700 4600
37917 [12602-23-2]	**Cobalt (II) Carbonate Basic pure, 45-50% Co** (Cobaltous Carbonate)	100 Gms 250 Gms 500 Gms	1600 3600 6600
99412 [7791-13-1]	**Cobalt (II) Chloride Hexahydrate pure, 97%** (Cobaltous Chloride Hexahydrate)	100 Gms 500 Gms	1300 6050
36018 [7791-13-1]	**Cobalt (II) Chloride Hexahydrate ACS, 98%** (Cobaltous Chloride Hexahydrate)	100 Gms 500 Gms	1410 6380
88282 [7791-13-1]	**Cobalt (II) Chloride Hexahydrate ExiPlus, Multi-Compendial, 98%** (Cobaltous Chloride Hexahydrate) meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	1520 6590
14386 [7440-48-4]	**Cobalt Metal Powder, 99.5%, -325mesh**	100 Gms 250 Gms	3100 6600
14941 [10026-22-9]	**Cobalt (II) Nitrate Hexahydrate extrapure, 97%** (Cobaltous Nitrate Hexahydrate)	100 Gms 500 Gms	900 3300
85942 [10026-22-9]	**Cobalt (II) Nitrate Hexahydrate ACS, 98%** (Cobaltous Nitrate Hexahydrate)	100 Gms 250 Gms 500 Gms	1100 2500 4100
21722 [1308-06-1]	**Cobalt (II,III) Oxide pure, 99%** (Tricobalt Tetraoxide)	100 Gms 250 Gms 500 Gms	2800 4000 7000
25710 [10026-24-1]	**Cobalt (II) Sulphate Heptahydrate pure, 98%** (Cobaltous Sulphate Heptahydrate)	100 Gms 500 Gms	800 3300
42752 [10026-24-1]	**Cobalt (II) Sulphate Heptahydrate extrapure AR, 99%** (Cobaltous Sulphate Heptahydrate)	100 Gms 250 Gms 500 Gms	900 2100 4100
36672	**Collodion Flexible pure**	500 ml	2020

Code	Name	Packing	Price
57380	**Copper AAS Standard Solution in 1N HNO3 (Traceable to NIST)**	100 ml 250 ml	1800 2900
72033 [7440-50-8]	**Copper Metal powder extrapure, 99.5%, -200 mesh**	500 Gms	3200
42577 [7440-50-8]	**Copper Metal Turnings, 99.5%, 40-100 mesh**	100 Gms 500 Gms	750 3100
28007 [95-48-7]	**o-Cresol pure, 98%**	500 ml 2500 ml	810 3220
77744 [108-39-4]	**m-Cresol pure, 98%** (3-Methylphenol)	500 ml 2500 ml	1560 6900
48620 [108-39-4]	**m-Cresol extrapure AR, 99%** (3-Methylphenol)	500 ml	1960
63997 [106-44-5]	**p-Cresol pure, 98.5%**	500 ml 2500 ml	1270 5640
89599 [1733-12-6]	**Cresol Red**	5 Gms 25 Gms 100 Gms	350 810 2880
21442 [10510-54-0]	**Cresyl Violet Acetate, 70%**	1 Gms 5 Gms	920 3800
89117 △ [123-73-9]	**Crotonaldehyde extrapure, 99%** (Crotonic Aldehyde, ß-Methacrolein)	500 ml 2500 ml	1270 5980
28376 [548-62-9]	**Crystal Violet (Gentian Violet)** (Methyl Violet 10B, C.I. No. 42555)	25 Gms 100 Gms 500 Gms	920 3220 15180
17610 [548-62-9]	**Crystal Violet (Gentian Violet) ExiPlus, Multi-Compendial** (Methyl Violet 10B, C.I. No. 42555) meets compendial specs of BP, Ph.Eur	25 Gms 100 Gms 500 Gms	1040 3680 17250
99062 [548-62-9]	**Crystal Violet (Gentian Violet) PrimeDx**	1 Kg 5 Kg	36800 172500
94666 ✔ [17455-13-9]	**18-Crown-6-Ether pure, 98%**	25 Gms 100 Gms 500 Gms	750 2400 11000
74472 [8014-13-9]	**Cumin Seed Oil extrapure, aldehyde 28%**	100 ml 500 ml	2070 7590
87431 [6046-93-1]	**Cupric Acetate Monohydrate pure, 98%** (Copper (II) Acetate Monohydrate)	100 Gms 500 Gms	600 2400
67912 [6046-93-1]	**Cupric Acetate Monohydrate ACS, 98%** (Copper (II) Acetate Monohydrate)	250 Gms	1800
57060 [6046-93-1]	**Cupric Acetate Monohydrate ExiPlus, Multi-Compendial, 99%** (Copper (II) Acetate Monohydrate) meets compendial specs of BP, Ph.Eur	250 Gms	2100
11204 [7789-45-9]	**Cupric Bromide pure, 99%** (Copper (II) Bromide)	100 Gms 500 Gms	1600 6500
13945 [12069-69-1]	**Cupric Carbonate Basic pure, 95%** (Copper (II) Carbonate Basic, Copper (II) Hydroxide Carbonate)	500 Gms 5 Kg	2200 21000
92315 [7447-39-4]	**Cupric Chloride pure, 98%** (Copper (II) Chloride)	500 Gms 1 Kg	1850 3500
94809 [10125-13-0]	**Cupric Chloride Dihydrate ACS, 99%**	500 Gms	1450
89797 [10125-13-0]	**Cupric Chloride Dihydrate ExiPlus, Multi-Compendial, 99%** (Copper (II) Chloride Dihydrate) meets compendial specs of USP, BP, Ph.Eur	500 Gms 1 Kg	1600 3000
92272 [10125-13-0]	**Cupric Chloride 0.3M Solution**	500 ml 1000 ml	700 1250
87367 [10031-43-3]	**Cupric Nitrate Trihydrate extrapure, 98%** (Copper (II) Nitrate Trihydrate)	250 Gms 500 Gms	800 1400
46704 [10031-43-3]	**Cupric Nitrate Trihydrate ACS, 98%** (Copper (II) Nitrate Trihydrate)	500 Gms	1600
51120 [1317-38-0]	**Cupric Oxide pure, 79% Cu** (Copper (II) Oxide)	100 Gms 500 Gms 5 Kg	700 2900 26400
35119 [7758-98-7]	**Cupric Sulphate Anhydrous extrapure, 99%** (Copper (II) Sulphate Anhydrous)	100 Gms 500 Gms	750 2840
61103 [7758-98-7]	**Cupric Sulphate Anhydrous extrapure AR, 99%** (Copper (II) Sulphate Anhydrous)	100 Gms 500 Gms	950 3110
58352 [7758-98-7]	**Cupric Sulphate Anhydrous ExiPlus, Multi-Compendial, 99%** (Copper (II) Sulphate Anhydrous) meets compendial specs of USP, Ph.Eur	100 Gms 500 Gms	1620 3510
38869 [7758-99-8]	**Cupric Sulphate Pentahydrate pure, 99%** (Copper (II) Sulphate Pentahydrate)	500 Gms 5 Kg	1100 10400
61298 [7758-99-8]	**Cupric Sulphate Pentahydrate ACS, 99.5%** (Copper (II) Sulphate Pentahydrate)	500 Gms 5 Kg	1200 11400
71244 [7758-99-8]	**Cupric Sulphate Pentahydrate ExiPlus, Multi-Compendial, 99%** (Copper (II) Sulphate Pentahydrate) meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	1300 12400
90445 [7758-99-8]	**Cupric Sulphate 0.1M Solution**	1000 ml	610
29497 [7787-70-4]	**Cuprous Bromide pure, 98.5%** (Copper (I) Bromide)	100 Gms 500 Gms	920 4140
26375 [7758-89-6]	**Cuprous Chloride pure, 97%** (Copper (I) Chloride , Copper Monochloride)	500 Gms 5 Kg	1800 17000
21136 [7758-89-6]	**Cuprous Chloride ACS, 97%** (Copper (I) Chloride, Copper Monochloride)	250 Gms	1400
93095 [7681-65-4]	**Cuprous Iodide extrapure AR, 99%** (Copper (I) Iodide)	25 Gms 100 Gms 500 Gms	1040 3450 16330
83086 [1317-39-1]	**Cuprous Oxide pure, 86% Cu** (Copper (I) Oxide)	100 Gms 500 Gms	750 3000
97461 [458-37-7]	**Curcumin (Turmeric Yellow)**	5 Gms 25 Gms 100 Gms	410 1130 3740

Code	Name	Packing	Price
20259 [110-82-7]	**Cyclohexane extrapure, 99%**	500 ml	300
		2500 ml	1130
		25 Ltr	10580
46616 [110-82-7]	**Cyclohexane ACS, 99.5%**	500 ml	350
		2500 ml	1250
23292 [110-82-7]	**Cyclohexane GC-HS, 99.9%**	250 ml	580
		1000 ml	1040
46300 [110-82-7]	**Cyclohexane for HPLC & UV Spectroscopy, 99.9%**	500 ml	520
		1000 ml	810
		2500 ml	1730
92604 [110-82-7]	**Cyclohexane Dried, 99.5%, water 0.005%**	500 ml	460
		1000 ml	870
34781 [110-82-7]	**Cyclohexane ACS, ExiPlus, Multi-Compendial, 99.5%**	500 ml	400
		2500 ml	1330
	meets compendial specs of USP, BP, Ph.Eur		
76191 [1735-17-7]	**Cyclohexane-d12 for NMR spectroscopy, 99.5 Atom %D** (Dodecadeuterocyclohexane)	5 ml	17000
31726 [108-94-1]	**Cyclohexanone extrapure, 99%**	500 ml	410
		2500 ml	1840
89635 [108-94-1]	**Cyclohexanone ACS, 99%**	500 ml	490
		2500 ml	2300
12183 [108-91-8]	**Cyclohexylamine pure, 99%** (Aminocyclohexane)	500 ml	460
		2500 ml	2070
45658 [108-91-8]	**Cyclohexylamine Reagent Grade, 99.5%** (Aminocyclohexane)	500 ml	500
		2500 ml	2420
85606 [91771-62-9]	**Cypriol Oil extrapure, 99%** (Cyperus Root Oil (Cyperus Scariosus), Nagarmotha Oil)	100 ml	8000
95322 [52-90-4]	**L-Cysteine extrapure CHR, 99%**	25 Gms	350
		100 Gms	1100
		500 Gms	4800
		5 Kg	44000
25961 [52-90-4]	**L-Cysteine for tissue culture, 99%**	25 Gms	450
		100 Gms	1250
		500 Gms	5200
69029 [7048-04-6]	**L-Cysteine Hydrochloride Monohydrate ExiPlus, Multi-Compendial, 99%**	100 Gms	950
		500 Gms	3570
		1 Kg	6250
		5 Kg	30000
	meets compendial specs of USP		
35751 [56-89-3]	**L-Cystine extrapure CHR, 99%**	25 Gms	450
		100 Gms	1100
		1 Kg	8500
		5 Kg	40000
83172 [56-89-3]	**L-Cystine for tissue culture, 99%**	25 Gms	520
		100 Gms	1190
		500 Gms	5340
38779 [56-89-3]	**L-Cystine ExiPlus, Multi-Compendial, 99%**	250 Gms	3000
		1 Kg	10200
		5 Kg	48000
	meets compendial specs of USP, BP, Ph.Eur		
10323 [71-30-7]	**Cytosine ExiPlus, Multi-Compendial, 99%**	25 Gms	700
		100 Gms	2500
		500 Gms	8000
	meets compendial specs of BP, Ph.Eur		
56096 [13419-61-9]	**Decane Sulphonic Acid Sodium Salt Anhydrous for HPLC, 99%**	25 Gms	1300
		100 Gms	4600
		250 Gms	11300
71801 [13419-61-9]	**Decane Sulphonic Acid Sodium Salt Monohydrate for HPLC, 99%**	25 Gms	1300
		100 Gms	4600
		250 Gms	11300
93995 [8049-11-4]	**Devardas Alloy (45% Al, 50% Cu, 5% Zn) extrapure**	100 Gms	590
		500 Gms	2360
37163 [8049-11-4]	**Devardas Alloy (45% Al, 50% Cu, 5% Zn) ACS**	100 Gms	790
		250 Gms	1690
61500 [9004-53-9]	**Dextrin (White) practical**	500 Gms	300
		5 Kg	2800
		25 Kg	9000
42738 [50-99-7]	**Dextrose extrapure** (D-Glucose)	500 Gms	250
		1 Kg	485
		5 Kg	2000
		25 Kg	8800
51758 [50-99-7]	**Dextrose ACS** (D-Glucose)	500 Gms	310
		1 Kg	600
		5 Kg	2850
		25 Kg	12370
16191 [14431-43-7]	**Dextrose Monohydrate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%** (D-Glucose Monohydrate)	500 Gms	270
		5 Kg	2490
		25 Kg	10890
	meets compendial specs of USP, BP		
40799 [3012-65-5]	**Diammonium Hydrogen Citrate Anhydrous (Ammonium Citrate Dibasic) ACS, 98%**	500 Gms	870
		5 Kg	8200
81434 [84-74-2]	**Dibutyl Phthalate pure, 98%**	500 ml	580
		2500 ml	2300
12088 [79-43-6]	**Dichloroacetic Acid pure, 99%**	100 ml	2720
		500 ml	9780
		2500 ml	46000
52153 [95-50-1]	**1,2-Dichlorobenzene pure, 99%** (o-Dichlorobenzene)	500 ml	460
		1000 ml	870
		2500 ml	2020
61012 [95-50-1]	**1,2-Dichlorobenzene Reagent Grade, 99%** (o-Dichlorobenzene)	250 ml	370
		2500 ml	2300
79510 [107-06-2]	**1,2-Dichloroethane extrapure, 99%** (Ethylene Dichloride)	500 ml	340
		2500 ml	1010
		25 Ltr	9500

Code	Name	Packing	Price
69566 [107-06-2]	**1,2-Dichloroethane ACS, 99%** (Ethylene Dichloride)	500 ml 2500 ml	370 1150
87988 [107-06-2]	**1,2-Dichloroethane for HPLC & UV Spectroscopy, 99.9%** (Ethylene Dichloride)	500 ml 1000 ml	700 1240
63462 [107-06-2]	**1,2-Dichloroethane Dried, 99.5%, water 0.005%** (Ethylene Dichloride)	1000 ml	1070
79126 [76-54-0]	**2,7-Dichlorofluorescein (DCF) ACS**	5 Gms 25 Gms	600 2500
56825 [75-09-2]	ⓘ **Dichloromethane (MDC) pure, 99%** (Methylene Chloride, DCM)	500 ml 2500 ml 25 Ltr	280 1150 10000
24532 [75-09-2]	ⓘ **Dichloromethane (MDC) ACS, 99.5%** (Methylene Chloride)	500 ml 2500 ml 25 Ltr	320 1250 11000
45368 [75-09-2]	ⓘ **Dichloromethane (MDC) ExiPlus, Multi-Compendial, 99%** (Methylene Chloride) meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	350 1350
29103 [75-09-2]	ⓘ **Dichloromethane (MDC) GC-HS, 99.9%** (Methylene Chloride)	250 ml 500 ml 1000 ml	430 700 1100
58628 [75-09-2]	ⓘ **Dichloromethane (MDC) for HPLC & UV Spectroscopy, 99.9%**	500 ml 1000 ml 2500 ml	410 600 1300
43551 [75-09-2]	ⓘ **Dichloromethane (MDC) Dried, 99.5%, water 0.005%** (Methylene Chloride)	500 ml 1000 ml	340 680
41512 [75-09-2]	ⓘ **Dichloromethane (MDC) for molecular biology, 99.9%** (Methylene Chloride)	250 ml 500 ml	470 800
93075 [1665-00-5]	**Dichloromethane-d2 (MDC-d2) for NMR spectroscopy, 99.5 Atom %D** (Methylene Chloride)	10 ml	21000
55287 [620-45-1]	**2,6-Dichlorophenol Indophenol Sodium ACS, 98%**	1 Gms 5 Gms 25 Gms	310 1200 4910
35551 [620-45-1]	**2,6-Dichlorophenol Indophenol Sodium Salt ExiPlus, Multi-Compendial, 98%** meets compendial specs of BP, Ph.Eur	1 Gms 5 Gms 25 Gms	330 1340 5240
31499 [94-75-7]	**2,4-Dichlorophenoxyacetic Acid (2,4-D) for tissue culture, 98%**	25 Gms 100 Gms 500 Gms	680 1250 5060
77833 [461-58-5]	**Dicyandiamide pure, 98%** (1-Cyanoguanidine, Dicyanodiamide)	500 Gms	750
61335 [77-73-6]	ⓘ **Dicyclopentadiene (Stabilized with BHT) extrapure, 93%**	100 ml 500 ml 2500 ml	240 700 3400
71164 [111-42-2]	**Diethanolamine pure, 98%**	500 ml 2500 ml	350 1400
66443 [111-42-2]	**Diethanolamine ACS, 98.5%**	500 ml 2500 ml	400 1600
47001 [111-42-2]	**Diethanolamine ExiPlus, Multi-Compendial, 98.5%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	450 1800
59161	**5M Diethanolamine Substrate Buffer(5M DEA w/ 2.5mM MgCl2) for phosphatase**	100 ml	2250
31255 [109-89-7]	**Diethylamine pure, 99%**	500 ml 2500 ml 25 Ltr	350 1500 13500
64912 [109-89-7]	**Diethylamine ACS, 99%**	500 ml 2500 ml	380 1800
28298 [109-89-7]	**Diethylamine for HPLC, 99.5%**	1000 ml	1600
77141 [100-37-8]	**2-Diethylaminoethanol (DEAE) extrapure, 99%**	500 ml 2500 ml	650 3140
59351 [91-66-7]	**N,N-Diethylaniline pure, 98%**	500 ml 2500 ml	900 4150
44369 [91-66-7]	**N,N-Diethylaniline Reagent Grade, 99%**	500 ml 2500 ml	1120 5160
36303 [105-58-8]	**Diethyl Carbonate extrapure, 99%**	500 ml 2500 ml	560 2580
83040 [111-46-6]	**Diethylene Glycol (Digol) pure, 98%**	500 ml 2500 ml 25 Ltr	300 1010 9190
68811 [111-46-6]	**Diethylene Glycol (Digol) Reagent Grade, 99%**	500 ml 2500 ml	360 1120
93616 [112-36-7]	**Diethylene Glycol Diethyl Ether (Ethyl Diglyme, DEGDEE) pure, 99%** (2-Ethoxy Ethyl Ether)	100 ml 500 ml	840 3590
44284 [111-96-6]	**Diethylene Glycol Dimethyl Ether (Diglyme, DEGDME) pure, 99%** (2-Methoxy Ethyl Ether)	500 ml 2500 ml 25 Ltr	1010 4480 42560
94361 [111-96-6]	**Diethylene Glycol Dimethyl Ether (Diglyme, DEGDME) Dried, 99%, water 0.05%** (2-Methoxy Ethyl Ether)	500 ml	1400
20882 [111-90-0]	**Diethylene Glycol Monoethyl Ether (DGMEE) pure, 98%** (Carbitol, Ethyldiglycol)	500 ml 2500 ml	500 2130
53422 [96-22-0]	**Diethyl Ketone (3-Pentanone) extrapure, 99%**	500 ml	1350
98389 [105-53-3]	**Diethyl Malonate extrapure, 99%**	500 ml 2500 ml	590 2690
39378 [105-53-3]	**Diethyl Malonate Reagent Grade, 99.5%**	500 ml 2500 ml	680 3140
61252 [95-92-1]	**Diethyl Oxalate pure, 99%**	500 ml 2500 ml	680 3030

Code	Name	Packing	Price
75995 [6283-63-2]	**N,N-Diethyl-p-Phenylenediamine Sulphate Salt (DEPPDA Sulphate) extrapure, 98%** (4-Amino-N,N-Diethylaniline Sulfate Salt)	25 Gms 100 Gms	2200 8400
16150 [84-66-2]	**Diethyl Phthalate pure, 99%**	500 ml 2500 ml	410 1910
78652 [64-67-5]	**Diethyl Sulphate extrapure, 99%**	500 ml 2500 ml	450 2020
15979 [75-11-6]	△ **Diiodomethane (Stabilized with Copper wire) extrapure, 98%** (Methylene iodide) For Separation of minerals	25 Gms 100 Gms 250 Gms	2800 6050 13440
77205 [108-18-9]	**Diisopropylamine pure, 98%**	500 ml 2500 ml 25 Ltr	1010 4710 44800
24465 [110-97-4]	ⓘ **Diisopropanolamine (DIPOA) extrapure AR, 98%** (Bis(2-hydroxypropyl)amine,1,1'-Iminodi-2-propanol, Iminodipropanol)	500 Gms 1 Kg	900 1550
67519 [108-20-3]	**Diisopropylether pure, 98%**	500 ml 2500 ml 25 Ltr	480 2020 18480
63282 [108-20-3]	**Diisopropylether ACS, 99%**	500 ml 2500 ml	520 2360
89299 [108-20-3]	**Diisopropylether Dried, 99%, water 0.0025%**	500 ml 1000 ml	900 1680
91837 [7087-68-5]	**N,N-Diisopropylethylamine (DIPEA) extrapure, 98%** (Ethylene Diisopropylamine, N,N-Ethyldiisopropylamine, DIEA)	100 ml 250 ml 1000 ml	400 900 3430
61125 [110-71-4]	**1,2-Dimethoxyethane (Monoglyme, EGDME) extrapure, 99%** (Ethylene Glycol Dimethyl Ether)	500 ml 2500 ml 25 Ltr	690 3190 30400
94378 [110-71-4]	**1,2-Dimethoxyethane (Monoglyme, EGDME) extrapure AR , 99.5%** (Ethylene Glycol Dimethyl Ether)	500 ml 2500 ml	780 3520
80675 [127-19-5]	**N,N-Dimethylacetamide (DMAc) extrapure, 99%**	500 ml 2500 ml 25 Ltr	380 1630 15000
60991 [127-19-5]	**N,N-Dimethylacetamide (DMAc) Reagent Grade, 99.5%**	500 ml 2500 ml	450 1880
79220 [127-19-5]	**N,N-Dimethylacetamide (DMAc) ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	500 2130
70926 [127-19-5]	**N,N-Dimethylacetamide (DMAc) GC-HS, 99.9%**	250 ml 500 ml 1000 ml	500 880 1630
98189 [127-19-5]	**N,N-Dimethylacetamide (DMAc) for HPLC, 99.8%**	500 ml 1000 ml	750 1250
73469 [116057-81-9]	**N,N-Dimethylacetamide-d9 (DMA-d9) for NMR spectroscopy, 99 Atom %D**	5 ml	70000
10652 [124-40-3]	**Dimethylamine 40% Aq. Solution**	500 ml 2500 ml	300 1150
14003 [100-10-7]	**p-Dimethylaminobenzaldehyde (DMAB) ACS, 98%** (Ehrlichs reagent)	100 Gms 1 Kg	1830 17060
77086 [100-10-7]	**p-Dimethylaminobenzaldehyde (DMAB) ExiPlus, Multi-Compendial, 98%** (Ehrlichs reagent) meets compendial specs of USP, BP, Ph.Eur	100 Gms 1 Kg	2020 18190
12566 [121-69-7]	**N,N-Dimethylaniline pure, 98%**	500 ml 2500 ml	610 2880
20364 [121-69-7]	**N,N-Dimethylaniline Reagent Grade, 99%**	500 ml 2500 ml	670 3110
16616 [68-12-2]	**N,N-Dimethylformamide (DMF) extrapure, 99%**	500 ml 2500 ml 25 Ltr	330 1580 15000
21974 [68-12-2]	**N,N-Dimethylformamide (DMF) ACS, 99.8%**	500 ml 2500 ml	370 1700
79121 [68-12-2]	**N,N-Dimethylformamide (DMF) ExiPlus, Multi-Compendial, 99.8%** meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	410 1850
50591 [68-12-2]	**N,N-Dimethylformamide (DMF) for HPLC & UV Spectroscopy, 99.9%**	500 ml 1000 ml 2.5 Ltr (Wide-Mouth)	510 870 4830
46165 [68-12-2]	**N,N-Dimethylformamide (DMF) GC-HS, 99.9%**	250 ml 500 ml 1000 ml	410 640 1130
93186 [68-12-2]	**N,N-Dimethylformamide (DMF) Dried, 99.5%, water 0.005%**	500 ml 1000 ml	580 1040
24017 [68-12-2]	**N,N-Dimethylformamide (DMF) for molecular biology, 99.9%**	100 ml 500 ml	520 980
14238 [4472-41-7]	**N,N-Dimethylformamide-d7 (DMF-d7) for NMR spectroscopy, 99.5 Atom %D** (Heptadeutero-N,N-Dimethylformamide, DMF-d)	1 ml 5 ml	16000 60000
52794 [4637-24-5]	**N,N-Dimethylformamide Dimethylacetal extrapure, 95%** (1,1-Dimethoxy-N,N-Dimethylmethylamine)	100 ml 500 ml	650 2880
81890 [95-45-4]	**Dimethylglyoxime ACS, 99%** (2,3-Butanedione Dioxime)	100 Gms 500 Gms	670 2880
19711 [95-45-4]	**Dimethylglyoxime ExiPlus, Multi-Compendial, 99%** (2,3-Butanedione Dioxime) meets compendial specs of BP, Ph.Eur	100 Gms 500 Gms	810 3220
92574 [75006-64-3]	**Dimethylglyoxime Disodium Salt Octahydrate extrapure AR, 99%**	250 Gms	1220
33158 [80-73-9]	**1,3-Dimethyl-2-Imidazolidinone (DMI) pure, 99%**	100 ml 500 ml 2500 ml	1370 5790 28380

Code	Name	Packing	Price
35161 [80-73-9]	**1,3-Dimethyl-2-Imidazolidinone (DMI) GC-HS, 99.9%**	100 ml 500 ml	1600 6100
60496 [80-73-9]	**1,3-Dimethyl-2-Imidazolidinone (DMI) ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	100 ml 500 ml	1820 6270
54272 [77-78-1]	**Dimethyl Sulphate pure, 99%** (Sulfuric Acid Dimethyl Ester)	500 ml	410
43404 [67-68-5]	**Dimethyl Sulphoxide (DMSO) extrapure, 99%**	500 ml 2500 ml 25 Ltr	1050 4800 45000
28580 [67-68-5]	**Dimethyl Sulphoxide (DMSO) ACS, 99.9%**	500 ml 2500 ml 25 Ltr	1150 5000 47000
66444 [67-68-5]	**Dimethyl Sulphoxide (DMSO) ACS, ExiPlus, Multi-Compendial, 99.9%** meets compendial specs of USP, Ph.Eur	500 ml 2500 ml	1300 5400
12767 [67-68-5]	**Dimethyl Sulphoxide (DMSO) GC-HS, 99.9%**	250 ml 500 ml 1000 ml	950 1800 3200
30239 [67-68-5]	**Dimethyl Sulphoxide (DMSO) for HPLC & UV Spectroscopy, 99.8%**	500 ml 1000 ml 2.5 Ltr (Wide-Mouth)	1700 3000 9000
71545 [67-68-5]	**Dimethyl Sulphoxide (DMSO) ULSI Semiconductor Grade, 99.9%**	1000 ml 2500 ml	9000 18000
66120 [67-68-5]	**Dimethyl Sulphoxide (DMSO) Dried, 99.5%, water 0.005%**	500 ml	1860
24075 [67-68-5]	**Dimethyl Sulphoxide (DMSO) for molecular biology, 99.8%**	100 ml 250 ml 500 ml	650 1300 2300
47809 [2206-27-1]	**Dimethyl Sulphoxide-d6 (DMSO-d6) for NMR spectroscopy, 99.8 Atom %D** (Methyl Sulfoxide-d6; Hexadeuterodimethyl Sulfoxide)	10 Gms 25 Gms	6300 12000
69545 [2206-27-1]	**Dimethyl Sulphoxide-d6 (with 0.03% TMS) (DMSO-d6 w/ 0.03% TMS) for NMR spectroscopy, 99.8 Atom %D** (Methyl Sulfoxide-d6; Hexadeuterodimethyl Sulfoxide)	10 Gms 25 Gms	6700 13200
36007 [2206-27-1]	**Dimethyl Sulphoxide-d6 (DMSO-d6) for NMR spectroscopy, 99.8 Atom %D** (Methyl Sulfoxide-d6, Hexadeuterodimethyl Sulfoxide)	10 ml 25 ml	7000 13000
42009 [96-31-1]	**N,N-Dimethylurea (1,3-Dimethylurea) extrapure, 99%**	500 Gms 5 Kg 25 Kg	580 5180 23000
19191 [96-31-1]	**N,N-Dimethylurea (1,3-Dimethylurea) (DMU) PrimeDx, 99%**	1 Kg 5 Kg 25 Kg	1730 6900 29900
43981 [117-81-7]	**Dioctyl Phthalate pure, 98%**	500 ml 2500 ml	530 2300
56353 [577-11-7]	**Dioctyl Sodium Sulfosuccinate (DOSS) ExiPlus, Multi-Compendial, 98%** (Aerosol OT, Ducusate Sodium) meets compendial specs of USP, BP	500 Gms 5 Kg	2190 20700
43324 [123-91-1]	**1,4-Dioxane pure, 99%**	500 ml 2500 ml 25 Ltr	490 2070 19550
88797 [123-91-1]	**1,4-Dioxane ACS, 99%**	500 ml 2500 ml	580 2300
95536 [123-91-1]	**1,4-Dioxane ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	690 2760
29469 [123-91-1]	**1,4-Dioxane Dried, 99.5%, water 0.005%**	500 ml	730
79218 [123-91-1]	**1,4-Dioxane GC-HS, 99.9%**	250 ml 1000 ml	710 2520
43294 [123-91-1]	**1,4-Dioxane for HPLC & UV Spectroscopy, 99.9%**	500 ml 1000 ml	1270 1940
82736 [646-06-0]	**1,3-Dioxolane extrapure AR, 99%**	500 ml 2500 ml	810 3450
98342 [122-39-4]	**Diphenylamine ACS, 99%**	100 Gms 500 Gms	300 920
41149 [122-39-4]	**Diphenylamine ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	100 Gms 500 Gms	370 1040
62567 [538-62-5]	**1,5-Diphenylcarbazone (mixture with Diphenylcarbazide) extrapure AR, ACS** (Phenylazoformic Acid-2-Phenylhydrazide)	5 Gms 25 Gms 100 Gms	640 2790 10090
75522 [92-71-7]	**2,5-Diphenyloxazole scintillation grade (PPO), 99%**	25 Gms 100 Gms 500 Gms	1960 7130 29900
94777 [108-65-6]	**Propylene Glycol Monomethyl Ether 2-Acetate (MPA) extrapure, 98%** (1-Methoxy-2-Propyl Acetate, PGMEA, 2-Acetoxy-1-Methoxypropane 1,2-Propanediol 1-Monomethyl Ether 2-	500 ml 2500 ml	1380 6670
43478 [34590-94-8]	**Dipropylene Glycol Monomethyl Ether (Mixture of Isomers) extrapure, 99%** (DPM, Dowanol)	500 ml 2500 ml	520 2300
88738 [6106-24-7]	**Disodium Tartrate Dihydrate extrapure, 99%**	500 Gms 5 Kg	920 8740
44700 [6106-24-7]	**Disodium Tartrate Dihydrate ACS, 99%**	500 Gms 5 Kg	1040 9660
16261 [6106-24-7]	**Disodium Tartrate Dihydrate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP	250 Gms 500 Gms	690 1150

Code	Name	Packing	Price
12852 [6106-24-7]	**Disodium Tartrate Dihydrate Calibration Standard (For Karl Fischer Titration), 99.5%**	100 Gms 250 Gms	6660 12650
88283 ✔ [112-40-3]	**Dodecane extrapure, 99%** (n-Dodecane)	100 ml 500 ml	1500 6330
25821 [2386-53-0]	**Dodecane Sulphonic Acid Sodium Salt Anhydrous for HPLC, 99%** (Sodium 1-Dodecanesulfonate)	5 Gms 25 Gms 100 Gms	1040 2530 8630
56085 [112-53-8]	**1-Dodecanol (Dodecyl Alcohol, Lauryl Alcohol) pure, 98%**	500 ml 2500 ml	780 3520
23715 [112-53-8]	**1-Dodecanol (Dodecyl Alcohol, Lauryl Alcohol) ACS, 98%**	500 ml 2500 ml	870 4030
93498 [112-53-8]	**1-Dodecanol (Dodecyl Alcohol, Lauryl Alcohol) ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	960 4370
88147	**DPX Mountant for histology**	250 ml 500 ml 1000 ml	440 830 1430
85675	**Dragendorff`s Reagent**	125 ml	1960
71867 [67762-27-0]	**Emulsifying Wax pure**	500 Gms 5 Kg	600 5520
56605 [548-24-3]	**Eosin Blue (C.I. No. 45400), 85%**	25 Gms 100 Gms	920 3400
29391 [17372-87-1]	**Eosin Yellow (water soluble) ACS**	25 Gms 100 Gms 500 Gms	350 990 4820
22115 [106-89-8]	**Epichlorohydrin pure, 98%** (1-Chloro-2,3-Epoxypropane)	500 ml 2500 ml 25 Ltr	520 2300 20700
73317 [106-89-8]	**Epichlorohydrin extrapure AR, 99%** (1-Chloro-2,3-Epoxypropane)	500 ml 2500 ml	600 2650
40876 [1787-61-7]	**Eriochrome Black T ACS**	25 Gms 100 Gms 500 Gms	410 1200 4600
76512 [3564-18-9]	**Eriochrome Cyanine R**	10 Gms 25 Gms 100 Gms	750 1700 6200
69059 [16423-68-0]	**Erythrosin B, 90%**	25 Gms 100 Gms	1040 3890
73961 [8007-09-8]	**Eschka`s Mixture extrapure AR, 38%**	250 Gms 1 Kg	6990 26500
40942 [141-43-5]	**Ethanolamine (MEA, Monoethanolamine, 2-Aminoethanol) extrapure, 99%**	500 ml 2500 ml	410 1840
57100 [141-43-5]	**Ethanolamine (MEA, Monoethanolamine, 2-Aminoethanol) ACS, 99%**	500 ml 2500 ml	460 2070
74683 [141-43-5]	**Ethanolamine (MEA, Monoethanolamine, 2-Aminoethanol) ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	560 2420
17220 [1239-45-8]	**Ethidium Bromide extrapure, 95%** (Homidium Bromide)	1 Gms 5 Gms	620 2690
93079 [1239-45-8]	**Ethidium Bromide for molecular biology, 95%** (Homidium Bromide)	1 Gms 5 Gms	670 2940
16201	**Ethidium Bromide Solution (10mg/ml)**	10 ml 50 ml	800 2000
89362 [141-78-6]	**Ethyl Acetate pure, 99%**	500 ml 2500 ml 25 Ltr	300 1340 11770
66183 [141-78-6]	**Ethyl Acetate ACS, 99.5%**	500 ml 2500 ml	370 1500
62952 [141-78-6]	**Ethyl Acetate ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	410 1720
72598 [141-78-6]	**Ethyl Acetate Dried, 99.5%, water 0.005%**	500 ml 1000 ml	430 830
73106 [141-78-6]	**Ethyl Acetate GC-HS, 99.9%**	500 ml 1000 ml	560 970
26614 [141-78-6]	**Ethyl Acetate for HPLC & UV Spectroscopy, 99.9%**	500 ml 1000 ml 2500 ml	480 700 1610
22636 [141-97-9]	**Ethyl Acetoacetate pure, 99%**	500 ml 2500 ml	490 2250
33123 [75-04-7]	**Ethylamine Solution 70% in water pure**	500 ml 2500 ml	540 2360
37810 ⓘ [74-96-4]	**Ethyl Bromide extrapure, 99%**	250 ml 500 ml	1240 2360
87095 [110-80-5]	**Ethyl Cellosolve pure, 99%** (Ethylene Glycol Monoethyl Ether)	500 ml 2500 ml	410 1610
48154 △ [541-41-3]	**Ethyl Chloroformate pure, 98%** (Ethyl Chlorocarbonate, Ethoxycarbonyl Chloride, Chloroformic Acid Ethyl Ester)	500 ml 2500 ml	600 2850
30166 [1678-91-7]	**Ethylcyclohexane extrapure, 99%**	500 ml 2500 ml	450 2140
65428 [107-15-3]	**Ethylene Diamine pure, 99%** (1,2-Diaminoethane)	500 ml 2500 ml 25 Ltr	390 1660 14980
70790 [107-15-3]	**Ethylene Diamine Reagent Grade, 99%** (1,2-Diaminoethane)	500 ml 2500 ml	450 1820
40648 [60-00-4]	**EDTA Free Acid extrapure, 98%** (Ethylenediaminetetraacetic Acid Free Acid)	500 Gms 5 Kg	620 5600

Code	Name	Packing	Price
18240 [60-00-4]	**EDTA Free Acid ACS, 99.4%** (Ethylenediaminetetraacetic Acid Free Acid)	100 Gms 500 Gms 5 Kg	280 650 6050
73816 [60-00-4]	**EDTA Free Acid ExiPlus, Multi-Compendial, 98%** (Ethylenediaminetetraacetic Acid Free Acid) meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	730 6500
44478 [25102-12-9]	**EDTA Dipotassium Salt pure, 98%** (Ethylenediaminetetraacetic Acid Dipotassium Salt)	100 Gms 500 Gms 5 Kg	300 650 6000
62196 [25102-12-9]	**EDTA Dipotassium Salt extrapure AR, 99%** (Ethylenediaminetetraacetic Acid Dipotassium Salt)	100 Gms 500 Gms 5 Kg	340 720 6300
12070 [6381-92-6]	**EDTA Disodium Salt Dihydrate pure, 98%** (Ethylenediaminetetraacetic Acid Disodium Salt Dihydrate)	500 Gms 5 Kg 25 Kg	600 5200 22000
40088 [6381-92-6]	**EDTA Disodium Salt Dihydrate ACS, 99%** (Ethylenediaminetetraacetic Acid Disodium Salt Dihydrate)	100 Gms 500 Gms 5 Kg	245 700 6800
50027 [6381-92-6]	**EDTA Disodium Salt Dihydrate ExiPlus, Multi-Compendial, 99%** (Ethylenediaminetetraacetic Acid Disodium Salt Dihydrate) meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms 5 Kg	270 720 7000
43272 [6381-92-6]	**EDTA Disodium Salt Dihydrate for molecular biology, 99.5%** (Ethylenediaminetetraacetic Acid Disodium Salt Dihydrate)	100 Gms 500 Gms 1 Kg 5 Kg	330 1000 1900 9000
64762 [6381-92-6]	**EDTA Disodium Salt Dihydrate PrimeDx, 98%**	1 Kg 5 Kg 25 Kg	1400 5800 26000
35108 [6381-92-6]	**EDTA Disodium Salt Dihydrate for tissue culture, 99.5%** (Ethylenediaminetetraacetic Acid Disodium Salt Dihydrate)	100 Gms 500 Gms	260 800
96540 [6381-92-6]	**EDTA Disodium Salt 0.1M Solution**	500 ml 1000 ml	370 670
48960 [6381-92-6]	**EDTA Disodium Salt 0.01M Solution**	500 ml 1000 ml	350 640
35888 [6381-92-6]	**EDTA Disodium Salt 0.5M Solution**	500 ml 1000 ml	500 950
73557 [6381-92-6]	**EDTA Disodium Salt 0.05M Solution**	500 ml 1000 ml	250 450
59389 [15708-41-5]	**EDTA Ferric Monosodium Salt extrapure, 12-14% Fe** (Ethylenediaminetetraacetic Acid Ferric Monosodium Salt)	500 Gms	640
41000 [14402-88-1]	**EDTA Magnesium Disodium Complex (EDTA Mg) pure** (Ethylenediaminetetraacetic Acid Magnesium Disodium Complex)	100 Gms 500 Gms 5 Kg	300 700 6500
15078 [10378-23-1]	**EDTA Tetrasodium Salt Dihydrate extrapure, 98%** (Ethylenediaminetetraacetic Acid Tetrasodium Salt Dihydrate)	100 Gms 500 Gms	300 700
78036 [10378-23-1]	**EDTA Tetrasodium Salt Dihydrate extrapure AR, 99%** (Ethylenediaminetetraacetic Acid Tetrasodium Salt Dihydrate)	100 Gms 500 Gms	400 1200
45247 [10378-23-1]	**EDTA Tetrasodium Salt Dihydrate for molecular biology, 99.5%** (Ethylenediaminetetraacetic Acid Tetrasodium Salt Dihydrate)	100 Gms 500 Gms	500 1500
24377 [13235-36-4]	**EDTA Tetrasodium Salt Tetrahydrate extrapure, 99%** (Ethylenediaminetetraacetic Acid Tetrasodium Salt Tetrahydrate)	100 Gms 500 Gms	400 900
19950 [85715-60-2]	**EDTA Trisodium Salt Hydrate extrapure, 95%** (Ethylenediaminetetraacetic Acid Trisodium Hydrate)	100 Gms 500 Gms 5 Kg	320 1300 12000
47471 [107-21-1]	**Ethylene Glycol (MEG, Monoethylene Glycol) pure, 98%** (Ethanediol)	500 ml 2500 ml 25 Ltr	340 1400 12000
95361 [107-21-1]	**Ethylene Glycol (MEG, Monoethylene Glycol) Reagent Grade, 99%**	500 ml 2500 ml	365 1600
27466 [122-99-6]	**Ethylene Glycol Monophenylether pure, 99%** (Phenyl cellosolve, 2-Phenoxyethanol)	500 ml 2500 ml	560 2800
82832 [2807-30-9]	**Propyl Cellosolve (Ethylene Glycol Monopropyl Ether) extrapure, 99%** (2-Propoxyethanol, Propyl glycol)	500 ml 2500 ml	1350 5380
25049 ✔ [60-29-7]	**Ethyl Ether pure, 98%** (Diethyl ether)	500 ml 2500 ml 25 Ltr	500 1980 16690
64665 ✔ [60-29-7]	**Ethyl Ether ACS, 98%** (Diethyl ether)	500 ml 2500 ml	560 2300
83416 ✔ [60-29-7]	**Ethyl Ether ExiPlus, Multi-Compendial, 98%** (Diethyl ether) meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	620 2380
83131 ✔ [60-29-7]	**Ethyl Ether Dried, 99.5%, water 0.005%** (Diethyl Ether)	500 ml 1000 ml 2500 ml	610 1160 2520
68329 ✔ [60-29-7]	**Ethyl Ether for HPLC & UV Spectroscopy, 99.9%**	500 ml 1000 ml	2580 4620
87796 ✔ [60-29-7]	**Ethyl Ether GC-HS, 99.9%**	500 ml 1000 ml	2590 4890
36979 ⓘ [109-94-4]	**Ethyl Formate pure, 98%**	500 ml 2500 ml 25 Ltr	560 2470 22400
50371 [2390-59-2]	**Ethyl Violet (C.I. No. 42600) for microscopy**	25 Gms	990
89474	**Fehling`s solution A**	500 ml	450
12603	**Fehling`s solution B**	500 ml	680

Code	Name	Packing	Price
25106 [13268-42-3]	**Ferric Ammonium Oxalate Trihydrate pure, 98%** (Ammonium Iron Oxalate Trihydrate, Ammonium Iron (III) Oxalate Trihydrate)	100 Gms 500 Gms	1010 3360
72287 [7705-08-0]	**Ferric Chloride Anhydrous pure, 98%** (Iron (III) Chloride)	500 Gms 5 Kg 25 Kg	320 2610 9800
66492 [7705-08-0]	**Ferric Chloride Anhydrous ExiPlus, Multi-Compendial, 98%** (Iron (III) Chloride) meets compendial specs of BP	500 Gms 5 Kg	460 3920
18831 [3522-50-7]	**Ferric Citrate Trihydrate pure, 18-20% Fe** (Iron (III) Citrate Trihydrate)	500 Gms	1100
92273 [3522-50-7]	**Ferric Citrate Trihydrate for tissue culture, 18-20% Fe** (Iron (III) Citrate Trihydrate)	500 Gms	1160
72655 [10025-77-1]	**Ferric Chloride Hexahydrate ACS, 97%** (Iron (III) chloride Hexahydrate)	500 Gms 5 Kg	770 7060
87144 [7782-61-8]	**Ferric Nitrate Nonahydrate ACS, 98%** (Iron (III) Nitrate Nonahydrate)	500 Gms	560
28454 ✔ [10045-86-0]	**Ferric Phosphate pure, 97%** (Iron(III) Phosphate)	500 Gms	1240
69325 [1309-37-1]	**Ferric Oxide Red extrapure, 95%** (Iron (III) Oxide Red)	500 Gms	500
39102 [15244-10-7]	**Ferric Sulphate Hydrate pure, 19-21% Fe** (Iron (III) Sulphate Hydrate)	500 Gms 5 Kg	320 3030
14621 [15244-10-7]	**Ferric Sulphate Hydrate extrapure AR, 19-21% Fe** (Iron (III) Sulphate Hydrate)	500 Gms 5 Kg	360 3360
74871 [15244-10-7]	**Ferric Sulphate Hydrate for tissue culture, 19-21% Fe** (Iron (III) Sulphate Hydrate)	500 Gms 5 Kg	480 4480
59609 [2944-68-5]	**Ferric Tartrate pure, Fe 18-23%** (Iron (III) Tartrate, Iron(II)Tartratemonohydrate)	100 Gms 500 Gms	2110 8400
57888 [14634-91-4]	**Ferroin solution AR**	100 ml	1060
37323 [14634-91-4]	**Ferroin Solution AR, ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	100 ml	1090
82984 [13478-10-9]	**Ferrous Chloride Tetrahydrate extrapure, 98%** (Iron (II) Chloride Tetrahydrate)	500 Gms 5 Kg	620 5600
94185 [1317-37-9]	**Ferrous Sulphide Fused Sticks, 23-27% S** (Iron(II) Sulfide, Ferrous Sulfide)	1 Kg 5 Kg	900 4320
93668 [7782-63-0]	**Ferrous Sulphate Heptahydrate for tissue culture, 99.5%** (Iron (II) Sulphate Heptahydrate)	500 Gms	460
31929 [7782-63-0]	**Ferrous Sulphate Heptahydrate pure, 98%** (Iron (II) Sulphate Heptahydrate)	500 Gms 5 Kg 25 Kg	260 1570 6720
97868 [7782-63-0]	**Ferrous Sulphate Heptahydrate ACS, 99%** (Iron (II) Sulphate Heptahydrate)	500 Gms 5 Kg	330 2940
53072 [75-12-7]	**Formamide pure, 98.5%** (Methanamide)	500 ml 2500 ml 25 Ltr	500 2360 21280
71714 [75-12-7]	**Formamide ACS, 99.5%** (Methanamide)	500 ml 2500 ml	650 2800
74841 [75-12-7]	**Formamide ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	730 3140
30349 [75-12-7]	**Formamide for molecular biology, 99.5%** (Methanamide)	500 ml 1000 ml	960 1570
83256 [75-12-7]	**Formamide for UV Spectroscopy, 99%** (Methanamide)	500 ml 1000 ml	1350 2360
62673 [64-18-6]	**Formic Acid Reagent Grade, 98%**	500 ml	560
59894 [64-18-6]	**Formic Acid for HPLC, 98%**	500 ml	2800
21547 △ [925-94-0]	**Formic Acid-d for NMR spectroscopy, 98 Atom%D, 95%w in D2O** (Formic-d Acid-d, Perdeuterioformic Acid)	5 ml	26000
42868 [57-48-7]	**D-Fructose extrapure, 98%** (D Fructose, D (-) Fructose)	500 Gms 5 Kg 25 Kg	550 5200 25000
68823 [57-48-7]	**D-Fructose for tissue culture, 98%** (D Fructose, D (-) Fructose)	500 Gms	950
65628 [57-48-7]	**D-Fructose ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	700 6000
86364 [110-17-8]	**Fumaric Acid pure, 99%**	500 Gms 5 Kg 25 Kg	400 3800 15000
79503 [110-17-8]	**Fumaric Acid extrapure AR, 99%**	500 Gms 5 Kg	600 5700
97328 [110-17-8]	**Fumaric Acid extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	500 Gms 5 Kg	800 7500
14277 [632-99-5]	**Fuchsin Basic** (Basic Fuchsin) (C.I. No. 42510)	25 Gms 100 Gms 500 Gms	280 770 3400
62967	**Fusion Mixture extrapure**	500 Gms	430
24016	**Fusion Mixture extrapure AR**	500 Gms	540
13142 [5995-86-8]	**Gallic Acid pure, 98%** (3,4,5-Trihydroxy Benzoic Acid)	100 Gms 500 Gms	850 3000
41291 [5995-86-8]	**Gallic Acid ACS, 98%** (3,4,5-Trihydroxy Benzoic Acid)	100 Gms 500 Gms	850 3000

Code	Name	Packing	Price
42306 [5995-86-8]	**Gallic Acid ExiPlus, Multi-Compendial, 98%** (3,4,5-Trihydroxy Benzoic Acid) meets compendial specs of BP, Ph.Eur	100 Gms 500 Gms	1100 3400
83740 [9000-70-8]	**Gelatin Powder ex. Porcine (Type B), 120g Bloom**	500 Gms 5 Kg	1500 13000
95110 [77-06-5]	**Gibberellic Acid (GA3), 90%**	1 Gms 10 Gms 100 Gms	200 1500 11600
24367 [468-44-0, 510-75-8]	**Gibberellic Acid (GA4+7), 90%**	1 Gms 5 Gms 25 Gms	900 2900 12100
94857 [77-06-5]	**Gibberellic Acid (GA3) for tissue culture, 90%**	1 Gms 10 Gms 100 Gms	200 1700 12500
45881 [51811-82-6]	**Giemsa Stain**	25 Gms 100 Gms 500 Gms	600 2290 11000
72402 ✔	**Giemsa Stain Solution** (Giemsa Azur-Eosin-Methylene Blue Solution)	100 ml 250 ml 500 ml	310 610 1100
20149 [56-86-0]	**L-Glutamic Acid extrapure CHR, 99%**	250 Gms 1 Kg 5 Kg	650 2100 9000
50277 [56-86-0]	**L-Glutamic Acid for tissue culture, 99%**	100 Gms 500 Gms	450 1500
18376 [56-86-0]	**L-Glutamic Acid ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP	100 Gms 250 Gms 1 Kg	380 750 2800
23229 [6106-04-3]	**L-Glutamic Acid Monosodium Salt Monohydrate (MSG) extrapure, 99%** (Monosodium Glutamate, Monosodium L-Glutamate)	500 Gms 5 Kg 25 Kg	470 4360 19620
69382 [6106-04-3]	**L-Glutamic Acid Monosodium Salt Monohydrate (MSG) PrimeDx, 99%**	5 Kg 25 Kg	6000 25000
42595 [56-81-5]	**Glycerol (Glycerine) Anhydrous pure, 99%**	500 ml 2500 ml 25 Ltr	450 2100 21000
77453 [56-81-5]	**Glycerol (Glycerine) Anhydrous ACS, 99.5%**	500 ml 2500 ml	500 2350
59991 [56-81-5]	**Glycerol (Glycerine) Anhydrous ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	570 2550
62417 [56-81-5]	**Glycerol (Glycerine) for molecular biology, 99.5%**	100 ml 500 ml	520 1790
85916 [60-01-5]	**Glycerol Tributyrate (Tributyrin) pure, 97%**	100 ml 500 ml	1620 6750

Code	Name	Packing	Price
25853 [56-40-6]	**Glycine pure, 99%**	500 Gms 5 Kg 25 Kg	730 6930 33600
66327 [56-40-6]	**Glycine ACS, 99.5%**	500 Gms 5 Kg 25 Kg	760 7250 35180
64072 [56-40-6]	**Glycine for molecular biology, 99.5%**	100 Gms 500 Gms 1 Kg 5 Kg	280 1220 2210 9050
72938 [56-40-6]	**Glycine for tissue culture, 99.5%**	500 Gms 5 Kg	890 8160
52574 [56-40-6]	**Glycine for electrophoresis, 99.5%**	100 Gms 500 Gms	1550 6620
69422 [56-40-6]	**Glycine ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	800 7500
40682 [56-40-6]	**Glycine IP**	1 Kg 5 Kg 25 Kg	1330 5520 24820
74072 [56-40-6]	**Glycine BP**	1 Kg 5 Kg 25 Kg	1330 5520 24820
37449 [56-40-6]	**Glycine USP**	1 Kg 5 Kg 25 Kg	1550 6900 33080
38063 [79-14-1]	**Glycolic Acid extrapure, 70% soln in water**	500 ml 2500 ml	1000 4520
27194	**Gold AAS Standard Solution in 0.5N HCl (Traceable to NIST)**	100 ml 250 ml	7600 18000
59611 🌡 [16961-25-4]	**Gold Chloride Hydrate (Tetrachloroauric Acid) ultrapure, 99.99%, ~49% Au**	500 Mg 1 Gms 5 Gms	24500 42000 205000
12023 🌡 [16961-25-4]	**Gold Chloride Trihydrate (Tetrachloroauric Acid) pure, ~49% Au**	500 Mg 1 Gms 5 Gms	23000 39000 190000
81173 🌡 [13874-02-7]	**Gold Chloride Sodium, 49% Au** (Sodium tetrachloroaurate(III) dihydrate)	1 Gms 5 Gms	43000 200000
84763 🌡 [16903-35-8]	**Gold Chloride Solution 30% in HCl, 17% Au**	5 ml 10 ml	62000 110000
93281 △ [13682-61-6]	**Gold Potassium Chloride (Potassium Tetrachloroaurate (III)) extrapure, 99.99%** (Gold Potassium Chloride)	500 Mg 1 Gms	25000 42000
10436	**Gram`s Iodine, 33% I**	5 Gms 25 Gms 100 Gms 500 Gms	500 1200 4500 21200

Code	Name	Packing	Price
59538 [7782-42-5]	**Graphite Synthetic (Type 1) pure, 98% C, 50mesh**	500 Gms	350
		2.5 Kg	1500
19123 [7782-42-5]	**Graphite Synthetic (Type 2) pure, 97% C, 100mesh**	500 Gms	400
		2.5 Kg	1600
26870 [7782-42-5]	**Graphite Synthetic (Type 3) pure, 96% C, 325mesh**	500 Gms	450
		2.5 Kg	1800
86601 [593-85-1]	**Guanidine Carbonate extrapure, 99%**	100 Gms	420
		500 Gms	1250
		5 Kg	11000
54542 [50-01-1]	**Guanidine Hydrochloride (GHC) pure, 99%**	100 Gms	350
		500 Gms	1200
		1 Kg	2200
		5 Kg	10000
45539 [50-01-1]	**Guanidine Hydrochloride (GHC) for molecular biology, 99.5%**	100 Gms	700
		500 Gms	2800
		1 Kg	4600
		5 Kg	22000
80272 [593-84-0]	**Guanidine Thiocyanate (GTC) for molecular biology, 99%**	100 Gms	800
		500 Gms	3200
		1 Kg	5400
		5 Kg	26000
68413 [50-01-1]	**Guanidine Hydrochloride (GHC) ExiPlus, Multi-Compendial, 99%**	250 Gms	1500
		1 Kg	4300
		5 Kg	20000
	meets compendial specs of USP, BP		
20222 [9000-30-0]	**Guar Gum extrapure**	500 Gms	420
		5 Kg	3500
90970 [9000-01-5]	**Acacia (Gum Acacia) extrapure** (Gum Arabic)	500 Gms	560
		5 Kg	4830
		25 Kg	23580
90938 [9000-01-5]	**Acacia (Gum Acacia) (enzyme free) extrapure AR** (Gum Arabic)	500 Gms	850
39758	**Hematoxylin solution (Delafield`s) for microscopy**	125 ml	490
78631	**Hematoxylin solution (Ehrlich) for microscopy**	125 ml	490
40362	**Hematoxylin solution (Harris) for microscopy**	125 ml	650
48441	**Hematoxylin solution (Mayer`s) for microscopy**	100 ml	270
19241 [142-82-5]	**n-Heptane extrapure, 99%**	500 ml	360
		2500 ml	1500
		25 Ltr	14000
94979 [142-82-5]	**n-Heptane Reagent Grade, 99%**	500 ml	430
		1000 ml	860
		2500 ml	1800
		25 Ltr	17000
96979 [142-82-5]	**n-Heptane ExiPlus, Multi-Compendial, 99%**	500 ml	460
		2500 ml	2100
	meets compendial specs of USP, BP, Ph.Eur		

Code	Name	Packing	Price
90731 [142-82-5]	**n-Heptane for HPLC & UV Spectroscopy, 99.5%**	500 ml	600
		1000 ml	1100
99295 [142-82-5]	**n-Heptane Dried, 99%, water 0.005%**	500 ml	580
94373 [22767-50-6]	**Heptane Sulphonic Acid Sodium Salt Anhydrous for HPLC, 99%**	25 Gms	1300
		100 Gms	4600
		250 Gms	11300
55272 [207300-90-1]	**Heptane Sulphonic Acid Sodium Salt Monohydrate for HPLC, 99%**	25 Gms	1300
		100 Gms	4600
		250 Gms	11300
55120 [111-70-6]	**1-Heptanol extrapure, 98%** (n-Heptanol, n-Heptyl Alcohol)	500 ml	1570
		2500 ml	7280
60805 [67-72-1]	**Hexachloroethane pure, 99%**	500 Gms	730
		5 Kg	6720
36027 [18497-13-7]	**Hexachloroplatinic (IV) Acid Hexahydrate ACS, 38-40% Pt** (Platinic Chloride Hexahydrate, Hydrogen Hexachloroplatinate (IV), Platinum Chloride Hexahydrate)	250 Mg	5000
		1 Gms	14000
		5 Gms	64000
		25 Gms	310000
37465 [629-73-2]	**1-Hexadecene (Cetene) extrapure, 95%**	250 ml	2850
		500 ml	4500
63539 [544-76-3]	**Hexadecane (Cetane) extrapure, 98%** (n-Hexadecane)	100 ml	1850
		500 ml	6000
81451 [100-97-0]	**Hexamine pure, 99%** (Hexamethylenetetramine, HMTA)	500 Gms	300
		5 Kg	2600
		25 Kg	11000
92391 [100-97-0]	**Hexamine ACS, 99%**	500 Gms	350
		5 Kg	3000
91454 [110-54-3]	**n-Hexane pure, 95%**	500 ml	500
		2500 ml	2300
		25 Ltr	19000
65764 [110-54-3]	**n-Hexane for HPLC, 95%**	500 ml	620
		1000 ml	1200
		2500 ml	2800
77045 [110-54-3]	**n-Hexane pure, 99%**	500 ml	600
		2500 ml	2600
		25 Ltr	24000
12534 [110-54-3]	**n-Hexane Reagent Grade, 99%**	500 ml	680
		2500 ml	2800
27226 [110-54-3]	**n-Hexane ExiPlus, Multi-Compendial, 99%**	500 ml	740
		2500 ml	3000
	meets compendial specs of USP, BP, Ph.Eur		
33823 [110-54-3]	**n-Hexane for HPLC, 99%**	500 ml	800
		1000 ml	1500
		2500 ml	3500
		2.5 Ltr (Wide-Mouth)	6400
		4 Ltr (Wide-Mouth)	9000

Code	Name	Packing	Price
40501 [110-54-3]	**n-Hexane GC-HS, 99.9%**	250 ml	820
		1000 ml	2700
29114 [110-54-3]	**n-Hexane Dried, 99%, water 0.005%**	500 ml	1050
		1000 ml	2000
49346 [2832-45-3]	**Hexane Sulphonic Acid Sodium Salt Anhydrous for HPLC, 99%**	25 Gms	1300
		100 Gms	4600
		250 Gms	11300
78597 [207300-91-2]	**Hexane Sulphonic Acid Sodium Salt Monohydrate for HPLC, 99%**	25 Gms	1300
		100 Gms	4600
		250 Gms	11300
70389 [107-41-5]	**Hexylene Glycol pure, 99%**	500 ml	680
		2500 ml	2920
		25 Ltr	26880
38885 [351-50-8]	**D-Histidine extrapure, 99%**	5 Gms	3030
		25 Gms	8000
94495 [71-00-1]	**L-Histidine (base) extrapure CHR, 99%**	25 Gms	800
		100 Gms	2800
		500 Gms	13000
		1 Kg	24500
		5 Kg	116000
42437 [71-00-1]	**L-Histidine ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	25 Gms	1000
		100 Gms	3390
		1 Kg	26000
28061 [5934-29-2]	**L-Histidine Hydrochloride Monohydrate extrapure CHR, 99%**	25 Gms	710
		100 Gms	2200
		500 Gms	8000
		1 Kg	14200
		5 Kg	68000
51565 [4998-57-6]	**DL-Histidine extrapure, 99%**	5 Gms	2000
		25 Gms	8000
		100 Gms	30000
40599 [12055-62-8]	**Holmium (III) Oxide ultrapure, 99.9%**	5 Gms	2000
		10 Gms	3500
		25 Gms	7500
59764 [68514-28-3]	**Humic Acid Potassium Salt (Potassium Humate) pure, 70% Humic**	500 Gms	380
		5 Kg	3500
53931 [7803-57-8]	**Hydrazine Hydrate extrapure, 80%** (Hydrazinium Hydroxide)	500 ml	790
		2500 ml	3140
		25 Ltr	28000
99313 [10034-93-2]	**Hydrazine Sulphate extrapure, 99%**	100 Gms	500
		500 Gms	2190
33243 [10034-93-2]	**Hydrazine Sulphate ACS, 99%**	100 Gms	620
		500 Gms	2580
78110 [10034-93-2]	**Hydrazine Sulphate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	800
		500 Gms	2930
21435 [10035-10-6]	**Hydrobromic Acid pure, 48-49% in water**	500 ml	1070
		2500 ml	4930
30106 [10035-10-6]	**Hydrobromic Acid, 31-34% in Glacial Acetic Acid**	500 ml	1180
		2500 ml	5830
99790 [7647-01-0]	**Hydrochloric Acid 0.1N Aq. Solution**	500 ml	170
		1000 ml	310
62889 [7647-01-0]	**Hydrochloric Acid 1N Aq. Solution**	500 ml	320
		1000 ml	580
30988 [7647-01-0]	**Hydrochloric Acid 0.5N Aq. Solution**	500 ml	190
		1000 ml	350
87039 [7647-01-0]	**Hydrochloric Acid 2N Aq. Solution**	500 ml	350
		1000 ml	660
80096 [7647-01-0]	**Hydrochloric Acid 3N Aq. Solution**	500 ml	380
		1000 ml	700
65955 [7647-01-0]	**Hydrochloric Acid 4N Aq. Solution**	500 ml	410
		1000 ml	770
72775 [7647-01-0]	**Hydrochloric Acid 4M Solution in Dioxane**	1000 ml	1470
		2500 ml	3640
34472 [7647-01-0]	**Hydrochloric Acid 5N Aq. Solution**	500 ml	450
		1000 ml	810
17560 [7647-01-0]	**Hydrochloric Acid 6N Aq. Solution**	500 ml	470
		1000 ml	860
79661 [7647-01-0]	**Hydrochloric Acid 0.05N Aq. Solution**	500 ml	210
		1000 ml	380
40918 [7647-01-0]	**Hydrochloric Acid 0.01M Aq. Solution**	500 ml	210
		1000 ml	380
56059 [5470-11-1]	**Hydroxylamine Hydrochloride pure, 95%**	100 Gms	250
		500 Gms	730
		5 Kg	6720
66164 [5470-11-1]	**Hydroxylamine Hydrochloride ACS, 96%**	100 Gms	310
		500 Gms	900
		5 Kg	7840
65508 [5470-11-1]	**Hydroxylamine Hydrochloride ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	340
		500 Gms	1010
		5 Kg	8960
47748 [7803-49-8]	**Hydroxylamine 50% Soln. in Water**	50 ml	2240
		100 ml	3590
		500 ml	12320
13877 [10039-54-0]	**Hydroxylamine Sulphate pure, 99%**	100 Gms	240
		500 Gms	620
83032 [10039-54-0]	**Hydroxylamine Sulphate ACS, 99%**	100 Gms	270
		500 Gms	730
44439 [622-40-2]	**4-(2-Hydroxyethyl)morpholine (4-HEM) extrapure, 99%** (2-Morpholinoethanol, 4-Morpholineethanol)	100 ml	440
		500 ml	1600

Code	Name	Packing	Price
91113 [148-24-3]	**8-Hydroxyquinoline (8-Quinolinol) extrapure, 99.5%**	100 Gms	1010
		500 Gms	3920
23264 [148-24-3]	**8-Hydroxyquinoline (8-Quinolinol) ACS, 99.5%**	100 Gms	1100
		500 Gms	4480
46528 [148-24-3]	**8-Hydroxyquinoline (8-Quinolinol) PrimeDx, 99.5%**	25 Gms	14560
		100 Gms	53760
26356 [6303-21-5]	**Hypophosphorous Acid pure, 30-32% in water** (Phosphinic Acid)	500 ml	1350
92379 [6303-21-5]	**Hypophosphorous Acid extrapure AR, 50% in water** (Phosphinic Acid)	500 ml	1570
64556 [8029-68-3]	**Ichthammol pure, 50%**	500 Gms	5600
		1 Kg	10080
87537 [8029-68-3]	**Ichthammol ExiPlus, Multi-Compendial, 50%** (Ammonium Ichthosulfonate, Sulfonated Bitumen, Ictiol, Ichthyol, Bituminol) meets compendial specs of USP, BP, Ph.Eur	500 Gms	8740
		1 Kg	16800
65712 [288-32-4]	**Imidazole pure, 99%**	100 Gms	450
		500 Gms	1240
		5 Kg	11500
		25 Kg	52500
32822 [288-32-4]	**Imidazole ACS, 99%**	100 Gms	1120
		500 Gms	4040
		5 Kg	36630
		25 Kg	163520
78496 [288-32-4]	**Imidazole ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	1240
		500 Gms	4040
61510 [288-32-4]	**Imidazole for molecular biology, 99.5%**	25 Gms	1010
		100 Gms	3140
		500 Gms	13440
10722 [288-32-4]	**Imidazole for tissue culture, 99.5%**	100 Gms	900
73017 [8000-27-9]	**Immersion Oil Synthetic for microscopy** (Cedarwood Oil Synthetic, Immersion Oil Thickened)	30 ml	350
		125 ml	900
		500 ml	1900
		5 Ltr	18500
63214	**Indium AAS Standard Solution in 1N HCl (Traceable to NIST)**	100 ml	1800
		250 ml	2900
88318 [87-51-4]	**Indole-3-Acetic Acid (IAA) pure, 98%**	5 Gms	470
		25 Gms	1510
		100 Gms	5300
96788 [87-51-4]	**Indole-3-Acetic Acid (IAA) for tissue culture, 98%**	5 Gms	540
		25 Gms	2160
		100 Gms	6480
18676 [771-51-7]	**Indole-3-Acetonitrile pure, 98%**	5 Gms	1840
		25 Gms	6810
40407 [133-32-4]	**Indole-3-Butyric Acid (IBA) pure, 98%**	5 Gms	400
		25 Gms	1200
		100 Gms	4200
60289 [133-32-4]	**Indole-3-Butyric Acid (IBA) for tissue culture, 98%**	5 Gms	500
		25 Gms	1400
		100 Gms	4500
35095 [771-50-6]	**Indole-3-Carboxylic Acid pure, 99%**	5 Gms	1370
		25 Gms	5440
67892 [7790-99-0]	**Iodine Monochloride extrapure, 98%**	25 Gms	900
		100 Gms	2800
		250 Gms	6200
18685 [7790-99-0]	**Iodine Monochloride ACS, 98%**	25 Gms	1200
		100 Gms	3200
		250 Gms	7000
66098 [7790-99-0]	**Iodine Monochloride Solution (1M in Acetic Acid) (Wij`s Solution) extrapure** (Chloroiodide Solution)	250 ml	2100
		500 ml	3900
40415 [12029-98-0]	**Iodine Pentoxide (Iodopentoxide) extrapure AR, 99%**	25 Gms	1950
		100 Gms	6380
52539 [7553-56-2]	**Iodine Resublimed ACS, 99.8%**	100 Gms	4000
		500 Gms	18000
37524 [7553-56-2]	**Iodine Resublimed ExiPlus, Multi-Compendial, 99.8%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	4500
		500 Gms	19000
68017 [7553-56-2]	**Iodine 0.01N Solution**	500 ml	550
		1000 ml	1030
59802 [7553-56-2]	**Iodine 1N Solution**	100 ml	1600
		500 ml	7500
12445 [7553-56-2]	**Iodine 0.1N Solution**	500 ml	1200
84644 [865-44-1]	**Iodine Trichloride extrapure, 97%**	25 Gms	3780
36388 [75-47-8]	**Iodoform pure, 99%**	25 Gms	1600
		100 Gms	5500
		250 Gms	12000
58593	**Iodophor, 1.60% I**	1000 ml	1520
		5 Ltr	6700
70031 [516-12-1]	**N-Iodosuccinimide extrapure, 99%** (1-Iodo-2,5-Pyrrolidinedione)	5 Gms	540
		25 Gms	1620
		100 Gms	6050
22157 [7439-88-5]	**Iridium Metal Powder extrapure, 99.9%**	100 Mg	12000
		500 Mg	48000
		1 Gms	90000
50716 [14996-61-3]	**Iridium Trichloride Hydrate extrapure, 99%**	250 Mg	18000
		1 Gms	62000
		5 Gms	300000
84425	**Iron AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml	1800
		250 ml	2900

Code	Name	Packing	Price
64857 [7439-89-6]	**Iron Metal Powder (Electrolytic), 99.5%, mesh ~200**	500 Gms 1 Kg 5 Kg	710 1100 4800
90354 [91-56-5]	**Isatin ACS, 98%** (2,3-Indolinedione)	25 Gms 100 Gms	900 3000
57020 [123-51-3]	**Isoamyl Alcohol (Isopentyl Alcohol) pure, 98.5%**	500 ml 2500 ml 25 Ltr	400 1840 16680
69931 [123-51-3]	**Isoamyl Alcohol (Isopentyl Alcohol) ACS, 98.5%**	500 ml 2500 ml	440 2020
88719 [123-51-3]	**Isoamyl Alcohol (Isopentyl Alcohol) ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	510 2300
69412 [78-83-1]	**Isobutanol pure, 99%**	500 ml 2500 ml 25 Ltr	420 1680 14170
98646 [78-83-1]	**Isobutanol ACS, 99%**	500 ml 2500 ml	520 2030
34159 [78-83-1]	**Isobutanol GC-HS, 99.9%**	250 ml 500 ml 1000 ml	650 1170 2040
52497 [110-19-0]	**Isobutyl Acetate extrapure, 98%** (Acetic Acid Isobutyl Ester)	500 ml 2500 ml 25 Ltr	640 2880 27600
41867 [78-84-2]	**Isobutyraldehyde extrapure, 99%** (2-Methylpropionaldehyde)	500 ml 2500 ml	510 2250
53502 [540-84-1]	**Isooctane extrapure, 99%**	500 ml 2500 ml 25 Ltr	1010 4740 42620
57075 [540-84-1]	**Isooctane ACS, 99%**	500 ml 2500 ml 25 Ltr	1250 5430 47020
81399 [540-84-1]	**Isooctane ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP	500 ml 2500 ml	1300 5860
68930 [540-84-1]	**Isooctane for HPLC & UV Spectroscopy, 99.5%**	500 ml 1000 ml 2500 ml	1580 2790 6100
67959 [540-84-1]	**Isooctane Fuel Grade for HPLC & UV Spectroscopy, 99.9%** (2,2,4 -Trimethylpentane)	500 ml 2500 ml	2420 10120
15485 [540-84-1]	**Isooctane GC-HS, 99.9%** (2,2,4 -Trimethylpentane)	500 ml	4370
98605 [540-84-1]	**Isooctane Dried, 99.5%, water 0.005%** (2,2,4 -Trimethylpentane)	1000 ml	2370
67800 [67-63-0]	**Isopropanol (IPA) pure, 99%** (Isopropyl alcohol, 2-Propanol)	500 ml 2500 ml 25 Ltr	360 1700 16000
62986 [67-63-0]	**Isopropanol (IPA) ACS, 99.5%** (Isopropyl alcohol, 2-Propanol)	500 ml 2500 ml	380 1850
11768 [67-63-0]	**Isopropanol (IPA) ACS, ExiPlus, Multi-Compendial, 99.5%** (Isopropyl alcohol, 2-Propanol) meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	395 2050
38445 [67-63-0]	**Isopropanol (IPA) for molecular biology, 99.8%** (Isopropyl alcohol, 2-Propanol)	100 ml 500 ml 1000 ml	290 810 1560
10140 [67-63-0]	**Isopropanol (IPA) GC-HS, 99.9%** (Isopropyl alcohol, 2-Propanol)	250 ml 500 ml 1000 ml	370 670 1270
16666 [67-63-0]	**Isopropanol (IPA) Dried, 99.5%, water 0.005%** (Isopropyl alcohol, 2-Propanol)	500 ml 1000 ml	750 1380
44036 [67-63-0]	**Isopropanol (IPA) for HPLC & UV Spectroscopy, 99.8%** (Isopropyl alcohol, 2-Propanol)	500 ml 1000 ml 2.5 Ltr (Wide-Mouth) 2500 ml 4 Ltr (Wide-Mouth)	560 980 3740 2070 5750
48906 [67-63-0]	**Isopropanol (IPA) electronic grade, 99.8%** (Isopropyl alcohol, 2-Propanol)	2500 ml	2760
62616 [67-63-0]	**Isopropanol (IPA) VLSI Semiconductor Grade, 99.8%** (Isopropyl alcohol, 2-Propanol)	1000 ml 2500 ml	9500 20000
88773 [22739-76-0]	**Isopropanol (IPA)-d8 for NMR spectroscopy, 99.5 Atom %D** (Isopropanol-d8, Octadeuteroisopropanol, 2-Propanol-d8)	5 ml	23000
13627 [108-22-5]	**Isopropenyl Acetate extrapure, 99%**	100 ml 500 ml	920 3800
35858 [75-31-0]	**Isopropylamine (MIPA) extrapure, 99%** (Monoisopropylamine)	500 ml 2500 ml 25 Ltr	400 1840 17830
73716 [110-27-0]	**Isopropyl Myristate pure, 90%**	500 ml 2500 ml	870 4140
32009 [110-27-0]	**Isopropyl Myristate ExiPlus, Multi-Compendial, 90%** meets compendial specs of USP, Ph.Eur	500 ml 2500 ml	980 4490
42822 [110-27-0]	**Isopropyl Myristate (0.2um filtered) extrapure AR, 98%**	500 ml 2500 ml	1500 6900
36494 [9002-13-5]	**Jack Bean Meal ex. Urease extrapure** (Crude Urease)	100 Gms 500 Gms	270 1220
20489 [2869-83-2]	**Janus Green B (JG-B), 65%** (Diazin Green, Union Green B)	5 Gms 25 Gms 100 Gms	600 2300 7000
86543 [62851-42-7]	**Jenners Stain (Eosin-methylene blue)**	25 Gms 100 Gms	600 2100

Code	Name	Packing	Price
57080 [1332-58-7]	**Kaolin Light extrapure**	500 Gms 5 Kg	200 750
34469	**Karl Fischer Reagent (Pyridine free, Single Solution)**	500 ml 2 x250 ml	4000 4000
62878 [79-33-4]	**L-Lactic Acid for tissue culture, 88%** (L-(+)-Lactic Acid)	100 ml	690
93762 [79-33-4]	**L-Lactic Acid pure, 88%** (L-(+)-Lactic Acid)	500 ml 2500 ml	760 3220
41242 [79-33-4]	**L-Lactic Acid Reagent Grade, 88%** (L-(+)-Lactic Acid)	500 ml 2500 ml	810 3450
21206 [79-33-4]	**L-Lactic Acid ExiPlus, Multi-Compendial, 88%** (L-(+)-Lactic Acid) meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	870 3800
54492 [79-33-4]	**L-Lactic Acid 0.1N Solution**	500 ml 1000 ml	290 550
41158 [79-33-4]	**L-Lactic Acid 1N Solution**	500 ml 1000 ml	490 900
97513	**Lacto Fuchsin Staining Solution**	100 ml	510
92598	**Lactophenol Cotton Blue Solution**	500 ml	350
20863	**Lactophenol Mountant Solution**	500 ml	380
25957 [10039-26-6, 64044-51-5]	**Lactose Monohydrate ACS**	500 Gms 5 Kg	750 7020
33525 [10039-26-6, 64044-51-5]	**Lactose Monohydrate for tissue culture**	500 Gms 5 Kg	810 7280
67091 [10039-26-6, 64044-51-5]	**Lactose Monohydrate ExiPlus, Multi-Compendial** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	870 7710
19459	**Lanthanum AAS Standard Solution in 2N HNO3 (Traceable to NIST)**	100 ml 250 ml	2400 3600
70514 [143-07-7]	**Lauric Acid pure, C12-99%** (Dodecanoic Acid)	500 Gms 5 Kg	640 5750
90346 [1335-32-6]	**Lead Subacetate extrapure, 33% Pb** (Lead (II) Acetate Basic, Lead Hydroxide Acetate, Horne`s Compound)	500 Gms 5 Kg	1070 10030
86215 ✔ [1335-32-6]	**Lead Subacetate ACS, 33% Pb** (Lead (II) Acetate Basic, Lead Hydroxide Acetate, Horne`s Compound)	500 Gms	1200
72518 [6080-56-4]	**Lead (II) Acetate Trihydrate pure, 99%** (Plumbous Acetate Trihydrate, Lead Diacetate)	500 Gms 5 Kg	580 5400
69586 [10099-74-8]	**Lead AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml 250 ml	1800 2900
69789 [7758-95-4]	**Lead (II) Chloride Anhydrous extrapure, 98%**	500 Gms	800
87373 [1309-60-0]	**Lead (IV) Oxide extrapure, 95%** (Lead Dioxide, Lead Peroxide)	100 Gms 500 Gms	490 2000
60843 [7783-46-2]	**Lead (II) Fluoride pure, 98%** (Plumbous Fluoride)	500 Gms	1200
61913 [51404-69-4]	**Lead (II) Hydroxide Acetate pure, 72-75% Pb**	500 Gms	950
64972 [10099-74-8]	**Lead (II) Nitrate extrapure, 99%** (Plumbous Nitrate, Lead Dinitrate)	500 Gms 5 Kg	600 5500
68030 [10099-74-8]	**Lead (II) Nitrate ACS, 99%** (Plumbous nitrate, lead dinitrate)	500 Gms 5 Kg	750 7200
50314 [1317-36-8]	**Lead (II) Oxide (Lead Monoxide, Litharge) pure, 99%**	500 Gms 1 Kg 5 Kg 25 Kg	600 1100 5000 23000
60795 [7446-14-2]	**Lead (II) Sulphate extrapure, 98%** (Anglesite)	500 Gms	700
18279 [1317-36-8]	**Lead (II) Oxide (Lead Monoxide, Litharge) ACS, 99%**	100 Gms 250 Gms 500 Gms	350 650 1200
68015	**Leishman Staining Solution**	250 ml	300
82184 [12627-53-1]	**Leishman Stain**	25 Gms 100 Gms 500 Gms	660 2000 9000
63007 [61-90-5]	**L-Leucine extrapure CHR, 99%**	100 Gms 500 Gms 1 Kg 5 Kg	650 3000 5500 25000
59178 [61-90-5]	**L-Leucine ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, Ph.Eur	500 Gms 1 Kg 5 Kg	3800 7000 31000
75333 [1305-62-0]	**Lime Water Saturated Solution**	500 ml	280
75105	**Lithium AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml 250 ml	3240 5220
31665 [7550-35-8]	**Lithium Bromide Anhydrous extrapure, 99%**	100 Gms 500 Gms 5 Kg	900 3960 36000
33662 [7550-35-8]	**Lithium Bromide Anhydrous ultrapure, 99.999% , -10mesh**	5 Gms 25 Gms	10370 34560
73154 [554-13-2]	**Lithium Carbonate pure, 98%**	100 Gms 250 Gms 500 Gms 5 Kg	1710 3240 6120 54000
19558 [554-13-2]	**Lithium Carbonate ACS, 99%**	100 Gms 250 Gms 1 Kg	1890 3960 10800
74884 [554-13-2]	**Lithium Carbonate ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP , BP, Ph.Eur	250 Gms	4320

Code	Name	Packing	Price
39692 [7447-41-8]	**Lithium Chloride Anhydrous pure, 99%**	100 Gms	1170
		250 Gms	2570
		500 Gms	4860
		5 Kg	45000
76359 [7447-41-8]	**Lithium Chloride Anhydrous ACS, 99%**	100 Gms	1260
		250 Gms	2700
		1 Kg	10080
		5 Kg	48600
26341 [7447-41-8]	**Lithium Chloride 1M Solution in Acetic Acid**	250 ml	2340
20003 [7447-41-8]	**Lithium Chloride 1M Solution in Ethanol**	250 ml	2700
81034 [7789-24-4]	**Lithium Fluoride pure, 98.5%**	25 Gms	700
		100 Gms	2000
		500 Gms	7500
44318 [1310-65-2]	**Lithium Hydroxide Anhydrous extrapure, 98%**	25 Gms	1800
		100 Gms	4500
		500 Gms	17100
52019 [1310-66-3]	**Lithium Hydroxide Monohydrate pure, 97%**	100 Gms	1260
		500 Gms	5040
		5 Kg	46800
48204 [1310-66-3]	**Lithium Hydroxide Monohydrate ACS, 98%**	100 Gms	1440
		500 Gms	5760
		5 Kg	50400
89966 [1310-66-3]	**Lithium Hydroxide Monohydrate ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, Ph.Eur	100 Gms	1620
		500 Gms	6300
82712 [7790-69-4]	**Lithium Nitrate (High Purity) extrapure AR, 99.5%**	100 Gms	1800
		500 Gms	8000
17864 [7790-69-4]	**Lithium Nitrate (High Purity) extrapure AR, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP	100 Gms	2000
		500 Gms	9000
88181 [10102-25-7]	**Lithium Sulphate Monohydrate ACS, 99%**	100 Gms	1240
		250 Gms	2380
		500 Gms	4180
39478 [10102-25-7]	**Lithium Sulphate Monohydrate ACS, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	1430
		250 Gms	2660
		500 Gms	4560
11373 [12007-60-2]	**Lithium Tetraborate ACS, 98%** (di-Lithium Tetraborate)	25 Gms	800
		100 Gms	2500
		500 Gms	8500
78724	**Lugols Iodine 1% Solution (Conc)** (Iodine-Potassium iodide Solution)	125 ml	580
		500 ml	1730
42603	**Lugols Iodine 5% Solution (Conc)**	125 ml	2070
		500 ml	6210
72349 [7274-88-6]	**D-Lysine Monohydrochloride extrapure, 99%**	5 Gms	3910
		25 Gms	14950

Code	Name	Packing	Price
30163 [56-87-1]	**L-Lysine (free base) Anhydrous extrapure, 98%**	25 Gms	2000
		100 Gms	7000
		500 Gms	30000
22527 [56-87-1]	**L-Lysine (free base) Anhydrous for tissue culture, 98%**	25 Gms	1800
15942 [56-87-1]	**L-Lysine (free base) Anhydrous ExiPlus, Multi-Compendial, 98.5%** meets compendial specs of USP	10 Gms	1600
		100 Gms	8200
85187 [39665-12-8]	**L-Lysine Monohydrate extrapure, 99%**	25 Gms	1970
		100 Gms	6000
		500 Gms	25000
76199 [39665-12-8]	**L-Lysine Monohydrate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP	25 Gms	2300
		100 Gms	6800
		500 Gms	27000
45976 [39665-12-8]	**L-Lysine Monohydrate for molecular biology, 99%**	25 Gms	3000
89361 [657-27-2]	**L-Lysine Monohydrochloride extrapure CHR, 99%**	100 Gms	250
		500 Gms	1070
		5 Kg	9000
77872 [657-27-2]	**L-Lysine Monohydrochloride ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	380
		500 Gms	1310
		5 Kg	11000
35443 [16674-78-5]	**Magnesium Acetate Tetrahydrate extrapure, 98%**	500 Gms	780
		5 Kg	6700
53888 [16674-78-5]	**Magnesium Acetate Tetrahydrate extrapure AR, 99%**	500 Gms	930
		5 Kg	8090
64328 [16674-78-5]	**Magnesium Acetate Tetrahydrate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	1040
50488 [16674-78-5]	**Magnesium Acetate Tetrahydrate for molecular biology, 99%**	100 Gms	470
		500 Gms	1390
52917	**Magnesium AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml	1890
		250 ml	3050
28935 [13446-53-2]	**Magnesium Bromide Hexahydrate pure, 98%**	500 Gms	2210
63818 [39409-82-0]	**Magnesium Carbonate Basic Light, 40-45% MgO** (Magnesium Carbonate Light)	500 Gms	460
		5 Kg	4310
31196 [7786-30-3]	**Magnesium Chloride Anhydrous extrapure, 98%**	100 Gms	900
		500 Gms	1400
		5 Kg	9400
39750 [7786-30-3]	**Magnesium Chloride Anhydrous for tissue culture, 98%**	100 Gms	1300
		500 Gms	3100
69396 [7791-18-6]	**Magnesium Chloride Hexahydrate ACS, 99%**	500 Gms	290
		5 Kg	2200

Code	Name	Packing	Price
54269 [7791-18-6]	**Magnesium Chloride Hexahydrate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP	500 Gms 5 Kg	340 2550
91417 [7791-18-6]	**Magnesium Chloride Hexahydrate for molecular biology, 99.5%**	250 Gms 500 Gms 5 Kg	760 1040 9240
16203 [7791-18-6]	**Magnesium Chloride Hexahydrate for tissue culture, 99.5%**	500 Gms 5 Kg	840 7510
28577 [7791-18-6]	**Magnesium Chloride 0.5M Solution**	500 ml 1000 ml	210 410
23772 [7791-18-6]	**Magnesium Chloride 1M Solution**	500 ml 1000 ml	240 450
58742 [153531-96-5]	**Magnesium Citrate Nonahydrate pure, 97%**	500 Gms	740
25349 [7783-40-6]	**Magnesium Fluoride Hydrate pure, 80%**	500 Gms	860
93838 [1309-42-8]	**Magnesium Hydroxide extrapure, 95%**	500 Gms 5 Kg	370 3420
16626 [1309-42-8]	**Magnesium Hydroxide ExiPlus, Multi-Compendial, 95%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	580 5520
45149 [7439-95-4]	**Magnesium Metal coarse powder, 99.8%, -36 mesh**	100 Gms 500 Gms	370 1260
19684 [7439-95-4]	**Magnesium Metal Ribbon extrapure AR, 99.5%**	25 Gms	360
94410 [7439-95-4]	**Magnesium Metal Turnings, 99.8%**	250 Gms 500 Gms 2.5 Kg	980 1530 7460
77697 [13446-18-9]	**Magnesium Nitrate Hexahydrate ACS, 98%**	500 Gms	450
50014 [7487-88-9]	**Magnesium Sulphate Dried Reagent Grade, 62-70% MgSO4, 99%**	500 Gms 2.5 Kg	600 2760
64187 [7487-88-9]	**Magnesium Sulphate Dried ExiPlus, Multi-Compendial, 62-70% MgSO4, 99%** meets compendial specs of USP, BP	500 Gms 2.5 Kg	650 2990
85611 [10034-99-8]	**Magnesium Sulphate Heptahydrate ACS, 98%**	500 Gms 1 Kg 5 Kg	280 500 2420
26129 [10034-99-8]	**Magnesium Sulphate Heptahydrate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	300 2630
74440 [10034-99-8]	**Magnesium Sulphate Heptahydrate for molecular biology, 99.5%**	500 Gms 1 Kg	950 1680
83890 [10034-99-8]	**Magnesium Sulphate Heptahydrate for tissue culture, 99.5%**	500 Gms 5 Kg	480 4520
90046 [10034-99-8]	**Magnesium Sulphate 0.1M Solution**	500 ml 1000 ml	290 520
10301 [10034-99-8]	**Magnesium Sulphate 0.5M Solution**	500 ml 1000 ml	310 560
29811 [10034-99-8]	**Magnesium Sulphate 1M Solution**	500 ml 1000 ml	330 590
78222 [557-04-0]	**Magnesium Stearate ExiPlus, Multi-Compendial** (Stearic Acid Magnesium Salt) meets compendial specs of USP, BP, Ph.Eur	500 Gms 1 Kg	480 790
57499 [1309-48-4]	**Magnesium Oxide Heavy Exiplus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	800 7500
85850 [1309-48-4]	**Magnesium Oxide Heavy extrapure, 98%**	500 Gms 5 Kg	650 6000
78382 [1309-48-4]	**Magnesium Oxide Light Exiplus, Multi-Compendial, 98%** meets compendial specs of BP, USP meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	900 8500
15547 [1309-48-4]	**Magnesium Oxide Light extrapure, 98%**	500 Gms 5 Kg	750 7000
64741 [2437-29-8]	**Malachite Green Oxalate**	100 Gms 500 Gms	610 2530
97890 [110-16-7]	**Maleic Acid extrapure AR, 99.5%**	100 Gms 500 Gms 2.5 Kg	270 640 2760
68589 [110-16-7]	**Maleic Acid extrapure AR, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur, IP	100 Gms 500 Gms	290 790
17289 [123-33-1]	**Maleic Hydrazide extrapure, 98%** (1,2-Dihydro-3,6-pyridazinedione, 3,6-Dihydroxypyridazine, Maleic Acid Hydrazide, 3,6-Pyridazinediol	100 Gms 1 Kg	1210 10350
24226 [6915-15-7]	**DL-Malic Acid pure, 99%**	500 Gms 5 Kg	600 5520
73638 ✔ [9050-36-6]	**Maltodextrin ExiPlus, Multi-Compendial, 20% De** meets compendial specs of BP	500 Gms 5 Kg 25 Kg	360 2880 12650
90827 ✔ [9050-36-6]	**Maltodextrin for tissue culture, 16.5-19.5% De**	500 Gms 5 Kg	260 2300
59338 [6363-53-7]	**D-Maltose Monohydrate for bacteriology and biochemistry, 92%**	100 Gms 500 Gms 5 Kg	320 1250 12080
42844 [6363-53-7]	**D-Maltose Monohydrate for tissue culture, 92%**	100 Gms 500 Gms 5 Kg	380 1560 14380
17451 [6363-53-7]	**D-Maltose Monohydrate ExiPlus, Multi-Compendial, 92%** meets compendial specs of USP	500 Gms 5 Kg	1500 13800

Code	Name	Packing	Price
97526 [6156-78-1]	**Manganese (II) Acetate Tetrahydrate pure, 97%** (Manganous Acetate Tetrahydrate)	500 Gms	500
78017 [6156-78-1]	**Manganese (II) Acetate Tetrahydrate extrapure AR, 99%** (Manganous Acetate Tetrahydrate)	500 Gms	580
94486	**Manganese AAS Standard Solution in 0.5N HCl (Traceable to NIST)**	100 ml 250 ml	2070 3340
45757 [598-62-9]	**Manganese (II) Carbonate pure, Mn content 42-46%** (Manganous Carbonate)	500 Gms	720
75113 [13446-34-9]	**Manganese (II) Chloride Tetrahydrate extrapure, 97%** (Manganous Chloride Tetrahydrate)	500 Gms 5 Kg	680 6210
12551 [13446-34-9]	**Manganese (II) Chloride Tetrahydrate ACS, 98%** (Manganous Chloride Tetrahydrate)	500 Gms 5 Kg	830 7590
73970 [7439-96-5]	**Manganese Metal Powder pure, 99%, ~325 mesh**	250 Gms	690
58011 [10034-96-5]	**Manganese (II) Sulphate Monohydrate pure, 98%** (Manganous Sulphate Monohydrate)	500 Gms 5 Kg	300 2800
12386 [10034-96-5]	**Manganese (II) Sulphate Monohydrate ACS, 99%** (Manganous Sulphate Monohydrate)	500 Gms 5 Kg	800 7400
78822 [10034-96-5]	**Manganese (II) Sulphate Monohydrate ExiPlus, Multi-Compendial, 99%** (Manganous Sulphate Monohydrate) meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	900 8000
38593 [1313-13-9]	**Manganese (IV) Oxide pure, 80%** (Manganese Dioxide)	500 Gms 5 Kg	500 4000
24248 [69-65-8]	**D-Mannitol extrapure, 98.5%**	500 Gms 1 Kg 5 Kg 25 Kg	940 1800 8700 42500
79887 [69-65-8]	**D-Mannitol ACS, 99%**	250 Gms 500 Gms 5 Kg	520 1020 9100
21304 [69-65-8]	**D-Mannitol ACS, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	1160 10920
23495 [108-78-1]	**Melamine extrapure, 98%** (2,4,6-Triamino-1,3,5-Triazine Monomer, sym-Triaminotriazine)	500 Gms 5 Kg	900 7500
64370 ✔ [2216-51-5]	**Menthol Synthetic pure**	100 Gms 250 Gms 1 Kg	950 2000 7500
69892 [60-24-2]	**2-Mercaptoethanol pure, 99%**	100 ml 500 ml	460 1610
83759 [60-24-2]	**2-Mercaptoethanol for molecular biology, 99%**	100 ml 250 ml 500 ml	810 1500 2070
79306 [1600-27-7]	**Mercuric Acetate pure, 98%** (Mercury (II) Acetate)	25 Gms 100 Gms 250 Gms	4100 15000 35000
80243 [1600-27-7]	**Mercuric Acetate ACS, ExiPlus, Multi-Compendial, 98%** (Mercury (II) Acetate) meets compendial specs of Ph.Eur	25 Gms 100 Gms 250 Gms	4500 17000 37000
58263 [7789-47-1]	**Mercuric Bromide ACS, 99%** (Mercury (II) Bromide)	25 Gms 100 Gms 250 Gms	4800 18000 42000
25699 [7487-94-7]	**Mercuric Chloride extrapure, 98%** (Mercury (II) Chloride)	25 Gms 100 Gms 250 Gms 500 Gms	4000 14000 32500 63000
29828 [7487-94-7]	**Mercuric Chloride extrapure AR, ACS, ExiPlus, Multi-Compendial, 98%** (Mercury (II) Chloride) meets compendial specs of USP	25 Gms 100 Gms 250 Gms 500 Gms	4300 15500 33500 66000
98926 [7774-29-0]	**Mercuric Iodide Red pure, 99%** (Mercury (II) Iodide)	25 Gms 100 Gms 250 Gms	3200 12000 26000
36881 [7774-29-0]	**Mercuric Iodide Red ACS, 99%** (Mercury (II) Iodide)	25 Gms 100 Gms 250 Gms	3400 13000 28000
93383 [7774-29-0]	**Mercuric Iodide Red ExiPlus, Multi-Compendial, 99%** (Mercury (II) Iodide) meets compendial specs of USP, BP, Ph.Eur	25 Gms 100 Gms 250 Gms	3700 14000 30000
68434 [7783-34-8]	**Mercuric Nitrate Monohydrate extrapure, 98%** (Mercury (II) Nitrate Monohydrate)	25 Gms 100 Gms 250 Gms	5100 19000 40000
89168 [7783-34-8]	**Mercuric Nitrate Monohydrate ACS, 98.5%** (Mercury (II) Nitrate Monohydrate)	25 Gms 100 Gms 250 Gms	5500 21000 42000
62847 [21908-53-2]	**Mercuric Oxide Red extrapure, 99%** (Mercury (II) Oxide)	25 Gms 100 Gms 250 Gms	4800 18000 42000
60178 [21908-53-2]	**Mercuric Oxide Red ACS, ExiPlus, Multi-Compendial, 99%** (Mercury (II) Oxide) meets compendial specs of Ph.Eur	25 Gms 100 Gms 250 Gms	5200 20000 45000
86666 [21908-53-2]	**Mercuric Oxide Yellow extrapure, 99%** (Mercury (II) Oxide)	25 Gms 100 Gms 250 Gms 500 Gms	4700 18000 42500 83000

Code	Name	Packing	Price
37053 [21908-53-2]	**Mercuric Oxide Yellow ACS, 99%** (Mercury (II) Oxide)	25 Gms	5200
		100 Gms	20000
		250 Gms	43500
		500 Gms	85000
32543 [7783-35-9]	**Mercuric Sulphate pure, 99%** (Mercury (II) Sulphate)	25 Gms	2900
		100 Gms	11000
		250 Gms	25000
		500 Gms	48000
20846 [7783-35-9]	**Mercuric Sulphate ACS, 99%** (Mercury (II) Sulphate)	25 Gms	3100
		100 Gms	12000
		250 Gms	27000
43302 [7783-35-9]	**Mercuric Sulphate ACS, ExiPlus, Multi-Compendial, 99%** (Mercury (II) Sulphate) meets compendial specs of USP, Ph.Eur	25 Gms	3500
		100 Gms	13000
		250 Gms	29000
48274 [592-85-8]	**Mercuric Thiocyanate extrapure, 98%** (Mercury (II) Thiocyanate)	25 Gms	3900
		100 Gms	14500
		250 Gms	34000
93119 [592-85-8]	**Mercuric Thiocyanate extrapure AR, ExiPlus, Multi-Compendial, 98%** (Mercury (II) Thiocyanate) meets compendial specs of Ph.Eur	25 Gms	4200
		100 Gms	15500
		250 Gms	36000
79301 [592-85-8]	**Mercuric Thiocyanate PrimeDx, 98%**	500 Gms	48000
		1 Kg	90000
47184 [10112-91-1]	**Mercurous Chloride pure, 98%** (Mercury (I) Chloride, Calomel)	25 Gms	4700
		100 Gms	18000
		250 Gms	42000
		500 Gms	80000
37921 [10112-91-1]	**Mercurous Chloride ACS, 99.5%** (Mercury (I) Chloride, Calomel)	25 Gms	5200
		100 Gms	20000
		250 Gms	44000
		500 Gms	82000
36004 [14836-60-3]	**Mercurous Nitrate Dihydrate extrapure, 97%** (Mercury (I) Nitrate Dihydrate)	25 Gms	4600
		100 Gms	17500
		250 Gms	41000
44682 [14836-60-3]	**Mercurous Nitrate Dihydrate extrapure AR, ACS, 98%** (Mercury (I) Nitrate Dihydrate)	25 Gms	4900
		100 Gms	19000
		250 Gms	43000
84229	**Mercury AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml	4900
		250 ml	7200
83568 [10124-48-8]	**Mercury (II) Amidochloride extrapure, 98%** (Ammoniated Mercury)	25 Gms	5000
		100 Gms	17000
		500 Gms	80000
43665 ⓘ [108-67-8]	**Mesitylene extrapure, 98%** (1,3,5-Trimethyl Benzene)	500 ml	1270
		1000 ml	2190
		2500 ml	4830
36193 [141-79-7]	**Mesityl Oxide (MOX) extrapure, 99%** (4-Methyl-3-Penten-2-One)	500 ml	2130
53383 ✔ [7439-97-6]	**Mercury Metal extrapure, 99%**	250 Gms	38000
		500 Gms	70000
86556 [79-41-4]	**Methacrylic Acid extrapure, 99%**	500 ml	810
		2500 ml	3450
87726 [75-75-2]	**Methanesulphonic Acid (MSA) pure, 99%**	500 ml	1700
		2500 ml	7600
		25 Ltr	72000
46989 [124-63-0]	**Methanesulphonyl Chloride (MSCl) pure, 98%**	500 ml	810
		2500 ml	3450
		25 Ltr	31050
46644 [124-63-0]	**Methanesulphonyl Chloride (MSCl) extrapure, 99%**	500 ml	870
		2500 ml	3570
		25 Ltr	33350
59029 [67-56-1]	**Methanol pure, 99%**	500 ml	290
		2500 ml	780
		25 Ltr	6620
65524 [67-56-1]	**Methanol ACS, 99.8%**	500 ml	320
		1000 ml	600
		2500 ml	840
37152 [67-56-1]	**Methanol ExiPlus, Multi-Compendial, 99.8%** meets compendial specs of USP, BP, Ph.Eur	500 ml	370
		1000 ml	630
		2500 ml	950
87655 [67-56-1]	**Methanol GC-HS, 99.9%**	500 ml	810
		1000 ml	1420
74022 [67-56-1]	**Methanol for LC-MS, 99.9%**	1 Ltr (Wide-Mouth)	4500
		2.5 Ltr (Wide-Mouth)	10000
		4 Ltr (Wide-Mouth)	14000
79345 [67-56-1]	**Methanol for HPLC & UV Spectroscopy, 99.8%**	500 ml	430
		2500 ml	1120
		2.5 Ltr (Wide-Mouth)	3740
		4 Ltr (Wide-Mouth)	5750
10036 [67-56-1]	**Methanol electronic grade, 99.8%**	2500 ml	1370
62456 [67-56-1]	**Methanol Dried extrapure AR, 99.8%, water 0.02%**	500 ml	420
		1000 ml	580
		2500 ml	1270
96446 [67-56-1]	**Methanol for molecular biology, 99.5%**	500 ml	1050
		1000 ml	1890
		2500 ml	4040
21944 🌡 [811-98-3]	**Methanol-d4 for NMR spectroscopy, 99.8 Atom %D**	10 Gms	12000
		50 Gms	53000

Code	Name	Packing	Price
19305 [63-68-3]	**L-Methionine extrapure CHR, 99%**	25 Gms	330
		100 Gms	700
		1 Kg	5800
		5 Kg	27000
65167 [109-86-4]	**Methyl Cellosolve extrapure, 99%**	500 ml	430
		2500 ml	2020
		25 Ltr	18480
38200 [109-86-4]	**Methyl Cellosolve ACS, 99.5%**	500 ml	480
		2500 ml	2240
60451 [108-87-2]	**Methyl Cyclohexane extrapure, 99%**	500 ml	440
		2500 ml	1800
46985 [845-10-3]	**Methyl Red Sodium Salt (water soluble) ACS, 95%** (C.I. No. 13020) (Acid Red 2)	25 Gms	300
		500 Gms	5100
38862 [110-26-9]	**N,N-Methylene Bisacrylamide (bis-Acrylamide) pure, 98%**	100 Gms	1040
		500 Gms	2910
		5 Kg	22000
38516 [110-26-9]	**N,N-Methylene Bisacrylamide (bis-Acrylamide) 3x cryst. extrapure AR, 99.5%**	25 Gms	540
		100 Gms	1700
		500 Gms	7020
67320 [110-26-9]	**N,N-Methylene Bisacrylamide 3x cryst. for molecular biology, 99.5%** (Bis-Acrylamide)	25 Gms	580
		100 Gms	1820
		250 Gms	4480
		500 Gms	8590
58779 [110-26-9]	**N,N`-Methylene Bisacrylamide 2% Aq. Solution for molecular biology**	25 ml	1800
		100 ml	5250
60363 [7220-79-3]	**Methylene Blue Trihydrate**	25 Gms	440
		100 Gms	1640
		500 Gms	8000
97866	**Methylene Blue 1% Aq. Solution**	125 ml	150
		250 ml	250
77248 [74-88-4]	**Methyl Iodide (Iodomethane) extrapure, 99%**	100 ml	7210
		250 ml	17280
64889 [108-11-2]	**Methyl Isobutyl Carbinol extrapure, 98%**	500 ml	1190
15926 [108-10-1]	**Methyl Isobutyl Ketone (MIBK) pure, 98%**	500 ml	360
		2500 ml	1520
		25 Ltr	14040
61929 [108-10-1]	**Methyl Isobutyl Ketone (MIBK) ACS, 98.5%**	500 ml	430
		2500 ml	1840
		25 Ltr	16200
10416 [108-10-1]	**Methyl Isobutyl Ketone (MIBK) for UV spectroscopy, 99.5%**	500 ml	1320
13077 [563-80-4]	**Methyl Isopropyl Ketone (MIPK) pure, 98%** (3-Methyl-2-Butanone, Isopropyl Methyl Ketone, Methyl Isopropyl Ketone)	250 ml	900
		1000 ml	3400
44191 [80-62-6]	**Methyl Methacrylate (Stabilized) extrapure, 99%** (Methacrylic Acid Methyl Ester, MMA)	500 ml	720
		2500 ml	2990

Code	Name	Packing	Price
72387 [109-02-4]	**4-Methyl Morpholine pure, 99%** (N-Methyl Morpholine)	250 ml	520
		500 ml	920
		1000 ml	1730
70280 [547-58-0]	**Methyl Orange**	25 Gms	210
		100 Gms	570
		1 Kg	4800
51369 [872-50-4]	**N-Methyl-2-Pyrrolidone (NMP) extrapure, 99%**	500 ml	640
		2500 ml	2880
		25 Ltr	26450
86263 [872-50-4]	**N-Methyl-2-Pyrrolidone (NMP) ACS, 99%**	500 ml	750
		2500 ml	3110
11118 [872-50-4]	**N-Methyl-2-Pyrrolidone (NMP) ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 ml	870
		2500 ml	3450
96644 [872-50-4]	**N-Methyl-2-Pyrrolidone (NMP) Dried, 99%, water 0.005%**	500 ml	1040
25903 [872-50-4]	**N-Methyl-2-Pyrrolidone (NMP) GC-HS, 99.9%**	250 ml	1040
		500 ml	1380
		1000 ml	2530
10331 [872-50-4]	**N-Methyl-2-Pyrrolidone (NMP) for HPLC, 99.5%**	1000 ml	2300
90214 [872-50-4]	**N-Methyl-2-Pyrrolidone (NMP) VLSI Semiconductor Grade, 99.5%**	1000 ml	14950
		2500 ml	32200
89975 [872-50-4]	**N-Methyl-2-Pyrrolidone (NMP) for molecular biology, 99.5%**	250 ml	880
44055 [99-76-3]	**Methyl-4-Hydroxybenzoate (Methyl Paraben) ExiPlus, Multi-Compendial, 98-102%** meets compendial specs of USP, BP	500 Gms	1610
		5 Kg	14490
60905 [5026-62-0]	**Methyl-4-Hydroxybenzoate Sodium Salt (Methyl Paraben Sodium) ExiPlus, Multi-Compendial, 95-102%** meets compendial specs of USP, BP	500 Gms	1500
		5 Kg	13920
49449 [493-52-7]	**Methyl Red extrapure AR**	25 Gms	300
		100 Gms	1100
		1 Kg	8300
40431 [493-52-7]	**Methyl Red ACS, ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	25 Gms	460
		100 Gms	1090
		1 Kg	8300
99983 [96-47-9]	**2-Methyl Tetrahydrofuran pure, 98%**	500 ml	1040
		2500 ml	4140
		25 Ltr	39100
14488 [96-47-9]	**2-Methyl Tetrahydrofuran Dried, 99%, water 0.003%**	100 ml	600
		500 ml	2000
37689 [8004-87-3]	**Methyl Violet**	25 Gms	350
		100 Gms	690
90121 [55-55-0]	**Metol extrapure, 99%** (4-Methylaminophenol Sulphate)	100 Gms	1130
		500 Gms	4830

Code	Name	Packing	Price
32198 [1319-77-3]	**Mixed Cresol pure, 85%** (Mixture of Ortho, Meta & Para Cresol)	500 ml	870
		2500 ml	4030
98254 [308080-99-1]	**Molecular Sieve 3A (0.3nm) (1.6-3mm Beads)**	250 Gms	550
		1 Kg	1530
		5 Kg	6110
72866 [70955-01-0]	**Molecular Sieve 4A (0.4nm) (1.6-3mm Beads)**	250 Gms	550
		1 Kg	1530
		5 Kg	6110
75690 [70955-01-0]	**Molecular Sieve 5A (0.5nm) (1.6-3mm Beads)**	250 Gms	600
		1 Kg	1640
		5 Kg	6540
27764 [63231-69-6]	**Molecular Sieve 13X (2-4mm Beads)**	250 Gms	550
		1 Kg	1530
		5 Kg	6110
38717	**Molybdenum AAS Standard Solution in 0.5N H_2SO_4 (Traceable to NIST)**	100 ml	1800
		250 ml	2900
31444 [1313-27-5]	**Molybdenum Trioxide extrapure, 99.5%**	100 Gms	3600
		250 Gms	8500
		500 Gms	16000
72615 [1313-27-5]	**Molybdenum Trioxide ACS, 99.5%**	100 Gms	3900
		250 Gms	9000
		500 Gms	17000
49664 [7782-91-4]	**Molybdic Acid pure, 85%**	100 Gms	2400
		250 Gms	5400
		500 Gms	10400
98466 [7782-91-4]	**Molybdic Acid extrapure AR, 85%**	100 Gms	2600
		250 Gms	5800
		500 Gms	11000
71772 [7782-91-4]	**Molybdic Acid extrapure AR, ACS, ExiPlus, Multi-Compendial, 85%**	100 Gms	2800
		250 Gms	6200
		500 Gms	12000
	meets compendial specs of USP		
44630 [110-91-8]	**Morpholine extrapure, 99%**	500 ml	560
		2500 ml	2530
		25 Ltr	23000
96740 [110-91-8]	**Morpholine ACS, 99%**	250 ml	370
		500 ml	640
11287 [110-91-8]	**Morpholine ACS, ExiPlus, Multi-Compendial, 99.5%**	500 ml	810
	meets compendial specs of BP, Ph.Eur		
98293 [3051-09-0]	**Murexide (Ammonium Purpurate) extrapure AR , 98%**	5 Gms	300
		25 Gms	1130
		100 Gms	4290
		500 Gms	18970
99571 [3051-09-0]	**Murexide (Ammonium Purpurate) extrapure AR, ACS, ExiPlus, Multi-Compendial, 98%**	5 Gms	360
		25 Gms	1260
		100 Gms	4500
		500 Gms	20790
	meets compendial specs of Ph.Eur		
13280 [91-20-3]	**Naphthalene extrapure AR, ExiPlus, Multi-Compendial, 99.5%**	500 Gms	490
		1 Kg	810
		5 Kg	3680
	meets compendial specs of USP, BP, Ph.Eur		
26320 [91-20-3]	**Naphthalene scintillation grade, 99%**	500 Gms	640
		1 Kg	1040
89740 [86-87-3]	**Naphthalene-1-Acetic Acid (NAA) extrapure, 99%** (1-Naphthylacetic Acid)	25 Gms	250
		100 Gms	740
		500 Gms	3260
75928 [86-87-3]	**Naphthalene-1-Acetic Acid (NAA) for tissue culture, 99%** (1-Naphthylacetic Acid)	25 Gms	300
		100 Gms	810
		500 Gms	3770
84226 [90-15-3]	**a-Naphthol pure, 99%**	100 Gms	350
		500 Gms	1270
36309 [90-15-3]	**a-Naphthol Reagent Grade, 99.5%**	25 Gms	250
		100 Gms	600
		500 Gms	2300
30585 [90-15-3]	**a-Naphthol ExiPlus, Multi-Compendial, 99.5%**	100 Gms	720
		500 Gms	2530
	meets compendial specs of USP, BP, Ph.Eur		
76018 [135-19-3]	**b-Naphthol pure, 99%**	500 Gms	580
74770 [135-19-3]	**b-Naphthol extrapure AR, 99%**	100 Gms	370
		500 Gms	1380
19261 [135-19-3]	**b-Naphthol extrapure AR, ExiPlus, Multi-Compendial, 99%**	100 Gms	430
		500 Gms	1610
	meets compendial specs of USP, BP, Ph.Eur		
62718 ⓘ [19381-50-1]	**Naphthol Green B (C.I. 10020) (water soluble) for microscopy** (Green Acid 1, Naphthol Green B (W.S.) (M.S.), Acid Green 1)	25 Gms	290
		100 Gms	900
		500 Gms	4030
77807 [145-50-6]	**1-Naphtholbenzein extrapure** (a-Naphtholbenzein, p-naphtholbenzein, 4,4'-(a-Hydroxybenzylidene)di-1-naphthol)	5 Gms	850
		25 Gms	3800
36059 🌡 [7681-93-8]	**Natamycin (NAT) extrapure, 90-102%** (Pimaricin)	100 Mg	2500
		1 Gms	6000
86079	**Neissers Stain A (Methylene Blue)**	125 ml	230
20379	**Neissers Stain B (Crystal Violet)**	125 ml	250
23058	**Neissers Stain C (Chrysoidine)**	125 ml	230
36248 [553-24-2]	**Neutral Red (High Purity), 85%**	10 Gms	1100
		25 Gms	2400
		100 Gms	8800
90691 [553-24-2]	**Neutral Red extrapure AR, 60%** (C.I. No. 50040)	10 Gms	500
		25 Gms	900
		100 Gms	3300
		500 Gms	12900

Code	Name	Packing	Price
18559 [553-24-2]	**Neutral Red (High Purity) ACS, ExiPlus, Multi-Compendial, 85%** meets compendial specs of BP	10 Gms 25 Gms	1150 2650
73342 [12635-27-7]	**Nickel Aluminium Alloy powder, 50% Ni, 50% Al**	100 Gms 500 Gms	870 3450
48615	**Nickel AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml 250 ml	1800 2900
42591 [6018-89-9]	**Nickel (II) Acetate Tetrahydrate pure, 98%** (Nickelous Acetate Tetrahydrate)	100 Gms 500 Gms	800 3300
66489 [958638-02-3]	**Nickel (II) Carbonate Basic Hydrate extrapure, 44-47% Ni** (Nickelous Carbonate Hydroxide Hydrate, Nickel Hydroxide Carbonate Hydrate)	500 Gms	5400
58481 [12607-70-4, 39430-27-8]	**Nickel (II) Carbonate Basic Tetrahydrate extrapure, 44-47% Ni** (Nickelous Carbonate Hydroxide Tetrahydrate, Nickel Hydroxide Carbonate Tetrahydrate)	500 Gms 5 Kg	3200 28900
84421 [7791-20-0]	**Nickel (II) Chloride Hexahydrate pure, 98%** (Nickelous Chloride Hexahydrate)	500 Gms 5 Kg	1700 16000
27303 [7791-20-0]	**Nickel (II) Chloride Hexahydrate ACS, 98%** (Nickelous Chloride Hexahydrate)	100 Gms 500 Gms	450 1900
10081 [7791-20-0]	**Nickel (II) Chloride Hexahydrate ExiPlus, Multi-Compendial, 98%** (Nickelous Chloride Hexahydrate) meets compendial specs of BP	100 Gms 500 Gms	500 2100
75546 [7440-02-0]	**Nickel Metal powder extrapure, 99.5%, ~100 mesh**	100 Gms 500 Gms	1400 5100
57709 [13478-00-7]	**Nickel (II) Nitrate Hexahydrate pure, 98%** (Nickelous Nitrate Hexahydrate)	500 Gms 5 Kg	1400 11600
30795 [13478-00-7]	**Nickel (II) Nitrate Hexahydrate extrapure AR, 99%** (Nickelous Nitrate Hexahydrate)	250 Gms 500 Gms	900 1600
33568 [13478-00-7]	**Nickel (II) Nitrate Hexahydrate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%** (Nickelous Nitrate Hexahydrate) meets compendial specs of Ph.Eur	250 Gms 500 Gms	1000 1800
49487 [1313-99-1]	**Nickel (II) Oxide pure, 77% Ni** (Nickel Monoxide)	100 Gms 500 Gms	2100 9100
12125 [10101-98-1]	**Nickel (II) Sulphate Heptahydrate pure, 98%** (Nickelous Sulphate Heptahydrate)	500 Gms 5 Kg	1150 11000
93377 [10101-98-1]	**Nickel (II) Sulphate Heptahydrate extrapure AR, 99%** (Nickelous Sulphate Heptahydrate)	500 Gms 5 Kg	1250 12000
48105 [10101-98-1]	**Nickel (II) Sulphate Heptahydrate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%** (Nickelous Sulphate Heptahydrate) meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	780 1350

Code	Name	Packing	Price
61397 [10101-97-0]	**Nickel (II) Sulphate Hexahydrate ACS, 99%** (Nickelous Sulphate Hexahydrate)	100 Gms 500 Gms 5 Kg	600 1800 15000
86097 [59-67-6]	**Nicotinic Acid (Pyridine-3-Carboxylic Acid) pure, 99%** (3-Picolinic Acid, Niacin, Vitamin B3)	100 Gms 500 Gms	490 1440
67066 [8005-03-6]	**Nigrosin (water soluble)**	100 Gms 500 Gms	600 2760
75529 [11099-03-9]	**Nigrosin (alcohol soluble)**	100 Gms 500 Gms	610 2880
71131	**Nigrosin 10% w/v Solution**	100 ml	230
36317 [2381-85-3]	**Nile Blue Chloride**	25 Gms 100 Gms	1270 5750
62887 [3625-57-8]	**Nile Blue Sulphate**	10 Gms 25 Gms 100 Gms	940 1890 6900
47353 [7385-67-3]	**Nile Red** (Nile Blue A Oxazone)	100 Mg 1 Gms 5 Gms	3000 16000 30000
55058 [485-47-2]	**Ninhydrin ACS**	10 Gms 25 Gms 100 Gms	700 1600 6200
93940 [485-47-2]	**Ninhydrin ExiPlus, Multi-Compendial** meets compendial specs of USP, BP, Ph.Eur	10 Gms 25 Gms 100 Gms	740 1700 6500
48081 [139-13-9]	**Nitrilotriacetic Acid (NTA) ACS, 98%** (N,N-Bis(carboxymethyl)glycine, Tris(carboxymethyl)amine)	100 Gms 500 Gms 5 Kg	300 800 7000
15319 [88-74-4]	**o-Nitroaniline extrapure, 99%**	250 Gms	860
34003 [99-09-2]	**m-Nitroaniline extrapure, 99%**	250 Gms	2020
99409 [100-01-6]	**p-Nitroaniline extrapure, 98%**	250 Gms	780
69114 [100-01-6]	**p-Nitroaniline (High Purity) extrapure AR, 99.5%**	25 Gms 100 Gms	400 1020
55410 [98-95-3]	**Nitrobenzene pure, 98%**	500 ml 2500 ml	400 1730
75010 [98-95-3]	**Nitrobenzene ACS, 99%**	500 ml 2500 ml	440 1960
40208 [98-95-3]	**Nitrobenzene ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	510 2190
79007 [75-52-5]	**Nitromethane pure, 97%** (Nitrocarbol, Nitroparaffin)	500 ml 2500 ml 25 Ltr	870 3800 34500

Code	Name	Packing	Price
20509 [75-52-5]	**Nitromethane ACS, 98%** (Nitrocarbol, Nitroparaffin)	500 ml	1220
		1000 ml	2160
		2500 ml	5020
50337 [100-02-7]	**p-Nitrophenol pure, 98%** (4-Nitrophenol)	100 Gms	340
		500 Gms	1100
25170 [100-02-7]	**p-Nitrophenol (High Purity) extrapure AR, 99.5%**	25 Gms	490
		100 Gms	1370
		500 Gms	6030
59687 [100-02-7]	**p-Nitrophenol (High Purity) extrapure AR, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of BP, Ph.Eur	25 Gms	560
		100 Gms	1510
49301 [100-02-7]	**p-Nitrophenol for Spectrophotometry, 99.5%** (4-Nitrophenol)	5 Gms	2160
		25 Gms	5460
54373 [35192-74-6]	**Nonane Sulphonic Acid Sodium Salt Anhydrous for HPLC, 98%** (Sodium 1-Nonanesulfonate)	5 Gms	2700
		25 Gms	11600
67956 [5324-84-5]	**Octane Sulphonic Acid Sodium Salt Anhydrous for HPLC, 99%**	25 Gms	1300
		100 Gms	4600
		250 Gms	11300
72554 [207596-29-0]	**Octane Sulphonic Acid Sodium Salt Monohydrate for HPLC, 99%**	25 Gms	1300
		100 Gms	4600
		250 Gms	11300
38732 [111-87-5]	**1-Octanol (1-Octyl Alcohol) pure, 99%** (n-Octanol, n-Octyl Alcohol)	500 ml	750
		2500 ml	3500
40334 [111-87-5]	**1-Octanol (1-Octyl Alcohol) for HPLC & UV Spectroscopy, 99%** (n-Octanol, n-Octyl Alcohol)	500 ml	2000
		1000 ml	3800
52761 [142-31-4]	**Octyl Sulphate Sodium Salt for HPLC, 99%** (Sodium Octyl Sulphate, SOS)	5 Gms	2000
		25 Gms	9000
87572 [112-80-1]	**Oleic Acid pure, 65%** (cis-9-Octadecenoic acid, Elainic acid)	500 ml	400
		2500 ml	1600
28460 ■ [112-80-1]	**Oleic Acid Reference Standard, 99%(GC)**	1 Gms	1200
		5 Gms	3000
		25 Gms	11000
55081 [112-80-1]	**Oleic Acid extrapure, 98%** (cis-9-Octadecenoic acid, Elainic acid)	25 ml	3000
		100 ml	6800
87519 [112-90-3]	**Oleylamine pure, 95%** (cis-9-Octadecenylamine)	250 ml	390
		1000 ml	1350
		2500 ml	3000
79576 [8001-25-0]	**Olive Oil**	100 ml	500
		250 ml	900
		500 ml	1700
13756 [1936-15-8]	**Orange G**	25 Gms	260
		100 Gms	810
		500 Gms	3600
66366 [504-15-4]	**Orcinol Anhydrous (3,5-Dihydroxytoluene) extrapure, 99%**	5 Gms	2300
		25 Gms	5640
		100 Gms	13920
49330 [6153-39-5]	**Orcinol Monohydrate (3,5-Dihydroxytoluene Monohydrate) extrapure, 99%**	10 Gms	1610
		25 Gms	3450
		100 Gms	12650
36916 [6153-39-5]	**Orcinol Monohydrate (3,5-Dihydroxytoluene Monohydrate) ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP	10 Gms	1840
		25 Gms	3680
		100 Gms	13800
36678 [7664-38-2]	**Orthophosphoric Acid extrapure, 85%** (Phosphoric Acid)	500 ml	800
		2500 ml	3600
		20 Ltr	27000
50056 [7664-38-2]	**Orthophosphoric Acid ACS, 85%** (Phosphoric Acid)	500 ml	900
		2500 ml	4000
78255 [7664-38-2]	**Orthophosphoric Acid ExiPlus, Multi-Compendial, 85%** (Phosphoric Acid) meets compendial specs of USP, BP, Ph.Eur	500 ml	980
		2500 ml	4300
38048 [7664-38-2]	**Orthophosphoric Acid for HPLC & UV Spectroscopy, 85%**	250 ml	1500
		500 ml	2600
		1000 ml	4700
83577 🌡 [20816-12-0]	**Osmic Acid extrapure, 99%**	1 Gms	8050
40510 [6153-56-6]	**Oxalic Acid extrapure AR, 99.5%**	500 Gms	380
		5 Kg	3500
72425 [6153-56-6]	**Oxalic Acid Dihydrate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	430
		5 Kg	4000
64056 [3375-31-3]	**Palladium (II) Acetate (Palladium Diacetate) extrapure, 99%, 47% Pd**	250 Mg	4000
		1 Gms	12000
		5 Gms	55000
		25 Gms	260000
59548 ✔ [14024-61-4]	**Palladium (II) Acetylacetonate (Pd(acac)2) extrapure, 34.7% Pd** (Bis(2,4-Pentanedionato) Palladium(II), Bis(Acetylacetonato) Palladium(II))	250 Mg	3200
		1 Gms	10000
		5 Gms	44000
69263	**Palladium AAS Standard Solution in 0.5N HCl (Traceable to NIST)**	100 ml	3600
		250 ml	5900
33080 [7647-10-1]	**Palladium (II) Chloride pure, 99%, 59-60% Pd** (Palladium Dichloride)	250 Mg	5000
		1 Gms	16000
		5 Gms	72000
		25 Gms	350000
96917 [7647-10-1]	**Palladium (II) Chloride ACS, ExiPlus, Multi-Compendial, 99%, 59-60% Pd** (Palladium Dichloride) meets compendial specs of USP, BP	250 Mg	5500
		1 Gms	17000
		5 Gms	75000
		25 Gms	365000
19694 [32916-07-7]	**Palladium (II) Nitrate Dihydrate extrapure, 99%, 40% Pd** (Palladium Dinitrate Dihydrate)	250 Mg	3500
		1 Gms	11000
		5 Gms	50000
92345 [7440-05-3]	**Palladium 5% On Asbestos**	1 Gms	1400
		10 Gms	12000

Code	Name	Packing	Price
57033 [7440-05-3]	**Palladium 10% On Asbestos**	1 Gms	2400
		10 Gms	23000
44002 [7440-05-3]	**Palladium 5% On Activated Carbon**	1 Gms	1400
		10 Gms	12000
		50 Gms	58000
45685 [7440-05-3]	**Palladium 5% On Activated Carbon (50% wet)**	1 Gms	700
		10 Gms	6000
		50 Gms	29000
52944 [7440-05-3]	**Palladium 10% On Activated Carbon**	1 Gms	2400
		10 Gms	23000
		50 Gms	112000
46494 [7440-05-3]	**Palladium 10% On Activated Carbon (50% wet)**	1 Gms	1400
		10 Gms	12000
		50 Gms	58000
40822 [7440-05-3]	**Palladium 10% on Calcium Carbonate**	1 Gms	2400
		5 Gms	11500
32049 [1314-08-5]	**Palladium (II) Oxide pure, 85% Pd** (Palladium Monoxide)	250 Mg	6000
		1 Gms	20000
81249 [13566-03-5]	**Palladium (II) Sulphate extrapure, 99%**	250 Mg	4200
		1 Gms	13000
18211 [8012-95-1, 8042-47-5]	**Paraffin Liquid Light (Mineral Oil Light) extrapure, 25-80mPas**	500 ml	270
		2500 ml	1150
		25 Ltr	10500
79167 [8012-95-1, 8042-47-5]	**Paraffin Liquid Heavy (Mineral Oil Heavy) extrapure, 110-230mPas**	500 ml	360
		2500 ml	1560
		25 Ltr	15000
86715 [8002-74-2]	**Paraffin Wax Cubes (Type 1 - 56-58) extrapure**	500 Gms	430
		1 Kg	790
47994 [8002-74-2]	**Paraffin Wax Pellets (Type 2 - 58-60) extrapure**	500 Gms	430
		1 Kg	790
84381 [8002-74-2]	**Paraffin Wax Pellets (Type 3 - 60-62) extrapure**	500 Gms	430
		1 Kg	790
74193 [8002-74-2]	**Paraffin Wax White Soft pure** (Petroleum Jelly White)	100 Gms	180
		500 Gms	440
66400 [8002-74-2]	**Paraffin Wax w/ Ceresin Blocks, 58-60°C**	500 Gms	800
		1 Kg	1400
		2 Kg	2700
60536 [8002-74-2]	**Paraffin Yellow pure**	500 Gms	720
94464 [30525-89-4]	**Paraformaldehyde (PFA) extrapure, 96%** (Polyoxymethylene, Paraform, p-Formaldehyde)	500 Gms	500
		5 Kg	4630
45057 [9000-69-5]	**Pectin pure** (Poly-D-Galacturonic Acid Methyl Ester)	100 Gms	950
		500 Gms	3800
		1 Kg	6440
		5 Kg	28750
96096 [109-66-0]	**n-Pentane extrapure, 99%**	500 ml	810
		2500 ml	3680
		25 Ltr	33350
41730 [109-66-0]	**n-Pentane Reagent Grade, 99%**	500 ml	920
		2500 ml	4030
94510 [109-66-0]	**n-Pentane ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 ml	1040
		2500 ml	4260
49135 [109-66-0]	**n-Pentane for HPLC & UV Spectroscopy, 99%**	500 ml	1380
		1000 ml	2530
		2500 ml	5750
47040 [5343-92-0]	**1,2-Pentanediol pure, 98%** (1,2-Dihydroxypentane)	25 ml	1610
		250 ml	9090
18486 [22767-49-3]	**Pentane Sulphonic Acid Sodium Salt Anhydrous for HPLC, 99%**	25 Gms	1300
		100 Gms	4600
		250 Gms	11300
90760 [207605-40-1]	**Pentane Sulphonic Acid Sodium Salt Monohydrate for HPLC, 99%**	25 Gms	1300
		100 Gms	4600
		250 Gms	11300
23283 [7601-90-3]	**Perchloric Acid ACS, for Diamond Industry, 70%** Suitable fo Diamond Industry	500 ml	2530
26312 [10450-60-9]	**Periodic Acid extrapure, 99%**	25 Gms	1270
		100 Gms	4030
		500 Gms	18400
11639 [10450-60-9]	**Periodic Acid ACS, 99%**	25 Gms	1440
		100 Gms	5180
11778 [8032-32-4]	**Petroleum Ether 35-60 extrapure**	500 ml	460
		2500 ml	1840
26440 [8032-32-4]	**Petroleum Ether 35-60 ACS**	500 ml	580
		2500 ml	2190
65643 [8032-32-4]	**Petroleum Ether 40-60 ExiPlus, Multi-Compendial** meets compendial specs of USP, BP, Ph.Eur	500 ml	610
		2500 ml	2420
58415 [8032-32-4]	**Petroleum Ether 60-80 pure**	500 ml	350
		2500 ml	1600
		25 Ltr	15000
36651 [8032-32-4]	**Petroleum Ether 60-80 Reagent Grade**	500 ml	380
		2500 ml	1750
15340 [8032-32-4]	**Petroleum Ether 60-80 for HPLC & UV Spectroscopy**	500 ml	640
		1000 ml	1150
		2500 ml	2300
15590 [8032-32-4]	**Petroleum Ether 80-100 pure**	500 ml	490
		2500 ml	1950
		25 Ltr	17830
52226 [8032-32-4]	**Petroleum Ether 100-120 extrapure**	500 ml	510
		2500 ml	1840
69031 [136112-00-0]	**Phenol:Chloroform:Isoamyl Alcohol (25:24:1) pH 8.0 for molecular biology**	100 ml	1610
		250 ml	3680
		500 ml	4950

Code	Name	Packing	Price
17080	**Phenol:Chloroform:Isoamyl Alcohol (125:24:1) pH 4.5 for molecular biology**	100 ml 500 ml	2990 10010
59606 [617-65-2]	**DL-Glutamic Acid extrapure, 98%** (DL-a-Aminoglutaric Acid)	25 Gms 100 Gms 500 Gms	1050 2590 8500
79030	**Phenol:Chloroform:Isoamyl Alcohol (49.5:49.5:1) pH 6.7 for molecular biology**	100 ml 500 ml	3110 10470
36976	**Phenol:Chloroform: Isoamyl Alcohol (49.5:49.5:1) pH 8.0 for molecular biology**	100 ml 500 ml	3220 11270
28104 [108-95-2]	**Phenol Crystalline pure, 99%** (Hydroxybenzene, Carbolic Acid)	500 Gms 2.5 Kg 25 Kg	730 3310 27010
14892 [108-95-2]	**Phenol Crystalline ACS, 99%** (Hydroxybenzene, Carbolic Acid)	500 Gms 2.5 Kg	820 3590
97745 [108-95-2]	**Phenol Crystalline ExiPlus, Multi-Compendial, 99.5%** (Hydroxybenzene, Carbolic Acid) meets compendial specs of USP, BP, Ph.Eur	500 Gms 2.5 Kg	900 3870
17286 [108-95-2]	**Phenol Crystalline for molecular biology, 99.5%** (Hydroxybenzene, Carbolic Acid)	100 Gms 500 Gms	580 2160
47484 [108-95-2]	**Phenol Equilibrated with 0.1M Citrate Buffer pH 4.5 for molecular biology w/ Stabilizer**	100 ml 500 ml	1350 5410
12692 [108-95-2]	**Phenol Tris Equilibrated for molecular biology w/ Stabilizer** (Hydroxybenzene, Carbolic Acid)	60 ml 100 ml 500 ml	1750 2370 6210
78348 [108-95-2]	**Phenol Tris Equilibrated for molecular biology w/o Stabilizer**	100 ml 500 ml	2370 6210
83275 [108-95-2]	**Phenol Saturated w/ 10% water for molecular biology, 90%** (Phenol Liquid w/ 10% water)	100 ml 500 ml	1130 3620
90179 [77-09-8]	**Phenolphthalein ACS, 98%**	50 Gms 100 Gms 500 Gms	430 640 3060
74785 [77-09-8]	**Phenolphthalein extrapure ExiPlus, Multi-Compendial, 98%** meets compendial specs of BP, Ph.Eur	100 Gms 500 Gms	770 3260
22757 [77-09-8]	**Phenolphthalein Indicator 0.5% Soln in Ethanol** (Phthalin;3,3-bis(4-Hydroxyphenyl) -2-Benzofuran-1-one)	125 ml 500 ml	180 460
79563 [77-09-8]	**Phenolphthalein Indicator 1% Soln in Ethanol**	125 ml 500 ml	220 520
90817 [143-74-8]	**Phenol Red**	25 Gms 100 Gms 500 Gms	670 2600 11880
68313 [143-74-8]	**Phenol Red ACS, ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	25 Gms 100 Gms 500 Gms	760 2810 12960
59462 [34487-61-1]	**Phenol Red Sodium Salt (water soluble) ACS, 95%** (Phenolsulfonephthalein Sodium Salt)	5 Gms 25 Gms 100 Gms	200 600 2200
40919 [673-06-3]	**D-Phenylalanine extrapure CHR, 99%**	5 Gms 25 Gms 100 Gms	440 1200 3000
67602 [150-30-1]	**DL-Phenylalanine extrapure CHR, 99%**	25 Gms 100 Gms	710 2410
85081 [63-91-2]	**L-Phenylalanine extrapure CHR, 99%**	100 Gms 500 Gms 1 Kg 5 Kg	500 2200 3800 16000
40170 [63-91-2]	**L-Phenylalanine ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP	100 Gms 500 Gms 1 Kg	650 2500 4600
72373 [106-50-3]	**p-Phenylenediamine pure, 98%** (1,4-Benzenediamine, 1,4-Diaminobenzene, 1,4-Phenylenediamine)	250 Gms 1 Kg	920 3220
58620 [60-12-8]	**2-Phenylethanol (PEA) extrapure, 99%** (Benzyl Carbinol, Phenethyl Alcohol, Phenylethyl Alcohol solution)	500 ml 2500 ml	2070 6900
64975 [100-63-0]	**Phenylhydrazine pure, 98%**	100 ml 500 ml 2500 ml	810 3220 16330
31804 [100-63-0]	**Phenylhydrazine Reagent Grade, 99%**	100 ml 500 ml	1380 6330
40555 [59-88-1]	**Phenylhydrazine Hydrochloride pure, 99%**	100 Gms 250 Gms	920 2070
12109 [59-88-1]	**Phenylhydrazine Hydrochloride extrapure AR, 99%**	100 Gms 250 Gms	1270 2990
26892 [59-88-1]	**Phenylhydrazine Hydrochloride extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 250 Gms	1610 3800
64732 [59-88-1]	**Phenylhydrazine Hydrochloride Superior extrapure AR, 99%**	5 Gms 100 Gms	1150 2880
28022 [138-24-9]	**Phenyltrimethyl Ammonium Chloride (PTMAC) pure, 97%**	100 Gms 500 Gms	810 3450
82322 [51429-74-4]	**Phosphomolybdic Acid ACS**	25 Gms 100 Gms 500 Gms	1300 4100 18000
78861 [13598-36-2]	**Phosphorous Acid Crystals (Phosphonic Acid) extrapure, 98%** (Orthophosphorous Acid)	500 Gms 2.5 Kg	600 2500
72656 [7723-14-0]	**Phosphorus Red pure 97%**	500 Gms	2200
80653 [12501-23-4]	**Phosphotungstic Acid Hydrate extrapure** (Tungstophosphoric Acid)	25 Gms 100 Gms 500 Gms	1500 5300 25000

Code	Name	Packing	Price
57855 [12501-23-4]	**Phosphotungstic Acid Hydrate extrapure AR**	25 Gms 100 Gms 500 Gms	1650 5700 27000
67397 [88-99-3]	**Phthalic Acid pure, 99%**	500 Gms	900
32624 [85-44-9]	**Phthalic Anhydride pure, 97%**	500 Gms 5 Kg	400 3000
93665 [85-44-9]	**Phthalic Anhydride ACS, 99%**	500 Gms 5 Kg	800 6900
44745 [17211-15-3]	**Phytic Acid Sodium Salt Hydrate (InsP6) extrapure, 70%** (Inositol hexaphosphoric acid dodecasodium salt, myo-Inositol hexakis(dihydrogen phosphate) dodecaso	5 Gms 25 Gms	3800 16700
49298 [83-86-3]	**Phytic Acid (InsP6) pure, 92%** (Inositol Hexaphosphoric Acid, Myo-Inositol Hexakisphosphate, Inositol Hexaphosphate)	5 Gms 25 Gms 100 Gms	2100 4000 14000
94847 [88-89-1]	**Picric Acid moistened with ~30% water extrapure, 99%**	100 Gms 500 Gms	440 1100
47998 [88-89-1]	**Picric Acid moistened with ~30% water ACS , 99.8%**	100 Gms 500 Gms	550 1320
45808 [88-89-1]	**Picric Acid moistened with ~30% water PrimeDx, 99.8%**	5 Kg 25 Kg	9500 45000
59828 [10025-65-7]	**Platinum (II) Chloride extrapure, 73% Pt** (Platinum Dichloride)	250 Mg 1 Gms	9000 30000
54489 [13454-96-1]	**Platinum (IV) Chloride extrapure, 99% Pt** (Platinum Tetrachloride, Platinic Chloride)	250 Mg 1 Gms	5300 15000
14184 [1314-15-4]	**Platinum (IV) Oxide Anhydrous extrapure, 84-86% Pt**	250 Mg 1 Gms 5 Gms	12000 35000 140000
16462 [52785-06-5]	**Platinum (IV) Oxide Hydrate extrapure, 76-80% Pt** Platinum Dioxide Hydrate, Adam's Catalyst	250 Mg 1 Gms	9000 28000
41175 [9000-55-9]	**Podophyllum Resin pure, 50%**	25 Gms 100 Gms	5640 19320
79133 [25322-68-3]	**Polyethylene Glycol 400 (PEG 400)** (Carbowax 400)	500 Gms 5 Kg	520 4600
87919 [25322-68-3]	**Polyethylene Glycol 400 (PEG 400) ExiPlus, Multi-Compendial** (Carbowax 400) meets compendial specs of BP	500 Gms 5 Kg	640 5750
16053 [25322-68-3]	**Polyethylene Glycol 1000 (PEG 1000)** (Poly(ethylene glycol), Polyethylene glycol 1000, Polyglycol, Polyethylene oxide, Polyoxyethylene)	500 Gms 5 Kg	670 6040
61213 [25322-68-3]	**Polyethylene Glycol 3350 Flakes (PEG 3350) ExiPlus, Multi-Compendial** (Carbowax 3350)	500 Gms 5 Kg	950 8740
48759 [25322-68-3]	**Polyethylene Glycol 3350 Powder (PEG 3350)** (Carbowax 3350)	500 Gms 5 Kg	830 7590
30970 [25322-68-3]	**Polyethylene Glycol 4000 Flakes (PEG 4000) Flakes** (Carbowax 4000)	500 Gms 5 Kg 25 Kg	580 5590 26500
76937 [25322-68-3]	**Polyethylene Glycol 4000 Flakes (PEG 4000) Flakes ExiPlus, Multi-Compendial** (Carbowax 4000) meets compendial specs of USP	500 Gms 5 Kg	690 6470
35207 [25322-68-3]	**Polyethylene Glycol 4000 Powder (PEG 4000)** (Carbowax 4000)	500 Gms 5 Kg	670 6210
67429 [25322-68-3]	**Polyethylene Glycol 6000 Flakes (PEG 6000)** (Carbowax 6000)	500 Gms 5 Kg	640 5520
26029 [25322-68-3]	**Polyethylene Glycol 6000 Flakes (PEG 6000) ExiPlus, Multi-Compendial** (Carbowax 6000) meets compendial specs of USP	500 Gms 5 Kg	750 7020
49194 [25322-68-3]	**Polyethylene Glycol 6000 Powder (PEG 6000)** (Carbowax 6000)	500 Gms 5 Kg	690 6330
54866 [25322-68-3]	**Polyethylene Glycol 8000 Flakes (PEG 8000) for molecular biology** (Carbowax 8000)	500 Gms 5 Kg	2020 18190
94701 [25322-68-3]	**Polyethylene Glycol 8000 Flakes (PEG 8000) ExiPlus, Multi-Compendial** (Carbowax 8000)	500 Gms 5 Kg	880 8200
44766 [25322-68-3]	**Polyethylene Glycol 8000 Powder (PEG 8000) for molecular biology** (Carbowax 8000)	500 Gms 5 Kg	1060 10300
64651 [25322-68-3]	**Polyethylene Glycol 9000 Flakes (PEG 9000)** (Carbowax 9000)	500 Gms 5 Kg	870 7770
46505 [25322-68-3]	**Polyethylene Glycol 9000 Flakes (PEG 9000) ExiPlus, Multi-Compendial** (Carbowax 9000)	500 Gms 5 Kg	1100 10350
22531 [25322-68-3]	**Polyethylene Glycol 9000 Powder (PEG 9000)** (Carbowax 9000)	500 Gms 5 Kg	920 8340
39569 [25322-68-3]	**Polyethylene Glycol 10000 (PEG 10000)** (Polyethylene oxide, Polyoxy ethylene, Carbowax 10000)	500 Gms 5 Kg	2760 25300
75727 [25322-68-3]	**Polyethylene Glycol 20000 Flakes (PEG 20000)** (Carbowax 20000)	100 Gms 500 Gms 1 Kg 5 Kg	920 4030 6900 20700
42248 [25322-68-3]	**Polyethylene Glycol 20000 Powder (PEG 20000)** (Carbowax 20000)	500 Gms 1 Kg	4370 7480
33817 [25322-68-3]	**Polyethylene Glycol 20000 Flakes (PEG 20000) ExiPlus, Multi-Compendial** (Carbowax 20000) meets compendial specs of USP, BP	100 Gms 500 Gms 1 Kg	1380 4600 8050

Code	Name	Packing	Price
41098 [8017-16-1]	**Polyphosphoric Acid extrapure, 83-87% P2O5**	500 ml 2500 ml	950 4000
28599 [9005-64-5]	**Polysorbate 20 (Tween 20) extrapure**	500 Gms 5 Kg 25 Kg	980 8050 38360
33712 [9005-64-5]	**Polysorbate 20 (Tween 20) for tissue culture**	500 Gms	1340
23610 [9005-64-5]	**Polysorbate 20 (Tween 20) ExiPlus, Multi-Compendial** meets compendial specs of BP	500 Gms 5 Kg	1200 10500
65296 [9005-64-5]	**Polysorbate 20 (Tween 20) for molecular biology**	100 Gms 500 Gms	1820 2940
54077 [9005-64-5]	**Polysorbate 20 (Tween 20) PrimeDx**	1 Kg 5 Kg	8050 35650
31453 [9005-66-7]	**Polysorbate 40 (Tween 40) extrapure**	250 Gms 5 Kg	710 10930
62539 [9005-66-7]	**Polysorbate 40 (Tween 40) ExiPlus, Multi-Compendial** meets compendial specs of USP, BP, Ph.Eur	250 Gms 5 Kg	1260 14380
12502 [9005-67-8]	**Polysorbate 60 (Tween 60) extrapure**	100 Gms 500 Gms 5 Kg	560 1220 9200
86878 [9005-67-8]	**Polysorbate 60 (Tween 60) ExiPlus, Multi-Compendial** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	1400 10350
28940 [9005-65-6]	**Polysorbate 80 (Tween 80) extrapure**	500 Gms 5 Kg 25 Kg	1050 7680 34870
95188 [9005-65-6]	**Polysorbate 80 (Tween 80) ExiPlus, Multi-Compendial** meets compendial specs of BP	500 Gms 5 Kg 25 Kg	1210 9080 43230
40606 [9005-65-6]	**Polysorbate 80 (Tween 80) PrimeDx**	5 Kg 25 Kg	10350 48300
81613 [9003-39-8]	**Polyvinylpyrrolidone K25 (Povidone, PVP K-25)**	100 Gms 500 Gms 5 Kg	750 2400 21000
41030 [9003-39-8]	**Polyvinylpyrrolidone K30 (Povidone, PVP K-30)**	100 Gms 500 Gms 5 Kg 25 Kg	550 1400 12000 55000
65155 [9003-39-8]	**Polyvinylpyrrolidone K30 (Povidone, PVP K-30) ExiPlus, Multi-Compendial** meets compendial specs of BP	100 Gms 500 Gms 5 Kg	650 1600 13500
24564 [9003-39-8]	**Polyvinylpyrrolidone K90 (Povidone, PVP K-90)**	100 Gms 500 Gms 5 Kg	1100 3300 23000
44333 [127-08-2]	**Potassium Acetate pure, 99%**	500 Gms 5 Kg	600 5500
79156 [127-08-2]	**Potassium Acetate ACS, 99%**	500 Gms 5 Kg	680 6400
40480 [127-08-2]	**Potassium Acetate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	730 6900
96248 [127-08-2]	**Potassium Acetate for molecular biology, 99.5%**	100 Gms 500 Gms	400 1200
81955 [127-08-2]	**Potassium Acetate for tissue culture, 99.5%**	100 Gms 500 Gms	300 900
99644 [7758-01-2]	**Potassium Bromate pure, 99.5%**	500 Gms 5 Kg	1350 12000
80823 [7758-02-3]	**Potassium Bromide pure, 98.5%**	500 Gms 5 Kg	1100 10200
66453 [7758-02-3]	**Potassium Bromide ACS, 99%**	500 Gms 5 Kg	1600 15000
25707 [7758-02-3]	**Potassium Bromide ExiPlus, Multi-Compendial, 98.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	1800 17000
13873 [7758-02-3]	**Potassium Bromide Powder for IR & FTIR spectroscopy, 99.5%**	25 Gms 100 Gms 500 Gms	1200 4270 18000
87005	**Potassium AAS Standard Solution in 0.5N HCl (Traceable to NIST)**	100 ml 250 ml	1800 2900
93649 [7758-02-3]	**Potassium Bromide 0.1M Solution**	500 ml 1000 ml	310 560
79982 [7758-02-3]	**Potassium Bromide 0.5M Solution**	500 ml 1000 ml	330 610
29988 [7758-02-3]	**Potassium Bromide 1M Solution**	500 ml 1000 ml	360 640
24794 [298-14-6]	**Potassium Bicarbonate extrapure, 99%** (Potassium Hydrogen Carbonate)	500 Gms 5 Kg	430 4100
22423 [298-14-6]	**Potassium Bicarbonate ACS, 99.7%** (Potassium Hydrogen Carbonate)	500 Gms 5 Kg	450 4200
57674 [298-14-6]	**Potassium Bicarbonate ExiPlus, Multi-Compendial, 99.7%** (Potassium Hydrogen Carbonate) meets compendial specs of USP, BP	500 Gms 5 Kg	470 4500
87272 [584-08-7]	**Potassium Carbonate anhydrous extrapure, 99%**	500 Gms 5 Kg	440 3810
83342 [584-08-7]	**Potassium Carbonate anhydrous extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP	500 Gms 5 Kg	510 4360
38630 [7447-40-7]	**Potassium Chloride ACS, 99.5%**	500 Gms 5 Kg 25 Kg	310 2290 10690

Code	Name	Packing	Price
50016 [7447-40-7]	**Potassium Chloride ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	340 2620
84984 [7447-40-7]	**Potassium Chloride for molecular biology, 99.5%**	500 Gms 1 Kg	900 1370
71088 [7447-40-7]	**Potassium Chloride for tissue culture, 99.5%**	250 Gms	770
84390 [7789-00-6]	**Potassium Chromate pure, 99%**	500 Gms 5 Kg	920 8400
93448 [7789-00-6]	**Potassium Chromate ACS, 99%**	500 Gms 5 Kg	1000 9000
31751 [7789-00-6]	**Potassium Chromate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	1110
49890 [6100-05-6]	**Potassium Citrate Tribasic Monohydrate extrapure AR, 99%**	500 Gms 5 Kg	550 5250
74176 [6100-05-6]	**Potassium Citrate Tribasic Monohydrate extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP	500 Gms 5 Kg	650 6250
13250 [7778-50-9]	**Potassium Dichromate pure, 99.5%**	500 Gms 5 Kg	950 8600
80233 [7778-50-9]	**Potassium Dichromate ACS, 99.9%**	500 Gms 5 Kg	1050 9100
55353 [7778-50-9]	**Potassium Dichromate ACS, ExiPlus, Multi-Compendial, 99.9%** meets compendial specs of BP	500 Gms 5 Kg	1200 9800
52403 [7778-77-0]	**Potassium Dihydrogen Orthophosphate pure, 99%**	500 Gms 5 Kg 25 Kg	630 5760 26400
50451 [7778-77-0]	**Potassium Dihydrogen Orthophosphate extrapure AR, 99%**	500 Gms 1 Kg 5 Kg	770 1440 6360
78297 [7778-77-0]	**Potassium Dihydrogen Orthophosphate ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 1 Kg 5 Kg	840 1620 6840
54358 [7778-77-0]	**Potassium Dihydrogen Orthophosphate for molecular biology, 99.5%** (Potassium Phosphate Monobasic)	250 Gms 500 Gms 5 Kg	720 1380 13200
30955 [7778-77-0]	**Potassium Dihydrogen Orthophosphate for HPLC & UV Spectroscopy, 99.5%**	500 Gms 5 Kg	1180 11400
32360 [7778-77-0]	**Potassium Dihydrogen Orthophosphate for tissue culture, 99.5%** (Potassium Phosphate Monobasic)	500 Gms 1 Kg 5 Kg	900 1740 7200
77012 [140-89-6]	**Potassium Ethyl Xanthate (PEX) extrapure, 98%** (O-Ethyl Xanthic Acid Potassium Salt)	100 Gms 500 Gms	1400 5800
15766 [13746-66-2]	**Potassium Ferricyanide extrapure, 98%**	100 Gms 500 Gms 5 Kg	450 1950 18000
59558 [13746-66-2]	**Potassium Ferricyanide ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP	250 Gms 500 Gms 5 Kg	1150 2130 19840
32294 [14459-95-1]	**Potassium Ferrocyanide ACS, 98.5%**	500 Gms 5 Kg	1040 9780
12914 [14459-95-1]	**Potassium Ferrocyanide ExiPlus, Multi-Compendial, 98.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	1100 10350
97535 [7789-23-3]	**Potassium Fluoride Anhydrous extrapure, 97%**	500 Gms 5 Kg	800 7000
91489 [16921-30-5]	**Potassium Hexachloroplatinate (IV) pure, 40% Pt** (Platinum Potassium Chloride, Potassium Chloroplatinate)	250 Mg 1 Gms 5 Gms 25 Gms	5000 14000 64000 310000
42414 [17084-13-8]	**Potassium Hexafluorophosphate extrapure, 99%**	25 Gms 100 Gms 500 Gms	810 1960 6900
40180 [16919-27-0]	**Potassium Hexafluorotitanate pure, 98%** (Potassium Fluotitanate, Potassium Titanium Fluoride, Dipotassium hexafluorotitanate, Titanium potas	100 Gms 500 Gms 5 Kg	580 1380 9200
37466 [16923-95-8]	**Potassium Hexafluorozirconate pure, 97%**	100 Gms 500 Gms	1240 2880
39535 [16905-14-9]	**Potassium Hexaiodoplatinate (IV) pure, 18-19% Pt**	250 Mg 1 Gms	3910 11500
80633 [877-24-7]	**Potassium Hydrogen Phthalate ACS, 99.95%**	500 Gms	830
29467 [877-24-7]	**Potassium Hydrogen Phthalate ExiPlus, Multi-Compendial, 99.95%** meets compendial specs of BP, Ph.Eur	500 Gms	950
49017 [877-24-7]	**Potassium Hydrogen Phthalate for HPLC & UV Spectroscopy, 99.9%**	500 Gms	1130
65031 [7646-93-7]	**Potassium Hydrogen Sulphate (Potassium Bisulphate) pure, 99%**	500 Gms 5 Kg	690 6670
10942 [7646-93-7]	**Potassium Hydrogen Sulphate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%** (Potassium Bisulphate) meets compendial specs of BP, Ph.Eur	500 Gms	830
76131 [868-14-4]	**Potassium Hydrogen Tartrate (Potassium Bitartrate) pure, 99%** (Potassium L-Tartrate Monobasic)	500 Gms 5 Kg	980 9090
84749 [1310-58-3]	**Potassium Hydroxide pellets extrapure, 85%**	500 Gms 1 Kg 5 Kg 25 Kg	560 980 4720 18400

Code	Name	Packing	Price
78269 [1310-58-3]	**Potassium Hydroxide Pellets Reagent Grade, 85%**	500 Gms	600
		1 Kg	1150
		5 Kg	5060
18543 [1310-58-3]	**Potassium Hydroxide Pellets ExiPlus, Multi-Compendial, 85%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	650
		5 Kg	5520
26035 [1310-58-3]	**Potassium Hydroxide Powder extrapure, 85%, 80-100 mesh**	500 Gms	580
		1 Kg	1100
		5 Kg	4830
67864 [7758-05-6]	**Potassium Iodate pure, 99%**	100 Gms	2360
		250 Gms	5600
		500 Gms	10760
11041 [7758-05-6]	**Potassium Iodate ACS, 99.4%**	100 Gms	2470
		250 Gms	5830
		500 Gms	11090
10809 [7758-05-6]	**Potassium Iodate ExiPlus, Multi-Compendial, 99.4%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	2690
		250 Gms	6280
		500 Gms	11880
78240 [7681-11-0]	**Potassium Iodide pure, 99%**	25 Gms	900
		100 Gms	3480
		250 Gms	8070
		500 Gms	15680
27874 [7681-11-0]	**Potassium Iodide ACS, 99.8%**	25 Gms	960
		100 Gms	3640
		250 Gms	8400
		500 Gms	16520
97028 [7681-11-0]	**Potassium Iodide ExiPlus, Multi-Compendial, 99.8%** meets compendial specs of USP, BP	25 Gms	1040
		100 Gms	3810
		250 Gms	8960
		500 Gms	17360
11677 [7681-11-0]	**Potassium Iodide 0.1M Solution**	500 ml	450
		1000 ml	680
74567 [7681-11-0]	**Potassium Iodide 1M Solution**	500 ml	2360
		1000 ml	4600
26645 [7681-11-0]	**Potassium Iodide 3M Solution**	500 ml	7060
		1000 ml	13560
93842 [4429-42-9]	**Potassium Metabisulphite pure, 95%**	500 Gms	580
		5 Kg	5410
29078 △ [7758-09-0]	**Potassium Nitrite ACS, 96%**	100 Gms	810
		500 Gms	2420
13413 [7790-21-8]	**Potassium Periodate ACS, 99.8%**	100 Gms	3110
		500 Gms	14610
86742 [7790-21-8]	**Potassium Periodate ExiPlus, Multi-Compendial, 99.8%** meets compendial specs of BP, Ph.Eur	100 Gms	3450
		500 Gms	16100
77407 [6487-48-5]	**Potassium Oxalate Monohydrate pure, 99%**	500 Gms	740
		5 Kg	6900
40707 [6487-48-5]	**Potassium Oxalate Monohydrate extrapure AR, ACS, 99.5%**	500 Gms	860
		5 Kg	8050
36903 [7722-64-7]	**Potassium Permanganate extrapure, 99%**	500 Gms	640
		5 Kg	4600
51470 [7722-64-7]	**Potassium Permanganate ACS, 99%**	500 Gms	690
		5 Kg	5060
26500 [7722-64-7]	**Potassium Permanganate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	750
		5 Kg	5750
10322 [7727-21-1]	**Potassium Persulphate extrapure, 97%**	500 Gms	600
		5 Kg	5750
52254 [7727-21-1]	**Potassium Persulphate extrapure AR, 98%**	500 Gms	690
		5 Kg	6330
48319 [7727-21-1]	**Potassium Persulphate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	760
		5 Kg	7130
90654 [7758-11-4]	**Potassium Phosphate Dibasic Anhydrous pure, 98%**	500 Gms	690
		5 Kg	5750
		25 Kg	27600
27387 [7758-11-4]	**Potassium Phosphate Dibasic Anhydrous ACS, 98%**	500 Gms	750
		5 Kg	6330
31964 [7758-11-4]	**Potassium Phosphate Dibasic Anhydrous ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	810
		5 Kg	6900
41295 [7758-11-4]	**Potassium Phosphate Dibasic Anhydrous for HPLC, 99.5%**	500 Gms	1380
57419 [7758-11-4]	**Potassium Phosphate Dibasic Anhydrous for tissue culture, 99%**	500 Gms	780
		5 Kg	6560
37859 [16788-57-1]	**Potassium Phosphate Dibasic Trihydrate extrapure, 99%**	500 Gms	1040
		5 Kg	9780
15914 [7778-53-2]	**Potassium Phosphate Tribasic Anhydrous extrapure AR, ACS, 98%** (Tripotassium Phosphate)	500 Gms	1380
37423 [12208-13-8]	**Potassium Pyroantimonate extrapure, 98%** (Potassium Hexahydroxoantimonate(V), Potassium Antimonate Trihydrate)	100 Gms	1840
		500 Gms	8050
94527 [12208-13-8]	**Potassium Pyroantimonate extrapure AR, 99%** (Potassium Hexahydroxoantimonate(V), Potassium Antimonate Trihydrate)	100 Gms	2070
		500 Gms	8510
50862 [7320-34-5]	**Potassium Pyrophosphate extrapure, 98%** Tetra Potassium Pyrophosphate, Potassium diphosphate tetrabasic, Tetrapotassium diphosphate, Potassi	500 Gms	1100

Code	Name	Packing	Price
53158 [7790-62-7]	**Potassium Pyrosulphate (Potassium Disulfate) extrapure, 97.5%**	100 Gms 500 Gms	2070 7480
86344 [7492-30-0]	**Potassium Ricinoleate extrapure** (Ricinoleic Acid Potassium Salt, Potassium Ricinolate)	500 Gms 5 Kg	1100 10350
44277 [7778-80-5]	**Potassium Sulphate extrapure, 98.5%**	500 Gms 5 Kg	490 3910
52182 [7778-80-5]	**Potassium Sulphate for tissue culture, 99%**	500 Gms	640
26253 [7778-80-5]	**Potassium Sulphate ACS 99%**	500 Gms 5 Kg	530 4830
73900 [7778-80-5]	**Potassium Sulphate ExiPlus, Multi-Compendial, 98.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	580 5180
66691 [7790-58-1]	**Potassium Tellurite Anhydrous extrapure, 98%** (Dipotassium Tellurite, Potassium Tellurate, Potassium Tellurium Oxide)	10 Gms 25 Gms 100 Gms	2070 4370 17830
35889 [123333-66-4]	**Potassium Tellurite Monohydrate extrapure, 90%** (Potassium Tellurate (IV) Hydrate, Dipotassium Tellurite Monohydrate)	10 Gms 25 Gms 100 Gms	1610 3860 14950
16273 [10025-99-7]	**Potassium Tetrachloroplatinate (II) pure, 47% Pt** (Potassium Platinum (II) Chloride)	250 Mg 1 Gms 5 Gms	6330 20700 97750
34570 [14484-69-6]	**Potassium Tetrafluoroaluminate pure, 98%** (Potassium Aluminium Fluoride, PAF)	500 Gms 5 Kg	1040 8280
99928 [14075-53-7]	**Potassium Tetrafluoroborate pure, 98%** (Potassium Fluoroborate, Potassium Fluoborate, Potassium Borofluoride)	250 Gms 1 Kg	290 1040
48218 [333-20-0]	**Potassium Thiocyanate pure, 97%**	500 Gms	900
41027 [333-20-0]	**Potassium Thiocyanate ACS, 98.5%**	500 Gms 5 Kg	1060 9500
59115 [333-20-0]	**Potassium Thiocyanate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	1330 11000
69285 [9005-25-8]	**Potato Starch pure**	500 Gms 5 Kg	360 3150
95365 [344-25-2]	**D-Proline extrapure CHR, 99%**	5 Gms 25 Gms	1500 4840
35139 [147-85-3]	**L-Proline extrapure CHR, 99%**	25 Gms 100 Gms 1 Kg 5 Kg	300 750 6500 30000
56073 [147-85-3]	**L-Proline for tissue culture, 99%**	25 Gms 100 Gms 500 Gms	350 900 3600
96877 [147-85-3]	**L-Proline ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	100 Gms 1 Kg 5 Kg	1260 8000 36420
61988 [147-85-3]	**Proline BP** (L-Proline BP)	1 Kg 5 Kg 25 Kg	8000 37500 175000
68128 [147-85-3]	**L-Proline USP**	1 Kg 5 Kg 25 Kg	9200 42500 200000
39236 [71-23-8]	**n-Propanol pure, 99%** (1-Propanol)	500 ml 2500 ml	300 1200
24501 [71-23-8]	**n-Propanol ACS, 99.5%** (1-Propanol)	500 ml 2500 ml	330 1300
53790 [71-23-8]	**n-Propanol ExiPlus, Multi-Compendial, 99.5%** (1-Propanol) meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	360 1500
42328 [71-23-8]	**n-Propanol GC-HS, 99.9%** (1-Propanol)	250 ml 1000 ml	420 1400
45909 [71-23-8]	**n-Propanol for HPLC, 99.8%** (1-Propanol)	1000 ml	1600
85994 [71-23-8]	**n-Propanol Dried, 99.5%, 0.005% water** (1-Propanol)	1000 ml	1200
52111 [123-38-6]	**Propionaldehyde (Propanal) pure, 97%**	500 ml 2500 ml	390 2000
12931 [79-09-4]	**Propionic Acid extrapure, C3-99% (GC)**	500 ml 2500 ml	340 1600
43883 [79-09-4]	**Propionic Acid ACS, C3-99.5% (GC)**	500 ml 2500 ml	380 1800
42263 [79-09-4]	**Propionic Acid ExiPlus, Multi-Compendial, C3-99.5% (GC)** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	410 2000
52810 [123-62-6]	**Propionic Anhydride extrapure, 98%**	500 ml 2500 ml	1040 4890
21587 [79-03-8]	**Propionyl Chloride extrapure, 99%**	100 ml 500 ml	330 950
85159 [79-03-8]	**Propionyl Chloride (Specially purified) Reagent Grade, 99.5%**	500 ml	1380
28849 [107-10-8]	**n-Propylamine pure, 99%**	500 ml	640
30754 [57-55-6]	**Propylene Glycol (1,2-Propanediol) pure, 99%**	500 ml 2500 ml 25 Ltr	480 2220 20400
53054 [57-55-6]	**Propylene Glycol (1,2-Propanediol) ACS, 99.5%**	500 ml 2500 ml	560 2520
95961 [57-55-6]	**Propylene Glycol (1,2-Propanediol), ExiPlus, Multi-Compendial , 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	600 2880

Code	Name	Packing	Price
40011 [121-79-9]	**Propyl Gallate extrapure, 98%**	100 Gms 500 Gms	1090 4260
47638 [121-79-9]	**Propyl Gallate IP**	1 Kg 5 Kg 25 Kg	8000 37500 175000
14086 [121-79-9]	**Propyl Gallate BP**	1 Kg 5 Kg 25 Kg	8000 37500 175000
16067 [121-79-9]	**Propyl Gallate NF**	1 Kg 5 Kg 25 Kg	9200 42500 200000
53898 [94-13-3]	**Propyl-4-Hydroxybenzoate (Propyl Paraben) ExiPlus, Multi-Compendial, 98-102%** meets compendial specs of USP, BP	500 Gms 5 Kg	1340 11890
52950 [35285-69-9]	**Propyl-4-Hydroxybenzoate Sodium Salt (Propyl Paraben Sodium) ExiPlus, Multi-Compendial, 94-102%** meets compendial specs of USP, BP	500 Gms 5 Kg	1400 14380
61808 [110-86-1]	**Pyridine extrapure, 99%**	250 ml 500 ml 2500 ml	520 980 4600
41249 [110-86-1]	**Pyridine ACS, 99%**	500 ml 2500 ml	1040 4780
51952 [110-86-1]	**Pyridine Dried, 99.5%, water 0.01%**	500 ml	1500
83329 [110-86-1]	**Pyridine for HPLC & UV Spectroscopy, 99.5%**	500 ml 1000 ml	1380 2510
75453 [110-86-1]	**Pyridine (For Karl Fischer Titration), 99.5%**	500 ml	1210
92665 [110-86-1]	**Pyridine ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	250 ml 500 ml 2500 ml	600 1120 4950
48918 [26299-14-9]	**Pyridinium Chlorochromate (PCC) extrapure, 98%**	25 Gms 100 Gms 500 Gms	580 1560 6890
13584 [20039-37-6]	**Pyridinium Dichromate (Cornforth Reagent, PDC) extrapure, 98%**	25 Gms 100 Gms 500 Gms	860 2330 10190
10012 [24057-28-1]	**Pyridium p-Toluenesulfonate (PPTS) extrapure, 98%** (p-Toluenesulfonic Acid Pyridine salt, PPTS)	100 Gms 500 Gms	2070 4600
71646 [87-66-1]	**Pyrogallol ACS, 99%**	100 Gms 500 Gms	1730 7480
30205 [87-66-1]	**Pyrogallol ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP	100 Gms 500 Gms	1900 7820
87445 △	**R.B.C. Diluting Fluid (Hayem`s Solution)** (RBC Diluting Reagent, Hayems Reagent)	500 ml	460
58050 [81-88-9]	**Rhodamine B** (Tetraethylrhodamine, Violet 10)	25 Gms 100 Gms 500 Gms	280 1010 4370
20707 [7791-11-9]	**Rubidium (I) Chloride ultrapure, 99.5%**	5 Gms 25 Gms	4140 19550
75420 [7440-18-8]	**Ruthenium Metal Powder, 99.9%**	1 Gms	14950
60633 [11103-72-3]	**Ruthenium Red Tetrahydrate**	1 Gms 5 Gms	7500 35000
75268 [14898-67-0]	**Ruthenium (III) Trichloride Hydrate pure, 38-40% Ru**	1 Gms 5 Gms	6000 28000
77862 [90-02-8]	**Salicylaldehyde pure, 99%** (2-Hydroxybenzaldehyde, o-Hydroxybenzaldehyde)	100 ml 250 ml 1000 ml	540 1100 4200
24922 [90-02-8]	**Salicylaldehyde Reagent Grade, 99%** (2-Hydroxybenzaldehyde, o-Hydroxybenzaldehyde)	100 ml 500 ml	700 2500
63504 [90-02-8]	**Salicylaldehyde ExiPlus, Multi-Compendial, 99%** (2-Hydroxybenzaldehyde, o-Hydroxybenzaldehyde) meets compendial specs of USP, BP, Ph.Eur	100 ml 250 ml 500 ml	750 1600 2800
11453 [69-72-7]	**Salicylic Acid pure, 99%**	500 Gms 5 Kg 25 Kg	840 8000 38000
62657 [69-72-7]	**Salicylic Acid ACS, 99%**	500 Gms 2.5 Kg	920 4400
16448 [69-72-7]	**Salicylic Acid ACS, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	980 9200
27769	**Selenium AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml 250 ml	1800 2900
52228	**Seralite SRA-400, Strong Basic Anion Exchange Resin** (equivalent to Amberlite IRA-400, Indion FFIP Anion Resin)	500 Gms 5 Kg	850 7800
14891	**Seralite SRC-120, Strong Acid Cation Exchange Resin** (equivalent to Amberlite IRC-120, Indion 225H Cation Resin)	500 Gms 5 Kg	550 5000
89830	**Seralite Resin WRC-50, Weak Acid Cation Exchange Resin** (equivalent to Amberlite IRC-50, Indion 236 Resin)	500 Gms 5 Kg	700 6500
15776	**SeraPreserve 10X Tissue Fixative Concentrate (Formalin Free)**	5 Ltr	5800
40148	**SeraPreserve 25X Tissue Fixative Concentrate (Formalin Free)**	5 Ltr	8760

Code	Name	Packing	Price
21569	**SeraPreserve 50X Tissue Fixative Concentrate (Formalin Free)**	5 Ltr	12500
53457 [312-84-5]	**D-Serine extrapure, 99%**	25 Gms	800
		100 Gms	3000
		500 Gms	10000
24425 [302-84-1]	**DL-Serine extrapure CHR, 99%**	25 Gms	600
		100 Gms	2000
		500 Gms	9000
84131 [56-45-1]	**L-Serine extrapure CHR, 99%**	25 Gms	500
		100 Gms	1500
		500 Gms	6500
63503 [56-45-1]	**L-Serine for tissue culture, 99%**	25 Gms	600
		100 Gms	1600
		500 Gms	7200
84494 [56-45-1]	**L-Serine ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP	25 Gms	820
		100 Gms	2400
		1 Kg	16000
85148 [112926-00-8]	**Silica Gel Blue (Self Indicating) (Coarse), 5-8 mesh**	500 Gms	560
		5 Kg	5240
84210 [112945-52-5]	**Silica Gel Fumed 200 (Hydrophilic)** (Silicon Dioxide amorphous , Aerosil 200)	100 Gms	400
		250 Gms	750
		2.5 Kg	6930
71548 [112926-00-8]	**Silica Gel Orange (Self Indicating), 2-5mm**	500 Gms	630
		5 Kg	6090
65856 [112926-00-8]	**Silica Gel 60-120 mesh**	500 Gms	500
		5 Kg	4600
		25 Kg	21610
40449 [112926-00-8]	**Silica Gel 60-200 mesh**	500 Gms	590
		5 Kg	5270
95178 [112926-00-8]	**Silica Gel 100-200 mesh**	500 Gms	500
		5 Kg	4600
		25 Kg	21610
63025 [112926-00-8]	**Silica Gel 200-400 mesh**	500 Gms	590
		5 Kg	5500
		25 Kg	26290
96671 [112926-00-8]	**Silica Gel 230-400 mesh**	500 Gms	660
		5 Kg	6200
		25 Kg	29210
36834 [112926-00-8]	**Silica Gel 400-700 mesh**	500 Gms	830
		1 Kg	1540
51849 [112926-00-8]	**Silica Gel G for TLC, 325mesh, w/ CaSO4**	500 Gms	570
		5 Kg	5150
		25 Kg	24540
52797 [112926-00-8]	**Silica Gel H for TLC, 325mesh, w/o CaSO4**	500 Gms	660
		5 Kg	5930
		25 Kg	28230
38062 [112926-00-8]	**Silica Gel GF254 for TLC, 325mesh, w/ Fluorescent Indicator w/ CaSO4**	500 Gms	4720
		5 Kg	45890
29774 [112926-00-8]	**Silica Gel HF254 for TLC, 325mesh, w/ Fluorescent Indicator w/o CaSO4**	500 Gms	4720
		5 Kg	45890
81618	**Silicon Antifoaming Agent 30% (Non-ionic)**	100 ml	430
		500 ml	1150
		1000 ml	2060
72243 [63148-62-9]	**Silicone Oil 100 CPS**	250 ml	1040
		1000 ml	3630
		2500 ml	8220
82045 [63148-62-9]	**Silicon Oil 330-380 CPS**	250 ml	910
		1000 ml	3140
		2500 ml	7490
85160 [563-63-3]	**Silver Acetate pure, 98%**	10 Gms	10500
		25 Gms	24000
		100 Gms	90000
32737	**Silver AAS Standard Solution in 0.5N HNO3 (Traceable to NIST)**	100 ml	3000
		250 ml	6500
77701 [7785-23-1]	**Silver Bromide extrapure, 98%**	10 Gms	9400
		25 Gms	22000
		100 Gms	85000
90097 [534-16-7]	**Silver Carbonate pure, 99%**	10 Gms	13000
		25 Gms	30000
		100 Gms	115000
80647 [7783-90-6]	**Silver Chloride extrapure AR, 99%**	10 Gms	13000
		25 Gms	30000
		100 Gms	115000
14417 [1470-61-7]	**Silver Diethyldithiocarbamate ACS**	5 Gms	4000
		25 Gms	19000
		100 Gms	70000
43800 [26042-64-8]	**Silver Hexafluoroantimonate(V) extrapure, 98%**	1 Gms	9500
		5 Gms	42000
86633 [7783-97-3]	**Silver Iodate extrapure, 99%**	10 Gms	6700
		25 Gms	16000
		100 Gms	60000
24607 [7783-96-2]	**Silver Iodide extrapure AR, 99%**	10 Gms	7800
		25 Gms	18500
		100 Gms	70000
32422 [7440-22-4]	**Silver Metal powder, 99.9%, -400 mesh**	10 Gms	16000
		25 Gms	37000
		100 Gms	140000
43491 [2386-52-9]	**Silver Methanesulfonate (Silver MSA) extrapure, 47.8%-58.4% Ag** (Methanesulfonic Acid Silver Salt, Silver Methylsulfonate, Silver Methylsulphonate, Methanesulphonic	1 Gms	3500
		5 Gms	14000
		25 Gms	60000
32814 [7761-88-8]	**Silver Nitrate extrapure, 99.5%**	10 Gms	9000
		25 Gms	21600
		100 Gms	83700
		500 Gms	414000
91233 [7761-88-8]	**Silver Nitrate for tissue culture, 99.5%**	10 Gms	9900
		25 Gms	23400

Code	Name	Packing	Price
94118 [7761-88-8]	**Silver Nitrate ACS, 99.9%**	10 Gms	9900
		25 Gms	22500
		100 Gms	85500
		500 Gms	423000
96045 [7761-88-8]	**Silver Nitrate ExiPlus, Multi-Compendial, 99.9%** meets compendial specs of USP, BP	25 Gms	23400
		100 Gms	87300
		500 Gms	432000
51535 [7761-88-8]	**Silver Nitrate 0.1N Solution**	500 ml	4700
		1000 ml	8600
		2500 ml	20700
54023 [7761-88-8]	**Silver Nitrate 0.02N Solution**	500 ml	2200
		1000 ml	4000
47262 [7761-88-8]	**Silver Nitrate 0.05N Solution**	500 ml	9900
		1000 ml	18000
49430 [7761-88-8]	**Silver Nitrate 1N Solution**	100 ml	35100
		500 ml	171000
69668 [20667-12-3]	**Silver Oxide pure, 99%**	10 Gms	15000
		25 Gms	35000
		100 Gms	135000
76610 [20667-12-3]	**Silver Oxide extrapure AR, 99%**	10 Gms	16000
		25 Gms	37000
63976 [7784-09-0]	**Silver Phosphate extrapure, 98%** (Silver Orthophosphate, Trisilver Phosphate)	5 Gms	4500
		25 Gms	21000
86887 [7784-09-0]	**Silver Phosphate extrapure AR, 99%** (Silver Orthophosphate, Trisilver Phosphate)	5 Gms	4800
		25 Gms	22000
46241 [9008-42-8]	**Silver Proteinate extrapure, 8% Ag**	5 Gms	2000
		25 Gms	7000
60112 [22199-08-2]	**Silver Sulfadiazine ExiPlus, Multi-Compendial, 98%** (Silver (I) Sulfadiazine, Sulfadiazine Silver) meets compendial specs of USP	5 Gms	3000
		25 Gms	12000
		100 Gms	46000
17929 [10294-26-5]	**Silver Sulphate pure, 98%**	10 Gms	9900
		25 Gms	23400
		100 Gms	90000
		500 Gms	432000
59940 [10294-26-5]	**Silver Sulphate ACS, 98%**	10 Gms	10800
		25 Gms	25200
		100 Gms	94500
		500 Gms	450000
27821 [2923-28-6]	**Silver Trifluoromethanesulfonate extrapure, 98%** (Silver Triflate, Trifluoromethanesulfonic Acid Silver Salt)	1 Gms	3000
		5 Gms	7000
		25 Gms	30000
69717 [13497-94-4]	**Silver Vanadate pure**	5 Gms	6000
38119 [8006-28-8]	**Soda Lime Granules with Indicator extrapure AR, 2-5mm, CO2 abs. capacity 27%**	500 Gms	690
		5 Kg	6560

Code	Name	Packing	Price
31378 [127-09-3]	**Sodium Acetate Anhydrous pure, 98%**	500 Gms	300
		5 Kg	2730
		25 Kg	13000
74537 [127-09-3]	**Sodium Acetate Anhydrous ACS, 99%**	500 Gms	440
		5 Kg	3470
		25 Kg	16500
51483 [127-09-3]	**Sodium Acetate Anhydrous ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	530
		5 Kg	4200
		25 Kg	18500
22342 [127-09-3]	**Sodium Acetate Anhydrous for molecular biology, 99%**	250 Gms	550
		500 Gms	950
70509 [127-09-3]	**Sodium Acetate Anhydrous for HPLC & UV Spectroscopy, 99%**	500 Gms	1380
54966 [6131-90-4]	**Sodium Acetate Trihydrate ACS, 99%**	500 Gms	310
		5 Kg	2730
		25 Kg	13000
44981 [6131-90-4]	**Sodium Acetate Trihydrate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	340
		5 Kg	3150
		25 Kg	14500
88035 [6131-90-4]	**Sodium Acetate Trihydrate for molecular biology, 99.5%**	250 Gms	360
		500 Gms	690
84658 [6131-90-4]	**Sodium Acetate Trihydrate for tissue culture, 99.5%**	500 Gms	280
56315 [6131-90-4]	**Sodium Acetate Trihydrate for HPLC, 99.5%**	500 Gms	1580
21636 [7758-16-9]	**Sodium Acid Pyrophosphate (SAPP) extrapure, 95%** (Sodium Pyrophosphate Dibasic, Disodium Diphosphate Dibasic, Disodium Dihydrogen Diphosphate)	500 Gms	490
80726 [9005-38-3]	**Sodium Alginate LV extrapure, 4-12cps, 30mesh**	500 Gms	3000
		5 Kg	28000
52775 [9005-38-3]	**Sodium Alginate HV extrapure, 1000cps, 100mesh** (Alginic Acid Sodium Salt)	500 Gms	3200
		5 Kg	31000
75063	**Sodium AAS Standard Solution in 0.5N HCl (Traceable to NIST)**	100 ml	2070
		250 ml	3340
45349 [15432-85-6]	**Sodium Antimonate pure, 99%**	100 Gms	1610
		250 Gms	3450
84070 [26628-22-8]	**Sodium Azide ACS, 99%**	100 Gms	640
		500 Gms	2530
		5 Kg	24150
		25 Kg	115000
17782 [26628-22-8]	**Sodium Azide for molecular biology, 99%**	100 Gms	830
		500 Gms	3570
11077 [26628-22-8]	**Sodium Azide PrimeDx, 99%**	5 Kg	34500
		25 Kg	161000

Code	Name	Packing	Price
65171 [532-32-1]	**Sodium Benzoate pure, 99%**	500 Gms	550
		5 Kg	4830
11114 [532-32-1]	**Sodium Benzoate Exiplus, Multicompendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	710
		5 Kg	5980
56398 [144-55-8]	**Sodium Bicarbonate extrapure, 99%**	500 Gms	280
		5 Kg	2190
		25 Kg	10350
45437 [144-55-8]	**Sodium Bicarbonate ACS, 99.7%**	500 Gms	330
		5 Kg	2530
		25 Kg	10930
89399 [144-55-8]	**Sodium Bicarbonate ACS, ExiPlus, Multi-Compendial, 99.7%**	500 Gms	360
		5 Kg	2760
36328 [144-55-8]	**Sodium Bicarbonate for molecular biology, 99.7%**	500 Gms	920
92957 [144-55-8]	**Sodium Bicarbonate for tissue culture, 99.7%**	500 Gms	810
		5 Kg	6900
93982 [16940-66-2]	**Sodium Borohydride (powder) ACS, 98%** (Sodium Tetrahydridoborate)	25 Gms	450
		100 Gms	950
		500 Gms	4200
		5 Kg	39000
43471 [16940-66-2]	**Sodium Borohydride (crystals) ACS, 98%** (Sodium Tetrahydridoborate)	25 Gms	500
		100 Gms	1050
		500 Gms	4600
		5 Kg	42000
32769 [7789-38-0]	**Sodium Bromate extrapure AR, 99%**	100 Gms	350
		500 Gms	1150
43327 [7647-15-6]	**Sodium Bromide extrapure, 99%**	500 Gms	920
		5 Kg	8510
84486 [6131-99-3]	**Sodium Cacodylate Trihydrate extrapure, 98%** (Dimethyl Arsenic Acid Sodium Salt Trihydrate)	10 Gms	4140
		25 Gms	10010
		100 Gms	38640
42215 [1984-06-1]	**Sodium Caprylate extrapure, 99%** (Caprylic Acid Sodium Salt)	100 Gms	700
		500 Gms	2000
50090 [1984-06-1]	**Sodium Caprylate BP, 99-101%**	1 Kg	7800
		5 Kg	36900
81340 [1984-06-1]	**Sodium Caprylate NF, 99-101%**	1 Kg	8000
		5 Kg	36400
89382 [497-19-8]	**Sodium Carbonate Anhydrous pure, 99.5%**	500 Gms	280
		5 Kg	2420
		25 Kg	11500
51938 [497-19-8]	**Sodium Carbonate Anhydrous ACS, 99%, 80-100mesh**	500 Gms	380
		5 Kg	3450
64079 [497-19-8]	**Sodium Carbonate Anhydrous ACS, 99.9%**	500 Gms	360
		5 Kg	3080
		25 Kg	14660
93857 [497-19-8]	**Sodium Carbonate Anhydrous ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	410
		5 Kg	3640
41721 [7647-14-5]	**Sodium Chloride ACS, 99.9%**	500 Gms	210
		1 Kg	400
		5 Kg	1300
		25 Kg	6010
76945 [7647-14-5]	**Sodium Chloride ExiPlus, Multi-Compendial, 99.9%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	220
		1 Kg	430
		5 Kg	1350
		25 Kg	6200
33205 [7647-14-5]	**Sodium Chloride for molecular biology, 99.9%**	500 Gms	660
		1 Kg	1200
		5 Kg	5560
26416 [7647-14-5]	**Sodium Chloride for tissue culture, 99.9%**	500 Gms	700
		1 Kg	1200
76154 [7647-14-5]	**Sodium Chloride for HPLC, 99.9%**	500 Gms	530
21522 [6132-05-4]	**Sodium Citrate Dibasic Sesquihydrate extrapure, 99%** (Citric Acid Disodium Salt Sesquihydrate, Disodium Hydrogen Citrate Sesquihydrate)	100 Gms	600
		500 Gms	1200
83407 [6132-04-3]	**Sodium Citrate Tribasic Dihydrate extrapure, 98%** (Trisodium Citrate Dihydrate)	500 Gms	300
		1 Kg	550
		5 Kg	2400
		25 Kg	9500
85919 [6132-04-3]	**Sodium Citrate Tribasic Dihydrate ACS, 99%**	500 Gms	320
		5 Kg	2600
		25 Kg	10200
94968 [6132-04-3]	**Sodium Citrate Tribasic Dihydrate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	450
		5 Kg	3500
67331 [6132-04-3]	**Sodium Citrate Tribasic Dihydrate for molecular biology, 99%**	250 Gms	550
		1 Kg	950
43819 [13600-98-1]	**Sodium Cobaltinitrite ACS** (Sodium Hexanitrocobaltate III)	25 Gms	1150
		100 Gms	2300
		500 Gms	7940
12680 [25895-60-7]	**Sodium Cyanoborohydride pure, 95%** (Sodium Cyanotrihydridoborate, Sodium Borocyanohydride)	5 Gms	1040
		25 Gms	3110
		100 Gms	9660
		500 Gms	41060
11446 [126-96-5]	**Sodium Diacetate pure** (Sodium Hydrogen Diacetate)	500 Gms	400
94963 [7789-12-0]	**Sodium Dichromate Dihydrate pure, 98%**	500 Gms	920
		5 Kg	7710
54750 [7789-12-0]	**Sodium Dichromate Dihydrate ACS, 99.5%**	500 Gms	1040
		5 Kg	8970

Code	Name	Packing	Price
11996 [7789-12-0]	**Sodium Dichromate Dihydrate ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP	500 Gms 5 Kg	1150 10350
11226 [20624-25-3]	**Sodium Diethyldithiocarbamate Trihydrate ACS,99 %**	100 Gms 500 Gms	710 2910
22607 [20624-25-3]	**Sodium Diethyldithiocarbamate Trihydrate ACS, ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP	100 Gms 500 Gms	810 3570
49653 [7775-14-6]	**Sodium Dithionite extrapure, 85%** (Sodium Hypodisulphite, Sodium Hydrosulfite)	500 Gms 5 Kg	620 6000
68792 [141-52-6]	**Sodium Ethoxide extrapure, 95%** (Sodium Ethylate)	500 Gms 1 Kg	1200 2000
10961 [14434-22-1]	**Sodium Ferrocyanide Decahydrate extrapure AR, 99%** (Sodium Hexacyanoferrate(II), Tetrasodium Hexacyanoferrate Decahydrate)	500 Gms 5 Kg	580 4830
57010 [7681-49-4]	**Sodium Fluoride pure, 98%**	500 Gms 5 Kg	1150 10240
29821 [7681-49-4]	**Sodium Fluoride ACS, 99%**	500 Gms 5 Kg	1150 10930
40544 [7681-49-4]	**Sodium Fluoride ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	1270 11620
86307 [141-53-7]	**Sodium Formate pure, 98%**	500 Gms 5 Kg	280 2420
39728 [141-53-7]	**Sodium Formate ACS, 99%**	500 Gms	520
84439 [149-44-0]	**Sodium Formaldehyde Sulfoxylate Hydrate (Rongalite) pure, 95%** (Monosodium Hydroxymethanesulfinate)	500 Gms	690
11407 [149-44-0]	**Sodium Formaldehyde Sulfoxylate Hydrate (Rongalite) ExiPlus, Multi-Compendial, 95%** (Monosodium Hydroxymethanesulfinate) meets compendial specs of USP	500 Gms	920
54762 [16893-85-9]	**Sodium Hexafluorosilicate pure, 99%** (Sodium Fluorosilicate, Sodium Silicofluoride)	1 Kg 5 Kg	690 2990
32537 [68915-31-1]	**Sodium Hexametaphosphate extrapure AR, 65-70%** (SHMP, Calgon, Sodium Polyphosphate)	500 Gms 5 Kg	810 7250
98060 [10034-88-5]	**Sodium Hydrogen Sulphate (Sodium Bisulphate Monohydrate) extrapure AR, 99%**	500 Gms	580
73082 [1310-73-2]	**Sodium Hydroxide Flakes extrapure, 97%**	500 Gms 1 Kg 5 Kg 25 Kg	300 540 1930 9080
96311 [1310-73-2]	**Sodium Hydroxide Pellets purified, 97%**	500 Gms 1 Kg 5 Kg 25 Kg	320 550 2010 9540
68151 [1310-73-2]	**Sodium Hydroxide Pellets Reagent Grade, 97%**	500 Gms 1 Kg 5 Kg 25 Kg	380 680 2690 12600
13913 [1310-73-2]	**Sodium Hydroxide Pellets, ExiPlus, Multi-Compendial, 97%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 1 Kg 5 Kg	370 660 2880
39700 [1310-73-2]	**Sodium Hydroxide Powder extrapure, 96%, 80-100 mesh**	500 Gms 1 Kg 5 Kg	350 620 2270
13994 [1310-73-2]	**Sodium Hydroxide 0.1N Solution**	500 ml 1000 ml	160 280
69978 [1310-73-2]	**Sodium Hydroxide 1N Solution**	500 ml 1000 ml 2500 ml	170 310 750
25366 [7681-52-9]	**Sodium Hypochlorite Solution, 4-7% Cl** (Soda Bleach, Liquid Bleach)	500 ml 1000 ml 5 Ltr	180 320 740
98760 [7681-55-2]	**Sodium Iodate extrapure, 99%**	25 Gms 100 Gms 500 Gms	950 3220 13800
71964 [7681-55-2]	**Sodium Iodate extrapure AR, 99.5%**	25 Gms 100 Gms 500 Gms	1100 3570 14950
77758 [7681-82-5]	**Sodium Iodide pure, 99%**	25 Gms 100 Gms 250 Gms 1 Kg	830 3170 7190 27600
21704	**Sodium Iodide 1M Solution**	100 ml 500 ml	870 2880
11004 [9004-82-4]	**Sodium Lauryl Ether Sulphate (SLES, AES), 28-30%** (Sodium Laureth Sulfate)	1000 ml 5 Ltr	950 2540
46821 [151-21-3]	**Sodium Lauryl Sulphate (SLS, SDS) (Needle-shaped) extrapure, 90%** (Sodium Dodecyl Sulphate)	500 Gms 5 Kg	860 8190
14374 [151-21-3]	**Sodium Lauryl Sulphate (SDS, SLS) ExiPlus, Multi-Compendial, 85%** (Sodium Dodecyl Sulphate) meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	810 7340
54468 [151-21-3]	**Sodium Lauryl Sulphate (SDS, SLS) ACS, 99%** (Sodium Dodecyl Sulphate)	100 Gms 250 Gms 500 Gms 1 Kg	1340 2920 5340 10310

Code	Name	Packing	Price
99671 [151-21-3]	**Sodium Lauryl Sulphate (SDS, SLS) for HPLC, 99%** (Sodium Dodecyl Sulphate)	25 Gms	1610
		100 Gms	4410
		250 Gms	9100
		1 Kg	22440
32096 [151-21-3]	**Sodium Lauryl Sulphate (SDS, SLS) for molecular biology, 99%** (Sodium Dodecyl Sulphate)	25 Gms	590
		100 Gms	1820
		500 Gms	6010
		1 Kg	11400
35825 [151-21-3]	**Sodium Lauryl Sulphate (SDS, SLS) High Purity, 99.5%** (Sodium Dodecyl Sulphate)	25 Gms	810
		100 Gms	2100
85369	**Sodium Lauryl Sulphate (SDS, SLS) 10% Solution** (10% SDS Solution)	100 ml	400
87547	**Sodium Lauryl Sulphate (SDS, SLS) 20% Solution** (20% SDS Solution)	250 ml	550
34378 [151-21-3]	**Sodium Lauryl Sulphate (SDS, SLS) BP, 85%** (Sodium Dodecyl Sulphate)	1 Kg	1800
		5 Kg	8200
		25 Kg	36000
21705 [7681-57-4]	**Sodium Metabisulphite extrapure AR, 98%**	500 Gms	350
		5 Kg	3280
		25 Kg	15530
52693 [7681-57-4]	**Sodium Metabisulphite extrapure AR, ACS, ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP	500 Gms	380
		5 Kg	3630
51153 [7790-28-5]	**Sodium Metaperiodate (Sodium Periodate) pure 99%**	100 Gms	2000
		500 Gms	9600
23621 [7790-28-5]	**Sodium Metaperiodate (Sodium Periodate) ACS, 99.8%**	100 Gms	2300
		500 Gms	10000
30037 [7790-28-5]	**Sodium Metaperiodate (Sodium Periodate) ExiPlus, Multi-Compendial, 99.8%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	2400
		500 Gms	10200
94601 [10213-79-3]	**Sodium Metasilicate Pentahydrate (Sodium Silicate Hydrated) extrapure**	500 Gms	350
		1 Kg	650
		5 Kg	3100
20660 [13517-24-3]	**Sodium Metasilicate Nonahydrate (Sodium Silicate Hydrated) extrapure**	500 Gms	300
		1 Kg	500
		5 Kg	2200
41670 ✔ [124-41-4]	**Sodium Methoxide Powder pure, 98%**	500 Gms	550
		5 Kg	5000
13547 [10102-40-6]	**Sodium Molybdate Dihydrate pure, 98%**	100 Gms	1380
		250 Gms	2750
		500 Gms	5280
18108 [10102-40-6]	**Sodium Molybdate Dihydrate ACS, 99.5%**	100 Gms	1540
		250 Gms	3250
		500 Gms	5720
70851 [10102-40-6]	**Sodium Molybdate Dihydrate ExiPlus, Multi-Compendial, 98%** meets compendial specs of BP, Ph.Eur	100 Gms	1650
		250 Gms	3520
		500 Gms	6050
44618 [7631-99-4]	**Sodium Nitrate ACS, 99%**	500 Gms	450
		5 Kg	4000
71910 □ [7632-00-0]	**Sodium Nitrite ACS, 97%**	500 Gms	580
		5 Kg	5470
28182 [13755-38-9]	**Sodium Nitroprusside Dihydrate pure, 98%** (Sodium Nitroferricyanide)	100 Gms	1270
		500 Gms	5750
71750 [13755-38-9]	**Sodium Nitroprusside Dihydrate ACS, 99%**	100 Gms	1380
		500 Gms	5980
58403 [13755-38-9]	**Sodium Nitroprusside Dihydrate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	1500
		500 Gms	6330
62381 [10101-89-0]	**tri-Sodium Orthophosphate Dodecahydrate pure, 98%** (Sodium Phosphate Tribasic Dodecahydrate (Trisodium Phosphate Dodecahydrate))	500 Gms	430
		5 Kg	3870
73670 [10101-89-0]	**tri-Sodium Orthophosphate Dodecahydrate ACS, 99%** (Sodium Phosphate Tribasic Dodecahydrate, Trisodium Phosphate Dodecahydrate)	500 Gms	520
		5 Kg	5000
69199 [10101-89-0]	**tri-Sodium Orthophosphate Dodecahydrate ACS, ExiPlus, Multi-Compendial, 99%** (Sodium Phosphate Tribasic Dodecahydrate, Trisodium Phosphate Dodecahydrate) meets compendial specs of USP	500 Gms	630
		5 Kg	5720
36176 [62-76-0]	**Sodium Oxalate pure, 99%**	500 Gms	560
		5 Kg	5290
57554 [62-76-0]	**Sodium Oxalate ACS, 99.5%**	500 Gms	700
		5 Kg	6300
36984 [62-76-0]	**Sodium Oxalate ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of BP, Ph.Eur	500 Gms	800
39815 [15630-89-4]	**Sodium Percarbonate 13% Active Oxygen** (Sodium carbonate hydrogen peroxide)	500 Gms	410
		5 Kg	3680
93122 [7791-07-3]	**Sodium Perchlorate Monohydrate extrapure, 98.5%**	100 Gms	2850
		500 Gms	6900
53388 [1313-60-6]	**Sodium Peroxide ACS, 93%** (Sodium Superoxide)	100 Gms	5060
		250 Gms	12080
		500 Gms	25300
97295 [7775-27-1]	**Sodium Persulphate pure, 98%** (Sodium Peroxodisulfate)	500 Gms	570
53046 [7558-79-4]	**Sodium Phosphate Dibasic Anhydrous extrapure, 99%**	500 Gms	600
		5 Kg	5520

Code	Name	Packing	Price
66273 [7558-79-4]	**Sodium Phosphate Dibasic Anhydrous ACS, 99%**	500 Gms 1 Kg 5 Kg	690 1270 6100
45932 [7558-79-4]	**Sodium Phosphate Dibasic Anhydrous ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 1 Kg 5 Kg	790 1380 6560
21669 [7558-79-4]	**Sodium Phosphate Dibasic Anhydrous for molecular biology, 99.5%**	250 Gms 500 Gms 5 Kg	960 1790 16100
32080 [10028-24-7]	**Sodium Phosphate Dibasic Dihydrate extrapure, 99%**	500 Gms 5 Kg 25 Kg	580 4600 21280
87258 [10028-24-7]	**Sodium Phosphate Dibasic Dihydrate extrapure AR, 99.5%**	500 Gms 5 Kg 25 Kg	690 5180 24730
52234 [10028-24-7]	**Sodium Phosphate Dibasic Dihydrate ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	1040 8050
87027 [10028-24-7]	**Sodium Phosphate Dibasic Dihydrate for molecular biology, 99.5%**	250 Gms 500 Gms 5 Kg	690 1150 10120
86858 [10028-24-7]	**Sodium Phosphate Dibasic Dihydrate for HPLC, 99.5%**	500 Gms	1380
57085 [10039-32-4]	**Sodium Phosphate Dibasic Dodecahydrate extrapure, 98.5%** (Sodium Hydrogen Phosphate Dodecahydrate)	500 Gms 5 Kg	520 4030
83417 [10039-32-4]	**Sodium Phosphate Dibasic Dodecahydrate extrapure AR, 99%** (Sodium Hydrogen Phosphate Dodecahydrate)	500 Gms 5 Kg	580 5290
18134 [7782-85-6]	**Sodium Phosphate Dibasic Heptahydrate ACS, 98-102%** (Disodium Hydrogen Phosphate Heptahydrate, Disodium Phosphate)	100 Gms 500 Gms 2.5 Kg	920 2530 11730
22249 [7558-80-7]	**Sodium Phosphate Monobasic Anhydrous Reagent Grade, 99%**	500 Gms 5 Kg	690 5750
98336 [7558-80-7]	**Sodium Phosphate Monobasic Anhydrous for HPLC, 99%** (Sodium Dihydrogen Orthophosphate Anhydrous)	500 Gms	1300
59443 [7558-80-7]	**Sodium Phosphate Monobasic Anhydrous for molecular biology, 99%**	250 Gms 500 Gms	1130 1610
91826 [7558-80-7]	**Sodium Phosphate Monobasic Anhydrous for tissue culture, 99%**	250 Gms 1 Kg	460 1500
27092 [13472-35-0]	**Sodium Phosphate Monobasic Dihydrate extrapure, 99%**	500 Gms 5 Kg	510 4600
40597 [13472-35-0]	**Sodium Phosphate Monobasic Dihydrate extrapure AR, 99%**	500 Gms 5 Kg	570 5060
79088 [13472-35-0]	**Sodium Phosphate Monobasic Dihydrate extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP	500 Gms 5 Kg	640 5520
93494 [13472-35-0]	**Sodium Phosphate Monobasic Dihydrate for molecular biology, 99%** (Sodium Hydrogen Phosphate Dihydrate)	500 Gms 5 Kg	950 8280
61707 [10049-21-5]	**Sodium Phosphate Monobasic Monohydrate ACS, 98%** (Sodium Dihydrogen Orthophosphate Monohydrate, Monosodium phosphate, Sodium dihydrogen phosphate mon	500 Gms 5 Kg	690 6160
97768 [10049-21-5]	**Sodium Phosphate Monobasic Monohydrate for HPLC, 99%** (Sodium Dihydrogen Orthophosphate Monohydrate, Monosodium phosphate, Sodium dihydrogen phosphate mon	500 Gms	1130
78535 [7601-54-9]	**Sodium Phosphate Tribasic Anhydrous extrapure, 96%** (Sodium Orthophosphate, Trisodium Phosphate)	500 Gms 5 Kg	750 6790
63059 [6381-59-5]	**Sodium Potassium Tartrate Tetrahydrate pure, 98%**	500 Gms 5 Kg	720 6900
18241 [6381-59-5]	**Sodium Potassium Tartrate Tetrahydrate ACS, 99%**	500 Gms 5 Kg	760 7250
88662 [6381-59-5]	**Sodium Potassium Tartrate Tetrahydrate ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	820 7590
33650 [137-40-6]	**Sodium Propionate pure, 99%**	500 Gms 5 Kg	920 8170
57003 [12507-68-5]	**Sodium Pyroantimonate extrapure, 60-65% Sb_2O_3**	100 Gms 250 Gms	2800 6000
58568 [5323-95-5]	**Sodium Ricinoleate extrapure, 45% solids** (Ricinoleic Acid Sodium Salt, Sodium Ricinolate)	500 Gms 1 Kg	400 720
43874 [54-21-7]	**Sodium Salicylate Reagent Grade, 99%** (Salicylic Acid Sodium Salt, 2-Hydrobenzoic Acid Sodium Salt, Sodium-2-Hydrobenzoate)	500 Gms	950
36590 [54-21-7]	**Sodium Salicylate ExiPlus, Multi-Compendial, 99.5%** (Salicylic Acid Sodium Salt, 2-Hydrobenzoic Acid Sodium Salt, Sodium-2-Hydrobenzoate) meets compendial specs of USP, BP, Ph.Eur	250 Gms 500 Gms	800 1400
62399 [13410-01-0]	**Sodium Selenate Anhydrous extrapure AR, 99%**	25 Gms 100 Gms 500 Gms	580 1380 5750

Code	Name	Packing	Price
83627 [9063-38-1]	**Sodium Starch Glycolate ex. Maize (SSG) ExiPlus, Multi-Compendial** meets compendial specs of USP	500 Gms 5 Kg	290 2530
36313 [6106-21-4]	**Sodium Succinate Hexahydrate Dibasic extrapure, 99%**	100 Gms 500 Gms	220 800
87578 [6106-21-4]	**Sodium Succinate Hexahydrate Dibasic extrapure AR, 99%**	100 Gms 500 Gms	280 1100
59977 [7757-82-6]	**Sodium Sulphate Anhydrous ACS, 99%**	500 Gms 1 Kg 5 Kg 25 Kg	270 450 1960 7590
74725 [7757-82-6]	**Sodium Sulphate Anhydrous ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	330 2300
86566 [7757-82-6]	**Sodium Sulphate Anhydrous PrimeDx, 99%**	5 Kg 25 Kg	3220 14380
37233 [1317-33-5]	**Molybdenum Disulfide pure, 98%**	100 Gms 250 Gms 500 Gms	3200 7600 14500
25365 [27610-45-3]	**Sodium Sulphide Flakes extrapure, 60%**	500 Gms 5 Kg	350 2990
32246 [7757-83-7]	**Sodium Sulphite Anhydrous ACS, 98%**	500 Gms 5 Kg	410 3220
29814 [7757-83-7]	**Sodium Sulphite Anhydrous ExiPlus, Multi-Compendial, 95%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	460 3800
16638 ✔ [10102-20-2]	**Sodium Tellurite extrapure, 98%** (Disodium Trioxotellurate, Sodium Tellurite IV, Tellurous Acid Disodium Salt)	25 Gms 100 Gms	6500 21000
56580 [1330-43-4]	**Sodium Tetraborate (Borax Anhydrous) extrapure AR, 99%** (Disodium Tetraborate)	100 Gms 500 Gms 2.5 Kg	980 2070 9550
38631 [1303-96-4]	**Sodium Tetraborate Decahydrate (Borax Decahydrate) extrapure, 99%** (Disodium Tetraborate Decahydrate)	500 Gms 5 Kg	440 3910
19486 [1303-96-4]	**Sodium Tetraborate Decahydrate (Borax Decahydrate) ACS, 99.5%** (Disodium Tetraborate Decahydrate)	500 Gms 5 Kg	590 4890
77668 [1303-96-4]	**Sodium Tetraborate Decahydrate (Borax Decahydrate) ExiPlus, Multi-Compendial, 99%** (Disodium Tetraborate Decahydrate) meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	640 5750
29701 [12179-04-3]	**Sodium Tetraborate Pentahydrate (Borax Pentahydrate) extrapure AR, 99.5%** (Disodium tetraborate pentahydrate)	500 Gms 5 Kg	490 4600
58581 [13755-29-8]	**Sodium Tetrafluoroborate pure, 98.5%** (Sodium Fluoroborate, Sodium Fluoborate)	500 Gms	1020
54285 [143-66-8]	**Sodium Tetraphenyl Borate (Kalignost) ACS 99.5%** (Sodium Tetraphenyl Boron)	10 Gms 25 Gms 100 Gms 1 Kg	1270 2760 10010 74520
30888 [143-66-8]	**Sodium Tetraphenyl Borate (Kalignost) ExiPlus, Multi-Compendial, 99.5%** (Sodium Tetraphenyl Boron) meets compendial specs of USP, BP, Ph.Eur	10 Gms 100 Gms 1 Kg	1380 10470 77050
94408 [540-72-7]	**Sodium Thiocyanate pure, 98%**	500 Gms	560
21444 [540-72-7]	**Sodium Thiocyanate ACS, 98%**	500 Gms	690
89728 [7772-98-7]	**Sodium Thiosulphate Anhydrous extrapure, 99%**	500 Gms 5 Kg	440 4000
52304 [10102-17-7]	**Sodium Thiosulphate Pentahydrate extrapure, 99%**	500 Gms 5 Kg 25 Kg	220 2000 9000
49479 [10102-17-7]	**Sodium Thiosulphate Pentahydrate ACS, 99.5%**	500 Gms 5 Kg 25 Kg	270 2500 12000
62148 [10102-17-7]	**Sodium Thiosulphate Pentahydrate ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	290 2750
33737 [7772-98-7]	**Sodium Thiosulphate 0.1N Solution**	500 ml 1000 ml	190 330
40023 [7785-84-4]	**Sodium Trimetaphosphate (STMP) extrapure AR, 95%** (Trisodium Trimetaphosphate, Sodium Cyclotriphosphate)	100 Gms 500 Gms	1380 3220
85124 [7758-29-4]	**Sodium Tripolyphosphate Anhydrous (STPP) extrapure, 57-59% as P2O5** (Pentasodium Tripolyphosphate, Sodium Triphosphate Pentabasic)	1 Kg 5 Kg	650 3100
81499 [10213-10-2]	**Sodium Tungstate Dihydrate extrapure, 98%**	100 Gms 250 Gms	6200 15000
49110 [10213-10-2]	**Sodium Tungstate Dihydrate ACS, 99%**	100 Gms 250 Gms	6800 16000
46215 [10213-10-2]	**Sodium Tungstate Dihydrate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 250 Gms	7200 17000
85446 [2538-85-4]	**Solochrom Dark Blue (Calcon), 50%**	5 Gms 25 Gms 100 Gms	400 1140 2940

Code	Name	Packing	Price
59985 [110-44-1]	**Sorbic Acid pure, 98%**	500 Gms	920
		1 Kg	1610
		5 Kg	7480
62016 [1338-39-2]	**Sorbitan Monolaurate (SPAN 20®) extrapure**	100 ml	2300
		500 ml	6900
		2500 ml	32200
78338 [1338-43-8]	**Sorbitan Monooleate (SPAN 80®) extrapure**	500 ml	1030
		2500 ml	4830
36179 [1338-41-6]	**Sorbitan Monostearate (SPAN 60®) extrapure**	500 Gms	870
71440 △ [8007-43-0]	**Sorbitan Sesquioleate (SPAN 83) extrapure** (Sorbitan Sesqui Oleate)	250 ml	2070
		1000 ml	6900
14281 [50-70-4]	**D-Sorbitol Powder extrapure, 98%**	500 Gms	610
		1 Kg	1200
		5 Kg	5750
33320 [50-70-4]	**D-Sorbitol Powder for tissue culture, 98%**	500 Gms	710
		1 Kg	1270
63443 [50-70-4]	**D-Sorbitol Powder ExiPlus, Multi-Compendial, 97%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	720
		1 Kg	1380
		5 Kg	6150
99103 [50-70-4]	**D-Sorbitol Powder for molecular biology, 98%**	500 Gms	1740
		1 Kg	3330
		5 Kg	11130
47886 [23647-14-5]	**SPADNS ACS**	1 Gms	580
		5 Gms	1500
		25 Gms	4600
30467 [7646-78-8]	**Stannic Chloride anhydrous extrapure, 99%** (Tin (IV) Chloride, Tin Tetrachloride)	100 ml	2990
		500 ml	13800
79862 ✔ [7646-78-8]	**Stannic Chloride 1M in Dichloromethane** (Tin (IV) Chloride/Tin Tetrachloride in 1M MDC)	100 ml	1430
		500 ml	5180
60006 🌡 [7646-78-8]	**Stannic Chloride 1M in Heptane** (Tin (IV) Chloride/Tin Tetrachloride in 1M Heptane)	100 ml	2980
98017 ✔ [7646-78-8]	**Stannic Chloride 1M in Toluene** (Tin (IV) Chloride/Tin Tetrachloride in 1M Toluene)	100 ml	2980
91068 [10025-69-1]	**Stannous Chloride Dihydrate pure, 97%** (Tin (II) Chloride Dihydrate)	100 Gms	1110
		250 Gms	2710
		500 Gms	5150
23329 [10025-69-1]	**Stannous Chloride Dihydrate ACS, 98%** (Tin (II) Chloride Dihydrate)	100 Gms	1240
		250 Gms	2830
		500 Gms	5400
89673 [10025-69-1]	**Stannous Chloride Dihydrate ExiPlus, Multi-Compendial, 98%** (Tin (II) Chloride Dihydrate) meets compendial specs of USP, BP, Ph.Eur	100 Gms	1330
		250 Gms	3020
		500 Gms	5520
58254 [9005-25-8]	**Starch Corn pure**	500 Gms	260
		5 Kg	2300
86106 [9005-84-9]	**Starch Hydrolysed extrapure**	100 Gms	1190
		500 Gms	5590
14418 [9005-84-9]	**Starch Soluble extrapure**	500 Gms	1330
		5 Kg	12420
64698 [9005-84-9]	**Starch Soluble ACS**	500 Gms	1440
		5 Kg	13570
76601 [9005-84-9]	**Starch Soluble ExiPlus, Multi-Compendial** meets compendial specs of USP, BP, Ph.Eur	500 Gms	1550
		5 Kg	14720
91944	**Strontium AAS Standard Solution in 0.5N HCl (Traceable to NIST)**	100 ml	2070
		250 ml	3340
66731 [10476-81-0]	**Strontium Bromide (Strontium Dibromide) extrapure, 99%**	5 Gms	2300
		25 Gms	5750
		100 Gms	13800
67911 [1633-05-2]	**Strontium Carbonate extrapure, 98%**	500 Gms	750
61832 [10025-70-4]	**Strontium Chloride Hexahydrate extrapure, 98%**	250 Gms	450
		500 Gms	580
		5 Kg	5520
60979 [1311-10-0]	**Strontium Hydroxide Octahydrate pure, 98%**	500 Gms	1040
		5 Kg	9780
29563 [10042-76-9]	**Strontium Nitrate Anhydrous extrapure, 98%**	500 Gms	690
96330 [110-15-6]	**Succinic Acid extrapure, 99%**	500 Gms	660
		5 Kg	5980
85752 [110-15-6]	**Succinic Acid ACS, 99%**	100 Gms	490
		500 Gms	920
		5 Kg	8630
62321 [110-15-6]	**Succinic Acid ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	530
		500 Gms	1100
10917 [123-56-8]	**Succinimide (Butanimide) extrapure AR, 99%**	100 Gms	560
		500 Gms	2020
12243 [56038-13-2]	**Sucralose extrapure, 98%** (1,6-Dichloro-1,6-Dideoxy-ß-D-Fructofuranosyl-4-Chloro-4-Deoxy-a-D-Galactopyranoside,Trichlorosucros	25 Gms	2000
		100 Gms	5250
84973 [57-50-1]	**Sucrose pure**	500 Gms	290
		5 Kg	1860
		25 Kg	8720
90701 [57-50-1]	**Sucrose ACS**	500 Gms	440
		1 Kg	820
		5 Kg	3930
		25 Kg	17440
27580 [57-50-1]	**Sucrose for molecular biology**	500 Gms	600
		1 Kg	1190
		5 Kg	4520

Code	Name	Packing	Price
92696 [57-50-1]	**Sucrose for tissue culture**	500 Gms 1 Kg 5 Kg	440 770 3270
82958 [57-50-1]	**Sucrose IP (Non-Parenteral Use)**	1 Kg 5 Kg 25 Kg	770 3750 18000
38142 [57-50-1]	**Sucrose BP (Non-Parenteral Use)**	1 Kg 5 Kg 25 Kg	770 3750 18000
56799 [57-50-1]	**Sucrose NF (USP) (Non-Parenteral Use)**	1 Kg 5 Kg 25 Kg	1200 5500 25000
57335 [85-83-6]	**Sudan IV for microscopy (C.I. No. 26105), 95%**	100 Gms 500 Gms	600 2400
47256 [4197-25-5]	**Sudan Black**	10 Gms 25 Gms 100 Gms	600 1200 4000
24073 [5329-14-6]	**Sulphamic Acid pure, 99%** (Amido Sulphonic Acid)	500 Gms 5 Kg	350 2870
42773 [5329-14-6]	**Sulphamic Acid ACS, 99.3%** (Amido Sulphonic Acid)	100 Gms 500 Gms	460 2070
81645 [5329-14-6]	**Sulphamic Acid ExiPlus, Multi-Compendial, 99.3%** (Amido Sulphonic Acid) meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	580 2530
89762 [63-74-1]	**Sulphanilamide pure, 99%** (p-Aminobenzenesulphonamide)	100 Gms 500 Gms 5 Kg	800 3000 28000
65404 [63-74-1]	**Sulphanilamide extrapure AR, 99%** (p-Aminobenzenesulphonamide)	100 Gms 500 Gms	900 4090
19689 [63-74-1]	**Sulphanilamide extrapure AR, ExiPlus, Multi-Compendial, 99%** (p-Aminobenzenesulphonamide) meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	1070 4470
15354 [121-57-3]	**Sulphanilic Acid pure, 99%** (p-Aminobenzenesulphonamide)	100 Gms 500 Gms	200 600
93360 [121-57-3]	**Sulphanilic Acid (High Purity) ACS, 98%**	100 Gms 500 Gms	400 1800
97979 [121-57-3]	**Sulphanilic Acid (High Purity) ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	600 2200
44445 [5965-83-3]	**5-Sulphosalicylic Acid Dihydrate extrapure, 99%**	250 Gms 500 Gms 5 Kg	560 950 9000
13100 [5965-83-3]	**5-Sulphosalicylic Acid Dihydrate ACS, 99%** (2-Hydroxy-5-Sulfobenzoic Acid Dihydrate)	100 Gms 500 Gms	430 1300
92635 [14807-96-6]	**Talc Powder Fine ExiPlus, Multi-Compendial, -325 mesh** (Hydrous Magnesium Silicate, Talcum Powder, Soapstone) meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	170 1310
92101 [1401-55-4]	**Tannic Acid ExiPlus, Multi-Compendial** (Gallotannin, Tannin) meets compendial specs of USP, BP, Ph.Eur	100 Gms 250 Gms 500 Gms 5 Kg	1210 2700 5100 48000
52931 [1934-21-0]	**Tartrazine, 87%**	100 Gms 1 Kg	380 3000
20428 [100-21-0]	**Terephthalic Acid pure, 98%**	500 Gms 5 Kg	380 3500
90444 [10534-59-5]	**Tetrabutylammonium Acetate (TBAA) extrapure, 98%** (N,N,N-Tributyl-1-Butanaminium Acetate, Tetra-n-Butylammonium Acetate)	25 Gms 100 Gms	2800 9000
20669 [1643-19-2]	**Tetrabutylammonium Bromide (TBAB) pure, 98%**	250 Gms 500 Gms 1 Kg	440 830 1540
57594 [1643-19-2]	**Tetrabutylammonium Bromide (TBAB) ACS, 98%**	25 Gms 100 Gms 500 Gms	550 1760 8250
54423 [1643-19-2]	**Tetrabutylammonium Bromide (TBAB) for HPLC, 99.5%**	25 Gms 100 Gms	1540 5720
42120 [1112-67-0]	**Tetrabutylammonium Chloride (TBAC) pure, 97%** (N,N,N-Tributyl-Butanaminium Chloride)	25 Gms 100 Gms	2200 5000
37256 [5574-97-0]	**Tetrabutylammonium Dihydrogen Phosphate (TBADP) extrapure AR, 98%** (Tetrabutylammonium Phosphate Monobasic)	5 Gms 25 Gms	3630 13200
78262 [5574-97-0]	**Tetrabutylammonium Dihydrogen Phosphate (TBADP) for HPLC, 99%** (Tetrabutylammonium phosphate monobasic)	5 Gms 25 Gms	5040 17150
84444 [5574-97-0]	**Tetrabutylammonium Dihydrogen Phosphate (TBADP) 0.1M in HPLC Water** (Tetrabutylammonium Phosphate Monobasic 0.1M in HPLC Water)	100 ml 500 ml	5720 27500
57968 [5574-97-0]	**Tetrabutylammonium Dihydrogen Phosphate (TBADP) 0.5M in HPLC Water** (Tetrabutylammonium Phosphate Monobasic 0.5M in HPLC Water)	10 ml 100 ml 500 ml	2200 8800 33000
53856 [5574-97-0]	**Tetrabutylammonium Dihydrogen Phosphate (TBADP) 1M in HPLC Water** (Tetrabutylammonium Phosphate Monobasic 1M in HPLC Water)	50 ml 100 ml 500 ml	8250 13200 36300

Code	Name	Packing	Price
50283 [87749-50-6]	**Tetrabutylammonium Fluoride Trihydrate (TBAF Trihydrate) extrapure, 98%**	25 Gms 100 Gms	1430 4180
30315 [429-41-4]	**Tetrabutylammonium Fluoride (TBAF) 1M Soln. in THF** (N,N,N-Tributyl-1-butanaminium Fluoride)	100 ml 500 ml	1600 5500
91958 [32503-27-8]	**Tetrabutylammonium Hydrogen Sulphate (TBAHS) pure, 98%**	100 Gms 250 Gms 500 Gms	690 1310 2310
83955 [32503-27-8]	**Tetrabutylammonium Hydrogen Sulphate (TBAHS) extrapure AR, 99%**	25 Gms 100 Gms 500 Gms	630 1810 8400
10778 [32503-27-8]	**Tetrabutylammonium Hydrogen Sulphate (TBAHS) extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 250 Gms	2170 5060
88383 [32503-27-8]	**Tetrabutylammonium Hydrogen Sulphate (TBAHS) for HPLC, 99.5%**	25 Gms 100 Gms 500 Gms	730 1980 9600
52162 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 0.1N Aq. Solution**	500 ml	1940
28342 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 0.1N in Isopropanol**	500 ml	2090
94836 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 1M Aq. Solution**	100 ml 500 ml	2530 11550
31916 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 1M in Methanol** (N,N,N-Tributyl-1-butanaminium hydroxide, TBAH, Tetra-n-butylammonium hydroxide)	100 ml 500 ml	2750 13200
41728 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 10% Aq. Solution**	100 ml 500 ml	1380 5720
15134 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 10% in Methanol**	100 ml 500 ml	1490 5500
41715 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 20% Aq. Solution**	100 ml 500 ml	1980 7700
54722 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 25% Aq. Solution**	100 ml 500 ml	2090 9020
68032 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 40% Aq. Solution**	100 ml 250 ml 500 ml 1000 ml	3410 8250 15400 26400
30064 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 54-56% Aq. Solution**	100 ml 500 ml	8250 33000
13656 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 25% in Methanol**	100 ml 250 ml 500 ml	2090 4730 9020
74212 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 40% in Methanol**	100 ml 250 ml 500 ml	3410 8250 15400
40878 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 0.1N in Methanol**	500 ml	1870
45205 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 0.1N in Methanol/Toluene**	500 ml	1680
44991 [2052-49-5]	**Tetrabutylammonium Hydroxide (TBAH) 20% in Methanol/Toluene**	100 ml 500 ml	2200 8500
28051 [311-28-4]	**Tetrabutylammonium Iodide (TBAI) extrapure AR, 99%**	25 Gms 100 Gms	1870 6600
62390 [311-28-4]	**Tetrabutylammonium Iodide (TBAI) extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	25 Gms 100 Gms	2200 6930
26158 [5593-70-4]	**Tetra-N-Butyl Orthotitanate (TNBT) extrapure, 23-24% TiO2** (Titanium(IV) n-Butoxide, Tetrabutyl Titanate, Orthotitanic Acid Tetrabutyester)	100 ml 500 ml 2500 ml	1040 2760 8050
69105 [127-18-4]	**Tetrachloroethylene pure, 99%** (Perchloroethylene)	500 ml 2500 ml 25 Ltr	680 2920 28000
25058 [127-18-4]	**Tetrachloroethylene ACS, 99%** (Perchloroethylene)	500 ml 2500 ml	860 3700
74438 [127-18-4]	**Tetrachloroethylene for UV spectroscopy, 99.5%** (Perchloroethylene)	500 ml 1000 ml	960 1740
48448 [127-18-4]	**Tetrachloroethylene GC-HS, 99.9%** (Perchloroethylene)	250 ml 1000 ml	790 2130
35793 [33685-54-0]	**1,1,2,2-Tetrachloroethane-d2 (TCE-d2) for NMR spectroscopy, 99.5 Atom %D** (1,2-Dideutero-1,1,2,2-Tetrachloroethane, TCE-d2)	1 Gms	8500
52435 [33685-54-0]	**1,1,2,2-Tetrachloroethane-d2 (TCE-d2) for NMR spectroscopy, 99.5 Atom %D** (1,2-Dideutero-1,1,2,2-Tetrachloroethane, TCE-d2)	5 ml	18000
70427 [14937-42-9]	**Tetradecylammonium Bromide (TDAB) extrapure AR, 99%** (Tetrakis(decyl) Ammonium Bromide)	5 Gms 25 Gms 100 Gms	3030 11430 38080
59256 [14937-42-9]	**Tetradecylammonium Bromide (TDAB) for HPLC, 99%** (Tetrakis(decyl) Ammonium Bromide)	5 Gms 25 Gms	3590 14000
67902 [71-91-0]	**Tetraethylammonium Bromide (TEAB) pure, 98%**	250 Gms 1 Kg 5 Kg	1240 4600 21840
22711 [56-34-8]	**Tetraethylammonium Chloride (TEAC) extrapure, 98%** (TEA Chloride)	100 Gms 500 Gms	2600 10500
64674 [16873-13-5]	**Tetraethylammonium Hydrogen Sulphate (TEAHS) for HPLC, 99%**	10 Gms 25 Gms 100 Gms	4260 8740 24640
93314 [78-10-4]	**Tetraethyl Orthosilicate (TEOS) extrapure, 98%, 28-28.8% SiO2** (Silicon Tetraethoxide, Tetraethoxysilane, Tetraethyl Silicate)	100 ml 500 ml 2500 ml	1010 1800 6950

Code	Name	Packing	Price
21957 [4368-51-8]	**Tetraheptyl Ammonium Bromide (THAB) extrapure, 99%**	5 Gms	2660
		25 Gms	9890
		100 Gms	35830
46986 [4368-51-8]	**Tetraheptyl Ammonium Bromide (THAB) for HPLC, 99%**	5 Gms	3810
		25 Gms	13440
		100 Gms	47040
82863 [109-99-9]	**Tetrahydrofuran (THF) (Stabilized) pure, 99%**	500 ml	640
		2500 ml	3100
		25 Ltr	29000
74497 [109-99-9]	**Tetrahydrofuran (THF) (Stabilized) ACS, 99%**	500 ml	700
		2500 ml	3300
83642 [109-99-9]	**Tetrahydrofuran (THF) (Stabilized) ExiPlus, Multi-Compendial, 99.8%** meets compendial specs of USP, BP, Ph.Eur	500 ml	730
		2500 ml	3500
32661 [109-99-9]	**Tetrahydrofuran (THF) GC-HS, 99.9%**	250 ml	600
		1000 ml	2000
94842 [109-99-9]	**Tetrahydrofuran (THF) for HPLC & UV Spectroscopy, 99.9%**	500 ml	1100
		1000 ml	1800
		2500 ml	4000
59167 [109-99-9]	**Tetrahydrofuran (THF) Dried, 99.5%, 0.005% water**	500 ml	1000
		1000 ml	1700
92562 [109-99-9]	**Tetrahydrofuran (THF) for molecular biology, 99.8%**	250 ml	600
		500 ml	1000
22698 [1693-74-9]	**Tetrahydrofuran-d8 (THF-d8) for NMR spectroscopy, 99.5 Atom %D**	5 ml	24000
		10 ml	42000
11088 [119-64-2]	**1,2,3,4-Tetrahydronaphthalene (Tetralin) pure, 98%**	500 ml	1400
		2500 ml	5830
98727 [80526-82-5]	**Tetramethylammonium Hydrogen Sulphate (TMAHS) for HPLC, 99%** (Tetramethylammonium Hydrogensulfate)	5 Gms	730
		25 Gms	2920
		100 Gms	10640
32391 [103812-00-6]	**Tetramethylammonium Hydrogen Sulphate Monohydrate (TMAHS Monohydrate) extrapure, 98%** (N,N,N-Trimethyl Methanaminium Sulfate Hydrate)	25 Gms	900
		100 Gms	2800
91496 [75-59-2]	**Tetramethylammonium Hydroxide (TMAH) 10% Aq. Solution**	50 ml	910
		250 ml	3480
94620 [75-59-2]	**Tetramethylammonium Hydroxide (TMAH) 25% Aq. Solution**	100 ml	3140
		250 ml	5940
		500 ml	10640
21371 [75-59-2]	**Tetramethylammonium Hydroxide (TMAH) 25% in Methanol**	100 ml	3360
		250 ml	5600
		500 ml	9520
38342 [75-59-2]	**Tetramethylammonium Hydroxide (TMAH) 25% Aq. Solution, VLSI Semiconductor Grade**	1000 ml	23520
		2500 ml	43120
42209 [75-59-2]	**Tetramethylammonium Hydroxide (TMAH) 0.1N in Isopropanol/Methanol**	100 ml	1570
		500 ml	3480

Code	Name	Packing	Price
56874 [75-59-2]	**Tetramethylammonium Hydroxide (TMAH) 1M Aq. Solution ACS**	100 ml	1240
		250 ml	2920
29740 [10424-65-4]	**Tetramethylammonium Hydroxide Pentahydrate (TMAH Pentahydrate) extrapure, 98%**	5 Gms	930
		25 Gms	1640
		100 Gms	5520
84666 [110-18-9]	**N,N,N,N-Tetramethyl Ethylenediamine (TEMED) extrapure AR, 99%**	100 ml	950
		250 ml	1560
		500 ml	2560
84198 [110-18-9]	**N,N,N,N-Tetramethyl Ethylenediamine (TEMED) extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of Ph.Eur	100 ml	1080
		500 ml	2830
52145 [110-18-9]	**N,N,N,N-Tetramethyl Ethylenediamine (TEMED) for molecular biology, 99.5%**	25 ml	620
		100 ml	1250
80909 [5810-42-4]	**Tetrapropylammonium Chloride (TPAC) extrapure, 98%**	10 Gms	3360
		50 Gms	14560
70177 [7722-88-5]	**Tetrasodium Pyrophosphate Anhydrous (TSPP anhy) extrapure, 98%** (Sodium Diphosphate Tetrabasic, Sodium Pyrophosphate Tetrabasic)	500 Gms	840
		5 Kg	7620
49473 [13472-36-1]	**Tetrasodium Pyrophosphate Decahydrate (TSPP Decahyd) ACS, 99%** (Sodium Diphosphate Tetrabasic, Sodium Pyrophosphate Tetrabasic)	500 Gms	2020
		2.5 Kg	9190
56135 [58-55-9]	**Theophylline anhydrous pure, 99%**	25 Gms	340
		100 Gms	1010
		500 Gms	1910
17942 [62-55-5]	**Thioacetamide ACS, 99%**	25 Gms	560
		100 Gms	1910
		500 Gms	8400
48979 [62-55-5]	**Thioacetamide ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	25 Gms	650
		100 Gms	2130
		500 Gms	9750
42174 [100-68-5]	**Thioanisole extrapure, 99%** (Methyl Phenyl Sulfide Methylthiobenzene)	25 ml	1070
		100 ml	2800
		500 ml	10200
75241 [504-17-6]	**2-Thiobarbituric Acid (TBA) extrapure AR, 99%** (4,6-Dihydroxy-2-Mercaptopyrimidine, 4,6-Dihydroxypyrimidine-2-Thiol)	25 Gms	1240
		100 Gms	4260
75493 [2390-54-7]	**Thioflavin T for microscopy (C.I. No. 49005), 90%** (Basic Yellow 1)	5 Gms	1680
		25 Gms	3670
35357 [68-11-1]	**Thioglycolic Acid pure, 80%**	500 ml	1500
		2500 ml	6500
74118 [68-11-1]	**Thioglycolic Acid ExiPlus, Multi-Compendial, 80%** meets compendial specs of USP	500 ml	1650
		2500 ml	7200

Code	Name	Packing	Price
80599 [78338-22-4]	**Thionin Acetate, 90%**	5 Gms 25 Gms	2680 11140
74681 [79-19-6]	**Thiosemicarbazide pure, 98%**	25 Gms 100 Gms 500 Gms	370 1290 5940
58658 [79-19-6]	**Thiosemicarbazide extrapure AR, 99%**	25 Gms 100 Gms	590 1680
88196 [62-56-6]	**Thiourea pure, 98%**	500 Gms 5 Kg	500 4000
16994 [62-56-6]	**Thiourea ACS, 99%**	100 Gms 250 Gms 500 Gms	610 1200 2300
20491 [62-56-6]	**Thiourea ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	2600 24000
27179 [76-61-9]	**Thymol Blue ACS**	5 Gms 25 Gms 100 Gms	300 1010 3800
79617 [76-61-9]	**Thymol Blue ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	5 Gms 25 Gms	330 1210
40817	**Thymol Blue Solution (0.1% in methanol)**	125 ml	220
24474 [125-20-2]	**Thymolphthalein ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	5 Gms 25 Gms 100 Gms	240 1000 3800
19710	**Tin AAS Standard Solution in 2N HCl (Traceable to NIST)**	100 ml 250 ml	2100 3300
56283 [7440-31-5]	**Tin Metal Granular extrapure, 99.5%**	100 Gms 250 Gms 500 Gms	3000 6800 13000
24737 [7440-31-5]	**Tin Metal Granular ACS, 99.5%**	100 Gms 250 Gms 500 Gms	3200 7300 14250
42780 [7440-31-5]	**Tin Metal Powder ACS, 99.5%, ~325 mesh**	100 Gms 250 Gms	3000 6500
60001 △ [149-45-1]	**Tiron Monohydrate extrapure, 98%** (Catechol-3,5-Disulfonic Acid Disodium Salt Monohydrate) Tiron is used as a metal indictor of Iron (Fe), and as a colorimetric reagent of Fe, Al, Ti, and other metal ions.	5 Gms 25 Gms 100 Gms	600 2200 8000
37139 🌡 [31011-57-1]	**Titanium (IV) Chloride THF Complex** (Tetrachlorobis(Tetrahydrofuran) Titanium)	5 Gms 25 Gms	3130 12100
55788 [13463-67-7]	**Titanium Dioxide extrapure, 98%** (Titanium (IV) Oxide)	500 Gms 1 Kg 5 Kg	610 1140 4750
91482 [3087-36-3]	**Titanium Ethoxide (ET) extrapure, 33-35% TiO2** (Tetraethyl Titanate, Titanium (IV) Ethoxide, Tetraethyl Orthotitanate)	100 ml 250 ml 1000 ml	1240 2020 6720
56993 [546-68-9]	**Titanium Isopropoxide (TTIP) extrapure, 27.8-28.6% TiO2** (Tetraisopropyl Orthotitanate, Titanium (IV) Isopropoxide, Titanium Tetraisopropoxide)	100 ml 500 ml 2500 ml	690 1330 6180
70472 [1829-00-1]	**Titan Yellow, 65%**	25 Gms 50 Gms 250 Gms	750 1400 5000
47556 [108-88-3]	**Toluene pure, 99%**	500 ml 2500 ml 25 Ltr	330 1270 10930
39984 [108-88-3]	**Toluene (Sulphur free) extrapure, 99.5%**	500 ml 2500 ml 25 Ltr	350 1380 11390
85577 [108-88-3]	**Toluene ACS, 99.5%**	500 ml 2500 ml 25 Ltr	410 1500 12420
95227 [108-88-3]	**Toluene ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of BP, Ph.Eur	500 ml 2500 ml	460 1730
17213 [108-88-3]	**Toluene Dried, 99.5%, 0.005% water**	500 ml 1000 ml 2500 ml	450 860 1610
80156 [108-88-3]	**Toluene Dried Special (Low Benzene, 25ppm), 99.5%, 0.005% water**	250 ml 1000 ml	520 1270
83662 [108-88-3]	**Toluene GC-HS, 99.9%**	500 ml 1000 ml	870 1560
29490 [108-88-3]	**Toluene for HPLC & UV Spectroscopy, 99.8%**	500 ml 1000 ml 2500 ml	550 950 1670
85590 [2037-26-5]	**Toluene-d8 for NMR spectroscopy, 99.5 Atom %D**	10 ml 25 ml	21000 45000
33635 [108-44-1]	**m-Toluidine extrapure, 99%** (3-Methylaniline, m-Aminotoluene)	500 ml 2500 ml	660 2590
82225 [70-55-3]	**p-Toluenesulfonamide (PTSamide) pure, 98%** (4-Toluenesulfonamide)	100 Gms 500 Gms	360 1380
29423 [6192-52-5]	**p-Toluenesulphonic Acid Monohydrate pure, 98%** (Toluene-4-Sulphonic Acid, PTSA)	500 Gms 5 Kg	600 4540
77948 [6192-52-5]	**p-Toluenesulphonic Acid Monohydrate ACS, 98.5%** (Toluene-4-Sulphonic Acid, PTSA)	100 Gms 500 Gms	370 1480
55580 [6192-52-5]	**p-Toluenesulphonic Acid Monohydrate (High Purity) Reagent Grade, 99.5%** (Toluene-4-Sulphonic Acid, PTSA) meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	720 2580

Code	Name	Packing	Price
86763 [657-84-1]	**p-Toluenesulphonic Acid Sodium Salt pure, 98%** (Sodium p-Toluenesulfonate, PTSA Sodium)	500 Gms	2130
73222 [102-76-1]	**Triacetin pure, 99%**	500 ml	940
		2500 ml	4520
21559 [102-82-9]	**Tributylamine extrapure, 99%**	250 ml	270
		1000 ml	920
74315 [126-73-8]	**Tributylphosphate extrapure, 99%**	500 ml	2070
		2500 ml	9660
91094 [6138-23-4]	**D-Trehalose Dihydrate extrapure CHR**	5 Gms	370
		25 Gms	1280
		100 Gms	4250
		1 Kg	12500
15886 [6138-23-4]	**D-Trehalose Dihydrate for tissue culture**	5 Gms	550
		25 Gms	1950
		100 Gms	4730
10952 [6138-23-4]	**D-Trehalose Dihydrate ExiPlus Multi-Compendial** meets compendial specs of USP, BP, Ph.Eur	25 Gms	1500
		100 Gms	4800
		1 Kg	17000
90544 ✔ [76-03-9]	**Trichloroacetic Acid extrapure, 99%**	100 Gms	230
		500 Gms	750
60677 ✔ [76-03-9]	**Trichloroacetic Acid ACS, 99%**	100 Gms	340
		500 Gms	830
31445 [76-03-9]	**Trichloroacetic Acid 10% solution**	100 ml	230
		500 ml	500
54290 [76-03-9]	**Trichloroacetic Acid 20% solution**	500 ml	520
15106 [79-01-6]	**Trichloroethylene extrapure, 99%**	500 ml	500
		2500 ml	1800
		25 Ltr	16100
92390 [79-01-6]	**Trichloroethylene ACS, 99.5%**	500 ml	500
		2500 ml	1900
		25 Ltr	17900
86976 [79-01-6]	**Trichloroethylene ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of BP, Ph.Eur	500 ml	600
		2500 ml	2100
76338 [79-01-6]	**Trichloroethylene electronic grade, 99.5%**	500 ml	900
		2500 ml	3300
97681 [79-01-6]	**Trichloroethylene for UV spectroscopy, 99.5%**	500 ml	900
		1000 ml	1500
42943 [79-01-6]	**Trichloroethylene Dried, 99.5%, 0.005% water**	500 ml	1200
20779 [102-71-6]	**Triethanolamine pure, 98%**	500 ml	640
		2500 ml	2990
29306 [121-44-8]	**Triethylamine pure, 98%**	500 ml	700
		2500 ml	3200
		25 Ltr	30000
75198 [121-44-8]	**Triethylamine Reagent Grade, 99.5%**	500 ml	800
		2500 ml	3400
95205 [121-44-8]	**Triethylamine ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	1000 ml	1800
		2500 ml	3700
67352 [121-44-8]	**Triethylamine for HPLC, 99.5%**	250 ml	650
		500 ml	1100
		1000 ml	2000
55221 [121-44-8]	**Triethylamine for molecular biology, 99.5%**	500 ml	900
12475 [77-93-0]	**Triethyl Citrate extrapure, 99%** (Citric Acid Triethyl Ester, Ethyl Citrate)	250 ml	600
		1000 ml	2200
		2500 ml	5200
37817 [112-27-6]	**Triethylene Glycol extrapure, 99%**	500 ml	450
		2500 ml	1910
96475 [112-27-6]	**Triethylene Glycol Reagent Grade, 99%**	500 ml	590
		2500 ml	2240
86570 [112-24-3]	**Triethylenetetramine (mix. of isomers) pure, 97%**	500 ml	2020
		2500 ml	8400
81085 [122-51-0]	**Triethyl Orthoformate (TEOF) extrapure, 99%**	500 ml	900
		2500 ml	4000
78984 [76-05-1]	**Trifluoroacetic Acid (TFA) pure, 99%**	100 ml	650
		500 ml	2500
		2500 ml	12000
65415 [76-05-1]	**Trifluoroacetic Acid (TFA) Reagent Grade, 99%**	100 ml	900
		250 ml	1800
		500 ml	2800
67869 [76-05-1]	**Trifluoroacetic Acid (TFA) ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	100 ml	1050
		500 ml	3200
82818 [76-05-1]	**Trifluoroacetic Acid (TFA) for HPLC & UV Spectroscopy, 99.9%**	25 ml	840
		100 ml	1400
		500 ml	4200
40801 [76-05-1]	**Trifluoroacetic Acid (TFA) for molecular biology, 99.9%**	100 ml	1700
		500 ml	5000
23947 [76-05-1]	**Trifluoroacetic Acid (TFA) for HPLC & Sequencing, 99.9%**	10 x1 ml	18000
78716 [599-00-8]	**Trifluoroacetic Acid-d (TFA-d1) for NMR spectroscopy, 99.5 Atom %D**	5 ml	6200
		10 ml	9500
82395 [4574-04-3]	**Trimethyltetradecylammonium Chloride (TMTDAC) extrapure, 99%**	25 Gms	1500
		100 Gms	4200
38476 [519-73-3]	**Triphenylmethane pure, 98%**	100 Gms	1210
		500 Gms	4960
14393 ✔ [512-56-1]	**Trimethyl Phosphate (TMeP) extrapure, 98%** (Phosphoric Acid Trimethyl Ester, TMPH, TMPO)	100 ml	900
		500 ml	1500
		2500 ml	7000
71033 [77-86-1]	**Tris Buffer ACS for molecular biology, 99.9%** (Tris (Hydroxymethyl) Aminomethane, Tromethamine)	100 Gms	880
		500 Gms	4150
		5 Kg	34880

Code	Name	Packing	Price
37969 [77-86-1]	**Tris Buffer Superior for molecular biology, 99.9%** (Tris (Hydroxymethyl) Aminomethane, Tromethamine)	100 Gms 500 Gms 5 Kg	1110 5240 44690
56995 [77-86-1]	**Tris Buffer for HPLC, 99.9%** (Tris (Hydroxymethyl) Aminomethane, Tromethamine)	100 Gms 500 Gms	1530 6540
16622 [77-86-1]	**Tris Buffer USP, 99-101% (Non Injectable Use)** (Tris (Hydroxymethyl) Aminomethane, Tromethamine)	1 Kg 5 Kg	8000 37500
56489 [77-86-1]	**Tris Buffer PrimeDx, 99.9%** (Tris (Hydroxymethyl) Aminomethane, Tromethamine)	1 Kg 5 Kg 25 Kg	10000 42000 180000
99438 [1185-53-1]	**Tris Hydrochloride (Tris HCl) extrapure AR, 99%** (Tris (Hydroxymethyl) Aminomethane Hydrochloride)	100 Gms 500 Gms 5 Kg	1090 5240 49050
89781 [1185-53-1]	**Tris Hydrochloride (Tris HCl) for molecular biology, 99%** (Tris (Hydroxymethyl) Aminomethane Hydrochloride)	100 Gms 500 Gms 1 Kg	1370 6000 9810
96922 [1185-53-1]	**Tris Hydrochloride (Tris HCl) PrimeDx, 99%** (Tris (Hydroxymethyl) Aminomethane Hydrochloride)	1 Kg 5 Kg	13000 55000
30190 [9002-93-1]	**Triton X-100 extrapure for scintillation**	500 ml 2500 ml	660 3170
64518 [9002-93-1]	**Triton X-100 for molecular biology**	100 ml 500 ml 2500 ml	880 1310 5500
93400 [9002-93-1]	**Triton X-100 ExiPlus, Multi-Compendial** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	1050 4950
39857 [9002-93-1]	**Triton X-100 PrimeDx**	2500 ml 25 Ltr	8000 72000
82724 [9002-93-1]	**Triton X-114 extrapure** (TX-114, Octylphenoxypolyethoxyethanol, tert-Octylphenol Ethoxylate)	500 ml 2500 ml	1200 5000
55219 [9036-19-5]	**Triton X-405 (70% solution) extrapure** (Octylphenolethoxylate,Polyethyleneglycol t-Octylphenylether,Polyoxyethylene(40)isooctylphenylether)	500 ml 1000 ml 2500 ml	2200 3600 8000
64315 [547-57-9]	**Tropaeolin O Sodium Salt** (Acid Orange 6, Chrysoin)	25 Gms 100 Gms	640 2290
10233 [54-12-6]	**DL-Tryptophan extrapure CHR, 99%**	25 Gms 100 Gms 500 Gms	1500 4800 21000
10263 [73-22-3]	**L-Tryptophan extrapure CHR, 99%**	100 Gms 1 Kg 5 Kg	750 6200 30000
49102 [73-22-3]	**L-Tryptophan for tissue culture, 99%**	25 Gms 100 Gms	800 2200
48288 [73-22-3]	**L-Tryptophan ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	100 Gms 1 Kg	900 6700
26579 [153-94-6]	**D-Tryptophan (base) extrapure, 99%**	10 Gms 25 Gms 100 Gms	2800 5800 15000
10726	**Tungsten AAS Standard Solution in H20 (Traceable to NIST)**	100 ml 250 ml	1000 2000
71290 [12027-43-9]	**Tungstosilicic Acid Hydrate extrapure AR** (Silicotungstic Acid Hydrate)	25 Gms 100 Gms	3600 13000
67134 [12027-43-9]	**Tungstosilicic Acid Hydrate extrapure AR, ExiPlus, Multi-Compendial** (Silicotungstic Acid Hydrate) meets compendial specs of BP, Ph.Eur	25 Gms 100 Gms	3800 14000
70600 [556-02-5]	**D-Tyrosine extrapure, 98%**	5 Gms 25 Gms	2000 7000
10342 [556-03-6]	**DL-Tyrosine extrapure CHR, 99%**	5 Gms 25 Gms 100 Gms	1500 3240 12000
16796 [60-18-4]	**L-Tyrosine extrapure CHR, 99%**	25 Gms 100 Gms 500 Gms 1 Kg	350 850 4000 7000
18917 [60-18-4]	**L-Tyrosine for tissue culture, 99%**	25 Gms 100 Gms 500 Gms	500 1400 5000
91772 [60-18-4]	**L-Tyrosine ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	25 Gms 100 Gms 500 Gms	400 1050 4300
62762 [57-13-6]	**Urea extrapure, 99%**	500 Gms 5 Kg	350 3350
69120 [57-13-6]	**Urea ACS, 99.5%**	500 Gms 5 Kg 25 Kg	620 5800 27000
86854 [57-13-6]	**Urea ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP	500 Gms 5 Kg	720 6200
21113 [57-13-6]	**Urea for molecular biology, 99.5%**	500 Gms 1 Kg 5 Kg	950 1700 7890
44255 [57-13-6]	**Urea for tissue culture, 99.5%**	500 Gms 1 Kg	600 1100
12859 [72-18-4]	**L-Valine extrapure CHR, 99%**	25 Gms 100 Gms 1 Kg	300 700 4300
55269 [72-18-4]	**L-Valine for tissue culture, 99%**	5 Gms 100 Gms 500 Gms	410 1250 4980

Code	Name	Packing	Price
42749	**Vanadium AAS Standard Solution in 2N H2SO4 (Traceable to NIST)**	100 ml 250 ml	1800 2900
76132 [1314-62-1]	**Vanadium Pentoxide extrapure AR, 99.5%**	25 Gms 100 Gms 500 Gms	350 1150 5400
83822 [1314-62-1]	**Vanadium Pentoxide extrapure AR, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP	25 Gms 100 Gms 500 Gms	400 1280 5650
83080 [121-34-6]	**Vanillic Acid extrapure, 99%** (4-Hydroxy-3-Methoxybenzoic Acid)	25 Gms 100 Gms	1500 4900
98842 [121-33-5]	**Vanillin pure, 99%** (4-hydroxy-3-Methoxybenzaldehyde)	100 Gms 500 Gms	500 1800
20563 [121-33-5]	**Vanillin extrapure AR, 99%** (4-hydroxy-3-Methoxybenzaldehyde)	25 Gms 100 Gms 500 Gms	300 550 2100
49456 [121-33-5]	**Vanillin extrapure AR, ExiPlus, Multi-Compendial, 99%** (4-hydroxy-3-Methoxybenzaldehyde) meets compendial specs of BP	25 Gms 100 Gms 500 Gms	340 650 2400
14553 ✔ [127-47-9]	**Vitamin A Acetate (Retinyl Acetate) for tissue culture** (Retinol Acetate)	5 Gms 25 Gms 100 Gms	650 2200 6000
76911 ✔ [79-81-2]	**Vitamin A Palmitate (Retinyl Palmitate) for tissue culture, 1.615-1.860MIU/g**	5 Gms 25 Gms	850 3640
60514 ■ [68-26-8]	**Vitamin A (Retinol), 90%** (Axerophthol, Vitamin A1)	5 Gms 25 Gms	2500 5000
92605 [7732-18-5]	**Water for HPLC**	1000 ml 2500 ml	400 900
21897 [7732-18-5]	**Water for LC-MS** (DI Water, Deionized Water, Distilled Water)	1 Ltr (Wide-Mouth) 2.5 Ltr (Wide-Mouth) 4 Ltr (Wide-Mouth)	1800 6000 12000
58348 [7732-18-5]	**Water Distilled**	5 Ltr	500
62104	**Winkler`s Solution 1**	100 ml 500 ml	980 4850
64236 [11138-66-2]	**Xanthan Gum ExiPlus, Multi-Compendial** meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	430 1500
23749 [95-47-6]	**o-Xylene pure, 98%** (1,2-Dimethylbenzene)	500 ml 2500 ml	420 1580
62231 [95-47-6]	**o-Xylene Reagent Grade, 99%** (1,2-Dimethylbenzene)	500 ml 2500 ml	650 3100
29280 [108-38-3]	**m-Xylene pure, 98%** (1,3-Dimethylbenzene)	500 ml	620
24923 [106-42-3]	**p-Xylene pure, 99%** (1,4-Dimethylbenzene)	500 ml 2500 ml	640 2790
54717 [1330-20-7]	**Xylene pure, 98%** (mixture of isomers)	500 ml 2500 ml 25 Ltr	360 1480 10850
63797 [2650-17-1]	**Xylene Cyanol FF** (Acid Blue 147, C.I. No. 42135)	5 Gms 25 Gms 100 Gms	1300 5200 17690
75122 [2650-17-1]	**Xylene Cyanol FF extrapure AR, ExiPlus, Multi-Compendial** (Acid Blue 147, C.I. No. 42135) meets compendial specs of USP	5 Gms 25 Gms 100 Gms	1500 5600 20410
49672 [125-31-5]	**p-Xylenol Blue extrapure AR, 95%** (p-Xylenosulphonaphthalein)	1 Gms 5 Gms	550 2550
49992 [3618-43-7]	**Xylenol Orange Tetrasodium Salt ACS**	5 Gms 10 Gms	520 900
31557 [3618-43-7]	**Xylenol Orange Tetrasodium Salt ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	10 Gms 25 Gms	1100 2600
96240 [1330-20-7]	**Xylene Sulphur Free extrapure, 98%** (mixture of isomers)	500 ml 2500 ml 25 Ltr	380 1520 11210
90998 [1330-20-7]	**Xylene ACS, 98.5%**	500 ml 2500 ml	450 1640
89159 [1330-20-7]	**Xylene ExiPlus, Multi-Compendial, 98.5%** meets compendial specs of USP, BP, Ph.Eur	500 ml 2500 ml	500 1880
45122 [1330-20-7]	**Xylene for molecular biology, 98.5%**	250 ml 500 ml 2500 ml	410 730 2180
85471 [1330-20-7]	**Xylene Dried, 99.5%, 0.005% water**	1000 ml	710
46561 [1314-36-9]	**Yttrium (III) Oxide extrapure AR, 99.9%**	10 Gms 25 Gms	950 2000
76192 [5970-45-6]	**Zinc Acetate Dihydrate extrapure, 98.5%**	500 Gms	750
76205 [5970-45-6]	**Zinc Acetate Dihydrate ACS, 99.5%**	500 Gms	870
53174 [5970-45-6]	**Zinc Acetate Dihydrate ACS, ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms	980
43758	**Zinc AAS Standard Solution in 1N HNO3 (Traceable to NIST)**	100 ml 250 ml	2070 3340
97765 [5263-02-5, 51839-25-9]	**Zinc Carbonate Basic pure, 58% Zn** (Zinc Hydroxide Carbonate, Carbonic Acid Zinc Salt Basic)	500 Gms	1150
99367 [5990-32-9]	**Zinc Citrate Dihydrate pure, 98%** (Citric Acid Zinc Salt)	500 Gms	1040
87288 [7646-85-7]	**Zinc Chloride pure, 97%**	500 Gms 5 Kg	580 4950

Code	Name	Packing	Price
80004 [7646-85-7]	**Zinc Chloride ACS, 98%**	500 Gms 5 Kg	1040 9660
67514 [7646-85-7]	**Zinc Chloride for tissue culture, 98%**	500 Gms 1 Kg	920 1610
86783 [7783-49-5]	**Zinc Fluoride pure, 98%**	100 Gms 500 Gms	460 1040
46621 [7440-66-6]	**Zinc Metal Dust pure, 95%, -325 mesh**	500 Gms 2.5 Kg	1040 4830
27008 [7440-66-6]	**Zinc Metal Dust pure, 99%, -325 mesh**	500 Gms 2.5 Kg	1250 6080
20050 [7440-66-6]	**Zinc Metal Dust ACS, 99.8%, -325 mesh**	500 Gms 2.5 Kg	1450 6840
73932 [7440-66-6]	**Zinc Metal Dust ExiPlus, Multi-Compendial, 99.5%, -325 mesh** meets compendial specs of BP, Ph.Eur	500 Gms 2.5 Kg	1590 7250
91229 [7440-66-6]	**Zinc Metal Granules, 99%, 3-12 mm**	500 Gms	1380
75001 ✔ [10196-18-6]	**Zinc Nitrate Hexahydrate extrapure, 98%**	500 Gms	690
24693 ✔ [10196-18-6]	**Zinc Nitrate Hexahydrate Reagent Grade, 98%**	500 Gms	810
87531 [1314-13-2]	**Zinc Oxide extrapure, 99%**	500 Gms 5 Kg	900 7920
67062 [1314-13-2]	**Zinc Oxide ACS, 99%**	500 Gms 5 Kg	930 8580
51508 [1314-13-2]	**Zinc Oxide ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	990 9240
13439 [557-05-1]	**Zinc Stearate extrapure, 10-15% ZnO** (Stearic acid zinc salt)	500 Gms 1 Kg	930 1590
67106 [7446-20-0]	**Zinc Sulphate Heptahydrate for tissue culture, 99%**	500 Gms 1 Kg	760 1340
76455 [7446-20-0]	**Zinc Sulphate Heptahydrate pure, 99%**	500 Gms 1 Kg 5 Kg	390 700 3370
78494 [13597-44-9]	**Zinc Sulphite Anhydrous extrapure, 80%**	500 Gms	1830
62294 [7488-52-0]	**Zinc Sulphite Dihydrate extrapure, 98%**	500 Gms	1760
47482 [1314-98-3]	**Zinc Sulphide extrapure, 98%**	100 Gms 500 Gms	960 3170
75738 [7446-20-0]	**Zinc Sulphate Heptahydrate ACS, 99.5%**	500 Gms 5 Kg	450 4100
35243 [7446-20-0]	**Zinc Sulphate Heptahydrate ExiPlus, Multi-Compendial, 99.5%** meets compendial specs of USP, BP, Ph.Eur	500 Gms 5 Kg	570 5150
12077 [13520-92-8]	**Zirconium Oxychloride Octahydrate extrapure, 99.5%**	100 Gms 250 Gms 500 Gms	500 900 1650
22047 [13520-92-8]	**Zirconium Oxychloride Octahydrate ACS, 99.5%**	100 Gms 250 Gms	1050 1800
38028 [14644-61-2]	**Zirconium Sulphate extrapure, 99.5%**	100 Gms 500 Gms	1100 4600
47985 [13826-66-9]	**Zirconyl Nitrate Hydrate ExiPlus, Multi-Compendial, 99.5%** (Zirconium Oxynitrate) meets compendial specs of USP, Ph.Eur	100 Gms 500 Gms	950 3600
26405 □ [116797-51-4]	**DL-Cysteine Hydrochloride Monohydrate extrapure, 98%** (H-DL-Cys-OH·HCl Monohydrate)	25 Gms 100 Gms	3200 5500
52647 △ [9016-45-9]	**Nonylphenol Ethoxylate (NPE) 9.5 Mole** (NP 9.5, Nonylphenyl, Nonidet P40, Octylphenyl Polyethylene Glycol)	100 ml 500 ml	700 2500

Part B

Bioreagents, Biochemicals & Speciality Fine Chemicals

PrimeDx™ An exclusive grade



PrimeDx™ is a specialized line of high-purity reagents crafted specifically for the diagnostic and molecular biology industries. Developed under stringent quality protocols, these reagents ensure maximum purity, consistency, and reliability, making them ideal for use in clinical diagnostics, molecular assays, and R&D applications.

This dedicated range of raw material chemicals is specially produced for diagnostic kit manufacturers, delivering application-specific performance and unmatched quality tailored to the demands of the diagnostic industry.

Key Features of PrimeDx™ Grade:

- Ensures high accuracy and reliability in diagnostic applications.
- Minimizes toxicological risks, especially crucial for sensitive diagnostic assays.
- Each batch meets stringent QC standards to maintain uniformity and reproducibility across diagnostic platforms.
- Designed to eliminate environmental and particulate contamination, preserving material integrity from production to use.

Applications:

- Clinical Chemistry
- Molecular Diagnostics
- ELISA & Immunoassays
- Genomic & Proteomic Analysis
- In-vitro Diagnostic Kit Manufacturing

PrimeDx™ Product Portfolio (Selected)

Prdt. Code	Product Name
69657	4-Aminoantipyrine (Ampyrone) PrimeDx, 99%
19575	Casein acc. to Hammarsten, PrimeDx, 95% Protein
53499	Dodecyltrimethylammonium Bromide (DOTAB) PrimeDx, 99%
61936	Dodecyltrimethylammonium Chloride (DOTAC) PrimeDx, 98%
64762	EDTA Disodium Salt Dihydrate PrimeDx, 98%
54479	L-Alanine PrimeDx, 99%
40495	L-Arginine (free base) PrimeDx, 99%
61310	MOPS Buffer PrimeDx, 99%
84868	p-Nitrophenylphosphate Disodium Salt Hexahydrate (pNPP) PrimeDx, 99%
86180	PIPES Buffer PrimeDx, 99%
26511	ß-Nicotinamide Adenine Dinucleotide Disodium Salt (Reduced) (ß-NADH.Na2, DPNH.Na2) PrimeDx, 98%
56489	Tris Buffer PrimeDx, 99.9%

Note: This is a selected listing. A comprehensive product range is available on request.

Packaging & Quality Assurance

Every PrimeDx™ batch comes with a Certificate of Analysis (CoA) to show it meets strict quality checks. It's made under SRL's high standards and packed in clean containers to keep it safe and pure all the way from the lab to where it's used.

Code	Name	Packing	Price
66308 [21293-29-8]	■ **(+)-Abscisic Acid Natural (S-ABA) for tissue culture, 98%** (Dormin, Abscissin II)	25 Mg 100 Mg 500 Mg	1470 4100 18900
31240 [14375-45-2]	■ **(+/-)-Abscisic Acid Synthetic (ABA) for tissue culture, 98%** (Dormin, Abscissin II)	25 Mg 100 Mg	1680 4310
40157 [30931-67-0]	● **ABTS (2,2-Azino-bis(3-Ethylbenzothiazoline-6-Sulfonic Acid) Diammonium Salt) extrapure AR , 98.5%**	100 Mg 500 Mg	2800 8500
28042 [30931-67-0]	● **ABTS (2,2-Azino-bis(3-Ethylbenzothiazoline-6-Sulfonic Acid) Diammonium Salt) extrapure, 98%**	500 Mg 1 Gms 5 Gms	4600 7000 26980
83088	**ABTS Peroxidase Stop Solution (5% SDS solution) suitable for molecular biology**	25 ml 100 ml	1630 4480
33108	🌡 **ABTS Substrate Solution (Single Solution)**	25 ml 100 ml	1630 3470
65457 [56180-94-0]	**Acarbose extrapure, 95%**	1 Gms 5 Gms	6000 20000
17965 [56180-94-0]	**Acarbose ExiPlus, Multi-Compendial, 95%** meets compendial specs of USP, BP, Ph.Eur	1 Gms 5 Gms	7500 25000
66822 [7365-82-4]	**ACES Buffer extrapure, 99%** (2-[(2-Amino-2-Oxoethylamino)Ethane] Sulphonic Acid)	25 Gms 100 Gms 250 Gms	1900 7400 18000
27907 [7365-82-4]	**ACES Buffer for molecular biology, 99.5%**	5 Gms 25 Gms 100 Gms	720 2700 9200
62316 [540-92-1]	**Acetone Sodium Bisulfite pure, 98%** (2-Hydroxy-2-Propanesulfonic Acid Monosodium Salt)	100 Gms 250 Gms	2000 3500
89561 [127-06-0]	**Acetoxime (Acetone Oxime) extrapure, 98%**	25 Gms 100 Gms 250 Gms	1400 3000 7200
89204 [30827-99-7]	■ **AEBSF Hydrochloride extrapure AR, 98%** (4-(2-Aminoethyl) benzenesulfonylfluoride Hydrochloride)	25 Mg 100 Mg 500 Mg	1500 3200 10000
59750 [132489-69-1]	■ **PUGNAc extrapure, 95%** (O-(2-Acetamido-2-Deoxy-D-Glucopyranosylidenamino) N-Phenylcarbamate)	5 Mg	16230
20140 [52574-08-0]	● **3-(Acetamidomethylthio)propanoic Acid (MPA(Acm)) extrapure, 98%**	1 Gms	5260
12052 [14691-89-5]	🌡 **4-Acetamido-2,2,6,6-Tetramethylpiperidine-1-Oxyl (Acetamido TEMPO) pure, 97%**	1 Gms 5 Gms 25 Gms 100 Gms	1950 3200 7000 20000
34437 [572-09-8]	■ **Acetobromo-D-Glucose Stab. 1% CaCO3 extrapure, 98%**	10 Gms 25 Gms	5200 12000
28590 [3068-32-4]	■ **Acetobromo-D-Galatose w/ 2% Calcium Carbonate extrapure, 98%**	5 Gms	6500

Code	Name	Packing	Price
43717 [941-98-0]	**1-Acetonaphthone pure, 98%** (Acetylnaphthalene, Naphthyl Ketone)	100 ml 250 ml	1370 2660
21009 [93-08-3]	**2-Acetonaphthone pure, 98%** (Acetylnaphthalene, Methyl-2-Naphthyl Ketone)	100 Gms 250 Gms	1300 3000
24168 [66-23-9]	■ **Acetylcholine Bromide extrapure, 98%**	5 Gms 25 Gms	2760 9200
17884 [97-69-8]	🌡 **N-Acetyl-L-Alanine extrapure, 99%**	5 Gms 25 Gms	1250 5040
70538 [997-55-7]	🌡 **N-Acetyl-L-Aspartic Acid extrapure, 98%**	5 Gms 25 Gms	2000 8500
23157 [6974-32-9]	🌡 **1-O-Acetyl-2,3,5-Tri-O-Benzoyl-b-D-Ribofuranose extrapure, 98%** (b-D-Ribofuranose 1-Acetate 2,3,5-Tribenzoate)	5 Gms 25 Gms	2500 8000
23598 [3080-30-6]	🌡 **1-O-Acetyl-2,3,5-Tri-O-Benzoyl-b-L-Ribofuranose extrapure, 98%** (b-L-Ribofuranose 1-Acetate 2,3,5-Tribenzoate)	1 Gms 5 Gms	12210 39540
41439 [14215-97-5]	🌡 **1-O-Acetyl-2,3,5-Tri-O-Benzoyl-D-Ribofuranose extrapure, 98%** (D-Ribofuranose 1-Acetate 2,3,5-Tribenzoate)	5 Gms 25 Gms	2450 11050
72981 [91110-24-6]	🌡 **1-O-Acetyl-2,3,5-Tri-O-Benzyl-b-D-Ribofuranose extrapure, 95%**	50 Mg 200 Mg	22680 50000
32993 [58381-23-0]	🌡 **1-O-Acetyl-2,3,5-Tri-O-Benzyl-D-Ribofuranose extrapure, 98%**	50 Mg 200 Mg	22680 50000
62929 [60-31-1]	🌡 **Acetylcholine Chloride extrapure, 99%**	25 Gms	6000
71332 [2260-50-6]	🌡 **Acetylcholine Iodide extrapure, 98%**	5 Gms 25 Gms	2900 10100
26904 [75520-41-1]	■ **Acetyl Coenzyme A Trilithium Salt extrapure, 90%**	5 Mg 25 Mg	7500 30000
47866 [616-91-1]	🌡 **N-Acetyl-L-Cysteine ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP	25 Gms 100 Gms 500 Gms	600 1600 6000
99328 [14215-68-0]	🌡 **N-Acetyl-D-Galactosamine extrapure, 98%**	100 Mg 1 Gms 5 Gms	2520 4000 12000
59012 [7512-17-6]	🌡 **N-Acetyl-D-Glucosamine extrapure AR, 99%**	10 Gms 100 Gms 500 Gms	600 2000 7000
50241 [1188-37-0]	🌡 **N-Acetyl-L-Glutamic Acid extrapure, 99%**	10 Gms 100 Gms	500 4280
68761 [2490-97-3]	🌡 **N-Acetyl-L-Glutamine extrapure, 99%**	25 Gms 100 Gms	1550 4880
56115 [543-24-8]	🌡 **N-Acetylglycine extrapure, 99%** (Acetamidoacetic Acid, Aceturic Acid)	25 Gms 100 Gms 500 Gms	1030 2800 6600

Code	Name	Packing	Price
76097 [1188-21-2]	**N-Acetyl-L-Leucine extrapure, 99%**	5 Gms	450
		25 Gms	1800
		100 Gms	5200
26124 [65-82-7]	**N-Acetyl-L-Methionine extrapure, 99%**	5 Gms	810
		25 Gms	3380
		100 Gms	11820
90430 [68-95-1]	**N-Acetyl-L-Proline extrapure, 98%**	5 Gms	1560
		25 Gms	5880
		100 Gms	15640
81977 [2466-76-4]	**N-Acetylimidazole extrapure AR, 98%**	25 Gms	1970
		100 Gms	6020
32414 [1115-47-5]	**N-Acetyl-DL-Methionine extrapure, 99%**	10 Gms	780
		100 Gms	1750
86419 [131-48-6]	**N-Acetylneuraminic Acid (NANA) extrapure, 98%**	100 Mg	1200
		1 Gms	3000
		5 Gms	5500
10197 [1122-62-9]	**2-Acetyl Pyridine extrapure, 99%**	25 ml	1850
		100 ml	6600
14507 [350-03-8]	**3-Acetyl Pyridine pure, 99%**	25 ml	1600
		100 ml	4000
38998 [1866-15-5]	**S-Acetylthiocholine Iodide extrapure AR, 99%**	1 Gms	2220
		5 Gms	5000
		25 Gms	14000
99703 [88-15-3]	**2-Acetyl Thiophene pure, 99%**	100 ml	1800
		500 ml	8100
11644 [1218-34-4]	**N-Acetyl-L-Tryptophan extrapure, 99%**	5 Gms	1420
		25 Gms	6120
48861 [1218-34-4]	**N-Acetyl-L-Tryptophan ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	5 Gms	1630
		25 Gms	6050
51031 [87-32-1]	**N-Acetyl-DL-Tryptophan extrapure, 99%**	25 Gms	2050
		100 Gms	7980
		1 Kg	75600
98070 [537-55-3]	**N-Acetyl-L-Tyrosine extrapure, 98%**	5 Gms	2080
		25 Gms	7940
16756 [36546-50-6]	**N-Acetyl-L-Tyrosine Ethyl Ester (ATEE) extrapure, 99%**	5 Gms	4200
		10 Gms	7000
		25 Gms	13000
61942 [96-81-1]	**N-Acetyl-L-Valine extrapure, 99%**	1 Gms	2700
97412 [10127-02-3]	**Acridine Orange hemi (Zinc Chloride) Salt** (C.I. No. 46005)	5 Gms	1100
		10 Gms	1850
		25 Gms	4190
37751 [10127-02-3]	**Acridine Orange hemi(Zinc Chloride) Salt for molecular biology**	5 Gms	1650
		10 Gms	2800
		25 Gms	5500
62255 [10127-02-3]	**Acridine Orange hemi(Zinc Chloride) Salt for tissue culture**	5 Gms	1580
		25 Gms	4800
82407 [65-61-2]	**Acridine Orange Hydrochloride Hydrate for molecular biology, 99%** (3,6-Bis(dimethylamino)Acridine Hydrochloride)	1 Gms	900
		5 Gms	2800
		25 Gms	12500
42070 [135-49-9]	**Acridine Yellow G** (Basic Yellow K, 3,6-Diamino-2,7-Dimethylacridine Hydrochloride)	5 Gms	1200
		25 Gms	4500
42428 [8048-52-0]	**Acriflavine Neutral for molecular biology, 98%** (3,6-Diamino-10-Methylacridinium Chloride, Acriflavine Base)	5 Gms	420
		25 Gms	1800
		100 Gms	6000
35188 [8063-24-9]	**Acriflavine Hydrochloride for molecular biology, 98.5%-105%** (3,6-Diamino-10-Methylacridinium Chloride Hydrochloride, Euflavine)	5 Gms	600
		25 Gms	2600
88834 [5165-97-9]	**2-Acrylamido-2-Methyl-1-Propanesulfonic Acid Sodium Salt Solution (Stabilized w/ MEHQ), 50% in Water** (NaATBS, Sodium AMPS 50% Solution, 2-Acrylamido-2-Methyl-1-Propanesulfonic Acid Sodium Salt Solution	250 ml	2000
		1000 ml	4500
39309 [51005-14-2]	**Actin extrapure, 99%**	1 Mg	34200
28269 [50-76-0]	**Actinomycin D (AMD) ex. Streptomyces Sp., 98%** (Actinomycin C1, Actinomycin IV, Dactinomycin)	2 Mg	5000
		5 Mg	8250
		10 Mg	14500
44597 [9012-37-7]	**Acylase 1 (ACY1) ex. Porcine Kidney, 1000U/mg**	50 Mg	2900
		250 Mg	10200
		1 Gms	31800
37690 [26239-55-4]	**ADA Buffer extrapure, 98%** (N-(2-Acetamido)Iminodiacetic Acid)	25 Gms	1200
		100 Gms	3000
		500 Gms	7400
83468 [7415-22-7]	**ADA Monosodium Salt Buffer extrapure, 98%**	5 Gms	4660
		25 Gms	9310
35609 [41689-31-0]	**ADA Disodium Salt Buffer extrapure, 98%**	25 Gms	3800
50300 [73-24-5]	**Adenine extrapure, 98%**	5 Gms	320
		25 Gms	1290
87583 [6509-19-9, 321-30-2]	**Adenine Hemisulphate Dihydrate extrapure, 99%** (Adenine Sulphate)	10 Gms	550
		25 Gms	1320
		100 Gms	4940
81951 [6509-19-9, 321-30-2]	**Adenine Hemisulphate Dihydrate for tissue culture, 99%** (Adenine Sulphate)	10 Gms	870
		25 Gms	2100
		100 Gms	7350
56603 [321-30-2]	**Adenine Hemisulphate Dihydrate ExiPlus, Multi-Compendial, 99%** (Adenine Sulphate) meets compendial specs of USP	25 Gms	2520
		100 Gms	8400
54865 [58-61-7]	**Adenosine extrapure, 99%**	5 Gms	350
		25 Gms	1340
		100 Gms	2500
94149 [58-61-7]	**Adenosine ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	25 Gms	1600
		100 Gms	3000

Code	Name	Packing	Price
84910 [37063-35-7] ■	**Adenosine-2',3' Cyclic Monophosphate Sodium Salt extrapure, 98%**	25 Mg 100 Mg	3510 10730
98333 [60-92-4] ■	**Adenosine-3',5'-Cyclic-Monophosphoric Acid extrapure, 98%**	100 Mg 500 Mg	2000 6000
43801 [16178-48-6] ■	**Adenosine-5-Diphosphate Disodium Salt (ADP-Na2) extrapure, 98%**	500 Mg 1 Gms 5 Gms 25 Gms	750 1400 5000 16000
99118 [72696-48-1] ■	**Adenosine-5-Diphosphate Monopotassium Salt Dihydrate (ADP-K) extrapure, 97%**	500 Mg 1 Gms 5 Gms	1760 3000 9000
90697 [4578-31-8] ■	**Adenosine-5-Monophosphate Disodium Salt (AMP.Na2) extrapure, 95%**	5 Gms 25 Gms 100 Gms 500 Gms	840 2000 7000 32000
84878 [34369-07-8] ■	**Adenosine-5'-Triphosphate Disodium Salt (ATP-Na2) extrapure vanadium free, 98%**	5 Gms 25 Gms 100 Gms	900 3200 12000
46600 [979-92-0] ■	**S-(5-Adenosyl)-L-Homocysteine extrapure, 99%** 5'-Deoxy-S-adenosyl-L-homocysteine	10 Mg	37210
85934 [97540-22-2] ■	**(S)Adenosyl-L-Methionine DisulphateTosylate, 98%**	500 Mg 1 Gms	1840 3470
40257 [1071-93-8] ■	**Adipic Acid Dihydrazide (ADH) extrapure, 98%** (Hexanedioic Acid Dihydrazide)	25 Gms 100 Gms 250 Gms	3000 4700 7200
75101 [488-81-3]	**Adonitol (Ribitol) extrapure, 99%**	5 Gms 25 Gms	1800 7550
37575 [51-42-3] 🌡	**Adrenaline Bitartarate (Epinephrine Bitartrate) extrapure, 98%** (L-Adrenaline Hydrogen Bitartrate, L-Epinephrine Hydrogen Bitartrate)	1 Gms 5 Gms	2300 8500
59783 [51-43-4] 🌡	**L-Adrenaline (L-Epinephrine) extrapure, 98%**	1 Gms 5 Gms	3000 10250
60645 [9012-36-6]	**Agarose High EEO for molecular biology**	25 Gms 100 Gms 500 Gms	1840 6610 31820
36601 [9012-36-6]	**Agarose Low EEO for molecular biology**	10 Gms 25 Gms 100 Gms 500 Gms 1 Kg	680 1350 4290 19580 38550
23287 [9012-36-6]	**Agarose Low EEO Superior Grade for molecular biology**	10 Gms 25 Gms 100 Gms 500 Gms 1 Kg	800 1720 5020 22640 42830
10423 [9012-36-6]	**Agarose Medium EEO for molecular biology**	10 Gms 25 Gms 100 Gms 500 Gms	900 2200 5800 25700
27813 [9012-36-6]	**Agarose UltraPhor for molecular biology**	10 Gms 25 Gms 100 Gms 500 Gms	6200 11000 40000 185000
32417 [9012-36-6]	**Agarose Low Melting for molecular biology**	5 Gms 10 Gms 25 Gms 100 Gms 250 Gms	3000 3500 8000 30000 70000
93419	**Agarose TAE mix 1.0%**	100 Gms	8530
67900	**Agarose TAE mix 2.0%**	100 Gms	8530
76171	**Agarose TBE mix 1.0%**	100 Gms	8530
12375	**Agarose TBE mix 2.0%**	100 Gms	8530
19339 [483-04-5] 🌡	**Ajmalicine (Raubasine) extrapure, 98%**	10 Mg 50 Mg 1 Gms	1500 6000 11200
52821 [4360-12-7] 🌡	**Ajmaline pure, 98%** (Ritmos)	50 Mg 1 Gms	8000 15800
74590 [33208-99-0]	**L-Alaninamide Hydrochloride extrapure, 97%** (L-Alanineamide Hydrochloride)	5 Gms 25 Gms	2400 9300
64422 [338-69-2]	**D-Alanine extrapure CHR, 99%**	25 Gms 100 Gms 250 Gms	1750 6000 12000
96244 [302-72-7]	**DL-Alanine extrapure CHR, 99%**	100 Gms 500 Gms 1 Kg	610 1650 3000
32344 [77471-41-1] ■	**D-Alanine-7-Amido-4-Methylcoumar in Free Base extrapure, 97%** (D-Alanine-7-Amido-4-Methylcoumarin)	10 Mg 25 Mg	11420 20760
11180 [201847-54-3] ■	**b-Alanine-7-Amido-4-Methylcoumar in Trifluoroacetate Salt extrapure, 99%** (ß-Alanine-7-Amido-4-Methylcoumarin Trifluoroacetate Salt, ß-Alanine-AMC.TFA salt)	10 Mg 25 Mg	12340 25220
28760 [201847-52-1] ■	**D-Alanine-7-Amido-4-Methylcoumar in Trifluoroacetate Salt extrapure, 99%** (D-Alanine-7-Amido-4-Methylcoumarin Trifluoroacetate Salt, D-Alanine-AMC.TFA salt)	10 Mg 25 Mg	8880 17260
25696 [96594-10-4] ■	**L-Alanine-7-Amido-4-Methylcoumar in Trifluoroacetate Salt extrapure, 98%** (L-Alanine-7-Amido-4-Methylcoumarin Trifluoroacetate Salt, L-Alanine-AMC.TFA salt)	25 Mg 100 Mg	5910 16870
94518 [42854-62-6] ■	**L-Alanine Benzyl Ester p-Toluene Sulfonate extrapure, 98%**	25 Gms	6230
17737 [5557-83-5] ■	**L-Alanine Benzyl Ester Hydrochloride extrapure, 98%**	25 Gms	5410

Code	Name	Packing	Price
64048 [3513-56-2]	■ **L-Alanine-ß-Naphthylamide Hydrobromide extrapure, 98%**	100 Mg	6360
94349 [1948-31-8]	🌡 **L-Alanyl-L-Alanine extrapure, 99%**	100 Mg 500 Mg	4030 19340
31622 [39537-23-0]	🌡 **L-Alanyl-L-Glutamine extrapure, 99%**	5 Gms 25 Gms	1200 3000
46620 [687-69-4]	🌡 **L-Alanyl Glycine extrapure, 99%**	1 Gms 5 Gms	6000 24000
73721 [3303-34-2]	🌡 **L-Alanyl-L-Leucine extrapure, 99%**	500 Mg	9320
31550	● **L-Alanyl-L-Leucine-L-Phenylalanine extrapure**	100 Mg	8000
97088 [2491-20-5]	■ **L-Alanyl Methyl Ester Hydrochloride extrapure, 99%**	5 Gms 25 Gms	1790 6830
41859 [3061-88-9]	■ **L-Alanyl-L-Tyrosine extrapure, 99%**	100 Mg	3180
38211 [1596-84-5]	🌡 **Alar extrapure AR, 99%**	5 Gms 25 Gms	1200 4500
97350 [9048-46-8]	● **Albumin Bovine (BSA) (pH 5.2) fraction V, 98%** (Bovine Serum Albumin)	5 Gms 50 Gms 100 Gms 250 Gms	2000 12100 23700 52500
83803 [9048-46-8]	● **Albumin Bovine (BSA) (pH 7) fraction V, 98%** (Bovine Serum Albumin)	5 Gms 25 Gms 100 Gms 250 Gms 1 Kg	900 4000 15300 35700 120000
85171 [9048-46-8]	● **Albumin Bovine (BSA) (pH 7) fraction V for molecular biology, 98%** (Bovine Serum Albumin)	5 Gms 25 Gms 100 Gms 500 Gms 1 Kg	1000 4200 16600 75600 140000
62560 [9048-46-8]	● **Albumin Bovine (BSA) (pH 7) fraction V PrimeDx, 98%** (Bovine Serum Albumin)	1 Kg 5 Kg	190000 850000
14438 [9048-46-8]	● **Albumin Bovine (BSA) (pH 7) fraction V Lipid Enhanced, 98%** (Bovine Serum Albumin)	25 Gms 100 Gms	8400 25500
32257 [9048-46-8]	● **Albumin Bovine (BSA) (pH 7) fraction V for tissue culture, 98%** (Bovine Serum Albumin)	5 Gms 25 Gms	2500 9300
65186 [9048-46-8]	● **Albumin Bovine (BSA) (pH 7) fraction V Reagent Grade, 99%** (Bovine Serum Albumin)	10 Gms 25 Gms 100 Gms	5400 9800 21200
93630 [9048-46-8]	● **Albumin Bovine (BSA) (pH 7) fraction V Cohn Analog, 96%** (Bovine Serum Albumin)	10 Gms 25 Gms 100 Gms	5400 9800 21200
19007 [9048-46-8]	🌡 **Albumin Bovine Solution (20mg/ml)**	1 ml 5 x1 ml	1300 2600
77761 [9048-46-8]	🌡 **Albumin Bovine Solution 10% Diluent Solution in PBS (BSA 10% Diluent/Blocking Solution in PBS)**	10 ml	3500
45017 [9006-59-1]	● **Albumin Egg Powder extrapure, 78%**	100 Gms 500 Gms	750 3500
56587 [9006-59-1]	● **Albumin Egg Flakes extrapure, 85%**	100 Gms 500 Gms	1700 7760
21840 [70024-90-7]	■ **Albumin Human Serum (rHSA) (Lipid Enhanced), 95%**	100 Mg 1 Gms 5 Gms	2000 5500 14000
83428 [70024-90-7]	■ **r-Albumin Human (rHSA) - Type LE ex. Yeast, 99%** (Human Serum Albumin)	100 Mg 1 Gms	3000 24000
48261 [33864-99-2]	**Alcian Blue** (C.I. No. 74240, Alcian Blue 8 GX)	5 Gms 10 Gms 100 Gms	1800 2800 16000
46822 [33864-99-2]	**Alcian Blue for tissue culture** (C.I. No. 74240, Alcian Blue 8 GX)	5 Gms 25 Gms 100 Gms	2200 5200 17000
90008 [9031-72-5]	■ **Alcohol Dehydrogenase (ALD) ex. Saccharomyces Cerevisiae, 300U/mg protein** (ex. Bakers Yeast)	7500 Units 15 K.Units	5400 10200
67482 [1715-30-6]	■ **Alexidine Dihydrochloride extrapure, 96%** (1,1-Hexamethylene-bis(5-[2-Ethylhexyl]biguanide))	250 Mg	5510
86220 [9005-32-7]	**Alginic Acid ex. Brown Algae, 19-25%**	100 Gms 250 Gms 500 Gms	400 940 1780
84554 [109-57-9]	**N-Allylthiourea extrapure, 98%** (1-Allyl-2-thiourea, Thiosinamine)	25 Gms 100 Gms	1350 4800
29186 [130-22-3]	**Alizarin Red S** (C.I. No. 58005)	25 Gms 100 Gms 500 Gms	750 2700 13000
60352 [130-22-3]	**Alizarin Red S for tissue culture** (C.I. No. 58005)	25 Gms 100 Gms	950 3200
95786 [584-42-9]	**Alizarin Yellow GG, 90%** (Metachrome Yellow, Mordant Yellow 1, Salicyl Yellow) (C.I. No. 14025)	5 Gms 25 Gms	660 1320
12151 [9001-78-9]	■ **Alkaline Phosphatase (ALP, AP, CIP) ex. Calf Intestine Mucosa, 20 DEA U/mg powder**	100 Mg 500 Mg 1 Gms	6490 12000 21000
97534	■ **Alkalophilic Proteinase ex. Streptomyces Sp., 20U/mg powder**	100 Units 200 Units	36050 67440
58020 [315-30-0]	**Allopurinol extrapure, 98%**	5 Gms 25 Gms	600 2200
46287 [7283-09-2]	■ **b-D-Allose extrapure, 98%**	100 Mg	7740
35401 [7635-11-2]	■ **L-Allose extrapure, 97%**	25 Mg 50 Mg	12670 21810
73583 [9046-28-0]	■ **L-Alpha Glycerophosphate Oxidase (GPO) ex. Microorganism, 15U/mg powder**	500 Units 1 K.Units 5 K.Units	1700 3000 7000

Code	Name	Packing	Price
61973 [1990-29-0]	**D-Altrose extrapure**	25 Mg	14710
		100 Mg	56280
		200 Mg	125470
74880 [26093-31-2]	**7-Amino-4-Methylcoumarin (7-AMC) extrapure, 99%**	1 Gms	3030
81396 [37517-28-5]	**Amikacin (AMK) free base, 98%**	50 Mg	6000
		250 Mg	25000
18036 [39831-55-5]	**Amikacin Sulphate (AMKS)**	1 Gms	2000
		5 Gms	6000
98062 [122-80-5]	**p-Aminoacetanilide pure, 98%**	100 Gms	510
		500 Gms	1820
40182 [99-92-3]	**p-Aminoacetophenone extrapure AR, 99%**	25 Gms	850
		100 Gms	3020
62360 [7620-28-2]	**D-a-Amino Adipic Acid extrapure, 98%**	250 Mg	2200
		1 Gms	6000
35418 [1118-90-7]	**L-a-Amino Adipic Acid extrapure, 98%**	250 Mg	2400
		1 Gms	7000
46832	**Amino Acids Kit**	1 Kit	3730
86936	**L-Amino Acids Kit**	1 Kit	4430
68384 [9000-88-8]	**D-Amino Acid Oxidase (DAA, DAO, DAMOX) ex. Porcine Kidney, 6000U/g powder**	50 Units	11260
		100 Units	19700
98437 [60-09-3]	**p-Aminoazobenzene pure, 98%**	100 Gms	300
		500 Gms	800
79793 [99-05-8]	**m-Aminobenzoic Acid pure, 98%**	100 Gms	700
		500 Gms	2600
74004 [99-05-8]	**m-Aminobenzoic Acid extrapure AR, 99%**	25 Gms	800
		100 Gms	2500
99839 [150-13-0]	**p-Aminobenzoic Acid (pABA) pure, 99%**	100 Gms	400
		250 Gms	800
		500 Gms	1500
		5 Kg	13000
19196 [150-13-0]	**p-Aminobenzoic Acid (pABA) ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	500
		500 Gms	1700
		1 Kg	3000
56960 [2835-77-0]	**2-Aminobenzophenone pure, 98%**	5 Gms	1800
		25 Gms	4640
26345 [56-12-2]	**4-Aminobutyric Acid extrapure CHR, 99%**	10 Gms	360
		25 Gms	890
65081 [2835-81-6]	**DL-a-Aminobutyric Acid extrapure, 99%** (DL-2-Aminobutyric Acid)	25 Gms	1960
		100 Gms	5850
10498 [719-59-5]	**2-Amino-5-Chlorobenzophenone pure, 98%**	100 Gms	4330
		500 Gms	21110
18125 [127946-77-4]	**1-Amino-1-Cyclopropanecarbonitrile Hydrochloride extrapure, 98%**	1 Gms	2500
		5 Gms	8000

Code	Name	Packing	Price
68678 [132-32-1]	**3-Amino-9-Ethyl Carbazole (AEC) extrapure, 95%**	1 Gms	800
		10 Gms	3500
		25 Gms	8500
44285 [1750-12-5]	**4-Amino-3-Hydrazino-5-Mercapto-1, 2,4-Triazole (AHMT), 98%** (equivalent to Purpald)	1 Gms	2200
		5 Gms	5100
		25 Gms	15000
94283 [16867-03-1]	**2-Amino-3-Hydroxypyridine pure, 98%**	25 Gms	1200
		100 Gms	4000
33864 [5451-09-2]	**5-Amino Levulinic Acid Hydrochloride extrapure, 97%**	100 Mg	1500
		1 Gms	3000
		5 Gms	9000
19902 [144-90-1]	**DL-3-Aminoisobutyric Acid extrapure, 98%** (2-Amino-2-Methyl-Propionic Acid)	1 Gms	4200
		5 Gms	12600
86898 [214139-20-5]	**DL-3-Aminoisobutyric Acid Hydrate extrapure, 98%**	25 Mg	1160
		100 Mg	2630
87907 [99-31-0]	**5-Aminoisophthalic Acid (5-AIPA) pure, 99%** (5-Aminobenzene-1,3-Dicarboxylic acid)	25 Gms	1300
		100 Gms	4000
36615 [23687-26-5]	**6-Aminoisoquinoline pure, 98%** (6-Isoquinolinylamine)	250 Mg	1900
		1 Gms	4500
50407 [124-68-5]	**2-Amino-2-Methyl-1-Propanol (AMP Buffer) pure, 95%**	500 ml	2200
		2500 ml	10600
51812 [115-69-5]	**2-Amino-2-Methyl-1,3-Propanediol (AMPD Buffer) extrapure, 99%** (Ammediol)	25 Gms	3800
		100 Gms	13500
91085 [84-89-9]	**5-Aminonaphthalene-1-Sulphonic Acid pure, 90%** (1-Naphthylamine-5-Sulfonic Acid, Laurents Acid)	25 Gms	2630
		100 Gms	9450
98372 [119-79-9]	**5-Aminonaphthalene-2-Sulphonic Acid pure, 98%**	25 Gms	2940
		100 Gms	9980
79911 [116-63-2]	**1-Amino-2-Naphthol-4-Sulphonic Acid (1,2,4-Acid) ACS, 90%** (1-Amino-2-Hydroxy-4-Naphthalenesulfonic acid)	25 Gms	520
		100 Gms	1740
		500 Gms	7200
21407 [1775-95-7]	**2-Amino-5-Nitrobenzophenone extrapure, 98%**	25 Gms	2940
		100 Gms	9670
58513 [2973-50-4]	**2-Aminophenyl-Acetonitrile pure, 97%**	250 Mg	5980
		1 Gms	15230
33993 [919-30-2]	**(3-Aminopropyl)triethoxysilane (APTES) extrapure, 98%** (3-(Triethoxysilyl)Propylamine)	100 ml	1200
		500 ml	3400
17081 [452-06-2]	**2-Aminopurine extrapure, 98%**	100 Mg	3220
		250 Mg	6910
		1 Gms	23300
65377 [504-29-0]	**2-Aminopyridine pure, 98%**	100 Gms	500
		250 Gms	1000
		500 Gms	1800
29786 [462-08-8]	**3-Aminopyridine pure, 98%**	25 Gms	900
		100 Gms	3000
		500 Gms	10500

Code	Name	Packing	Price
16104 [504-24-5]	**4-Aminopyridine pure, 98%**	10 Gms 25 Gms 100 Gms	470 850 2800
12920 [584-13-4]	**4-Amino-1,2,4-Triazole pure, 98%**	25 Gms 100 Gms	900 3000
40057 [54723-94-3]	**Ammonium Phosphomolybdate Hydrate extrapure, 61% Mo** (Ammonium-12-Molybdophosphate,Molybdophosphoric Acid Ammonium,Phosphomolybdic Acid Ammonium Hydrate)	25 Gms 100 Gms	1300 3200
31438 [13820-40-1]	**Ammonium Tetrachloropalladate (II) extrapure, 98%** (Palladium (II) Ammonium Chloride)	250 Mg 1 Gms	4000 10000
32285 [13820-41-2]	**Ammonium Tetrachloroplatinate (II) extrapure, 99%** (Platinum (II) Ammonium Chloride)	250 Mg 1 Gms	4200 13000
24645 [61336-70-7]	**Amoxicillin Trihydrate (AMOT), 97%**	1 Gms 5 Gms 10 Gms	300 1290 2330
54713 [1397-89-3]	**Amphotericin B (AMT)**	250 Mg 1 Gms 5 Gms	1290 4910 23950
61314 [69-52-3]	**Ampicillin Sodium Salt (AMP-Na)**	1 Gms 5 Gms 25 Gms	380 1730 6290
26647 [69-52-3]	**Ampicillin Sodium Salt (AMP-Na) for tissue culture**	1 Gms 5 Gms 25 Gms	460 1710 6720
82904 [7177-48-2]	**Ampicillin Trihydrate (AMP), 99%**	1 Gms 5 Gms 10 Gms	330 1420 2570
96459 [68399-79-1]	**AMPSO Buffer extrapure, 99%** (N-(1,1-Dimethyl-2-Hydroxyethyl) -3-Amino-2-Hydroxypropane Sulphonic Acid)	5 Gms 25 Gms 100 Gms	3000 8000 17000
85441 [29883-15-6]	**D-Amygdalin extrapure, 97%**	1 Gms 10 Gms 25 Gms	1440 4410 8000
95114 [9000-90-2]	**a-Amylase ex. Porcine Pancreas, 5U/mg powder**	1 Gms 5 Gms	1200 5000
28588 [9000-90-2]	**a-Amylase ex. Porcine Pancreas, 10-25U/mg powder**	100 Mg	4200
46504 [9000-92-4]	**a-Amylase ex. Malt (Fungal Diastase ex. Aspergillus Oryzae) 1:2000, (2000U/g powder)**	100 Gms 500 Gms	700 2600
40321 [206659-00-9]	**8-Anilino-1-Naphthalene Sulphonic Acid Ammonium Salt Hydrate (ANS Hydrate) extrapure, 98%** (8-Anilino-1-Naphthalenesulfonic Acid Ammonium Salt Hydrate)	5 Gms 25 Gms 100 Gms	2900 8300 31900
42275 [82-76-8]	**8-Anilino-1-Naphthalenesulphonic Acid (ANS) extrapure, 98%** (8-Anilino-1-Naphthalenesulfonic Acid, N-Phenyl Peri Acid)	5 Gms 25 Gms 100 Gms	3100 8900 19900
99610 [18108-68-4]	**8-Anilino-1-Napthalenesulfonic Acid Hemimagnesium Salt Hydrate extrapure, 98%** (ANSA Magnesium Salt, Magnesium(II) 8-Anilino-1-Naphthalenesulfonate, 1,8-ANS Magnesium Salt, N-Phen	5 Gms 25 Gms	3500 12000
89119 [123-11-5]	**Anisaldehyde extrapure, 99%** (p-Methoxybenzaldehyde)	100 ml 250 ml 500 ml	600 1100 2200
90820 [100-09-4]	**p-Anisic Acid extrapure, 98%** (p-Methoxybenzoic Acid)	100 Gms 500 Gms	600 2200
70940 [100-66-3]	**Anisole extrapure, 99%**	500 ml	800
64447 [100-66-3]	**Anisole Reagent Grade, 99%**	500 ml	1300
59553 [90-44-8]	**Anthrone pure, 95%**	25 Gms 100 Gms	500 1600
96476 [90-44-8]	**Anthrone ACS, 97%**	25 Gms 100 Gms	700 2000
63714 [90-44-8]	**Anthrone ExiPlus, Multi-Compendial, 97%** meets compendial specs of BP	25 Gms 100 Gms	800 2400
44089 [85251-63-4]	**Apple Seed (PyrusMalus) Oil**	250 ml 500 ml	2200 4200
62179 [9087-70-1]	**Aprotinin ex. Bovine Lung, 3-8 TIU/mg powder** (BPTI, Trypsin Inhibitor)	5 Mg 25 Mg 100 Mg 250 Mg	2000 7100 18900 44100
79815 [9036-66-2]	**Arabinogalactan ex. Larch Wood, 95%** ((+)-Arabinogalactan, Polyarabinogalactan)	25 Gms 100 Gms	5200 15000
23848 [10323-20-3]	**D-Arabinose extrapure** (D-(-)-Arabinose)	5 Gms 25 Gms 100 Gms	400 1400 5300
52392 [87-72-9, 5328-37-0]	**L-Arabinose extrapure** (L-(+)-Arabinose)	100 Gms 500 Gms	1000 4700
52384 [87-72-9, 5328-37-0]	**L-Arabinose ExiPlus, Multi-Compendial** (L-(+)-Arabinose) meets compendial specs of BP, Ph.Eur	25 Gms 500 Gms	900 6300
74926 [488-82-4]	**D(+)Arabitol (D-Arabinitol) extrapure 99%**	1 Gms 5 Gms 25 Gms	800 3200 10900
91121 [7643-75-6]	**L-(-) Arabitol (L-Arabinitol) extrapure 99%**	1 Gms 10 Gms	1300 5500
89657 [506-30-9]	**Arachidic Acid extrapure, 99%** (Eicosanoic Acid, Arachic Acid, Arachidinic Acid)	1 Gms 5 Gms 25 Gms	2700 6300 16800
20975 [506-32-1]	**Arachidonic Acid (non-animal) extrapure, 99%**	10 Mg 25 Mg 100 Mg	2700 4800 11600

Code	Name	Packing	Price
13492 [497-76-7]	**p-Arbutin extrapure, 98%** (4-Hydroxyphenyl-B-D-glucopyranoside, Hydroquinone-B-D-glucopyranoside)	10 Gms 25 Gms 100 Gms	3200 6300 12600
89184 [157-06-2]	**D-Arginine (free base) extrapure CHR, 99%** ((R)-2-Amino-5-Guanidinopentanoic Acid)	1 Gms 5 Gms	2500 5300
34537 [36589-29-4]	**L-Arginine Ethyl Ester Dihydrochloride extrapure, 98%**	5 Gms 25 Gms 100 Gms	1400 5200 16900
11249 [26340-89-6]	**L-Arginine Methyl Ester Dihydrochloride extrapure, 98%**	25 Gms 100 Gms	2900 8400
47014 [1668-00-4]	**Arsenazo III Free Acid (High Purity) extrapure AR, 85%** (2,7-Bis(2-Arsonophenylazo) -1,8-Dihydroxynaphthaline-3,6-Disulfonic Acid)	1 Gms 5 Gms	7900 36800
57656 [9029-44-1]	■ **L-Ascorbate Oxidase (Type 1) ex. Cucurbita Sp., 150U/mg powder**	2500 Units 10000 Units	5800 16300
60402 [9029-44-1]	■ **L-Ascorbate Oxidase (Type 2) ex. Acremonium Sp., 350U/mg powder**	1000 Units	7500
27999 [464-92-6]	● **Asiatic Acid pure, 97%** (2a,23-Dihydroxyursolic Acid, Dammarolic acid, NSC 166063)	50 Mg 200 Mg	4300 15100
11479 [5794-24-1]	**D-Asparagine Monohydrate extrapure, 99%**	5 Gms 25 Gms	800 2400
86883 [5794-13-8]	**L-Asparagine Monohydrate extrapure CHR, 99%**	100 Gms 1 Kg 5 Kg	800 6100 28400
92468 [5794-13-8]	**L-Asparagine Monohydrate for tissue culture, 99%**	25 Gms 100 Gms 500 Gms	600 1200 4800
44444 [5794-13-8]	**L-Asparagine Monohydrate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 1 Kg	1000 6700
90601 [1783-96-6]	**D-Aspartic Acid extrapure CHR, 99%**	5 Gms 25 Gms	900 2700
71353 [16115-68-7]	**L-Aspartic Acid Diethyl Ester Hydrochloride extrapure, 98%**	5 Gms 25 Gms	2100 6300
50144 [32213-95-9]	**L-Aspartic Acid Dimethyl Ester Hydrochloride extrapure, 98%**	1 Gms 5 Gms 25 Gms 100 Gms	1200 4200 12800 43200
57358 [5908-99-6]	△ **Atropine Sulfate Monohydrate extrapure, 98%** (a-(Hydroxymethyl)Benzeneacetic Acid 8-Methyl-8-Azabicyclo[3.2.1]Oct-3-yl Ester, Tropine Tropate)	5 Gms 25 Gms	1600 7200
21505 [1405-69-2]	■ **Avidin ex. Egg White for molecular biology, 10U/mg Protein**	10 Mg 25 Mg	5000 10000
34476 [320-67-2]	■ **5-Azacytidine extrapure, 98%**	100 Mg 1 Gms	1900 3200
59070 [76801-85-9]	**Azaerythromycin (AZAE)**	1 Gms 5 Gms	2400 5700
81947 [932-53-6]	**6-Azathymine extrapure, 98%**	1 Gms	6200
61898 [30516-87-1]	● **3-Azido-3-Deoxythymidine extrapure, 99%** (AZT, Azidothymidine, ZDV, Zidovudine)	25 Mg 1 Gms 5 Gms	900 3200 8000
17156 [117772-70-0]	**Azithromycin Dihydrate (AZM), 96%**	1 Gms 5 Gms 10 Gms	800 3500 6300
47138 [2997-92-4]	**2,2-Azobis(2-Methylpropionamidine) Dihydrochloride (AAPH) extrapure, 98%** a,a'-Azodiisobutyramidine dihydrochloride	5 Gms 25 Gms 100 Gms	2700 4500 7200
11029 [78-67-1]	**2,2-Azobisisobutyronitrile (AIBN) pure, 98%** (2,2-Azobis(2-Methylpropionitrile))	100 Gms 500 Gms 5 Kg	600 2100 20000
95426 [102110-73-6]	■ **Azoalbumin (Bovine) prepared from Bovine fraction V powder extrapure**	250 Mg 500 Mg 1 Gms	3600 6900 13100
21070 [102110-74-7]	● **Azocasein extrapure**	250 Mg 500 Mg 1 Gms	3360 6160 11540
37020 [206752-32-1]	**Azomethine H Monosodium Salt Hydrate, 97%**	1 Gms 5 Gms 25 Gms	790 3480 12320
78130 [5941-07-1]	**Azomethine H Monosodium Salt, 97%**	1 Gms 5 Gms 25 Gms	1240 5040 13440
13741 [78110-38-0]	**Aztreonam (AZN), 96%**	500 Mg 1 Gms	3590 6720
53953 [1192-88-7]	**1-Cyclohexene-1-Carboxaldehyde extrapure, 97%**	1 Gms 5 Gms	8400 40320
76432 [21668-81-5]	**3-(Carbamimidoylthio) -1-Propanesulfonic Acid extrapure, 98%**	5 Gms 25 Gms 100 Gms	1200 2500 5000
11730 [531-55-5]	**Azure B (Azure I)** (C.I. No. 52010) (N,N,N-Trimethylthionin, Methylene Azure A)	25 Gms 100 Gms 500 Gms	700 2600 11600
62151 [37247-10-2]	**Azur II** (C.I. 52010/52015) (3-Chloro-6-Hydroxy-2-Pyridinecarboxylic Acid, Methylene Azure II)	10 Gms 25 Gms 100 Gms	400 700 2400
25003 [531-53-3]	**Azure A Chloride** (C.I. No. 52005) (3-Amino-7-(Dimethylamino)phenothiazin-5-ium Chloride)	25 Gms 100 Gms 500 Gms	700 2600 11600
17327 [1405-87-4]	● **Bacitracin (BCT) ex. Bacillus Licheniformis**	50000 Units 250000 Units	1600 3900

Code	Name	Packing	Price
33496 ■ [126150-97-8]	**BAPTA AM extrapure, 98%** Tetrakis (Acetoxymethyl) 1,2-Bis(2-Aminophenoxy)Ethane-N,N,N',N'-Tetraacetate, 1,2-Bis(2-Aminophenox	10 Mg 25 Mg	5500 12000
34501 🌡 [85233-19-8]	**BAPTA Free Acid extrapure, 98%** (1,2-Bis(2-Aminophenoxy) Ethane-N,N,N',N'-Tetraacetic Acid, 2,2'-(Ethylenedioxy)Dianiline-N,N,N',N'-	100 Mg 1 Gms 5 Gms	1900 4000 18000
52571	**Barfoed`s Reagent**	250 ml	700
29090 [6211-24-1]	**Barium Diphenylamine-4-Sulfonate ACS** (4-(Phenylamino)benzenesulfonic Acid Barium Salt, Diphenylamine-4-Sulfonic Acid Barium Salt)	5 Gms 25 Gms 100 Gms	310 960 3400
27297 [4733-39-5]	**Bathocuproine (BCP) ACS, 99%**	500 Mg 1 Gms 5 Gms	2500 4000 18000
28242 [52698-84-7]	**Bathocuproinedisulphonate Disodium Salt extrapure AR, 98%**	100 Mg 250 Mg 1 Gms 5 Gms	1800 4200 12000 50000
50221 [1662-01-7]	**Bathophenanthroline extrapure AR, 99%** (4,7-Diphenyl-1,10-Phenanthroline, BPhen)	1 Gms 5 Gms	2800 12000
98680 [52746-49-3, 98645-86-4, 53744-42-6]	**Bathophenanthroline Disulphonate Disodium Salt Hydrate extrapure, 99%**	250 Mg 1 Gms 5 Gms	3700 12100 52500
36333 [21205-91-4]	**9-BBN Triflate Solution (0.5M in hexanes) pure** (9-[{9-(Trifluoromethyl)sulfonyl} oxy]9-Borabicyclo(3.3.1)Nonane)	500 ml 1000 ml	8300 13400
64733 ✔ [112-85-6]	**Behenic Acid pure, 85%** (Docosanoic acid, 1-Docosanoic acid, n-Docosanoic acid, Docosoic acid)	25 Gms 250 Gms	500 900
35591 ✔ [112-85-6]	**Behenic Acid extrapure, 99%** (Docosanoic acid, 1-Docosanoic acid, n-Docosanoic acid, Docosoic acid)	5 Gms 25 Gms	4200 7400
74710 🌡 [661-19-8]	**1-Docosanol (Behenyl Alcohol) extrapure, 98%**	25 Gms 500 Gms	800 5800
40270 🌡 [2798-25-6]	**Bellidifolin (BLF) extrapure, 99%** (2,4-Dichlorobenzyl Alcohol)	10 Mg	14000
82142 [94-41-7]	**Benzalacetophenone (Chalcone) pure, 98%** (1,3-Diphenyl-2-Propenone)	25 Gms 100 Gms 500 Gms	810 2190 8510
93014 🌡 [206752-36-5]	**Benzamidine Hydrochloride Hydrate extrapure, 98%**	5 Gms 25 Gms 100 Gms	1150 2990 9090
13740 [93-98-1]	**Benzanilide pure, 98%**	100 Gms 500 Gms	810 2880
42689 [98-09-9]	**Benzenesulphonyl Chloride pure, 98%**	500 ml	1270
41035 [98-09-9]	**Benzenesulphonyl Chloride ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP	500 ml	1500
57730 [91-00-9]	**Benzhydrylamine pure, 97%** (a-Aminodiphenylmethane)	25 ml 100 ml 500 ml	810 1500 6100
17851 🌡 [495-18-1]	**Benzohydroxamic Acid extrapure, 99%**	5 Gms 25 Gms	3220 6850
17192 [134-81-6]	**Benzil extrapure, 98%**	250 Gms	1040
93528 [76-93-7]	**Benzilic Acid extrapure, 99%** (Diphenylglycolic Acid)	100 Gms 250 Gms	1610 3800
35034 [51-17-2]	**Benzimidazole extrapure, 99%**	25 Gms 100 Gms 500 Gms	350 920 4260
31398 [441-38-3]	**a-Benzoinoxime extrapure AR, 99%** (Cupron)	25 Gms 100 Gms	580 1610
44371 🌡 [5350-57-2]	**Benzophenone Hydrazone pure, 99%**	100 Gms 500 Gms	3680 14610
47554 [93-91-4]	**1-Benzoylacetone pure, 98%**	10 Gms 50 Gms	2760 10240
76356 ● [2645-08-1]	**N-a-Benzoyl-L-Arginine Ethyl Ester Hydrochloride (BAEE) extrapure AR, 99%**	1 Gms 5 Gms 25 Gms	1000 4700 22100
96249 ● [911-77-3]	**N-a-Benzoyl-DL-Arginine-4-Nitroanilide Hydrochloride (BANI, L-BAPNA) extrapure, 98%**	100 Mg 250 Mg 1 Gms	1900 3200 9300
51003 ● [913-04-2]	**N-a-Benzoyl-DL-Arginine-ß-Naphthylamide Hydrochloride (BANA) extrapure, 97%**	250 Mg 1 Gms	3100 9700
46170 ● [23405-15-4]	**N-(Benzoyloxy)Succinimide (Bz-OSu) extrapure, 99%** (Benzoic Acid N-Succinimidyl Ester, N-Succinimidyl Benzoate)	1 Gms 5 Gms	4200 8700
46651 ● [1214-39-7]	**6-Benzyladenine (6-BAP) extrapure AR, 99%** (6-Benzylaminopurine)	5 Gms 25 Gms 100 Gms 500 Gms	640 940 3330 15000
54395 ● [1214-39-7]	**6-Benzyladenine (6-BAP) for tissue culture, 99%** (6-BAP, BA, N6-Benzyladenine, 6-Benzylaminopurine, N-benzyl-7H-purin-6-amine, 6-Benzyl adenine)	5 Gms 25 Gms 100 Gms 500 Gms	1050 5100 8800 20000
38728 [100-46-9]	**Benzylamine pure, 99%**	500 ml 2500 ml 25 Ltr	920 4140 38640
64554 [100-46-9]	**Benzylamine Reagent Grade, 99%**	250 ml 500 ml	690 1150
91695 [146631-00-7]	**4-Benzyloxy-Phenylboronic Acid extrapure, 97%**	5 Gms 25 Gms	4500 8500
47102 [3054-01-0]	**S-Benzyl-L-(+)-Cysteine extrapure, 98%**	25 Gms	18400

Code	Name	Packing	Price
26862 [472696-9]	**O-Benzyl-L-Serine extrapure, 99%**	1 Gms	1800
		5 Gms	6400
		10 Gms	12300
		25 Gms	28000
50562 [2899-66-3]	**3-(Benzylthio)-Propionic Acid extrapure, 98%**	1 Gms	3800
		5 Gms	10000
75036 [21871-47-6]	**5-(Benzylthio)-1H-Tetrazole (BTT Activator) extrapure, 99%**	5 Gms	1300
		25 Gms	4000
29161 [56-37-1]	**Benzyltriethylammonium Chloride pure, 98%**	500 Gms	690
		5 Kg	6210
98924 [56-39-9]	**Benzyltrimethylammonium Chloride pure, 98%**	100 Gms	460
		500 Gms	2070
		5 Kg	18400
60168 [56-39-9]	**Benzyltrimethylammonium Chloride ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	580
		500 Gms	2420
30307 [100-85-6]	**Benzyltrimethyl Ammonium Hydroxide 40% in Methanol extrapure** (Triton B)	100 ml	1610
		500 ml	4950
18437 ✔ [8007-75-8]	**Bergamot Oil extrapure**	500 ml	6000
		1000 ml	10700
34489 [10191-18-1]	**BES Buffer extrapure, 99%** (N,N-bis(2-Hydroxyethyl)2-Aminoethan esulphonic Acid)	25 Gms	900
		100 Gms	3000
		500 Gms	12000
		1 Kg	22000
59299 [66992-27-6]	**BES Sodium Salt Buffer extrapure, 99%**	25 Gms	3000
		100 Gms	7000
		250 Gms	14000
83710	**Bial's Reagent**	250 ml	700
84867 🌡 [979-88-4]	**Bicinchoninic Acid Disodium Salt (BCA) extrapure, 98%** (2,2-Biquinoline-4,4-Dicarboxylic Acid Disodium Salt)	1 Gms	3000
		5 Gms	8900
39849 [150-25-4]	**BICINE Buffer extrapure, 99%** (N,N-bis(2-Hydroxyethyl) Glycine)	25 Gms	630
		100 Gms	1900
		500 Gms	8500
66466 ● [86303-22-2]	**Big CHAPS, 95%** N,N'-Bis(3-D-gluconamidopropyl) cholamide	500 Mg	10400
		1 Gms	18200
		5 Gms	58300
64839 🌡 [635-65-4]	**Bilirubin extrapure AR, 99%**	100 Mg	1000
		1 Gms	6800
		5 Gms	27300
66915 ■ [80619-01-8]	**Bilirubin Oxidase ex. Myrothecium, 1.2U/mg powder**	25 Units	5800
18888 🌡 [58-85-5]	**D-Biotin extrapure, 99%** (Vitamin H)	1 Gms	600
		10 Gms	4600
		25 Gms	8500
		100 Gms	28000
89932 🌡 [58-85-5]	**D-Biotin for tissue culture, 99%** (Vitamin H)	1 Gms	660
		10 Gms	5500
		25 Gms	10000
		100 Gms	34000
92196 [92-52-4]	**Biphenyl extrapure, 99%**	100 Gms	470
		1 Kg	1500
		5 Kg	6000
79363 [366-18-7]	**2,2-Bipyridyl ACS, 99.5%**	5 Gms	800
		25 Gms	3000
		100 Gms	8300
37025 [14592-56-4]	**Bis(Acetonitrile)Dichloropalladium (II) extrapure, 99%, ~41% Pd** (Palladium (II) Chloride Diacetonitrile Complex)	500 Mg	6000
		1 Gms	11000
		5 Gms	53000
98369 [10526-07-5]	**1,3-Bis(3-Aminophenoxy) Benzene (APB) pure, 98%** (Resorcinol bis(3-Aminophenyl) Ether)	5 Gms	6790
		25 Gms	25190
47063 [13080-89-2]	**4,4-Bis(4-Aminophenoxy)Diphenyl Sulfone pure, 98%**	5 Gms	3600
		25 Gms	13500
11866 ✔ [14220-64-5]	**Bis(benzonitrile)palladium(II) Dichloride extrapure, 27.4% Pd**	500 Mg	4200
		1 Gms	7900
		5 Gms	37000
24637 🌡 [542-88-1]	**Bis(chloromethyl) Ether pure, 96%**	25 Gms	5800
		100 Gms	16800
96246 [32005-36-0]	**Bis(Dibenzylideneacetone) Palladium extrapure, 16.5%-20.5% Pd** (Palladium(0) bis(Dibenzylideneacetone), Pd(dba)2)	250 Mg	2500
		1 Gms	6000
		5 Gms	25000
52119 🌡 [140-95-4]	**N,N-Bis(hydroxymethyl)urea (DMU) pure, 95%** (1,3-Bis(hydroxymethyl)urea, Dimethylolurea)	25 Gms	1400
		100 Gms	3500
		500 Gms	15800
98932 [20734-58-1]	**1,8-Bis(dimethylamino) Naphthalene (DMAN) extrapure, 99%** (N,N,N',N'-Tetramethyl-1,8-Diaminonaphthalene)	5 Gms	1800
		25 Gms	5500
16660 [72287-26-4]	**1,1`-Bis(Diphenylphosphino)Ferrocene Palladium(II) Dichloride extrapure, 14.5% Pd** (Dichloro [(1,1-Bis-Diphenylphosphino Ferrocene)] Palladium (II) (14.5% Pd))	250 Mg	3000
		1 Gms	5800
		5 Gms	24000
87281 [95464-05-04]	**1,1`-Bis(Diphenylphosphino)Ferrocene-Palladium(II)Dichloride-Dichloromethane Complex extrapure, 98%**	1 Gms	5000
		5 Gms	18000
		25 Gms	80000
49463 ■ [29934-17-6]	**Bis(Tricyclohexylphosphine)Palladium(II) Dichloride extrapure, 98%** (Dichlorobis(tricyclohexylphosphine)Palladium(II))	250 Mg	3000
		1 Gms	5800
		5 Gms	24000
35994 🌡 [13965-03-2]	**Bis(Triphenylphosphine) Palladium(II) Dichloride extrapure, 14%-15.50% Pd**	1 Gms	5000
		5 Gms	22000
		25 Gms	100000
26477 ✔ [10416-59-8]	**N,O-Bis-(Trimethylsilyl) Acetamide (BSA) extrapure, 95%**	25 ml	900
		100 ml	2500
		250 ml	5300

Code	Name	Packing	Price
57677 ✔ [25561-30-2]	**N,O-Bis-(Trimethylsilyl) Trifluoro Acetamide (BSTFA) extrapure, 99%**	5 ml 25 ml 100 ml 500 ml	500 1500 5100 23100
52353 🌡	**BSA+TMCS(3-5%), 93-97%** (N,O-Bis(trimethylsilyl) Acetamide, Trimethylchlorosilane)	10 ml 25 ml	3100 6200
22555 ✔	**BSA+TMCS+TMSI, (3:2:3)** N,O-Bis(trimethylsilyl)Acetamide, Trimethylchlorosilane,Trimethylsilylimidazole (GC Derivatization)	10 x1 ml 25 ml	5800 8900
46976 ● [25561-30-2]	**N,O-Bis(trimethylsilyl) trifluoro acetamide, Trimethylchlorosilane (BSTFA+TMCS) (99:1)** for GC Derivatization	10 x1 ml 25 ml 100 ml	5800 2600 9000
31817 [80-05-7]	**Bisphenol A pure, 97%** (4,4-Isopropylidenediphenol)	500 Gms	500
25417 [6976-37-0]	**BIS-TRIS Buffer extrapure, 99%** (bis [2-Hydroxyethyl] Aminotris [Hydroxymethyl] Methane (bis-Tris) (free base)	25 Gms 100 Gms 500 Gms	1310 4790 22360
15893 [124763-51-5]	**BIS-TRIS Hydrochloride Buffer extrapure, 98%**	10 Gms 25 Gms 100 Gms	3110 5820 14860
49084 [64431-96-5]	**BIS-TRIS Propane Buffer extrapure, 99%**	25 Gms 100 Gms 500 Gms	4000 13000 52500
77031	**Biuret Reagent Solution ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	250 ml	300
20899 🌡 [9041-93-4]	**Bleomycin Sulfate (BLM-S), 90%**	10 Mg 50 Mg	6300 23100
22649 ● [1871-22-3]	**Blue Tetrazolium Chloride (BTC, Tetrazolium Blue Chloride) extrapure AR**	1 Gms 5 Gms 25 Gms	1200 5500 25000
86204 ● [7764-95-6]	**BOC-D-Alanine extrapure, 99%**	10 Gms 25 Gms	3000 5000
42849 ● [15761-38-3]	**BOC-L-Alanine extrapure, 99%**	10 Gms 25 Gms	750 1780
50483 ● [28875-17-4]	**BOC-L-Alanine Methyl Ester extrapure, 98%**	10 Gms 25 Gms	8340 14590
75597 ● [4248-19-5]	**BOC Amide pure, 98%** (tert-Butyl Carbamate)	5 Gms 25 Gms 100 Gms	880 1980 6600
23875 ●	**BOC-L-Amino Acids Kit**	1 Kit	31620
50647 ● [73874-95-0]	**N-BOC-4-Aminopiperidine pure, 97%**	5 Gms 25 Gms	3600 13170
33534 ● [23680-31-1]	**N-BOC-O-Benzyl-L-Serine extrapure, 98%** (BOC-L-Ser(Bzl)-OH)	10 Gms 25 Gms	27350 47850
12320 ● [2130-96-3]	**N-BOC-O-Benzyl-L-Tyrosine extrapure, 98%** (BOC-L-Tyr(Bzl)-OH)	10 Gms 25 Gms	11490 20100
61900 ● [35897-34-8]	**BOC-L-Arginine Hydrochloride extrapure, 98.5%**	25 Gms 100 Gms	4830 13580
12202 ● [7536-58-5]	**BOC-L-Aspartic Acid-4-Benzylester extrapure, 98%**	10 Gms 25 Gms	2600 6500
19362 ● [75647-01-07]	**BOC-D-Aspargine extrapure, 99%**	10 Gms 25 Gms	4710 10980
67012 ● [7536-55-2]	**BOC-L-Aspargine extrapure, 99%**	10 Gms 25 Gms 100 Gms	1070 2380 9020
38230 ● [62396-48-9]	**BOC-D-Aspartic Acid extrapure, 99%**	10 Gms 25 Gms	3910 8800
30394 ● [13726-67-5]	**BOC-L-Aspartic Acid extrapure, 98%**	10 Gms	3000
35516 ● [19746-37-3]	**BOC-L-Acetamidomethyl-L-Cysteine extrapure, 98%**	10 Gms 25 Gms	17400 30440
44054 ● [5068-28-0]	**BOC-L-Benzyl-L-Cysteine extrapure, 98%**	10 Gms 25 Gms	11960 20930
53340 ● [132844-48-5]	**BOC-Cyclopropylamine extrapure, 98%**	10 Gms 25 Gms	10510 18390
14930 ● [10389-65-8]	**BOC-L-Cystine extrapure, 99%**	10 Gms 25 Gms	1060 2350
92584 ● [34404-28-9]	**BOC-D-Glutamic Acid extrapure, 98%**	10 Gms 25 Gms	2500 5000
69376 ● [2419-94-5]	**BOC-L-Glutamic acid extrapure, 99%**	10 Gms 25 Gms 100 Gms	1070 2380 9020
82284 ● [13574-13-5]	**BOC-L-Glutamic Acid-5-Benzyl Ester extrapure, 98%** (BOC-l-Glu(OBzl)-OH)	10 Gms 25 Gms	12330 21560
45737 ● [61348-28-5]	**BOC-D-Glutamine extrapure, 99%**	10 Gms 25 Gms	6500 14620
46975 ● [13726-85-7]	**BOC-L-Glutamine extrapure, 99%**	25 Gms 100 Gms	2350 8820
71943 ● [4530-20-5]	**BOC-Glycine extrapure, 99%**	25 Gms 100 Gms	2080 6260
44986 ● [31954-27-5]	**BOC-Glycine Methyl Ester extrapure, 98%** (BOC-Gly-OMe)	10 ml 25 ml	1900 4500
10733 ● [17791-52-5]	**BOC-L-Histidine extrapure, 98.5%**	10 Gms 25 Gms	4450 8610
38544 ● [13726-69-7]	**BOC-L-Hydroxyproline extrapure, 98%** (BOC-Hyp-OH)	5 Gms 25 Gms	5240 11860
74956 ● [13139-16-7]	**BOC-L-Isoleucine extrapure, 99%**	10 Gms 25 Gms 100 Gms	1190 2670 9000
55665 ● [84358-13-4]	**BOC-Isonipecotic Acid extrapure, 99%** (1-BOC-Piperidine-4-Carboxylic Acid)	10 Gms 25 Gms	6390 11170

Code	Name	Packing	Price
39190	● **BOC-Isonipecotic Acid Methyl Ester extrapure, 98%** (N-Boc-Piperidine-4-Carboxylic Acid Methyl Ester)	10 Gms 25 Gms	10230 17890
22237 [16937-99-8]	● **BOC-D-Leucine extrapure, 98%**	10 Gms 25 Gms	3800 8000
25833 [13139-15-6]	● **BOC-L-Leucine Monohydrate extrapure, 99%**	25 Gms 100 Gms	1700 6000
88985 [2389-45-9]	● **Na-BOC-Ne-Z-L-Lysine (BOC-Lys(Z)-OH) extrapure, 98%** (BOC-Lys(Z)-OH)	5 Gms 25 Gms	4200 11980
57147 [15098-69-8]	● **N-a-e-T-BOC-L-Lysine-Dicyclohexyl ammonium Salt extrapure, 99%** (BOC-L-Lys(BOC)-OH.DCHA)	10 Gms 25 Gms	5820 10180
69681 [2418-95-3]	● **BOC-L-Lysine extrapure, 98%**	10 Gms 25 Gms	3500 7000
81350 [2488-15-5]	● **BOC-L-Methionine extrapure, 98.5%**	25 Gms 100 Gms	2400 5700
91418 [88495-54-9]	● **BOC-Nipecotic Acid extrapure, 99%**	10 Gms 25 Gms	6390 11170
82933 [1487631-41-1]	● **BOC-Nipecotic Acid Methyl Ester extrapure, 98%**	10 Gms 25 Gms	11240 19660
53717 [2188-18-3]	● **BOC-Nitro-L-Arginine extrapure, 99%** (Boc-Arg(NO2)-OH)	10 Gms 25 Gms	5630 10580
14910 [33305-77-0]	● **BOC-4-Nitro-L-Phenylalanine extrapure, 98%** (BOC-L-Phe(4 No2)-OH)	10 Gms 25 Gms	6560 11480
76600 [7535-56-0]	● **BOC-4-Nitro-L-Phenylalanine Ethyl Ester extrapure, 98%** (BOC-Phe-Onp)	10 Gms 25 Gms	11640 20360
40474 [18942-49-9]	● **BOC-D-Phenylalanine extrapure, 99%**	10 Gms 25 Gms	5270 11830
23293 [13734-34-4]	● **BOC-L-Phenylalanine extrapure, 99%**	10 Gms 25 Gms 100 Gms	1250 2870 9230
25694 [35909-92-3]	● **BOC-D-Phenylalanine Methyl Ester extrapure, 98%**	10 Gms 25 Gms	13050 22830
60090 [51987-73-6]	● **BOC-L-Phenylalanine Methyl Ester extrapure, 98%**	10 Gms 25 Gms	10150 17760
46437 [33125-05-2]	● **BOC-D-Phenylglycine extrapure, 99%**	1 Gms	1510
40931 [2900-27-8]	● **BOC-L-Phenylglycine extrapure, 99%**	10 Gms 25 Gms	5440 9520
53665 [26250-84-0]	● **BOC-Pipecolic Acid extrapure, 99%**	10 Gms 25 Gms	5880 10280
41311 [132910-79-3]	● **BOC-Pipecolic Acid Methyl Ester extrapure, 98%**	10 Gms 25 Gms	13120 22960
12098 [57260-71-6]	● **1-BOC-Piperazine extrapure, 99%**	25 Gms 100 Gms	1500 4500
43977 [35150-07-3]	● **BOC-L-Prolinamide extrapure, 98%**	1 Gms 5 Gms 10 Gms	4240 12690 18320
93481 [37784-17-1]	● **BOC-D-Proline extrapure, 99%**	10 Gms 25 Gms	3570 8030
91894 [15761-39-4]	● **BOC-L-Proline extrapure, 99%** (BOC-Pro-OH)	25 Gms 100 Gms 250 Gms	1800 5900 14000
71026 [733323-65-6]	● **BOC-D-Proline Methyl Ester extrapure, 98%**	10 Gms 25 Gms	16740 29300
31641 [59936-29-7]	● **BOC-L-Proline Methyl Ester extrapure, 98%**	10 Gms 25 Gms	15220 26630
91810 [6368-20-3]	● **BOC-D-Serine extrapure, 98%**	25 Gms 100 Gms	2500 4800
77874 [3262-72-4]	● **BOC-L-Serine extrapure, 99%**	10 Gms 25 Gms 100 Gms	2300 5160 16020
20256 [2766-43-0]	● **BOC-L-Serine Methyl Ester extrapure, 98%** (BOC-Ser-OMe)	10 ml 25 ml	4200 8000
73138 [55674-67-4]	● **BOC-D-Threonine extrapure, 98%**	10 Gms 25 Gms	4710 10980
77340 [2592-18-9]	● **BOC-L-Threonine extrapure, 99%**	10 Gms 25 Gms 100 Gms	1190 2670 9000
61017 [5241-64-5]	● **BOC-D-Tryptophan extrapure, 99%** (BOC-D-Trp-OH)	10 Gms 25 Gms	4000 9000
86659 [13139-14-5]	● **BOC-L-Tryptophan extrapure, 99%** (BOC-L-Trp-PH, BOC-Trp-OH)	25 Gms 100 Gms	2150 8000
61475 [3978-80-1]	● **BOC-L-Tyrosine extrapure, 99%**	10 Gms 25 Gms 100 Gms	2540 5340 16990
77170 [4326-36-7]	● **BOC-L-Tyrosine Methyl Ester extrapure, 98%**	10 Gms 25 Gms	10150 17760
35049 [22838-58-0]	● **BOC-D-Valine extrapure, 99%**	25 Gms	3800
80577 [13734-41-3]	● **BOC-L-Valine extrapure, 99%**	10 Gms 25 Gms 100 Gms	690 1380 4990
97810 [58561-04-9]	● **BOC-L-Valine Methyl Ester extrapure, 98%**	10 Gms 25 Gms	7330 12820
59946 [56602-33-6]	🌡 **BOP Reagent extrapure AR, 98%** (Castros Reagent)	5 Gms 25 Gms 100 Gms	750 3500 13500
24790 [2802-68-8]	🌡 **Boron Trifluoride in Methanol (50%) Complex**	50 ml 250 ml 1000 ml	4370 18500 40000
48847 [373-57-9]	🌡 **Boron Trifluoride in Methanol (14%) Solution**	25 ml 100 ml 500 ml	4850 6840 12850
68774 [9007-83-4]	■ **Bovine Immunoglobulin G (Bovine IgG) Reagent Grade, 90% IgG Total Protein**	1 Gms 5 Gms 25 Gms	6000 18000 60000

Code	Name	Packing	Price
19219	**Bradford Reagent for Proteins**	100 ml	600
		500 ml	1850
51195 [76962-43-7]	**Brassinolide (BR) synthetic technical grade, 90%**	25 Mg	3500
		100 Mg	9000
		500 Mg	40000
64222 [6104-58-1]	**Brilliant Blue G-250 for molecular biology** (Coomassie Brilliant Blue G250)	5 Gms	620
		25 Gms	2310
93473 [6104-59-2]	**Brilliant Blue R-250 for molecular biology** (Coomassie Brilliant Blue R250)	5 Gms	580
		25 Gms	2200
84778	**Brilliant Blue R Staining Solution (Coomassie Brilliant Blue R-250 solution)**	100 ml	3140
75824 [81029-05-2]	**Brilliant Cresyl Blue** (C.I. No. 51010 Brilliant Blue C, Brilliant Blue ALD)	5 Gms	800
		25 Gms	3000
		100 Gms	8000
60343 [9001-00-7]	**Bromelain ex. Pineapple Stem, 2400 GDU/g powder**	25 Gms	450
		100 Gms	1200
		500 Gms	3500
36953 [7252-83-7]	**Bromoacetaldehyde Dimethyl Acetal (Stabilized) pure, 97%** (2-Bromo-1,1-Dimethoxyethane)	25 ml	900
		100 ml	3000
		500 ml	12000
96141 [683-57-8]	**2-Bromoacetamide pure, 97%**	5 Gms	2300
		25 Gms	7500
25038 [103-88-8]	**p-Bromoacetanilide pure, 98%** (4-Bromoactanilide)	100 Gms	1600
		500 Gms	6800
39903 [79-08-3]	**Bromoacetic Acid pure, 99%**	100 Gms	800
		500 Gms	3400
15180 [99-90-1]	**4-Bromoacetophenone pure, 98%** (p-Bromoacetophenone)	25 Gms	600
		100 Gms	1500
		500 Gms	7000
30930 [2142-63-4]	**3-Bromoacetophenone extrapure, 99%** (1-Acetyl-3-Bromobenzene, m-Bromoacetophenone)	25 Gms	1100
		100 Gms	3200
47085 [591-19-5]	**m-Bromoaniline pure, 97%**	25 ml	2100
		100 ml	7300
22912 [106-40-1]	**p-Bromoaniline pure, 98%** (4-Bromoaniline)	100 Gms	1500
		500 Gms	6000
47405 [2398-37-0]	**3-Bromoanisole pure, 98%** (m-Bromoanisole)	50 ml	1100
		100 ml	2000
		500 ml	8700
54364 [104-92-7]	**4-Bromoanisole pure, 98%** (p-Bromoanisole)	100 Gms	1100
		500 Gms	3500
27879 [108-86-1]	**Bromobenzene extrapure, 99%**	250 ml	800
		1000 ml	3400
61714 [108-86-1]	**Bromobenzene Reagent Grade, 99%**	250 ml	1000
20107 [585-76-2]	**3-Bromobenzoic Acid extrapure, 99%**	10 Gms	1000
		25 Gms	2800
17833 [90-90-4]	**p-Bromobenzophenone pure, 98%**	5 Gms	1300
		25 Gms	4200
		100 Gms	13200
69530 [589-15-1]	**p-Bromobenzyl Bromide pure, 98%**	100 Gms	7300
79118 [536-38-9]	**4-Chlorophenacyl Bromide pure, 98%** (2-Bromo-4-Chloro Acetophenone, o-Bromo-4-Chloro Acetophenone)	25 Gms	1840
		100 Gms	4600
		500 Gms	17250
88685 [106-39-8]	**1-Bromo-4-Chloro Benzene pure, 98%**	100 Gms	1040
		500 Gms	3910
11447 [107-04-0]	**1-Bromo-2-Chloroethane pure, 97%**	100 Gms	2200
		500 Gms	9800
83788 [3252-36-6]	**5-Bromo-4-Chloro-3-Indolyl Acetate (X-3-Acetate, X-Acetate) extrapure, 98%**	25 Mg	3000
		250 Mg	12000
49199 [129541-42-0]	**5-Bromo-4-Chloro-3-Indolyl Caprylate (X-Caprylate) extrapure, 98%** (5-Bromo-4-chloro-3-Indolyl Octanoate)	25 Mg	2000
		250 Mg	8000
82860 [209347-94-4]	**5-Bromo-6-Chloro-3-Indolyl-Caprylate (Magenta Caprylate) extrapure, 97%**	25 Mg	1250
		100 Mg	2600
		500 Mg	9000
93775 [107021-38-5]	**5-Bromo-4-Chloro-3-Indolyl-a-D-Galactopyranoside (X-a-Gal) for molecular biology, 98%** (X-a-D-Galactoside)	25 Mg	1200
		100 Mg	3000
		1 Gms	12000
74021 [108789-36-2]	**5-Bromo-4-Chloro-3-Indolyl-a-D-Glucopyranoside (X-a-Glc) for molecular biology, 98%** (X-a-D-Glucoside)	50 Mg	3200
		1 Gms	18000
90349 [177966-52-8]	**5-Bromo-4-Chloro-3-Indolyl-b-D-Cellobioside (X-Cellobioside, X-B-D-Cel) extrapure, 99%**	25 Mg	2300
		100 Mg	5000
45904 [7240-90-6]	**5-Bromo-4-Chloro-3-Indolyl-b-D-Galactopyranoside (X-Gal) for tissue culture, 98%**	100 Mg	1300
		500 Mg	5000
		1 Gms	9000
10513 [7240-90-6]	**5-Bromo-4-Chloro-3-Indolyl-ß-D-Galactopyranoside (X-Gal) for molecular biology, 98%**	100 Mg	1000
		500 Mg	3400
		1 Gms	6000
19185 [93863-89-9]	**5-Bromo-6-Chloro-3-Indolyl b-D-Glucopyranoside (Magenta-Glc) for molecular biology, 98%** (Magenta Glucoside)	250 Mg	7000
		1 Gms	20000
16651 [129541-41-9]	**5-Bromo-4-Chloro-3-Indolyl-ß-D-Glucuronide Sodium Anhydrous (X-Gluc Sodium) extrapure, 98%** (X-GlcA Sodium Salt, X-ß-D-GLcA-Na , X-Glucuronide, X-Glucurono Sodium Salt, BC-Indicator) (X-GlcA	25 Mg	2200
		100 Mg	4800
		250 Mg	9800
		1 Gms	32000
55702 [370100-64-4]	**5-Bromo-4-Chloro-3-Indolyl-ß-D-Glucuronide Sodium Salt Trihydrate extrapure, 98%**	25 Mg	1200
		100 Mg	2400
		1 Gms	12000

Code	Name	Packing	Price
20402 [114162-64-0]	● **5-Bromo-4-Chloro-3-Indolyl-B-D-Glucuronide Cyclohexylammonium Salt (X-Gluc CHX) extrapure, 99%** (X-beta-Glucuronide CHX Salt)	100 Mg 250 Mg 1 Gms	2500 4000 8000
20077 [114162-64-0]	● **5-Bromo-4-Chloro-3-Indolyl-B-D-Glucuronide Cyclohexylammonium Salt (X-Gluc CHX) for tissue culture,** (X-beta-Glucuronide CHX Salt)	100 Mg 500 Mg	3000 5200
59056 [144110-43-0]	● **5-Bromo-6-Chloro-3-Indolyl-ß-D-Glucuronide Cyclohexylammonium Salt extrapure, 98%** (Magenta-GlcA, Magenta Glucuronide)	50 Mg 250 Mg	3000 6000
82224 [93863-88-8]	● **5-Bromo-6-Chloro-3-Indoxyl-B-D-Galactopyranoside (Magenta-Gal) extrapure, 98%** (Magenta-beta-D-Gal)	50 Mg 250 Mg 1 Gms	2500 6000 12000
62949 [15548-60-4]	● **5-Bromo-4-Chloro-3-Indolyl-b-D-Glucopyranoside (X-Glu, X-Glc) for tissue culture, 98%**	25 Mg 100 Mg	800 2000
51683 [15548-60-4]	● **5-Bromo-4-Chloro-3-Indolyl-ß-D-Glucopyranoside (X-Glu, X-Glc) for molecular biology, 98%** (X-Glucoside)	100 Mg 500 Mg 1 Gms	1500 3500 6000
24714 [102185-33-1]	● **5-Bromo-4-Chloro-3-Indolyl Phosphate Disodium Salt (BCIP) extrapure, 98%** (X-Phosphate Disodium Salt)	100 Mg 500 Mg 1 Gms 5 Gms	2000 6000 10500 46000
66620 [102185-33-1]	● **5-Bromo-4-Chloro-3-Indolyl Phosphate Disodium Salt (BCIP) for tissue culture, 98%**	100 Mg 500 Mg	3000 7500
14847 [6769-80-8]	● **5-Bromo-6-Chloro-3-Indolylphosphate-p-Toluidine Salt (BCIP Red) extrapure, 95%** (Magenta Phosphate p-Toludine Salt)	25 Mg 100 Mg 500 Mg	2000 4000 8000
45366 [6578-06-9]	● **5-Bromo-4-Chloro-3-Indolyl Phosphate p-Toluidine Salt (BCIP p-toluidine salt) for molecular biology,** (X-phosphate p-Toluidine Salt)	25 Mg 100 Mg 500 Mg	2450 4500 9000
75531	● **BCIP Red/NBT Solution A (BCIP Red Solution)**	100 ml	17750
73654	● **BCIP Red/NBT Solution B (NBT Solution)**	100 ml	2270
71007 [518033-33-5]	● **5-Bromo-4-Chloro-3-Indolyl-ß-D-Ribofuranoside extrapure, 95%**	25 Mg 100 Mg	5000 14000
45924 [212515-11-2]	■ **5-Bromo-4-Chloro-3-Indoxyl-myo-Inositol-1-Phosphate Ammonium Salt (X-Phos-Inositol, X-IP) extrapure,** (myo-Inositol 1-(5-Bromo-4-Chloro-3-Indoxylphosphate) Ammonium Salt)	20 Mg 100 Mg	5150 19580
62861 [160369-85-7]	● **5-Bromo-4-Chloro-3-Indoxyl-a-D-N-Acetylneuraminic Acid Sodium Salt extrapure (X-NANA.Na), 98%**	2 Mg 5 Mg	16910 28170
37869 [109-70-6]	**1-Bromo-3-Chloropropane pure, 98%** (Trimethylene Bromochloride)	100 ml 500 ml	1040 4030
78388 [112-29-8]	**1-Bromodecane pure, 98%** (Decyl Bromide)	100 ml 250 ml	1400 3100

Code	Name	Packing	Price
14433 [3832-48-2]	■ **Bromodifluoroacetyl Chloride extrapure, 97%**	1 Gms	10000
38209 [348-57-2]	✔ **1-Bromo-2,4-Difluorobenzene pure, 98%**	25 ml 100 ml	3300 11400
28503 [143-15-7]	**1-Bromododecane pure, 98%** (Dodecyl Bromide, Lauryl Bromide)	250 ml	3500
57418 [59-14-3]	● **5-Bromo-2-Deoxyuridine extrapure, 98%**	100 Mg 500 Mg 1 Gms 5 Gms	1200 3700 6300 26100
11849 [59-14-3]	● **5-Bromo-2-Deoxyuridine ExiPlus, Multi-Compendial, 98%** meets compendial specs of BP, Ph.Eur	1 Gms 5 Gms	8000 27500
59866 [2576-47-8]	✔ **2-Bromoethylamine Hydrobromide extrapure AR, 99%**	100 Gms 500 Gms 5 Kg	1200 5000 21000
94471 [592-55-2]	**2-Bromo Ethyl Ether pure, 98%**	250 ml	19700
71170 [41604-19-7]	✔ **4-Bromo-2-Fluorobiphenyl pure, 98%**	10 Gms 25 Gms 100 Gms	4400 7300 15200
58421 [629-04-9]	**1-Bromoheptane pure, 98%** (Heptyl Bromide)	100 ml 250 ml	1300 2800
52240 [112-82-3]	**1-Bromohexadecane pure, 96%** (Hexadecyl Bromide)	100 ml 250 ml 500 ml	2760 6560 8700
26459 [111-25-1]	**1-Bromohexane pure, 99%** (Hexyl Bromide)	100 ml 500 ml	1000 3100
25036 [52488-36-5]	✔ **4-Bromoindole pure, 96%**	1 Gms 5 Gms	2100 8100
54635 [10075-50-0]	✔ **5-Bromoindole pure, 99%**	5 Gms 25 Gms	500 1700
38092 [18656-96-7]	● **5-Bromo-3-Indolyl-B-D-Glucuronide Cyclohexylammonium Salt (X-GlcA CHX salt) extrapure, 97%** (Blue Glucuronide CHX, X-beta-D-GlcA.CHX)	25 Mg	18100
27348 [97753-82-7]	■ **5-Bromo-3-Indolyl-B-D-Galactopyranoside (Bluo-Gal) for molecular biology, 99%**	100 Mg 1 Gms	5400 16200
43675 [17357-14-1]	✔ **5-Bromoindoxylacetate extrapure AR, 99%**	100 Mg 500 Mg	2200 8700
85999 [80008-69-1]	● **5-Bromo-3-Indoxyl Phosphate p-Toluidine Salt extrapure, 98%** (Lapis-phos.p-tol, Lapis-phos PT)	250 Mg	12000
82645 [107-82-4]	**1-Bromo-3-Methylbutane (Stabilized w/ Ag) pure, 98.5%** (Isopentyl Bromide, Isoamyl Bromide)	100 ml 250 ml 500 ml	1700 2900 4600
89943 [90-11-9]	**1-Bromonaphthalene pure, 98%**	100 ml 250 ml	2100 5000
97167 [580-13-2]	**2-Bromonaphthalene pure, 98%**	5 Gms 25 Gms	1700 5700

Code	Name	Packing	Price
95479 ✔ [5773-80-8]	**6-Bromo-2-Naphthoic Acid pure, 98%**	5 Gms 25 Gms	2700 9600
99876 ✔ [15231-91-1]	**6-Bromo-2-Naphthol pure, 98%**	100 Gms	13200
22753 ■ [15572-30-2]	**6-Bromo-2-Naphthyl-ß-D-Galactopyranoside extrapure, 98%**	1 Gms	7200
86429 ■ [15548-61-5]	**6-Bromo-2-Naphthyl-ß-D-Glucopyranoside extrapure, 98%**	100 Mg 500 Mg	3200 14200
79945 ■ [15548-61-5]	**6-Bromo-2-Naphthyl-a-D-Mannopyranoside extrapure, 98%**	25 Mg 100 Mg	13200 48600
51879 ✔ [30007-47-7]	**5-Bromo-5-Nitro-1,3-Dioxane (Bronidox) extrapure, 99%**	1 Gms 5 Gms	1600 4000
93052 [52-51-7]	**2-Bromo-2-Nitro-1,3-Propanediol (Bronopol, BNPD) ExiPlus, Multi-Compendial, 99-101%** meets compendial specs of BP	100 Gms 500 Gms	600 2200
31868 [111-83-1]	**1-Bromooctane pure, 98%** (Octyl Bromide)	100 ml 250 ml	1000 2100
64493 [112-89-0]	**1-Bromooctadecane pure, 97%** (Octadecyl Bromide)	100 ml 500 ml	2100 9900
10459 [110-53-2]	**1-Bromopentane pure, 98%** (Pentyl Bromide, Amyl Bromide)	100 ml 250 ml 500 ml	1200 1800 3400
59505 [1119-51-3]	**5-Bromo-1-Pentene extrapure, 99%** (1-Bromo-4-Pentene, 4-Pentenyl Bromide)	5 Gms 25 Gms	2200 8500
78478 [591-20-8]	**3-Bromophenol pure, 98%**	10 Gms 25 Gms 100 Gms	1600 3300 8100
45609 [106-41-2]	**p-Bromophenol pure, 99%** (4-Bromophenol)	100 Gms 500 Gms	1800 7400
62771 ✔ [2800-80-8]	**Bromophenol Red (BPR), 90%** (3,3-Dibromophenolsulfonphthalein)	5 Gms 25 Gms	800 3000
35675 ✔ [5467-74-3]	**4-Bromophenylboronic acid extrapure, 97%**	5 Gms 25 Gms	2600 11100
88783 [86-90-8]	**4-Bromophthalic Anhydride pure, 97%**	5 Gms 25 Gms	2500 8300
98418 🌡 [223463-14-7]	**2-Bromopyridine-5-boronic acid extrapure, 95%** (Boronic Acid)	1 Gms 5 Gms	3700 14500
93829 [128-08-5]	**N-Bromosuccinimide pure, 98%**	100 Gms 500 Gms 5 Kg	500 1550 14000
51059 [128-08-5]	**N-Bromosuccinimide extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP	100 Gms 500 Gms 5 Kg	530 1600 14000
22628 [112-71-0]	**1-Bromotetradecane pure, 98%** (Myristyl Bromide)	250 ml 1000 ml	2500 9000
55570 [1003-09-4]	**2-Bromothiophene pure, 98%**	50 ml 100 ml	1800 3200
46714 [872-31-1]	**3-Bromothiophene pure, 98%**	10 ml 25 ml 100 ml	1200 2200 8100
11521 [95-46-5]	**2-Bromotoluene pure, 98%**	25 ml 100 ml	900 2200
92125 [591-17-3]	**3-Bromotoluene pure, 98%**	10 ml 25 ml	1100 2600
52250 [106-38-7]	**4-Bromotoluene extrapure, 99%**	100 Gms 500 Gms	1600 6800
27964 ● [51-20-7]	**5-Bromouracil extrapure, 99%**	5 Gms 25 Gms	1500 4900
12211 ● [957-75-5]	**5-Bromouridine extrapure, 99%**	100 Mg 500 Mg	1000 1700
42338 ✔	**Borax Carmine (Grenacher)**	5 Gms 25 Gms	600 2400
27691 [1638332-73-6]	**Brucine Hydrate Reagent Grade, 98%** (10,11-Dimethoxystrychnine, 2,3-Dimethoxystrychnidin-10-one)	5 Gms 25 Gms 100 Gms	600 2300 9000
32202 ✔ [7128-64-5]	**2,5-bis(5-tert-Butylbenzoxazol-2-yl) Thiophene (BBOT) scintillation grade, 99%**	5 Gms 25 Gms	2300 5060
34186 [109-65-9]	**n-Butyl Bromide pure, 98%** (1-bromobutane)	250 ml 500 ml 1000 ml	750 1380 2300
58390 [78-76-2]	**sec-Butyl Bromide pure, 98%** (2-bromobutane)	100 ml 500 ml	580 2300
33781 [507-19-7]	**2-Bromo-2-Methylpropane (tert-Butyl Bromide) (Stabilized w/ Ag) pure, 98%**	100 ml 500 ml	2300 8970
88194 ✔ [98-29-3]	**p-tert-Butyl Catechol pure, 98%** (1,2-Benzenediol)	100 Gms 500 Gms	460 1730
91666 ✔ [870-46-2]	**tert-Butyl Carbazate extrapure, 98%** (BOC Hydrazide)	25 Gms 100 Gms	1840 6210
22942 [507-20-0]	**tert-Butyl Chloride extrapure, 99%** (2-Chloro-2-Methylpropane)	500 ml 1000 ml	920 1380
32431 🌡 [58479-61-1]	**tert-Butylchlorodiphenylsilane (TBDPSCl) Reagent Grade, 98%** (tert-Butyldiphenylchlorosilane, tert-Butyl(chloro)diphenylsilane)	10 ml 25 ml 100 ml	1500 3450 12080
90611 [54925-64-3]	**N-tert-Butyldimethylsilyl Imidazole (TBDMSIM) pure, 97%**	1 ml 5 ml	2350 10990
39227 [18162-48-6]	**tert-Butyldimethylsilyl Chloride (TBDMS.Cl) extrapure, 98%**	25 Gms 100 Gms 500 Gms	610 2190 9430
52573 🌡 [69739-34-0]	**tert-Butyldimethylsilyl Trifluoromethanesulfonate (TBDMS triflate) extrapure, 98%** (Trifluoromethanesulfonic Acid tert-Butyldimethylsilylester)	5 ml 25 ml 100 ml	1500 6330 17250
29972 [542-69-8]	**n-Butyl Iodide extrapure, 98%**	25 ml 100 ml 500 ml	1840 6330 28410

Code	Name	Packing	Price
48792 [513-48-4]	**sec-Butyl Iodide pure, 98%** (2-Iodobutane)	25 ml 100 ml	2880 10350
28791 ✔ [79917-90-1]	**1-Butyl-3-Methylimidazolium Chloride (BMIM Cl) extrapure, 98%**	25 Gms 250 Gms	2200 10500
50807 ✔ [448245-52-1]	**1-Butyl-3-Methylimidazolium Dicyanamide (BMIM.DC) pure, 98%**	5 Gms 25 Gms	3700 12600
15748 ✔ [174501-64-5]	**1-Butyl-3-Methylimidazolium Hexafluorophosphate (BMIM PF6) extrapure, 98%**	5 Gms 25 Gms 100 Gms	1300 4200 13700
84903 ✔ [80432-09-3]	**1-Butyl-3-Methylimidazolium Tetrachloroaluminate (BMIM AlCl4) extrapure, 95%**	10 Gms 25 Gms 50 Gms	3200 5300 8400
92656 ✔ [174501-65-6]	**1-Butyl-3-Methylimidazolium Tetrafluoroborate (BMIMBF4) extrapure, 98%**	5 Gms 25 Gms	2400 5300
61005 ✔ [174899-66-2]	**1-Butyl-3-Methylimidazolium Trifluoromethanesulfonate (BMIM Otf) extrapure, 95%** (1-Butyl-3-Methylimidazolium Triflate)	1 Gms 5 Gms 25 Gms	1100 2000 7500
46178 ✔ [223437-11-4]	**1-Butyl-1-Methylpyrrolidinium Bis(trifluoromethylsulfonyl)imide (BMP TFSI) extrapure, 98%**	5 Gms 25 Gms	1800 7400
82426 ✔ [479500-35-1]	**1-Butyl-1-Methylpyrrolidinium Chloride extrapure, 98%**	5 Gms 25 Gms	2300 10000
52735 ✔ [330671-29-9]	**1-Butyl-1-Methylpyrrolidinium Hexafluorophosphate (BMP FP6) extrapure, 97.5%**	5 Gms 25 Gms	2100 8400
88656 ✔ [345984-11-4]	**1-Butyl-1-Methylpyrrolidinium Tetrafluoroborate (BMP BF4) extrapure, 97%**	1 Gms 5 Gms	900 3200
34317 ✔ [98-54-4]	**p-tert-Butyl Phenol pure, 98%**	500 Gms	810
80049 ● [18822-58-7]	**O-tert-Butyl-L-Serine extrapure, 97%**	1 Gms 5 Gms	2100 6300
12749 [106-31-0]	**n-Butyric Anhydride pure, 98%**	500 ml	2400
16490 [106-31-0]	**n-Butyric Anhydride extrapure AR, 99%**	100 ml 500 ml	800 2900
51971 [96-48-0]	**Y-Butyrolactone extrapure, 99%**	500 ml 2500 ml 25 Ltr	750 3500 31000
29850 [141-75-3]	**Butyryl Chloride extrapure, 98%**	100 ml 500 ml	350 1200
60400 ✔ [2963-78-2]	**Butyrylcholine Chloride extrapure, 96%**	1 Gms 5 Gms	2370 8200
82366 🌡 [1866-16-6]	**S-Butyrylthiocholine Iodide extrapure, 98%**	1 Gms 5 Gms 25 Gms	2500 5400 17500
56522 [561-20-6]	**Cacotheline Monohydrate extrapure, 95%** (Redox Indicator)	5 Gms 25 Gms	1000 4200
95546 ✔ [462-94-2]	**Cadaverine Free Base (1,5-Diaminopentane) extrapure, 97%** (Pentamethylenediamine, 1,5-Pentanediamine)	1 ml 5 ml 25 ml	1500 4250 8000
35547 [1476-39-7]	**Cadaverin Dihydrochloride extrapure, 98%** (1,5-Diaminopentane Dihydrochloride, Pentamethylenediamine Hydrochloride)	1 Gms 5 Gms 25 Gms	2700 4800 10300
50375 ✔ [1306-25-8]	**Cadmium Telluride ultrapure, 99.99%**	5 Gms 25 Gms	3800 11700
42062 ✔ [331-39-5]	**Caffeic Acid extrapure, 98%** (3,4-Dihydroxybenzeneacrylic Acid, 3,4-Dihydroxycinnamic Acid, 3-(3,4-Dihydroxyphenyl)-2-Propenoic A	5 Gms 25 Gms	2400 6300
65201 ✔ [1461-15-0]	**Calcein Reagent (Fluorescein Complexone) extrapure AR**	1 Gms 5 Gms 25 Gms	300 900 3400
99100 ✔ [50-14-6]	**Calciferol (Vitamin D2), 97%**	500 Mg 1 Gms 5 Gms	2000 3500 17000
47700 [9005-43-0]	**Calcium Caseinate extrapure, 90% Protein** (Casein Calcium Salt)	100 Gms 500 Gms	1000 3200
31269 [137-08-6]	**Calcium-D-Pantothenate pure, 98%** (D-Pantothenic Acid Calcium Salt, Vitamin B5)	25 Gms 100 Gms 500 Gms	500 1200 5500
65326 [3737-95-9]	**Calconcarboxylic Acid extrapure AR** (Patton & Reeders Indicator for Calcium)	5 Gms 25 Gms	300 900
66406 [3737-95-9]	**Calconcarboxylic Acid extrapure AR, ExiPlus, Multi-Compendial** (Patton & Reeders Indicator for Calcium) meets compendial specs of BP, Ph.Eur	25 Gms 100 Gms	1000 3300
42532 ✔ [3147-14-6]	**Calmagite ACS** (Calver II)	5 Gms 25 Gms 100 Gms	400 1400 4500
94340 [3144-16-9]	**D-10 Camphor Sulphonic Acid pure, 99%** (D-Camphor Sulphonic Acid)	25 Gms 100 Gms 500 Gms	1000 2500 11600
25421 [3144-16-9]	**D-10 Camphor Sulphonic Acid ExiPlus, Multi-Compendial, 99%** (D-Camphor Sulphonic Acid) meets compendial specs of BP, Ph.Eur	100 Gms 500 Gms	2700 12600
20940 [35963-20-3]	**L-10 Camphor Sulphonic Acid pure, 99%** (L-(-)-Camphor Sulphonic Acid)	25 Gms 100 Gms 500 Gms	650 2000 8000
74616 [8007-47-4]	**Canada Balsam Natural for microscopy**	100 ml 500 ml	400 1500
64210 [1135-40-6]	**CAPS Buffer extrapure, 99%** (3-[Cyclohexylamino]-1-Propane Sulphonic Acid)	25 Gms 100 Gms 500 Gms	1130 3490 13670

Code	Name	Packing	Price
39271 [73463-39-5]	**CAPSO Buffer extrapure, 99%** (3-(Cyclohexylamino) -2-Hydroxypropane Sulphonic Acid)	25 Gms 100 Gms 500 Gms	1400 3500 15000
71222 [102601-34-3]	**CAPSO Sodium Salt Buffer extrapure, 99%**	25 Gms 100 Gms 500 Gms	2200 6500 25000
54338 [86-74-8]	**Carbazole pure, 95%** (9-Azafluorene, 9H-Carbazole, Dibenzopyrrole, Diphenylenimine)	100 Gms 500 Gms	1500 6000
66964 [4800-94-6]	**Carbenicillin Disodium Salt for tissue culture, 90%**	250 Mg 1 Gms 5 Gms	1960 2860 6430
10566 [4800-94-6]	**Carbenicillin Disodium Salt, 90%**	500 Mg 1 Gms 5 Gms	880 1450 5800
32943 [638-23-3]	**Carbocysteine extrapure, 98%** (s-Carboxymethyl-l-Cysteine)	1 Gms 5 Gms	1100 3300
98536	**Carbohydrates Kit Type I**	1 Kit	4400
88515	**Carbohydrates Kit Type II**	1 Kit	8500
67852 [530-62-1]	**N,N-Carbonyldiimidazole (CDI) extrapure AR, 98%**	100 Gms 500 Gms 5 Kg	1000 4000 30000
87379 [25487-66-5]	**3-Carboxybenzeneboronic Acid extrapure, 97%** (3-Carboxyphenylboronic Acid)	5 Gms 25 Gms 100 Gms	1000 3000 9500
51171 [14047-29-1]	**4-Carboxyphenylboronic Acid extrapure, 97%**	1 Gms 5 Gms	1500 4800
35982 [23945-44-0]	**5-Carboxyuracil extrapure, 98%**	5 Gms 10 Gms	3700 7200
24352 [541-15-1]	**L-Carnitine extrapure, 99%**	25 Gms 100 Gms 500 Gms	2100 6300 21000
85641 [7235-40-7]	**b-Carotene extrapure, 90%**	1 Gms 5 Gms 25 Gms	1700 3600 8400
95727 [99-48-9]	**Carveol pure, 96%** ((-)-Carveol, Mixture of Isomers, p-Mentha-6,8-Dien-2-ol)	10 Gms 100 Gms 250 Gms	1850 17500 41000
16730 [2244-16-8]	**D-Carvone pure, 98%** ((S)-(+)-Carvone, p-Mentha-6,8-Dien-2-One)	25 Gms 100 Gms	2500 7000
71953 [9000-71-9]	**Casein Protein Rich Alkali Soluble extrapure, 80% Protein**	500 Gms 5 Kg	1600 13400
27012 [9000-71-9]	**Casein Fat Free Ultrapurified, 0.2% Fat**	500 Gms	2900
23845 [9000-71-9]	**Casein acc. to Hammarsten, 95% Protein**	100 Gms 500 Gms 1 Kg	2800 12200 21700
19575 [9000-71-9]	**Casein acc. to Hammarsten, PrimeDx, 95% Protein**	1 Kg 5 Kg	28000 125000
42168 [9001-05-2]	**Catalase (CAT) (Type 1) ex. Bovine Liver, 11700U/mg powder**	25 Mg 100 Mg	3700 8000
78695 [9001-05-2]	**Catalase (CAT) (Type 2) ex. Bovine Liver, 2000-5000U/mg protein**	1 Gms	7300
92063 [53994-73-3]	**Cefaclor (CFC), 95%**	1 Gms 5 Gms	1000 3200
15689 [66592-87-8]	**Cefadroxil Monohydrate (CFD), 95%**	1 Gms 5 Gms	1000 4200
85254 [62893-20-3]	**Cefoperazone Sodium (CPZ, Cefaperazone Sodium)**	1 Gms 5 Gms	1300 5200
26732 [74356-00-6]	**Cefotetan Disodium Salt extrapure, 85%**	25 Mg 100 Mg	6300 18900
14299 [27164-46-1]	**Cefazolin Sodium Salt (CZ), 89-110%**	100 Mg 500 Mg	4700 8000
13910 [125110-14-7]	**Cefixime Trihydrate (CFM), 95-101%**	1 Gms 5 Gms	1000 2100
13329 [80210-62-4]	**Cefpodoxime (CFPD) free acid extrapure, 96%**	25 Mg 100 Mg	3200 9500
58597 [33564-30-6]	**Cefoxitin Sodium Salt (CTX), 94%**	100 Mg 1 Gms	1300 3000
80223 [52152-93-9]	**Cefsulodin Sodium Salt (CFS) extrapure, 97%**	10 Mg 250 Mg	1000 7900
79364 [104376-79-6]	**Ceftriaxone Disodium Salt Hemiheptahydrate (CFTZ)** (Ceftriazone Sodium Salt)	1 Gms 5 Gms	1500 4200
13081 [104376-79-6]	**Ceftriaxone Disodium Salt Hemiheptahydrate (CFTZ) ExiPlus, Multi-Compendial** (Ceftriazone Sodium Salt)	1 Gms 5 Gms	3000 7250
84755 [56238-63-2]	**Cefuroxime Sodium Salt, 86-100%**	2 Gms 10 Gms 25 Gms	1400 3000 6000
85294 [528-50-7]	**D-Cellobiose extrapure**	5 Gms 25 Gms 100 Gms	1200 5100 19800
95382 [9012-54-8]	**Cellulase (Meicellase) ex. Aspergillus Niger, 0.3 U/mg powder** Cellulase (Meicellase) ex. Aspergillus Niger, 0.3U/mg powder	100 Mg 1 Gms 5 Gms 25 Gms	1600 3400 10000 42000
24801 [9012-54-8]	**Cellulase Onozuka R-10 ex. Trichoderma Viride, 10000U/g powder** (Meicellase)	1 Gms 5 Gms	8400 39900
65480 [9012-54-8]	**Cellulase Onozuka RS ex. Trichoderma Viride, 16000U/g powder** (Meicellase)	1 Gms	12700
45662 [9004-38-0]	**Cellulose Acetate Phthalate extrapure**	100 Gms 500 Gms	800 3200

Code	Name	Packing	Price
33340 [9004-34-6]	**Cellulose microcrystalline pure for TLC, 98%**	500 Gms	610
		5 Kg	5200
		25 Kg	21000
56338 [15686-71-2]	**Cephalexin Hydrate (CFL), 95-103%** (Cefalexin Hydrate)	1 Gms	1500
		5 Gms	4000
		25 Gms	11000
36413 [58-71-9]	**Cephalothin Sodium Salt (CF), 97%**	1 Gms	3400
		5 Gms	8600
19660 [64485-93-4]	**Cefotaxime Sodium Salt (CTX), 916-964ug/mg**	1 Gms	1500
		5 Gms	3200
40111 [76089-77-5]	**Cerium (III) Trifluoromethanesulfonate extrapure, 98%** (Cerium triflate, Trifluoromethanesulfonic Acid Cerium Salt)	5 Gms	9000
		25 Gms	30000
54784 [7787-69-1]	**Cesium Bromide extrapure AR, 99.5%**	10 Gms	1600
		50 Gms	6500
95125 [534-17-8]	**Cesium Carbonate extrapure, 99%**	25 Gms	2600
		100 Gms	9400
		500 Gms	43100
22966 [7647-17-8]	**Cesium Chloride extrapure AR for molecular biology, 99.9%**	25 Gms	3700
		100 Gms	14100
		500 Gms	64800
32355 [7647-17-8]	**Cesium Chloride ultrapure for molecular biology, 99.9%**	25 Gms	3900
		100 Gms	15200
		500 Gms	68100
36117 [13400-13-0]	**Cesium Fluoride extrapure AR, 99%**	25 Gms	3700
		100 Gms	13000
22941 [7789-17-5]	**Cesium Iodide extrapure, 99.5%**	10 Gms	2000
		50 Gms	9800
73254 [7789-17-5]	**Cesium Iodide ultrapure, 99.999%**	5 Gms	1900
		10 Gms	3500
		25 Gms	6400
		100 Gms	19500
52757 [1306-38-3]	**Cerium Metal Ultrapure, 200 mesh, 99.9%**	5 Gms	8500
		25 Gms	20000
88733 [7789-18-6]	**Cesium Nitrate extrapure, 99%** (Nitric Acid Cesium Salt)	5 Gms	1300
		10 Gms	1900
		50 Gms	8000
57455 [10294-54-9]	**Cesium Sulphate for molecular biology, 99.9%**	10 Gms	1700
		25 Gms	4000
		100 Gms	15500
40649 [122-18-9]	**Benzyldimethylhexadecylammonium Chloride (BDHDAC) pure, 97%** (Cetyldimethylbenzylammonium Chloride, Cetalkonium Chloride)	100 Gms	880
		1 Kg	5750
77040 [140-72-7]	**Cetylpyridinium Bromide (CPB) pure, 98%** (Hexadecylpyridinium Bromide)	100 Gms	1380
		500 Gms	5750
59088 [6004-24-6]	**Cetylpyridinium Chloride (CPC) extrapure, 98%**	100 Gms	810
		500 Gms	2880

Code	Name	Packing	Price
99698 [6004-24-6]	**Cetylpyridinium Chloride (CPC) ExiPlus, Multi-Compendial, 98%** meets compendial specs of BP, Ph.Eur	100 Gms	920
		500 Gms	3450
66302 [57-09-0]	**Cetyltrimethyl Ammonium Bromide (CTAB) extrapure AR, 99%** (Hexadecyltrimethyl Ammonium Bromide)	100 Gms	570
		500 Gms	2390
		5 Kg	21280
59771 [57-09-0]	**Cetyltrimethyl Ammonium Bromide (CTAB) extrapure AR, ExiPlus, Multi-Compendial, 99%** (Hexadecyltrimethyl Ammonium Bromide) meets compendial specs of BP, Ph.Eur	100 Gms	630
		500 Gms	2810
		5 Kg	25300
12779 [57-09-0]	**Cetyltrimethyl Ammonium Bromide (CTAB) for molecular biology, 99%** (Hexadecyltrimethyl Ammonium Bromide)	100 Gms	650
		500 Gms	3040
		5 Kg	27030
23595 [112-02-7]	**Cetyltrimethyl Ammonium Chloride (CTAC) extrapure AR, 99%** (Cetrimonium Chloride, Hexadecyltrimethyl Ammonium Chloride)	25 Gms	2330
		100 Gms	6840
		500 Gms	10250
68623 [112-02-7]	**Cetyltrimethyl Ammonium Chloride (CTAC) for molecular biology, 99%** (Cetrimonium Chloride, Hexadecyltrimethyl Ammonium Chloride)	25 Gms	3420
		100 Gms	9300
87868 [75621-03-3, 331717-45-4]	**CHAPS Buffer extrapure, 99%** (3-[(3-Cholamidopropyl) Dimethylaminonio] -1-Propanesulphonate)	1 Gms	900
		5 Gms	3200
		25 Gms	14000
		100 Gms	52000
21420 [75621-03-3, 331717-45-4]	**CHAPS Buffer for molecular biology, 99.5%**	1 Gms	1200
		5 Gms	4400
		25 Gms	18000
66903 [82473-24-3]	**CHAPSO Buffer extrapure, 99%** (3-[(3-Cholamidopropyl) Dimethylammonio]-2-Hydroxy 1-Propanesulphonate)	1 Gms	5540
		5 Gms	21990
18081 [103-47-9]	**CHES Buffer extrapure, 99%** (2-(N-Cyclohexylamino) Ethanesulphonic Acid)	25 Gms	2500
		100 Gms	5500
		500 Gms	25000
98507 [1398-61-4]	**Chitin (Poly-(b1-4)-N-acetyl glucosamine) extrapure**	25 Gms	3430
		100 Gms	10140
27391 [9012-76-4]	**Chitosan (Low MW) extrapure, 10-150m.Pas, 90% DA** (Deacetylated Chitin)	25 Gms	1200
		100 Gms	3000
		500 Gms	9500
18824 [9012-76-4]	**Chitosan (Medium MW) extrapure, 150-500m.Pas, 90% DA** (Deacetylated Chitin)	10 Gms	620
		25 Gms	1300
		100 Gms	3200
		500 Gms	11000
96266 [9012-76-4]	**Chitosan (High MW) extrapure, 800m.Pas, 90% DA** (Deacetylated Chitin)	25 Gms	1600
		100 Gms	5500
		500 Gms	21000
44097 [148411-57-8]	**Chitosan Oligosaccharide (water soluble) extrapure, 90% DA**	10 Gms	3500
		25 Gms	6000
		100 Gms	11000

Code	Name	Packing	Price
59811 [577-76-4]	🌡 **Chitosan Dimer Dihydrochloride extrapure, 98%** (Chitobiose Dihydrochloride)	5 Mg	15000
19363 [41708-93-4]	■ **Chitosan Trimer Trihydrochloride extrapure, 98%** (Chitotriose Trihydrochloride)	5 Mg	50000
62506 [117399-50-5]	■ **Chitosan Tetramer Tetrahydrochloride extrapure, 98%** (Chitotetraose Tetrahydrochloride)	5 Mg	17000
80517 [35061-50-7]	■ **Chitosan Pentamer Pentahydrochloride extrapure, 98%** (Chitopentaose Pentahydrochloride)	5 Mg	15520
36824 [41708-95-6]	■ **Chitosan Hexamer Hexahydrochloride extrapure, 95%** (Chitohexaose Hexahydrochloride)	5 Mg	30000
11975 [15879-93-3]	✔ **a-Chloralose pure, 98%** (Anhydro-D-Glucochloral, a-D-Glucochloralose)	25 Gms 100 Gms	2100 8000
77441 [87-88-7]	□ **Chloranilic Acid extrapure, 98%** (2,5-Dichloro-3,6-Dihydroxy-2,5-Cyclohexadiene-1,4-Dione, 2,5-Dichloro-3,6-Dihydroxy-p-Benzoquinone)	25 Gms 100 Gms	2400 8400
50146 [2829-43-8]	✔ **Chloranitine Fast Red BB** (C.I. No. 25380)	5 Gms 25 Gms	900 3000
97686 [56-75-7]	● **Chloramphenicol (CMP), 98-102%**	5 Gms 25 Gms 100 Gms	1000 2000 5000
28812 [56-75-7]	● **Chloramphenicol (CMP) for tissue culture, 98-102%**	5 Gms 25 Gms	1200 2250
57956 [56-75-7]	● **Chloramphenicol (CMP) ExiPlus, Multi-Compendial, 98.0-102%** meets compendial specs of USP, BP	25 Gms 100 Gms	2750 6000
37281 [539-03-7]	**p-Chloroacetanilide extrapure AR, 99%**	25 Gms 100 Gms	730 2220
77863 [539-03-7]	**p-Chloroacetanilide extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms	2360
70179 [532-27-4]	**2-Chloroactophenone pure, 97%** (Phenacyl Chloride)	25 Gms 100 Gms	2650 8400
40585 [99-91-2]	**p-Chloroacetophenone extrapure, 98%**	100 ml 500 ml 1000 ml	1640 6630 12520
92804 [108-42-9]	**m-Chloroaniline (MCA) pure, 99%**	500 ml 2500 ml	1610 7020
43265 [89-98-5]	**o-Chlorobenzaldehyde pure, 98%** (2-Chlorobenzaldehyde)	100 ml 500 ml 2500 ml	350 870 3800
17927 [104-88-1]	✔ **p-Chlorobenzaldehyde pure, 97%** (4-Chlorobenzaldehyde)	100 Gms 500 Gms 5 Kg	520 1210 10930
91019 [119-56-2]	✔ **4-Chlorobenzhydrol pure, 98%**	50 Gms 100 Gms	2760 4830
69159 [118-91-2]	**o-Chlorobenzoic Acid pure, 98%**	100 Gms 500 Gms	350 1380
30519 [118-91-2]	**o-Chlorobenzoic Acid extrapure AR, 99.5%**	100 Gms 500 Gms	460 1960
34366 [535-80-8]	**m-Chlorobenzoic Acid extrapure, 99%**	100 Gms 500 Gms	1200 4000
24711 [74-11-3]	**p-Chlorobenzoic Acid pure, 99%**	100 Gms 500 Gms	380 1380
97199 [74-11-3]	**p-Chlorobenzoic Acid extrapure AR, 99.7%**	100 Gms	690
53378 [134-85-0]	**p-Chlorobenzophenone extrapure, 99%**	100 Gms 500 Gms	690 2530
98001 [89-97-4]	**2-Chlorobenzylamine pure, 98%** (o-Chlorobenzylamine)	25 ml 100 ml	1380 4830
95074 [628-20-6]	**4-Chlorobutyronitrile pure, 98%**	100 ml 500 ml	3480 15120
29440 [4635-59-0]	**4-Chlorobutyryl Chloride pure, 98%**	250 ml 500 ml 1000 ml	1270 2300 4140
23843 [999-81-5]	**Chlorocholine Chloride (CCC) extrapure, 98%** (2-Chloroethyl)trimethylammonium chloride, Chlormequat chloride, Choline dichloride, Cycocel)	5 Gms 25 Gms	1500 4600
74430 [12092-47-6]	**Chloro(1,5-cyclooctadiene) Rhodium (I) Dimer extrapure, 98%**	100 Mg 1 Gms	15000 130000
10231 [3140-73-6]	✔ **2-Chloro-4,6 Dimethoxy-1,3,5 Triazine extrapure, 99%**	5 Gms 25 Gms 100 Gms	1040 4030 13800
69297 [97-00-7]	**1-Chloro-2,4-Dinitrobenzene Reagent Grade, 99%**	100 Gms 500 Gms	350 920
93937 [97-00-7]	**1-Chloro-2,4-Dinitrobenzene ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP	100 Gms 500 Gms	520 1730
17691 [2578-45-2]	✔ **2-Chloro-3,5-Dinitropyridine extrapure AR, 99%**	1 Gms 5 Gms	5890 18080
38825 [1679-18-1]	**4-Chlorophenylboronic Acid extrapure, 97%**	5 Gms 25 Gms	2070 8550
22654 [16672-87-0]	**2-Chloroethylphosphonic Acid (ETHREL, Ethephon) 40% soln.**	100 ml 250 ml	920 1730
10455 [367-21-5]	✔ **3-Chloro-4-Fluoroaniline pure, 98%**	25 Gms 100 Gms 500 Gms	810 1610 6900
17703 [2252-51-9]	**2-Chloro-4-Fluorobenzoic Acid pure, 98%**	5 Gms 25 Gms	1900 6150
36393 [2252-50-8]	**2-Chloro-5-Fluorobenzoic Acid pure, 98%**	10 Gms 100 Gms	3390 17600
84538 [434-75-3]	**2-Chloro-6-Fluorobenzoic Acid pure, 98%**	5 Gms 25 Gms	3180 7420

Code	Name	Packing	Price
33676 [161957-55-7]	**3-Chloro-2-Fluorobenzoic Acid pure, 98%**	5 Gms 10 Gms	5160 8840
14208 [1035373-85-3]	□ **2-(2-((6-Chlorohexyl)oxy)Ethoxy) Ethylamine Hydrochloride extrapure, 97%**	100 Mg 1 Gms	4000 15000
69932 [157381-11-8]	● **2-Chloro-4-Nitrophenyl 4-O-b-D-Galactopyranosylmaltosid e (Gal-G2-a-CNP, Gal CNPG2) extrapure, 90%** (a-(2-Chloro-4-Nitrophenyl) -ß-1,4-Galactopyranosylmaltoside)	250 Mg 1 Gms 5 Gms	7000 24000 110000
85563	**Chloroform : Isoamyl Alcohol (24:1) for molecular biology**	100 ml 500 ml	1300 2730
41117 [936-59-4]	✔ **3-Chloro-1-Phenyl Propanone pure, 98%** (3-Chloropropiophenone, B-Chloropropiophenone)	25 Gms 100 Gms 500 Gms	2070 4830 19550
72928 [1914-99-4]	**Chlorophosphonazo III for Spectrophotometry** (Bis(4-chloro-2-Phosphonobenzolazo) Chromotropic acid)	100 Mg 500 Mg	4180 17600
51562 [6001-76-9]	🌡 **Chlorogenic Acid Hemihydrate extrapure, 98%**	100 Mg 1 Gms	480 2300
78193 [10310-21-1]	● **6-Chloroguanine extrapure, 98%** (2-Amino-6-Chloropurine)	1 Gms 5 Gms	1900 7580
99721 [159954-35-5]	■ **6-Chloro-3-Indolyl-Caprylate (Salmon Caprylate) for molecular biology, 97%** (6-Chloro-3-Indoxyl-Octanoate)	50 Mg	6000
79621 [135313-63-2]	● **4-Chloro-3-Indolyl-ß-Galactopyrano side extrapure, 98%**	100 Mg	6050
98551 [138182-21-5]	● **6-Chloro-3-Indolyl ß-D-Galactopyranoside (Salmon Gal) extrapure, 98%** (Red Gal, Rose Gal)	100 Mg 500 Mg 1 Gms	2000 3200 6000
82947 [138182-20-4]	● **6-Chloro-3-Indolyl-ß-D-Glucuronide Cyclohexylammonium Salt extrapure, 97%** (RED-GLUC, Rose-Glucuronide CHA Salt, Salmon-Glucuronide CHA Salt)	50 Mg 1 Gms	4000 18000
89432 [467214-46-6]	■ **6-Chloro-3-Indoxyl-a-D-Glucopyran oside (Salmon-a-D-Glc) extrapure, 98%** (Salmon-alpha-D-Glucoside)	100 Mg 500 Mg	2500 9500
71848 [130-26-7]	**5-Chloro-7-Iodo-8-Quinolinol (Clioquinol) ExiPlus, Multi-Compendial, 97-103%** (5-Chloro-8-Hydroxy-7-Iodoquinoline, Iodochlorohydroxyquinoline) meets compendial specs of USP, BP, Ph.Eur	25 Gms 100 Gms 500 Gms	2500 6500 24700
25898 [14110-97-5]	■ **4-Chloromercuribenzenesulphonic Acid Monosodium Salt (pCMBS) extrapure, 95%** (p-Chloromercuriphenylsulfonic Acid Sodium Salt)	10 Mg	30600
26004 [138-85-2]	✔ **p-Chloromercuribenzoic Acid Sodium Salt extrapure, 98%**	1 Gms 5 Gms	3390 12260
81929 [604-44-4]	● **4-Chloro-1-Naphthol Substrate Grade extrapure, 99%**	2.5 Gms 5 Gms 25 Gms	3030 4580 22150
85267 [121-73-3]	**1-Chloro-3-Nitrobenzene pure, 98%**	250 Gms 500 Gms	630 1200
77459 [2516-96-3]	**2-Chloro-5-Nitrobenzoic Acid pure, 99%**	25 Gms 100 Gms 500 Gms	420 1210 4600
21259 [157843-41-9]	● **2-Chloro-4-Nitrophenyl a-L-Fucopyranoside (CNPF) extrapure, 98%** (AFU)	25 Mg	9720
48958 [118291-90-0]	● **2-Chloro-4-Nitrophenyl a-D-Maltotrioside (CNPG3) extrapure, 95%**	25 Mg 100 Mg	3920 9000
66684 [3970-40-9]	**2-Chloro-3-Nitrotoluene extrapure, 98%**	1 Gms 5 Gms	1730 6900
41107 [108-43-0]	**3-Chlorophenol pure, 98%** (m-Chlorophenol)	25 Gms 100 Gms 500 Gms	1130 4030 19550
94280 [95-57-8]	**o-Chlorophenol pure, 98%**	500 ml 1000 ml 2500 ml	1020 1990 3850
74230 [106-48-9]	**p-Chlorophenol pure, 98%** (4-Chlorophenol)	500 Gms	1270
61937 [106-48-9]	**p-Chlorophenol PrimeDx, 98%** (4-Chlorophenol)	1 Kg	4370
18344 [4430-20-0]	**Chlorophenol Red**	5 Gms 25 Gms	720 2800
87976 [99792-79-7]	● **Chlorophenol Red-ß-D-Galactopyranoside (CPRG) extrapure, 90%**	25 Mg	6000
21524 [122-88-3]	✔ **Chlorophenoxyacetic Acid (4-CPA) for tissue culture, 98%** (Parachlorophenoxyacetic acid)	25 Gms 100 Gms 500 Gms	800 2800 12000
91145 [101-10-0]	✔ **2-(3-Chlorophenoxy)-Propionic Acid (3-CPA) pure, 99%**	10 Gms 25 Gms	2000 4500
28272 [1878-66-6]	✔ **4-Chlorophenylacetic Acid pure, 98%**	100 Gms 500 Gms	2000 9000
40881 [327-97-9]	✔ **Chlorogenic Acid extrapure, 98%** (1,3,4,5-Tetrahydroxycyclohexanecarb oxylic acid 3-(3,4-dihydroxycinnamate), 3-(3,4-Dihydroxycinnamoy	100 Mg 1 Gms 5 Gms	500 2500 11750
23727 [598-78-7]	**2-Chloropropionic Acid pure, 97%**	100 ml 500 ml	1560 6770
50217 [29617-66-1]	**S-2-Chloropropionic Acid pure, 98%** (S-(-)-2-Chloropropionic Acid)	5 ml 25 ml 100 ml	1130 3520 12680
70774 [34841-35-5]	**3-Chloropropiophenone pure, 98%**	25 Gms 100 Gms 500 Gms	3250 11780 51500
53313 [6285-05-8]	✔ **4-Chloropropiophenone pure, 98%**	100 Gms 500 Gms	3260 14690

Code	Name	Packing	Price
62020 [87-42-3]	● **6-Chloropurine extrapure, 99%**	1 Gms	880
		5 Gms	3110
		25 Gms	12300
57430 [2004-06-0]	● **6-Chloropurine Riboside extrapure**	250 Mg	2160
		1 Gms	7310
18741 [95-88-5]	**4-Chlororesorcinol pure, 98%** (1-Chloro-2,4-Dihydroxybenzene, 1,3-Dihydroxy-4-Chlorobenzene)	100 Gms	2300
		500 Gms	9200
26176 [321-14-2]	**5-Chlorosalicylic Acid pure, 98%**	100 Gms	660
		500 Gms	1900
59838 [106-43-4]	**p-Chlorotoluene pure, 98%**	500 ml	1510
		1000 ml	2550
97389 [67-97-0]	🌡 **Cholecalciferol (Vitamin D3) for tissue culture, 97%** (7-Dehydrocholesterol Activated)	1 Gms	1400
		5 Gms	5700
		25 Gms	27000
54181 [57-88-5]	● **Cholesterol ex. Sheep Wool (Lanolin) extrapure AR, 99%**	5 Gms	850
		25 Gms	2000
		100 Gms	7500
		500 Gms	33000
38615 [57-88-5]	● **Cholesterol ex. Sheep Wool (Lanolin) Reference Standard Grade for molecular biology, 99%**	5 Gms	1100
		25 Gms	3200
		100 Gms	11000
51626 [9026-00-0]	■ **Cholesterol Esterase (CE) ex. Porcine Pancreas, 35U/mg powder**	10 Mg	9000
		50 Mg	40000
15787 [9028-76-6]	■ **Cholesterol Oxidase ex. Microorganism, 10U/mg powder**	250 Units	3000
		1 K.Units	5000
		2.5 K.Units	7500
		5 K.Units	12000
		25 K.Units	52000
77645 [604-35-3]	● **Cholesteryl Acetate extrapure, 99%**	5 Gms	860
		25 Gms	3800
54908 [81-25-4]	✔ **Cholic Acid extrapure, 98%**	25 Gms	3000
		100 Gms	10000
		500 Gms	48000
53645 [361-09-1, 206986-87-0]	✔ **Cholic Acid Sodium Salt Hydrate pure, 99%** (Sodium Cholate Hydrate)	25 Gms	3200
		100 Gms	11000
		500 Gms	50000
93319 [4197-07-3]	✔ **Chromotrope 2R**	10 Gms	700
		25 Gms	1400
31723 [5808-22-0]	**Chromotropic Acid Disodium Salt Dihydrate ACS, 99%**	10 Gms	650
		25 Gms	1300
		100 Gms	3600
		500 Gms	14000
72580 [5808-22-0]	**Chromotropic Acid Disodium Salt Dihydrate ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	10 Gms	800
		25 Gms	1500
		100 Gms	4200
		500 Gms	16000
35085 [9004-07-3]	■ **a-Chymotrypsin ex. Bovine Pancreas, 40 BTEE U/mg protein**	25 Mg	2700
		250 Mg	5000
		1 Gms	8000
		5 Gms	25000
13225 [9004-07-3]	■ **a-Chymotrypsin ex. Porcine, 40 BTEE U/mg protein**	250 Mg	2500
		1 Gms	4500
		5 Gms	21000
87063 [9035-75-0]	■ **a-Chymotrypsinogen A (5x cryst) ex. Bovine Pancreas, 45U/mg protein**	250 Mg	7000
88213 [1885-38-7]	**Cinnamonitrile extrapure, 97%** (3-Phenylacrylonitrile)	25 ml	880
		100 ml	2620
26591 [102-92-1]	🌡 **Cinnamoyl Chloride pure, 98%**	100 ml	2370
		500 ml	9380
78079 [86393-32-0]	● **Ciprofloxacin Hydrochloride Hydrate (CIP), 98-102%**	1 Gms	500
		5 Gms	1200
29965 [5392-40-5]	**Citral pure, 95%** (3,7-Dimethyl-2,6-Octadienal)	100 ml	580
		500 ml	2400
13294 [8000-29-1]	✔ **Citronella Oil extrapure, 34%**	250 ml	2750
		500 ml	5040
35417 [372-75-8]	**L-Citrulline extrapure CHR, 99%**	100 Gms	600
		500 Gms	2800
94708 [372-75-8]	**L-Citrulline ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP	100 Gms	750
		500 Gms	3200
41694 [99-11-6]	**Citrazinic Acid extrapure, 97%** (2,6-Dihydroxypyridine-4-Carboxylic acid, Dihydroxyisonicotinic Acid)	25 Gms	2400
		100 Gms	5300
		500 Gms	16300
78746 [8000-34-8]	**Clove Oil extrapure, 80%**	100 ml	800
		500 ml	3000
68320 [128595-07-3]	● **CP-OSu (N-(Cyclopentyloxycarbonyloxy)Succinimide) extrapure, 99%**	25 Gms	4500
		100 Gms	13000
68636 [7081-44-9]	● **Cloxacillin Sodium Monohydrate (CSX), 95-102%**	1 Gms	800
		5 Gms	3200
31104 [55779-48-1]	■ **Coelenterazine Native (CLZN) ex. Aequorea Sp., 95%**	250 mcg	15100
		1 Mg	24800
40804 [85-61-0]	■ **Coenzyme A Free Acid extrapure, 90%**	10 Mg	3000
		25 Mg	5000
		100 Mg	15000
95225 [55672-92-9 (anhy)]	■ **Coenzyme A Sodium Salt Hydrate (CoA Na2) for molecular biology, 85%**	10 Mg	7500
		25 Mg	18000
64147 [18439-24-2]	■ **Coenzyme A Trilithium Salt ex. Yeast extrapure, 90% CoA**	25 Mg	5600
15039 [303-98-0]	■ **Coenzyme Q10 extrapure, 99%**	100 Mg	800
		500 Mg	2600
		1 Gms	4800
		5 Gms	11500

Code	Name	Packing	Price
17701 ✔ [64-86-8]	**Colchicine extrapure, 97%**	1 Gms	3600
		5 Gms	17000
		10 Gms	30000
93156 ✔ [64-86-8]	**Colchicine ExiPlus, Multi-Compendial, 97%** meets compendial specs of USP, BP, Ph.Eur	1 Gms	4500
		5 Gms	20000
45830 ✔ [64-86-8]	**Colchicine for tissue culture, 97%**	1 Gms	4000
		5 Gms	18500
51681 🌡 [1264-72-8]	**Colistin Sulphate (CL)**	500 Mg	1800
		1 Gms	3500
		5 Gms	10400
90443 ✔ [9007-34-5]	**Collagen ex. Marine Fish extrapure, 95%**	5 Gms	2400
		25 Gms	7000
		100 Gms	22000
51581 ✔	**Colloidal Protein Staining Solution**	100 ml	3600
94278 [108-75-8]	**2,4,6-Collidine pure, 98%** (2,4,6-Trimethylpyridine)	100 ml	1500
		500 ml	6600
38496 ■ [1391-06-6]	**Conalbumin (Ovotransferrin) ex. Chicken Egg White, 80% Protein**	100 Mg	5000
		500 Mg	20000
24275 [573-58-0]	**Congo Red ACS, 75%**	25 Gms	400
		100 Gms	1400
		500 Gms	6500
10408 ● [8064-90-2]	**Cotrimoxazole (CTX), 95%**	100 Mg	2700
		500 Mg	12600
64305 [501-98-4]	**p-Coumaric Acid pure, 98%** (Trans-4-Hydroxycinnamic Acid)	5 Gms	2000
		25 Gms	3200
		100 Gms	8500
42074 ✔ [91-64-5]	**Coumarin pure, 98%** (1-Benzopyran-2-one)	100 Gms	600
		500 Gms	2200
76072 [6020-87-7]	**Creatine Monohydrate extrapure, 99%**	25 Gms	400
		100 Gms	1000
		1 Kg	9500
25758 ● [71519-72-7]	**Creatine Phosphate Disodium Salt Tetrahydrate extrapure, 99%** (Sodium Creatine Phosphate Dibasic Tetrahydrate, Phosphocreatine Disodium Salt Tetrahydrate)	5 Gms	3000
		25 Gms	11000
		100 Gms	36000
25074 ■ [9025-13-2]	**Creatininase ex. Recombinant E.Coli, 500U/mg powder**	2 K.Units	4200
		5 K.Units	8400
65142 ✔ [60-27-5]	**Creatinine extrapure CHR, 99%**	25 Gms	900
		100 Gms	3200
		500 Gms	12600
19428 ✔ [60-27-5]	**Creatinine extrapure AR, 99%**	25 Gms	1000
		100 Gms	3600
		500 Gms	14700
62000 ■ [37289-15-9]	**Creatinine Deiminase ex. Microorganism, 10U/mg powder**	250 Units	10300

Code	Name	Packing	Price
29121 [596-27-0]	**o-Cresolphthalein AR, 95%**	5 Gms	200
		25 Gms	700
		100 Gms	2400
64982 [2411-89-4]	**o-Cresolphthalein Complexone (Phthalein Purple) ACS, 80%**	1 Gms	600
		5 Gms	2000
		25 Gms	6300
72278 [94442-10-1]	**o-Cresolphthalein Complexone Disodium Salt extrapure**	1 Gms	1400
		5 Gms	5800
31393 [2303-01-7]	**m-Cresol Purple extrapure AR**	1 Gms	300
		5 Gms	1200
60837 [107-93-7]	**Crotonic Acid pure, 99%**	100 Gms	700
		250 Gms	1200
		500 Gms	2200
51665 ■	**r-CRP ex. E.Coli (Recombinant C-reactive protein) extrapure**	1 Mg	12400
12295 [3572-06-3]	**Cuelure pure, 95%** (4-[4-(Acetyloxy)phenyl]-2-Butanone)	25 Gms	800
		100 Gms	2700
		500 Gms	12300
		1 Kg	23400
		2.5 Kg	51800
37870 ✔ [126747-14-6]	**4-Cyanophenylboronic Acid extrapure, 95%**	1 Gms	1600
		5 Gms	6900
		25 Gms	20000
86320 ✔ [150255-96-2]	**3-Cyanophenylboronic Acid extrapure, 97%**	1 Gms	1500
		5 Gms	3600
		25 Gms	14700
57446 🌡 [108-77-0]	**Cyanuric Chloride extrapure, 99%** (2,4,6-Trichloro-1,3,5-Triazine)	100 Gms	400
		250 Gms	700
		500 Gms	1200
51458 ● [59016-56-7]	**1-Cyano-4-Dimethylaminopyridiniu m Tetrafluoroborate (CDAP) extrapure, 97.5%**	100 Mg	7200
		250 Mg	16400
		1 Gms	63300
57902 ✔ [879866-74-7]	**1-(3-Cyanopropyl) -3-Methylimidazolium Dicyanamide (CPMIM.DC) pure, 98%** (1-Methyl-3-(3-Cyanopropyl)imidazolium Dicyanamide)	1 Gms	14000
		5 Gms	60000
65291 ● [506-68-3]	**Cyanogen Bromide pure, 98%**	25 Gms	800
		100 Gms	2000
31787 ● [506-68-3]	**Cyanogen Bromide ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, Ph.Eur	25 Gms	1000
		100 Gms	2200
81778 🌡 [506-68-3]	**Cyanogen Bromide Solution (5.0M in Acetonitrile)**	100 ml	6300
94952 🌡 [5457-28-3]	**3-Cyanoindole pure, 98%**	5 Gms	3100
		25 Gms	12400
33202 🌡 [15861-24-2]	**5-Cyanoindole pure, 98%**	10 Gms	27800
		25 Gms	68900
83244 🌡 [7585-39-9]	**b-Cyclodextrin, base (BCD) for tissue culture, 98%**	100 Gms	700
		500 Gms	2100

Code	Name	Packing	Price
76496 [7585-39-9]	**ß-Cyclodextrin, base (BCD), 98%**	25 Gms 100 Gms 500 Gms	700 2000 4200
98395 [2043-61-0]	**Cyclohexanecarboxaldehyde pure, 97%** (Hexahydrobenzaldehyde)	25 ml 100 ml	3700 14300
37594 [637-88-7]	**1,4-Cyclohexanedione pure, 99%**	25 Gms 100 Gms 500 Gms	1300 4400 21000
55799 [108-93-0]	**Cyclohexanol extrapure, 99%**	500 ml 2500 ml	700 2800
82626 [636-82-8]	**1-Cyclohexene-1-Carboxylic Acid extrapure, 97%** (3,4,5,6-Tetrahydrobenzoic Acid)	1 Gms 5 Gms	2400 10500
69160 [100-50-5]	**3-Cyclohexene-1-Carboxaldehyde extrapure, 98%** (1,2,3,6-Tetrahydrobenzaldehyde, 4-Formyl-1-cyclohexene)	25 ml 100 ml	2100 5500
31668 [110-83-8]	**Cyclohexene (Tetrahydrobenzene) pure, 99%**	100 ml 500 ml 2500 ml	400 1800 8200
86620 [66-81-9]	**Cycloheximide extrapure, 94%** (Actidione)	100 Mg 1 Gms 5 Gms	500 2400 9900
95581 [66-81-9]	**Cycloheximide for tissue culture, 94%** (Actidione)	100 Mg 1 Gms 5 Gms	800 2800 11400
38772 [120-92-3]	**Cyclopentanone (Keto Cyclopentane) pure, 99%** (Keto Cyclopentane)	100 ml 500 ml 2500 ml	370 1210 5520
88410 [5614-37-9]	**Cyclopentyl Methyl Ether (CPME) Reagent Grade, 99.5%**	250 ml 500 ml 1000 ml	1600 3000 5800
74083 [1759-53-1]	**Cyclopropanecarboxylic Acid (CPC-Acid) pure, 98%**	25 Gms 100 Gms 250 Gms	1600 3200 7400
62160 [411235-57-9]	**Cyclopropylboronic Acid Monohydrate extrapure, 97%**	1 Gms 5 Gms	900 3400
41269 [68-41-7]	**D-Cycloserine (D-CSR), 98%** (D-4-Amino-3-Isoxazolidone)	1 Gms 5 Gms	1500 5000
73383 [921-01-7]	**D-Cysteine (base) extrapure, 99%** ((S)-2-Amino-3-Mercaptopropionic Acid)	1 Gms 5 Gms 25 Gms	1800 3700 14000
83610 [207121-46-8]	**D-Cysteine Hydrochloride Monohydrate extrapure, 98%**	5 Gms 25 Gms	2000 4000
98973 [65-46-3]	**Cytidine extrapure, 99%**	25 Gms 100 Gms 500 Gms	850 2800 12000
87600 [63-37-6]	**Cytidine-5-Monophosphate (5-CMP) extrapure, 99%**	1 Gms 5 Gms	620 2450
48666 [81012-87-5]	**Cytidine-5-Triphosphate Disodium Salt Dihydrate (5-CTP.Na2) extrapure, 97%**	100 Mg 500 Mg 1 Gms	1250 4500 7000
81551 [9007-43-6]	**Cytochrome C (Oxidized) (Type 1) ex. Horse Heart extrapure, 90%**	25 Mg 100 Mg 500 Mg 1 Gms	3300 10500 43000 80000
90777 [9007-43-6]	**Cytochrome C (Oxidized) (Type 2) ex. Horse Heart, 95%**	25 Mg 100 Mg 500 Mg 1 Gms	4000 13000 60000 115000
62682 [147-94-4]	**Cytosine-B-D-Arabinofuranoside (Ara-C, Cytarabine) extrapure, 98%** (Cytosine Arabinoside)	100 Mg 1 Gms 5 Gms	2000 3000 7000
75543 [146998-31-4]	**DABCYL-N-Succinimidyl Ester extrapure, 98%** (4-[4-(Dimethylamino)phenylazo]Benzoic Acid N-Succinimidyl Ester)	50 Mg 250 Mg	3000 10000
31920 [552-66-9]	**Daidzin extrapure, 98%** (Daidzein-7-O-ß-D-Glucopyranoside)	25 Mg 100 Mg	5000 12000
90441 [119478-55-6]	**Danofloxacin Mesylate (DM), 98%**	100 Mg	2500
64018 [605-65-2]	**Dansyl Chloride (DNSCl) extrapure, 99%** (5-Dimethylaminonaphthalene-1-Sulfonyl Chloride)	1 Gms 5 Gms 25 Gms	3600 12000 35000
18668 [28718-90-3]	**DAPI Dihydrochloride (4,6-Diamidino-2-Phenylindole Dihydrochloride) for molecular biology, 95%**	10 Mg 25 Mg	3000 6000
28488 [28718-90-3]	**DAPI Dihydrochloride (4,6-Diamidino-2-Phenylindole Dihydrochloride) for tissue culture, 95%**	100 Mg	4000
80145 [28718-90-3]	**DAPI Dihydrochloride (4,6-Diamidino-2-Phenylindole Dihydrochloride) Solution (1mg/ml in H2O)**	1 ml 5 x1 ml	750 2850
22499 [38078-09-0]	**DAST (Diethylamino Sulfurtrifluoride) extrapure AR, 95%**	1 Gms 5 Gms 25 Gms 100 Gms	675 3000 10000 37000
11915 [22059-21-8]	**DBCO Acid extrapure, 97%**	250 Mg 1 Gms	3000 7000
40396 [14933-08-5]	**DDAPS (N-Dodecyl-N,N-dimethyl-3-ammonio-1-propanesulfonate, Lauryl sulfobetaine) extrapure, 98%**	2 Gms 10 Gms 50 Gms	990 2960 14500
47478 [9013-34-7]	**DEAE Cellulose 11 analytical grade for molecular biology**	5 Gms 25 Gms	2450 10040
57136 [9013-34-7]	**DEAE Cellulose 23 analytical grade for molecular biology**	5 Gms 25 Gms	2760 11630
92321 [9013-34-7]	**DEAE Cellulose 32 analytical grade for molecular biology**	5 Gms 25 Gms	3430 14440

Code	Name	Packing	Price
10529 [9013-34-7]	**DEAE Cellulose 52 analytical grade for molecular biology**	5 Gms	2170
		25 Gms	10060
		100 Gms	35000
54665 [91-17-8]	**Decalin (Decahydronaphthalene) pure, 98%**	500 ml	1900
		2500 ml	9000
47112 [112-47-0]	**1,10-Decanediol pure, 98%** (Decamethylene Glycol, 1,10-Dihydroxydecane)	100 Gms	1500
		500 Gms	6500
72007 [143-10-2]	**1-Decanethiol pure, 95%**	100 ml	3200
		500 ml	13000
91106 [112-13-0]	**Decanoyl Chloride pure, 98%** (Capric Acid Chloride)	100 ml	2000
		500 ml	7600
48467 [112-30-1]	**Decyl Alcohol pure, 99%** (1-Decanol)	100 ml	280
		500 ml	870
		2500 ml	4000
43940 [2082-84-0]	**Decyltrimethylammonium Bromide (DTMAB) extrapure, 98%**	25 Gms	2500
		100 Gms	8000
83158 [16373-93-6]	● **2-Deoxyadenosine Monohydrate extrapure, 98%**	1 Gms	1270
		5 Gms	5660
39192 [2922-74-9]	● **2'-Deoxyadenosine-5'-Monophosphate Disodium Salt (dAMP-Na2) extrapure, 98%**	100 Mg	2000
49092 [951-77-9]	● **2-Deoxycytidine extrapure, 98%**	1 Gms	2500
55193 [13085-50-2]	● **2`-Deoxycytidine -5`-Monophosphate Disodium Salt (dCMP-Na2) extrapure, 98%**	500 Mg	10000
84476 [1949-89-9]	● **2-Deoxy-D-Galactose extrapure, 98%** (2-Deoxygalactose)	1 Gms	6500
		5 Gms	28000
43421 [154-17-6]	● **2-Deoxy-D-Glucose (2-DG) extrapure AR, 99.9%**	1 Gms	900
		5 Gms	4100
		25 Gms	20000
22350 [312693-72-4]	● **2-Deoxyguanosine Monohydrate extrapure, 99%**	1 Gms	6000
25695 [52558-16-4, 33430-61-4]	● **2'-Deoxyguanosine -5'-Monophosphate Sodium Salt (dGMP) extrapure, 98%**	1 Gms	17000
16516 [890-38-0]	● **2'-Deoxyinosine extrapure, 98%**	1 Gms	2800
61824 [9003-98-9]	■ **Deoxyribonuclease I (DNase I) - Type 1A ex. Bovine Pancreas for MB, 2000Kunitz U/mg, 80% Protein**	10 K.Units	2420
		50 K.Units	7700
		200 K.Units	25300
15409 [9003-98-9]	■ **Deoxyribonuclease I (DNase I) - Type 1B ex. Bovine Pancreas for MB, 2000Kunitz U/mg, 80% Protein**	10 Mg	3910
		50 Mg	17600
		100 Mg	31900
14658 [9003-98-9]	■ **Deoxyribonuclease I - Type 2A (DNase I) ex. Bovine Pancreas for MB, 180Kunitz U/mg, 80% Protein**	5000 Units	1100
		50000 Units	5600
		100000 Units	10560

Code	Name	Packing	Price
52942 [9003-98-9]	■ **Deoxyribonuclease I - Type 2B (DNase I) ex. Bovine Pancreas for MB,180Kunitz U/mg protein,85%Protein**	10 Mg	620
		100 Mg	5880
		1 Gms	22030
74322 [100403-24-5]	■ **Deoxyribonucleic Acid Free Acid (DNA) ex. Fish Sperm**	5 Gms	2100
		10 Gms	4000
		25 Gms	8000
30295 [9007-49-2, 68938-01-2]	■ **Deoxyribonucleic Acid Sodium Salt (DNA Sodium Salt) ex. Salmon Milt** (SS-DNA, sm-DNA)	5 Gms	2250
		25 Gms	8500
68451 [73049-39-5]	■ **Deoxyribonucleic Acid Sodium Salt (DNA Sodium Salt Highly Polymerized) ex. Calf Thymus extrapure**	25 Mg	2750
		100 Mg	5000
		250 Mg	17000
84384 [533-67-5]	● **2-Deoxy-D-Ribose extrapure**	1 Gms	800
		5 Gms	1990
		25 Gms	6000
		100 Gms	18000
92627 [951-78-0]	● **2'-Deoxyuridine extrapure, 98%**	1 Gms	1800
		5 Gms	6000
41602 [42155-08-8]	● **2`-Deoxyuridine-5`-Monophosphate Disodium Salt (dUMP-Na2) extrapure, 98%**	100 Mg	3000
		1 Gms	25000
41143 [83-44-3]	● **Deoxycholic Acid extrapure, 99%**	25 Gms	3600
		100 Gms	12800
24417 [302-95-4]	● **Deoxycholic Acid Sodium Salt (Sodium Deoxycholate) Bacto grade, 99%**	25 Gms	2700
		100 Gms	9700
		500 Gms	47300
96876 [145224-92-6]	● **Deoxycholic Acid Sodium Salt Monohydrate (Sodium Deoxycholate Monohydrate) Bacto grade, 99%**	25 Gms	3000
		100 Gms	10000
		500 Gms	49700
43767 [87413-09-0]	🌡 **Dess Martin Periodinane extrapure, 97%**	5 Gms	750
		25 Gms	3400
		100 Gms	12500
47601 [9004-54-0]	🌡 **Dextran ex. Leuconostoc Sp. - Tech 10 (9-11)**	5 Gms	3000
		25 Gms	4000
		100 Gms	7800
		500 Gms	26000
41570 [9004-54-0]	🌡 **Dextran ex. Leuconostoc Sp. - Tech 40 (35-45)**	5 Gms	3000
		25 Gms	4000
		100 Gms	7500
		500 Gms	25000
12891 [9004-54-0]	🌡 **Dextran ex. Leuconostoc Sp. - Tech 70 (65-75)**	5 Gms	3000
		25 Gms	4000
		100 Gms	7500
		500 Gms	25000
35270 [9004-54-0]	🌡 **Dextran ex. Leuconostoc Sp. - Tech 110 (100-120)**	5 Gms	3000
		25 Gms	4000
		100 Gms	7500
		500 Gms	25000

Code	Name	Packing	Price
62753 [9004-54-0]	🌡 **Dextran ex. Leuconostoc Sp. - Tech 250 (225-275)**	5 Gms	3000
		25 Gms	4000
		100 Gms	7500
		500 Gms	25000
19301 [9004-54-0]	🌡 **Dextran ex. Leuconostoc Sp. - Tech 500 (450-550)**	10 Gms	3000
		25 Gms	4000
		100 Gms	7500
		500 Gms	25000
28354 [87915-38-6]	■ **Dextran Blue 5 ex. Leuconostoc Sp.**	500 Mg	15000
82575 [87915-38-6]	■ **Dextran Blue 40 ex. Leuconostoc Sp.**	1 Gms	25000
57356 [87915-38-6]	■ **Dextran Blue 110 ex. Leuconostoc Sp.**	1 Gms	15000
14420 [87915-38-6]	■ **Dextran Blue 500 ex. Leuconostoc Sp.**	1 Gms	15000
67045 [87915-38-6]	■ **Dextran Blue 2000 ex. Leuconostoc Sp.**	1 Gms	11000
		5 Gms	50000
83826 [9015-73-0]	■ **DEAE-Dextran 500 ex. Leuconostoc Sp.**	1 Gms	9680
		5 Gms	24890
33434 [9064-91-9]	■ **DEAE-Dextran Hydrochloride 500 ex. Leuconostoc Sp.**	5 Gms	2000
		25 Gms	4500
		100 Gms	16000
62250 [60842-46-8]	🌡 **FITC-Dextran 4000 ex. Leuconostoc Sp.** (Fluorescein Isothiocyanate–Dextran)	100 Mg	12000
		250 Mg	24000
36454 [60842-46-8]	🌡 **FITC-Dextran 40000 ex. Leuconostoc Sp.** (Fluorescein Isothiocyanate–Dextran)	100 Mg	12000
		250 Mg	24000
92901 [60842-46-8]	🌡 **FITC-Dextran 70000 ex. Leuconostoc Sp.** (Fluorescein Isothiocyanate–Dextran)	100 Mg	12000
		250 Mg	24000
82221 [60842-46-8]	🌡 **FITC-Dextran 2000000 ex. Leuconostoc Sp.** (Fluorescein Isothiocyanate–Dextran)	100 Mg	12000
		250 Mg	24000
65163	● **Dextran Phenyl 40 (30-50)**	1 Gms	7310
		5 Gms	23970
76203 [9011-18-1]	🌡 **Dextran Sulphate Sodium Salt 30-50(DSS 35-50) ex. Leuconostoc Sp.Colitis Grade for molecular biology**	10 Gms	3500
		25 Gms	7000
		100 Gms	24000
		250 Gms	42000
99629 [9011-18-1]	🌡 **Dextran Sulphate Sodium Salt 500 (DSS 500) ex. Leuconostoc Sp. for molecular biology**	10 Gms	3800
		25 Gms	7500
		100 Gms	26000
		250 Gms	45000
17379 [3749-36-8]	🌡 **DHP Linker (3,4-Dihydro-2H-Pyran-2-Methanol) extrapure, 98%**	5 ml	4500
		25 ml	20000
27515 [325146-81-4]	🌡 **Diacetato [(R)-(+) -2,2`-bis(Diphenylphosphino) -1,1`-Binaphthyl]Ruthenium (II)) extrapure, 98%** ((R)-Ru(OAc)2(BINAP))	250 Mg	5000
		1 Gms	18000
91379 [123-42-2]	**Diacetone Alcohol extrapure, 98%** (4-Hydroxy-4-Methyl-2-Pentanone)	500 ml	400
		2500 ml	1800
34224 [582-52-5]	✔ **Diacetone-D-Glucose (DAG) extrapure, 98%** (1,2:5,6-Di-O-Isopropylidene-a-D-Glucofuranose, D-Glucose Diacetonide)	25 Gms	1500
		100 Gms	4500
		500 Gms	20000
72352 [57-71-6]	**Diacetylmonoxime extrapure AR, 99%** (2,3-Butanedione Monoxime)	25 Gms	420
		100 Gms	1300
		500 Gms	5500
10006 [6342-17-2]	✔ **Piperazine Diacrylamide extrapure, 99%** (1,4-(Diacryloyl) Piperazine)	1 Gms	2000
		5 Gms	3500
		25 Gms	12000
62736 [1482-98-0]	✔ **L-2,4 Diaminobutyric Acid Monohydrochloride extrapure, 98%**	1 Gms	2490
		5 Gms	11800
		25 Gms	45340
94524 [91-95-2]	🌡 **3,3-Diaminobenzidine (DAB) pure, 98%**	5 Gms	1300
		25 Gms	5200
		100 Gms	18000
17076 [868272-85-9]	🌡 **3,3'-Diaminobenzidine Tetrahydrochloride Hydrate (DAB.4HCl.xH2O) extrapure AR, 98%**	1 Gms	1310
		5 Gms	4500
		25 Gms	16000
98465	🌡 **3,3-Diaminobenzidine Tetrahydrochloride (DAB.4HCl) Buffer Substrate Solution for Peroxidase**	100 ml	7460
17758 [535-87-5]	**3,5-Diaminobenzoic Acid extrapure, 98%**	100 Gms	920
		500 Gms	4170
20723 [618-56-4]	**3,5-Diaminobenzoic Acid Dihydrochloride extrapure, 98%**	25 Gms	5510
16368 [125572-95-4]	**trans-1,2-Diaminocyclohexane-N,N,N,N Tetraacetic Acid Monohydrate (CDTA) ACS, 97%**	10 Gms	1410
		25 Gms	2600
		100 Gms	9500
		500 Gms	43000
85122 [137-09-7]	**2,4-Diaminophenol Dihydrochloride (Amidol) extrapure, 98%**	5 Gms	3850
		25 Gms	9320
34124 [1904-98-9]	✔ **2,6-Diaminopurine extrapure, 98%**	1 Gms	800
		5 Gms	3800
		25 Gms	15500
13261 [20325-40-0]	🌡 **o-Dianisidine Dihydrochloride extrapure, 98%** (3,3-Dimethoxybenzidine Dihydrochloride)	1 Gms	700
		5 Gms	2200
		25 Gms	10500
94912 [82-22-4]	✔ **1,1-Dianthrimide extrapure AR, 98%** (Di(1-anthraquinoyl)amine) (1,1-Iminodianthraquinone) (for determination of Boron)	1 Gms	3080
		5 Gms	13200
90405 [9001-18-7]	■ **Diaphorase ex. Clostridium Sp., 30U/mg powder**	500 Units	5340
99956 [50978-11-5]	**Diatrizoic Acid Dihydrate extrapure, 98%** (Sodium Diatrizoate Hydrate, Amidotrizoic Acid Dihydrate)	5 Gms	2000
		25 Gms	6000
61048 [50978-11-5]	**Diatrizoic Acid Dihydrate ExiPlus, Multi-Compendial, 98%** (Sodium Diatrizoate Hydrate, Amidotrizoic Acid Dihydrate)	1 Gms	1250
		5 Gms	3500
		25 Gms	8000

Code	Name	Packing	Price
23144 [103-49-1]	**Dibenzylamine pure, 98%**	100 ml 500 ml	720 3500
96661 [103-50-4]	**Dibenzyl Ether pure, 97%** (Benzyl Ether)	500 ml 2500 ml	1800 8000
88749 ✔ [106-37-6]	**1,4-Dibromobenzene pure, 98%**	500 Gms	1700
43589 [110-52-1]	**1,4-Dibromobutane pure, 98%**	250 ml	1200
74600 ✔ [77-48-5]	**1,3-Dibromo-5,5-Dimethylhydantoin (DBDMH) extrapure, 98%** (Dibromantin)	500 Gms	1200
33384 [106-93-4]	**1,2-Dibromoethane pure, 98%** (Ethylene Dibromide)	250 ml 1000 ml	1050 4050
31575 [626-41-5]	**3,5-Dibromophenol pure, 98%**	1 Gms 5 Gms	3080 13950
34890 [109-64-8]	**1,3-Dibromopropane pure, 99%** (Trimethylene Dibromide)	100 ml 500 ml	1130 4600
96763 🌡 [24424-99-5]	**Di-tert-Butyldicarbonate (BOC Anhydride, DiBOC) extrapure, 98.5 -101.5%**	100 Gms 500 Gms 5 Kg	640 2070 17250
69246 🌡 [16980-89-5]	**N6,2-O-Dibutyryl-Adenosine-3,5-Cyclic-Monophosphate Sodium Salt, 97%**	25 Mg 100 Mg	5600 13670
26414 🌡 [1918-00-9]	**DICAMBA for tissue culture, 95%** (3,6-Dichloro-2-Methoxybenzoic Acid)	250 Mg 1 Gms 5 Gms	1500 4000 9000
62578 🌡 [74124-79-1]	**Di-(N-Succinimidyl) Carbonate extrapure, 95%**	5 Gms 25 Gms 100 Gms	850 2800 7200
29154 [2234-16-4]	**2,4-Dichloroacetophenone pure, 98%**	100 Gms 500 Gms	700 2200
94453 [541-73-1]	**1,3-Dichlorobenzene extrapure, 99%** (m-Dichlorobenzene)	250 ml 1000 ml 2500 ml	700 2600 5800
91427 [50-45-3]	**2,3-Dichlorobenzoic Acid pure, 98%**	25 Gms 100 Gms	1200 4200
62122 [50-84-0]	**2,4-Dichlorobenzoic Acid extrapure, 98%**	100 Gms 500 Gms	350 1600
66481 [51-44-5]	**3,4-Dichlorobenzoic Acid pure, 98%**	25 Gms 100 Gms	2300 7000
40583 ■ [95408-45-0]	**Dichloro[1,1`-bis(di-t-Butylphosphino)Ferrocene]Palladium(II) extrapure, 98%** (1,1`-Bis(di-tert-Butylphosphino)Ferrocene)Dichloropalladium(II))	250 Mg 1 Gms 5 Gms	2500 6000 22000
33604 [19978-61-1]	**Dichloro(1,2-bis(Diphenylphosphino)Ethane)Palladium (II) extrapure, 98%** ([1,2-Bis(Diphenylphosphino)Ethane]Dichloropalladium(II))	250 Mg 1 Gms 5 Gms	3000 5800 24000
97111 [59831-02-6]	**Dichloro(1,3-bis(Diphenylphosphino)Propane)Palladium(II) extrapure, 98%** ([1,2-Bis(Diphenylphosphino)Ethane]Dichloropalladium(II))	250 Mg 1 Gms 5 Gms	3000 5800 24000
99494 [29964-62-3]	**Dichloro[1,4-bis(Diphenylphosphino)Butane]Palladium(II) extrapure, 98%**	250 Mg 1 Gms 5 Gms	3000 5800 24000
70264 [927706-57-8]	**Dichloro-[1,3-bis(Diisopropylphenyl)-2-Imidazolidinylidene]-(3-Chloropyridyl) Palladium(II), 98%** ((1,3-Bis(2,6-Diisopropylphenyl)Imidazolidene)(3-Chloropyridyl) Palladium(II) Dichloride)	250 Mg 1 Gms 5 Gms	2200 5000 18000
18602 🌡 [821-10-3]	**1,4-Dichloro-2-Butyne pure, 97%**	5 Gms 25 Gms	3200 11100
21942 🌡 [13965-03-02]	**Dichlorobis(Triphenylphosphine)Palladium(II) extrapure, 98%** (Pd(PPh2)2C)	1 Gms 5 Gms 25 Gms	5500 21000 95000
74660 □ [50982-12-2]	**Dichloro (Cycloocta-1,5-Diene) Ruthenium (II) pure, 35.50-36.50% Ru** ((1,5-Cyclooctadiene) Ruthenium (II) Chloride Polymer)	1 Gms 5 Gms	10000 50000
85048 ■ [4091-99-0]	**2,7-Dichlorofluorescein Diacetate (DCFH) extrapure, 97%** (DCFH-DA, 2,7-Dichlorodihydrofluorescein Diacetate)	25 Mg 100 Mg	2600 7200
86428 🌡 [84-58-2]	**2,3-Dichloro-5,6-Dicyano-1,4-Benzoquinone (DDQ) pure, 98%**	25 Gms 100 Gms 500 Gms	1510 6050 30000
52979 ✔ [4885-02-3]	**Dichloromethyl Methyl Ether pure, 97%**	5 ml 25 ml 100 ml	2200 5500 21000
88912 [89-61-2]	**2,5-Dichloronitrobenzene pure, 99%** (1,4-Dichloro-2-Nitrobenzene)	250 Gms	250
52039 ✔ [120-83-2]	**2,4-Dichlorophenol pure, 98%**	25 Gms 250 Gms	500 1500
87188 [87-65-0]	**2,6-Dichlorophenol extrapure, 99%**	25 Gms 100 Gms 500 Gms	400 1200 4000
28181 ■ [350-67-4]	**3,5-Dichlorophenyldiazonium Tetrafluroborate pure, 98%** (3,5-Dichlorobenzenediazonium Tetrafluoroborate)	5 Gms 25 Gms	12000 35000
42180 [68716-47-2]	**2,4-Dichlorophenylboronic Acid extrapure, 98%**	1 Gms 5 Gms 25 Gms	1070 2970 10380
93611 [151169-75-4]	**3,4-Dichlorophenylboronic Acid extrapure, 97%**	1 Gms 5 Gms 25 Gms	1000 2500 8900

Code	Name	Packing	Price
33017 [101-38-2]	**2,6-Dichloroquinone-4-Chlorimide (Gibbs Reagent) extrapure AR, 99%** (N,2,6-Trichloro-p-Benzoquinoneimide)	5 Gms 25 Gms 100 Gms	1610 6530 21350
38762 [2831-66-5]	**2,4-Dichloro-1,3,5-Triazine pure, 98%**	1 Gms 5 Gms	17650 58810
57262 [101-83-7]	**Dicyclohexylamine pure, 99%**	500 ml 2500 ml	830 3800
79596 [538-75-0]	**N,N-Dicyclohexyl Carbodiimide (DCC) pure, 98%**	100 Gms 500 Gms 5 Kg	500 1800 16000
95525 [4097-22-7]	**2,3-Dideoxyadenosine (ddA) extrapure, 98%**	5 Mg	2940
47248 [7481-89-2]	**2,3-Dideoxycytidine (ddC) extrapure, 98%**	25 Mg	3150
73273 [69655-05-6]	**2,3-Dideoxyinosine (ddI, ddIno) extrapure, 98%**	25 Mg	5200
22328 [3416-05-5]	**2,3-Dideoxythymidine (3-Deoxythymidine, ddT) extrapure, 98%**	25 Mg	6300
17412 [5983-09-5]	**2,3-Dideoxyuridine (ddU) extrapure, 98%**	25 Mg 100 Mg	8400 25200
16312 [165534-43-0]	**DEPBT (3-(Diethoxyphosphoryloxy) -1,2,3-benzotriazin-4(3H)-one) extrapure, 98%**	5 Gms 25 Gms	3680 10240
86283 [120-21-8]	**p-Diethylaminobenzaldhyde pure, 99%**	100 Gms	2610
19885 [869-24-9]	**2-Diethylaminoethyl Chloride Hydrochloride extrapure, 99%** (2-Chlorotriethylamine Hydrochloride)	25 Gms 100 Gms 500 Gms	800 1600 6000
33684 [10369-83-2]	**Diethyl Aminoethyl Hexanoate (DA-6) pure, 98%**	5 Gms 25 Gms	2200 7710
12915 [20605-01-0]	**Diethyl bis(hydroxymethyl)malonate pure, 97%** (Diethyl 2,2-bis(hydroxymethyl)malonate)	5 Gms 25 Gms 100 Gms	2400 6000 22000
79628 [814-49-3]	**Diethyl Chlorophosphate pure, 97%** (Chlorophosphoric Acid Diethyl Ester, Diethyl Phosphoryl Chloride)	25 Gms 100 Gms	3800 12000
56701 [298-06-6]	**O,O-Diethyl Dithiophosphate extrapure, 85%**	5 Gms 25 Gms 100 Gms	1800 4200 15000
64699 [67-43-6]	**Diethylenetriaminepentaacetic Acid (DTPA) extrapure AR, 99%** (Pentetic Acid)	100 Gms 500 Gms	490 1250
81346 [140-01-2]	**Diethylenetriaminepentaacetic Acid Pentasodium Salt (DTPA) extrapure, 40% (Liquid)**	500 ml 1000 ml	500 950
17598 [111-40-0]	**Diethylenetriamine pure, 98%**	500 ml 2500 ml	1100 5000
46791 [1609-47-8]	**Diethyl Pyrocarbonate (DEPC) for molecular biology, 99%**	5 ml 25 ml 100 ml	700 2200 8000
84254 [91-67-8]	**N,N-Diethyl-m-Toluidine pure, 98%**	100 ml 500 ml	1120 3470
99823 [367-25-9]	**2,4-Difluoroaniline extrapure, 99%**	25 Gms 100 Gms	1480 5270
29553 [372-18-9]	**1,3-Difluorobenzene (MDFB) extrapure, 99%** (m-Diflourobenzene)	25 ml 100 ml 500 ml	1100 2330 10930
72081 [1583-58-0]	**2,4-Difluorobenzoic Acid pure, 98%**	5 Gms 25 Gms	2490 8490
83488 [327-92-4]	**1,5-Difluoro-2,4-Dinitrobenzene (DFDNB) pure, 97%**	5 Gms 25 Gms	3080 12100
10670 [446-35-5]	**2,4-Difluoronitrobenzene extrapure, 99%**	25 Gms 100 Gms 500 Gms	2190 5210 23920
90309 [144025-03-6]	**2,4-Difluorophenylboronic Acid extrapure, 95%**	5 Gms 10 Gms	6010 10440
53755 [11024-24-1]	**Digitonin (Digitin) Reagent Grade**	100 Mg 1 Gms	3400 30000
48712 [20830-75-5]	**Digoxin extrapure, 95%**	50 Mg 1 Gms 5 Gms	600 3500 15000
91374 [28230-32-2]	**3,4-Dihydro-3-Hydroxy-4-Oxo-1,2,3-Benzotriazine (DHOBT, DHBT) extrapure, 98%**	5 Gms 25 Gms	1890 6510
23615 [5524-05-0]	**Dihydrocarvone pure, 97%** ((+)-Dihydrocarvone)	25 Gms 100 Gms	4800 12000
85408 [20777-49-5]	**Dihydrocarvyl Acetate (mix of Isomers) extrapure, 95%** ((-)-Dihydrocarvyl Acetate)	25 Gms 100 Gms	6000 16500
86580 [110-87-2]	**3,4-Dihydro-2H-pyran (3,4-DHP) pure, 97%** (2,3-Dihydropyran)	25 ml 100 ml 500 ml	450 1200 4500
65437 [89-84-9]	**2,4-Dihydroxyacetophenone extrapure, 99%** (Resacetophenone)	25 Gms 100 Gms	770 2450
41384 [490-78-8]	**2,5-Dihydroxyacetophenone extrapure, 99%** (2-Acetylhydroquinone, Quinacetophenone)	5 Gms 25 Gms	1380 5500
58494 [699-83-2]	**2,6-Dihydroxyacetophenone extrapure, 99%** (2-Acetylresorcin)	5 Gms 10 Gms 25 Gms 100 Gms	1570 2820 6180 21530
15395 [89-86-1]	**2,4-Dihydroxybenzoic Acid pure, 98%** (ß-Resocylic Acid)	25 Gms 100 Gms 500 Gms	1050 2060 7890
69911 [132-86-5]	**1,3-Dihydroxynaphthalene (Naphthoresorcinol) extrapure AR, 99%** (1,3-Naphthalenediol)	1 Gms 5 Gms 25 Gms	2750 11000 41800
41110 [83-56-7]	**1,5-Dihydroxynaphthalene pure, 98%**	100 Gms 500 Gms	4180 16500
65796 [575-38-2]	**1,7-Dihydroxynaphthalene pure, 97%** (1,7-Naphthalenediol)	5 Gms 25 Gms	6540 23490

Code	Name	Packing	Price
33580 [582-17-2]	2,7-Dihydroxynaphthalene pure, 98% (2,7-Naphthalenediol)	5 Gms 25 Gms 100 Gms	1430 2860 6320
63370 [116313-85-0]	3,4-Dihydroxy-5-Nitro Benzaldehyde pure, 98%	25 Gms 100 Gms	4290 10720
16517 [14894-18-9]	DL-Dihydrozeatin (DHZ) pure, 98% (6-(4-Hydroxy-3-Methylbutylamino)Purine)	10 Mg 25 Mg	4000 7500
19022 [83-73-8]	5,7-Diiodo-8-Quinolinol (Diiodohydroxyquinoline, Iodoquinol) ExiPlus, 96.5% (5,7-Diiodo-8-Hydroxyquinoline)	25 Gms 100 Gms 500 Gms	2750 7150 27170
43913 [693-13-0]	N,N-Diisopropyl Carbodiimide (DIPC) extrapure, 99%	25 Gms 100 Gms 500 Gms	900 3000 13000
95423	1,2,3,4-Di-O-Isopropylidene L-Arabinopyranoside extrapure	10 Gms	7610
94784 [4064-06-6]	1,2,3,4-Di-O-Isopropylidene a-D-Galactopyranoside extrapure, 98%	10 Gms	7610
10232 [304-55-2]	2,3-Dimercaptosuccinic Acid (DMSA, Succimer) extrapure, 98% (meso-2,3-Dimercaptosuccinic Acid)	1 Gms 5 Gms 25 Gms	1800 4000 12000
98526 [1072-71-5]	Bismuthiol (2,5-Dimercapto-1,3,4-Thiadiazole, DMTD) pure, 98% (1,3,4-Thiadiazole-2,5-Dithiol)	25 Gms 100 Gms 500 Gms	300 700 2500
76620 [1201-38-3]	2,5-Dimethoxyacetophenone pure, 99%	100 Gms 500 Gms	8790 34230
40824 [93-02-7]	2,5-Dimethoxybenzaldehyde (2,5-DMB) pure, 98%	25 Gms 100 Gms	1510 5780
24116 [613-45-6]	2,4-Dimethoxybenzaldehyde extrapure AR, 98%	25 Gms 100 Gms	2610 4920
82410 [150-78-7]	1,4-Dimethoxybenzene pure, 99% (Hydroquinone Dimethyl Ether)	250 Gms	2060
66053 [77-76-9]	2,2-Dimethoxypropane pure, 98%	500 ml 1000 ml 2500 ml	550 1080 2640
84774 [58-15-1]	4-Dimethylaminoantipyrine (Aminopyrine) extrapure, 99% (4-Dimethylamino-2,3-Dimethyl-1-Phenyl-3-Pyrazolin-5-one)	25 Gms 100 Gms	1850 6930
42442 [536-17-4]	5-(4-Dimethylaminobenzylidene) Rhodanine pure, 98%	1 Gms 5 Gms 25 Gms	1210 2640 7920
58451 [6203-18-5]	p-Dimethylaminocinnamaldehyde (DMACA) extrapure, 98%	5 Gms 25 Gms	2990 12710
48690 [6203-18-5]	p-Dimethylaminocinnamaldehyde (DMACA) Solution, 1% in HCl Used for flavonoid, proanthocyanidin or polyphenol detection	10 ml 50 ml	770 1650
58048 [108-01-0]	2-Dimethylaminoethanol extrapure, 98% (N,N-Dimethylethanolamine)	500 ml 2500 ml	640 2880
18814 [4584-46-7]	2-Dimethylaminoethyl Chloride Hydrochloride pure, 98% (2-Chloro-N,N-Dimethylethylamine Hydrochloride)	100 Gms 500 Gms	470 1740
49235 [25952-53-8]	1-(3-Dimethylaminopropyl)-3-Ethyl Carbodiimide Hydrochloride (EDC.HCl, EDAC.HCl) extrapure, 99% (N-Ethyl-N'-(3-Dimethylaminopropyl Ethylcarbodimide Hydrochloride)	25 Gms 100 Gms 500 Gms	750 2400 10000
93088 [1122-58-3]	4-Dimethylaminopyridine (DMAP) extrapure, 98%	25 Gms 100 Gms 500 Gms	400 900 3400
50593 [95-68-1]	2,4-Dimethylaniline pure, 98% (2,4-Xylidine)	250 ml 2500 ml	700 6050
33029 [87-62-7]	2,6-Dimethylaniline pure, 99% (2,6-Xylidine)	250 ml 1000 ml	940 3520
49579 [611-01-8]	2,4-Dimethylbenzoic Acid pure, 98% (m-Xylene-4-Carboxylic Acid, asym. -m-Xylylic Acid)	10 Gms 50 Gms	2310 9900
86330 [616-38-6]	Dimethyl Carbonate pure, 99%	500 ml 2500 ml 25 Ltr	280 1270 11500
21800 [75-78-5]	Dimethyl Dichlorosilane (DMDCS) pure, 98% (Dichlorodimethylsilane)	250 ml 500 ml 1000 ml	1150 2200 4000
10192 [1643-20-5]	N,N-Dimethyldodecylamine N-Oxide Solution (DDAO) pure, 30% (Lauryldimethylamine N-Oxide, LDAO)	100 ml 500 ml	1200 3000
58764 [1643-20-5]	N,N-Dimethyldodecylamine N-Oxide (DDAO) extrapure, 95% (Lauryldimethylamine N-Oxide, LDAO)	5 Gms 25 Gms	4300 18000
70976 [2491-06-7]	N,N-Dimethyl Glycine Hydrochloride pure, 98%	10 Gms 50 Gms	1600 6000
14157 [77-71-4]	5,5-Dimethylhydantoin (DMH) extrapure, 99% (4,4-Dimethyl-2,5-Dioxoimidazolidine, Fennosurf 300)	100 Gms 500 Gms	320 950
83177 [108-59-8]	Dimethyl Malonate pure, 98%	250 ml 500 ml	600 660
22247 [13138-48-2]	3,3-Dimethylnaphthidine extrapure AR, 98% (4,4'-Diamino-3,3'dimethyl-1,1'-binaphthyl)	200 Mg 1 Gms	4050 11090
35047 [13290-96-5]	Dimethyl 5-Nitroisophthalate pure, 99%	25 Gms 100 Gms	1500 3000
20923 [1123-49-5]	3,5-Dimethyl-4-Nitroisoxazole extrapure, 98%	5 Gms 25 Gms	3040 9710
98950 [95-87-4]	2,5-Dimethyl Phenol pure, 98% (2,5-Xylenol)	100 Gms	2190
83834 [17469-89-5]	N,N-Dimethyl-L-Phenylalanine extrapure, 98%	1 Gms 5 Gms	3000 9000

Code	Name	Packing	Price
93287 [99-98-9]	✔ **N,N-Dimethyl-p-Phenylenediamine (DMPPDA) pure, 98%** (4-(Dimethylamino)Aniline, 4-Amino-N,N-Dimethylaniline, N,N-Dimethyl-1,4-Phenylenediamine, DMPD)	5 Gms 25 Gms 100 Gms	1200 3250 11000
67318 [120-61-6]	✔ **Dimethyl Terephthalate pure, 99%**	500 Gms	800
91273 [518-67-2]	✔ **Dimidium Bromide extrapure, 98%** (3,8-Diamino-5-Methyl-6-Phenylphenan thridinium Bromide, Trypadine)	500 Mg 1 Gms 5 Gms	1500 2700 11000
84846 [97-02-9]	✔ **2,4-Dinitroaniline pure, 98%**	100 Gms 500 Gms	1250 5210
58586 [99-34-3]	**3,5-Dinitrobenzoic Acid extrapure, 98%**	100 Gms 500 Gms	620 2440
66764 [99-34-3]	**3,5-Dinitrobenzoic Acid extrapure AR, 99%**	100 Gms 500 Gms	800 3300
87464 [99-34-3]	**3,5-Dinitrobenzoic Acid extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	100 Gms 500 Gms	930 3580
56569 [14323-43-4]	**Dinitrodiaminepalladium (II) extrapure, 45-50% Pd**	250 Mg 1 Gms	5000 13000
88639 [119-26-6]	**2,4-Dinitrophenylhydrazine (DNPH) Reagent Grade, 99%** (Brady`s Reagent)	25 Gms 100 Gms 500 Gms	300 800 3800
85387 [119-26-6]	**2,4-Dinitrophenylhydrazine (DNPH) ExiPlus, Multi-Compendial, 99%** (Brady`s Reagent) meets compendial specs of BP, Ph.Eur	25 Gms 100 Gms	300 1000
75477	**2,4-Dinitrophenylhydrazine (DNPH) Reagent Solution** (2,4-DNP Reagent)	100 ml	350
10253 [609-99-4]	**3,5-Dinitrosalicylic Acid pure, 98%**	25 Gms 100 Gms 500 Gms	800 2000 8000
38707 [609-99-4]	**3,5-Dinitrosalicylic Acid extrapure AR, 99%**	25 Gms 100 Gms	1200 3000
13313 [609-99-4]	**3,5-Dinitrosalicylic Acid Reagent Solution**	100 ml	800
41131 [512-04-9]	**Diosgenin extrapure, 95%** ((25R)-5-Spirosten-3ß-ol, 3ß-Hydroxy-5-Spirostene)	1 Gms 5 Gms	1800 5300
37187 [102-04-5]	● **1,3-Diphenylacetone extrapure, 98%**	25 Gms 100 Gms 500 Gms	2750 7370 30360
54164 [86-29-3]	✔ **Diphenylacetonitrile pure, 98%**	100 Gms	1100
43511 [140-22-7]	● **1,5-Diphenylcarbazide ACS** (DPC, 1,5-Diphenylcarbohydrazide, sym.-Diphenylcarbazide)	25 Gms 100 Gms	2200 8000
57719 [91-01-0]	✔ **Diphenylcarbinol pure, 99%** (Benzhydrol)	100 Gms 500 Gms	550 1980
23957 [101-84-8]	**Diphenyl Ether pure, 98%** (Diphenyl Oxide)	500 ml 1000 ml	880 1650
90574 [536-46-9]	✔ **N,N-Dimethyl-p-Phenylenediamine Dihydrochloride (DMPPDA.2HCl) pure, 98%** (4-(Dimethylamino)aniline dihydrohchloride, 4-Amino-N,N-dimethylaniline dihydrochloride, N,N-Dimethy	10 Gms 50 Gms	4400 19250
66525 [1720-32-7]	🌡 **1,6-Diphenyl-1,3,5-Hexatriene (DPH) pure, 98%** (Dicinnamyl, 1,6-Diphenylhexatriene)	1 Gms 5 Gms	11000 50000
34747 [57-41-0]	✔ **5,5-Diphenylhydantoin (Phenytoin) pure, 98%**	100 Gms	1760
78190 [530-47-2]	✔ **Diphenylhydrazine Hydrochloride, 97%**	5 Gms 25 Gms	1980 5610
41855 [101-81-5]	🌡 **Diphenylmethane pure, 98%**	500 Gms	900
29128 [1898-66-4]	🌡 **2,2-Diphenyl-1-Picrylhydrazyl (DPPH) extrapure, 95%**	250 Mg 1 Gms 5 Gms	3500 6900 25800
38201 [120-46-7]	✔ **1,3-Diphenyl-1,3-Propanedione (Dibenzoylmethane) extrapure, 98%**	25 Gms 100 Gms	3520 7920
52196 [102-08-9]	✔ **Diphenyl Thiourea pure, 98%**	100 Gms 500 Gms	990 2860
12781 [68399-80-4]	**DIPSO Buffer extrapure, 99%** (2,2-bis(2-Hydroxyethyl) -3-Amino-2-Hydroxypropane Sulphonic Acid)	25 Gms 100 Gms	3800 13000
16722 [102783-62-0]	**DIPSO Buffer Sodium Salt extrapure, 98%**	10 Gms 50 Gms	3000 6000
93634 [119-91-5]	✔ **2,2-Biquinoline (2,2-Diquinolyl, Cuproine) extrapure AR, 99%** (2,2'-Diquinoline)	1 Gms 5 Gms	900 4000
32363 [69-78-3]	● **5,5'-Dithiobis(2-Nitro Benzoic Acid) extrapure (DTNB), 98%** (Ellmans Reagent)	1 Gms 5 Gms 25 Gms 100 Gms	1100 3200 13000 48000
27956 [6892-68-8]	● **Dithioerythritol (DTE) extrapure, 99%** (Clelands Other reagent)	1 Gms 5 Gms 25 Gms	2200 7150 28600
53384 [6892-68-8]	● **Dithioerythritol (DTE) for molecular biology, 99%** (Clelands Other reagent)	1 Gms 5 Gms	2750 8250
48086 [79-40-3]	● **Dithiooxamide (Rubeanic Acid) extrapure AR, 99%**	1 Gms 5 Gms 25 Gms	1100 3850 16500
84834 [3483-12-3]	● **DL-Dithiothreitol (DTT) extrapure, 98%** (Clelands Reagent)	1 Gms 5 Gms 25 Gms 100 Gms	580 2050 8000 30000
46220 [3483-12-3]	● **DL-Dithiothreitol (DTT) ExiPlus, Multi-Compendial, 98%** meets compendial specs of BP	5 Gms	2260

Code	Name	Packing	Price
17315 [3483-12-3]	● **DL-Dithiothreitol (DTT) for molecular biology, 98%** (Clelands Reagent)	1 Gms 5 Gms 25 Gms 100 Gms	690 2500 9000 33000
15248	🌡 **DL-Dithiothreitol 0.1M Solution (DTT 0.1M Solution)**	10 ml 100 ml	700 1600
16927	🌡 **DL-Dithiothreitol 1M Solution (DTT 1M Solution)**	10 ml 50 ml	1200 3000
97980 [60-10-6]	**Dithizone pure, 70%** (Diphenyl Thiocarbazone, Dithiozone)	5 Gms	650
89257 [60-10-6]	**Dithizone ACS, 85%** (Diphenyl Thiocarbazone, Dithiozone)	5 Gms 25 Gms 100 Gms	780 3400 12500
46141 [60-10-6]	**Dithizone ExiPlus, Multi-Compendial, 85%** (Diphenyl Thiocarbazone, Dithiozone) meets compendial specs of BP, Ph.Eur	5 Gms 25 Gms	850 3900
56266	**Dithizone Reagent**	100 ml	1200
37237 [77-77-0]	**Divinyl Sulfone (DVS) pure, 96%** (Vinyl Sulfone)	5 ml 25 ml	2100 6100
80704 [112-16-3]	ⓘ **Lauroyl Chloride extrapure, 99%** (Dodecanoyl Chloride)	100 Gms 250 Gms 500 Gms	1200 2600 4500
93089 [124-22-1]	✔ **Dodecylamine (Laurylamine) pure, 98%**	100 Gms 500 Gms	440 1650
85985 [219947-93-0]	**1-Dodecyl-3-Methylimidazolium Hexafluorophosphate (DMIM.PF6, C12MIMPF6) extrapure, 98%** (1-Lauryl-3-Methylimidazolium Hexafluorophosphate)	5 Gms 25 Gms	5170 19580
85755 [69227-93-6]	🌡 **n-Dodecyl-b-D-Maltoside (DDM) extrapure, 98%** (Dodecyl-Maltoside, n-Dodecyl-b-D-Maltopyranoside)	500 Mg 1 Gms 5 Gms	4620 8140 34650
24940 [1119-94-4]	**Dodecyltrimethylammonium Bromide (DOTAB) pure, 99%** (Lauryl Trimethylammonium Bromide)	100 Gms 500 Gms 5 Kg	900 3000 27000
53499 [1119-94-4]	**Dodecyltrimethylammonium Bromide (DOTAB) PrimeDx, 99%**	5 Kg 25 Kg	24000 110000
42020 [112-00-5]	**Dodecyltrimethylammonium Chloride (DOTAC) pure, 98%** (Lauryl Trimethylammonium Chloride)	100 Gms 500 Gms 5 Kg	900 3000 27000
61936 [112-00-5]	**Dodecyltrimethylammonium Chloride (DOTAC) PrimeDx, 98%**	5 Kg 25 Kg	24000 110000
77201 [59-92-7]	**L-Dopa extrapure, 99%** (L-Hydroxy Phenylalanine)	5 Gms 25 Gms	1000 3500
87943 [63-84-3]	**DL-Dopa extrapure, 99%** (DL-Dihydroxy Phenylalanine)	1 Gms 5 Gms 25 Gms	400 800 3200
45462 [62-31-7]	🌡 **Dopamine Hydrochloride (Dopamine HCl) extrapure, 98%**	5 Gms 25 Gms 100 Gms	1050 3680 13500
65094 [24390-14-5]	● **Doxycycline Hyclate (DXH), 800U/mg**	1 Gms 5 Gms 10 Gms	2400 3200 4500
78985 [7688-25-7]	✔ **DPPB (1,4-Bis(Diphenylphosphino)Butane) extrapure, 98%** (1,4-Butanediylbis(Diphenylphosphine) , Tetramethylenebis(Diphenylphosphine))	5 Gms 25 Gms	1900 7500
48237 [1663-45-2]	✔ **DPPE (1,2-Bis(Diphenylphosphino)Ethane) extrapure, 98%** (Diphos, Ethylenebis(Diphenylphosphine))	5 Gms 25 Gms	1100 5000
23660 [12150-46-8]	✔ **DPPF (1,1-Bis(Diphenylphosphino)Ferrocene) extrapure, 98%** (1,1-Ferrocenediyl-bis(Diphenylphosphine))	1 Gms 5 Gms 25 Gms	1100 3800 12900
34103 [6737-42-4]	✔ **DPPP (1,3-Bis(Diphenylphosphino)Propane) extrapure, 98%**	1 Gms 5 Gms 25 Gms	800 2300 10000
75470 [608-66-2]	🌡 **Dulcitol (D-Galactitol) extrapure, 99%** (D-Galactitol)	25 Gms 100 Gms	1600 5300
62168 [1308-87-8]	**Dysprosium Oxide extrapure, 99.99%**	5 Gms 25 Gms	8000 32000
62813 [16357-59-8]	🌡 **EEDQ extrapure, 99%** (N-Ethoxycarbonyl-2-Ethoxy-1,2-di-Hydroxyquinoline)	10 Gms 25 Gms 100 Gms	1200 2200 7000
17398	**Ehrlich`s Aldehyde Reagent**	100 ml	300
99640	**Ehrlich`s Amino Sugar Reagent**	100 ml	300
81802 [112-79-8]	■ **Elaidic Acid Reference Standard, 99%(GC)** (trans-9-Octadecenoic Acid, trans-Oleic Acid)	1 Gms 5 Gms	4000 16000
11814 [39445-21-1]	■ **Elastase (ELA, PPE) ex. Porcine Pancreas, 10U/mg Protein**	10 Mg	6500
19764 [155569-91-8]	**Emamectin Benzoate extrapure, 90%**	100 Mg 1 Gms	1800 4000
73363 [518-82-1]	● **Emodin ex. Polygonum Cuspidatum, 90%** (Emodin ex. Fallopia Japonica, 1,3,8-Trihydroxy-6-Methyl-9,10-Anthracenedione, 1,3,8-Tri-hydroxy-6-M	25 Mg 500 Mg 1 Gms 5 Gms	1500 3200 5000 12000
90310	**Emulgen B-66 (Eumulgin B66)** (Polyoxyethylene Alkyl Ether)	100 ml 500 ml	4300 16500
20790 [93106-60-6]	● **Enrofloxacin (EFX), 98%**	5 Gms	1800

Code	Name	Packing	Price
59508 [21087-74-1]	● **Epiprogoitrin Potassium Salt Reference Standard Grade, 80%** ((S)-2-Hydroxy-3-Butenylglucosinolate Potassium Salt)	10 Mg	27700
97857 [56390-09-1]	● **Epirubicin Hydrochloride (EPR.HCl), 97%**	10 Mg 25 Mg	3800 7400
49470 [16052-06-5]	**EPPS Buffer extrapure, 99.5%**	10 Gms 50 Gms 100 Gms	1860 5920 11720
50062 [55750-63-5]	■ **N-(epsilon-Maleimidocaproyloxy) Succinimide extrapure, 98%** (EMCS)	10 Mg 50 Mg	5400 9900
41120 [10025-75-9]	**Erbium (III) Chloride Hexahydrate ultrapure, 99.99%**	5 Gms 25 Gms	8000 18000
54788 [13476-05-6]	**Erbium (III) Nitrate Hexahydrate extrapure, 99.9%**	5 Gms 25 Gms	5000 12000
29831 [12061-16-4]	**Erbium (III) Oxide ultrapure, 99.99%**	5 Gms 25 Gms	5000 12000
12485 [139177-64-3]	**Erbium (III) Trifluoromethanesulfonate extrapure, 97%** (Erbium Triflate, Trifluoromethanesulfonic Acid Erbium Salt)	5 Gms 25 Gms	10000 40000
56488 [149-32-6]	✔ **Erythritol extrapure, 99%** (Erythrite, meso-Erythritol)	25 Gms 100 Gms 500 Gms 5 Kg	400 700 3000 28000
48535 [114-07-8]	🌡 **Erythromycin (E) for tissue culture**	1 Gms 5 Gms	400 1800
51114 [66778-17-4, 531-75-9]	🌡 **Esculin Sesquihydrate extrapure, 98%** (Aesculin)	5 Gms 25 Gms	2600 12000
77688 [1070-11-7]	✔ **Ethambutol Dihydrochloride (ETB.2HCl), 98%**	1 Gms 5 Gms	800 2100
89737 [540-63-6]	**1,2-Ethanedithiol extrapure, 98.5%**	25 ml 100 ml 500 ml	2100 7000 31000
16844 [31005-02-4]	✔ **7-Ethoxycoumarin extrapure, 98%**	25 Mg 100 Mg	1500 3500
23545 [578-54-1]	✔ **2-Ethylaniline pure, 98%**	100 ml 500 ml	1500 2700
27542 [667-27-6]	**Ethyl Bromodifluoroacetate extrapure, 98%** (Ethyl 2-Bromo-2,2-Difluoroacetate, Ethyl 2-Bromo-2,2-Difluoroethaneperoxoate, Ethyl Difluorobromoac	25 Gms 100 Gms	3200 8000
96669 [105-39-5]	ⓘ **Ethyl Chloroacetate pure, 98%**	500 ml 2500 ml	750 3500
95953 [535-13-7]	**Ethyl 2-Chloropropionate practical grade, 95%**	100 ml	1400
92052 [105-56-6]	**Ethyl Cyanoacetate extrapure, 98%**	500 ml 2500 ml	800 3500
69674 [126213-50-1]	🌡 **3,4-Ethylenedioxythiophene (EDOT) pure, 98%**	5 Gms 25 Gms	2200 4400
62858 [67-42-5]	🌡 **Ethylene Glycol-O-O-bis-(2-Aminoethyl) N,N,N,N-Tetraacetic Acid (EGTA) extrapure, 98%**	10 Gms 25 Gms 100 Gms	2400 4200 10500
78214 [149-57-5]	**2-Ethylhexanoic Acid extrapure, 99%**	250 ml 1000 ml	300 1000
32021 [104-76-7]	**2-Ethyl-1-Hexanol pure, 99%**	500 ml 2500 ml	400 1500
69647 [124-03-8]	**Ethylhexadecyldimethylammonium Bromide (EHDDMAB) extrapure, 98%** (Cetyldimethylethylammonium Bromide, Dimethylethylhexadecylammonium Bromide)	25 Gms 100 Gms 500 Gms	1200 3500 11000
71407 [103-11-7]	**2-Ethylhexyl Acrylate extrapure, 99%** (Acrylic acid 2-ethylhexyl Ester)	500 ml 2500 ml	770 2970
30484 [88122-99-0]	**Ethylhexyl Triazone extrapure, 98%** (2,4,6-Trianilino(p-Carbo-2-Ethylhexyl-1 -oxy)-1,3,5-Triazine)	100 Mg 1 Gms 5 Gms	2310 7920 11000
91018 [13224-99-2]	✔ **Ethylidene Glucose extrapure, 90%**	1 Gms 5 Gms	4730 12210
10674 [75-03-6]	**Ethyl Iodide extrapure, 99%**	25 ml 100 ml 250 ml	2640 9460 22220
78503 [128-53-0]	🌡 **N-Ethylmaleimide extrapure, 99%**	5 Gms 25 Gms 100 Gms	2310 10450 41580
81176 [128-53-0]	🌡 **N-Ethylmaleimide ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	5 Gms 25 Gms 100 Gms	2640 11550 46200
85108 [62-50-0]	✔ **Ethyl Methanesulphonate (EMS) extrapure, 98%**	10 Gms 100 Gms	1050 9350
27035 [62-50-0]	✔ **Ethyl Methanesulphonate (EMS) for tissue culture, 98%**	10 Gms	2370
38303 [623-53-0]	**Ethyl Methyl Carbonate (EMC) extrapure, 98%** (Carbonic Acid Ethyl Methyl Ester, Methyl Ethyl Carbonate)	10 ml 50 ml	2200 5500
15371 [65039-09-0]	✔ **1-Ethyl-3-Methylimidazolium Chloride (EMIM Cl) extrapure, 98%**	5 Gms 25 Gms 100 Gms	2400 4200 8400
10956 [342573-75-5]	✔ **1-Ethyl-3-Methylimidazolium Ethyl Sulfate (EMIM EtOSO3) extrapure, 99%**	5 Gms 25 Gms	2300 7200
34543 [155371-19-0]	✔ **1-Ethyl-3-Methylimidazolium Hexafluorophosphate (EMIM PF6) extrapure, 97%**	5 Gms 25 Gms	1600 6600
96150 [80432-05-9]	✔ **1-Ethyl-3-Methylimidazolium Tetrachloroaluminate (EMIM AlCl4) extrapure, 95%**	5 Gms 25 Gms	2600 9500

Code	Name	Packing	Price
60887 [143314-16-3]	✔ **1-Ethyl-3-Methylimidazolium Tetrafluoroborate (EMIM BF4) extrapure, 97%**	5 Gms 25 Gms	5000 6500
90329 [145022-44-2]	✔ **1-Ethyl-3-Methylimidazolium Trifluoromethanesulfonate (EMIM Otf), 98%** (1-Ethyl-3-Methylimidazolium Triflate)	1 Gms 5 Gms 25 Gms	1400 4200 12900
49544 [104-90-5]	**5-Ethyl-2-Methylpyridine pure, 97%**	100 ml 250 ml	950 1680
64770 [223436-99-5]	✔ **1-Ethyl-1-Methylpyrrolidinium Bis(trifluoromethylsulfonyl)imide (EMP TFSI) extrapure, 99%**	1 Gms 5 Gms 25 Gms	2200 3520 14300
70945 [111-62-6]	🌡 **Ethyl Oleate pure**	500 ml 2500 ml	720 3410
24489 [90002-36-1]	**2-Ethylphenylboronic Acid pure, 95%** (2-Ethylbenzeneboronic Acid, o-Ethylphenylboronic Acid)	1 Gms 5 Gms	2440 10400
42676 [5308-25-8]	✔ **1-Ethylpiperazine pure, 99%** (N-Ethylpiperazine)	50 Gms 250 Gms	880 2200
72424 [80-40-0]	**Ethyl p-Toluenesulfonate pure, 98%** (Ethyl Tosylate, p-Toluenesulfonic Acid Ethyl Ester)	25 ml 500 ml	850 7120
59863 [89797-68-2]	**5-(Ethylthio)-1H-Tetrazole (ETT Activator) extrapure, 99%**	5 Gms 25 Gms 100 Gms	2000 5000 12000
65770 [1530-32-1]	**Ethyltriphenylphosphonium Bromide pure, 97%**	25 Gms 100 Gms	590 1460
57056 [121-32-4]	**Ethyl Vanillin extrapure, 98%**	100 Gms 500 Gms	920 4120
24820 [109-92-2]	✔ **Ethyl Vinyl Ether (Stabilized) pure, 98%**	100 ml 500 ml 2500 ml	680 2130 10290
93050 [8000-48-4]	**Eucalyptus Oil extrapure, 60%**	250 ml	1340
69553 [97-53-0]	🌡 **Eugenol extrapure, 99%** (2-Methoxy-4-(2-Propenyl)phenol, 4-Allyl-2-Methoxyphenol, 4-Allylguaiacol)	100 ml 500 ml	1340 4840
98345 [13759-92-7]	**Europium (III) Chloride Hexahydrate ultrapure, 99.99%**	1 Gms 5 Gms 25 Gms	9600 24800 71600
97625 [10031-53-5]	**Europium (III) Nitrate Hexahydrate ultrapure, 99.99%**	5 Gms 25 Gms	26000 74800
44088 [1308-96-9]	**Europium (III) Oxide extrapure, 99.9%**	1 Gms 5 Gms	6000 21200
91883 [52093-25-1]	**Europium (III) Trifluoromethanesulfonate extrapure, 98%** (Europium Triflate, Trifluoromethanesulfonic Acid Europium Salt)	1 Gms 5 Gms	8800 22200
46650 [314-13-6]	**Evans Blue**	10 Gms 25 Gms 100 Gms	1200 2500 8000
33138	**Exton`s Reagent**	100 ml	260
66977 [14263-94-6]	🌡 **Fast Blue B Salt (DBB)** (o-Dianisidine bis(Diazotized) Zinc Double Salt, Azoic Diazo No. 48, Diazo Blue B)	25 Gms 100 Gms 500 Gms	700 2500 11000
94315 [5486-84-0]	🌡 **Fast Blue BB Salt Hemi Zinc Chloride Salt, 99%** (4-Amino-2,5-Diethoxybenzanilide Diazotated Zinc Double Salt, 4-Benzoylamino-2,5-Diethoxybenzenediaz	5 Gms 10 Gms 25 Gms	2410 4010 6010
62003 [14726-29-5]	🌡 **Fast Blue RR Hemi(Zinc Chloride) Salt** 4-Benzoylamino-2,5-Dimethoxybenzenediazonium Chloride Hemi(Zinc Chloride) salt, Azoic Diazo No. 24S	1 Gms 5 Gms	4000 18000
60339 [2353-45-9]	**Fast Green FCF (Food Green 3), 85%** (C.I. No. 42053)	5 Gms 25 Gms 100 Gms	440 1680 6000
79950 [89453-69-0]	**Fast Red TR Salt Hemi Zinc Chloride Salt** (4-Chloro-2-Methylbenzenediazonium Salt, Azoic Diazo No.11) (C.I. No. 37085)	100 Mg 500 Mg	2940 9000
77504 [79551-14-7]	✔ **Ferene Disodium Salt extrapure, 98%** (3-(2-Pyridyl)-5,6-di(2-Furyl) -1,2,4-Triazine-5,5-Disulfonic Acid Disodium Salt)	1 Gms 10 Gms	5500 33200
14977 [102-54-5]	✔ **Ferrocene extrapure, 98%**	25 Gms 100 Gms 500 Gms	550 1700 7000
38082 [69898-45-9]	✔ **Ferrozine Monosodium extrapure, 97%** (5,6-Diphenyl-3-(2-pyridyl) -1,2,4-triazine-4',4"-disulfonic acid sodium salt, PDT disulfonate)	1 Gms 5 Gms 25 Gms	2300 9900 46700
33920 [69898-45-9]	✔ **Ferrozine Monosodium PrimeDx, 97%**	100 Gms 250 Gms	195000 450000
60096 [1135-24-6]	🌡 **Ferulic Acid pure, 98%** (4-Hydroxy-3-Methoxycinnamic Acid)	5 Gms 25 Gms 100 Gms	800 2300 7000
23108 [9001-32-5]	🌡 **Fibrinogen ex. Bovine Plasma, 65% Protein** (Coagulation Factor I)	1 Gms 5 Gms	5000 10500
45460 [26873-85-8]	● **Ficoll 400® for molecular biology** (Polysucrose 400)	5 Gms 25 Gms 100 Gms	2000 8500 30000
87939 [146-14-5, 84366-81-4]	● **Flavin Adenine Dinucleotide Sodium Salt Hydrate (FAD Hydrate) extrapure, 95%**	25 Mg 100 Mg 250 Mg 1 Gms	2500 8000 13000 42000
69510 [525-82-6]	**Flavone extrapure**	25 Gms	1800
89465 [1343-88-0]	**Florisil 60-100 mesh for chromatography**	100 Gms 250 Gms 1 Kg	2750 5500 16000

Code	Name	Packing	Price
22239 [86386-73-4]	**Fluconazole (FLC), 98%**	1 Gms 5 Gms	800 1600
54083 ● [28920-43-6]	**9-Fluorenylmethyl Chloroformate (FMOC Chloride, FMOC-Cl) extrapure, 99%**	5 Gms 25 Gms 100 Gms 1 Kg	550 1800 5000 38000
29175 ● [24324-17-2]	**9-Fluorenylmethanol (FMOC-OH) extrapure, 98%** (9-Fluorenemethanol, 9-(Hydroxymethyl)fluorene)	25 Gms 100 Gms 1 Kg	1500 4000 30000
46068 ● [38183-12-9]	**Fluorescamine extrapure AR, 99%**	100 Mg 500 Mg 1 Gms	5000 10500 17000
85445 ● [38183-12-9]	**Fluorescamine extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	25 Mg 100 Mg 500 Mg	3000 5500 12000
50389 [2321-07-5]	**Fluorescein Free Acid (C.I. No. 45350:1), 96%** (Solvent Yellow 94)	25 Gms 100 Gms	270 900
69943 [596-09-8]	**Fluorescein Diacetate pure, 98%** (3,6-Diacetoxyfluoran, Di-O-Acetylfluorescein)	5 Gms 25 Gms 100 Gms	4800 12000 42000
73698 ● [3326-32-7]	**Fluorescein Isothiocyanate Isomer I (FITC) extrapure, 95%**	25 Mg 100 Mg 250 Mg 1 Gms	1800 2700 6000 18000
55091 [518-47-8]	**Fluorescein Sodium Salt (C.I. No. 45350)** (Acid Yellow 73)	25 Gms 100 Gms 1 Kg	320 1200 10500
92458 [403-42-9]	**4-Fluoroacetophenone pure, 98%**	25 Gms 100 Gms	1430 4500
57482 [348-54-9]	**2-Fluoroaniline pure, 98%**	100 ml 500 ml	1700 8000
61820 [372-19-0]	**3-Fluoroaniline pure, 98%**	25 ml 100 ml	3000 11500
62395 [371-40-4]	**4-Fluoroaniline pure, 98%**	100 ml 500 ml	1800 8400
55482 [321-28-8]	**2-Fluoroanisole pure, 98%**	25 ml 100 ml	3150 9980
45810 [456-49-5]	**3-Fluoroanisole pure, 98%**	25 Gms 100 Gms	2630 9450
51009 [459-60-9]	**4-Fluoroanisole extrapure, 99%**	25 Gms 100 Gms	1200 4500
90335 [459-57-4]	**4-Fluorobenzaldehyde pure, 98%**	100 ml 500 ml	3500 14000
88879 [462-06-6]	**Fluorobenzene pure, 99%**	100 ml 500 ml 1000 ml	1000 4000 7600
68427 [456-22-4]	**4-Fluorobenzoic Acid pure, 98%**	25 Gms 100 Gms	840 2310
94304 [1194-02-1]	**4-Fluorobenzonitrile pure, 98%**	25 ml 100 ml	2430 8450
13212 [140-75-0]	**4-Fluorobenzylamine extrapure, 99%**	5 Gms 25 Gms	840 2210
11821 🌡 [2341-22-2]	**5-Fluorocytidine (5-FCytd) extrapure, 99%**	1 Gms	4200
24245 🌡 [2022-85-7]	**5-Fluorocytosine (5-FCtys) extrapure, 98%** (5-FA, 4-Amino-5-Fluoro-2-Hydroxypyrimidine)	1 Gms 5 Gms	2100 6830
81015 🌡 [50-91-9]	**5-Fluoro-2-Deoxyuridine extrapure, 98%**	25 Mg 100 Mg 1 Gms	1680 5260 32440
67782 [70-34-8]	**1-Fluoro-2,4-Dinitrobenzene (FDNB) extrapure AR, 99%**	25 ml 100 ml 250 ml	1400 4800 10250
42199 [70-34-8]	**1-Fluoro-2,4-Dinitrobenzene (FDNB) for HPLC, 99%** (Sanger Reagent)	5 ml	8000
68632 [179897-94-0]	**5-Fluoro-2-Methoxyphenyl Boronic Acid extrapure, 95%**	5 Gms 10 Gms	5000 8000
78447 [350-46-9]	**1-Fluoro-4-Nitrobenzene extrapure, 99%**	100 Gms 500 Gms	1100 4300
27492 🌡 [3013-27-2]	**2-Fluoro-3-Nitrotoluene extrapure, 98%** (3-Fluoro-2-nitrotoluene)	1 Gms	4000
78510 ■ [703-95-7]	**5-Fluoroorotic Acid (5-FOA) for molecular biology, 98%**	250 Mg 1 Gms	3970 13050
50958 ■ [207291-81-4]	**5-Fluoroorotic Acid Hydrate (5-FOA Hydrate) for molecular biology, 98%**	250 Mg 1 Gms	4200 13610
11933 ■ [220141-70-8]	**5-Fluoroorotic Acid Monohydrate (5-FOA Monohydrate) for molecular biology, 98%**	250 Mg 1 Gms	4090 13390
66079 [367-12-4]	**2-Fluorophenol pure, 98%**	10 ml 100 ml	1910 9810
47982 [372-20-3]	**3-Fluorophenol pure, 98%**	10 ml 50 ml	3760 12710
41593 [371-41-5]	**4-Fluorophenol pure, 98%**	100 Gms 500 Gms	3760 14280
10573 [405-50-5]	**4-Fluorophenylacetic Acid pure, 98%**	5 Gms 25 Gms 100 Gms	1060 3110 9530
86567 🌡 [2629-55-2]	**o-Fluoro-DL-Phenylalanine extrapure, 99%**	1 Gms	7420
53933 🌡 [2629-54-1]	**m-Fluoro-DL-Phenylalanine extrapure, 98%**	500 Mg	4410
96044 🌡 [51-65-0]	**p-Fluoro-DL-Phenylalanine extrapure, 99%**	1 Gms 5 Gms	2800 13000

Code	Name	Packing	Price
11280 [1993-03-9]	**2-Fluorophenylboronic Acid extrapure, 97%**	1 Gms	980
		5 Gms	1990
		25 Gms	8710
54271 [768-35-4]	**3-Fluorophenylboronic Acid extrapure, 98%**	1 Gms	790
		5 Gms	1720
		25 Gms	7330
68733 [1765-93-1]	**4-Fluorophenylboronic Acid extrapure, 98%**	5 Gms	2000
		25 Gms	5000
		100 Gms	17000
37042 [352-32-9]	**4-Fluorotoluene pure, 99%**	100 ml	3140
21015 [446-11-7]	**4-Fluoro-3-Nitrotoluene extrapure, 98%**	5 Gms	2200
		25 Gms	7700
15088 [51-21-8]	**5-Fluorouracil (5-FU) extrapure, 99%**	1 Gms	530
		5 Gms	1680
		25 Gms	5670
59333 [51-21-8]	**5-Fluorouracil (5-FU) Exiplus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	5 Gms	3300
		25 Gms	9000
28329 [316-46-1]	**5-Fluorouridine (5-FUrd) extrapure, 98%**	1 Gms	4200
		5 Gms	9980
53061	● **FMOC-L-Amino Acids Kit**	1 Kit	27080
97306 [86060-81-3]	● **FMOC-S-Acetamido Methyl-L-Cysteine (Fmoc-Cys(Acm) -OH) extrapure, 98%**	5 Gms	2000
		25 Gms	9000
77535 [79990-15-1]	● **FMOC-D-Alanine extrapure, 99%**	5 Gms	2470
		25 Gms	8000
33213 [35661-39-3]	● **FMOC-L-Alanine extrapure, 99%**	5 Gms	520
		25 Gms	2300
		100 Gms	8450
81433 [135112-27-5]	● **FMOC-2-Aminobutric Acid extrapure, 99%**	1 Gms	3650
		5 Gms	16610
66725 [130752-32-8]	● **FMOC-D-Arginine extrapure, 99%**	1 Gms	2180
		5 Gms	9670
		25 Gms	43500
88459 [91000-69-0]	● **FMOC-L-Arginine extrapure, 98%**	1 Gms	430
		5 Gms	1470
		25 Gms	6180
		100 Gms	22200
23508 [108321-39-7]	● **FMOC-D-Asparagine extrapure, 99%**	1 Gms	720
		5 Gms	3000
		25 Gms	12500
66815 [71989-16-7]	● **FMOC-L-Asparagine extrapure, 99%**	5 Gms	520
		25 Gms	2300
		100 Gms	8450
49323 [132388-59-1]	● **FMOC-N-Trityl-L-Asparagine (FMOC-Asn(Trt)-OH) extrapure, 99%**	5 Gms	3510
		100 Gms	22200
30074 [119062-05-4]	● **FMOC-L-Aspartic Acid extrapure, 99%**	1 Gms	6850
		5 Gms	12330
59022 [86060-84-6]	● **FMOC-L-Aspartic Acid-4- Benzy Ester extrapure, 98%**	10 Gms	11410
		25 Gms	19950
24407 [53298-33-2]	● **FMOC-S-Benzyl-L-Cysteine extrapure, 98%**	10 Gms	12150
		25 Gms	21260
52208 [83792-48-7]	● **FMOC-O-Benzyl-L-Serine extrapure, 98%** (FMOC-Ser(Bzl)-OH)	10 Gms	9000
		25 Gms	17000
31407 [71989-40-7]	● **FMOC-O-Benzyl-L-Tyrosine extrapure, 98%**	10 Gms	11410
		25 Gms	19970
22303 [200344-33-8]	● **FMOC-D-Citrulline extrapure, 99%**	100 Mg	3060
		1 Gms	4550
61458 [133174-15-9]	● **FMOC-L-Citrulline extrapure, 99%** (FMOC-CIT-OH)	1 Gms	1790
		5 Gms	6300
65668 [135248-89-4]	● **FMOC-L-Cysteine extrapure, 95%**	25 Mg	4490
		100 Mg	14940
86254 [121343-82-6]	● **FMOC-L-Glutamic Acid extrapure, 99%**	10 Gms	4990
		25 Gms	9840
68902 [123639-61-2]	● **FMOC -L-Glutamic Acid g-Benzyl Ester extrapure, 98%**	10 Gms	11540
		25 Gms	20190
90673 [71989-18-9]	● **FMOC-L-Glutamic Acid 5-Tert-Butyl Ester (FMOC-Glu(OtBu)-OH) extrapure, 99%**	10 Gms	4740
		25 Gms	11830
		100 Gms	22410
21323 [104091-08-9]	● **FMOC-D-Glutamic Acid 5-Tert-Butyl Ester (FMOC-D-Glu(OtBu)-OH) extrapure, 99%**	5 Gms	3000
		25 Gms	12000
59874 [112898-00-7]	● **FMOC-D-Glutamine extrapure, 99%**	500 Mg	630
		2.5 Gms	2520
		25 Gms	22360
27631 [71989-20-3]	● **FMOC-L-Glutamine extrapure, 99%**	25 Gms	1800
		100 Gms	6000
21282 [132327-80-1]	● **Na-FMOC-No-Trityl-L-Glutamine (FMOC-Gln(TRT)-OH) extrapure, 98%**	5 Gms	2500
		100 Gms	18000
18680 [29022-11-5]	● **FMOC-Glycine extrapure, 99%**	5 Gms	840
		25 Gms	2730
17849 [116611-64-4]	● **FMOC-L-Histidine extrapure, 98%**	5 Gms	3860
90688 [143688-83-9]	● **FMOC-D-Isoleucine extrapure, 99%**	500 Mg	17160
22834 [71989-23-6]	● **FMOC-L-Isoleucine extrapure, 99%**	5 Gms	630
		25 Gms	2610
		100 Gms	9200
16515 [114360-54-2]	● **FMOC-D-Leucine extrapure, 99%**	1 Gms	630
		5 Gms	2520
		25 Gms	10490
99284 [35661-60-0]	● **FMOC-L-Leucine extrapure, 99%**	25 Gms	2300
		100 Gms	6490

Code	Name	Packing	Price
20340 [71989-26-9]	● **Na-FMOC-Ne-BOC-L-Lysine (FMOC-Lys(BOC)-OH) extrapure, 98%**	10 Gms 25 Gms 100 Gms	3100 7350 13130
67377 [86060-82-4]	● **FMOC-N-e-Z-L-Lysine extrapure, 98%**	10 Gms 25 Gms	14260 24240
77712 [105047-45-8]	● **FMOC-L-Lysine extrapure, 99%**	1 Gms 5 Gms	1410 5380
23694 [112883-40-6]	● **FMOC-D-Methionine extrapure, 99%**	1 Gms 5 Gms 25 Gms	1040 4330 18030
50606 [71989-28-1]	● **FMOC-L-Methionine extrapure, 99%**	5 Gms 25 Gms 100 Gms	520 2300 8450
61215 [167393-62-6]	🌡 **N-a-Fmoc-N-e-4-Methyltrityl-L-Lysine (FMOC-Lys(MTT)-OH) extrapure, 98%**	1 Gms 5 Gms	3000 11000
95562 [112883-41-7]	● **FMOC-D-Norleucine extrapure, 99%**	1 Gms 5 Gms	3810 17330
35941 [77284-32-3]	● **FMOC-L-Norleucine extrapure, 99%**	1 Gms 5 Gms	2700 12260
73648 [144701-24-6]	● **FMOC-D-Norvaline extrapure, 99%**	1 Gms 5 Gms	3810 17330
31631 [135112-28-6]	● **FMOC-L-Norvaline extrapure, 99%**	1 Gms 5 Gms	3020 13730
24626 [95753-55-2]	● **FMOC-4-Nitro-L-Phenylalanine extrapure (FMOC-L-Phe(4NO2) -OH), 98%**	10 Gms 25 Gms	8550 14950
47450 [86123-10-6]	● **FMOC-D-Phenylalanine extrapure, 99%**	1 Gms 5 Gms 25 Gms 100 Gms	630 2560 8950 32200
39061 [35661-40-6]	● **FMOC-L-Phenylalanine extrapure, 99%**	5 Gms 25 Gms	600 2000
54836 [111524-95-9]	● **FMOC-D-Phenylglycine extrapure, 98%**	25 Gms	5000
33014 [102410-65-1]	● **FMOC-L-Phenylglycine extrapure, 98%**	10 Gms 25 Gms	1800 3000
53412 [101555-62-8]	● **FMOC-D-Proline extrapure, 99%**	5 Gms	4330
12047 [71989-31-6]	● **FMOC-L-Proline extrapure, 99%**	5 Gms 25 Gms 100 Gms	630 2610 9200
53534 [116861-26-8]	● **FMOC-D-Serine extrapure, 99%**	25 Gms	5000
18126 [73724-45-5]	● **FMOC-L-Serine extrapure, 99%**	1 Gms 5 Gms	1490 3230
13334 [71989-33-8]	● **FMOC-O-Tert-Butyl-L-Serine (FMOC-Ser(tBu)-OH) extrapure, 99%**	10 Gms 100 Gms	3000 16000
56859 [157355-81-2]	● **FMOC-D-Threonine extrapure, 98%**	1 Gms 5 Gms	3500 13000
64712 [73731-37-0]	● **FMOC-L-Threonine extrapure, 99%**	5 Gms 25 Gms 100 Gms	960 4030 14390
36877 [109425-51-6]	● **Na-FMOC-N(im)-Trityl-L-Histidine (FMOC-His(Trt)-OH) extrapure, 99%**	10 Gms 100 Gms	3000 16000
38217 [138797-71-4]	● **FMOC-O-Tert-Butyl-D-Threonine (FMOC-D-Thr(tBu)-OH) extrapure, 98%**	1 Gms 5 Gms	1500 6000
42366 [86123-11-7]	● **FMOC-D-Tryptophan extrapure, 99%**	5 Gms 25 Gms	2200 8000
70309 [35737-15-6]	● **FMOC-L-Tryptophan extrapure, 99%**	5 Gms 25 Gms	1000 2630
96532 [163619-04-3]	● **Na-FMOC-N(in)-BOC-D-Tryptophan (FMOC-D-Trp(BOC)-OH) extrapure, 98%**	1 Gms 5 Gms	2200 8000
54034 [112883-29-1]	● **FMOC-D-Tyrosine extrapure, 99%**	1 Gms 5 Gms 25 Gms	1040 4330 18030
79199 [92954-90-0]	● **FMOC-L-Tyrosine extrapure, 99%**	5 Gms 25 Gms 100 Gms	1430 6000 21160
22292 [71989-38-3]	● **FMOC-O-Tert-Butyl-L-Tyrosine (FMOC-Tyr(tBu)-OH) extrapure, 98%**	5 Gms 100 Gms	2110 16000
15917 [84624-17-9]	● **FMOC-D-Valine extrapure, 99%**	1 Gms 5 Gms 25 Gms	1360 5620 23410
47834 [68858-20-8]	● **FMOC-L-Valine extrapure, 99%**	5 Gms 25 Gms 100 Gms	410 1870 7010
79235 [82911-69-1]	● **FMOC-OSu extrapure AR, 99%**	5 Gms 25 Gms 100 Gms 500 Gms	650 2400 7250 33000
33866 [5025-82-1/54323-50-1]	✔ **Folcisteine (ATCA, NATCA) pure, 99%** (3-Acetyl-1,3-Thiazolidine-4-Carboxylic acid, R-Aminofol)	5 Gms	2000
43164 [59-30-3]	✔ **Folic Acid pure, 98%** (Vitamin B9, Vitamin M, Pteroyl-L-Glutamic acid)	25 Gms 100 Gms 500 Gms	800 2800 12400
24791 [59-30-3]	✔ **Folic Acid ExiPlus, Multi-Compendial, 98%** (Vitamin B9, Vitamin M, Pteroyl-L-Glutamic acid) meets compendial specs of BP	10 Gms 25 Gms	500 1100
17989 [59-30-3]	✔ **Folic Acid for tissue culture, 98%** (Vitamin B9, Vitamin M, Pteroyl-L-Glutamic acid)	25 Gms 100 Gms 500 Gms	1100 3500 15200
39520	**Folin & Ciocalteus Phenol (FCP) Reagent AR**	100 ml 250 ml 500 ml	600 1400 2500

Code	Name	Packing	Price
64684	**Folin's Uric Acid Reagent**	250 ml	590
49244 ✔ [68157-60-8]	**Forchlorfenuron (KT-30, CPPU, N-(2-Chloro-4-pyridyl)-N-phenylurea) extrapure, 99%**	1 Gms 5 Gms	1500 5000
53870 [87199-17-5]	**4-Formylphenylboronic Acid extrapure, 95%**	5 Gms 25 Gms	1550 6200
85925	**Fouchet`s Reagent**	100 ml	250
54704 ● [81028-91-3]	**D-Fructose-1,6-Diphosphate Octahydrate Trisodium Salt extrapure, 99%**	250 Mg 1 Gms	500 1600
25760 ● [26177-86-6]	**D-Fructose-6-Phosphate Disodium Salt extrapure, 95%**	100 Mg 500 Mg	3170 13000
73212 🌡 [2438-80-4]	**L-Fucose extrapure, 99%** (L-Fuc, 6-Deoxy-L-Galactose)	1 Gms 5 Gms 25 Gms	2000 7000 30000
76007 🌡 [3615-37-0]	**D-(+)-Fucose extrapure** (6-Deoxy-L-Galactose, Rhodeose)	500 Mg 1 Gms	7200 12500
47462 ⓘ [98-01-1]	**Furfural pure, 98%**	500 ml 2500 ml	700 2800
86133 ⓘ [98-01-1]	**Furfural ACS, 98%**	500 ml	1130
11049 ⓘ [617-89-0]	**Furfurylamine pure, 98%**	500 ml 2500 ml	1480 6700
18641 [54-31-9]	**Furosemide (Frusemide), 98.5%**	5 Gms 25 Gms	3560 9930
81309 ✔ [88-14-2]	**2-Furoic Acid pure, 98%** (Furan-2-Carboxylic Acid)	100 Gms 500 Gms	1000 3800
81294 ✔ [488-93-7]	**3-Furoic Acid pure, 98%** (Furan-3-Carboxylic Acid)	5 Gms 25 Gms	2500 9500
99124 ✔ [527-69-5]	**2-Furoyl Chloride pure, 98%** (Furan-2-Carbonyl Chloride)	20 ml 100 ml	1500 5000
90068 ✔ [26214-65-3]	**3-Furoyl Chloride pure, 97%** (Furan-3-Carbonyl Chloride)	1 Gms 5 Gms	4000 14500
79108 ✔ [539-47-9]	**Furylacrylic Acid pure, 99%** (3-(2-Furyl)acrylic Acid, Furan-2-Acrylic Acid, Furfurylideneacetic Acid)	25 Gms	2020
58327 🌡 [108321-42-2]	**G-418 Sulphate (Geneticin) extrapure**	100 Mg 1 Gms 5 Gms	1260 6300 29400
87858 🌡 [108321-42-2]	**G-418 Sulphate (Geneticin) for tissue culture**	100 Mg 1 Gms 5 Gms	2220 8060 36510
45734	**Gaby-Hadley Reagent A**	100 ml	300
94536	**Gaby-Hadley Reagent B**	100 ml	430
45781 [13450-84-5]	**Gadolinium (III) Chloride Hexahydrate ultrapure, 99.99%**	5 Gms 25 Gms	5000 12000
10438 [19598-90-4]	**Gadolinium (III) Nitrate Hexahydrate ultrapure, 99.99%**	5 Gms 25 Gms	10400 21200
31353 [12064-62-9]	**Gadolinium (III) Oxide extrapure, 99.9%**	5 Gms 25 Gms 100 Gms	920 3360 11900
52480 [52093-29-5]	**Gadolinium (III) Trifluoromethanesulfonate extrapure, 98%** (Gadolinium(III) Triflate, Trifluoromethanesulfonic Acid Gadolinium(III) Salt)	1 Gms 5 Gms	6000 12000
14725 [1772-03-8]	**D-Galactosamine Hydrochloride extrapure, 98%**	1 Gms 5 Gms 25 Gms	800 3000 12000
88176 [59-23-4]	**D-Galactose extrapure, 98%**	25 Gms 100 Gms 500 Gms	500 1500 6200
65528 [59-23-4]	**D-Galactose ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP	100 Gms 500 Gms	2100 8000
27146 ■ [9031-11-2]	**B-Galactosidase (BGAL) ex. E. Coli, 500U/mg powder** (Lactase, B-D-Galactoside Galactohydrolase)	5 K.Units	32000
20096 [4163-60-4]	**b-D-Galactose Pentaacetate extrapure, 99%**	25 Gms	8800
83055 ■ [1953-04-4]	**Galantamine Hydrobromide extrapure, 98%**	250 Mg	7700
14101 [69365-72-6]	**Gallium (III) Nitrate Hydrate extrapure, 99.9%**	5 Gms 25 Gms	4250 16160
93956 🌡 [63699-78-5]	**L-Gamma-Glutamyl-3-Carboxy-4-Nitroanilide (Glupa) extrapure, 99%** (L-Glutamic acid 5-(3-Carboxy-4-Nitroanilide) Ammonium Salt)	1 Gms 5 Gms	8800 37400
23813 🌡 [30562-34-6]	**Geldanamycin (GLD) ex. Streptomyces Hygroscopicus, 98%**	5 Mg 25 Mg	6600 24200
72394 [71010-52-1]	**Gelrite Gellan Gum (Low Acyl)**	100 Gms 250 Gms 500 Gms 1 Kg 5 Kg	2010 4200 6000 11000 52000
59431 [71010-52-1]	**Gelrite Gellan Gum (High Acyl)**	100 Gms 500 Gms 1 Kg 5 Kg	1800 5000 10000 48000
53409 🌡 [6902-77-8]	**Genipin extrapure, 98%**	25 Mg 100 Mg	3000 8000
37636 🌡 [1405-41-0]	**Gentamicin Sulphate (GEN)**	1 Gms 5 Gms 25 Gms	410 1640 6080
50827 🌡 [1405-41-0]	**Gentamicin Sulphate (GEN) for tissue culture**	1 Gms 5 Gms 25 Gms	960 2400 7890

Code	Name	Packing	Price
95677 [554-91-6]	**b-D-Gentiobiose extrapure, 95%** (Amygdalose)	100 Mg 500 Mg	3500 8000
56397 [4562-36-1]	**Gitoxin extrapure, 98%**	25 Mg	12000
38490 [5996-06-5]	**Glaucine Hydrobromide, 97%**	10 Gms	6000
76052 [9074-98-0]	**b-Glucanase ex. Aspergillus Niger, 10U/mg powder**	5 Gms 25 Gms	900 4000
40193 [21087-78-5]	**Glucobarbarin Potassium Salt Reference Standard Grade, 90%**	10 Mg	47250
63878 [15592-36-6]	**Glucocheirolin Potassium Salt Reference Standard Grade, 95%** (3-Methylsulfonylpropyl Glucosinolate Potassium Salt)	10 Mg	57750
61404 [15592-37-7]	**Glucoerucin Potassium Salt Reference Standard Grade, 95%** (4-Methylthiobutyl Glucosinolate Potassium Salt)	10 Mg	58800
20378 [15592-34-4]	**Glucoiberin Potassium Salt Reference Standard, 85%** (3-Methylsulfinylpropyl Glucosinolate Potassium Salt)	10 Mg	35700
84553 [245550-57-6]	**Gluconapin Potassium Salt Reference Standard, 90%** (3-Butenyl Glucosinolate Potassium Salt)	10 Mg	35700
27429 [18425-76-8]	**Gluconasturtiin Potassium Salt Reference Standard, 95%** (2-Phenylethyl Glucosinolate Potassium Salt)	10 Mg	44100
45063 [1198-69-2]	**D-Glucono-1,4-Lactone pure, 85%**	50 Mg	7140
57400 [52153-09-0]	**L-Glucono-1,5-Lactone pure, 85%** (3,4,5-Trihydroxy-6-(Hydroxymethyl)Tetrahydro-2H-Pyran-2 -one, L-Gluconolactone)	25 Mg	6510
55195 [90-80-2]	**D-Glucono-1,5-Lactone extrapure, 99%** (D-(+)-Gluconic Acid d-Lactone, D-Gluconolactone, Glucono-Delta-Lactone)	100 Gms 500 Gms	1580 2420
74587 [21414-41-5]	**Glucoraphanin Potassium Salt Reference Standard Grade, 90%** (4-Methylsulfinylbutyl Glucosinolate Potassium Salt)	10 Mg	33600
15546 [108844-81-1]	**Glucoraphenin Potassium Salt Reference Standard Grade, 95%** (4-Methylsulfinyl 3-Butenylglucosinolate Potassium Salt)	10 Mg	36750
73360 [32449-92-6]	**D-Glucorono-3,6-Lactone extrapure, 99%**	25 Gms 100 Gms 500 Gms	2310 3460 10290
28453 [9032-08-0]	**Glucoamylase ex. Rhizopus Sp., 30U/mg powder**	1000 Units 5000 Units	7050 30000
45648 [66-84-2]	**Glucosamine Hydrochloride extrapure, 99%**	100 Gms 500 Gms 5 Kg	900 3200 25000
49531 [66-84-2]	**Glucosamine Hydrochloride Exiplus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms 5 Kg	1300 3800 32000
14810 [9028-53-9]	**Glucose Dehydrogenase ex. Microorganism, 250U/mg powder**	5 K.Units	15500
61788 [9001-37-0]	**Glucose Oxidase (GOD) ex. Aspergillus Niger, 100U/mg powder**	10 K.Units 25 K.Units	2170 4870
12419 [604-68-2]	**a-D-Glucose Pentaacetate extrapure, 99%**	25 Gms	1260
65776 [604-69-3]	**b-D-Glucose Pentaacetate extrapure, 99%**	25 Gms 100 Gms	1080 3570
23208 [6736-77-2 (anhy)]	**D-Glucose-1-Phosphate Dipotassium Salt Hydrate (G1P.K2 Hydrate) extrapure, 98%** (Cori Ester Dipotassium Salt Hydrate)	1 Gms 5 Gms	1800 8200
62591 [150399-99-8]	**D-Glucose-1-Phosphate Disodium Salt Tetrahydrate (G1P.Na2 Tetrahydrate) extrapure, 98%** (Cori Ester Disodium Salt Tetrahydrate)	1 Gms 5 Gms 25 Gms	700 2600 11700
32994 [3671-99-6]	**D-Glucose-6-Phosphate Disodium Salt Dihydrate (G6P.Na2 Dihydrate) extrapure, 99%** (Robinson Ester)	1 Gms 5 Gms 25 Gms	1300 5300 25200
71037 [3671-99-6]	**D-Glucose-6-Phosphate Disodium Salt Dihydrate (G6P.Na2 Dihydrate) PrimeDx, 99%**	100 Gms	105600
11997 [9001-40-5]	**r-Glucose-6-Phosphate Dehydrogenase (rG6PD) ex. Leuconostoc Mesenteroides, 400U/mg powder**	5 K.Units 25 K.Units	5500 20000
75551 [9001-42-7]	**a-Glucosidase (Maltase) ex. Yeast, 100U/mg powder**	10 Mg 100 Mg	3210 29050
13275 [9001-22-3]	**B-Glucosidase ex. Sweet Almonds (Type 2A), 10U/mg powder** (Gentiobiase, Cellobiase)	250 Units 1000 Units	4000 9000
41369 [9001-22-3]	**B-Glucosidase ex. Sweet Almonds (Type 2B), 10U/mg powder** (Gentiobiase, Cellobiase)	25 Mg 100 Mg 1 Gms	2100 6300 23100
70077 [5115-71-9]	**Glucotropaeolin Potassium Salt Reference Standard Grade, 95%**	10 Mg	18900
88229 [6556-12-3]	**D-Glucuronic Acid extrapure, 98%**	1 Gms 5 Gms 25 Gms	1400 3300 11000
74117 [9029-12-3]	**Glutamate Dehydrogenase (L-GLDH) ex. Bovine Liver, 10U/mg powder**	100 Mg	8600
61126 [6893-26-1]	**D-Glutamic Acid extrapure CHR, 99%**	25 Gms 100 Gms 250 Gms	1700 4500 9000
20655 [72669-53-5]	**g-L-Glutamic Acid 7-Amido-4-Methylcoumarin extrapure, 98%** (y-L-Glutamic Acid-AMC)	10 Mg 25 Mg 100 Mg	7200 13900 39700

Code	Name	Packing	Price
76228 [1676-73-9]	● **L-Glutamic Acid r-Benzyl Ester extrapure, 99%** (5-Benzyl-L-Glutamate)	5 Gms 25 Gms 100 Gms	1300 2700 8400
98087 [9000-97-9]	■ **Glutamic Oxaloacetic Transaminase (GOT) ex. Porcine Heart, 50-300U/mg powder**	1 K.Units	14800
17632 [5959-95-5]	**D-Glutamine extrapure, 99%**	1 Gms 5 Gms	3200 5800
12381 [56-85-9]	**L-Glutamine extrapure CHR, 99%**	100 Gms 500 Gms 1 Kg	600 1800 3400
86612 [56-85-9]	**L-Glutamine ExiPlus, Multi-Compendial, 99%** ((S)-2,5-Diamino-5-Oxopentanoic acid, L-Glutamic acid 5-amide, Levoglutamide) meets compendial specs of USP, BP	100 Gms 500 Gms 1 Kg 5 Kg	750 2200 4000 16000
31498 [67953-08-6]	🌡 **L-r-Glutamyl-p-Nitroanilide Hydrochloride extrapure, 99%**	500 Mg 1 Gms	2900 4900
92577 [111-30-8]	✔ **Glutardialdehyde 25% Aq. Solution** (Glutaric Dialdehyde Solution, Pentane-1,5-dial, Glutaraldehyde)	500 ml 2500 ml	800 3300
23475 [111-30-8]	✔ **Glutardialdehyde 50% Aq. Solution** (Glutaric Dialdehyde Solution, Pentane-1,5-dial)	500 ml 2500 ml	1200 4000
37394 [110-94-1]	✔ **Glutaric Acid pure, 99%**	100 Gms 500 Gms	1100 4700
87175 [108-55-4]	✔ **Glutaric Anhydride pure, 98%**	100 Gms 500 Gms	1600 6300
94980 [27025-41-8]	● **Glutathione Oxidized (GSSG) extrapure, 99%**	250 Mg 1 Gms 5 Gms 25 Gms 100 Gms	1500 3200 9000 30000 75000
41422 [27025-41-8]	● **Glutathione Oxidized (GSSG) ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	1 Gms 25 Gms 100 Gms	4000 36000 85000
22151 [27025-41-8]	● **Glutathione Oxidized (GSSG) for molecular biology, 99.5%**	250 Mg 1 Gms 5 Gms	2500 4500 18000
48938 [70-18-8]	● **Glutathione Reduced (GSH) extrapure, 98%**	5 Gms 25 Gms 100 Gms 500 Gms	900 3000 7000 24000
16111 [70-18-8]	● **Glutathione Reduced (GSH) for tissue culture, 99.5%**	5 Gms 25 Gms 100 Gms	1050 3300 7300
77285 [70-18-8]	● **Glutathione Reduced (GSH) ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP, BP, Ph.Eur	5 Gms 25 Gms 100 Gms	1300 3600 7600
11514 [70-18-8]	● **Glutathione Reduced (GSH) for molecular biology, 99.5%**	1 Gms 5 Gms 25 Gms 100 Gms	800 2000 5400 12000
14210 [453-17-8]	✔ **D-(+)-Glyceraldehyde extrapure, 85%**	1 Gms	15800
62423 [9030-66-4]	■ **Glycerol Kinase (GK) ex. Microorganism, 30U/mg powder**	1 K.Units	10000
36451 [106-91-2]	**Glycidyl Methacrylate pure, 97%**	100 ml 500 ml	880 2860
52160 [3033-77-0]	🌡 **Glycidyl Trimethyl Ammonium Chloride (GTMAC) pure, 90%** ((2,3-Epoxypropyl)trimethylammonium Chloride, EPTMAC)	100 Gms 500 Gms	880 3300
63956 [113728-13-5]	🌡 **Glycine-7-Amido-4-Methylcoumarin Hydrobromide Salt extrapure, 98%** (Glycine-AMC.HBr sal)	10 Mg 25 Mg	10500 20400
56204 [2462-31-9]	■ **Glycine Benzyl Ester Hydrochloride extrapure, 98%**	5 Gms 25 Gms	3200 12200
67536 [1738-76-7]	🌡 **Glycine Benzyl Ester p-Toluene Sulfonate extrapure, 98%**	25 Gms	6400
85055 [107-43-7]	🌡 **Betaine (Glycine Betaine) pure, 98%**	25 Gms 100 Gms 500 Gms	900 1900 5000
99993 [623-33-6]	🌡 **Glycine Ethyl Ester Hydrochloride extrapure, 99%**	100 Gms 500 Gms	600 2300
18308 [7758-33-0]	🌡 **Glycyl Histidine Glycine extrapure, 95%**	100 Mg	9600
47055 [72957-37-0]	🌡 **Glycine-L-Histidine-L-Lysine extrapure, 98%**	100 Mg	12900
98803 [31692-85-0]	△ **Glycofurol extrapure, 99%** (Tetraglycol, Tetrahydrofurfurylalcohol Polyethylene Glycolether)	100 ml 500 ml	4600 18000
74078 [9005-79-2]	🌡 **Glycogen ex. Oyster extrapure, 85%**	5 Gms 25 Gms	7000 22000
49740 [9005-79-2]	🌡 **Glycogen ex. Oyster for molecular biology, 85%**	1 Gms 5 Gms 25 Gms	2200 8000 26000
50006 [23147-58-2]	✔ **Glycolaldehyde Dimer extrapure,95%**	500 Mg	4200
50391 [64480-66-6]	🌡 **Glycoursodeoxycholic Acid (GUDCA) extrapure, 97%** (Glycylursodeoxycholic Acid, Ursodeoxycholylglycine)	1 Gms 5 Gms	5000 18000
97291 [13115-71-4]	🌡 **Glycyl-L-Glutamine extrapure, 99%**	1 Gms 5 Gms	2100 8300
11530 [556-50-3]	🌡 **Glycyl Glycine extrapure, 99%**	25 Gms 100 Gms 500 Gms	1700 5800 23100
72998 [102601-58-1]	🌡 **Z-Glycyl-Glycyl-L-Arginine 7-Amido-4-Methylcoumarin Hydrochloride Salt extrapure, 95%** (Z-Glycine-Glycine-Arginine-AMC.HCl Salt)	5 Mg 25 Mg	9900 40900

Code	Name	Packing	Price
15746 [65022-15-3]	■ **Z-Glycyl-LProline-4-Nitroanilide (Z-Gly-Pro-4-Nitroanilide) extrapure, 98%**	25 Mg 100 Mg	4200 8400
41675 [556-33-2]	**Glycyl-Glycyl-L-Glycine (Triglycine) extrapure, 99%**	1 Gms 5 Gms 25 Gms	1800 7700 28300
27057 [69242-40-6]	**Glycyl-Glycyl-L- Isoleucine extrapure, 99%**	1 Gms 5 Gms	4900 20200
70039 [14857-82-0]	**Glycyl-Glycyl-L-Leucine extrapure, 99%**	1 Gms 5 Gms	4700 19900
29764 [20274-89-9]	**Glycyl-Glycyl-L-Valine extrapure, 99%**	25 Mg 100 Mg	7500 22500
82725 [19461-38-2]	**Glycyl-L-Isoleucine extrapure, 99%**	1 Gms 5 Gms	2100 8000
96010 [869-19-2]	**Glycyl-L-Leucine extrapure, 99%**	1 Gms 5 Gms	1800 7900
27587 [103213-38-3]	**Glycyl Leucine Phenylalanine extrapure, 99%**	100 Mg	17200
25952 [115035-46-6]	**Glycyl-L-Proline-7-Amido-4-Methylc oumarin Hydrobromide Salt extrapure, 98%** (H-Gly-Pro-AMC-HBr, GP-AMC.HBr Salt)	25 Mg	35700
74329 [658-79-7]	**Glycyl-L-Tyrosine extrapure, 99%**	1 Gms 5 Gms	2100 8300
34601 [1963-21-9]	**Glycyl-L-Valine extrapure, 99%**	1 Gms 5 Gms	2600 10800
92245 [107-22-2]	**Glyoxal 40% Aq. Solution pure**	500 ml 2500 ml	580 2070
35657	**Griess Reagent for Nitrite extrapure**	100 ml	4260
53613 [90-05-1]	**Guaiacol extrapure, 99%**	250 Gms 1 Kg	1150 4140
44260 [73-40-5]	**Guanine extrapure, 99%**	25 Gms 100 Gms 1 Kg	800 2500 20000
37567 [635-39-2]	**Guanine Hydrochloride extrapure, 98%**	5 Gms 25 Gms	700 3000
48306 [10333-92-3]	**Guanine Sulphate extrapure, 98%**	1 Gms	1700
57892 [118-00-3]	**Guanosine extrapure, 99%**	25 Gms 100 Gms 500 Gms	600 1800 7000
95494 [5550-12-9]	■ **Guanosine-5-Monophosphate Disodium Salt (5-GMP-Na2) extrapure, 98%**	25 Gms 100 Gms	3700 8400
70564 [56001-37-7]	■ **Guanosine-5-Triphosphate Disodium Salt (5-GTP-Na2) extrapure, 90%**	100 Mg 500 Mg	3200 5300
11945 [65202-07-5]	**Guayule (DCPTA) for tissue culture, 98%** (2-(3,4-Dichlorophenoxy) -N,N-diethylethanamine)	25 Gms 100 Gms	2500 8300
11777 [6322-07-2]	■ **D-Gulono-1,4-Lactone extrapure, 98%**	1 Gms 5 Gms	2500 11800
97514 [1128-23-0]	■ **L-Gulono-1,4-Lactone extrapure, 98%**	10 Gms 25 Gms	5800 10400
37421	**Gunzburg Reagent**	100 ml	700
95068 [9008-02-0]	**Haemoglobin Powder ex. Bovine Erythrocytes**	100 Gms 500 Gms	2500 12000
52293 [304-21-2]	**Harmaline pure, 98%**	1 Gms	5200
99132 [148893-10-1]	**HATU extrapure, 98%**	5 Gms 25 Gms 100 Gms 500 Gms	1000 2800 9500 42000
44424 [94790-37-1]	**HBTU extrapure, 98%**	25 Gms 100 Gms 500 Gms	1200 4200 18000
73322 [330645-87-9]	**HCTU extrapure, 98%** (O-(6-Chlorobenzotriazol-1-yl) -N,N,N',N'-Tetramethyluronium Hexafluorophosphate)	5 Gms 25 Gms 100 Gms	900 2600 9000
23315 [150-39-0]	**HEDTA Buffer (N-(2-Hydroxyethyl) Ethylenediaminetriacetic Acid) extrapure, 99%**	25 Gms 100 Gms 500 Gms	1300 3200 14000
28394 [517-28-2]	**Hematoxylin Monohydrate** (Natural Black I)	5 Gms 25 Gms 100 Gms	2300 8500 30000
39573 [517-28-2]	**Hematoxylin Monohydrate PrimeDx**	250 Gms 1 Kg	80000 280000
78372 [16009-13-5]	**Hemin Chloride cryst. ex. Bovine extrapure, 98%** (Hematin Hydrochloride, Chlorohemin)	1 Gms 5 Gms	1600 5600
98457 [16009-13-5]	**Hemin Chloride cryst. ex. Porcine extrapure, 98%** (Hematin Hydrochloride, Chlorohemin)	1 Gms 5 Gms	1600 5600
58389 [9041-08-1]	**Heparin Sodium Salt ex. Bovine Intestinal Mucosa, 150 IU/mg**	20 K.Units 100 K.Units	1200 3600
50661 [9041-08-1]	**Heparin Sodium Salt ex. Bovine Intestinal Mucosa, 150 IU/mg**	100 Mg 1 Gms	1300 5300
29768 [161308-36-7]	**HEPBS Buffer extrapure, 99%**	10 Gms 25 Gms 100 Gms	7410 9870 35760
63732 [7365-45-9]	**HEPES Buffer extrapure, 99%** (N-[2-Hydroxyethyl] Piperazine-N[2-Ethane Sulphonic Acid])	25 Gms 100 Gms 500 Gms 1 Kg	530 1680 7350 13650
48489 [7365-45-9]	**HEPES Buffer for tissue culture, 99%** (N-[2-Hydroxyethyl] Piperazine-N[2-Ethane Sulphonic Acid])	25 Gms 100 Gms	790 2100

Code	Name	Packing	Price
16826 [7365-45-9]	**HEPES Buffer for molecular biology, 99.5%** (N-(2-Hydroxyethyl) Piperazine-N(2-Ethane Sulphonic Acid)	25 Gms 100 Gms 500 Gms	840 2310 8400
78805	**HEPES Buffer 1M Solution pH-7.3 for molecular biology**	100 ml	3360
68744 [75277-39-3]	**HEPES Sodium Salt extrapure, 99%** (N-2-Hydroxyethyl piperazine N-2-Ethane Sulphonic Acid Sodium Salt)	25 Gms 100 Gms 500 Gms	1000 2900 9800
32225 [75277-39-3]	**HEPES Sodium Salt for molecular biology, 99%** (N-2-Hydroxyethyl piperazine N-2-Ethane Sulphonic Acid Sodium Salt)	25 Gms 100 Gms 500 Gms	1700 4000 13300
39764 [75277-39-3]	**HEPES Sodium Salt for tissue culture, 99.5%** (N-2-Hydroxyethyl piperazine N-2-Ethane Sulphonic Acid Sodium Salt)	25 Gms 100 Gms	1600 3600
86017 [110-43-0]	**2-Heptanone (Methyl-n-Amyl Ketone, MAK) extrapure 99%** (Methyl Pentyl Ketone)	100 ml 500 ml	500 2000
18338 [2528-61-2]	**Heptanoyl Chloride pure, 98%** (Enanthic Chloride, Oenanthic Chloride)	25 Gms 100 Gms	650 2000
29681 [10534-89-1]	**Hexaammine Cobalt (III) Chloride extrapure, 99%** (Cobalt Hexammine Trichloride, Hexaamminecobalt Trichloride, Luteo-cobaltic Chloride)	5 Gms 25 Gms	3200 9500
76570 [16941-12-1]	**Hexachloroplatinic Acid extrapure, 37-40% Pt** (Hydrogen Hexachloroplatinate (IV))	250 Mg 1 Gms 5 Gms	5300 15000 68000
28373 [143-27-1]	**Hexadecylamine (Cetylamine) pure, 98%**	5 Gms 25 Gms 100 Gms	1040 2190 4140
21437 [920-66-1]	**1,1,1,3,3,3-Hexafluoro-2-Propanol (HFIP) extrapure, 99%** (Hexafluoroisopropanol)	25 ml 100 ml 500 ml	2990 8740 41060
98573 [920-66-1]	**1,1,1,3,3,3-Hexafluoro-2-Propanol (HFIP) for UV spectroscopy, 99%** (Hexafluoroisopropanol)	25 ml 100 ml	4830 12080
28437 [999-97-3]	**1,1,1,3,3,3-Hexamethyldisilazane (HMDS) pure, 98%**	100 ml 500 ml 2500 ml	330 1400 5500
83721 [124-09-4]	**Hexamethylenediamine (HMDA) pure, 98%** (1,6-Diaminohexane, 1,6-Hexanediamine)	100 Gms 500 Gms	1200 3000
79404 [34338-96-0]	**(2S,5S)-2,5-Hexanediol extrapure, 99%**	1 Gms 5 Gms	3000 9000
56464 [111-27-3]	**n-Hexanol pure, 98%**	500 ml 2500 ml	1100 4830
98669 [111-27-3]	**n-Hexanol Reagent Grade, 99.5%**	250 ml 500 ml 2500 ml	840 1380 5750
68497 [142-61-0]	**Hexanoyl Chloride pure, 98%**	100 ml 500 ml	1100 4200
69502 [9001-51-8]	**Hexokinase ex. Saccharomyces Sp. 150U/mg powder**	5000 Units	6900
46639 [136-77-6]	**4-Hexylresorcinol extrapure, 99%** (Antascarin, Ascarinol, 4-Hexyl-1,3-Dihydroxybenzene)	5 Gms 25 Gms 100 Gms	800 2000 6500
35696 [495-69-2]	**Hippuric Acid pure, 99%** (Benzoylglycine, Benzoylaminoacetic Acid)	100 Gms 500 Gms	500 1380
97843 [495-69-2]	**Hippuric Acid Reference Standard, 100%**	10 Gms	5180
98808 [744-59-2]	**Hippuryl-L-Phenylalanine extrapure, 99%**	1 Gms	7000
83597 [56-92-8]	**Histamine Dihydrochloride (HSM), 98%**	1 Gms 5 Gms 25 Gms	1070 2970 8000
81240 [39968-33-7]	**HOAT (1-Hydroxy-7-azobenzotriazole) pure, 98%**	1 Gms 5 Gms 25 Gms	780 2000 8500
38653 [2592-95-2]	**HOBT Anhydrous extrapure, 99%**	100 Gms 500 Gms 1 Kg	800 3000 5200
27645 [123333-53-9]	**HOBT Monohydrate pure, 99%**	100 Gms 500 Gms 1 Kg	700 2400 4100
93864 [121080-96-4]	**D-Homocitrulline extrapure CHR, 98%**	1 Gms 5 Gms	4000 12000
94194 [82795-51-5]	**D-Homophenylalanine extrapure, 97%**	500 Mg 5 Gms	2500 10000
78523 [943-73-7]	**L-Homophenylalanine extrapure, 97%**	500 Mg 1 Gms 5 Gms 25 Gms	1900 3200 7500 23700
58931 [89-51-0]	**Homophthalic Acid pure, 98%**	25 Gms 100 Gms	1600 5400
69080 [306-08-1]	**Homovanillic Acid extrapure, 99%**	250 Mg 1 Gms 10 Gms	1300 4100 23700
27698 [9067-32-7]	**Hyaluronic Acid Sodium Salt - Type 1 (600-1000 KD) extrapure, 90%** (HA, Hyaluronan, Hyaluronate, Hyasis 850P, Sodium hyaluronate, Poly(b-glucuronic acid-[1?3]-b-N-acet	100 Mg 1 Gms 5 Gms	4000 12000 35000
75810 [9067-32-7]	**Hyaluronic Acid Sodium Salt - Type 2 (1000-2000 kD) extrapure, 90%** (HA, Hyaluronan, Hyaluronate, Hyasis 850P, Sodium hyaluronate, Poly(b-glucuronic acid-[1?3]-b-N-acet	100 Mg 1 Gms 5 Gms	4000 12000 35000
39678 [9067-32-7]	**Hyaluronic Acid Sodium Salt - Type 3 (80-100 kD) extrapure, 90%** (HA, Hyaluronan, Hyaluronate, Sodium hyaluronate, Poly(b-glucuronic acid-[1->3]-b-N-acetylglucosamin	100 Mg 1 Gms 5 Gms	4000 12000 35000

Code	Name	Packing	Price
89002 [9067-32-7]	■ **Hyaluronic Acid Sodium Salt - Type 4 (40-50 kD) extrapure, 90%** (HA, Hyaluronan, Hyaluronate, Sodium hyaluronate, Poly(b-glucuronic acid-[1?3]-b-N-acetylglucosamine	100 Mg 1 Gms 5 Gms	4000 12000 35000
43384 [9067-32-7]	■ **Hyaluronic Acid Sodium Salt - Type 5 (8-15 kD) extrapure, 90%** (HA, Hyaluronan, Hyaluronate, Sodium hyaluronate, Poly(b-glucuronic acid-[1?3]-b-N-acetylglucosamine	100 Mg 1 Gms 5 Gms	4000 12000 35000
98512 [5341-61-7]	**Hydrazine Dihydrochloride pure, 98%**	100 Gms 500 Gms	340 1400
36681 [619-67-0]	**4-Hydrazinobenzoic Acid extrapure, 95%**	25 Gms 100 Gms	3300 9900
40935 [122-66-7]	**Hydrazobenzene pure, 97%** (1,2-Diphenylhydrazine)	25 Gms 100 Gms	600 2100
13542 [50-23-7]	**Hydrocortisone pure, 97%**	1 Gms	500
28449 [123-31-9]	**Hydroquinone pure, 98%** (1,4-Dihydroxybenzene, Quinol)	100 Gms 500 Gms	280 1200
25311 [123-31-9]	**Hydroquinone extrapure AR, 99%**	100 Gms 500 Gms	330 1500
17732 [123-31-9]	**Hydroquinone extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	400 1700
91893 [118-93-4]	**o-Hydroxyacetophenone pure, 95%**	100 ml 500 ml	1100 4400
58245 [121-71-1]	**m-Hydroxyacetophenone pure, 98%**	25 Gms 100 Gms	600 1500
54564 [99-93-4]	**p-Hydroxyacetophenone extrapure, 98%**	100 Gms 500 Gms	600 2300
66616 [100-83-4]	**m-Hydroxybenzaldehyde pure, 98%**	25 Gms 100 Gms 500 Gms	600 1300 5800
85692 [123-08-0]	**p-Hydroxybenzaldehyde pure, 99%**	100 Gms 500 Gms 5 Kg	600 2200 18900
78957 [99-96-7]	**p-Hydroxybenzoic Acid pure, 99%**	100 Gms 500 Gms	400 1000
11635 [114-63-6]	□ **p-Hydroxybenzoic Acid Sodium Salt (HBA Sodium) extrapure, 98%** (4-Hydroxybenzoic Acid Sodium Salt, Sodium p-Hydroxybenzoate)	25 Gms 100 Gms 500 Gms	1950 3800 18000
59645 [1137-42-4]	**4-Hydroxybenzophenone pure, 98%**	100 Gms 500 Gms	3900 13100
29893 [623-05-2]	**4-Hydroxybenzyl Alcohol (4-HBA) pure, 98%**	25 Gms 100 Gms 500 Gms	1400 4900 9500
12197 [9028-38-0]	■ **D-3-Hydroxybutyrate Dehydrogenase (3-HBDH) ex. Pseudomonas Sp., 100U/mg powder** ((b-Hydroxybutyrate Dehydrogenase, (R)-3-Hydroxybutanoate: NAD+ oxidoreductase)	1000 Units	28700
82739 [54970-72-8]	**2-Hydroxy-3,5,- Dichlorobenzene Sulphonic Acid Sodium Salt (DCHBS), 98%** (DHBS, 3,5-Dichloro-2-hydroxybenzenesulfonic acid, Sodium 3,5-dichloro-2-hydroxybenzenesulfonate, Tr	5 Gms 25 Gms 100 Gms	3410 8500 32000
81073 [3006-96-0]	● **HMBA Linker (4-(Hydroxymethyl)benzoic Acid) extrapure, 99%**	5 Gms 25 Gms	2200 7500
84818 [68858-21-9]	● **HMPA Linker (4-(Hydroxymethyl)phenoxyacetic Acid) extrapure, 98%**	1 Gms 5 Gms 25 Gms	2000 3500 14000
56043 [128446-35-5]	✔ **Hydroxypropyl-B-Cyclodextrin (HPBCD) ExiPlus, Multi-Compendial, 98%** meets compendial specs of USP	25 Gms 100 Gms	3200 8800
53133 [122-37-2]	**4-Hydroxydiphenylamine extrapure, 98%** (4-(Phenylamino)phenol)	25 Gms 100 Gms	4000 7000
58240 [13161-30-3]	**2-Hydroxypyridine-N-Oxide (HOPO) pure, 98%**	25 Gms 100 Gms	5200 19000
44923 [1306-06-5]	🌡 **Hydroxylapatite ~ 25% solids suspension (10g solids/40ml) in 0.001M**	40 ml 100 ml	2100 4000
62953 [67-47-0]	🌡 **5-Hydroxymethylfurfural (HMF) extrapure, 98%**	1 Gms 5 Gms 25 Gms 100 Gms	2070 6670 30250 70040
77491 [67-47-0]	🌡 **5-Hydroxymethylfurfural (HMF) ExiPlus, Multi-Compendial, 98%** meets compendial specs of BP, Ph.Eur	1 Gms 5 Gms 25 Gms	2650 9890 33810
79897 [87199-15-3]	**3-(Hydroxymethyl) Phenylboronic acid extrapure, 95%**	1 Gms 5 Gms 10 Gms	2700 12400 23100
93740 [63451-35-4]	**Hydroxynaphthol Blue**	5 Gms 25 Gms 100 Gms	350 1300 4700
69399 [524-38-9]	**N-Hydroxyphthalimide pure, 98%**	25 Gms 100 Gms 500 Gms	500 1600 5300
26921 [51-35-4]	**L-Hydroxyproline extrapure CHR, 99%**	25 Gms 100 Gms	700 2500
57803 [142-08-5]	**2-Hydroxypyridine pure, 98%**	100 Gms 500 Gms	3000 11100

Code	Name	Packing	Price
64675 [109-00-2]	**3-Hydroxypyridine pure, 98%**	25 Gms	800
		100 Gms	2600
		500 Gms	11600
87228 [626-64-2]	**4-Hydroxypyridine practical grade, 95%**	25 Gms	1600
		100 Gms	5400
47913 [6066-82-6]	**N-Hydroxysuccinimide pure, 97%** (NHS)	25 Gms	800
		100 Gms	2700
		500 Gms	12000
88997 ● [14348-40-4]	**3-Hydroxy-2,4,6-Tribromo-Benzoic Acid (TBHBA), 97%**	1 Gms	800
		5 Gms	3500
84688 ● [53279-72-4]	**3-Hydroxy-2,4,6-Triiodobenzoic Acid (HTBA), 97%**	1 Gms	3000
		5 Gms	9000
67941 [4350-09-8]	**5-Hydroxy-L-Tryptophan extrapure, 98%** (Serotonin precursor)	1 Gms	1900
		5 Gms	8400
11056 🌡 [31282-04-9]	**Hygromycin B (HGR) for tissue culture, 90%**	100 Mg	2600
		250 Mg	6300
67317 🌡 [31282-04-9]	**Hygromycin B (HGR), 90%**	100 Mg	1500
		250 Mg	3600
		1 Gms	9000
94663 🌡 [31282-04-9]	**Hygromycin B (HGR Solution), 50mg/ml in Aq. Solution**	20 ml	10300
		100 ml	40500
63573 ✔ [68-94-0]	**Hypoxanthine extrapure, 99%**	5 Gms	350
		25 Gms	1040
		100 Gms	3910
36693 [3034-50-2]	**4-Imidazolecarboxaldehyde pure, 98%** (4-Formylimidazole, 1H-Imidazole-4-Carbaldehyde, 5-Imidazolecarboxaldehyde)	1 Gms	2990
		5 Gms	6670
52106 [606-23-5]	**1,3-Indanedione pure, 96%** (1,3-Diketohydrindene, 1,3-Dioxoindan)	5 Gms	2990
		25 Gms	12080
32328 🌡 [615-13-4]	**2-Indanone pure, 98%**	10 Gms	2420
		25 Gms	4830
		100 Gms	16910
46486 [860-22-0]	**Indigo Carmine ACS**	25 Gms	410
		100 Gms	1210
		500 Gms	4950
31016 [22519-64-8]	**Indium (III) Chloride Tetrahydrate ultrapure, 99.99%** (Indium Chloride)	5 Gms	6200
		25 Gms	13000
52859 [20661-21-6]	**Indium (III) Hydroxide ultrapure, 99.99%**	5 Gms	4000
		25 Gms	10000
84697 [7440-74-6]	**Indium Metal Ingots, 99.9%**	10 Gms	3200
95867 [207398-97-8]	**Indium (III) Nitrate Hydrate extrapure, 99.9%**	5 Gms	4500
		10 Gms	8000
		50 Gms	24000
43095 [1312-43-2]	**Indium (III) Oxide extrapure, 99.9%** (Diindium Trioxide, Indium Sesquioxide)	5 Gms	3000
		25 Gms	6900
		100 Gms	20000
64331 [13464-82-9]	**Indium (III) Sulphate Anhydrous extrapure, 99%**	10 Gms	2400
		25 Gms	5300
75920 [13464-82-9]	**Indium (III) Sulphate Anhydrous ultrapure, 99.99%**	10 Gms	4800
		50 Gms	22100
96535 [7440-74-6]	**Indium Wire Ø 0.5 mm, 99.9%**	5 Gms	3000
		10 Gms	5200
19228 [120-72-9]	**Indole crystalline extrapure AR, 99%**	25 Gms	600
		100 Gms	1600
		500 Gms	7200
69522 🌡 [1912-33-0]	**Indole-3-Acetic Acid Methyl Ester extrapure, 98%** (Methyl-3-Indolylacetate)	1 Gms	1800
		10 Gms	9900
63661 🌡 [57105-50-7]	**Indole-3-Acetyl-L-Phenylalanine extrapure, 98%** (N-(3-Indolylacetyl)-L-Phenylalanine, IAA-L-Phe)	1 Gms	36300
		5 Gms	108700
13550 🌡 [57105-42-7]	**Indole-3-Acetyl-L-Valine extrapure, 99%**	25 Mg	2500
60704 ✔ [1204-06-4]	**Indole-3-Acrylic Acid extrapure, 99.5%** (3-Indoleacrylic Acid)	1 Gms	3100
		5 Gms	15200
97637 ✔ [830-96-6]	**Indole-3-Propionic Acid pure, 98%**	5 Gms	1200
		25 Gms	3500
69447 🌡 [487-89-8]	**Indole-3-Carboxyaldehyde extrapure, 99%** (3-Indolylaldehyde)	10 Gms	1400
		25 Gms	2000
		100 Gms	6300
96553 🌡 [126787-65-3]	**Indoxyl-ß-D-Galactopyranoside extrapure, 98%** (Indoxyl-Gal)	25 Mg	2500
39627 🌡 [3318-43-2]	**3-Indoxyl Phosphate Disodium Salt extrapure, 98%**	25 Mg	8000
		100 Mg	12000
		500 Mg	32000
25670 ■ [2642-37-7]	**3-Indoxyl Sulfate Potassium Salt extrapure, 98%**	50 Mg	4000
		250 Mg	9500
51113 🌡 [58-63-9]	**Inosine extrapure, 99%**	25 Gms	800
		100 Gms	2800
		500 Gms	12000
19594 🌡 [352195-40-5]	**Inosine-5-Monophosphate Disodium Salt Hydrate (5-IMP-Na2 Hydrate) extrapure, 99%** (Disodium Inosate)	5 Gms	2000
77985 🌡 [36051-67-9]	**Inosine-5-Triphosphate Disodium Salt (ITP-Na2) extrapure, 97%**	25 Mg	3200
		100 Mg	6300
98211 [87-89-8]	**Inositol pure, 99%** (Vitamin B8)	100 Gms	500
		500 Gms	1900
		1 Kg	3500
		5 Kg	15000
87114 [87-89-8]	**Inositol for tissue culture, 99%** (Vitamin B8)	100 Gms	550
		500 Gms	2100
		1 Kg	3800
		5 Kg	17000
53860 [9005-80-5]	**Inulin extrapure**	25 Gms	1200
		100 Gms	3700

Code	Name	Packing	Price
61568 [9001-57-4]	■ **Invertase (Saccharase) ex. Candida Sp., 100U/mg powder** (Invertase-ß-D-Fructofuranoside Fructohydrolase)	2 K.Units	12900
51611 [88-67-5]	🌡 **2-Iodobenzoic Acid pure, 98%**	25 Gms	1400
		100 Gms	4100
55129 [626-62-0]	✔ **Iodocyclohexane (Stabilized w/ Cu) extrapure, 98%** (Cyclohexyl Iodide)	10 Gms	3870
		25 Gms	8280
		100 Gms	12420
82104 [2050-77-3]	**1-Iododecane (Stabilized with Cu) extrapure, 98%** (Decyl Iodide, Decyliodide)	25 ml	2420
		100 ml	7710
83441 [54-42-2]	● **5-Iodo-2-Deoxyuridine extrapure, 98%**	500 Mg	1800
		5 Gms	11800
40686 [146-68-9]	● **2(4-Iodophenyl)-3 (4-Nitrophenyl) -5-Phenyl-2H-Tetrazolium Chloride (INT) extrapure AR, 98%**	1 Gms	1800
		5 Gms	8000
99918 [696-07-1]	● **5-Iodouracil extrapure, 98%**	5 Gms	2300
		10 Gms	3900
63088 [538-93-2]	**Isobutyl Benzene pure, 98%**	500 ml	690
		2500 ml	2880
34516 [84110-40-7]	**Isobutylboronic Acid extrapure, 98%** (2-Methylpropylboronic acid)	1 Gms	1700
		5 Gms	3700
		25 Gms	16400
78521 [78-77-3]	**Isobutyl Bromide (Stabilized w/ Ag) pure, 98%**	100 ml	1040
		500 ml	4830
30671 [513-38-2]	**Isobutyl Iodide (1-Iodo-2-Methylpropane) (Stabilized w/Ag) pure, 97%**	25 ml	5180
		100 ml	14260
		250 ml	32200
39634 [79-31-2]	**Isobutyric Acid extrapure, 99%**	500 ml	1270
		2500 ml	5750
35853 [79-30-1]	**Isobutyryl Chloride pure, 98.5%**	250 ml	1100
		1000 ml	3800
70170 [79-30-1]	**Isobutyryl Chloride specially purified, 98.5%**	250 ml	1500
		1000 ml	5500
42615 [1637-73-6]	**DL-Isocitrate Trisodium Salt, 94%**	1 Gms	2200
		5 Gms	7900
90189 [25339-17-7]	**Isodecanol pure, 98%**	500 ml	1200
60672 [97-54-1]	🌡 **Isoeugenol (mixture of cis and trans) pure, 98%**	100 ml	1300
		500 ml	5300
33401 [319-78-8]	**D-Isoleucine extrapure, 98%**	250 Mg	9000
		1 Gms	18000
		5 Gms	72000
26731 [73-32-5]	**L-Isoleucine extrapure CHR, 99%**	100 Gms	700
		500 Gms	3300
		1 Kg	6000
94082 [499-40-1]	✔ **Isomaltose extrapure, 98%**	25 Mg	4500
		100 Mg	8800
		1 Gms	31000
42021 [3371-50-4]	🌡 **Isomaltotriose extrapure, 97%** (6-O-[6-O-(a-D-Glucopyranosyl) -a-D-Glucopyranosyl]-D-Glucose)	25 Mg	5000
88109 [35997-20-7]	✔ **Isomaltotetraose extrapure, 95%**	5 Mg	3500
45433 [6082-32-2]	✔ **Isomaltopentaose extrapure, 95%**	1 Mg	51580
94885 [54-85-3]	✔ **Isoniazid extrapure, 99%** (Isonicotinic Acid Hydrazide, 4-Pyridinecarboxylic Acid Hydrazide)	5 Gms	300
		100 Gms	800
		500 Gms	3600
68866 [2365-40-4]	🌡 **N6-[2-Isopentenyl] Adenine extrapure, 98%**	250 Mg	2900
		500 Mg	4800
		1 Gms	8400
18018 [121-91-5]	**Isophthalic Acid pure, 99%**	500 Gms	460
67742 [99-63-8]	🌡 **Isophthaloyl Chloride pure, 98%**	100 Gms	1960
		500 Gms	6100
36321 [108-21-4]	**Isopropyl Acetate pure, 99%**	500 ml	580
		2500 ml	2070
61386 [7724-76-7]	🌡 **N6-(2-Isopentenyl) Adenosine pure, 98%** (6-(?,?-Dimethylallylamino)purine riboside, N6-(3,3-Dimethylallylamino)purine riboside, 6-(3,3-Dimet	50 Mg	2100
		200 Mg	6000
58571 [2210-25-5]	🌡 **N-Isopropylacrylamide (NIPAM) extrapure, 99%**	10 Gms	1600
		25 Gms	3200
		100 Gms	12000
45273 [80041-89-0]	**Isopropylboronic Acid extrapure, 97%**	5 Gms	3500
		25 Gms	10000
16390 [75-26-3]	**Isopropyl Bromide (Stab w/ Ag) pure, 99%**	100 ml	460
		500 ml	1730
18346 [105-48-6]	**Isopropyl Chloroacetate pure, 98%**	500 ml	1040
38525 [75-30-9]	**Isopropyl Iodide pure, 98%**	25 ml	2420
		100 ml	8050
		250 ml	19320
52469 [142-91-6]	**Isopropyl Palmitate pure, 90%**	500 ml	920
		2500 ml	3910
54110 [367-93-1]	🌡 **Isopropyl-B-D-Thiogalactopyranoside (IPTG) (dioxan free), 99%**	5 Gms	1500
		25 Gms	7000
		100 Gms	22000
		1 Kg	180000
67208 [367-93-1]	🌡 **Isopropyl-B-D-Thiogalactopyranoside (IPTG) (dioxan-free) for MB, 99%**	5 Gms	2000
		25 Gms	8000
		100 Gms	26000
		1 Kg	210000
82202	🌡 **IPTG Solution (20mg/ml)**	1 ml	1100
		5 x1 ml	2500
86577 [590-86-3]	🌡 **Isovaleraldehyde pure, 98%**	500 ml	750
		1000 ml	1300

Code	Name	Packing	Price
51656 [108-12-3]	**Isovaleryl Chloride pure, 98%** (Isopentanoyl Chloride, 3-Methylbutyryl Chloride)	100 Gms 500 Gms	1400 5300
35798 [97-65-4]	**Itaconic Acid pure, 99%**	250 Gms	320
18443 [70288-86-7]	**Ivermectin (IVM), 95-102%**	1 Gms 5 Gms 25 Gms	800 2000 6000
79238 [77026-92-7]	**Jasmonic Acid for molecular biology, 95%**	25 Mg	16120
98314 [61789-91-1]	**Jojoba oil extrapure**	250 ml 500 ml	2200 3500
70878 [64013-70-3]	**Kanamycin Acid Sulphate (KS) extrapure, 670U/mg** (Kanamycin Disulfate, Kanamycin A, Kanamycin Acid Sulfate)	1 Gms 5 Gms 25 Gms	500 2000 7500
89346 [64013-70-3]	**Kanamycin Acid Sulphate (KS) for tissue culture, 670U/mg** (Kanamycin Disulfate, Kanamycin A, Kanamycin Acid Sulfate)	1 Gms 5 Gms 25 Gms	700 3000 9500
99311 [25389-94-0]	**Kanamycin Monosulphate (KM) extrapure, 750mcg/mg**	1 Gms 5 Gms 25 Gms	800 2000 7000
76061 [25389-94-0]	**Kanamycin Monosulphate (KM) for tissue culture, 750mcg/mg**	1 Gms 5 Gms 25 Gms	950 2500 8500
14259 [600-18-0]	**2-Ketobutyric Acid pure, 97%** (a-Ketobutyric acid, Propionylformic Acid, 2-Oxobutyric Acid)	5 Gms 25 Gms	6500 30000
93929 [328-50-7]	**a-Ketoglutaric Acid extrapure, 99%**	25 Gms 100 Gms 500 Gms	700 2400 9900
28386 [328-50-7]	**a-Ketoglutaric Acid (High Purity) extrapure AR, 99.5%**	25 Gms 100 Gms 500 Gms	900 3300 15100
22826 [305-72-6]	**a-Ketoglutaric Acid Disodium Salt Anhydrous PrimeDx, 95%**	25 Gms 100 Gms 500 Gms	3400 11300 51600
18223 [305-72-6]	**a-Ketoglutaric Acid Disodium Dihydrate extrapure, 98%** (2-Oxoglutaric Acid Sodium Salt)	25 Gms 100 Gms 500 Gms	2900 9500 44800
90814 [22202-68-2]	**a-Ketoglutaric Acid Monosodium Salt extrapure, 98%** (2-Oxoglutaric acid monosodium salt, 2-Oxopentanedioic acid monosodium salt, Sodium 2-oxoglutarate m	25 Gms 100 Gms	5600 21300
61341 [525-79-1]	**Kinetin extrapure AR, 99%**	5 Gms 25 Gms 100 Gms	600 2000 7000
34154 [525-79-1]	**Kinetin for tissue culture, 99%**	1 Gms 5 Gms 25 Gms	400 1200 4000
83296 [4338-47-0]	**Kinetin Riboside extrapure, 98%**	100 Mg 1 Gms	6000 35000
19744 [501-30-4]	**Kojic acid extrapure, 99%** (5-Hydroxy-2-Hydroxymethyl-4-Pyranone)	5 Gms 25 Gms 100 Gms	210 550 2050
98158 [2140-29-6]	**Kojibiose extrapure, 99%**	10 Mg 50 Mg	3180 10780
81662	**Kovac`s Indole Reagent for microbiology**	100 ml	330
69068	**Kraut`s Reagent**	100 ml	680
53110 [9028-36-8]	**D-Lactate Dehydrogenase (D-LDH) ex. Microorganism, 400U/mg powder**	2500 Units 5000 Units	5000 9500
82334 [9001-60-9]	**r-D-Lactate Dehydrogenase (r-D-LDH) ex. E. Coli (Thermostable), 200U/mg powder**	5000 Units 10000 Units	7010 12620
15932 [8006-54-0]	**Lanolin anhydrous extrapure** (Wool Wax, Wool Grease)	500 Gms 5 Kg	1700 15000
56839 [79-63-0]	**Lanosterol pure, 50%**	1 Gms 10 Gms	2000 6000
20351 [10025-84-0]	**Lanthanum Chloride Heptahydrate extrapure AR, ACS, 99%**	100 Gms 500 Gms	1200 4000
73274 [7439-91-0]	**Lanthanum Metal Ultrapure, 200 mesh, 99.9%**	5 Gms 25 Gms	8500 14000
34032 [10277-43-7]	**Lanthanum Nitrate Hexahydrate Reagent Grade, 99%**	100 Gms 500 Gms	800 3000
91288 [10277-43-7]	**Lanthanum Nitrate 0.1M**	500 ml 1000 ml	1400 2500
36917 [312696-10-9]	**Lanthanum Oxalate Hydrate extrapure AR, 99.5%**	100 Gms 500 Gms	1600 7000
59415 [1312-81-8]	**Lanthanum Oxide extrapure, 99.9%**	100 Gms 500 Gms	800 3200
47835 [52093-26-2]	**Lanthanum (III) Trifluoromethanesulfonate extrapure, 98%** (Lanthanum Triflate, Trifluoromethanesulfonic Acid Lanthanum Salt)	1 Gms 5 Gms	1050 3250
92044 [111-82-0]	**Lauric Acid Methyl Ester Reference Standard, 99.5%(GC)** (Methyl Laurate)	5 Gms 25 Gms	1800 6500
87444 [142-18-7]	**L-Lauroyl-Rac-Glycerol extrapure (Imwitor 312), 99%** (1-Monododecanoyl-Rac-Glycerol, 1-Monolaurin, 1-Monolauroyl-Rac-Glycerol, Rrac-Glycerol 1-Laurate, D	1 Gms 5 Gms	3870 15570
74090 [137-16-6]	**N-Lauroylsarcosine Sodium Salt (Sarkosyl Sodium) pure, 98%**	25 Gms 100 Gms 500 Gms 2.5 Kg	1520 4000 9000 41000
15448 [137-16-6]	**N-Lauroylsarcosine Sodium Salt (Sarkosyl Sodium) for molecular biology, 98%** (Sodium Lauroyl Sarcosine)	25 Gms 100 Gms 250 Gms 1 Kg	2400 4800 7000 22000

Code	Name	Packing	Price
67983 [137-16-6]	**N-Lauroylsarcosine Sodium Salt (Sarkosyl Sodium Solution) 30% Aq. Solution**	500 ml 2500 ml	900 4000
63645 [8000-28-0]	**Lavender oil extrapure, 30-60% LA**	25 ml 100 ml 500 ml	900 2200 6000
71707 [7758-97-6]	**Lead (II) Chromate pure, 98%** (Plumbous Chromate)	500 Gms	1050
21854 [7439-92-1]	**Lead Metal Lumps, 99%**	500 Gms	680
62603 [7439-92-1]	**Lead Metal Powder extrapure AR, 99.5%**	500 Gms	1000
59557 [8002-43-5]	**Lecithin ex. Egg, 60%** (L-a-Phosphatidylcholine)	25 Gms 100 Gms	5470 19440
59561 [8002-43-5]	**Lecithin ex. Soya, 30%** (L-a-Phosphatidylcholine)	100 Gms 500 Gms	1770 3890
64613 [10466-61-2]	**L-Leucinamide Hydrochloride extrapure, 99%**	5 Gms	2000
61304 [328-38-1]	**D-Leucine extrapure, 99%**	5 Gms 25 Gms	1890 5000
99664 [33105-81-6]	**tert-DL-Leucine extrapure, 98%**	500 Mg	3000
47736 [62480-44-8]	**L-Leucine 7-Amido-4-Methylcoumarin Hydrochloride Salt extrapure, 98%** (L-Leucine-AMC.HCl Salt)	10 Mg 25 Mg	5500 11000
40528 [9082-71-7]	**Leucine Dehydrogenase (LeuDH) ex. Bacillus Sp., 20U/mg powder**	100 Units	14230
27029 [893-36-7]	**L-Leucine-ß-Naphthylamide Hydrochloride extrapure, 97%**	100 Mg	2630
11127 [7298-84-2]	**L-Leucyl-L-Alanine Hydrate extrapure, 98%**	1 Gms	4000
37520 [997-05-7]	**D-Leucyl-Glycine extrapure, 98%**	100 Mg 500 Mg	2520 11900
23467 [100986-85-4, 138199-71-0]	**Levofloxacin (LVX), 98%**	1 Gms 5 Gms 25 Gms	2800 5500 14000
60011 [5141-20-8]	**Light Green SF Yellowish, 65%** (C.I. No. 42095)	25 Gms 100 Gms 500 Gms	900 3250 15000
98099 [60-33-3]	**Linoleic Acid (Free Acid) extrapure, 98%**	5 ml 25 ml 100 ml	3000 5500 8000
13562 [60-33-3]	**Linoleic Acid (Free Acid) Reference Standard, 99%(GC)**	1 Gms 5 Gms	1640 6830
49462 [112-63-0]	**Linoleic Acid Methyl Ester Reference Standard, 99%(GC)** (Methyl Linoleate)	1 Gms 5 Gms 25 Gms	2630 4730 17850
60770 [9001-62-1]	**Lipase (Steapsin) ex. Aspergillus Sp., 40U/mg solids**	25 Gms 100 Gms 500 Gms	520 1500 6000
15471 [9004-02-8]	**Lipoprotein Lipase (LPL) ex. Pseudomonas Sp., 20U/mg powder**	10 Mg 50 Mg 250 Mg	2460 4000 12000
92661 [9029-60-1]	**Lipoxidase ex. Glycine Max (Soyabean), 50000U/mg powder**	100 Mg 250 Mg 1 Gms	6300 13130 39900
68457 [6108-17-4]	**Lithium Acetate Dihydrate ultrapure, 99%**	100 Gms 250 Gms 500 Gms	5580 13140 25740
86213 [6108-17-4]	**Lithium Acetate Dihydrate PrimeDx, 99%**	1 Kg	64800
86701 [16853-85-3]	**Lithium Aluminium Hydride (LAH) extrapure, 95%**	100 Gms	43200
24966 [6080-58-6]	**Lithium Citrate Tribasic Tetrahydrate (High Purity) extrapure AR, 99.5%**	100 Gms 250 Gms 1 Kg	1260 2520 9360
91740 [17023-24-4]	**Lithium Iodide Monohydrate pure, 98%**	25 Gms 100 Gms	5220 19800
43904 [2044-56-6]	**Lithium Lauryl Sulphate extrapure, 99%** (Lithium Dodecyl Sulphate)	5 Gms 25 Gms 100 Gms 500 Gms	3420 10080 29880 132480
89141 [867-55-0]	**Lithium Lactate (specially purified) extrapure, 99%** (Lactic Acid Lithium Salt)	100 Gms 500 Gms	5040 23220
69375 [7439-93-2]	**Lithium Metal, 99%**	100 Gms	38160
23410 [553-91-3]	**Lithium Oxalate extrapure, 99%**	25 Gms 100 Gms	3780 7020
16350 [14283-07-9]	**Lithium Tetrafluoroborate extrapure, 98%** (Lithium Borofluoride, Lithium Fluoroborate)	5 Gms 25 Gms	2880 8640
29288 [434-13-9]	**Lithocholic Acid (LCA) pure, 98%**	10 Gms 25 Gms	7500 18000
52737 [98079-52-8]	**Lomefloxacin Hydrochloride (LMF), 98.5%**	1 Gms 5 Gms	1260 3260
50320 [2591-17-5]	**D-Luciferin ex. Firefly, 99%** ((S)-4,5-Dihydro-2-(6-Hydroxy-2-Benzothiazolyl) -4-Thiazolecarboxylic Acid)	10 Mg 50 Mg 500 Mg	2500 5000 20000
58891 [115144-35-9]	**D-Luciferin Potassium Salt ex. Firefly, 99%**	100 Mg 500 Mg 1 Gms	4000 13000 20000
62851 [103404-75-7]	**D-Luciferin Sodium Salt ex. Firefly, 99%**	25 Mg 100 Mg 500 Mg	3200 5500 22000
49982 [521-31-3]	**Luminol extrapure, 98%**	5 Gms 25 Gms	1070 4160
60445	**Luminol - Peroxide Chemiluminisence Solution for ELISA**	100 ml	1900

Code	Name	Packing	Price
67273 [127-40-2]	■ **Lutein extrapure, 80% Lutein Ester** (Xanthophyll, a-Carotene-3,3`-diol)	25 Mg	9000
70795 [201853-23-8]	🌡 **L-Lysine 7-Amido-4-Methylcoumarin Acetate Salt extrapure, 98%** (L-Lysine-AMC.Acetate salt)	10 Mg 25 Mg	8870 18150
72890 [97281-40-8]	■ **L-a-Lysophosphatidylethanolamine (Egg Lyso PE) ex. Egg, 99%** (L-a-Lysocephalin, 3-sn-Lysophosphatidylethanolamine)	10 Mg 25 Mg	14500 24000
45822 [12650-88-3]	■ **Lysozyme (3x cryst) ex. Egg white (Muramidase) for molecular biology, 15000U/mg powder**	1 Gms 5 Gms 25 Gms	840 3090 14230
35355 [1114-34-7]	**D-Lyxose extrapure**	1 Gms 5 Gms	1890 5670
38894 [1949-78-6]	**L-Lyxose extrapure**	100 Mg 1 Gms	1470 8400
81889 [9032-75-1]	■ **Macerozyme "R-10" ex. Rhizopus Sp., 3000U/g powder** (Pectinase)	500 Mg 1 Gms 5 Gms	5000 9000 43000
66638 [9032-75-1]	■ **Macerozyme "R-10" ex. Trichoderma & Aspergillus Niger, 10000U/g powder** (Pectinase)	500 Mg 1 Gms	3250 5500
40169 [9001-64-3]	■ **Malate Dehydrogenase (MDH) ex. Microorganism, 40U/mg powder**	1 K.Units 5 K.Units	6500 28000
87448	■ **r-Malate Dehydrogenase (r-MDH) (Thermostable), 500U/mg**	5000 Units 10000 Units	7600 14020
84225 [108-31-6]	**Maleic Anhydride extrapure AR, 99%**	500 Gms 5 Kg	450 4200
18808 [13676-54-5]	✔ **Bis-Maleimide pure, 95%**	100 Gms	6850
55114 [80307-12-6]	● **4-Maleimidobutyric Acid N-Hydroxysuccinimide Ester (GMBS) extrapure, 99%**	100 Mg 1 Gms 5 Gms	3000 6000 13000
57096 [108-13-4]	**Malonamide extrapure, 98%**	25 Gms 100 Gms 500 Gms	1400 2100 8500
43541 [141-82-2]	**Malonic Acid pure, 99%**	250 Gms 500 Gms	550 950
26504 [141-82-2]	**Malonic Acid extrapure AR, 99%**	100 Gms 250 Gms 500 Gms	330 650 1200
32261 [109-77-3]	🌡 **Malononitrile extrapure, 98%** (Malonic Dinitrile, Dicyanomethane)	100 Gms 500 Gms	700 2900
43786 [1109-28-0]	**Maltotriose (DP3) extrapure, 90%** (a-1,4-Glucotriose)	250 Mg 1 Gms 5 Gms	3200 6900 33000
28083 [34612-38-9]	■ **Maltotetraose (DP4) extrapure, 97%** (a-D-Glucopyranosyl-1,4-O-a-D-Glucopyranosyl-1,4-O-a-D-Glucopyranosyl-1,4-D-Glucopyranose)	10 Mg 25 Mg	6000 10000
10450 [34620-76-3]	🌡 **Maltopentaose (DP5) extrapure, 95%**	10 Mg 25 Mg	12000 26900
27607 [34620-77-4]	🌡 **Maltohexaose (DP6) extrapure, 65%**	50 Mg	12000
12894 [34620-78-5]	🌡 **Maltoheptaose (DP7) extrapure, 80%** (Amyloheptaose)	25 Mg	7900
58199 [611-71-2]	✔ **R-(-)-Mandelic Acid extrapure, 99%** (D-(-)-Mandelic Acid)	25 Gms 100 Gms	920 2820
30727 [17199-29-0]	✔ **S-(+)-Mandelic Acid extrapure, 99%** (L-(+)-Mandelic Acid, S-Z-Hydroxyphenyl Acetic Acid)	25 Gms 100 Gms	920 2490
27072 [90-64-2]	**DL-Mandelic Acid pure, 99%**	100 Gms 500 Gms	1200 5500
94551 [5505-63-5]	🌡 **D-Mannosamine Hydrochloride extrapure, 98%**	100 Mg 500 Mg	2020 9610
13735 [3458-28-4]	**D-Mannose extrapure AR, 99%**	25 Gms 100 Gms 500 Gms 5 Kg	600 1100 4500 40000
58476 [3458-28-4]	**D-Mannose for tissue culture, 99%**	25 Gms 100 Gms 500 Gms	800 1500 5500
87875 [3458-28-4]	**D-Mannose extrapure AR, ExiPlus, Multi-Compendial, 99%** meets compendial specs of BP, Ph.Eur	25 Gms 100 Gms 1 Kg	700 1300 12000
72002 [104872-94-8]	🌡 **D-Mannose-6-Phosphate Barium Salt Hydrate extrapure, 98%**	25 Mg	25080
43804 [95713-52-3]	🌡 **Marfeys Reagent extrapure AR, 98%** (FDAA, FDNP-ALA-NH2, Na-(2,4-Dinitro-5-Fluorophenyl) -L-Alaninamide)	50 Mg 200 Mg 1 Gms	5000 9000 38000
70902 [73-31-4]	■ **Melatonin extrapure, 99%**	1 Gms 5 Gms 25 Gms	600 1400 5000
10030 [7057-57-0]	**Meldolas Blue extrapure** (8-Dimethylamino-2,3-Benzophenoxazine Hemi(zinc Chloride) Salt, Naphthol Blue R, Basic Blue 6)	5 Gms	2000
71405 [2033-24-1]	🌡 **Meldrum`s Acid extrapure, 98%** (2,2-Dimethyl-1,3-dioxane-4,6-dione cycl-Isopropylidene Malonate Malonic Acid Cyclic Isopropylidene	25 Gms 100 Gms 500 Gms	700 2000 7000
77393 [66009-10-7,585-99-9]	**D(+) Melibiose Monohydrate extrapure, 98%**	1 Gms 5 Gms 25 Gms	1600 3200 9000
83709 [10030-67-8]	**D-Melezitose Monohydrate extrapure, 99%**	25 Gms	15000
61495 [58-27-5]	**Menadione O/S (Vitamin K3) pure, 98%**	25 Gms 100 Gms	950 2100

Code	Name	Packing	Price
44111 [6147-37-1] ✔	**Menadione Sodium Bisulphite Trihydrate W/S (Vitamin K3 Sodium Bisulphite Trihydrate) pure, 95%** (2-Methyl-1,4-Napthoquinone Sodium Bisulfite)	25 Gms 100 Gms 500 Gms	750 2000 8200
26531 [96-27-5]	**1-Thioglycerol (3-Mercapto-1,2-Propanediol) pure, 90%** (Monothioglycerol)	25 ml 100 ml 500 ml	1300 3800 18300
54474 [32887-01-7] 🌡	**Mecillinam (MCL) extrapure, 95%**	25 Mg 100 Mg	3000 9000
56868 [107-96-0]	**3-Mercaptopropionic Acid pure, 98%**	100 ml 500 ml	1150 5000
18382 [6112-76-1] 🌡	**6-Mercaptopurine Monohydrate extrapure, 99%**	1 Gms 5 Gms	950 4600
42253 [2637-34-5] 🌡	**2-Mercaptopyridine extrapure, 98%** (2-Pyridinethiol, 2-Pyridyl Mercaptan)	25 Gms 100 Gms	2100 6500
81439 [4432-31-9]	**MES Free Acid Buffer extrapure, 99%** (4-Morpholineethanesulfonic Acid)	25 Gms 100 Gms 500 Gms	1050 3780 13650
31808 [145224-94-8]	**MES Monohydrate Buffer extrapure, 99%** (2-Morpholineethane Sulphonic Acid Monohydrate)	25 Gms 100 Gms 500 Gms 1 Kg	950 3000 9000 17000
86751 [145224-94-8]	**MES Monohydrate Buffer for tissue culture, 99%** (2-Morpholineethane Sulphonic Acid Monohydrate)	25 Gms 100 Gms	1050 3200
26854 [145224-94-8]	**MES Monohydrate Buffer for molecular biology, 99%** (2-Morpholineethane Sulphonic Acid Monohydrate)	25 Gms 100 Gms 500 Gms	1200 3500 11000
43396 [117961-21-4]	**MES Hemisodium Salt Buffer extrapure, 98%**	25 Gms 100 Gms	5430 18500
69408 [71119-23-8]	**MES Sodium Salt Buffer extrapure, 99%**	10 Gms 25 Gms 100 Gms	2290 3270 10120
72841 [104206-82-8]	**Mesotrione extrapure, 98%** (2-(4-Mesyl-2-Nitrobenzoyl) -1,3-Cyclohexanedione)(2-(4-Mesyl-2-nitrobenzoyl) -1,3-cyclohexanedione)	100 Mg 1 Gms	800 5000
47062 ✔	**MES1: Spectrometric Oil Standard D12 (12 elements in 10µg/g Aviation Reference Oil)** Ag, Al, Cr, Cu, Fe, Mg, Na, Ni, Pb, Si, Sn, Ti @ 10 µg/g Aviation Reference Oil	100 Gms	38000
19278 ✔	**MES2: Metals in Biodiesel Standard (5 elements in 10µg/g Biodiesel)** Ca, K, Mg, Na, P @ 10 µg/g	100 Gms	96000
55049 ✔	**MES3: ICP Alkali`s & Alkaline Earths (10 elements in 5% HNO3)** Ba, Be, Ca, Cs, K, Li, Mg, Na, Rb, Sr @ 100 µg/ml 5% HNO3	100 ml	33000
25562 ✔	**MES4: ICP Noble Metals (8 elements in 20% HCl)** Au, Ir, Os, Pd, Pt, Re, Rh, Ru @ 100 µg/ml 20% HCl	100 ml	42000
21551 ✔	**MES5: ICP Rare Earth & Geo Elements (21 elements in 5% HNO3)** Ba, Ce, Dy, Er, Eu, Gd, Ho, La, Lu, Nd, Pr, Rb, Sc, Sm, Sr, Tb, Th, Tm, U, Y, Yb @ 100 µg/ml 5% HNO3	100 ml	64000
82634 ✔	**MES6: ICP Non-Metals (7 elements in 5% HCl tr. HF)** (trace Hydrofluoric Acid (HF)) As, B, P, S, Se, Si, Te @ 100 µg/ml 5% HCl, tr. HF	100 ml	30000
18263 ✔	**MES7: ICP Common Elements Mix 1 (9 elements in 5% HNO3)** Cd, Co, Cr, Cu, Fe, Mn, Ni, V, Zn @ 100 µg/ml 5% HNO3	100 ml	35000
63924 ✔	**MES8: ICP Common Elements Mix 2 (20 elements in 5% HNO3, tr. HF)** Ag, Al, B, Ca, Co, Cr, Cu, Fe, K, Mg, Mn, Na, Ni, P, Pb, Si, Sn, Ti, V, Zn @ 100 µg/ml 5% HNO3, tr. HF	100 ml	54000
58498 ✔	**MES9: ICP Comprehensive Mix A (34 elements in 40% Aqua Regia)** Ag, Al, As, Ba, Bi, Ca, Cd, Ce, Dy, Er, Eu, Ga, Gd, Ho, La, Lu, Mg, Na, Nd, P, Pb, Pr, Rb, Sc, Se, Sm, Sr, Tb, Th, Tl, Tm, U, Y, Yb @ 10 µg/ml 40% Aqua Regia	100 ml	64000
36719 ✔	**MES10: ICP Comprehensive Mix B (32 elements in 40% Aqua Regia, tr. HF)** Au, B, Be, Co, Cr, Cu, Fe, Ge, Hf, Ir, K, Li, Mn, Mo, Nb, Ni, Os, Pd, Pt, Re, Rh, Ru, Sb, Si, Sn, Ta, Te, Ti, V, W, Zn, Zr @ 10 µg/ml 40% Aqua Regia, tr. HF	100 ml	85000
76134 ✔	**MES11: ICP 68 Element Multi Standard 1 (48 elements in 5% HNO3)** Al, As, B, Ba, Be, Bi, Ca, Cd, Ce, Co, Cr, Cs, Cu, Dy, Er, Eu, Fe, Ga, Gd, Ho, In, K, La, Li, Lu, Mg, Mn, Na, Nd, Ni, P, Pb, Pr, Rb, Re, Sc, Se, Sm, Sr, Tb, Th, Tl, Tm, U, V, Y, Yb, Zn @100 µg/ml (Standard 1 contains 48 elements) 5% HNO3	25 ml 100 ml	30000 78000
97799 ✔	**MES12: ICP 68 Element Multi Standard 2 (12 elements in 5% HNO3, tr. HF)** (trace Hydrofluoric Acid (HF)) Ag, Ge, Hf, Mo, Nb, Sb, Si, Sn, Ta, Ti, W, Zr @ 100µg/ml (Standard 2 contains 12 elements) 5% HNO3, tr. HF	100 ml	36000
45688 ✔	**MES13: ICP 68 Element Multi Standard 3 (8 elements in 10% HCl)** Au, Ir, Os, Pd, Pt, Rh, Ru, Te @ 100 µg/ml (Standard 3 contains 8 elements) 10% HCl	100 ml	45000

Code	Name	Packing	Price
53816	✔ **MES14: ICP Detector Calibration Multi- Element Mix (11 elements in 5% HNO3)** 7Li, Be, Mg, Sc, Co, Y, In, Ce, Tb, Tl, U (various conc.) 5% HNO3	100 ml	58000
91353	✔ **MES15: ICP Internal Standard Multi-Element Mix 1 (7 elements in 5% HNO3)** 6Li, Sc, Ga, Y, In, Tb, Bi @ 100 µg/ml 5% HNO3	100 ml	58000
51587	✔ **MES16: ICP Internal Standard Multi-Element Mix 2 (7 elements in 2% HNO3)** 6Li, Sc, Ga, Y, In, Tb, Bi (various conc.) 2% HNO3	100 ml	58000
67084	✔ **MES17: ICP Internal Standard Multi-Element Mix 3 (7 elements in 5% HNO3, tr. HF)** (trace Hydrofluoric Acid (HF)) 6Li, Sc, Ge, In, Tb, Lu, Bi @ 100 µg/ml 5% HNO3, tr. HF	100 ml	58000
31805	✔ **MES18: ICP Internal Standard Multi-Element Mix 4 (7 elements in 5% HNO3, tr. HF)** (trace Hydrofluoric Acid (HF)) 6Li, Sc, Ge, Te, In, Tb, Bi (various conc.) 5% HNO3, tr. HF	100 ml	58000
86382	✔ **MES19: ICP Internal Standard Multi-Element Mix 1 (25 elements in 75cSt Hydrocarbon Oil)**	100 Gms	50000
26027	✔ **MES20: ICP Calibration Standard for EPA 6020 (22 elements in 2% HNO3, tr. Tartaric Acid)** (trace Tartaric Acid)	100 ml	41000
83069	✔ **MES21: ICP 21 Calibration Standard (21 elements in 5% HNO3, tr. HF, tr. Tartaric Acid)** (trace Hydrofluoric Acid (HF), trace Tartaric Acid)	100 ml	41000
73584	✔ **MES22: ICP 19 Calibration Standard (18 elements in 5% HNO3)**	100 ml	41000
59878	✔ **MES23: ICP 29 Calibration Standard (29 elements in 5% HNO3)**	25 ml	43000
50138	✔ **MES24: ICP Environmental Calibration Standard A (24 elements in 5% HNO3, tr. HF, tr. Tartaric Acid)** (trace Hydrofluoric Acid (HF), trace Tartaric Acid)	100 ml	58000
73844	✔ **MES25: ICP Calibration Standard for EPA 6010 (18 elements in 5% HNO3, tr. Tartaric Acid)** (trace Tartaric Acid)	25 ml	43000
57183	✔ **MES26: ICP Water Pollution Standard 1 (15 elements in 5% HNO3)**	25 ml	35000
60648	✔ **MES27: ICP Common & Transition Elements Standard (1000µg/ml) (27 elements in 5% HNO3, 0.2% HF)** (0.2% Hydrofluoric acid (HF))	100 ml	58000

Code	Name	Packing	Price
76332 [121-47-1]	**Metanilic Acid pure, 98%**	250 Gms 1 Kg	200 720
20302 ✔ [7143-01-3]	**Methane Sulphonic Anhydride pure, 97%** (Mesic Anhydride)	25 Gms 100 Gms	1000 3000
54322 [348-67-4]	**D-Methionine extrapure, 99%**	5 Gms 25 Gms	1000 3000
65309 [59-51-8]	**DL-Methionine extrapure CHR, 99%**	100 Gms 500 Gms 1 Kg	460 1750 3200
98183 ✔ [15985-39-4]	**L-Methionine Sulfoximine (MSX) extrapure, 95%** (L-S-(3-Amino-3-Carboxypropyl) -S-Methylsulfoximine, L-Methionine Sulphoximine)	250 Mg 500 Mg 1 Gms 5 Gms	9400 16000 26000 105000
40043 🌡 [7314-32-1]	**L-Methionine Sulphone extrapure, 99%**	1 Gms	3430
19485 🌡 [3226-65-1]	**L-Methionine Sulphoxide extrapure, 99%**	1 Gms 5 Gms 25 Gms	2500 7500 28000
86596 🌡 [133073-73-1]	**Methotrexate Hydrate (MTR), 98%**	100 Mg 500 Mg 1 Gms	3000 9000 16000
91842 ✔ [100-06-1]	**p-Methoxyacetophenone pure, 99%** (4-Acetylanisole)	100 Gms 500 Gms	610 2000
81142 [2396-60-3]	**4-Methoxyazobenzene (4-MAB) extrapure AR, 98%** (p-Methyoxyazobenzene, 4-Phenylazoanisole)	1 Gms 5 Gms 25 Gms	2000 3500 8500
75797 □ [105-13-5]	**4-Methoxybenzyl Alcohol (p-Anisyl Alcohol) pure, 98%** (p-Methoxybenzyl Alcohol)	100 Gms 500 Gms	940 3200
13573 ● [93102-05-7]	**N-Methoxymethyl-N-Trimethylsilylm ethyl-Benzylamine extrapure, 98%**	5 Gms 25 Gms	11670 43730
51234 [5720-06-9]	**2-Methoxyphenylboronic Acid extrapure, 95%** (Boronic Acid, 2-Methoxyphenylboronic Acid)	1 Gms 5 Gms	1380 6170
40632 [5720-07-0]	**4-Methoxyphenylboronic Acid extrapure, 98%**	1 Gms 5 Gms 25 Gms	860 1670 8070
23438 [163105-89-3]	**2-Methoxypyridine-5-Boronic Acid extrapure, 98%** (6-Methoxy-3-Pyridinylboronic Acid)	1 Gms 5 Gms	1570 4800
24252 ■ [54344- 92-2]	**4-Methoxy-2,3,6-Trimethyl Benzaldehyde extrapure, 98%** (2,3,6-Trimethyl-p-Anisaldehyde)	1 Gms 5 Gms	3450 10350
88381 [105-45-3]	**Methyl Acetoacetate pure, 98%**	500 ml 2500 ml 25 Ltr	520 2300 20700
21392 [105-45-3]	**Methyl Acetoacetate Reagent Grade, 99%**	500 ml 2500 ml	1160 5150
64330 [122-00-9]	**p-Methylacetophenone pure, 98%**	100 Gms 250 Gms	400 800

Code	Name	Packing	Price
93149 [100-61-8]	**N-Methylaniline pure, 98%**	500 ml 2500 ml	750 3570
29879 ● [38894-11-0]	**3-Methyl-2-Benzothiazolinone Hydrazone HCl Monohydrate (MBTH) extrapure AR (High Purity), 99%**	1 Gms 5 Gms 25 Gms	1420 4580 20430
23343 [98139-72-1]	**3-Methylbutylboronic acid extrapure, 95%**	1 Gms 5 Gms	3610 17550
39876 ⓘ [96-34-4]	**Methyl Chloroacetate extrapure, 99%**	500 ml	680
64426 ■ [3100-04-7]	**1-Methylcyclopropene (1-MCP), 98%**	100 Mg 500 Mg	2800 8000
34428 [116-54-1]	**Methyl Dichloroacetate extrapure, 99%**	500 ml	1600
45880 [41798-81-6]	**3-Methyl-1,5-Di-p-Tolyl-1,4-Pentaza diene pure, 97%** (N,N-bis[(4-Methylphenyl)diazenyl]Methanamine)	1 Gms 5 Gms	4160 19570
31535 [1934-16-3]	**New Methylene Blue N extrapure, 70%** (Basic Blue 24, Methylene Blue N)	1 Gms 5 Gms 25 Gms	1400 3700 12000
28710 [6586-05-6]	**New Methylene Blue N Zinc Chloride Double Salt, 90%**	5 Gms 25 Gms	1540 6000
73603 [6586-05-6]	**New Methylene Blue N Zinc Chloride Double Salt for tissue culture, 90%**	5 Gms	3510
84753 [93-15-2]	**Methyl Eugenol extrapure, 99%**	100 ml 500 ml 2500 ml	1310 5580 25700
28409 [93-61-8]	**N-Methylformanilide extrapure, 99%**	100 ml 500 ml	1810 6000
77007 [107-31-3]	**Methyl Formate pure, 99%**	500 ml 2500 ml	700 3200
49281 [97-30-3]	**Methyl-a-D-Glucopyranoside extrapure, 99%** (a-Methyl-D-Glucoside)	100 Gms 500 Gms	1600 6500
37685 [603-76-9]	**1-Methylindole pure, 97%**	5 ml 25 ml 100 ml	1000 2500 8000
96595 [7114-03-6]	**Methyl Green Zinc Chloride Double Salt**	5 Gms 25 Gms	860 4000
74916 [7114-03-6]	**Methyl Green Zinc Chloride Double Salt ExiPlus, Multi-Compendial** meets compendial specs of BP, Ph.Eur	5 Gms 25 Gms	940 4140
85544 🌡 [134253-42-2]	**1-O-Methyl-beta-D-Glucuronic Acid Sodium Salt extrapure, 99%** (1-OMe-b-D-GlcA.Na)	50 Mg	3800
86720 🌡 [693-98-1]	**2-Methylimidazole pure, 99%**	100 Gms 500 Gms 1 Kg	420 1650 3020
15913 [93-16-3]	**Methylisoeugenol (mixture of cis and trans) (1,2-Dimethoxy-4-propenylbenzene) pure, 98%** (4-(1-Propenyl) -1,2-Dimethoxybenzene, Isoeugenyl Methyl Ether)	25 ml 100 ml 500 ml	660 1980 6600
73282 🌡 [2682-20-4]	**2-Methyl-4-Isothiazolin-3-One (MIT) pure, 95%** (2-Methyl-3(2H)-Isothiazolone, 2-Methyl-1,2-Thiazol-3-One, 2-Methylisothiazol-3(2H)-One)	5 Gms	7600
74384 [66-27-3]	**Methyl Methanesulphonate (MMS) extrapure, 98%**	25 ml 100 ml	2530 9350
42067 [66-27-3]	**Methyl Methanesulphonate (MMS) for tissue culture, 98%**	10 ml	1740
71512 🌡 [50-66-8]	**6-Methylmercaptopurine extrapure, 95%** (6-Methyl(thio)purine)	100 Mg 1 Gms	1800 5000
66960 [93-04-9]	**Methyl 2-Naphthyl Ether pure, 99%** (Yara Yara)	250 Gms 1 Kg	510 1400
67915 [109-01-3]	**1-Methylpiperazine extrapure, 99%** (N-Methylpiperazine)	100 ml 500 ml 2500 ml	340 1460 6840
40610 [16419-60-6]	**2-Methylphenylboronic Acid extrapure, 97%**	1 Gms 5 Gms 25 Gms	1770 2800 11630
20622 [17933-03-8]	**3-Methylphenylboronic Acid extrapure (m-Tolylboronic Acid) extrapure, 97%**	1 Gms 5 Gms 25 Gms	950 2670 12880
43322 [107-87-9]	**Methyl-n-Propyl Ketone pure, 99%**	500 ml 2500 ml	1240 5160
87964 [1340-02-9]	**Methyl Purple Indicator 0.1% Solution extrapure AR**	250 ml	910
33878 🌡 [608-25-3]	**2-Methylresorcinol pure, 98%** (2,6-Dihydroxytoluene)	25 Gms 100 Gms 250 Gms	1460 4350 10670
29783 [119-36-8]	**Methyl Salicylate pure, 99%**	500 ml 2500 ml	840 3700
41652 [80-48-8]	**Methyl p-Toluenesulfonate pure, 98%** (Methyl p-Tosylate, p-Toluenesulfonic Acid Methyl Ester)	100 ml 500 ml	1010 2940
44243 [1779-49-3]	**Methyltriphenylphosphonium Bromide pure, 98%**	100 Gms 500 Gms	1200 4710
18040 🌡 [90-33-5]	**b-Methylumbelliferone extrapure AR, 98%**	10 Gms 25 Gms 100 Gms	1240 2240 6720
13050 🌡 [37067-30-4]	**4-Methylumbelliferyl-N-Acetyl-ß-D-Glucosaminide extrapure, 99%**	25 Mg 50 Mg 100 Mg	6050 10080 19160
15208 🌡 [20671-66-3]	**4-Methylumbelliferyl Caprylate extrapure, 99%**	25 Mg 50 Mg 100 Mg	2240 3920 6720

Code	Name	Packing	Price
21896 [72626-61-0]	**4-Methylumbelliferyl-ß-D-Cellobiosi de extrapure, 98%**	25 Mg 50 Mg	4160 7270
34391 [6160-78-7]	**4-Methylumbelliferyl-ß-D-Galactopyranoside extrapure**	100 Mg 250 Mg	1200 1800
48838 [18997-57-4]	**4-Methylumbelliferyl-ß-D-Glucopyranoside extrapure, 98%** B-Glucosidase Substrate	500 Mg 1 Gms 5 Gms	5830 10080 44800
92009 [24589-78-4]	**N-Methyl-N-Trimethylsilyl Trifluoroacetamide (MSTFA) extrapure, 98%**	5 ml 25 ml 100 ml	3140 11200 33600
83294 [24589-78-4]	**MSTFA+TMCS (99:1)** (N-Methyl-N-(trimethylsilyl)trifluoroacetamide, Trimethylchlorosilane)	10 x1 ml	6720
86338 [66966-09-4]	**4-Methylumbelliferyl-a-L-Iduronide Free Acid extrapure, 97%**	2 Mg	53760
49959	**4-Methylumbelliferyl Myo-Inositol-1-Phosphate, N-Methyl-Morpholine Salt, Biosynth patent(WO99/48899)**	25 Mg	17010
21400 [15220-11-8]	**4-Methylumbelliferyl Sulfate Potassium Salt extrapure, 99%**	100 Mg 500 Mg	5180 14090
51504 [199329-67-4, 881005-91-0]	**4-Methylumbelliferyl-b-D-Glucuroni de Trihydrate (MUG) extrapure, 99%**	100 Mg 1 Gms 5 Gms	2120 3590 13440
77337 [1463-10-1]	**5-Methyluridine extrapure, 98%** (Ribothymidine)	5 Gms 25 Gms	2370 7460
25828 [22832-87-7]	**Miconazole Nitrate (MCN)extrapure, 98%**	1 Gms 5 Gms 25 Gms	970 2800 5600
67424 [72432-03-2]	**Miglitol extrapure, 98%** (N-(b-Hydroxyethyl) -1-Deoxynojirimycin, Glyset)	1 Gms 5 Gms	6160 13440
15761	**Molisch Reagent**	100 ml	460
64563 [500-44-7]	**Mimosine ex. Koa Hoale Seeds, 98%**	25 Mg 100 Mg	4560 13000
55690 [22373-78-0]	**Monensin Sodium Salt (MSN), 97%**	1 Gms 5 Gms	3330 8640
41781 [13939-06-5]	**Molybdenum Hexacarbonyl extrapure, 98%** (Hexacarbonylmolybdenum, Molybdenum carbonyl)	5 Gms 25 Gms	6720 26880
75556 [7439-98-7]	**Molybdenum Metal Sheet 99.9%**	1 piece	28790
68078	**10X MOPS Buffer for molecular biology**	100 ml 500 ml	680 2200
69824 [1132-61-2]	**MOPS Buffer extrapure, 99%** (3-Morpholinopropane Sulphonic Acid)	25 Gms 100 Gms 500 Gms 1 Kg	1000 2310 9450 17250
61310 [1132-61-2]	**MOPS Buffer PrimeDx, 99%** (3-Morpholinopropane Sulphonic Acid)	1 Kg 5 Kg	20000 95000
85123 [117961-20-3]	**MOPS Hemisodium Salt Buffer extrapure, 99%**	25 Gms 100 Gms	5190 16650
57851 [71119-22-7]	**MOPS Sodium Salt Buffer extrapure, 99.5%**	50 Gms 100 Gms 500 Gms	3210 6050 25000
66043 [1132-61-2]	**MOPS Buffer for molecular biology, 99%** (3-Morpholinopropane Sulphonic Acid)	25 Gms 100 Gms 500 Gms	1370 3360 12000
38602 [68399-77-9]	**MOPSO Buffer extrapure, 99%** (3-Morpholino-2-Hydroxy Propane Sulphonic Acid)	25 Gms 100 Gms 500 Gms	1200 3500 16000
82595 [79803-73-9]	**MOPSO Sodium Salt Buffer extrapure, 99%**	50 Gms 250 Gms 1 Kg	3000 10000 35000
24829 [480-16-0, 6472-38-4]	**Morin Hydrate (Aluminium Ionophore I, C.I. 75660) extrapure AR, 98%** (2',3,4',5,7-Pentahydroxyflavone)	1 Gms 5 Gms	1140 4200
77810 [64953-12-4]	**Moxalactam Sodium Salt (LTX) extrapure, 85%** (Latamoxef)	100 Mg 1 Gms	7500 24000
14189 [526-99-8]	**Mucic Acid extrapure, 99%**	25 Gms 100 Gms 500 Gms	690 2210 7880
93027 [84082-64-4]	**Mucin ex. Porcine Stomach, Type III, 90%**	10 Gms	6200
15358 [529-44-2]	**Myricetin cryst pure, 98%**	10 Mg 25 Mg 100 Mg 1 Gms	1800 3000 5000 8000
45945 [544-63-8]	**Myristic Acid pure, C14-95% (GC)**	500 Gms	590
36459 [544-63-8]	**Myristic Acid extrapure AR, C14-99% (GC)**	500 Gms 5 Kg	950 8500
86442 [544-63-8]	**Myristic Acid (High Purity) extrapure, C14-99.5% (GC)**	5 Gms 25 Gms	300 1070
70718 [124-10-7]	**Myristic Acid Methyl Ester Reference Standard, 99%(GC)**	5 Gms 25 Gms 100 Gms	1600 3200 9000
28136 [112-72-1]	**Myristyl Alcohol (1-Tetradecanol) pure, 98%**	100 ml 500 ml	210 840
20346 [389-08-2]	**Nalidixic Acid Free Acid (NA) extrapure, 99%** (8-Naphthyridine-3-Carboxylic Acid)	1 Gms 5 Gms	300 1320
96676 [120-18-3]	**Naphthalene-2-Sulphonic Acid Monohydrate extrapure, 98%**	25 Gms 100 Gms	1250 3980
55335 [123409-01-8]	**1,3,6-Naphthalenetrisulfonic Acid Trisodium Salt Hydrate (NTSA Sodium Hydrate) extrapure, 95%**	250 Gms 1 Kg	7000 16000
57955 [604-59-1]	**a-Naphthoflavone extrapure AR, 98%**	1 Gms 5 Gms 25 Gms	730 2650 12040

Code	Name	Packing	Price
93051 [1163-67-3]	**Naphthol-AS-Acetate extrapure, 99%**	100 Mg	1430
		500 Mg	5220
56061 [1237-75-8]	**Naphthol AS-BI extrapure, 98%**	5 Gms	720
70430 [37-87-6]	**Naphthol-AS-BI-b-D- Glucuronide extrapure, 99%**	25 Mg	4500
		100 Mg	10000
87408 [13989-98-5]	**Naphthol-AS-Phosphate extrapure, 95%**	250 Mg	11000
		500 Mg	18000
11401 [1919-91-1]	**Naphthol-AS-BI- Phosphate extrapure, 98%**	100 Mg	4160
		250 Mg	9490
87791 [1064-48-8]	**Naphthol Blue Black B**	25 Gms	530
		100 Gms	1700
		500 Gms	8000
45838 [521-24-4]	**1,2-Naphthoquinone-4- Sulphonic Acid Sodium Salt extrapure AR, 99%**	5 Gms	600
		25 Gms	2400
10150 [86-86-2]	**1-Naphthyl Acetamide (NAD) technical grade, 98%**	5 Gms	3490
		25 Gms	13700
42383 [207569-06-0, 2183-17-7]	**1-Naphthyl Phosphate Disodium Salt extrapure, 97%** (a-Naphthyl Phosphate Disodium Salt, Disodium 1-Naphthyl Phosphate)	250 Mg	3200
		1 Gms	6400
		5 Gms	24000
59523 [120-23-0]	**2-Naphthoxyacetic Acid (BNOA) extrapure, 98%**	25 Gms	570
		100 Gms	1950
43036 [830-81-9]	**a-Naphthyl Acetate extrapure AR, 99%**	5 Gms	1000
		10 Gms	1810
		25 Gms	3660
		100 Gms	11360
79388 [1523-11-1]	**b-Naphthyl Acetate extrapure AR, 99%**	5 Gms	1000
		10 Gms	1850
		25 Gms	3170
		100 Gms	10730
11939 [132-75-2]	**1-Naphthylacetonitrile pure, 97%**	25 Gms	4900
		100 Gms	18980
14722 [2243-83-6]	**2-Naphthoyl Chloride pure, 98%** (2-Naphthalenecarbonyl Chloride)	5 Gms	2000
		25 Gms	6000
		100 Gms	18000
61166 [1465-25-4]	**N-(1-Naphthyl)Ethylenediamine Dihydrochloride (NEDA) extrapure AR, ACS, 98%** (Bratton-Marshall Reagent, Marshall's Reagent)	5 Gms	600
		10 Gms	900
		25 Gms	2000
		100 Gms	7000
84942 [1465-25-4]	**N-(1-Naphthyl)Ethylenediamine Dihydrochloride (NEDA) extrapure AR,ACS,ExiPlus, Multi-Compendial, 98%** (Bratton-Marshall Reagent, Marshall's Reagent) meets compendial specs of BP	5 Gms	800
		10 Gms	1200
		25 Gms	2300
27734 [10236-47-2]	**Naringin pure, 97%**	5 Gms	900
		25 Gms	2000
		100 Gms	6000
95956 [67604-48-2]	**Naringenin pure, 95%**	1 Gms	1100
		5 Gms	2000
		25 Gms	5000
93007 [484-11-7, 654054-57-6]	**Neocuproine ACS, 99%** (2,9-Dimethyl-1,10-Phenanthroline)	1 Gms	1310
		5 Gms	5200
		25 Gms	20000
35449 [7296-20-0]	**Neocuproine Hydrochloride Monohydrate extrapure AR,99%** (2,9 Dimethyl-1-10-Phenanthroline Hydrochloride)	1 Gms	1640
		5 Gms	7600
23049 [10024-93-8]	**Neodymium (III) Chloride Anhydrous, 99.9%**	5 Gms	3600
		25 Gms	12000
78442 [13477-89-9]	**Neodymium (III) Chloride Hexahydrate extrapure, 99.9%**	25 Gms	4500
		100 Gms	10000
37499 [1313-97-9]	**Neodymium (III) Oxide extrapure AR, 99.9%**	10 Gms	1400
		25 Gms	3200
		100 Gms	12000
57369 [7440-00-8]	**Neodymium Metal Ultrapure, 200 mesh, 99.9%**	5 Gms	34000
		25 Gms	56000
86879 [16454-60-7]	**Neodymium (III) Nitrate Hexahydrate extrapure, 99.9%**	25 Gms	10600
		100 Gms	28800
99080 [34622-08-7]	**Neodymium (III) Trifluoromethanesulfonate extrapure, 98%** (Neodymium Triflate, Neodymium Perfluoromethanesulfonate)	1 Gms	6400
		5 Gms	16000
35282 [1405-10-3]	**Neomycin Sulphate (NMS), 600U/mg**	25 Gms	650
		100 Gms	1500
27309 [126-30-7]	**Neopentyl Glycol pure, 98%**	500 Gms	540
13607 [15501-33-4]	**Neopentyl Iodide extrapure, 98%** (2,2-Dimethyl-1-Iodopropane, 1-Iodo-2,2-Mimethylpropane)	5 ml	6180
		25 ml	19780
14678	**Nessler`s Reagent**	100 ml	600
60609 [56391-57-2]	**Netilmicin Sulphate (NET), 595ug/mg**	100 Mg	1800
49814 [3264-82-2]	**Nickel (II) Acetylacetonate Anhydrous extrapure, 95%**	25 Gms	4000
		100 Gms	10500
		500 Gms	38000
72860 [98-92-0]	**Nicotinamide extrapure, 99.5%** (Niacinamide, Pyridine-3-Carboxylic Acid Amide)	100 Gms	320
		500 Gms	1200
26759 [53-84-9]	**ß-Nicotinamide Adenine Dinucleotide (Oxidized) (ß-NAD, DPN) extrapure, 98%**	1 Gms	900
		5 Gms	3800
		25 Gms	17000
		100 Gms	60000
19628 [53-84-9]	**ß-Nicotinamide Adenine Dinucleotide (Oxidized) (ß-NAD, DPN) ExiPlus, Multi-Compendial, 98%** meets compendial specs of BP, Ph.Eur	1 Gms	970
		5 Gms	4100
		25 Gms	19000
		100 Gms	70000

Code	Name	Packing	Price
87253 [53-84-9] ■	**ß-Nicotinamide Adenine Dinucleotide (Oxidized) (ß-NAD, DPN) for molecular biology, 98%**	1 Gms 5 Gms	1410 5100
86561 [53-84-9] ■	**ß-Nicotinamide Adenine Dinucleotide (Oxidized) (ß-NAD, DPN) PrimeDx, 98%**	100 Gms	120000
77268 [606-68-8] ■	**ß-Nicotinamide Adenine Dinucleotide Disodium Salt (Reduced) (ß-NADH.Na2, DPNH.Na2) extrapure, 98%**	500 Mg 1 Gms 5 Gms 25 Gms	1250 2000 8000 27000
54941 [606-68-8] ■	**ß-Nicotinamide Adenine Dinucleotide Disodium Salt (Reduced) (ß-NADH.Na2, DPNH.Na2) for MB, 98%**	100 Mg 1 Gms 5 Gms	700 2500 9800
26511 [606-68-8] ■	**ß-Nicotinamide Adenine Dinucleotide Disodium Salt (Reduced) (ß-NADH.Na2, DPNH.Na2) PrimeDx, 98%**	100 Gms	220000
44566 [64417-72-7] ■	**ß-Nicotinamide Adenine Dinucleotide Lithium Salt (Oxidized) (ß-NAD.Li) extrapure, 95%**	500 Mg 1 Gms 5 Gms	2280 4320 11380
82819 [1184-16-3] ■	**ß-Nicotinamide Adenine Dinucleotide Phosphate Monosodium Salt (ß-NADP.Na, TPN.Na) extrapure, 98%**	100 Mg 1 Gms 5 Gms 25 Gms	820 2000 8000 31000
72970 [1184-16-3] ■	**ß-Nicotinamide Adenine Dinucleotide Phosphate Monosodium Salt (ß-NADP.Na) for MB, 98%**	100 Mg 1 Gms	1200 3200
57451 [24292-60-2] ■	**ß-Nicotinamide Adenine Dinucleotide Phosphate Disodium Salt (ß-NADP.Na2, TPN.Na2) extrapure, 95%**	100 Mg 500 Mg 1 Gms 5 Gms	850 3280 5610 18110
56640 [24292-60-2] ■	**ß-Nicotinamide Adenine Dinucleotide Phosphate Disodium Salt (ß-NADP.Na2, TPN.Na2) for MB, 95%**	100 Mg 1 Gms	1230 10280
99197 [2646-71-1] ■	**ß-Nicotinamide Adenine Dinucleotide Phosphate Tetrasodium Salt (Reduced) (ß-NADPH.Na4) extrapure, 98%** (TPNH2.Na4)	25 Mg 100 Mg 1 Gms	940 2690 22780
57647 [2646-71-1] ■	**ß-Nicotinamide Adenine Dinucleotide Phosphate Tetrasodium Salt (Reduced) (ß-NADPH.Na4) for MB, 98%** (TPNH2.Na4)	25 Mg 100 Mg	2580 5610
55615 [68141-45-7] ■	**ß-Nicotinamide Adenine Dinucleotide Phosphate Monopotassium Salt (Oxidized) (ß-NADP.K)extrapure, 95%**	500 Mg 1 Gms	3500 6000
88991 [20111-18-6] ■	**ß-Nicotinamide Adenine Dinucleotide Sodium Salt (ß-NAD.Na) extrapure, 95%**	1 Gms 5 Gms 25 Gms	2000 8500 35000
95130 [1094-61-7] ■	**ß-Nicotinamide Mononucleotide (ß-NMN) extrapure, 98%** (ß-Nicotinamide ribose monophosphate, NMN, Nicotinamide ribotide, Nicotinamide-1-ium-1-ß-D-ribofuran	1 Gms 5 Gms 25 Gms	1200 5000 20000
79601 [108347-85-9] 🌡	**ß-Nicotinamide Mononucleotide Disodium Salt (Reduced) (ß-NMNH) extrapure, 98%**	100 Mg 1 Gms 5 Gms	800 4000 12000
84999 [23111-00-4] 🌡	**Nicotinamide Riboside Chloride extrapure, 97%**	5 Gms 25 Gms	3500 8000
70035 [2415659-01-5]	**Nicotinamide Riboside Malate extrapure, 98%** (Nicotinamide Riboside (Malate))	5 Gms 25 Gms	3500 8000
26478 [10026-12-7]	**Niobium Chloride ultrapure, 99.9%** (Niobium Pentachloride)	5 Gms 25 Gms	2240 6720
73387 [1313-96-8]	**Niobium Pentoxide extrapure AR, 99.9%**	10 Gms 25 Gms 100 Gms 500 Gms	710 1440 4720 22640
31162 [150824-47-8] 🌡	**Nitenpyram extrapure, 98%**	1 Gms 5 Gms	2590 9410
39180 [121-89-1]	**m-Nitroacetophenone pure, 98%**	100 Gms 500 Gms	710 3310
21254 [552-89-6]	**2-Nitrobenzaldehyde extrapure AR, 99%** (o-Nitrobenzaldehyde)	25 Gms 100 Gms 500 Gms	650 2200 10000
17735 [555-16-8]	**p-Nitrobenzaldehyde extrapure, 99%**	25 Gms 100 Gms 500 Gms	600 1500 6500
17644 [99-61-6]	**m-Nitrobenzaldehyde pure, 98%**	100 Gms 500 Gms	550 2000
17003 [127-68-4]	**m-Nitrobenzene Sulphonic Acid Sodium Salt pure, 98%** (Sodium 3-Nitrobenzenesulfonate, Resist Salt)	1 Kg 5 Kg	680 2800
96144 [121-92-6]	**m-Nitrobenzoic Acid pure, 99%**	100 Gms 500 Gms	560 2500
81431 [121-92-6]	**m-Nitrobenzoic Acid extrapure AR, 99%**	25 Gms 100 Gms 250 Gms	490 1440 3140
79303 [62-23-7]	**p-Nitrobenzoic Acid pure, 99%**	100 Gms 500 Gms	440 2000
16148 [552-16-9] □	**o-Nitrobenzoic acid pure, 95%** (2-Nitrobenzoic acid, Ortho-Nitrobenzoic acid)	25 Gms 100 Gms 500 Gms	1100 2690 12320
27500 [122-04-3]	**4-Nitrobenzoyl Chloride (PNBC) pure, 98%** (p-Nitrobenzoic Acid Chloride)	25 Gms 100 Gms 500 Gms	830 1650 4810
39996 [619-73-8]	**p-Nitrobenzyl Alcohol pure, 97%**	25 Gms 100 Gms	1230 3840

Code	Name	Packing	Price
12790 [100-11-8]	**p-Nitrobenzyl Bromide (PNBBR) pure, 99%**	25 Gms	580
		100 Gms	1020
		500 Gms	4830
97695 [1083-48-3]	**4-(4-Nitrobenzyl)Pyridine extrapure AR, 99%**	5 Gms	6390
		25 Gms	24640
11207 [298-83-9]	■ **Nitro Blue Tetrazolium Chloride (Nitro BT) (NBT) extrapure AR, 99%**	100 Mg	610
		250 Mg	1100
		1 Gms	3300
		5 Gms	15400
48898 [298-83-9]	■ **Nitro Blue Tetrazolium Chloride (Nitro BT) (NBT) for molecular biology, 99%**	100 Mg	880
		250 Mg	1540
		1 Gms	4840
13001 [298-83-9]	■ **Nitro Blue Tetrazolium Chloride (Nitro BT) (NBT) for tissue culture, 99%**	100 Mg	1760
15527 [88-73-3]	✔ **o-Nitrochlorobenzene (ONCB) pure, 99%**	500 Gms	530
		5 Kg	4260
76773 [4883-67-4]	**2-Nitrocyclohexanone extrapure, 99%**	5 Gms	22960
84258 [67-20-9]	🌡 **Nitrofurantoin extrapure, 98-102%** (N-(5-Nitro-2-Furfurylidene) -1-Aminohydantoin, Furadoxyl)	5 Gms	1480
		25 Gms	3330
		100 Gms	6660
90728 [88-75-5]	✔ **o-Nitrophenol extrapure AR, 99%**	25 Gms	440
		100 Gms	1330
72892 [51-28-5]	**2,4-Dinitrophenol (2,4-DNP) moistened with water, 99%**	25 Gms	280
		100 Gms	790
21059 [577-71-9]	**3,4-Dinitrophenol pure, 97%**	1 Gms	8740
		5 Gms	26880
18432 [830-03-5]	**p-Nitrophenyl Acetate extrapure, 99%**	1 Gms	500
		5 Gms	1680
		25 Gms	4880
13463 [3740-52-1]	**2-Nitrophenylacetic Acid pure, 98%**	5 Gms	790
		25 Gms	2730
		100 Gms	9490
68146 [104-03-0]	**p-Nitrophenylacetic Acid pure, 98%**	25 Gms	1390
		100 Gms	4820
72565 [3459-18-5]	■ **4-Nitrophenyl-N-Acetyl-B-D-Glucosaminide extrapure, 99%**	100 Mg	3590
		1 Gms	26880
72367 [207591-86-4]	🌡 **4-Nitro-L-Phenylalanine Monohydrate extrapure, 99%**	5 Gms	4480
		25 Gms	8960
83602 [3482-57-3]	🌡 **p-Nitrophenyl ß-D-Cellobioside extrapure, 98%** (4-Nitrophenyl-beta-D-Cellobioside, PNP-beta-D-cellobioside)	100 Mg	13220
		500 Mg	34100
38325 [96597-16-9]	● **4-Nitrophenyl-4,6-Ethylidene-a-D-Maltoheptaoside (pNP-G7, EPS), EPS 90%** (Ethylidene-4-Nitrophenyl-a-D-Maltoheptaoside, EPS-PNP-G7)	25 Mg	6180
		100 Mg	13030
66398 [369-07-3]	● **o-Nitrophenyl-ß-D-Galacto Pyranoside (ONPG) extrapure, 99%**	1 Gms	1010
		5 Gms	4260
		25 Gms	11760
64685 [7493-95-0]	● **p-Nitrophenyl a-D-Galactopyranoside extrapure, 98%**	250 Mg	4180
		1 Gms	8740
		5 Gms	42000
96226 [3150-24-1]	● **p-Nitrophenyl-ß-D-Galactopyranoside (PNPG) extrapure, 98%**	250 Mg	1220
		500 Mg	2200
		1 Gms	3950
		5 Gms	7000
92734 [2492-87-7]	● **p-Nitrophenyl-ß-D-Glucopyranoside extrapure, 98%**	100 Mg	920
		1 Gms	4400
		5 Gms	14500
12735 [3767-28-0]	● **p-Nitrophenyl-a-D-Glucopyranoside extrapure, 98%**	100 Mg	1170
		1 Gms	3300
		5 Gms	10200
92822 [10344-94-2]	● **p-Nitrophenyl-ß-D-Glucuronide (PNPG (Gus)) extrapure, 99%**	100 Mg	2500
		500 Mg	7000
		1 Gms	12000
27166 [13331-27-6]	**3-Nitrophenylboronic Acid extrapure, 97%**	1 Gms	1130
		5 Gms	2240
		25 Gms	10250
82582 [100-16-3]	**p-Nitrophenylhydrazine extrapure AR, 99%**	10 Gms	930
		25 Gms	1560
88485 [333338-18-4/ 4264-83-9]	■ **p-Nitrophenylphosphate Disodium Salt Hexahydrate (pNPP) Reagent Grade, 99%**	5 Gms	800
		25 Gms	2800
		250 Gms	25000
90341 [333338-18-4/ 4264-83-9]	■ **p-Nitrophenylphosphate Disodium Salt Hexahydrate (pNPP) for molecular biology, 99%**	5 Gms	1000
		25 Gms	3200
84868 [333338-18-4/ 4264-83-9]	■ **p-Nitrophenylphosphate Disodium Salt Hexahydrate (pNPP) PrimeDx, 99%**	250 Gms	30000
		500 Gms	57000
93413 [68189-42-4]	■ **p-Nitrophenylphosphate Ditris Salt crystalline (pNPP-d) extrapure, 99%**	1 Gms	700
		5 Gms	2000
		25 Gms	7100
76230 [68189-42-4]	■ **p-Nitrophenylphosphate Ditris Salt crystalline (PNPP-d) PrimeDx, 99%**	50 Gms	22000
53456 [398128-60-4]	**4-Nitrophenyl trans-Ferulate extrapure, 98%** (3-Methoxy-4-Hydroxycinnamic Acid 4-Nitrophenyl Ester, PNP-trans-Ferulate)	25 Mg	3000
		100 Mg	6000
		1 Gms	24000
85173 [603-11-2]	**3-Nitrophthalic Acid pure, 98%**	100 Gms	450
		500 Gms	1500
46384 [610-27-5]	**4-Nitrophthalic Acid pure, 98%**	100 Gms	450
		500 Gms	1500
43713 [525-05-3]	**Nitroso R-Salt extrapure AR**	25 Gms	450
		100 Gms	1500
		500 Gms	6500
21872 [99-08-1]	**3-Nitrotoluene extrapure, 99%**	250 ml	340
19394 [121-14-2]	ⓘ **2,4-Dinitrotoluene pure, 98%** (1,3-Dinitro-4-Methylbenzene)	25 Gms	2250
		100 Gms	3500
		500 Gms	16500

Code	Name	Packing	Price
60944 [611-08-5]	5-Nitrouracil extrapure, 98%	10 Gms 25 Gms	2700 5510
81787 [6635-20-7]	5-Nitrovanillin pure, 98% (3-Methoxy-4-Hydroxy-5-Nitrobenzaldehyde)	5 Gms 25 Gms	2800 5600
64659 [821-55-6]	2-Nonanone pure, 98% (n-Heptyl Methyl Ketone)	100 ml 500 ml	1300 3200
56986 [51-98-9]	Norethisterone Acetate extrapure, 99% (Norethindrone Acetate)	100 Mg 1 Gms 5 Gms	1500 4200 8000
43000 [70458-96-7]	Norfloxacin (NFX), 98%	1 Gms 5 Gms	1700 4000
82075 [327-57-1]	L-Norleucine pure, 97% ((S)-(+)-2-Aminohexanoic acid, (S)-2-Aminocaproic Acid)	1 Gms 5 Gms	1800 7000
49021 [616-06-8]	DL-Norleucine extrapure, 99%	5 Gms 25 Gms	800 3000
27434 [1476-53-5]	Novobiocin Sodium Salt (NVB), 95%	100 Mg 1 Gms 5 Gms	1110 7000 30000
56221 [6600-40-4]	L-Norvaline extrapure, 99%	1 Gms 5 Gms	1780 3680
68547 [6409-77-4]	Nuclear Fast Red (C.I. No. 60760)	1 Gms 5 Gms	1500 7400
89249	Nylander`s Reagent	50 ml	430
56083	Obermayor`s Reagent	100 ml	240
51038 [124-30-1]	Octadecylamine pure, 98%	250 Gms 1 Kg 5 Kg	1100 4200 20000
27703 [111-64-8]	Octanoyl Chloride pure, 99%	100 ml 500 ml	650 2700
48804 [111-86-4]	N-Octylamine pure, 99% (1-Aminooctane, n-Octylamine, Caprylamine)	250 ml 1000 ml	1350 3750
60507 [149342-80-3]	Octyl a-D-Galactopyranoside extrapure, 98%	50 Mg	2800
63777 [40427-75-6]	Octyl ß-D-Galactopyranoside extrapure, 98%	50 Mg	2800
25253 [29781-80-4]	Octyl a-D-Glucopyranoside extrapure, 98%	100 Mg	14000
33134 [29836-26-8]	Octyl ß-D-Glucopyranoside extrapure, 98%	1 Gms 5 Gms	5500 25000
33910 [9002-93-1]	Octylphenyl Polyethylene Glycol (IGEPAL CA-630®) for molecular biology ((Octylphenoxy)polyethoxyethanol)	100 ml 500 ml	1040 4600
62799 [82419-36-1]	Ofloxacin (OFX), 98%	1 Gms 5 Gms	2200 5000
23576 [1320-06-5]	Oil Red O for molecular biology, 75% (Solvent Red 27, Sudan Red 5B) (C. I. No. 26125)	25 Gms 100 Gms	710 2220
23534 [7060-74-4]	Oleandomycin Phosphate extrapure, 100U/mg (Matromycin)	100 Mg	2000
70855 [8008-57-9]	Orange Oil extrapure	100 ml 500 ml	840 3890
30774 [1400-62-0]	Orcein	5 Gms 10 Gms 25 Gms	1370 2400 5000
41950 [8007-11-2]	Oregano Oil extrapure, 60-80% Carvacrol (Lippia Graveolens Oil, Origanum, Oil of Origanum)	100 ml 500 ml	2000 5000
86848 [3184-13-2]	L-Ornithine Monohydrochloride extrapure CHR, 99%	25 Gms 100 Gms 500 Gms	550 990 4400
98499 [7240-38-2]	Oxacillin Sodium Salt Monohydrate (OXC.Na.H2O), 815ug/mg	5 Gms	4330
62704 [328-42-7]	Oxalacetic Acid extrapure, 98% (Oxaloacetic Acid, Ketosuccinic Acid, Oxobutanedioic Acid, 2-Oxosuccinic Acid)	1 Gms 5 Gms 25 Gms	1700 8200 40000
76207	Oxidase Reagent (Gordon-Macleod Reagent)	100 ml	500
91088 [14698-29-4]	Oxolinic Acid (OXO), 99%	5 Gms	5000
29779 [59587-05-8]	Oxolinic Acid Sodium (OXO-Na), 98%	1 Gms 5 Gms	2690 9140
55540 [6153-64-6]	Oxytetracycline Dihydrate (OTC.2H2O) extrapure, 98%	1 Gms	410
12455 [2058-46-0]	Oxytetracycline Hydrochloride (OTC), 95%	1 Gms 5 Gms 25 Gms	600 2300 4000
61799 [9006-59-1]	Ovalbumin ex. Chicken Egg White for molecular biology, 90%	1 Gms 5 Gms 25 Gms	3220 14380 69000
57590 [76738-62-0]	Paclobutrazol pure, 95% (a-tert-Butyl-ß-(4-Chlorobenzyl)-1H-1,2,4-Triazole-1-Ethanol, PACZEE)	5 Gms 25 Gms	900 3500
42883 [33069-62-4]	Paclitaxel (Taxol, PCXL), 99%	10 Mg 50 Mg 250 Mg	1260 4000 7500
45967 [7440-05-3]	Palladium Sponge, 99.9%	250 Mg 1 Gms	9000 30000
42095 [42196-31-6]	Palladium(II) Trifluoroacetate extrapure, 98%	1 Gms 5 Gms	9500 41000
15028 [57-10-3]	Palmitic Acid pure, C16-98%	500 Gms 5 Kg	550 5000
34262 [57-10-3]	Palmitic Acid Reference Standard, C16-99% (GC)	5 Gms 25 Gms 100 Gms	580 2200 7600
77065 [112-39-0]	Palmitic Acid Methyl Ester Reference Standard, 99%(GC) (Methyl Palmitate)	5 Gms 25 Gms	2000 5880

Code	Name	Packing	Price
62199 [112-67-4]	✔ **Palmitoyl Chloride pure, 98%** (Hexadecanoyl Chloride, Palmitic Acid Chloride)	25 ml 100 ml 500 ml	1800 4000 8000
83162 [130-85-8]	**Pamoic Acid pure, 98%**	25 Gms 100 Gms	1200 4000
74635 [85-85-8]	**PAN Indicator (1(2-Pyridylazo)-2-Naphthol) ACS, 99%** (1(2-Pyridylazo)-2- Naphthol)	1 Gms 5 Gms 25 Gms	1000 4000 18000
37740 [8049-47-6]	■ **Pancreatin 3x ex. Porcine Pancreas, 75U/mg powder**	100 Gms 500 Gms	1010 4910
41421 [8049-47-6]	■ **Pancreatin 4NF ex. Porcine Pancreas, 100U/mg**	100 Gms 500 Gms	2900 11950
95907 [9001-73-4]	🌡 **Papain 1x USP ex. Papaya Latex, 6000 USP U/mg powder**	100 Gms 500 Gms 5 Kg	400 1730 14000
80400 [9001-73-4]	🌡 **Papain 5x USP ex. Papaya Latex, 30000 USP U/mg powder** (Papainase, Carica Papaya)	25 Gms 100 Gms 500 Gms	750 2200 8500
59196 [16593-81-0]	**PAR extrapure, 99%** (4-(2-Pyridylazo)-Resorcinol Monosodium Salt Hydrate)	1 Gms 5 Gms	860 3680
27436 [1263-89-4]	🌡 **Paromomycin Sulfate (PRM)**	1 Gms 5 Gms 25 Gms	3150 8900 34270
13684 [9046-40-6]	**Pectic Acid pure**	100 Gms	3700
90464 [9032-75-1]	■ **Pectinase ex. Aspergillus Niger, 3U/mg powder** (Macerozyme)	1 K.Units 5 K.Units 10 K.Units	1500 5000 8000
27295 [9032-75-1]	■ **Pectinase ex. Aspergillus Niger, 3000 IU/g powder** (Macerozyme)	100 Gms 500 Gms	2800 8000
38215 [1002-84-2]	**Pentadecanoic Acid (Pentadecylic Acid) extrapure, 98%** (Pentadecyclate, n-Pentadecylate)	5 Gms 25 Gms 100 Gms	2500 6500 18900
41991 [602-94-8]	**Pentafluorobenzoic Acid pure, 98%** (2,3,4,5,6-Pentafluorobenzoic Acid)	5 Gms 25 Gms 100 Gms	940 3280 7540
40309 [69-57-8]	✔ **Penicillin G Sodium Salt (PNL), 98%**	1 Gms 5 Gms	2500 6500
13585 [771-61-9]	✔ **Pentafluorophenol (PFP) extrapure, 99%** (2,3,4,5,6-Pentafluorophenol)	5 Gms 25 Gms	1420 5880
20895 [42075-32-1]	🌡 **(2R,4R)-2,4-Pentanediol extrapure, 99%**	1 Gms 2.5 Gms	33240 80550
10610 [72345-23-4]	🌡 **(2S,4S)-2,4-Pentanediol extrapure, 99%**	1 Gms 2.5 Gms	27700 67160
16165 [54-95-5]	🌡 **Pentylenetetrazole (Pentetrazole) extrapure, 99%** (1,5-Pentamethylenetetrazole) (Pentylenetetrazole,1,5-Pentamethylenetetrazole)	5 Gms 25 Gms	2000 6000
32212 [8006-90-4]	**Peppermint oil extrapure, 50%**	100 ml 500 ml	600 2200
12647 [9001-75-6]	🌡 **Pepsin 1:3000 ex. Porcine Stomach Mucosa, 0.8Anson U/mg**	25 Gms 100 Gms 500 Gms	1050 2550 11550
48062 [9001-75-6]	🌡 **Pepsin 1:10000 ex. Porcine Stomach Mucosa, 2.5 Anson U/mg powder**	25 Gms 100 Gms 500 Gms	2080 7400 33500
11787 [9001-75-6]	🌡 **Pepsin 1:10000 ex. Porcine Stomach Mucosa extrapure AR, 2.5 Anson U/mg powder**	25 Gms 100 Gms 500 Gms	2310 8090 37550
63364 [9001-75-6]	■ **Pepsin 2x cryst. ex. Porcine Stomach Mucosa, 2500U/mg powder**	1 Gms 5 Gms 25 Gms 100 Gms	3300 10450 38500 82500
28498 [9003-99-0]	■ **Peroxidase ex. Horseradish RZ >3.0, w/o Stab., Salt free (HRP Type 1A), 250U/mg Powder**	5 K.Units 25 K.Units 100 K.Units	4400 9900 30800
17589 [9003-99-0]	■ **Peroxidase ex. Horseradish RZ >3.0, w/o Stab., Salt free (HRP Type 1B), 250U/mg Powder**	100 Mg 500 Mg	9400 35000
73292 [9003-99-0]	■ **Peroxidase ex. Horseradish RZ >2.0, w/ 30% Stab., Salt free (HRP Type 2), 110U/mg powder**	5 K.Units 25 K.Units 50 K.Units	2200 6880 12380
87913 [9003-99-0]	■ **Peroxidase ex. Horseradish RZ >2.0, w/o Stab., Salt free (HRP Type 3), 180U/mg powder**	5 K.Units 25 K.Units 50 K.Units	2750 7700 13200
30709 [66-71-7]	**1,10-Phenanthroline Anhydrous extrapure, 99%, water 2%**	5 Gms 25 Gms 100 Gms	800 2800 8000
65898 [5144-89-8]	**1,10-Phenanthroline Monohydrate ACS, 99.5%**	5 Gms 25 Gms 100 Gms 500 Gms	560 2000 5800 26000
31888 [18851-33-7]	**1,10-Phenanthroline Hydrochloride Monohydrate extrapure AR, 99.5%**	5 Gms 25 Gms	900 4000
55782 [299-11-6]	● **Phenazonium Methosulphate (PMS) extrapure, 99%** (N-Methylphenazonium Methyl Sulphate, N-Methylphenazonium Methosulphate, Phenazine Methosulphate)	1 Gms 5 Gms 25 Gms	810 3280 15600
13920 [299-11-6]	● **Phenazonium Methosulphate (PMS) for tissue culture, 99%** (N-Methylphenazonium Methyl Sulphate, N-Methylphenazonium Methosulphate, Phenazine Methosulphate)	1 Gms 5 Gms 25 Gms	1060 4910 20450
48914 [68807-90-9]	**Phenolphthalein Diphosphate Tetrasodium Salt extrapure AR, 95%** (Phenolphthalein Bisphosphate Tetrasodium Salt)	1 Gms 5 Gms 25 Gms	1800 4340 16650

Code	Name	Packing	Price
95444 [122-59-8]	**Phenoxyacetic Acid extrapure, 99%** (Glycolic Phenyl Ether)	100 Gms	750
		500 Gms	2800
19630 [108321-84-2]	**L-Phenylalanine-7-Amido-4-Methylcoumarin Trifluoroacetate Salt extrapure, 98%** (L-Phenylalanine-AMC.TFA Salt)	25 Mg	11650
		50 Mg	16940
79455 [2462-32-0]	● **L-Phenylalanine Benzyl Ester Hydrochloride extrapure, 99%**	25 Gms	11710
34351 [3182-93-2]	● **L-Phenylalanine Ethyl Ester Hydrochloride extrapure, 99%**	1 Gms	1000
		5 Gms	4110
		25 Gms	12210
78586 [7524-50-7]	● **L-Phenylalanine Methyl Ester Hydrochloride extrapure, 99%**	1 Gms	880
		5 Gms	2500
		25 Gms	7700
89311 [65147-22-0]	● **Z-Phenylalanylarginine 7-Amido-4-Methylcoumarin Hydrochloride Salt extrapure, 99%** (Z-Phe-Arg-AMC.HCl Salt)	25 Mg	35050
14059 [91-40-7]	**N-Phenylanthranilic Acid (2-Anilinobenzoic Acid) extrapure AR, 99%**	25 Gms	1250
		100 Gms	4200
53993 [1689-82-3]	**4-Phenylazophenol (Solvent Yellow 7, 4-HAB) extrapure, 98%** (p-Hydroxyazobenzene)(C.I. No. 11800)	25 Gms	1200
		100 Gms	3400
92890 [98-80-6]	**Phenylboronic Acid extrapure, 95%**	25 Gms	1900
		100 Gms	7000
70291 [1821-12-1]	**4-Phenyl Butyric Acid pure, 99%**	25 Gms	2400
		100 Gms	6500
27644 [95-54-5]	🌡 **o-Phenylenediamine free base (OPD) extrapure AR, 99%**	5 Gms	2400
		25 Gms	8000
		100 Gms	18000
73229 [615-28-1]	🌡 **o-Phenylenediamine Dihydrochloride (OPD.2HCl) extrapure AR, 99%**	5 Gms	3800
		25 Gms	10790
59576 [64-04-0]	**2-Phenyl Ethylamine pure, 99%**	250 ml	800
		500 ml	1400
27861 [2871-15-0]	🌡 **Phenyl-a-D-Galactopyranoside extrapure, 98%**	250 Mg	10950
77334 [2818-58-8]	🌡 **Phenyl-ß-D-Galactopyranoside extrapure, 97%**	1 Gms	6830
45754 [122-60-1]	🌡 **Phenyl Glycidyl Ether (PGE) extrapure, 98%** (2,3-Epoxypropyl Phenyl Ether, 1,2-Epoxy-3-phenoxypropane)	100 ml	650
		500 ml	2400
23116 [2935-35-5]	🌡 **L-Phenyl Glycine extrapure, 99%**	25 Gms	720
		100 Gms	2800
		500 Gms	12000
17083 [103-72-0]	🌡 **Phenyl Isothiocyanate (PITC) pure, 98%**	25 ml	1500
		100 ml	4000
		500 ml	17500
18962 [103-72-0]	🌡 **Phenyl Isothiocyanate (PITC) Reagent Grade, 99%**	25 ml	2800
		100 ml	6020
15127 [103-72-0]	🌡 **Phenyl Isothiocyanate (PITC) for HPLC & Sequencing, 99%**	5 x1 ml	2000
		5 ml	1800
		5 x5 ml	8200
98152 [941-69-5]	🌡 **N-Phenylmaleimide extrapure, 98%**	5 Gms	800
		25 Gms	2200
		100 Gms	3500
89793 [62-38-4]	**Phenylmercuric Acetate (PMA) pure, 98%** (Mercury Phenyl Acetate)	5 Gms	2000
		25 Gms	8500
		100 Gms	32000
		250 Gms	75000
23358 [8003-05-2]	**Phenylmercuric Nitrate (Basic) (PMN) pure, 98%** (Mercury Phenyl Nitrate)	5 Gms	3200
		25 Gms	14500
		100 Gms	52000
		250 Gms	120000
87606 [329-98-6]	● **Phenylmethane Sulphonyl Fluoride (PMSF) extrapure, 99%**	1 Gms	750
		5 Gms	2400
		25 Gms	8500
		100 Gms	28000
84375 [329-98-6]	● **Phenylmethane Sulphonyl Fluoride (PMSF) for molecular biology, 99%**	1 Gms	950
		5 Gms	3000
		25 Gms	12000
90365 [66778-08-3]	🌡 **Phenylphosphate Disodium Salt Dihydrate (High Purity) extrapure AR, 98%**	25 Gms	2140
		100 Gms	6410
		250 Gms	16010
50478 [118-55-8]	**Phenyl Salicylate (Salol) extrapure, 99%**	100 Gms	2000
		500 Gms	8200
16795 [103-85-5]	**N-Phenylthiourea extrapure, 97%**	10 Gms	2140
		25 Gms	3800
		100 Gms	14000
30746	■ **Phleomycin (PLM) Solution (20mg/mL in 20mM HEPES)**	1 ml	7500
69353 [108-73-6]	**Phloroglucinol Anhydrous extrapure AR, 99%**	25 Gms	1300
		100 Gms	4500
		500 Gms	19000
99529 [18472-87-2]	🌡 **Phloxin B**	25 Gms	400
		100 Gms	1400
		500 Gms	6000
95131	**1X Phosphate Buffered Saline (PBS) for molecular biology**	100 ml	530
		500 ml	1290
78529	**10X Phosphate Buffered Saline (PBS) for molecular biology**	100 ml	660
		500 ml	1520
13435	🌡 **10X Phosphate Buffered Saline Tween-20 (PBST) for molecular biology**	200 ml	600
		1000 ml	2500
40083 [10526-80-4]	● **Phosphoenolpyruvate Monocyclohexylammonium Salt (PEP MCHA Salt) extrapure, 98%**	1 Gms	3000
		5 Gms	14000
81145 [10526-80-4]	● **Phosphoenolpyruvate Monocyclohexylammonium Salt (PEP MCHA Salt) PrimeDx, 98%**	50 Gms	150000

Code	Name	Packing	Price
31569 [35556-70-8]	● **Phosphoenolpyruvate Tricyclohexylammonium Salt (PEP TCHA Salt) extrapure, 98%**	250 Mg 1 Gms 5 Gms	3000 7000 29000
36257 [35556-70-8]	● **Phosphoenolpyruvate Tricyclohexylammonium Salt (PEP TCHA Salt) PrimeDx, 98%**	25 Gms	150000
84426 [9067-77-0]	■ **Phosphoenolpyruvate Carboxylase (PPC) ex. Microorganism, 5U/mg solid**	1 K.Units	20100
69400 [9001-87-0]	■ **Phospholipase D (Lecithinase D) ex. Streptomyces chromofuscus, 40U/mg powder**	200 Units 400 Units	17960 28730
84122 [26016-99-9]	🌡 **Phosphomycin Disodium Salt (Fosfomycin Disodium Salt) extrapure, 99%** ((-)-(1R,2S)-(1,2-Epoxypropyl)Phosphonic Acid, Phosphonomycin)	1 Gms 5 Gms	3570 7140
65312 [9012-69-5]	■ **Phosphorylase B ex. Rabbit Muscle, 50U/mg powder**	10 Mg	6000
27329 [643-79-8]	🌡 **o-Phthalaldehyde (OPA) extrapure AR, 99%** (Phthaldialdehyde)	5 Gms 25 Gms 100 Gms	1000 4000 9000
89886 [88-96-0]	**Phthalamide (Phthalic Acid Diamide) pure, 98%**	25 Gms 100 Gms	1260 2310
62747 [85-41-6]	**Phthalimide extrapure, 99%** (Phthalic Dicarboximide, 1,3-Dihydro-1,3-Dioxoisoindole, 1,3-Dioxoisoindoline)	500 Gms 5 Kg	580 5570
45501 [87-41-2]	**Phthalide pure, 98%** (1-Iso-Benzofuranone)	100 Gms 500 Gms	400 1800
26603 [4192-28-3]	✔ **N-Phthaloyl-L-Alanine extrapure, 98%**	1 Gms 5 Gms	1080 4890
44729 [340-90-9]	🌡 **N-Phthaloyl-L-Glutamic Acid extrapure, 98%**	5 Gms	4450
19971	🌡 **N-Phthaloyl-L-Methionine extrapure, 98%**	1 Gms 5 Gms	1830 7690
15187 [5123-55-7]	🌡 **N-Phthaloyl-L- Phenylalanine extrapure, 99%**	1 Gms 5 Gms	910 3170
64073	🌡 **N-Phthaloyl-L-Valine extrapure, 98%**	1 Gms 5 Gms	1670 7690
92486 [1445-69-8]	**Phthalyl Hydrazide pure, 98%** (Phthalic Hydrazide)	25 Gms 100 Gms 500 Gms	1200 3000 13500
58924 [554-62-1]	■ **Phytosphingosine Hydrochloride extrapure, 98%**	1 Gms 5 Gms	6800 18000
67432 [1918-02-1]	**Picloram extrapure, 95%** (4-Amino-3,5,6-Trichloropicolinic Acid, 4-Amino-3,5,6-Trichloropyridine-2-Carb oxylic Acid)	1 Gms 5 Gms 25 Gms	1600 4000 12000
41238 [98-98-6]	**a-Picolinic Acid pure, 98%**	25 Gms 100 Gms 500 Gms	650 2000 8500
93418 [124-87-8]	**Picrotoxin pure, 98%**	1 Gms 5 Gms	2200 6000
90670 [111-16-0]	**Pimelic Acid (Heptanedioic acid) extrapure, 98%**	25 Gms 100 Gms	3290 11650
99487 [76-09-5]	✔ **Pinacol pure, 98%** (2,3-Dimethyl-2,3-Butanediol, Tetramethylethylene Glycol)	25 Gms 100 Gms 500 Gms	800 2500 11000
31826 [2620-50-0]	**Piperonylamine pure, 97%** (1,3-Benzodioxole-5-Methylamine, 3,4-(Methylenedioxy)Benzylamine)	5 ml	7000
70375 [5625-37-6]	**PIPES Buffer extrapure, 99%** (Piperazine-1,4-bis(2-Ethanesulphonic, Acid))	25 Gms 100 Gms 500 Gms 1 Kg	850 2500 10250 18000
49159 [5625-37-6]	**PIPES Buffer for molecular biology, 99%** (Piperazine-1,4-bis(2-Ethanesulphonic, Acid))	100 Gms 500 Gms	3000 11000
48112 [5625-37-6]	**PIPES Buffer for tissue culture, 99%** (Piperazine-1,4-bis(2-Ethanesulphonic, Acid))	25 Gms 100 Gms	1200 3500
86180 [5625-37-6]	**PIPES Buffer PrimeDx, 99%**	1 Kg 5 Kg	24000 110000
58842 [108321-27-3]	**PIPES Dipotassium Salt Buffer extrapure, 99%**	25 Gms 100 Gms	4000 9000
55025 [76836-02-7]	**PIPES Disodium Salt Buffer extrapure, 99%**	25 Gms 100 Gms 500 Gms	3000 7000 25000
44033 [100037-69-2]	**PIPES Sesquisodium Salt Buffer extrapure, 99%**	25 Gms 100 Gms 250 Gms	3200 11000 24000
94348 [10010-67-0]	**PIPES Sodium Salt Buffer extrapure, 99%**	25 Gms 100 Gms 250 Gms	2000 6500 14000
72328 [110-85-0]	**Piperazine Anhydrous pure, 98%**	250 Gms 500 Gms 5 Kg	900 1600 13500
15861 [3282-30-2]	**Pivaloyl Chloride extrapure, 99%**	100 ml 250 ml 500 ml	480 850 1650
11599 [7440-06-4]	△ **Platinum Sponge, 99.9%**	250 Mg 1 Gms	12000 35000
20924 [26763-19-9]	■ **Polyadenylic Acid Potassium Salt (POLY A Potassium) extrapure**	100 Mg 1 Gms	3500 25000
44196 [26936-40-3]	■ **Polycytidylic Acid Potassium Salt (POLY C Potassium) extrapure**	100 Mg	4000
62995 [26936-41-4]	■ **Polyinosinic Acid Potassium Salt (POLY I Potassium) extrapure**	100 Mg	4000
54422 [25988-63-0]	■ **Poly-L-Lysine Hydrobromide (PLL HBr), 70-150kDa** (L-Lysine Homopolymer Hydrobromide)	25 Mg 100 Mg 250 Mg	9500 33000 75000

Code	Name	Packing	Price
27007 [27964-99-4]	■ **Poly-D-Lysine Hydrobromide (PDL HBr), 70-150KDa** (D-Lysine Homopolymer Hydrobromide)	25 Mg 100 Mg	11000 36000
22107 [28211-04-3]	■ **Poly-e-L-Lysine Hydrochloride (PLL HCl), 3500-4500Da** (e-Poly-L-Lysine HCl, N6-[(Phenylmethoxy)carbonyl]-L-Lysine Homopolymer)	25 Gms 100 Gms	8000 25000
90347 [25513-46-6]	■ **r-Polyglutamic Acid technical grade, 20%**	1 Gms 5 Gms	800 3000
67258 [36562-70-6]	🌡 **Polyguluronic Acid extrapure, 80%**	1 Gms 5 Gms	14500 62000
83011 [1405-20-5]	■ **Polymyxin B Sulphate (PB, Aerosporin), 6000U/mg**	250 Mg 1 Gms	1400 4000
22810 [68412-54-4]	✔ **Polyoxyethylene (5) Nonylphenylether Branched (IGEPAL CO-520) for molecular biology**	100 Gms 500 Gms	1100 3700
20659 [68412-54-4]	✔ **Polyoxyethylene (9) Nonylphenylether Branched (IGEPAL CO-630) for molecular biology**	100 Gms 500 Gms	1200 4000
47380 [61788-85-0]	**Polyoxyl 40 Hydrogenated Castor Oil extrapure** (Macr(Macrogolglycerol Hydroxystearate, PEG-40 hydrogenated castor oil)	500 Gms	2520
14749 [27416-86-0]	■ **Polyuridylic Acid Potassium Salt (POLY U Potassium) extrapure**	25 Mg	3000
38610 [6226-79-5]	**Ponceau S**	25 Gms 100 Gms	790 3050
14092	**Ponceau S Staining Solution**	100 ml	1750
41186 [68189-43-5]	**POPSO Buffer extrapure, 99%** (Piperazine-1,4-bis (2-Hydroxypropanesulphonic Acid)	25 Gms 100 Gms 500 Gms	2000 5000 20000
81800 [108321-08-0]	**POPSO Sesquisodium Salt Buffer extrapure, 98%**	25 Gms	6000
42486 [2720-73-2]	**Potassium Amylxanthate (PAX) pure, 97%** (Amylxanthic Acid Potassium Salt, Potassium O-Amyl Dithiocarbonate)	100 Gms 500 Gms	1450 4000
20883 [299-27-4]	**Potassium Gluconate extrapure, 99%** (D-Gluconic Acid Potassium Salt, Potassium-D-Gluconate)	100 Gms 500 Gms	300 1260
48320 [865-47-4]	**Potassium tert-Butoxide pure, 98%**	100 Gms 500 Gms	750 3200
95056 [24634-61-5]	**Potassium Sorbate pure, 99%**	500 Gms 1 Kg	900 1600
14884 [12142-33-5]	**Potassium Stannate pure, 37% Sn**	100 Gms 500 Gms	3000 9500
28008 [19423-77-9]	**Praseodymium (III) Chloride Hydrate extrapure, 99.9%**	25 Gms 100 Gms	6000 14000
83218 [15878-77-0]	**Praseodymium (III) Nitrate Hexahydrate extrapure, 99.9%** (Praseodymium Trinitrate Hexahydrate)	25 Gms 100 Gms	5000 16000
99350 [52093-27-3]	**Praseodymium (III) Trifluoromethanesulfonate extrapure, 98%** (Praseodymium Triflate, Trifluoromethanesulfonic Acid Praseodymium Salt)	5 Gms	12000
63639 [12037-29-5]	**Praseodymium (III, IV) Oxide extrapure, 99.9%**	10 Gms 50 Gms	1400 6000
10693 [57-83-0]	● **Progesterone extrapure, 97-103%** (4-Pregnene-3,20-Dione)	5 Gms 25 Gms 100 Gms	1310 4000 14000
54236 [21087-77-4]	**Progoitrin Potassium Salt Reference Standard Grade, 75%** (R)-2-Hydroxy-3-Butenylglucosinolate Potassium Salt)	10 Mg	30000
92825 [127277-53-6]	● **Prohexadione Calcium for tissue culture, 95%**	100 Mg 1 Gms 5 Gms	2000 5000 14000
31651 [72556-74-2]	**Phosphocholine Chloride Calcium Salt Tetrahydrate pure, 97%** (Calcium phosphorylcholine chloride, Choline chloride phosphate calcium salt, Choline phosphate chlo	5 Gms 25 Gms 100 Gms	2530 5420 15650
16110 [16652-71-4]	🌡 **L-Proline Benzyl Ester Hydrochloride extrapure, 99%**	5 Gms 25 Gms	2200 8000
16680 [72162-84-6]	■ **Proline Specific Endopeptidase ex. Flavobacterium Sp., 5U/mg powder**	50 Units	9800
31798 [107-19-7]	**Propargyl Alcohol pure, 99%**	250 ml 500 ml	1125 2100
90557 [62937-45-5]	**D-Prolinamide extrapure, 98%** ((R)-Pyrrolidine-2-Carboxamide)	5 Gms 25 Gms	4750 17790
35981 [7531-52-4]	**L-Prolinamide extrapure, 98%** ((2S)-2-Pyrrolidinecarboxamide, (-) -Prolinamide, (2S) -2-Carbamoylpyrrolidine)	5 Gms 25 Gms 100 Gms 500 Gms	900 2500 9000 42000
57732 [68832-13-3]	🌡 **D-Prolinol extrapure, 98%** ((R)-(-)-2-(Hydroxymethyl)pyrrolidine, (R)-(-)-2-Pyrrolidinemethanol)	1 ml 5 ml 25 ml	2050 8010 13790
46457 [23356-96-9]	🌡 **L-Prolinol extrapure, 98%** ((S)-(+)-2-(Hydroxymethyl)pyrrolidine, (S)-(+)-2-Pyrrolidinemethanol)	5 ml 25 ml	2000 5500
41622 [93-55-0]	**Propiophenone pure, 99%**	100 ml 500 ml	600 2000
84215 [106-94-5]	**n-Propyl Bromide pure, 98%**	100 ml 500 ml 2500 ml	480 1600 6200
45799 [108-32-7]	**Propylene Carbonate pure, 99%**	500 ml 2500 ml	600 2660
56082 [107-08-4]	**n-Propyl Iodide extrapure, 99%**	25 ml 100 ml 250 ml	2300 8400 20200

Code	Name	Packing	Price
25034 [106-36-5]	**Propyl Propionate extrapure, 99%**	250 ml 1000 ml	1130 4240
11195 [25535-16-4]	**Propidium Iodide (PI) for molecular biology, 95%** (3,8-Diamino-5-[3-(diethylmethylammonio)propyl] -6-Phenylphenanthridinium Diiodide)	10 Mg 25 Mg 100 Mg	2400 4200 12500
14887 [25535-16-4]	**Propidium Iodide Solution (1mg/ml in H2O)**	10 ml	3000
78349 [9009-65-8]	**Protamine Sulfate for molecular biology, 90-110%** (Protamine Sulphate)	1 Gms 5 Gms 10 Gms	3000 11000 20000
79546	**Protein A soluble**	2 Mg	12190
84589 [9036-06-0]	**Protease ex. AspergillusOryzae (Fungal) (Type 1 - Acidic), 500kHUT/g powder** (Subtilisin A)	25 Gms 100 Gms	4000 12000
38207 [9036-06-0]	**Protease ex. Bacillus Licheniformis (Type 4 - Neutral), 10U/mg powder** (Subtilisin A)	25 Mg 100 Mg	3890 12500
81525 [9036-06-0]	**Protease ex. Bacillus Sp. (Type 5 - Neutral, granules), 8kNPU/g** (Subtilisin A)	100 Gms 500 Gms	600 2500
66789 [9036-06-0]	**Protease ex. Bacillus Subtilis (Type 2 - Alkaline, granules), 100kAPU/g** (Subtilisin A)	250 Gms 1 Kg	5200 19000
23434 [9036-06-0]	**Protease ex. Bacillus Subtilis (Type 3 - Alkaline, granules), 1000kAPU/g** (Subtilisin A)	500 Gms	2190
84244 [9036-06-0]	**Protease ex. Bacillus Subtilis (Type 6 - Neutral), 100kNPU/g powder** (Subtilisin A)	100 Gms 500 Gms	8000 35000
49936 [39450-01-6]	**Proteinase K ex. Tritirachium Album (Type A) for molecular biology, 30U/mg protein** (Endopeptidase K)	25 Mg 100 Mg 500 Mg 1 Gms 5 Gms	1100 3600 13000 20000 85000
32753 [39450-01-6]	**Proteinase K ex. PichiaPastoris (Type D - Recombinant) for molecular biology & PCR, 30U/mg protein** (Endopeptidase K)	25 Mg 100 Mg 500 Mg 1 Gms 5 Gms	1300 4200 14200 22000 92000
45073 [39450-01-6]	**Proteinase K ex. PichiaPastoris (Type E - Recombinant) for molecular biology & PCR, 40U/mg protein** (Endopeptidase K)	25 Mg 100 Mg	1600 5000
36331	**Proteinase K Solution (20 mg/ml)**	1 ml 5 x1 ml 5 x5 ml	2400 5000 17000
21496 [9029-47-4]	**Protocatechuate-3,4-Dioxygenase (PCD) ex. Pseudomonas Sp., 3U/mg solid** (Protocatechuate-3,4-Dioxygenase)	5 Units 25 Units	4000 8000
66893 [9057-02-7]	**Pullulan (Polymaltotriose, PULL) extrapure**	5 Gms 25 Gms 100 Gms	2000 4000 13000
93909 [58-58-2]	**Puromycin Dihydrochloride (PRM), 98%**	10 Mg 25 Mg 100 Mg	3280 6400 18000
85492 [110-60-1]	**Putrescine free base extrapure, 99%** (1,4-Diaminobutane, 1,4-Butanediamine, Tetramethylenediamine)	25 Gms 100 Gms 500 Gms	1700 5300 25000
81343 [333-93-7]	**Putrescine Dihydrochloride extrapure, 99%** (1,4-Butanediamine Dihydrochloride, 1,4-Diaminobutane Dihydrochloride, Tetramethylenediamine Dihydro	5 Gms 25 Gms 100 Gms 500 Gms	2000 4100 12000 30000
49270 [333-93-7]	**Putrescine Dihydrochloride for tissue culture, 99%**	1 Gms 5 Gms 25 Gms	940 3280 8650
57322 [98-96-4]	**Pyrazinamide extrapure, 98%** (Pyrazinecarboxamide, Pyrazinamide, Pyrazinoic Acid Amide)	5 Gms 25 Gms	1500 3200
31412 [288-13-1]	**Pyrazole extrapure, 98%** (1,2-Diazole)	5 Gms 25 Gms 100 Gms	1000 2100 5600
32357	**Pyrene Actin (10%) ex. Rabbit Skeletal Muscle extrapure, 99%**	1 Mg	88920
39433 [67000-89-9]	**1-Pyrenebutanol extrapure, 99%**	500 Mg 1 Gms	5000 9400
81284 [19694-02-1]	**1-Pyrenecarboxylic Acid pure, 97%**	500 Mg 1 Gms 5 Gms	4390 7830 34390
59503 [93324-65-3]	**1-Pyrenemethylamine Hydrochloride (PyNH2) pure, 95%**	500 Mg 1 Gms	4090 7480
67172 [289-80-5]	**Pyridazine pure, 98%**	5 Gms 25 Gms	2730 12600
31129 [7291-22-7]	**Pyridine-d5 For NMR Spectroscopy, 99.5 Atom %D**	10 ml	24000
73512 [1692-25-7]	**3-Pyridineboronic Acid extrapure, 98%**	1 Gms 5 Gms 25 Gms	2040 8480 37300
16539 [499-83-2]	**Pyridine-2,6-Dicarboxylic Acid extrapure, 99%** (DPA, Dipicolinic Acid)	25 Gms 100 Gms 500 Gms	850 2000 8000
12644 [636-73-7]	**Pyridine-3-Sulphonic Acid (P3SA) extrapure, 98%** (3-Pyridinesulfonic Acid)	25 Gms	2000
52928 [41468-25-1]	**Pyridoxal-5-Phosphate Monohydrate extrapure, 98%**	1 Gms 5 Gms 25 Gms	1020 3000 9000
88205 [58-56-0]	**Pyridoxine Hydrochloride (Vitamin B6) pure, 98%**	25 Gms 100 Gms 500 Gms	600 1600 6000

Code	Name	Packing	Price
47330 ✔ [1046-56-6]	**3(2-Pyridyl) -5,6-Diphenyl-1,2,4-Triazine (PDT) extrapure AR, 99%**	1 Gms 5 Gms	2400 11000
16315 🌡 [68181-17-9]	**3-(2-Pyridyldithio)propionic Acid N-Hydroxysuccinimide Ester (SPDP Crosslinker) extrapure, 95%** (N-Succinimidyl 3-(2-Pyridyldithio)Propionate, 3-(2-Pyridyldithio)-Propionic Acid-OSu)	100 Mg 1 Gms 5 Gms	1300 4200 16000
49970 [4042-36-8]	**D-Pyroglutamic Acid extrapure, 98%** (D-Glutamic Acid Lactam, (R)-(+) -2-Pyrrolidinone-5-Carboxylic Acid)	5 Gms 25 Gms 100 Gms	3000 5000 9000
77577 🌡 [98-79-3]	**L-Pyroglutamic Acid extrapure, 99%**	100 Gms 500 Gms	800 2000
36306 🌡 [66642-36-2]	**L-Pyroglutamic Acid-7-Amido-4-Methylcoumarin extrapure, 98%** (L-Pyroglutamic Acid-AMC)	10 Mg 25 Mg	10010 19450
97079 [89-32-7]	**Pyromellitic Dianhydride (PMDA) extrapure, 98%** (Benzene-1,2,4,5-Tetracarboxylic Dianhydride)	25 Gms 100 Gms 500 Gms	1370 2210 7040
75612 🌡 [92-32-0]	**Pyronin G (Pyronin Y), 50%** (C.I. No. 45005)	1 Gms 5 Gms 25 Gms	940 3350 16000
71692 [109-97-7]	**Pyrrole extrapure, 99%**	25 ml 100 ml 500 ml	700 2200 9000
91402 [123-75-1]	**Pyrrolidine pure, 99%**	25 ml 100 ml 500 ml	280 610 2750
88543 [123-75-1]	**Pyrrolidine Reagent Grade, 99.5%**	25 ml 100 ml 250 ml	360 880 1760
73694 [616-45-5]	**2-Pyrrolidone pure, 99%**	250 ml 500 ml 2500 ml 25 Ltr	550 940 3960 35200
13558 [72909-34-3]	**Pyrroloquinoline Quinone (PQQ) extrapure, 98%** (4,5-Dihydro-4,5-Dioxo-1H-Pyrrolo[2,3-f]quinoline-2,7,9-Tricarboxylic Acid, Methoxatin)	100 Mg 1 Gms 5 Gms	2260 5060 11550
56826 🌡 [122628-50-6]	**Pyrroloquinoline Quinone Disodium Salt (PQQ Disodium) extrapure, 98%** (4,5-Dihydro-4,5-Dioxo-1H-Pyrrolo[2,3-f]quinoline-2,7,9-Tricarboxylic Acid Disodium, Methoxatin Diso	100 Mg 1 Gms 5 Gms	3300 6600 13200
63884 ■ [9014-20-4]	**Pyruvate Dehydrogenase ex. Microorganism, 2.0U/mg powder**	10 Units	11360
29986 ■ [9001-59-6]	**Pyruvate Kinase ex. Rabbit Muscle, 300U/mg powder**	2000 Units 4000 Units	6600 12650
59150 🌡 [127-17-3]	**Pyruvic Acid pure, 98%**	25 ml 100 ml	990 2640

Code	Name	Packing	Price
71923 🌡 [6151-25-3, 849061-97-8]	**Quercetin Dihydrate extrapure, 98%**	25 Gms 100 Gms 500 Gms	2200 5500 24200
97899 🌡 [117-92-0]	**Quinaldine Red**	1 Gms 5 Gms	1580 7320
81076 🌡 [77-95-2]	**Quinic Acid extrapure, 98%**	5 Gms 25 Gms 100 Gms	1250 3920 13180
32648 [91-22-5]	**Quinoline pure, 97%**	500 ml 2500 ml	2030 9680
51649 [17629-30-0]	**D-Raffinose Pentahydrate extrapure**	25 Gms 100 Gms	1980 7480
21719 ● [212783-75-0]	**Ramage Linker (FMOC-Suberol) extrapure, 98%**	1 Gms 5 Gms	4400 14300
24527 🌡	**Rapeseed Glucosinolate Mixture Reference Standard Grade, 90%**	10 Mg	10400
36474 [2580-78-1]	**Reactive Blue 19 (Remazol Brilliant Blue R)** (Reactive Blue 19)	25 Gms 100 Gms	2490 8500
95238 [12236-86-1]	**Reactive Blue 21 (Remazol Turquoise Blue G-133)**	25 Gms 100 Gms	800 2500
12176 [12225-83-1]	**Reactive Orange 16 (Remazol Brilliant Orange 3R)** (C.I. No. 17757)	25 Gms 100 Gms	2000 4000
65637 [12220-08-5]	**Reactive Orange 107 (Remazol Golden Yellow RNL)**	25 Gms 100 Gms 500 Gms	800 3000 8000
15034 [11099-79-9]	**Reactive Red 21 (Remazol Red BB)**	25 Gms 100 Gms	1200 3200
35832 [12226-12-9]	**Reactive Red 35 (Remazol Brilliant Red 5B)**	25 Gms 100 Gms	600 2000
79740 [61951-82-4]	**Reactive Red 120** (Reactive Brilliant Red HE3B) (C.I. 292775)	25 Gms 100 Gms	1400 4000
23176 [98114-32-0]	**Reactive Red 180 (Remazol Brilliant Red F3B)**	5 Gms 25 Gms 100 Gms	700 1400 4000
80539 [12226-38-9]	**Reactive Violet 5 (Remazol Brilliant Violet 5R)** (C.I. No. 18097)	25 Gms 100 Gms	1400 4000
67074 [20317-19-5]	**Reactive Yellow 17 (Remazol Golden Yellow G)** (C.I. No. 18852)	25 Gms 100 Gms	2000 6000
55200 [62758-13-8]	**Resazurin Sodium extrapure AR, 20-30%**	1 Gms 5 Gms	1200 5500
42650 [62758-13-8]	**Resazurin Sodium extrapure AR, 85%**	1 Gms 5 Gms 25 Gms	1400 6200 29000
42817 [62758-13-8]	**Resazurin Sodium extrapure AR, ExiPlus, Multi-Compendial, 85%** meets compendial specs of USP	1 Gms 5 Gms 25 Gms	1500 6500 31000

Code	Name	Packing	Price
46861 ✔ [50-55-5]	**Reserpine extrapure, 99%**	1 Gms 5 Gms	3800 12650
39798 [108-46-3]	**Resorcinol ACS, 99%**	100 Gms 500 Gms	550 1980
22173 [108-46-3]	**Resorcinol ExiPlus, Multi-Compendial, 99%** meets compendial specs of USP, BP, Ph.Eur	100 Gms 500 Gms	610 2200
69040 [10030-85-0]	**L-Rhamnose Monohydrate extrapure**	5 Gms 25 Gms 100 Gms	650 1500 5000
59652 [7440-15-5]	**Rhenium Powder extrapure, 99.99%, -325mesh**	1 Gms 5 Gms	1000 3500
59774 □ [15956-28-2]	**Rhodium (II) Acetate Dimer extrapure, 45-47% Rh** (Dirhodium Tetraacetate, Tetrakis (Acetato) Dirhodium (II))	50 Mg 250 Mg	6000 21000
45165 [20765-98-4]	**Rhodium (III) Chloride Trihydrate pure, 98%** (Rhodium Trichloride Hydrate)	250 Mg 1 Gms	11000 43000
46729 ✔ [73482-96-9]	**Rhodium (II) Octanoate Dimer extrapure, 97%**	250 Mg 1 Gms	7500 28000
34392 [83-88-5]	**Riboflavin pure, 98%** (Vitamin B2)	25 Gms 100 Gms	600 2000
57443 🌡 [130-40-5]	**Riboflavine-5-Monophosphate Sodium Salt Hydrate extrapure, 73-79%**	25 Gms 100 Gms	740 2690
98915 ■ [9001-99-4]	**Ribonuclease A (RNase A) ex. Bovine Pancreas for molecular biology, 50Kunitz/mg**	25 Mg 100 Mg 500 Mg	1200 3500 15000
62910 ■ [63231-63-0]	**Ribonucleic Acid (RNA) Free Acid ex. Yeast, 85%**	25 Gms 100 Gms	1400 4000
61866 ■ [73049-77-1]	**Ribonucleic Acid Sodium Salt (RNA Sodium, Sodium Nucleinate) ex. Yeast for molecular biology, 85%**	10 Gms 25 Gms 100 Gms	950 1800 5200
48259 🌡 [50-69-1]	**D-Ribose extrapure**	25 Gms 100 Gms 1 Kg	550 900 7000
76538 🌡 [50-69-1]	**D-Ribose for tissue culture**	25 Gms 100 Gms	600 1100
72377 [15673-79-7]	**D-Ribose-5-Phosphate Barium Salt Hexahydrate extrapure, 98%**	100 Mg	4600
12310 🌡 [13292-46-1]	**Rifampicin (RFP) (Rifampin), 97-102%**	1 Gms 5 Gms	720 2400
81319 ■ [80621-81-4]	**Rifaximin (RFX) extrapure, 98%**	1 Gms 5 Gms	1800 7000
23311 🌡 [145069-56-3]	**Rink Amide Linker Resin (Knorr Linker) extrapure, 99%** (4-[(2,4-Dimethoxyphenyl) (Fmoc-amino)Methyl]phenoxyacetic Acid)	1 Gms 5 Gms 25 Gms	1400 4200 15000
34140 🌡 [431041-83-7]	**Rink Amide MBHA Resin (100-200 mesh), Loading 0.3-0.8mmol/g** (4-(2,4-Dimethoxyphenyl-FMOC-Aminomethyl) Phenoxyacetamido-Norleucyl-MBHA Resin)	1 Gms 5 Gms 25 Gms	2000 7000 30000
90780 [569-61-9]	**p-Rosaniline Hydrochloride** (Pararosaniline Hydrochloride, Pararosaniline Chloride, Basic Parafuchsin) (C.I. No. 42500)	25 Gms 100 Gms 500 Gms	900 3000 13000
60549 [632-69-9]	**Rose Bengal** (Acid Red 94) (C.I. No. 45440)	10 Gms 25 Gms 100 Gms	1100 2200 8200
10641 [8007-01-0]	**Rose Oil extrapure**	10 ml 25 ml 100 ml	1400 2600 8000
66412 [603-45-2]	**p-Rosolic Acid** (Aurin)	25 Gms 100 Gms 500 Gms	900 3000 13000
23916 [14284-93-6]	**Ruthenium AcetylAcetonate extrapure, 24% Ru** (Ruthenium(III) 2,4-Pentanedionate, Tris(2,4-Pentanedionato)Ruthenium(III))	1 Gms 5 Gms	7800 36000
55156 □ [10049-08-8]	**Ruthenium (III) Chloride Anhydrous extrapure, 45-55% Ru** (Ruthenium Trichloride)	1 Gms 5 Gms	14000 60000
18325 [250249-75-3]	**Rutin Trihydrate pure, 95%** (Vitamin P Trihydrate, Rutoside, Quercetin-3-Rutinoside)	25 Gms 100 Gms 500 Gms	900 1800 7000
16281	**Safranine 0.5 w/v Aq. Solution**	125 ml 500 ml	220 700
16593 [477-73-6]	**Safranine O** (Basic red 2) (C.I. No. 50240)	25 Gms 100 Gms 500 Gms	650 2400 11000
72335 [477-73-6]	**Safranine O PrimeDx**	1 Kg	23000
10673 [138-52-3]	**D-(-)-Salicin extrapure**	5 Gms 25 Gms	2020 5380
57467 ✔ [959-36-4]	**Salicylaldehyde Azine extrapure, 97%** (2-Hydroxybenzaldehyde azine)	1 Gms 5 Gms 25 Gms	950 2000 7500
77386 [65-45-2]	**Salicylamide pure, 98%**	500 Gms 5 Kg	1550 14000
45549	**Salt Mixture Bernhart Tommarelli** (Dietary Mineral Mix acc to Bernhart-Tommarelli)	500 Gms 1 Kg 5 Kg	1400 2690 12880
13464	**Salt Mixture H.M.W. (acc. to Hubble Mendal and Wakeman)** (Dietary Mineral Mix acc to HMW)	500 Gms 1 Kg 5 Kg	1400 2690 12880
71637	**Salt Mixture as per USP XIV** (Dietary Mineral Mix acc to USP XIV)	500 Gms 1 Kg 5 Kg	1680 3200 15120

Code	Name	Packing	Price
93538 [13465-55-9]	**Samarium (III) Chloride Hexahydrate pure, 99.9%** (Samarium trichloride hydrate)	10 Gms 50 Gms	5000 18000
92082 [13759-83-6]	**Samarium (III) Nitrate Hexahydrate ultrapure, 99.95%** (Trinitric Acid Samarium (III) Hexahydrate)	5 Gms 25 Gms	5000 18000
31372 [12060-58-1]	**Samarium (III) Oxide ultrapure, 99.95%** (Samaria)	5 Gms 25 Gms	5000 10000
54043 [52093-28-4]	**Samarium (III) Trifluoromethanesulfonate extrapure, 98%** (Samarium (III) Triflate, Trifluoromethanesulfonic Acid Samarium (III) Salt)	1 Gms 5 Gms	4000 10000
68580 [144026-79-9]	**Scandium (III) Trifluoromethanesulfonate extrapure, 98%** (Scandium Triflate, Trifluoromethanesulfonic Acid Scandium Salt)	1 Gms 5 Gms	4000 10000
80967 [8047-15-2]	**Saponin ex. Quillaja for molecular biology, 9% Sapogenin**	25 Gms 100 Gms	6000 11000
99474 [8047-15-2]	**Saponin ex. Gypsophila Roots for molecular biology, 25% Sapogenin**	5 Gms 25 Gms 100 Gms	1800 4500 12000
10147 ■ [9029-22-5]	**Sarcosine Oxidase (SOX) ex. Recombinant E.coli, 20U/mg powder**	500 Units 1 K.Units	7000 13000
21032 [12060-08-1]	**Scandium Oxide, 99.99%**	1 Gms 5 Gms	2630 7000
17386	**Schiff's Reagent**	500 ml	530
32380 [92-61-5]	**Scopoletin extrapure, 98%**	25 Mg 100 Mg 1 Gms	3000 7000 22000
96843 ✔ [111-19-3]	**Sebacoyl Chloride pure, 95%** (Decanedioyl Chloride Sebacic Acid Dichloride, Sebacyl chloride)	100 ml 500 ml	4000 12000
62330 [111-20-6]	**Sebacic Acid pure, 98%** (Decanedioic acid)	500 Gms	1080
19182 [563-41-7]	**Semicarbazide Hydrochloride extrapure, 98%**	100 Gms 500 Gms	670 3100
44779 [7446-08-4]	**Selenium Dioxide pure, 98%**	100 Gms 500 Gms	2000 8000
90945 [7782-49-2]	**Selenium Metal Powder extrapure, 99.9%, 100 mesh**	25 Gms 100 Gms 500 Gms	700 2500 7500
51200 [7782-49-2]	**Selenium Metal Powder ultrapure, 99.999%, 100 mesh**	25 Gms 100 Gms	3570 11000
77855 [7782-49-2]	**Selenium Metal Granular ultrapure, 99.999%, 1-6mm**	25 Gms 100 Gms	2000 5000
29815	**Seliwanoff`s Reagent**	100 ml	290
64088 🌡	**Seralose® 2B (2% Beaded Agarose Susp., 60-200micron)**	100 ml 500 ml	4500 18000
54540 🌡	**Seralose® 4B (4% Beaded Agarose Susp., 45-165micron)**	100 ml 500 ml	4500 18000
12585 🌡	**Seralose® 6B (6% Beaded Agarose Susp., 45-165micron)**	100 ml 500 ml	4500 18000
81862 🌡	**Seralose® CL-2B (2% Beaded Agarose Cross-Linked Susp., 60-200micron)**	100 ml 500 ml	6000 24000
11974 🌡	**Seralose® CL-4B (4% Beaded Agarose Cross-Linked Susp., 45-165micron)**	100 ml 500 ml	6000 24000
22264 🌡	**Seralose® CL-6B (6% Beaded Agarose Cross-Linked Susp., 45-165micron)**	100 ml 500 ml	6000 24000
85165 🌡 [55819-71-1]	**DL-Serine Hydrazide Hydrochloride pure, 98%** (DL-Serinohydrazide Hydrochloride)	5 Gms	15690
17017 🌡 [5680-80-8]	**L-Serine Methyl Ester Hydrochloride extrapure, 99%**	25 Gms 100 Gms 500 Gms	1800 4000 16000
18472 🌡 [3722-51-8]	**Sieber Linker (3-Hydroxyxanthone) extrapure, 99%**	1 Gms 5 Gms 25 Gms	2000 5000 12000
77996 🌡	**Silica Cysteine-Functionalized extrapure** (Particle Size — 40-63µm, Pore Size - 60A)	5 Gms 10 Gms	19970 33270
46400 🌡	**Silica DMT-Functionalized extrapure**	5 Gms 25 Gms	11230 46740
31156 🌡	**Silica Imidazole-Functionalized extrapure**	5 Gms	3000
39343 🌡	**Silica Thiol-Functionalized (STF) extrapure**	5 Gms 25 Gms	8030 33410
45570 🌡	**Silica TEMPO Catalyst extrapure**	5 Gms 10 Gms	18510 30850
78339 🌡	**Silica C18 (17%) Bonded monomeric extrapure**	5 Gms 25 Gms	5260 15190
43139 🌡	**Silica Carbonate Bonded extrapure**	5 Gms 25 Gms	9260 38570
35228 🌡	**Silica Cyano Bonded extrapure**	5 Gms 10 Gms	5870 11450
30753 🌡	**Silica Diol Bonded extrapure**	5 Gms 25 Gms	5780 24040
92868 🌡	**Silica TMA Acetate Bonded non-endcapped extrapure**	5 Gms 25 Gms	8240 34300
79427 [7440-21-3]	**Silicon Metal powder, 98.5%, -4 mesh**	500 Gms	1300
73081 [7440-21-3]	**Silicon Metal powder, 98.5%, -200 mesh**	100 Gms 500 Gms 5 Kg	320 980 9450

Code	Name	Packing	Price
90677 [16411-05-5]	**Sinalbin Potassium Salt Reference Standard Grade, 98%** (p-Hydroxybenzylglucosinolate Potassium Salt, Glucosinalbin Potassium Salt)	10 Mg	33600
28945 [3952-98-5, 64550-88-5]	**Sinigrin Potassium Salt Hydrate Reference Standard Grade, 97%** (2-Propenyl Glucosinolate Potassium Salt, Sinigroside Potassium Salt)	25 Mg 100 Mg	6300 10290
94846 [122-34-9]	**Simazine extrapure,98%** (2,4-Bis(ethylamino) -6-chloro-1,3,5-triazine)	5 Gms 25 Gms	2500 4000
56226 [83-46-5]	**b-Sitosterol ex. Rice Bran, 95%** (Beta Sitosterol)	5 Gms 25 Gms	2000 9400
13728 [83-46-5]	**b-Sitosterol ex. Glycine Max (Soyabean), 95%** (Beta Sitosterol)	25 Gms 100 Gms	4440 17410
89658 [1333-83-1]	**Sodium Bifluoride pure, 98%**	500 Gms	1430
53713 [7782-82-3]	**Sodium Biselenite pure, 98%**	100 Gms	4330
22308 [156-54-7]	**Sodium Butyrate Reagent Grade, 98.5%**	25 Gms 100 Gms 500 Gms	1100 1800 6800
41059 [9005-46-3]	**Sodium Caseinate extrapure** (Casein Sodium Salt)	100 Gms 500 Gms	1200 5400
43624 [7775-11-3]	**Sodium Chromate pure, 98%**	500 Gms	640
17296 [51580-86-0]	**Sodium Dichloroisocyanurate Dihydrate extrapure, 98%** (Dichloroisocyanuric Acid Sodium Salt Dihydrate, NaDCC)	100 Gms 500 Gms	600 1500
60947 [22722-98-1]	**Sodium Dihydro-bis-Methoxyethoxy Aluminate (VITRIDE) pure, 68-71%**	100 Gms 500 Gms	2000 9250
46288 [6152-67-6]	**Sodium Diphenylamine-4-Sulfonate ACS**	10 Gms 25 Gms 100 Gms	420 850 3000
25818 [5137-46-2]	**Sodium Diphenyl Reagent AR** (Sodium Biphenyl Complex)	15 ml 10 x15 ml	960 8900
67516 [5137-46-2]	**Sodium Diphenyl Reagent ExiPlus, Multi-Compendial** (Sodium Biphenyl Complex)	15 ml 10 x15 ml	1430 10440
55083 [25155-30-0]	**Sodium Dodecylbenzenesulfonate (SDBS) pure, 80%** (Dodecylbenzenesulfonic Acid Sodium Salt)	100 Gms 500 Gms 5 Kg	450 825 7200
23819 [126-92-1]	**Sodium 2-Ethylhexyl Sulfate 40% aq. Solution** (2-Ethylhexyl Sulfate Sodium Salt)	500 ml 1000 ml	1500 2700
36511 [751-94-0]	**Sodium Fucidate (Na-FC), 97.5%** (Fucidic Acid Sodium Salt)	1 Gms 5 Gms	3000 6000
55655 [17013-01-3]	**Sodium Fumarate Dibasic extrapure AR, 99.5%** (Disodium Fumarate)	100 Gms 500 Gms	900 4000
15331 [55073-41-1, 154804-51-0]	**Sodium-ß-Glycerophosphate Hydrate extrapure AR, 99%** (Sodium-ß-Glycerol Phosphate Hydrate, ß-GP, ß-Glycerolphosphate, Disodium ß-glycerol phosphate)	5 Gms 25 Gms 100 Gms	2000 4500 15000
97971 [16564-43-5]	**Sodium Glycochenodeoxycholate (SGCDC) extrapure, 97%** (Glycochenodeoxycholic Acid Sodium Salt)	1 Gms 5 Gms 25 Gms 100 Gms	8000 25000 105000 280000
31287 [863-57-0]	**Sodium Glycocholate Dried (SGC Dried) extrapure, 98%**	500 Mg 1 Gms 5 Gms 25 Gms	3690 5500 17000 80000
15556 [338950-81-5]	**Sodium Glycocholate Hydrate (SGC Hydrate) extrapure, 98%**	1 Gms 5 Gms 25 Gms 100 Gms	4000 15000 68000 180000
43453 [527-07-1]	**Sodium Gluconate pure, 99%**	500 Gms 5 Kg	260 2400
73281 [140-93-2]	**Sodium Isopropyl Xanthate (SIPX) pure, 90%** (Sodium Isopropoxy(Thioxo)Methanethiolate Proxan Sodium)	100 Gms 500 Gms	1150 3700
54565 [141-95-7]	**Sodium Malonate extrapure, 99%**	100 Gms 500 Gms	900 4200
81020 [114-21-6]	**Sodium DL-Mandelate extrapure, 99%** (DL-Mandelic Acid Sodium Salt)	5 Gms 25 Gms	900 4200
18444 [13718-26-8]	**Sodium Metavanadate extrapure, 99%** (Sodium Trioxovanadate, Sodium Vanadium Oxide, Sodium Vanadium Trioxide)	25 Gms 100 Gms 500 Gms	950 2400 10000
85247 [31681-98-8]	**Sodium-2-Naphthyl Phosphate Disodium Salt extrapure AR, 99%**	1 Gms 5 Gms	980 4740
86766 [13721-39-6]	**Sodium Orthovanadate pure, 45%** (Sodium Vanadium Oxide)	10 Gms 50 Gms 250 Gms	2000 8000 18000
14215 [55963-78-5]	**Polyanetholesulphonic Acid Sodium Salt (Sodium Polyanetholesulphonate, SPS) extrapure** (Sodium PolyanetholSulphonate)	5 Gms 25 Gms	3000 8000
23569 [113-24-6]	**Sodium Pyruvate (Pyruvic Acid Sodium Salt) extrapure AR, 99%**	25 Gms 100 Gms 500 Gms	400 800 2600
19099 [523-21-7]	**Sodium Rhodizonate Dibasic extrapure AR, 99%** Rhodizonic Acid Sodium Salt, 3,4,5,6-Tetraoxocyclohexene-1,2-diol disodium salt, Rhodizonic acid dis	1 Gms 5 Gms 25 Gms	580 1380 6340
95335 [10102-18-8]	**Sodium Selenite Anhydrous extrapure, 98%**	25 Gms 100 Gms 500 Gms	780 2610 10210

Code	Name	Packing	Price
80451 [10102-18-8]	**Sodium Selenite Anhydrous extrapure AR, 99%**	25 Gms	820
		100 Gms	2950
		500 Gms	11120
17169 [10102-18-8]	**Sodium Selenite Anhydrous for tissue culture, 99%**	25 Gms	960
		100 Gms	2840
		500 Gms	14070
58529 [6009-98-9]	**Sodium Taurochenodeoxycholate (STCDC) extrapure, 95%** (Taurochenodeoxycholic Acid Sodium Salt)	250 Mg	5000
		1 Gms	12000
63703 [345909-26-4]	**Sodium Taurocholate Hydrate (STC Hydrate) extrapure, 98%** (Taurocholic Acid Sodium Salt Hydrate)	1 Gms	2370
		5 Gms	6200
		25 Gms	27000
		100 Gms	85000
57861 [1180-95-6]	**Sodium Taurodeoxycholate Dried (STDC Dried) extrapure, 95%**	1 Gms	3500
		5 Gms	10000
		25 Gms	42000
40391 [207737-97-1]	**Sodium Taurodeoxycholate Hydrate (STDC Hydrate) extrapure, 95%**	1 Gms	2500
		5 Gms	8000
		25 Gms	34000
		100 Gms	120000
31102 [35807-85-3]	**Sodium Tauroursodeoxycholate (TUDCA Sodium) extrapure, 98%** (Tauroursodeoxycholic Acid Sodium Salt)	1 Gms	3000
		5 Gms	12000
78786 [367-51-1]	**Sodium Thioglycolate Reagent Grade, 99%** (Sodium Mercaptoacetate)	100 Gms	1400
		250 Gms	3100
		500 Gms	6000
48950 [56553-60-7]	**Sodium Triacetoxyborohydride (STAB) pure, 98%**	25 Gms	600
		100 Gms	1800
		500 Gms	8000
69128 [2926-29-6]	**Sodium Trifluoromethanesulfinate (Langlois reagent) extrapure, 95%** (Sodium Trifluoromethylsulfinate, Trifluoromethanesulfinic Acid Sodium Salt, Sodium Triflinate)	5 Gms	3800
		25 Gms	12000
99957 [20311-51-7]	**a-Solamargine extrapure, 98%** (Solasodine-L-Rhamnosido-D-Glucoside)	10 Mg	4500
		25 Mg	9000
53771 [126-17-0]	**Solasodine extrapure, 98%** (Purapuridine, Solancarpidine)	10 Mg	4500
		25 Mg	9000
89725 [19121-58-5]	**Solasonine extrapure, 98%** (Solasodamine)	25 Mg	9000
54281 [38916-34-6]	**Somatostatin extrapure, 97%**	5 Mg	2000
		25 Mg	8000
62468 [87-79-6]	**L-(-)-Sorbose extrapure, 99%**	5 Gms	2000
		25 Gms	4500
		100 Gms	16000
51719 [22189-32-8]	**Spectinomycin Dihydrochloride Pentahydrate extrapure, 98%**	1 Gms	1200
		5 Gms	4000
		25 Gms	13000
17030 [124-20-9]	**Spermidine extrapure, 99%**	1 Gms	4500
		5 Gms	13000
		25 Gms	55000
91710 [71-44-3]	**Spermine Free Base extrapure, 99%**	250 Mg	2000
		1 Gms	3200
		5 Gms	12000
22310 [306-67-2]	**Spermine Tetrahydrochloride extrapure, 99%**	1 Gms	3460
		5 Gms	8980
		10 Gms	17270
		25 Gms	31080
26413 [85187-10-6]	**Sphingomyelin (Type I) ex. Bovine Brain Sphingolipids extrapure, 99%**	10 Mg	6000
		100 Mg	25000
95852 [85187-10-6]	**Sphingomyelin (Type II) ex. Bovine Heart Sphingolipids extrapure, 99%**	10 Mg	6500
		100 Mg	28000
77061 [657408-07-6]	**S-Phos (2-Dicyclohexylphosphino-2,6-Dimethoxybiphenyl) extrapure, 98%**	1 Gms	1250
		5 Gms	2400
		25 Gms	10000
63056 [122341-56-4]	**3-(2`-Spiroadamantane)-4-Methoxy-4-(3```-Phosphoryloxy) Phenyl-1,2-Dioxetane (AMPPD) extrapure, 98%**	5 Mg	20920
67715 [52-01-7]	**Spironolactone extrapure, 97%**	1 Gms	3600
		5 Gms	9800
12590	**20X SSC Buffer pH 6.9-7.1 for molecular biology**	200 ml	330
		500 ml	740
78305	**20X SSPE Buffer pH-7.4 suitable for molecular biology**	200 ml	900
		1000 ml	2500
80268 [10094-58-3]	**StachyoseTetrahydrate extrapure, 98%**	25 Mg	3150
		100 Mg	12450
		250 Mg	28690
31384 [57-11-4]	**Stearic Acid pure, C18-98%**	500 Gms	1700
		5 Kg	15000
73878 [57-11-4]	**Stearic Acid extrapure, C18-99%**	5 Gms	600
		25 Gms	1400
		100 Gms	5000
73439 [57-11-4]	**Stearic Acid ExiPlus, Multi-Compendial, C18-99%** meets compendial specs of Ph.Eur	5 Gms	510
		25 Gms	2200
		100 Gms	7000
61609 [112-61-8]	**Stearic Acid Methyl Ester Reference Standard, 99% (GC)** (Methyl Stearate)	5 Gms	3500
		25 Gms	9800
61600	**STET Lysis Buffer**	500 ml	3000
39895 [83-48-7]	**Stigmasterol extrapure, 95%**	1 Gms	1170
		5 Gms	3500
		25 Gms	12000
87610 [9013-20-1]	**Streptavidin (SA) ex. Streptomyces Avidinii for molecular biology, 15U/mg Protein**	5 Mg	8120
		10 Mg	14930

Code	Name	Packing	Price
91014 [3810-74-0]	**Streptomycin Sulphate (STM)**	5 Gms	410
		25 Gms	1350
		100 Gms	2810
62463 [3810-74-0]	**Streptomycin Sulphate for tissue culture**	5 Gms	930
		25 Gms	1860
		100 Gms	3380
14653 [18883-66-4]	**Streptozotocin (STZ) extrapure, 98%**	100 Mg	1850
		500 Mg	8500
		1 Gms	16000
		5 Gms	65000
28065 [100-42-5]	**Styrene Monomer (Stab w/ tBC) extapure, 99.5%** (Phenylethylene, Vinylbenzene)	500 ml	950
		2500 ml	4400
48142 [108-30-5]	**Succinic Anhydride extrapure, 99.5%**	100 Gms	380
		500 Gms	1500
72368 [115223-09-1]	**N-Succinimidyl Ferrocenecarboxylate pure, 97%** (Ferrocenecarboxylic Acid N-Succinimidyl Ester)	1 Gms	4000
		5 Gms	18000
24088 [64987-85-5]	**Succinimidyl-4-(N-(N-maleimidomethyl) cyclohexane-1-carboxylate (SMCC), 98%**	250 Mg	3200
		1 Gms	8000
		5 Gms	20000
20721 [386229-69-2]	**Sucrose Heptasulphate Potassium Salt pure, 80%** (1,3,4,6-Tetra-O-Sulfo-b-D-Fructofuranosyl-a-DGlucopyranoside-2,3,4-Tris(Hydrogen Sulfate)Potassium)	5 Mg	24000
		25 Mg	85000
47973 [54182-58-0]	**Sucrose Octasulphate Aluminium Salt (Sucralfate) extrapure, 30-36%** (Sucrose Octasulfate-Aluminum Complex)	5 Gms	3450
		25 Gms	9200
37654 [73264-44-5]	**Sucrose Octasulphate Potassium Salt extrapure, 95%** (Sucrose Octasulfate-Potassium Complex, Potassium Sucrose Octasulfate)	5 Gms	8000
		25 Gms	37500
69538 [26446-38-8]	**Sucrose Palmitate extrapure, 90%** (b-D-Fructofuranosyl-a-D-Glucopyranoside Monohexadecanoate)	25 Gms	2200
		100 Gms	6700
67787 [85-86-9]	**Sudan III**	100 Gms	510
		500 Gms	2200
81329 [6209-17-2]	**Sulfacetamide Sodium Monohydrate extrapure, 99%-101%** (N-Acetylsulfanilamide Sodium Salt Monohydrate)	5 Gms	2200
		50 Gms	4000
16775 [68-35-9]	**Sulfadiazine (SFD) extrapure, 99-101%** (4-Amino-N-(2-Pyrimidinyl)Benzenesulfonamide)	50 Gms	2630
		100 Gms	4990
28277 [547-32-0]	**Sulfadiazine Sodium Salt (SFD Sodium) extrapure, 99%** (4-Amino-N-(2-Pyrimidinyl)Benzenesulfonamide Sodium Salt)	25 Gms	2900
		100 Gms	9100
69181 [7704-34-9]	**Sulfur Precipitated pure, 98%**	500 Gms	230
		5 Kg	2600
86684 [723-46-6]	**Sulphamethoxazole (SMX), 99%**	5 Gms	2350
		25 Gms	6410
72701 [126-33-0]	**Sulpholane Anhydrous pure, 99%** (Sulfolane Anhydrous)	500 ml	800
		2500 ml	3700
13674 [126-33-0]	**Sulpholane Anhydrous extrapure AR, 99%** (Sulfolane Anhydrous)	500 ml	1025
		2500 ml	5050
31249 [9054-89-1]	**Superoxide Dismutase (SOD) ex. Bovine Erythrocytes, 3000U/mg powder**	20 K.Units	12450
		100 K.Units	40310
47406 [14414-32-5]	**Syringaldazine pure, 95%**	1 Gms	4500
41719 [134-96-3]	**Syringaldehyde extrapure, 99%**	5 Gms	550
		25 Gms	1500
		100 Gms	5500
64614 [530-57-4]	**Syringic Acid pure, 96%** (3,5-Dimethoxy-4-Hydroxybenzoic Acid, Gallic Acid 3,5-Dimethyl Ether)	5 Gms	1600
		25 Gms	3000
		100 Gms	9000
71989	**10X TAE Buffer pH-8.3 for molecular biology**	100 ml	590
		500 ml	1040
83170	**10X TBE Buffer pH-8.3 for molecular biology**	200 ml	660
		500 ml	1130
51782	**10X TE Buffer pH-8.0 for molecular biology**	100 ml	1260
		500 ml	4200
89528 [7440-25-7]	**Tantalum Metal Sheet (Type 1), 99.9%**	1 piece	5000
10964 [29915-38-6]	**TAPS Buffer extrapure, 99%** (N-Tris (hydroxymethyl) Methyl-3-Aminopropane Sulphonic Acid)	25 Gms	1370
		100 Gms	3620
44194 [68399-81-5]	**TAPSO Buffer extrapure, 99%** (3-[Tris-(hydroxymethyl) Methylamino] -2- Hydroxypropane Sulphonic Acid)	25 Gms	2810
		100 Gms	9290
28682 [147-71-7]	**D-(-)-Tartaric Acid extrapure, 99%**	100 Gms	630
		250 Gms	1510
60301 [87-69-4]	**L-(+)-Tartaric Acid ExiPlus, Multi-Compendial, 99.7%** (L-Threaric Acid) meets compendial specs of USP, BP, Ph.Eur	500 Gms	690
		5 Kg	6210
63904 [133-37-9]	**DL-Tartaric Acid extrapure, 99%**	500 Gms	610
		5 Kg	5750
15711 [107-35-7]	**Taurine extrapure CHR, 99%** (2-Aminoethanesulfonic Acid)	100 Gms	700
		1 Kg	5000
		5 Kg	20000
47064 [125700-67-6]	**TBTU extrapure, 98%**	25 Gms	2000
		100 Gms	6000
		500 Gms	24000
92481 [7446-07-3]	**Tellurium Dioxide ultrapure, 99.999%**	5 Gms	5000
		25 Gms	12000
59479 [13494-80-9]	**Tellurium Metal powder, 99.9%, -200 mesh**	10 Gms	2200
		100 Gms	17000

Code	Name	Packing	Price
67813 [13798-24-8]	**Terbium (III) Chloride Hexahydrate extrapure, 99.9%**	5 Gms 25 Gms	10000 40000
52633 [13451-19-9]	**Terbium (III) Nitrate Hexahydrate ultrapure, 99.99%**	2 Gms 10 Gms	9000 34000
52313 [57584-27-7]	**Terbium (III) Nitrate Pentahydrate extrapure, 99.9%**	5 Gms 25 Gms	10000 40000
99186 [148980-31-8]	**Terbium (III) Trifluoromethanesulfonate extrapure, 98%** (Terbium(III) Triflate, Trifluoromethanesulfonic Acid Terbium(III) Salt)	1 Gms 5 Gms	8000 20000
62369 [12037-01-3]	**Terbium (III,IV) Oxide ultrapure, 99.99%**	1 Gms 5 Gms	10000 40000
96104 [100-20-9]	**Terephthaloyl Chloride pure, 99%**	100 Gms 500 Gms	1300 5400
14676 [7365-44-8]	**TES Buffer extrapure, 99%** (N-Tris(Hydroxymethyl) Methyl-2-Aminoethane Sulphonic Acid)	5 Gms 25 Gms 100 Gms	800 2800 10500
51318 [70331-82-7]	**TES Buffer Sodium Salt extrapure, 99%**	25 Gms 100 Gms	5900 17930
19914 [58-22-0]	● **Testosterone (free base) extrapure, 99%** (17ß-Hydroxy-4-Androsten-3-one)	1 Gms 5 Gms 25 Gms	1800 3500 12000
83651 [315-37-7]	● **Testosterone Enanthate extrapure, 99%** (17ß-Hydroxy-4-Androsten-3-one 17-Enanthate)	1 Gms 5 Gms 25 Gms	2000 5000 13000
71344 [57-85-2]	● **Testosterone Propionate extrapure, 98%** (17ß-Hydroxy-4-Androsten-3-one 17-Propionate)	5 Gms 25 Gms	1670 7420
93831 [13035-61-5]	✔ **1,2,3,5-Tetra-O-Acetyl-b-D-Ribofuranose extrapure, 98%** (beta-D-Ribofuranose-1,2,3,5-Tetraacetate)	5 Gms 25 Gms	2620 9490
76333 [79-34-5]	ⓘ **1,1,2,2-Tetrachloroethane Reagent Grade, 99%**	500 ml 1000 ml	1750 3430
82175 [64-75-5]	■ **Tetracycline Hydrochloride (TE) for tissue culture, 95%**	5 Gms 25 Gms 100 Gms	1050 4000 11100
38614 [64-75-5]	■ **Tetracycline Hydrochloride (TE), 95%**	5 Gms 25 Gms 100 Gms	440 2170 7950
58128 [1191-50-0]	**Tetradecyl Sulphate Sodium Salt (Tergitol-4), 95%** (Sodium Myristyl Sulfate, Sodium Tetradecyl Sulfate, STDS)	100 Mg 500 Mg	1500 5800
64090 [68-05-3]	**Tetraethylammonium Iodide (TEAI) extrapure, 99%**	25 Gms 100 Gms 500 Gms	1700 5100 23600

Code	Name	Packing	Price
23231 [74163-81-8]	✔ **(S)-(-)-1,2,3,4-Tetrahydroisoquinoline-3-Carboxylic Acid pure, 98%** ((S) -1,2,3,4-Tetrahydro-3-Isoquinolinecarboxylic Acid)	1 Gms 5 Gms 25 Gms	900 2000 6250
55633 [14221-01-3]	🌡 **Tetrakis(Triphenylphosphine) Palladium (0) extrapure, 97%**	1 Gms 5 Gms 25 Gms	4000 17000 75000
76265 [529-34-0]	**1-Tetralone pure, 98%**	100 ml 500 ml	1300 5250
98500 [25103-58-6]	**Tert-Dodecylmercaptan (mix. of isomers) extrapure, 98.5%** (tert-Dodecanethiol, tert-Laurylmercaptan)	100 ml 500 ml	1200 2400
19028 [75-57-0]	**Tetramethylammonium Chloride (TMAC) pure, 98%**	100 Gms 500 Gms	520 2200
82774 [75-58-1]	**Tetramethylammonium Iodide (TMAI) extrapure, 99%**	100 Gms	4100
84821 [54827-17-7]	🌡 **3,3,5,5-Tetramethyl Benzidine (TMB) pure, 98%**	1 Gms 5 Gms 25 Gms	1950 5000 18000
69860 [54827-17-7]	🌡 **3,3,5,5-Tetramethyl Benzidine (TMB) extrapure AR, 98%**	1 Gms 5 Gms 25 Gms	2400 6400 24000
68738 [54827-17-7]	🌡 **3,3,5,5-Tetramethylbenzidine (TMB) extrapure AR, ExiPlus, Multi-Compendial, 98%** meets compendial specs of BP, Ph.Eur	1 Gms 5 Gms	3200 8000
55597 [64285-73-0]	🌡 **3,3,5,5-Tetramethylbenzidine Dihydrochloride anhydrous (TMB.diHCl anhy) extrapure AR, 99%**	100 Mg 1 Gms 5 Gms	2730 5820 18980
36746 [207738-08-7]	🌡 **3,3,5,5-Tetramethyl Benzidine Dihydrochloride Hydrate (TMB.2HCl.xH2O) extrapure, 98%** (TMB DiHCl)	500 Mg 1 Gms 5 Gms	3000 5760 19610
64407	🌡 **TMB Solution A (1% TMB in DMSO)**	100 ml	1560
38783	🌡 **TMB Solution B (30% H2O2)**	100 ml	910
99986 [366-29-0]	✔ **N,N,N,N-Tetramethylbenzidine extrapure, 95%**	100 Mg 1 Gms 5 Gms	1520 3000 8800
13618 [3010-96-6]	🌡 **2,2,4,4-Tetramethyl-1,3-Cyclobutanediol (CBDO) (mixture of isomers) extrapure, 99%**	5 Gms 25 Gms	5250 15500
88001 [3277-26-7]	🌡 **1,1,3,3-Tetramethyldisiloxane (TMDS) extrapure, 98%**	25 ml 100 ml	1230 3740
53180 [80-70-6]	ⓘ **1,1,3,3-Tetramethylguanidine extrapure AR, 99%** N,N,N',N'-Tetramethylguanidine	100 ml 500 ml	2100 8190
52300 [637-01-4]	🌡 **N,N,N,N-Tetramethyl-p-phenylenediamine Dihydrochloride (Wursters Reagent, TMPPD) extrapure AR, 99%**	5 Gms 25 Gms	1240 5600

Code	Name	Packing	Price
91902 [2564-83-2]	2,2,6,6-Tetramethylpiperidine-1-Oxyl (TEMPO) pure, 98%	5 Gms	500
		25 Gms	2200
		100 Gms	6000
33379 [75-76-3]	■ Tetramethyl Silane (TMS) for NMR spectroscopy, 99.8%	25 ml	4910
		100 ml	14020
77204 [1941-30-6]	Tetrapropylammonium Bromide extrapure, 99%	100 Gms	700
		500 Gms	2500
84990 [4499-86-9]	Tetrapropylammonium Hydroxide (TPAH) 10% Aq. Solution	100 ml	1100
		500 ml	3200
98078 [4499-86-9]	Tetrapropylammonium Hydroxide (TPAH) 20% Aq. Solution	100 ml	1500
		500 ml	6000
58146 [4499-86-9]	Tetrapropylammonium Hydroxide (TPAH) 40% Aq. Solution	100 ml	3000
		500 ml	13000
75311 [288-94-8]	1-H-Tetrazole sublimed extrapure, 99%	1 Gms	300
		5 Gms	1350
		25 Gms	5090
92739 [1719-71-7]	Tetrazolium Violet extrapure, 95%	1 Gms	2950
62183 [1719-71-7]	Tetrazolium Violet (high purity) extrapure AR, 98%	100 Mg	1300
		1 Gms	4110
		5 Gms	16010
44826 [3081-61-6]	L-Theanine extrapure, 99%	25 Gms	850
		100 Gms	1800
		500 Gms	7000
83155 [326-91-0]	3(2-Thenoyl)-1,1,1-Trifluoroacetone (TTFA) extrapure AR, 99% (2-Thenoyltrifluoroacetone, TTA (extractant), TTB)	10 Gms	1260
		100 Gms	8000
71081 [67-03-8]	Thiamine Hydrochloride pure, 98.5% (Vitamin B1)	25 Gms	550
		100 Gms	2200
		500 Gms	8250
82500 [67-03-8]	Thiamine Hydrochloride for tissue culture, 98.5% (Vitamin B1)	25 Gms	800
		100 Gms	2500
39427 [67-03-8]	Thiamine Hydrochloride ExiPlus, Multi-Compendial, 98.5% meets compendial specs of USP, BP, Ph.Eur	25 Gms	1250
		100 Gms	3850
		500 Gms	9500
79903 [15318-45-3]	Thiamphenicol (THP), 98%	1 Gms	1200
		5 Gms	3000
33611 [298-93-1]	Thiazolyl Blue Tetrazolium Bromide (MTT) extrapure AR, 98%	100 Mg	480
		500 Mg	1600
		1 Gms	3000
		5 Gms	11000
58945 [298-93-1]	Thiazolyl Blue Tetrazolium Bromide (MTT) for tissue culture, 98%	500 Mg	2000
		1 Gms	3600
58910 [57360-69-7]	Thiazolyl Blue Formazan (MTT Formazan) extrapure AR, 98% (1-(4,5-Dimethylthiazol-2-yl) -3,5-Diphenylformazan)	250 Mg	4600
		1 Gms	11300
11652 [51707-55-2]	Thidiazuron (TDZ) extrapure, 97%	1 Gms	3500
		5 Gms	6500
		25 Gms	21000
88098 [51707-55-2]	Thidiazuron (TDZ) for tissue culture, 97%	250 Mg	2500
		1 Gms	5000
		5 Gms	10500
10083 [507-09-5]	Thioacetic Acid pure, 97%	100 ml	900
		500 ml	3400
94066 [602-41-5]	Thiocolchicoside extrapure, 98% (2-Demethoxy-2-Glucosidoxythiocolchicine)	50 Mg	1500
		250 Mg	2500
		1 Gms	6000
46913 [1077-28-7]	● DL-6,8-Thioctic Acid (a-Lipoic Acid) extrapure, 98%	5 Gms	2540
		25 Gms	6000
		250 Gms	9000
45487 [70-49-5]	Thiomalic Acid (Mercaptosuccinic Acid) pure, 99%	25 Gms	840
		100 Gms	1370
		500 Gms	6090
85090 [54-64-8]	Thiomersal (Thimerosal) ExiPlus, Multi-Compendial, 97% meets compendial specs of USP, BP, Ph.Eur	5 Gms	4300
		10 Gms	8250
		25 Gms	20000
		100 Gms	75000
27662 [54-64-8]	Thiomersal (Thimerosal) PrimeDx, 98%	10 Gms	9000
		25 Gms	22000
		100 Gms	82000
40413 [6165-68-0]	Thiophene-2-Boronic Acid extrapure, 98% (2-Thienylboronic Acid, Thien-5-ylboronic Acid)	1 Gms	810
		5 Gms	2220
		25 Gms	9980
57589 [527-72-0]	Thiophene-2-Carboxylic Acid pure, 98%	5 Gms	1640
		25 Gms	5430
		100 Gms	16160
39416 [4090-29-3]	■ Thionicotinamide Adenine Dinucleotide (Thio-NAD) extrapure, 90%	25 Mg	3470
		100 Mg	12980
35020 [1921-48-8]	■ Thionicotinamide Adenine Dinucleotide Disodium Salt Reduced (Thio-NADH) extrapure, 93%	25 Mg	2400
		100 Mg	6500
29808 [1918-77-0]	Thiophene-2-Acetic Acid pure, 98%	5 Gms	930
		25 Gms	3000
		100 Gms	10130
14055 [98-03-3]	2-Thiophenecarboxaldehyde pure, 98% (Thiophene-2-Carboxaldehyde, Thenaldehyde, Thiophene-2-Aldehyde)	25 ml	900
		100 ml	3400
		250 ml	8000
28033 [34592-47-7]	L-Thioproline extrapure, 98% (L-4-Thiazolidinecarboxylic Acid)	10 Gms	600
		25 Gms	1200
83203 [147-93-3]	Thiosalicylic Acid pure, 98% (2-Mercaptobenzoic Acid)	25 Gms	840
		100 Gms	2210
96768 [141-90-2]	2-Thiouracil extrapure, 98% (2-Mercaptopyrimidine-4-ol)	25 Gms	1050
		100 Gms	3890

Code	Name	Packing	Price
41862 [3688-92-4]	**Thorin extrapure AR**	1 Gms 5 Gms	3200 15000
50444 [632-20-2]	**D-Threonine (allo-free) extrapure, 99%**	5 Gms 25 Gms 100 Gms	1800 6580 23910
75164 [80-68-2]	**DL-Threonine extrapure CHR, 99%**	25 Gms 100 Gms	3400 12000
88355 [72-19-5]	**L-Threonine extrapure CHR, 99%**	100 Gms 500 Gms 1 Kg 5 Kg	650 2600 4500 20000
62667 [72-19-5]	**L-Threonine for tissue culture, 99%**	25 Gms 100 Gms 500 Gms	500 1800 7500
24815 [44520-55-0]	**D-Threoninol extrapure, 98%**	500 Mg 1 Gms	6130 10550
20724 [50-89-5]	**Thymidine extrapure, 99%** (2-Deoxythymidine, B-Thymidine, Thymine Deoxyriboside, 5-Methyldeoxyuridine)	5 Gms 25 Gms	900 3200
45764 [33430-62-5]	**Thymidine-5-Monophosphate Disodium Salt (TMP-Na2) extrapure, 98%**	100 Mg 250 Mg	700 2000
57046 [65-71-4]	**Thymine extrapure, 98%**	25 Gms 100 Gms	700 2000
47110 [89-83-8]	**Thymol Crystal pure, 99%**	100 Gms 500 Gms 5 Kg	550 1600 14000
77675 [1913-93-5]	**Thymolphthelexone**	1 Gms 5 Gms	3050 11620
86637 [7512-38-1]	**Thymol Violet**	5 Gms	470
90196 [80-59-1]	**Tiglic Acid extrapure, 98%**	25 Gms 100 Gms	1310 4620
77082 [411207-31-3]	**Tiotropium Bromide Monohydrate (TioBr) extrapure, 98%**	50 Mg 1 Gms	3250 35000
87347 [7550-45-0]	**Titanium Tetrachloride 1M in Dichloromethane** (Titanium (IV) Chloride in MDC)	100 ml 500 ml	1080 4990
83692 [7550-45-0]	**Titanium Tetrachloride 1M in Toluene** (Titanium (IV) Chloride in Toluene)	100 ml 500 ml	1080 4990
84161 [125700-73-4]	**TNTU extrapure, 99%**	25 Gms 100 Gms	4000 14000
56183 [32986-56-4]	**Tobramycin (TBM)**	50 Mg 500 Mg	1200 4000
59844 [7695-91-2]	**DL-a-Tocopherol Acetate (Vitamin E Acetate) for tissue culture, 98%**	5 Gms 25 Gms 100 Gms	360 1000 3500
52885 [887144-97-0]	**Togni`s Reagent (3,3-Dimethyl-1-(trifluoromethyl)-1,2-benziodoxole) extrapure, 97%** (1-Trifluoromethyl-3,3-Dimethyl-1,2-Benziodoxole)	250 Mg 1 Gms 5 Gms	6100 15700 26900
67593 [119-93-7]	**o-Tolidine extrapure AR, 98%** (3,3- Dimethyl Benzidine)	25 Gms 100 Gms 500 Gms	690 2100 3780
65954 [119-93-7]	**o-Tolidine extrapure AR, ExiPlus, Multi-Compendial, 98%** (3,3- Dimethyl Benzidine) meets compendial specs of USP, BP, Ph.Eur	25 Gms 100 Gms 500 Gms	790 2520 4200
82897 [118-90-1]	**o-Toluic Acid pure, 99.5%**	250 Gms 500 Gms	500 890
70097 [99-04-7]	**m-Toluic Acid pure, 98%**	250 Gms 500 Gms	550 900
61421 [90-94-5]	**p-Toluic Acid pure, 98%**	250 Gms 500 Gms	480 850
12358 [106-49-0]	**p-Toluidine pure, 98%**	500 Gms	1300
86003 [106-49-0]	**p-Toluidine specially purified extrapure, 99%**	500 Gms	1600
22134 [92-31-9]	**Toluidine Blue O, 90%** (Basic Blue 17, Blutene Chloride, Methylene Blue T50, Tolonium Chloride)	10 Gms 25 Gms 100 Gms	340 770 2660
37699 [92-31-9]	**Toluidine Blue O for tissue culture, 90%** (Basic Blue 17, Blutene Chloride, Methylene Blue T50, Tolonium Chloride)	10 Gms 25 Gms	910 2050
20611 [1576-35-8]	**p-Toluenesulfonyl Hydrazide pure, 98%** (p-Toluenesulfonic Acid Hydrazide, p-Toluenesulphonyl Hydrazide, Tosylhydrazide)	25 Gms 100 Gms 500 Gms	900 3400 11000
68913 [75737-38-1]	**meta-Topolin (mT) for molecular biology, 98%** (6-(3-Hydroxybenzylamino)purine)	25 Mg 100 Mg 1 Gms	1500 3000 12000
56521 [82692-93-1]	**TOOS extrapure, 98%** (Sodium 3-[Ethyl(m-Tolyl)amino] -2-Hydroxy-1-Propanesulfonate, (3-(N-Ethyl-3-methylanilino)-2-Hydroxy	1 Gms 5 Gms 25 Gms	4500 12500 55000
95595 [136849-72-4]	**TOTU extrapure, 99%**	5 Gms 25 Gms 100 Gms	2100 6000 22000
25793 [3682-35-7]	**TPTZ (2,4,6-Tri-(2-Pyridyl) -5-Triazine) extrapure AR , 99%**	1 Gms 5 Gms 25 Gms	1200 4500 16000
49011 [593-50-0]	**Triacontanol (TRIA) extrapure, 90%**	500 Mg 1 Gms	8000 14000
68084 [3588-17-8]	**trans,trans-Muconic acid (ttMA) extrapure, 98%** (trans,trans-2,4-Hexadienedioic Acid, trans,trans-1,3-Butadiene-1,4-Dicarboxylic Acid)	1 Gms 5 Gms	2380 11150

Code	Name	Packing	Price
42618 [288-88-0]	**1,2,4-Triazole pure, 98%**	100 Gms	290
		500 Gms	1300
93572 [41253-21-8]	**1,2,4-Triazole Sodium extrapure, 98%** (Sodium-1,2,4-Triazolide)	25 Gms	760
		100 Gms	2800
		500 Gms	5000
78411 [626-39-1]	**1,3,5-Tribromobenzene pure, 98%**	25 Gms	3700
76074 [90131-68-3]	**Tribulus Terrestris Extract, 45% FS**	5 Gms	2600
		25 Gms	5000
86118 [120-82-1]	**1,2,4-Trichlorobenzene (1,2,4-TCB) pure, 99%**	500 ml	500
		1000 ml	1050
		2500 ml	2000
99486 [108-70-3]	**1,3,5-Trichlorobenzene (1,3,5-TCB) pure, 99%**	100 Gms	2800
		500 Gms	9500
94076 [115-20-8]	**2,2,2-Trichloroethanol pure, 98%** (Trichloroethyl alcohol, ß,ß,ß-Trichloroethyl alcohol, 2,2,2-Trichloroethan-1-ol)	25 ml	2200
		100 ml	7500
		250 ml	12000
27740 [5704-04-1]	**TRICINE Buffer Reagent Grade, 99%** (N-Tris-(hydroxymethyl) Methylglycine)	25 Gms	1100
		100 Gms	4400
		500 Gms	20000
67490 [2622-14-2]	**Tricyclohexylphosphine (Cy3P) extrapure, 94%** (P(Cy)3, PCy3)	5 Gms	1500
		25 Gms	6500
		100 Gms	25000
63437 [58656-04-5]	**Tricyclohexylphosphine Tetrafluoroborate extrapure, 98%**	1 Gms	2170
		5 Gms	6490
80184 [637-39-8]	**Triethanolamine Hydrochloride extrapure AR, 99%**	100 Gms	500
		500 Gms	1730
80058 [75-89-8]	**2,2,2-Trifluoroethanol pure, 99%**	25 ml	490
		100 ml	1210
		500 ml	5060
14428 [10009-20-8]	**N-6-Trifluoroacetyl-L-Lysine extrapure, 98%** (H-Lys(Tfa)-OH)	5 Gms	1890
		25 Gms	3800
12631 [1423-26-3]	**3-Trifluoromethyl Phenylboronic Acid extrapure, 97%**	1 Gms	1550
		5 Gms	2880
		25 Gms	14320
79009 [70-00-8]	**5-Trifluorothymidine extrapure, 99%**	100 Mg	4830
		1 Gms	36800
37758 [88-82-4]	**2,3,5-Triiodobenzoic Acid (TIBA) extrapure AR, 98%**	5 Gms	690
		10 Gms	1380
		25 Gms	2880
35241 [6893-02-3]	**3,3,5-Triiodo-L-Thyronine extrapure, 95%** (3,3,5-Triiodothyronin)	100 Mg	4500
		1 Gms	26000
14219 [538-24-9]	**Trilaurin pure, 98%(GC)** (Glycerol Trilaurate, Glyceryl Trilaurate, Glyceryl Tridodecanoate)	5 Gms	1500
		25 Gms	3700
		100 Gms	8600
84217 [538-24-9]	**Trilaurin Reference Standard, 99% (GC)** (Glycerol Trilaurate, Glyceryl Trilaurate, Glyceryl Tridodecanoate)	100 Mg	2200
		1 Gms	6000
		5 Gms	13000
57355 [554-95-0]	**Trimesic Acid extrapure, 98%** (Benzene-1,3,5-Tricarboxylic acid, 1,3,5-Benzenetricarboxylic Acid)	5 Gms	740
		25 Gms	2630
		100 Gms	4200
20139 [75-77-4]	**Trimethyl Chlorosilane (TMCS) extrapure, 98%**	250 ml	700
		1000 ml	2400
		2500 ml	5400
		25 Ltr	52000
66889 [142-25-6]	**N,N,N-Trimethylethylenediamine pure, 97%**	5 Gms	3200
		25 Gms	11000
18879 [738-70-5]	**Trimethoprim (TMP), 98.5%**	5 Gms	750
		25 Gms	2400
66069 [86-81-7]	**3,4,5-Trimethoxy Benzaldehyde pure, 98%**	25 Gms	1040
		100 Gms	3450
41569 [621-23-8]	**1,3,5-Trimethoxybenzene pure, 99%**	25 Gms	3020
40017 [118-41-2]	**3,4,5-Trimethoxybenzoic Acid extrapure, 99%** (Gallic Acid Trimethyl Ether, Trimethylgallic Acid)	25 Gms	950
		100 Gms	3150
49601 [35141-30-1]	**N-(3-Trimethoxysilylpropyl) Diethylenetriamine (TMSPD) pure, 95%** (2-[2-(3-Trimethoxysilylpropylamino)ethylamino]ethylamine, 3-[2-(2-Aminoethylamino)ethylamino]propyl	100 ml	5000
		500 ml	20000
90265 [75-77-4]	**Trimethyl Chlorosilane (TMCS) extrapure, 99%**	250 ml	800
		1000 ml	2900
		2500 ml	6500
67883 [75-77-4]	**Trimethyl Chlorosilane (TMCS) 1M Solution in THF**	100 ml	3500
41020 [10108-86-8]	**Trimethyloctylammonium Chloride (TMOACl) extrapure, 98%** (Octyltrimethylammonium Chloride)	25 Gms	2000
		100 Gms	6500
68989 [149-73-5]	**Trimethyl Orthoformate (TMOF) pure, 98%**	500 ml	1350
		2500 ml	6000
		25 Ltr	53000
49901 [996-50-9]	**N-(Trimethylsilyl)diethylamine (TMSDEA) extrapure, 97%** ((Diethylamino)trimethylsilane, N,N-Diethyltrimethylsilylamine)	25 ml	5110
		100 ml	14580
45238 [18156-74-6]	**N-Trimethylsilylimidazole (TMSIM) extrapure, 98%** ((1-(Trimethylsilyl)imidazole)N-Trimethylsilylimidazole, TSIM)	5 ml	700
		25 ml	2100
		100 ml	6800
97039 [5683-30-7]	**3-Trimethylsilylpropionic Acid (TMSPA) extrapure, 95%**	1 ml	5470
		5 ml	17620
39671 [1774-47-6]	**Trimethylsulfoxonium Iodide (TMSI) extrapure, 99%**	25 Gms	1000
		100 Gms	2300
		500 Gms	9900
18779 [555-45-3]	**Trimyristin pure, 95%(GC)** (Glyceryl Trimyristate, Glycerol Trimyristate)	5 Gms	2000
		25 Gms	4600

Code	Name	Packing	Price
83676 [555-45-3]	**Trimyristin Reference Standard, 99%(GC)** (Glyceryl Trimyristate, Glycerol Trimyristate)	1 Gms 5 Gms 25 Gms	4000 12000 21000
84050 [1116-76-3]	**Trioctylamine (Tri-n-Octylamine) pure, 95%**	100 ml 500 ml	700 2000
37958 [122-32-7]	**Triolein Reference Standard, 99% (GC)**	500 Mg 1 Gms 5 Gms	3600 6000 28000
59708 [555-44-2]	**Tripalmitin pure, 85% (GC)** (Glycerol Tripalmitate, Glyceryl Tripalmitate)	25 Gms 100 Gms	5000 12000
75700 [555-44-2]	**Tripalmitin Reference Standard, 99% (GC)** (Glycerol Tripalmitate, Glyceryl Tripalmitate)	100 Mg 500 Mg 1 Gms	1200 3500 6000
38015 [5893-05-0]	**N-Triphenyl Methyl Glycine extrapure, 98%**	5 Gms 25 Gms	2000 9000
24651 [603-35-0]	**Triphenylphosphine pure, 98%**	100 Gms 500 Gms 1 Kg	650 2600 4500
43068 [14243-64-2]	**Triphenylphosphinegold(I) Chloride ((Ph3P)AuCl), 98%** (Chloro(triphenylphosphine)Gold(I))	100 Mg 500 Mg 1 Gms	6000 24000 42000
95474 [6399-81-1]	**Triphenylphosphine Hydrobromide pure, 98%**	5 Gms 25 Gms	800 1500
32098 [791-28-6]	**Triphenylphosphine Oxide (TPPO) pure, 98.5%** (Triphenyl Phosphorus Oxide)	25 Gms 100 Gms 500 Gms	350 1200 3000
70874 [1034-39-5]	**Triphenylphosphonium Dibromide pure, 98%**	5 Gms 25 Gms 100 Gms	1600 4000 10000
53730 [531-52-2]	**1,3,5-Triphenyltetrazolium Formazan (TTC Formazan) extrapure, 95%**	1 Gms 5 Gms 25 Gms	1600 4700 8500
65599 [298-96-4]	**2,3,5-Triphenyltetrazolium Chloride (TTC) extrapure, 99%**	10 Gms 25 Gms 50 Gms 100 Gms	2000 4000 7000 13000
98635	**10X Tris Acetate-SDS Buffer pH-8.3 for molecular biology** (TA-SDS Running Buffer)	500 ml 1000 ml	800 1500
47804	**20X Tris Acetate-SDS Buffer pH-8.3 for molecular biology** (TA-SDS Running Buffer)	500 ml 1000 ml	1500 2850
65421 [60804-74-2]	**Tris(2,2-bipyridine)ruthenium(II) Hexafluorophosphate ((Ru(BPY)3(PF6)2), 97%** (Ruthenium-tris(2,2'-bipyridyl) dihexafluorophosphate)	250 Mg 1 Gms	4200 15000
26055 [50525-27-4]	**Tris(2,2-bipyridyl)dichlororutheniu m (II) Hexahydrate (Ru(BPY)3), 98%** (Ruthenium-tris(2,2'-bipyridyl) dichloride, Tris(2,2'-bipyridyl)ruthenium(II) chloride hexahydrate)	250 Mg 1 Gms	3800 7000
55664 [51805-45-9]	**Tris(2-carboxyethyl) Phosphine Hydrochloride (TCEP) extrapure AR, 98%**	1 Gms 5 Gms 25 Gms	1200 2800 8000
79175 [51364-51-3]	**Tris(dibenzylideneacetone)dipallad ium(0) (Pd2(dba)3) extrapure, 22-25% Pd**	500 Mg 1 Gms 5 Gms 25 Gms	5000 8000 32000 150000
96455 [52522-40-4]	**Tris(Dibenzylideneacetone)Dipalla dium-Chloroform Adduct extrapure, 20% Pd** (Dipalladium-tris(Dibenzylideneacetone)Chloroform Complex, Pd2(dba)3.CHCl3)	250 Mg 1 Gms 5 Gms	2500 7000 25000
87326	**10X Tris-Glycine pH-8.3 Tank Buffer for molecular biology**	200 ml 1000 ml	840 3150
57806	**10X Tris-Glycine-SDS Buffer for molecular biology** (TG-SDS Running Buffer)	200 ml 1000 ml	480 1070
83471	**10X Tris Buffered Saline (TBS) for molecular biology** (10X TBS Buffer)	500 ml	1960
63825	**10X Tris Buffered Saline Tween-20 (TBST) for molecular biology**	500 ml	520
79091	**1M Tris Hydrochloride Buffer (1M Tris HCl) pH-7.2 for molecular biology** (Tris(hydroxymethyl)aminomethane, Hydrochloric Acid)	100 ml 500 ml	1110 4680
55904	**1M Tris Hydrochloride Buffer (1M Tris HCl) pH-8.0 for molecular biology** (Tris(hydroxymethyl)aminomethane, Hydrochloric Acid)	100 ml 500 ml	1640 4910
80743 [1185-53-1]	**2.6M Tris Hydrochloride Buffer (2.6M Tris HCl) pH-8 for molecular biology** (Tris(hydroxymethyl)aminomethane, Hydrochloric Acid)	100 ml 500 ml	1900 5400
37852	**10X Tris-Tricine-SDS Buffer for molecular biology** (TT-SDS Running Buffer)	500 ml	1210
62850 [6850-28-8]	**Tris Acetate extrapure AR, 99%** (Tris(Hydroxymethyl)aminomethane Acetate)	100 Gms 500 Gms	1400 6530
34969 [6850-28-8]	**Tris Acetate for molecular biology, 99%** (Tris(Hydroxymethyl) Aminomethane Acetate)	100 Gms 500 Gms	1900 8900
23017 [68123-29-5]	**Tris Carbonate extrapure, 99%** (Tris(Hydroxymethyl)aminomethane Carbonate)	100 Gms	2070
32119 [72200-76-1]	**Tris Maleate extrapure, 99.5%** (Tris(Hydroxymethyl)aminomethane Maleate)	100 Gms 500 Gms	1550 7000

Code	Name	Packing	Price
25913 [41521-38-4]	**Tris Nitrate extrapure, 99%** Tris(Hydroxymethyl)aminomethane Nitrate	100 Gms	2070
55088 [108321-13-7]	**Tris Oxalate extrapure, 99%** (Tris(Hydroxymethyl)aminomethane Oxalate)	100 Gms	1450
68353 [108321-11-5]	**Tris Phosphate Dibasic extrapure, 98%** (Tris(Hydroxymethyl)aminomethane Phosphate Dibasic)	100 Gms	1640
43235 [6992-39-8]	**Tris(Hydroxymethyl)-Aminomethane Phosphate Monobasic extrapure (Tris Phosphate Monobasic), 98%**	100 Gms	1640
63291 [85169-32-0]	**Tris Succinate extrapure, 99%** (Tris(Hydroxymethyl)aminomethane Succinate)	100 Gms	1640
47952 [6992-38-7]	**Tris Sulphate extrapure, 98%** (Tris(Hydroxymethyl)aminomethane Sulphate)	100 Gms	1090
59666 [14694-95-2]	**Tris(Triphenylphosphine) Rhodium (I) Chloride extrapure, 99.9%** (Rhodium(I) Tris(Triphenylphosphine) Chloride, Wilkinson's Catalyst)	250 Mg 1 Gms 5 Gms	9000 33000 150000
62837 [15529-49-4]	**Tris(Triphenylphosphine) Ruthenium (II) Dichloride pure, 10-11% Ru** (Tris(Triphenylphosphine) Dichlororuthenium)	1 Gms 5 Gms	3500 14000
50615 [555-43-1]	**Tristearin pure, 90%(GC)** (Glycerol Tristearate, Glyceryl Tristearate)	25 Gms 100 Gms	2380 8660
11085 [555-43-1]	**Tristearin Reference Standard, 99% (GC)** (Glycerol Tristearate, Glyceryl Tristearate)	500 Mg 1 Gms 5 Gms	2500 4000 16000
48445 [2799-07-7]	**S-Trityl-L-Cysteine extrapure, 98%**	1 Gms	4510
58026 [72-57-1]	**Trypan Blue** (Direct blue 14) (C.I. No. 23850)	25 Gms 100 Gms	1200 4500
31591 [72-57-1]	**Trypan Blue for tissue culture** (Direct blue 14) (C.I. No. 23850)	25 Gms	1670
12788 [72-57-1]	**Trypan Blue 0.4% Solution**	50 ml 125 ml	500 1100
60484 [9002-07-7]	■ **Trypsin (3x cryst.) ex Bovine Pancreas, 7500 BAEE U/mg (2500NF U/mg) powder**	250 Mg 1 Gms 5 Gms	1180 2380 10600
42657 [9035-81-8]	■ **Trypsin Inhibitor ex. Soyabean, 10000 BAEE U/mg**	25 Mg 100 Mg 250 Mg 500 Mg 1 Gms	3500 9000 18000 28000 51000
16406 [9035-81-8]	■ **Trypsin Inhibitor (Ovomucoid) ex. Chicken Egg White** (Trypsin Inhibitor ex. Gallus Gallus Domesticus)	100 Mg 1 Gms	5000 14000
12094 [343-94-2]	**Tryptamine Hydrochloride extrapure, 98%**	1 Gms 5 Gms 25 Gms	610 2000 9000

Code	Name	Packing	Price
67654 [7524-52-9]	**L-Tryptophan Methyl Ester Hydrochloride extrapure, 99%**	5 Gms 25 Gms	2140 9250
93900 [105832-38-0]	**TSTU extrapure, 98%** (N,N,N`,N`-Tetramethyl-O-(N-succinimidyl)uronium tetrafluoroborate, O-(N-Succinimidyl) -N,N,N',N'-tet	5 Gms 25 Gms 100 Gms	1500 3200 11000
51335 [7440-33-7]	**Tungsten Metal Powder, 99.9%, -325 mesh**	100 Gms 500 Gms	10000 45000
17472 [7783-03-1]	**Tungstic Acid pure, 98%**	100 Gms 250 Gms	7500 17000
55401 [7783-03-1]	**Tungstic Acid extrapure AR, 99%**	100 Gms 500 Gms	8500 38000
98328 [8002-64-2]	**Turpentine Oil extrapure**	500 ml 2500 ml 25 Ltr	350 1600 15000
33541 [60-19-5]	**Tyramine Hydrochloride extrapure, 98%**	25 Gms 100 Gms 500 Gms	1600 4000 9500
94870 [9002-10-2]	■ **Tyrosinase (Polyphenol Oxidase) ex. Mushroom, 600U/mg powder** (Catechol Oxidase, Monophenol Monooxygenase, Dihydroxyphenylalanine:Oxygen Oxidoreductase)	25 K.Units 100 K.Units	12000 42000
54194 [25301-02-4]	**Tyloxapol** Ethoxylated p-Tert-Octylphenol Formaldehyde Polymer	10 Gms 25 Gms 50 Gms	10200 24500 46000
46121 [93-35-6]	**Umbelliferone (7-Hydroxycoumarin) extrapure, 98%**	10 Gms 25 Gms	2100 4500
46487 [112-38-9]	**Undecylenic Acid extrapure,C11-99% (GC)**	100 Gms 500 Gms	650 2810
23622 [66-22-8]	**Uracil extrapure, 99%**	25 Gms 100 Gms 1 Kg 5 Kg	430 1550 12000 50000
81405 [6159-44-0]	✔ **Uranyl Acetate Dihydrate ACS, 98%**	1 Gms 5 Gms 25 Gms	3600 13500 72000
70780 [20596-93-4]	**Uranyl Magnesium Acetate extrapure, 98%**	1 Gms 5 Gms	4200 19000
74495 [13520-83-7]	✔ **Uranyl Nitrate Hexahydrate extrapure, 98%**	1 Gms 5 Gms 25 Gms	4000 17000 82000
64807 [9002-13-5]	■ **Urease ex. Jack Beans, 200U/mg powder**	3000 Units 10000 Units 25 K.Units 100 K.Units	1400 3140 7320 12630
64632	■ **r-Urease (Thermostable) ex. E. Coli, 150U/mg powder**	10000 Units	14720

Code	Name	Packing	Price
90739 [69-93-2]	**Uric Acid extrapure AR, 99%**	25 Gms	4200
		100 Gms	15680
		500 Gms	67200
95410 [9002-12-4]	**r-Uricase (Thermostable) ex. Yeast, 3U/mg powder**	500 Units	6490
62205 [58-96-8]	**Uridine extrapure, 99%**	5 Gms	800
		25 Gms	2400
		100 Gms	7000
20964 [28053-08-9, 117756-22-6]	**Uridine 5`-Diphosphoglucose Disodium Salt Hydrate extrapure (UDPG.2Na Hydrate), 98%**	25 Mg	1400
		100 Mg	3200
		1 Gms	12000
48301 [3387-36-8]	**Uridine 5-Monophosphate Disodium Salt Hydrate (5-UMP-Na2) extrapure, 98.5%**	5 Gms	900
		25 Gms	2500
56613 [19817-92-6]	**Uridine-5-Triphosphate Trisodium Salt Dihydrate (UTP-Na3) extrapure, 90%**	100 Mg	1500
		500 Mg	4000
		1 Gms	7000
50178 [128-13-2]	**Ursodeoxycholic Acid (UDCS) extrapure, 99%**	1 Gms	950
		5 Gms	1500
		25 Gms	6000
83198 [1143-70-0]	**Urolithin A extrapure, 97%** (3,8-Dihydroxy-6H-benzo[c]chromen-6 -one,	1 Gms	1800
		5 Gms	3000
		25 Gms	8500
24436 [7562-61-0]	**(+)- Usnic Acid, 98%**	1 Gms	1500
		5 Gms	3000
		25 Gms	9000
44762 [110-62-3]	**Valeraldehyde pure (Pentanal), 98%** (Pentanal)	100 ml	470
		500 ml	1200
31477 [109-52-4]	**n-Valeric Acid pure, C5-99% (GC)**	500 ml	1060
		1000 ml	1990
22128 [110-59-8]	**Valeronitrile (Butyl Cyanide) pure , 99%** (Butyl Cyanide)	100 ml	400
		500 ml	1350
80634 [1009-14-9]	**Valerophenone pure, 98%**	100 ml	500
		500 ml	1800
		2500 ml	8000
95091 [638-29-9]	**Valeroyl Chloride pure, 98%**	100 ml	700
		500 ml	2800
21562 [640-68-6]	**D-Valine extrapure CHR, 99%**	5 Gms	600
		25 Gms	1900
		100 Gms	6500
14159 [516-06-3]	**DL-Valine extrapure CHR, 99%**	25 Gms	800
		100 Gms	2900
		500 Gms	13500
		1 Kg	25000
83185 [2462-34-2]	**L-Valine Benzyl Ester Hydrochloride extrapure, 98%**	25 Gms	10520
38052 [6306-52-1]	**L-Valine Methyl Ester Hydrochloride extrapure, 99%**	5 Gms	960
		25 Gms	2500
		100 Gms	5000

Code	Name	Packing	Price
40652 [2001-95-8]	**Valinomycin, 95%**	10 Mg	4500
45344 [7440-62-2]	**Vanadium Metal Powder, 99.7%**	5 Gms	1700
		25 Gms	7000
61078 [1404-93-9]	**Vancomycin Hydrochloride (VA), 90%**	1 Gms	3000
		5 Gms	7000
72241 [55-10-7,2394 -20-9]	**DL-4-Hydroxy-3-Methoxymandelic Acid (VMA) extrapure, 99%** (Vanillylmandelic Acid, 4-Hydroxy-3-Methoxy-DL-Mandelic Acid)	250 Mg	1300
		1 Gms	3200
		5 Gms	9500
11547 [120-14-9]	**n-Veratraldehyde pure, 98%**	100 Gms	1500
		250 Gms	3000
96011 [91-16-7]	**Veratrole pure, 98%** (1,2-Dimethoxy Benzene)	100 ml	600
		500 ml	2200
77472 [68-19-9]	**Vitamin B12 for molecular biology, 97%** (Cyanocobalamin)	250 Mg	770
		1 Gms	1800
		5 Gms	7000
		25 Gms	24000
56437 [68-19-9]	**Vitamin B12 for tissue culture, 97%** (Cyanocobalamin)	250 Mg	1100
		1 Gms	2400
		5 Gms	8800
79446 [13422-55-4]	**Vitamin B12a (B12 Analog) ExiPlus, Multi-Compendial, 98%-102%** (Methylcobalamin, Methyl-5,6-Dimethylbenzimidazolylcob alamin, Methyl Vitamin B12) meets compendial specs of USP	1 Gms	3200
		5 Gms	5600
		25 Gms	28000
54911 [84-80-0]	**Vitamin K1 extrapure, 98%** (Phytonadione, Phyllaquinone, 2-Methyl-3-phytyl-1,4-Napthoquinone, 3-Phytylmenadione)	1 Gms	1580
		5 Gms	7140
50342 [3153-26-2]	**Vanadyl Acetylacetonate pure, 98%** (Vanadyl Acetoacetate, Vanadium(IV) oxide acetylacetonate, VO(acac)2, Vanadium(IV)-oxy acetylacetona	25 Gms	1050
		100 Gms	3600
23761 [123334-20-3]	**Vanadyl Sulphate Hydrate pure, 97%** (Vanadium (IV) Oxide Sulfate Hydrate, Oxysulphato Vanadium (IV) Hydrate, Vanadium Oxysulphate)	25 Gms	700
		100 Gms	2300
		500 Gms	10500
40996 [83480-29-9]	**Voglibose extrapure, 98%** (Basen, Glustat, N-(1,3-Dihydroxy-2-propyl)valiolamine)	100 Mg	3000
		500 Mg	9000
66532 [201058-08-4]	**Wang Resin, 1% DVB, 0.4-2mmol/g, 100-200 mesh** (p-Benzyloxybenzyl Alcohol Resin, DVB cross linked)	5 Gms	3210
		25 Gms	12600
		100 Gms	28350
76374 [247021-90-5]	**Weinreb Linker extrapure, 98%** (N-FMOC-N-Methoxy-3-Aminopropioni c Acid)	1 Gms	8310
13111 [68988-92-1]	**Wright Stain** (Eosine Methylene Blue solution according to Wright)	25 Gms	600
		100 Gms	2340
		500 Gms	11000
72465 [193149-74-5]	**WST-8, 95%** (2-(2-Methoxy-4-Nitrophenyl)-3-(4-Nitrophenyl)-5-(2,4-Disulfophenyl) -2H-Tetrazolium Sodium Salt)	50 Mg	8000

Code	Name	Packing	Price
52714	🌡 **X-Gal Solution (20mg/ml)**	1 ml	2140
20748 [69-89-6]	✔ **Xanthine extrapure CHR, 99%**	5 Gms 25 Gms	820 3280
39607 [9002-17-9]	🌡 **Xanthine Oxidase Suspension ex. Butter Milk, 1-4 U/mg protein**	5 Units 25 Units	4500 9000
62965 [9002-17-9]	■ **Xanthine Oxidase ex. Microorganism, 15U/mg powder**	250 Units	14610
94623 [6754-58-1]	🌡 **Xanthohumol ex. Hop (Humulus Lupulus), 94%**	5 Mg	18690
62630 [59-00-7]	✔ **Xanthurenic Acid extrapure, 90%**	500 Mg 1 Gms	3200 5500
55962 [161265-03-8]	● **XantPhos (4,5-Bis(Diphenylphosphino) -9,9-Dimethylxanthine) extrapure, 98%**	1 Gms 5 Gms 25 Gms 100 Gms	1200 3500 13500 50000
27540 [564483-18-7]	● **X-Phos (2-Dicyclohexylphosphino-2,4,6-Trii sopropylbiphenyl) extrapure, 98%**	1 Gms 5 Gms 25 Gms	1200 2600 9000
91186 [9014-63-5]	✔ **Xylan from Beechwood**	5 Gms 10 Gms 25 Gms 100 Gms	2700 5000 9000 40000
30565 [9014-63-5]	✔ **Xylan from Corn Cob, 95%**	10 Gms 25 Gms 100 Gms 500 Gms	1200 2500 7000 24000
13814 [37278-89-0]	● **Xylanase ex. Aspergillus Niger, 60MU/g powder**	25 Gms 100 Gms 500 Gms	1230 4000 18000
79074 [3761-53-3]	**Xylidine Ponceau 2R, 60%** (Acid Red 26, Ponceau 2 R, Xylidine Ponceau)	25 Gms 100 Gms	650 2400
64949 [87-99-0]	✔ **Xylitol extrapure**	100 Gms 500 Gms 1 Kg	300 900 1600
84974 [58-86-6]	**D-(+)-Xylose extrapure**	100 Gms 500 Gms	420 1100
29426 [58-86-6]	**D-(+)-Xylose ExiPlus, Multi-Compendial** meets compendial specs of USP, BP	100 Gms 500 Gms	510 1300
98468 [609-06-3]	**L-(-)-Xylose extrapure, 99%**	5 Gms 25 Gms	1500 6000
31832 [65-19-0]	■ **Yohimbine Hydrochloride extrapure, 99%**	1 Gms 5 Gms 25 Gms	2000 4200 12500
36744 [10361-91-8]	**Ytterbium (III) Chloride anhydrous, 99.9%**	5 Gms	8000
24099 [10035-01-5]	**Ytterbium (III) Chloride Hexahydrate extrapure, 99.9%** (Ytterbium Trichloride Hexahydrate)	10 Gms 50 Gms	5000 12000
20435 [35725-34-9]	**Ytterbium (III) Nitrate Pentahydrate extrapure, 99.9%**	10 Gms 50 Gms	6000 18000
58985 [1314-37-0]	**Ytterbium (III) Oxide extrapure, 99.9%**	10 Gms 50 Gms	7000 20000
61813 [252976-51-5]	**Ytterbium (III) Trifluoromethanesulfonate Hydrate extrapure, 98%** (Ytterbium triflate, Trifluoromethanesulfonic Acid Ytterbium Salt)	5 Gms 25 Gms	9000 25000
63190 [10361-92-9]	**Yttrium (III) Chloride anhydrous ultrapure, 99.99%** (Yttrium Trichloride)	10 Gms 50 Gms	7000 18000
63390 [10025-94-2]	**Yttrium (III) Chloride Hexahydrate extrapure AR, 99.9%**	5 Gms 25 Gms 100 Gms	2000 4000 9000
44716 [7440-65-5]	**Yttrium Metal Ultrapure, 200 mesh, 99.9%**	5 Gms 25 Gms	12000 35000
51515 [13494-98-9]	**Yttrium (III) Nitrate Hexahydrate extrapure AR, 99.9%**	25 Gms 100 Gms	3500 9000
21412 [52093-30-8]	**Yttrium (III) Trifluoromethanesulfonate extrapure, 98%** (Yttrium triflate, Trifluoromethanesulfonic Acid Yttrium Salt)	5 Gms 25 Gms	9000 32000
37482	● **Z-L-Amino Acids Kit**	1 Kit	18000
68280 [2304-94-1]	● **Z-b-Alanine extrapure, 99%**	5 Gms 25 Gms	1500 3500
15701 [26607-51-2]	● **Z-D-Alanine extrapure, 99%**	5 Gms 25 Gms	2500 6500
24368 [4132-86-9]	● **Z-DL-Alanine extrapure, 99%**	5 Gms	2410
46096 [1142-20-7]	● **Z-L-Alanine extrapure, 99%**	5 Gms 25 Gms 100 Gms	540 2220 7970
30211 [42918-86-5]	● **Z-2-Aminobutyric Acid extrapure, 99%**	1 Gms 5 Gms 25 Gms	950 4200 20330
94475 [6382-93-0]	● **Z-D-Arginine extrapure, 99%**	1 Gms 5 Gms 25 Gms	6080 20240 73570
60020 [1234-35-1]	● **Z-L-Arginine extrapure, 99%**	25 Gms	4210
55237 [56672-63-0]	● **Z-L-Arginine HCl extrapure, 99%**	25 Gms 100 Gms	5230 19600
92873 [4474-86-6]	● **Z-D-Asparagine extrapure, 99%**	1 Gms 5 Gms 25 Gms	1090 4490 18710
60164 [29880-22-6]	● **Z-DL-Asparagine extrapure, 98%**	1 Gms 5 Gms	1520 4980

Code	Name	Packing	Price
13328 [2304-96-3]	● **Z-L-Asparagine extrapure, 99%**	1 Gms	310
		5 Gms	990
		25 Gms	2590
		100 Gms	8630
32899 [78663-07-7]	● **Z-D-Aspartic Acid extrapure, 99%**	1 Gms	640
		5 Gms	2660
		25 Gms	11080
84275 [4515-21-3]	● **Z-DL-Aspartic Acid extrapure, 99%**	5 Gms	700
		25 Gms	2400
63143 [1152-61-0]	● **Z-L-Aspartic Acid extrapure, 99%**	5 Gms	460
		10 Gms	860
		25 Gms	2060
78064 [3479-47-8]	● **Z-L-Aspartic Acid-4-Benzyl Ester extrapure, 98%** (Z-L-Asp(OBzl)-OH)	10 Gms	13560
		25 Gms	23730
97300 [5680-86-4]	● **Z-L-Glutamic Acid-5-Benzyl Ester extrapure, 98%**	10 Gms	14130
		25 Gms	24720
73971 [63648-73-7]	● **Z-D-Glutamic Acid extrapure, 98%**	1 Gms	760
		5 Gms	2660
67271 [1155-62-0]	● **Z-L-Glutamic Acid extrapure, 99%**	5 Gms	360
		25 Gms	1460
		100 Gms	5220
83309 [13139-52-1]	● **Z-D-Glutamine extrapure, 99%**	1 Gms	1230
		5 Gms	5110
		25 Gms	21270
44096 [2650-64-8]	● **Z-L-Glutamine extrapure, 99%**	5 Gms	750
		25 Gms	2590
		100 Gms	8630
86632 [1138-80-3]	● **Z-Glycine extrapure, 99%**	5 Gms	700
		25 Gms	1420
		100 Gms	4320
21793 [19728-57-5]	● **Z-DL-Histidine extrapure, 99%**	1 Gms	3480
		5 Gms	16560
57156 [14997-58-1]	● **Z-L-Histidine extrapure, 99%**	1 Gms	1200
		5 Gms	5430
		25 Gms	21290
14088 [13504-85-3]	● **Z-L-Hydroxyproline extrapure, 99%**	1 Gms	1400
		5 Gms	6330
		25 Gms	24720
34913 [3160-59-6]	● **Z-L-Isoleucine extrapure, 95%**	5 Gms	2010
		25 Gms	6280
46817 [28862-79-5]	● **Z-D-Leucine extrapure, 99%**	1 Gms	1090
		5 Gms	4510
		25 Gms	18800
76353 [3588-60-1]	● **Z-DL-Leucine extrapure, 99%**	5 Gms	1830
		25 Gms	8700
85612 [2018-66-8]	● **Z-L-Leucine extrapure, 98%**	5 Gms	1520
		25 Gms	4330

Code	Name	Packing	Price
39218 [70671-54-4]	● **Z-D-Lysine extrapure, 99%**	1 Gms	5000
		5 Gms	12000
54420 [2212-75-1]	● **Z-L-Lysine extrapure, 99%**	1 Gms	1280
		5 Gms	5900
71691 [2212-76-2]	● **Z-N-e-Lysine Dicyclohexyl Ammonium Salt extrapure, 98%** (Z-Lysine(BOC)-OH.Dicyclohexylamine Salt)	5 Gms	5950
91642 [28862-80-8]	● **Z-D-Methionine extrapure, 99%**	1 Gms	1500
		5 Gms	3000
11751 [4434-61-1]	● **Z-DL-Methionine extrapure, 99%**	1 Gms	740
		5 Gms	3480
55768 [1152-62-1]	● **Z-L-Methionine extrapure, 99%**	1 Gms	690
		5 Gms	2850
		25 Gms	9350
21200	● **Z-4-Nitro-L-Phenylalanine extrapure, 99%** ((S)-2-(((Benzyloxy)Carbonyl)Amino)-3-(4-Nitrophenyl) Propanoic Acid, Cbz-L-4-nitrophenylalanine)	10 Gms	10840
		25 Gms	18970
66829 [15027-14-2]	● **Z-D-Norleucine extrapure, 99%**	500 Mg	4610
		2.5 Gms	20700
12185 [15027-13-1]	● **Z-DL-Norleucine extrapure, 99%**	1 Gms	3150
90175 [39608-30-5]	● **Z-L-Norleucine extrapure, 99%**	1 Gms	4130
		5 Gms	18790
34201	● **Z-D-Norvaline extrapure, 99%**	1 Gms	2790
		5 Gms	12660
21693 [21691-43-0]	● **Z-DL-Norvaline extrapure, 99%**	1 Gms	2310
		5 Gms	10480
75427 [21691-44-1]	● **Z-L-Norvaline extrapure, 99%**	1 Gms	2790
		5 Gms	12660
72850 [3304-51-6]	● **Z-L-Ornithine extrapure, 98%**	5 Gms	3000
		25 Gms	6000
24317 [2448-45-5]	● **Z-D-Phenylalanine extrapure, 99%**	1 Gms	640
		5 Gms	2660
		25 Gms	11040
32900 [3588-57-6]	● **Z-DL-Phenylalanine extrapure, 99%**	5 Gms	2310
		25 Gms	10480
42182 [1161-13-3]	● **Z-L-Phenylalanine extrapure, 99%**	25 Gms	2240
		100 Gms	7910
27486 [35909-92-3]	● **Z-L-Phenyalanine Methyl Ester extrapure, 98%**	10 Gms	12240
		25 Gms	21410
41739 [19647-71-3]	● **Z-Phenylalanine-p-Nitroanilide extrapure**	1 Gms	7590
16881 [2578-84-9]	● **Z-L-Phenylalanine-4-Nitrophenyl Ester extrapure, 98%**	1 Gms	3000
		10 Gms	11020
12399 [17609-52-8]	● **Z-D-Phenylglycine extrapure, 98%**	10 Gms	7400
		25 Gms	12940
74157 [53990-33-3]	● **Z-L-Phenylglycine extrapure, 99%**	10 Gms	6490
		25 Gms	11360

Code	Name	Packing	Price
24037 [6404-31-5]	● **Z-D-Proline extrapure, 99%**	500 Mg	930
		1 Gms	1570
		5 Gms	6040
		25 Gms	30180
31688 [1148-11-4]	● **Z-L-Proline extrapure, 99%**	5 Gms	1260
		25 Gms	3780
63273 [6081-61-4]	● **Z-D-Serine extrapure, 99%**	500 Mg	2180
		1 Gms	3360
		5 Gms	14030
75675 [1145-80-8]	● **Z-L-Serine extrapure, 99%**	1 Gms	420
		5 Gms	1100
		25 Gms	5000
52287 [80384-27-6]	● **Z-D-Threonine extrapure, 99%**	1 Gms	1300
		5 Gms	5200
11228 [19728-63-3]	● **Z-L-Threonine extrapure, 99%**	5 Gms	1300
39868 [2279-15-4]	● **Z-D-Tryptophan extrapure, 99%**	1 Gms	1800
		5 Gms	5500
57591 [7432-21-5]	● **Z-L-Tryptophan extrapure, 99%**	25 Gms	2100
59268 [64205-12-5]	● **Z-D-Tyrosine extrapure, 99%**	1 Gms	3700
		5 Gms	10000
31078 [1164-16-5]	● **Z-L-Tyrosine extrapure, 99%**	5 Gms	1950
		25 Gms	4980
39463 [16677-29-5]	● **Z-O-Benzyl-L-Tyrosine extrapure, 98%** (Z-L-Tyr(Bzl)-OH)	10 Gms	7580
74556 [1685-33-2]	● **Z-D-Valine extrapure, 99%**	1 Gms	1320
		5 Gms	3000
		25 Gms	10000
17520 [3588-63-4]	● **Z-DL-Valine extrapure, 99%**	5 Gms	2790
		25 Gms	13250
73743 [1149-26-4]	● **Z-L-Valine extrapure, 99%**	25 Gms	1800
		100 Gms	4680
92186 [10314-98-4]	● **Z-Isonipecotic Acid extrapure, 99%** (1-Carbobenzoxy-4-Piperidinecarboxylic Acid)	10 Gms	4000
		25 Gms	6000
88291 [78190-11-1]	● **Z-Nipecotic Acid, 95%**	10 Gms	6470
70443 [13139-17-8]	● **Z-OSu (Z-ONSu) extrapure, 99%** (Z-ONSu)	25 Gms	1200
		100 Gms	4500
81055 [13114-27-7]	🌡 **Zeatin mixed Isomers for tissue culture, ~95%**	25 Mg	2200
		100 Mg	6000
67299 [13114-27-7]	🌡 **Zeatin mixed Isomers, ~95%**	25 Mg	1800
		100 Mg	5000
		250 Mg	10500
80440 [6025-53-2]	🌡 **Zeatin Riboside extrapure, 99%**	10 Mg	2000
		25 Mg	3200
		100 Mg	8200
		500 Mg	36000


Code	Name	Packing	Price
25910 [7699-45-8]	**Zinc Bromide Anhydrous, 99%** (Zinc Dibromide)	100 Gms	1100
		500 Gms	2420
86849 [62625-22-3]	**Zincon Sodium Salt ACS, 98%**	1 Gms	450
		5 Gms	1400
52049 [1314-23-4]	**Zirconium Dioxide pure, 99%** (Zirconia, Zirconium (IV) Oxide)	100 Gms	840
		500 Gms	3000
94059 [7704-99-6]	**Zirconium Hydride extrapure,10-30µm, 99.9%** (Zirconium Dihydride)	10 Gms	8000
		50 Gms	30000

Part B1
Animal Cell Culture

SeraPreserve™

Formalin-Free Tissue Fixative Concentrate





A safe and advanced alternative. SeraPreserve™ is a revolutionary formalin-free tissue fixative designed to deliver superior preservation, hardening, and staining enhancement of biological tissues. It is the ideal choice for histology, cytology, pathology, and immunohistochemistry, ensuring clear, detailed tissue analysis with a focus on lab safety and diagnostic accuracy.

Why Choose SeraPreserve™?

- **Formalin-Free Safety:** Non-carcinogenic formulation – contains no formalin or formaldehyde
- **Optimal Preservation:** Glycols and pH-adjusted organics maintain tissue integrity and clarity
- **Reliable Stability:** Consistent performance at room temperature and neutral pH
- **Staining Enhancement:** Improves dye uptake and cellular detail visibility
- **Easy Handling & Extended Storage:** Facilitates sample preparation, slicing, and long-term storage
- **Versatile Use:** Suitable for IHC, histopathology, diagnostics, forensics, biobanking & research

Applications:

- Histology & Pathology
- Immunohistochemistry (IHC)
- Cytology (smear preparations)
- Clinical & Research Laboratories
- Biobanks & Forensic Labs



Available In:

Concentration	Prdt. Code	Pack Size
10×	15776	5 Ltr
25×	40148	5 Ltr
50×	21569	5 Ltr

Code	Name	Packing	Price
68161 [7365-82-4]	**ACES Buffer for cell culture, 99.5%, Endotoxin (BET) 1EU/mg**	25 Gms 100 Gms	5200 12250
23568 [87-32-1]	**N-Acetyl-DL-Tryptophan for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	5 Gms 25 Gms	3000 9000
52937 [3244-88-0]	**Acid Fuchsin for cell culture, Endotoxin (BET) 0.05EU/mg** (Acid Violet 19, Fuchsin S, Fuchsin acid, Rubine S) (C.I. No. 42685)	25 Gms	1600
75250 [56-41-7]	**L-Alanine for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	25 Gms 100 Gms	1050 2490
96409 [9048-46-8]	**Albumin Bovine (pH 7) fraction V for cell culture (Bovine Serum Albumin,BSA), 98%, (BET) 3EU/mg** (Bovine Serum Albumin, BSA)	5 Gms 25 Gms 100 Gms	2000 5000 17500
31539 [33864-99-2]	**Alcian Blue for cell culture, Endotoxin (BET) 0.05EU/mg**	5 Gms 25 Gms	4000 8000
68860 [7784-24-9]	**Aluminium Potassium Sulphate Dodecahydrate for cell culture, 99.5%, Endotoxin (BET) 0.05EU/mg**	500 Gms 1 Kg	1560 2880
15875 [1397-89-3]	**Amphotericin B (AMT) Solution (250ug/ml) for cell culture, Endotoxin (BET) 1EU/ml**	20 ml 100 ml	1860 5360
20969 [74-79-3]	**L-Arginine (free base) for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	25 Gms 100 Gms	1150 3550
33468 [1405-87-4]	**Bacitracin (BCT) ex. Bacillus Licheniformis for cell culture, Endotoxin (BET) 0.05EU/mg**	50000 Units 250000 Units	1400 4500
36884 [10191-18-1]	**BES Buffer for cell culture, 99%, Endotoxin (BET) 0.05EU/mg** (N,N-bis(2-Hydroxyethyl)2-Aminoethane Sulphonic Acid)	25 Gms 100 Gms	1650 5400
92370 [3513-03-9]	**Blasticidin S Hydrochloride for cell culture, 98%, Endotoxin (BET) 0.05EU/mg**	10 Mg 25 Mg	6120 14470
55594 [3513-03-9]	**Blasticidin S HCl Solution in HEPES Buffer (10mg/ml) Sterile for cell culture, Endotoxin (BET)1EU/ml**	1 ml 5 x1 ml	7000 25000
37936 [104376-79-6]	**Ceftriaxone Disodium Salt Hemiheptahydrate (CFTZ) for cell culture, Endotoxin (BET) 0.2EU/mg**	1 Gms 5 Gms	4500 14000
81642 [64485-93-4]	**Cephotaxime Sodium Salt (CFT) for cell culture, 916-964ug/mg, Endotoxin (BET) 0.20EU/mg**	100 Mg 1 Gms	2500 6000
84483 [57-09-0]	**Cetyltrimethyl Ammonium Bromide (CTAB) for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	25 Gms 100 Gms	2500 6000
17796 [56-75-7]	**Chloramphenicol cell culture, 99%, Endotoxin (BET) 0.2 EU/mg**	5 Gms 25 Gms 100 Gms	2600 7000 18000
97900 [57-88-5]	**Cholesterol ex. Sheep Wool (Lanolin) for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	5 Gms 25 Gms 100 Gms	2400 10500 30000
70553 [86393-32-0]	**Ciprofloxacin Hydrochloride Hydrate (CPFX) for cell culture, 98-102%, Endotoxin (BET) 0.05EU/mg**	1 Gms 5 Gms	690 2300
81395 [11028-71-0]	**Concanavalin A Lectin ex. Canavalia Ensiformis Type 1 for cell culture, 64ug/mg, (BET) 0.05EU/mg** (Con A) (Jack bean)	25 Mg 40 Mg	3680 5460
68722 [11028-71-0]	**Concanavalin A Lectin ex. Canavalia Ensiformis Type 2 for cell culture, 20 ug/ml, (BET) 0.05EU/mg**	25 Mg 100 Mg	2670 6010
15353 [60-27-5]	**Creatinine for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	25 Gms 100 Gms 1 Kg	1370 4350 42000
67133 [10510-54-0]	**Cresyl Violet Acetate for cell culture, Endotoxin (BET) 0.5EU/mg**	5 Gms	10000
29751 [548-62-9]	**Crystal Violet for cell culture, Endotoxin (BET) 1EU/mg** (Methyl Violet 10B) (C.I. No.42555)	25 Gms 100 Gms	1750 5800
70134 [7048-04-6]	**L-Cysteine Hydrochloride Monohydrate for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	25 Gms 100 Gms	1150 3250
54676 [50-99-7]	**Dextrose for cell culture, Endotoxin (BET) 0.05EU/mg** (D-Glucose Anhydrous)	100 Gms 500 Gms	500 1450
28525 [3483-12-3]	**DL-Dithiothreitol (DTT) for cell culture, 98%, Endotoxin (BET) 1EU/mg** (Clelands reagent)	1 Gms 5 Gms	2000 7000
86058	**Dulbecco`s Modified Eagles Medium 1X w/ 4500mg/lit Glucose, L-Glutamine, Sodium Pyruvate, Sodium Bic** (DMEM High Glucose, 1X DMEM)	500 ml	3500
55654	**Dulbecco`s Phosphate Buffered Saline 1X w/o Calcium & Magnesium** (1X DPBS w/o Ca, Mg)	500 ml	2500
70520	**Dulbecco`s Phosphate Buffered Saline 10X w/o Calcium and Magnesium** (10X DPBS w/o Ca, Mg)	500 ml	4000
81614 [17372-87-1]	**Eosin Yellow (water soluble) for cell culture, Endotoxin (BET) 1EU/mg**	25 Gms 100 Gms	2250 4500
74516 [122-99-6]	**Ethylene Glycol Monophenyl Ether for cell culture, 99%, Endotoxin (BET) 0.05EU/ml** (Phenyl Cellosolve, 2-Phenoxyethanol)	100 ml 500 ml	2400 4500
11324 [7782-63-0]	**Ferrous Sulphate Heptahydrate for cell culture, 99.5%, Endotoxin (BET) 2.0 EU/mg** (Iron (II) Sulphate Heptahydrate)	500 Gms	1200
98131 [518-47-8]	**Fluorescein Sodium Salt (C.I. No. 45350) for cell culture, Endotoxin (BET) 0.05EU/mg** (Acid Yellow 73)	100 Gms	2500
76611 [59-30-3]	**Folic Acid for cell culture, 98%, Endotoxin (BET) 0.05EU/mg** (Vitamin B9, Vitamin M, Pteroyl-L-Glutamic acid)	5 Gms 25 Gms	1200 5000
20596 [632-99-5]	**Fuchsin Basic for cell culture, Endotoxin (BET) 0.05EU/mg** (Basic Fuchsin) (C.I. No. 42510)	25 Gms 100 Gms	2200 8000

Code	Name	Packing	Price
77469 [51811-82-6]	**Giemsa Stain for cell culture, Endotoxin (BET) 0.05EU/mg**	5 Gms 25 Gms	4000 15500
75931 [6106-04-3]	**L-Glutamic Acid Monosodium Salt Monohydrate (MSG) for cell culture, 99%, Endotoxin (BET) 0.05EU/mg** (Monosodium Glutamate, Monosodium L-Glutamate)	100 Gms 500 Gms	900 2750
80952 [56-85-9]	**L-Glutamine for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	10 Gms 25 Gms	950 1800
74693 [27025-41-8] ●	**Glutathione Oxidized (GSSG) for cell culture, 99%, Endotoxin (BET) 1EU/mg**	100 Mg 250 Mg 1 Gms 5 Gms	2000 3700 7000 32500
65542 [56-81-5]	**Glycerol (Glycerine) Anhydrous for cell culture, 99.5%, Endotoxin (BET) 0.05EU/mg**	500 ml	1300
34268 [56-40-6]	**Glycine for cell culture, 99.5%, Endotoxin (BET) 1EU/mg** (Aminoacetic Acid)	100 Gms	1700
39149 [7365-45-9]	**HEPES Buffer for cell culture, 99.5%, Endotoxin (BET) 0.05EU/mg** (N-(2-Hydroxyethyl) piperazine-N(2-Ethane Sulphonic Acid)	25 Gms 100 Gms 500 Gms	1950 6400 26500
16509 [71-00-1]	**L-Histidine (base) for cell culture, 99%, Endotoxin (BET) 1EU/mg**	25 Gms 100 Gms	1250 3500
86877 [5934-29-2]	**L-Histidine Hydrochloride Monohydrate for cell culture, 99%, Endotoxin (BET) 1EU/mg**	25 Gms 100 Gms	1100 3500
54661 [3258-02-4] ■	**N-Hydroxycytidine (EIDD-1931) for cell culture, 98%, Endotoxin (BET) 0.05EU/mg**	10 Mg 25 Mg	7000 14000
75054 [99011-02-6] ✔	**Imiquimod for cell culture, 98%, Endotoxin (BET) 1EU/mg**	250 Mg	13500
77379 [87-89-8]	**Inositol for cell culture, 99%, Endotoxin (BET) 0.05EU/mg** (Vitamin B8)	25 Gms 100 Gms	1000 3200
72828 [11061-68-0] ■	**r-Insulin ex. Human for cell culture, 25 USP U/mg, Endotoxin (BET) 1.0EU/mg** (Recombinant Insulin)	50 Mg 100 Mg	7000 13000
11894 [73-32-5]	**L-Isoleucine for cell culture, 99%, Endotoxin (BET) 1EU/mg**	25 Gms 100 Gms	1300 4300
95793 [62851-42-7]	**Jenners Stain (Eosin-methylene blue) for cell culture, Endotoxin (BET) 0.05EU/mg**	25 Gms 100 Gms	3100 9800
95419 [61-90-5]	**L-Leucine for cell culture, 99%, Endotoxin (BET) 1EU/mg**	100 Gms 500 Gms	3000 8000
95173 [56-87-1] 🌡	**L-Lysine (free base) Anhydrous for cell culture, 98%, Endotoxin (BET) 0.05EU/mg**	25 Gms 100 Gms	4000 12000
29011 [12650-88-3] ■	**Lysozyme (3x cryst) ex.Egg white (Muramidase) for cell culture, 15000U/mg, Endotoxin (BET) 0.05EU/mg**	1 Gms 5 Gms 25 Gms	3000 13500 40000
22352 [7786-30-3]	**Magnesium Chloride Anhydrous for cell culture, 98%, Endotoxin (BET) 0.05EU/mg**	100 Gms	1800
13546 [7791-18-6]	**Magnesium Chloride Hexahydrate for cell culture, 99.5%, Endotoxin (BET) 0.05EU/mg**	500 Gms	3000
36502 [7487-88-9]	**Magnesium Sulphate Dried for cell culture, 62-70% MgSO4, 99%, Endotoxin (BET) 0.05EU/mg**	500 Gms	1750
12493 [6363-53-7]	**D-Maltose Monohydrate for cell culture, 92%, Endotoxin (BET) 0.05EU/mg**	100 Gms 500 Gms	1000 4000
20490 [60-24-2]	**2-Mercaptoethanol for cell culture, 99%, Endotoxin (BET) 0.05EU/ml**	100 ml 250 ml	1300 2650
55311 [145224-94-8]	**MES Monohydrate Buffer for cell culture, 99%, Endotoxin (BET) 0.05EU/mg** (2-Morpholineethane Sulphonic Acid Hydrate)	25 Gms 100 Gms 500 Gms	3800 8500 36000
42125 [63-68-3]	**L-Methionine for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	25 Gms 100 Gms	1100 4000
13908 [1132-61-2]	**MOPS Buffer for cell culture, 99%, Endotoxin (BET) 0.05EU/mg** (3-Morpholinopropane Sulphonic Acid)	25 Gms 100 Gms 500 Gms	3200 7800 17000
69599 [1405-10-3] 🌡	**Neomycin Sulphate (NMS) for cell culture, Endotoxin (BET) 0.05EU/mg**	5 Gms 25 Gms	800 1900
64257 [53-84-9] ■	**ß-Nicotinamide Adenine Dinucleotide (ß-NAD,DPN) for cell culture, 98%, Endotoxin (BET) 0.05EU/mg**	100 Mg 1 Gms 5 Gms	1010 4680 17250
10581 [606-68-8] ■	**ß-Nicotinamide Adenine Dinucleotide Disodium Salt(Reduced) (ß-NADH.Na2,DPNH.Na2) for cellculture,98%**	100 Mg 1 Gms 5 Gms	2010 10290 37830
88722 [24292-60-2] ■	**ß-Nicotinamide Adenine Dinucleotide Phosphate Disodium Salt for cell culture,95%, (BET) 0.05EU/mg**	100 Mg 1 Gms	3680 21140
28771 [1184-16-3] ■	**ß-Nicotinamide Adenine Dinucleotide Phosphate Monosodium Salt for cell culture, 98%, (BET) 0.05EU/mg**	100 Mg 1 Gms	3560 13350
93517 [98-92-0]	**Nicotinamide for cell culture, 99.5%, Endotoxin (BET) 0.05EU/mg** (Niacinamide, Pyridine-3-Carboxylic Acid Amide)	100 Gms 500 Gms	1250 3750
63509 [1094-61-7] ■	**ß-Nicotinamide Mononucleotide (ß-NMN) for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	25 Mg 1 Gms 5 Gms	500 2500 9000
50275 [59-67-6]	**Nicotinic Acid (Pyridine-3-Carboxylic Acid) for cell culture, 99%, Endotoxin (BET) 0.05EU/mg** (3-Picolinic Acid, Niacin, Vitamin B3)	100 Gms 500 Gms	950 3000
16828 [1400-61-9] 🌡	**Nystatin (Mycostatin) for cell culture, 85%, Endotoxin (BET) 0.05EU/mg** (Fungicidin)	250 Mg 1 Gms 5 Gms	670 2340 8900
17565 [1400-61-9] ■	**Nystatin (Mycostatin) Solution (10000U/ml) for cell culture, Endotoxin (BET) 10.0EU/ml** (Fungicidin)	20 ml	1240


Code	Name	Packing	Price
56341 [3184-13-2]	**L-Ornithine Monohydrochloride for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	25 Gms 100 Gms	1600 4500
30917 [328-42-7]	**Oxalacetic Acid for cell culture, 98%, Endotoxin (BET) 0.05EU/mg** (Oxaloacetic Acid, Ketosuccinic Acid, Oxobutanedioic Acid, 2-Oxosuccinic Acid)	1 Gms 5 Gms	3000 11000
79311 [34487-61-1]	**Phenol Red Sodium Salt (water soluble) for cell culture, 95%, Endotoxin (BET) 0.05EU/mg** (Phenolsulfonephthalein Sodium Salt)	5 Gms 25 Gms	1100 4000
68793 [63-91-2]	**L-Phenylalanine for cell culture, 99%, Endotoxin (BET) 1EU/mg**	25 Gms 100 Gms	1100 3000
77839 [329-98-6]	**Phenylmethane Sulphonyl Fluoride (PMSF) for cell culture, 99%, Endotoxin (BET) 0.05EU/mg**	1 Gms 5 Gms	1300 4400
62374 [11006-33-0]	**Phleomycin (PLM) ex. Streptomyces verticillus for cell culture, Endotoxin (BET) 1.0EU/mg** (Zeocin)	5 Mg 25 Mg 100 Mg	3700 9000 24000
52819 [1405-20-5]	**Polymyxin B Sulphate (PB, Aerosporin), 6000U/mg for cell culture, Endotoxin (BET) 1.0EU/mg**	250 Mg 1 Gms	2200 4870
67209 [9005-65-6]	**Polysorbate 80 (Tween 80) for cell culture, Endotoxin (BET) 0.05EU/mg** (Polyoxyethylene Sorbitan Monoleate)	500 Gms	3000
47425 [7447-40-7]	**Potassium Chloride for cell culture, 99.5%, Endotoxin (BET) 0.05EU/mg**	500 Gms	3450
65897 [147-85-3]	**L-Proline for cell culture, 99%, Endotoxin (BET) 1.0EU/mg**	25 Gms 100 Gms	1050 3350
82361 [58-58-2]	**Puromycin Dihydrochloride (PRM) for cell culture, Endotoxin (BET) 10EU/mg**	10 Mg 25 Mg	5000 9000
14050 [58-56-0]	**Pyridoxine Hydrochloride (Vitamin B6) for cell culture, 99%, Endotoxin (BET) 1EU/mg**	10 Gms 25 Gms	1500 3500
77387 [83-88-5]	**Riboflavin for cell culture, 98%, Endotoxin (BET) 0.05EU/mg** (Vitamin B2)	25 Gms 100 Gms	1150 3500
46143 [13292-46-1]	**Rifampicin (RFP) (Rifampin) for cell culture, 97-102%, Endotoxin (BET) 0.05EU/mg**	1 Gms 5 Gms	4200 15000
84130	**RPMI 1640 Medium 1X**	500 ml	2500
51024	**RPMI 1640 Medium w/ L-Glutamine w/ 25mM HEPES**	500 ml	4000
31904 [477-73-6]	**Safranine O for cell culture, Endotoxin (BET) 0.05EU/mg** (Basic Red 2)	25 Gms 100 Gms	5300 9500
59260 [56-45-1]	**L-Serine for cell culture, 99%, Endotoxin (BET) 1EU/mg**	25 Gms 100 Gms	2500 5300
16670 [156-54-7]	**Sodium Butyrate for cell culture, 99%, Endotoxin (BET) 1.0EU/mg** (Butyric Acid Sodium Salt)	5 Gms 100 Gms	3000 13000
89360 [7647-14-5]	**Sodium Chloride for cell culture, 99.9%, Endotoxin (BET) 0.05EU/mg**	500 Gms 1 Kg 5 Kg	2000 3900 12750
66003 [7558-79-4]	**Sodium Phosphate Dibasic Anhydrous for cell culture, 99.5%, Endotoxin (BET) 0.05EU/mg** (Disodium Hydrogen Phosphate Anhydrous)	500 Gms 5 Kg	3500 32000
16120 [10028-24-7]	**Sodium Phosphate Dibasic Dihydrate for cell culture, 99.5%, Endotoxin (BET) 0.05EU/mg** (Disodium hydrogen phosphate Dihydrate)	500 Gms 5 Kg	1700 10000
27621 [7558-80-7]	**Sodium Phosphate Monobasic Anhydrous for cell culture, 99%, Endotoxin (BET) 0.05EU/mg** (Sodium Dihydrogen Phosphate Anhydrous)	100 Gms 500 Gms	2000 5500
68350 [18883-66-4]	**Streptozotocin (STZ) for cell culture, 98%, Endotoxin (BET) 0.05EU/mg**	100 Mg 1 Gms	2800 19000
68250 [57-50-1]	**Sucrose for cell culture, Endotoxin (BET) 0.05EU/mg**	500 Gms 1 Kg 5 Kg	3000 5500 22000
88726 [72-19-5]	**L-Threonine for cell culture, 99%, Endotoxin (BET) 1EU/mg**	25 Gms 100 Gms	2500 7500
53073 [4697-14-7]	**Ticarcillin Disodium Salt extrapure, 85%, Endotoxin (BET) 0.05EU/mg**	1 Gms 5 Gms	1950 8860
24000 [92-31-9]	**Toluidine Blue O for cell culture, 90%, Endotoxin (BET) 0.05EU/mg** (Basic Blue 17, Blutene Chloride, Methylene Blue T50, Tolonium Chloride)	10 Gms 25 Gms	3300 6200
30766 [5704-04-1]	**TRICINE Buffer for cell culture, 99%, Endotoxin (BET) 0.1EU/mg** (N-Tris(Hydroxymethyl) Methylglycine)	25 Gms 100 Gms	2400 5500
35416 [38748-32-2]	**Triptolide ex. Tripterigium Wilfordii for cell culture, 98%** (PG 490)	5 Mg	15000
12695 [77-86-1]	**Tris Buffer for cell culture, 99.9%, Endotoxin (BET) 0.05 EU/mg** (Tris(hydroxymethyl)aminomethane, Tromethamine)	100 Gms 500 Gms	950 4000
65972 [9002-07-7]	**r-Trypsin - Type Standard ACF for cell culture, 800 USP U/mg powder, Endotoxin (BET) 0.05EU/mg** (Recombinant Trypsin)	100 Mg 1 Gms 5 Gms	5090 28000 125000
74055 [9002-07-7]	**Trypsin 1:250 ex. Bovine Pancreas for cell culture, 1000BAEE U/mg, Endotoxin (BET) 10.0EU/mg**	25 Gms 100 Gms 500 Gms	2400 7500 36000
23560	**r-Trypsin 1x w/ EDTA Soln in Phosphate Buffer for cell culture (AOF),150U/ml,Endotoxin (BET) 1EU/ml**	100 ml	4600
52682 [9002-07-7]	**Trypsin ex. Bovine Pancreas for cell culture 0.2 Anson U/g (2000U/g), Endotoxin (BET) 10.0EU/mg**	100 Gms 500 Gms 1 Kg	950 4640 8550
20374 [9002-07-7]	**Trypsin ex. Porcine Gastric Mucosa for cell culture, 0.2 Anson U/g (2000U/g),(BET) 0.05EU/mg powder**	100 Gms 500 Gms 1 Kg	990 4580 8930

Code	Name	Packing	Price
80149 [60-18-4]	**L-Tyrosine for cell culture, 99%, Endotoxin (BET) 1EU/mg**	25 Gms 100 Gms	1300 4800
78168 [57-13-6]	**Urea for cell culture, 99.5%, Endotoxin (BET) 0.05EU/mg**	500 Gms 5 Kg	2000 9500
20451 [72-18-4]	**L-Valine for cell culture, 99%, Endotoxin (BET) 1EU/mg**	25 Gms 100 Gms	850 3000
12540 [68-19-9]	**Vitamin B12 for cell culture, 97%, Endotoxin (BET) 0.05EU/mg** (Cyanocobalamin)	250 Mg 1 Gms	1570 4100
18681 [111072-31-2]	**XTT Sodium for cell culture, 90%, Endotoxin (BET) 5.0EU/mg**	25 Mg 100 Mg	7500 23000
30108 [7646-85-7]	**Zinc Chloride for cell culture, 98%, Endotoxin (BET) 0.05EU/mg**	500 Gms	2700
92647 [7446-20-0]	**Zinc Sulphate Heptahydrate for cell culture, 99.5%, Endotoxin (BET) 0.05EU/mg**	100 Gms 500 Gms	1600 3400

Part C
Nanomaterials and Renewable Energy

(Carbon Nanotubes, Nanodispersants and Nanodispersions, Nanowires, Nanofilms and Nanopowders)

Linkers & Cross-Linkers for Protein Chemistry



Cross-linking is the process of chemically joining two or more molecules by a covalent bond. **SRL** introduced linking and cross-linking reagents are important reagents used in protein chemistry applications. Typically, they are commonly used reagents to modify nucleic acids, drugs and solid surfaces. Areas of research where these products find application are bio-conjugation, immobilization, surface modification and bio-molecule labeling. These reagents contain reactive ends to specific functional groups, such as primary amines and sulfhydryls, on proteins or other molecules. The availability of several chemical groups in proteins and peptides make them targets for conjugation and for study using cross-linking methods. Some common applications of cross-linkers is for the identification of near-neighbor protein relationships and for ligand-receptor interactions.

- 3-(Acetamidomethylthio)propanoic Acid (MPA(Acm))
- 3-(Benzylthio)-Propionic Acid
- DABCYL-N-Succinimidyl Ester
- DHP Linker (3,4-Dihydro-2H-Pyran-2-Methanol)
- 5-(Ethylthio)-1H-Tetrazole (ETT)
- HMBA Linker (4-(Hydroxymethyl)benzoic Acid)
- HMPA Linker
- Ramage Linker
- Rink Amide Linker (Knorr Linker)
- Sieber Linker
- Weinreb Linker

Compounds for Peptide Synthesis

Peptide synthesis is characterized as the formation of a peptide bond between two amino acids. While the definition of a peptide is not definitive, it usually refers to chains of up to 30-50 amino acids. Peptide synthesis today is a common approach in even high-throughput biological research and product and drug development, although it is a technique about 100 years old. We offer various types of peptide reagents that can be used for biological production and research work.

Amino Acids and Derivatives
- BOC
- FMOC
- L/D/DL
- Z/CBZ

N,N'-Dicyclohexylcarbodimide DCC

Diisopropylcarbodimide DIPC

1-(3 Dimethylaminopropyl)-3-Carbodimide Hydrochloride (EDC.HCl), EDAC.HCl)

1,4 Dithioerythritol (DTE)

DL-Dithiothreitol (DTT)

EEDQ

1,2 Ethanedithiol (EDT)

9-Fluorenylmethanol

9-Fluorenylmethyl-N-succinimidyl Carbonate

1-Fluoro-2,4 Dinitrobenzene

HABU

HOAT

HOBT

Hydrobromic Acid

2-Hydroxypyridine-N-oxide (1-Hydroxy-2-pyridone) (HOPO)

N-Hydroxysuccinamide

Di-(n-succinimidyl) Carbonate

TBTU

TNTU

Trifluoroacetic Acid

Wang Resin (100-200 mesh)

Z-OSU

For more details, please feel free to mail us at **marketing@srlchem.com**

Code	Name	Packing	Price
75964 [1344-28-1]	**Alumina (Alpha) Nanopowder, 20-30nm, 99.9%**	25 Gms 100 Gms 250 Gms	2760 6790 12310
25536 [1344-28-1]	**Alumina (Gamma) Nanopowder, 15nm, 99.5%** (Hydrophilic)	25 Gms 100 Gms	11160 16680
96202 [7429-90-5]	**Aluminium Nanopowder (Type 1), 40nm, 99%**	1 Gms 5 Gms	14610 38990
18412 [7429-90-5]	**Aluminium Nanopowder (Type 2), 70 nm, 99.9%**	1 Gms 5 Gms	14610 38990
41892 [1318-23-6]	**Aluminium Oxide (Boehmite) Nanodispersion (50nm), 20% in water**	25 Gms 100 Gms	3450 10350
29603 [12004-39-6]	**Aluminium Titanate Nanopowder, 25nm, 98.5%**	2 Gms	16100
19158 [7782-42-5]	**Aminated Graphene Octadecylamine covalently linked (GnAm18d) Nanopowder** (Reduced Graphene Oxide-NH-(CH2)17-CH3)	10 Mg 50 Mg	5200 16700
67801 [1309-64-4]	**Antimony Trioxide Nanopowder, 80-200nm, 99.9%** (Diantimony Trioxide, Antimony (III) Oxide)	5 Gms 25 Gms	1400 3000
10842 [12047-27-7]	**Barium Titanate Nanopowder, 100nm, 99.5%**	25 Gms 100 Gms	3300 8800
78768 [1403674-49-6]	**Bis(2`,4`-Bis(Hexyloxy)Biphenyl-4-yl)Amine extrapure, 95%** (Z-907 Modified Red Dye)	25 Mg	13400
14683 [1304-76-3]	**Bismuth (III) Oxide Nanopowder, 200nm, 99.8%** (Dibismuth Trioxide)	25 Gms 100 Gms	6600 18600
21772 [1318-23-6]	**Boehmite Nanopowder, 20nm, 99.9%**	5 Gms 25 Gms	3000 11200
62554 [12069-32-8]	**Boron Carbide Nanopowder, 50nm, 99%**	10 Gms 25 Gms	8400 17700
58903 ⓘ [7440-42-8]	**Boron Micronpowder, 45micron, 96%**	10 Gms 25 Gms	8200 18000
50943 [10043-11-5]	**Boron Nitride (hBN) Aerosol Spray (13Oz/369g)** (hexagonal Boron Nitride)	1 piece	13300
75454 [10043-11-5]	**Boron Nitride Nanopowder, 70nm, 99%**	25 Gms 100 Gms	3300 8800
45882 [10043-11-5]	**Boron Nitride Nanotubes (BNNT-C) Cylindrical, 30-50nm**	10 Mg 50 Mg	30000 95000
74899 [471-34-1]	**Calcium Carbonate Nanopowder, 80nm, 98%**	25 Gms 100 Gms 500 Gms	1100 4000 18800
72893 [12049-50-2]	**Calcium Titanate Nanopowder, 90nm, 99.9%**	1 Gms	15100
12918 [12013-47-7]	**Calcium Zirconium Oxide Nanopowder, 50nm, 99.5%**	1 Gms 5 Gms	14900 47500
69126 🌡	**CdSe/ZnS Core Shell Quantum Dots, Fluorescence 560nm, (Hydrophobic)** (CdSe/ZnS560) (Hydrophobic)	10 Mg	22400
33618 🌡	**CdSe/ZnS Core Shell Quantum Dots, Fluorescence 600nm, (Hydrophobic)** (CdSe/ZnS600) (Hydrophobic)	10 Mg	22400
15377 🌡	**CdSe/ZnS Core Shell Quantum Dots, Fluorescence 630nm, (Hydrophobic)** (CdSe/ZnS639) (Hydrophobic)	10 Mg	22400
50439 🌡	**CdTe Core Type COOH Functionalized Quantum Dots, Fluorescence 510 nm (Hydrophilic)**	10 Mg	23000
37909 🌡	**CdTe Core Type COOH Functionalized Quantum Dots, Fluorescence 610 nm (Hydrophilic)**	10 Mg	23000
19960 🌡	**CdTe Core Type COOH Functionalized Quantum Dots, Fluorescence 710 nm (Hydrophilic)**	10 Mg	23000
37194 [7440-44-0]	**Carbon Porous Nanopowder, 5000-8000nm**	1 Gms 10 Gms 50 Gms	1300 4200 18800
28186	**CNT Dispersant AC**	20 ml	12000
35210	**CNT Dispersant AQ**	20 ml	12000
59800 [308068-56-6]	**TCNT Type 1- Carbon Nanotubes Thin Walled, Length 5-30µm, OD 2-5nm, 90%**	25 Mg 100 Mg	12900 33200
73953 [308068-56-6]	**SWCNT Type 1 Carbon Nanotubes Single- Walled, Length 5-30µm, OD 1-2nm, 95%**	250 Mg 1 Gms	18000 60000
20770 [308068-56-6]	**SWCNT Type 2 Carbon Nanotubes Single- Walled (Short), Length 1-3µm, OD 1-2nm, 90%**	100 Mg 1 Gms	7500 70000
29129 [308068-56-6]	**SWCNT Type 3 Carbon Nanotubes Single- Walled COOH Functionalized, Length 0.5-2.0µm, OD 1-2nm, 90%**	25 Mg 100 Mg	5200 9400
18989 [308068-56-6]	**SWCNT Type 4 Carbon Nanotubes Single-Walled, Length 5-30µm, OD 1-2nm, 90%** (-COOH Functionalized)	50 Mg 250 Mg	6300 25000
55024 [308068-56-6]	**SWCNT Type 5 Carbon Nanotubes Single-Walled** (-OH Functionalized)	50 Mg 250 Mg	11200 30400
27404 [308068-56-6]	**SWCNT Type 6 Carbon Nanotubes Single-Walled, Length 5-30µm, OD 1-2nm, 99%** (-NH2 Functionalized)	50 Mg 250 Mg	16900 65200
95425 [308068-56-6]	**SWCNT Type 7 Carbon Nanotubes Single-Walled (Short), 1mg/100ml in Water, Length 1-3µm, OD 1-2nm, 90%** (Semiconducting)	5 ml 20 ml	12000 40000
52710 [308068-56-6]	**SWCNT Type 8 Carbon Nanotubes Single-Walled, 1wt% in Water, Length 5-30µm, OD 1-2nm, 90%** (Semiconducting)	5 ml 20 ml	10000 38000
36078 [308068-56-6]	**DWCNT Type 1 Carbon Nanotubes Double Walled, Length 50µm, OD 2-4nm, 60%**	100 Mg 500 Mg	22500 56200

Code	Name	Packing	Price
85143 [308068-56-6]	**DWCNT Type 2 Carbon Nanotubes Double Walled (Short), Length 0.5-2µm, OD 2-4nm, 60%**	100 Mg	23600
45565 [308068-56-6]	**DWCNT Type 3 Carbon Nanotubes Double-Walled** (-COOH Functionalized)	50 Mg	12400
69072 [308068-56-6]	**DWCNT Type 4 Carbon Nanotubes Double-Walled** (-OH Functionalized)	50 Mg 250 Mg	12400 38600
34528 [308068-56-6]	**HMWCNT Type1- Carbon Nanotubes Multi Walled Helical, Length 10µm, OD 100-200nm,(C) 90%**	250 Mg 1 Gms	12900 17400
57743 [308068-56-6]	**MWCNT Type 1 Carbon Nanotubes Multi Walled, Length 10-30µm, OD 8nm, 95%**	1 Gms 5 Gms 25 Gms	3800 12900 37200
88440 [308068-56-6]	**MWCNT Type 2 Carbon Nanotubes Multi Walled, Length 10-30µm, OD 8-15nm, 95%**	5 Gms 25 Gms	5200 13600
68465 [308068-56-6]	**MWCNT Type 3 Carbon Nanotubes Multi Walled, Length 10-30µm, OD 10-20nm, 95%**	1 Gms 5 Gms 25 Gms	2100 5200 20800
41067 [308068-56-6]	**MWCNT Type 4 Carbon Nanotubes Multi Walled, Length 10-30µm, OD 20-30nm, 95%**	5 Gms 25 Gms	3700 10400
41335 [308068-56-6]	**MWCNT Type 5 Carbon Nanotubes Multi Walled, Length 10-30µm, OD 30-50nm, 95%**	5 Gms 25 Gms 100 Gms	3200 6300 22900
96390 [308068-56-6]	**MWCNT Type 6 Carbon Nanotubes Multi Walled, Length 10-30µm, OD 50nm, 95%**	5 Gms 25 Gms 100 Gms	3400 6800 24000
35203 [308068-56-6]	**MWCNT Type 7 Carbon Nanotube Multi Walled (Graphitized), 99.9%**	1 Gms 5 Gms	8400 27800
73235 [308068-56-6]	**MWCNT Type 8 Carbon Nanotubes Multi-Walled, Length 20µm, OD 10-20nm, 99.9%** (Graphitized, -COOH Functionalized)	500 Mg 1 Gms	4000 6300
46918 [308068-56-6]	**MWCNT Type 9 Carbon Nanotubes Multi-Walled, Length 20µm, OD 10-20nm, 99.9%** (Graphitized -OH Functionalized)	500 Mg 1 Gms	4000 6300
58255 [308068-56-6]	**MWCNT Type 10 Carbon Nanotubes Multi-Walled, Length 10-30µm, OD 10-20nm, 95%** (-COOH Functionalized)	1 Gms 5 Gms	2500 5200
29466 [308068-56-6]	**MWCNT Type 11 Carbon Nanotubes Multi-Walled, Length 10-30µm, OD 20nm, 95%** (-NH2 Functionalized)	50 Mg 100 Mg 1 Gms 5 Gms	3400 5100 16900 50600
82903 [308068-56-6]	**MWCNT Type 12 Carbon Nanotubes Multi-Walled, Length 10-30µm, OD 10-20nm, 95%** (-OH Functionalized)	500 Mg 1 Gms 5 Gms	2900 5100 10200
28658 [308068-56-6]	**MWCNT Type 13 Carbon Nanotubes Multi-Walled, Length 10-30µm, OD 30-50nm, 95%** (-COOH Functionalized)	500 Mg 1 Gms 5 Gms 25 Gms	2500 4500 9000 29300
65875 [308068-56-6]	**MWCNT Type 14 Carbon Nanotubes Multi-Walled, Length 10-30µm, OD 30-50nm, 95%** (-NH2 Functionalized)	50 Mg 100 Mg	11300 16900
24968 [308068-56-6]	**MWCNT Type 15 Carbon Nanotubes Multi-Walled, Length~ 10-30µm, OD 30-35nm, 95%** (-OH Functionalized)	500 Mg 1 Gms 5 Gms	2500 3200 7300
33670 [308068-56-6]	**MWCNT Type 16 Carbon Nanotubes Multi-Walled Short, Length 0.5-2µm, OD 8nm, 95%**	500 Mg 1 Gms 5 Gms	4700 7300 28100
28174 [308068-56-6]	**MWCNT Type 17 Carbon Nanotubes Multi-Walled Short, Length 0.5-2µm, OD 10-20nm, 95%**	500 Mg 1 Gms 5 Gms	3200 5500 20800
90797 [308068-56-6]	**MWCNT Type 18 Carbon Nanotubes Multi-Walled Short, Length 0.5-2.0µm, OD 30-50nm, 95%**	1 Gms 5 Gms	5700 19100
42927 [308068-56-6]	**MWCNT Type 19 - Carbon Nanotubes Multi Walled (Conducting), Length 3-6micron,OD 10nm**	1 Gms 5 Gms	20300 67400
58908 [308068-56-6]	**MWCNT-Ni Type 1 Carbon Nanotubes Multi Walled, Nickel coated, Length 10-30 micron, OD 5-20nm, 98%**	1 Gms 5 Gms	9400 35400
53105 [308068-56-6]	**MWCNT-Ni Type 2 Carbon Nanotubes Multi Walled , Nickel coated**	1 Gms 5 Gms	11800 43800
90500 [308068-56-6]	**MWCNT-Ni Type 3 Carbon Nanotubes Multi Walled , Nickel coated**	1 Gms 5 Gms	11800 43800
75288 [308068-56-6]	**MWCNT-Ni Type 4 Carbon Nanotubes Multi Walled , Nickel coated**	1 Gms 5 Gms	11800 43800
65468 [9012-76-4]	**Chitosan Nanopowder, 80-100nm, 99%**	5 Gms 25 Gms	3600 10700
44316 [1308-38-9]	**Chromium (III) Oxide Nanopowder, 90nm, 99%** (Chromia)	10 Gms	14600
64835 ■ [1152310-69-4]	**C106 Dye, 90%** (Cis-Bis(Isothiocyanato) (2,2Bipyr-4,4-Dicarboxylato) (4,4bis(5(Hexylthio)thiophen2yl)2,2-Bipy)Ru(II))	25 Mg	9600
35347 ✔ [847665-45-6]	**K19 Dye (Sensitizer), 85%** (cis-Bis(Thiocyanato) (2,2-Bipyr-4,4-Dicarboxylic Acid)- (4,4-bis(Phexyloxystyryl) -2,2-Bipyr)Ru(II))	25 Mg	9700
74373 [1332-58-7]	**Clay Nanopowder, 80-150nm, 99%**	5 Gms 25 Gms	7100 10700
31570 [12052-28-7]	**Cobalt Iron Oxide Nanopowder, 50nm, 99.9%** (Cobalt Ferrite)	1 Gms 5 Gms	1300 4200
74613 [12052-28-7]	**Cobalt Zinc Iron Oxide Nanopowder, 50nm, 99.9%**	5 Gms 25 Gms	6300 16700
74343 [7440-50-8]	**Copper Nanopowder, 50nm, 99.9%**	5 Gms	30800

Code	Name	Packing	Price
73612 [7440-50-8]	**Copper Nanowires**	10 Mg	36400
28954 [1317-38-0]	**Cupric Oxide Nanopowder, 40nm, 99%** (Copper (II) Oxide)	5 Gms 25 Gms 100 Gms	2000 5700 18000
19151 ✔ [7782-40-3]	**Diamond Nanopowder, 4-7nm**	1 Gms	22600
85476 🌡 [1141011-53-1]	**Diisobutyl 2,2-Bipyridine-4,4-Dicarboxylate extrapure, 95%**	100 Mg	13400
91287 [1309-37-1]	**Ferric Oxide (Alpha) Nanopowder, 30nm, 99%** (Iron (III) Oxide (Alpha))	5 Gms 25 Gms 100 Gms 250 Gms	2000 5400 17200 24800
51744 [1309-37-1]	**Ferric Oxide (Gamma) Nanopowder, 20-50nm, 99.9%** (Iron (III) Oxide (Gamma))	5 Gms 25 Gms 100 Gms 250 Gms	1700 4300 12700 22600
16708 [1317-61-9]	**Ferrous Ferric Oxide Nanopowder (Type 1), 20-50nm, 99.9%** (Iron (II, III) Oxide)	5 Gms 25 Gms 100 Gms	1800 3300 11000
85993 [146958-06-7]	**Formamidinium Bromide extrapure, 99%** (Methanimidamide Bromide, Formamidine Hydrobromide, Iminomethylamine Hydrobromide)	1 Gms 5 Gms	4500 10500
67537 [879643-71-7]	**Formamidinium Iodide extrapure, 99%** (Formamidine Hydroiodide, Methanimidamide Hydroiodide)	1 Gms 5 Gms	2600 8400
37461	**Fullerene C60, 99.5%**	250 Mg 1 Gms	16000 40000
44170	**Fullerene C70, 99%**	100 Mg 500 Mg	12500 41600
84500 [7440-57-5]	**Gold Nanopowder, 90nm, 99.99%**	1 Gms	97400
73225 [7440-57-5]	**Gold Nanoparticles Dispersion (Spherical) (AU05), 0.05mg/ml Citrate, 0.1mM in PBS**	10 ml 25 ml	3200 6000
36848 [7440-57-5]	**Gold Nanoparticles Dispersion (Spherical) (AU10), 0.05mg/ml Citrate, 0.1mM in PBS**	10 ml 25 ml	3200 6000
72123 [7440-57-5]	**Gold Nanoparticles Dispersion (Spherical) (AU20), 0.05mg/ml Citrate, 0.1mM in PBS**	10 ml 25 ml	3200 6000
60974 [7440-57-5]	**Gold Nanoparticles Dispersion (Spherical) (AU40), 0.05mg/ml Citrate, 0.1mM in PBS**	10 ml 25 ml	3200 6000
91082 [7440-57-5]	**Gold Nanoparticles Dispersion (Spherical) (AU60), 0.05mg/ml Citrate, 0.1mM in PBS**	10 ml 25 ml	3200 6000
27944 [7782-42-5]	**Single Layer Graphene (SLG) Nanopowder, 98%**	10 Mg 100 Mg	6300 18800
73949 [7782-42-5]	**Single Layer Graphene Factory (SLGF) Nanopowder, 98%**	100 Mg 250 Mg	5500 13700
89922 [7782-42-5]	**Graphene Carboxyl (GCOOH) Nanopowder, 1000-5000nm, 99%**	25 Mg	12500
60703 [7782-42-5]	**Graphene Carboxyl (GCOOHW) Nanodispersion, 5mg/ml in water**	20 ml	45800
88398 [7782-42-5]	**Graphene Film (Super Paper) (GFSP), 20000nm, 98%**	1 piece	156000
87185 [7782-42-5]	**Graphene Nitrogen-doped (GNdp) Nanopowder, 500-5000nm**	50 Mg 250 Mg	6400 29200
90251 [7782-42-5]	**Graphene Oxide Film (Super Paper) (GOFSP),20000nm, 98%**	1 piece	116700
55093 [7782-42-5]	**Graphene Platelet Nanopowder (GPN Type 1), 6-8nm, 99.5%**	5 Gms 25 Gms	5800 20800
36529 [7782-42-5]	**Graphene Platelet Nanopowder (GPN Type 2), 11-15 nm, 99.5%**	1 Gms 5 Gms	1400 4400
98585 [7782-42-5]	**Graphene Platelet Nanopowder (GPN Type 3), 2-10nm, 99%**	5 Gms 25 Gms	3200 9400
53927 [7440-40-0]	**Graphene Quantum Dots (1mg/ml) (GQD), 15nm, 80%**	5 ml 25 ml	5500 20000
49888 [7782-42-5]	**Single Layer Graphene Oxide (SLGO) Nanoflakes, 0.8-1.2nm, 99%**	25 Mg 100 Mg 500 Mg	4700 11500 43700
58531 [7782-42-5]	**Single Layer Graphene Oxide (SLGOE) Nanodispersion, 5mg/ml in ethanol**	10 ml	21700
72536 [7782-42-5]	**Single Layer Graphene Oxide (SLGOW) Nanodispersion, 5mg/ml in water**	10 ml 25 ml	10800 18200
68908 [7782-42-5]	**Graphite Nanopowder (Type 2), 1-18nm, 87%**	1 Gms 5 Gms	10300 38600
85833 [7782-42-5]	**Graphite Nanopowder (Type 1), 400nm, 98%**	5 Gms 25 Gms	2600 9400
69869 [7782-42-5]	**Graphite Nanopowder Lubricant Grade, 50nm, 99.5%**	5 Gms 25 Gms	6800 30000
38058 [51311-17-2]	**Graphite Fluoride (Carbon Monofluoride) Nanopowder, 1000-10000nm**	1 Gms 5 Gms	2800 7800
39525 [7782-42-5]	**Graphite Oxide Nanopowder, Dia 500-1500nm**	100 Mg 500 Mg 1 Gms	2800 9900 17700
11934 [7782-42-5]	**Graphite Pyrolytic Micronpowder, 50micron**	5 Gms	8800
69895 [12167-74-7]	**Hydroxyapatite Micronpowder, 12micron, 99%** (Calcium Hydroxyapatite, Calcium Phosphate Tribasic)	5 Gms 25 Gms	6100 12700
13616 [12167-74-7]	**Hydroxyapatite Nanopowder, 20nm, 99%** (Calcium Hydroxyapatite, Calcium Phosphate Tribasic)	5 Gms 25 Gms	6300 15900
66614 [1312-43-2]	**Indium (III) Oxide Nanopowder ultrapure, 50nm, 99.99%**	1 Gms	34400
45974 [1312-81-8]	**Lanthanum Oxide Nanopowder ultrapure, 90nm, 99.99%**	5 Gms 25 Gms	3900 13700

Code	Name	Packing	Price
52726 ✔ [15365-14-7]	**Lithium Iron (II) Phosphate Micronpowder, 2micron, 98%** (Ferrous Lithium Phosphate, Iron Lithium Phosphate, LFP, Triphylite)	5 Gms 25 Gms	8100 16200
53768 [12031-95-7]	**Lithium Titanate (LTO) (Spinel) Nanopowder, 200nm, 99%** (Tetralithium Titanate)	25 Gms 100 Gms	14400 25200
69667 [12068-51-8]	**Magnesium Aluminate (Spinel) Nanopowder, 50nm, 99%**	1 Gms 2 Gms	4950 8630
71829 [1309-48-4]	**Magnesium Oxide Nanopowder, 50nm, 99.5%**	5 Gms 25 Gms	2600 7800
56176 [73905-81-4]	**Magnetic Iron Oxide Nanocrystals Powder, 20nm, 99%**	20 Mg	124500
85918 [1317-34-6]	**Manganese (III) Oxide Nanopowder, 30nm,99%** (Manganese Sesquioxide)	10 Gms	11500
39543 [6876-37-5]	**Methylammonium Bromide extrapure, 99%** (Methanamine Hydrobromide, Methylamine Hydrobromide)	1 Gms 5 Gms	6100 15600
25334 [14965-49-2]	**Methylammonium Iodide extrapure, 98%** (Methanamine Hydriodide)	1 Gms 5 Gms 25 Gms	2400 4300 11500
31562 [247058-06-6]	**Methyl 4,4``-Dimethyl-[2,2`:6`,2``-Terpyridine] -4``-Carboxylate extrapure, 95%**	50 Mg 100 Mg	12200 18900
83005 [1317-33-5]	**Molybdenum Disulfide Nanopowder Lubricant Grade, 90nm, 99%**	5 Gms 10 Gms	11700 22100
95822 [12316-56-2]	**MXene Nanoflakes, Multilayer, 100-200nm**	100 Mg 500 Mg	18800 62400
22701 [12363-89-2]	**MXene Nanoflakes, 1-10 layers** (few-layer MXene Ti3C2Tx nanoflakes)	25 Mg 100 Mg	22000 70000
16806 ✔ [141460-19-7]	**N3 Red Dye, 95%** (N-3 Red Dye, Bis(Isothiocyanato)Bis(2,2-Bipyridyl-4,4 -Dicarboxylato) Ruthenium(II))	100 Mg	8200
30786 [207347-46-4]	**N719 Dye, 95%** Di-Tetrabutylammonium cis-Bis(Isothiocyanato)Bis(2,2-Bipyridyl -4,4-Dicarboxylato) Ruthenium(II), N-7	50 Mg 250 Mg 1 Gms	6100 12100 35000
52716 🌡 [359415-47-7]	**N749 Dye (Black Dye), 95%** (Triisothiocyanato-(2,2:6,6-Terpyridyl-4,4,4-Tricarboxylato) Ruthenium(II) Di(Tetra-Butylammonium))	50 Mg	21700
35001 [7440-02-0]	**Nickel Nanopowder, 20nm, 99.9%**	5 Gms	30100
32293 [1313-99-1]	**Nickel (II) Oxide Nanopowder, 50nm, 99.5%** (Nickel Monoxide)	25 Gms	7300
73985	**Nickel Titanium Alloy Nanopowder, 30-120nm**	250 Mg 1 Gms	19300 57800
69691 [12645-50-0]	**Nickel Zinc Iron Oxide (NZFO) Nanopowder, APS 10-30nm, 98.5%** (Nickel Zinc Ferrite)	5 Gms 25 Gms	3200 13600
74380 🌡	**Perovskite Quantum Dots (pQD 450), Fluorescence 450nm (Hydrophobic)**	50 Mg	19700
15360 🌡	**Perovskite Quantum Dots (pQD 480), Fluorescence 480nm (Hydrophobic)**	50 Mg	19700
77345 🌡	**Perovskite Quantum Dots (pQD 510), Fluorescence 510nm (Hydrophobic)**	50 Mg	19700
23300 🌡	**Perovskite Quantum Dots (pQD 550), Fluorescence 550nm (Hydrophobic)**	50 Mg	19700
22116 [7440-05-3]	**Palladium Nanopowder (Type I), 30nm, 99.9%**	100 Mg 500 Mg	20000 75000
14387 [7440-06-4]	**Platinum Nanopowder, 50nm, 99.9%**	250 Mg	30000
65469 [9002-88-4]	**Polyethylene Nanospheres - Type 1, 740nm-4990nm**	25 Mg 100 Mg	6800 13000
50349 [9002-88-4]	**Polyethylene Nanospheres - Type 2, 200nm-9900nm**	25 Mg 100 Mg	7000 16300
64798 [9011-14-7]	**Poly(Methyl Methacrylate) (pMMAC) Nanospheres, 1500nm-11500nm**	100 Mg 500 Mg	4700 18900
28699 [7782-49-2]	**Selenium Nanopowder in aq. media, 0.1mg/ml, 50nm**	2 ml	123600
26112 [7440-21-3]	**Silicon Nanopowder, 40nm, 98%**	1 Gms 5 Gms	5500 25000
43362 [409-21-2]	**Silicon Carbide Nanopowder, 50nm, 98%**	25 Gms 100 Gms	5800 20300
69294 [7631-86-9]	**Silicon Dioxide Nanopowder, 15nm, 99.5%**	25 Gms 100 Gms	2400 7400
22806 [7631-86-9]	**Silicon Dioxide (Silica) Nanodispersion Type A (20nm), 40% in water**	25 Gms 100 Gms	2100 6200
38475 [10097-28-6]	**Silicon Monoxide Nanopowder, 400nm, 99.99%**	1 Gms 5 Gms	8900 35000
63338 [12033-89-5]	**Silicon Nitride Nanopowder, 20nm, 99%**	5 Gms 25 Gms	5200 12500
44416 [7440-22-4]	**Silver Nanopowder (Type I), 100nm, 99.9%**	1 Gms 5 Gms 10 Gms	4400 13600 25000
36214 [7440-22-4]	**Silver Nanopowder (Type II), 20nm, 99.9%**	1 Gms 5 Gms 10 Gms	3200 11500 18800
53941 [7440-22-4]	**Silver Nanowires**	10 Mg 50 Mg	6100 24000
46912 🌡 [7440-22-4]	**Silver Nanoparticles Dispersion (Spherical) (AG20), 0.02mg/ml in Citrate (5mM)**	10 ml 25 ml	3200 6000
49241 🌡 [7440-22-4]	**Silver Nanoparticles Dispersion (Spherical) (AG40), 0.02mg/ml in Citrate (2mM)**	10 ml 25 ml	3200 6000
58322 🌡 [7440-22-4]	**Silver Nanoparticles Dispersion (Spherical) (AG60), 0.02mg/ml in Citrate (2mM)**	10 ml 25 ml	3200 6000
15669 [18282-10-5]	**Stannic Oxide Nanopowder, 60nm, 99.5%** (Tin (IV) Oxide)	5 Gms 25 Gms	4400 12600

Code	Name	Packing	Price
27873 [7440-25-7]	**Tantalum Nanopowder, 70nm, 99.9%**	250 Mg 1 Gms	12500 33300
74708 ■	**Tectomer 2-Tailed (C8H16(CH2NHGLY4)2.2HCl), 90%**	5 Mg	13700
73330 ■	**Tectomer 3-Tailed (CH3C(CH2-NH-Gly7)3.3CF3CO2H) , 90%**	5 Mg	13700
27570 ■ [318286-59-8]	**Tectomer 4-Tailed (C(CH2-NH-Gly7)4.4HCl), 90%**	5 Mg	13700
72951 [7440-32-6]	**Titanium Nanopowder, 50nm, 99.9%**	1 Gms	10000
74629 [13463-67-7]	**Titanium Dioxide Photocatalyst Nanopowder, 7nm, 95%** (Titanium (IV) Oxide)	25 Gms 100 Gms	4300 11000
90885 [13463-67-7]	**Titanium Dioxide Catalyst Nanopowder, 50nm, 98%** (Titanium (IV) Oxide)	25 Gms 100 Gms	3600 10700
94632 [13463-67-7]	**Titanium Dioxide AnataseNanodispersion (15nm), 15% in water** (Titanium (IV) Oxide)	25 Gms	9900
40566 [13463-67-7]	**Titanium Dioxide Ultrapure Nanopowder, 250nm, 99.97%** (Titanium (IV) Oxide)	25 Gms 100 Gms	4200 11700
35299 [13463-67-7]	**Titanium Dioxide Rutile Nanopowder, 20nm, 92%** (Titanium (IV) Oxide)	25 Gms	4100
10161 [13463-67-7]	**Titanium Dioxide Rutile Nanopowder, 50nm, 99%**	25 Gms	4200
86729 [13463-67-7]	**Titanium Dioxide Rutile Nanopowder (Hydrophilic), 50nm, 99%** (Titanium (IV) Oxide)	25 Gms 100 Gms	4100 12000
64394 [13463-67-7]	**Titanium Dioxide Rutile Nanodispersion (30nm), 40% in water** (Titanium (IV) Oxide)	25 Gms	8200
16547 [13463-67-7]	**Titanium Dioxide Nanowires (Type 1)** (Titanium (IV) Oxide)	50 Mg	13000
22091 [13463-67-7]	**Titanium Dioxide Nanowires (Type 2)** (Titanium (IV) Oxide)	50 Mg	13000
52163 [25583-20-4]	**Titanium Nitride Nanopowder, 20nm, 97%** (Titanium (III) Nitride)	5 Gms 10 Gms	6000 11500
75513 ✔ [1392221-69-0]	**FK 102 Co(II) PF6 Salt Complex, 98%** (Tris(2-(1H-Pyrazol-1-yl)Pyridine)Cobalt(II) Di[hexafluorophosphate])	1 Gms 5 Gms	3900 13200
51632 ✔ [1346416-71-4]	**FK 102 Co(III) PF6 Salt Complex, 98%** (Tris(2-(1H-Pyrazol-1-yl)pyridine)cobalt(III) Tri[Hexafluorophosphate])	1 Gms 5 Gms	5000 16400
66142 ✔	**FK 102 Co(II) TFSI Salt Complex, 98%** (Tris(2-(1H-Pyrazol-1-yl)pyridine)cobalt(II) Di[bis(Trifluoromethane)sulfonimide])	1 Gms 5 Gms	3900 13200
84076 ✔ [2057441-94-6]	**FK 102 Co(III) TFSI Salt Complex, 98%** (Tris(2-(1H-pyrazol-1-yl)pyridine)cobalt(III) Tri[bis(Trifluoromethane)sulfonimide)	1 Gms 5 Gms	2800 9300
33139 🌡 [1197995-33-7]	**Triethyl [2,2`:6`,2``-Terpyridine] -4,4`,4``-Tricarboxylate extrapure, 95%**	50 Mg	11100
88917 [7440-33-7]	**Tungsten Nanopowder, 50nm, 99.9%**	1 Gms 5 Gms	4700 14600
51456 [12070-12-1]	**Tungsten Carbide Nanopowder, 200nm, 99%**	5 Gms 25 Gms 100 Gms	3500 9500 35000
23882 [12138-09-9]	**Tungsten Disulfide (WS2) Aerosol Spray (13Oz/369g)**	1 piece	16000
53032 [12138-09-9]	**Tungsten Disulfide Micronpowder, 0.6micron, 99.9%** (Tungsten (IV) Sulphide)	25 Gms 100 Gms	3700 12000
78519 [12138-09-9]	**Tungsten Disulfide Nanopowder, 90nm, 99%**	1 Gms 5 Gms	3300 14000
23851 [1314-35-8]	**Tungsten Trioxide Nanopowder, 60-90nm, 99.5%** (Tungsten (VI) Oxide, Tungstic Anhydride)	5 Gms 25 Gms	3000 12000
82874 ✔ [502693-09-6]	**Z907 Dye, 95%** cis-Bis(Isothiocyanato) (2,2-Bipyridyl-4,4-Dicarboxylato) (4,4-Di-Nonyl-2-Bipyridyl) Ruthenium(II), Mo	50 Mg 250 Mg	6100 17000
19366	**Zeolite - Mesoporous Carbon Nanopowder (CMK-3 Type),Dia. 3-6nm**	100 Mg 250 Mg	13700 20800
31755	**Zeolite - Mesoporous Silica Molecular Sieve KIT-5**	1 Gms	13200
65276	**Zeolite - Mesoporous Silica Molecular Sieve Kit-6**	1 Gms	13200
83881 [7631-86-9]	**Zeolite - Mesoporous Silica Nanopowder (SBA-15 Type), Dia. 6-11nm**	250 Mg 1 Gms	3300 9900
97621 [7631-86-9]	**Zeolite - Mesoporous Silica Nanopowder (1D-Hexagonal SBA-41 Type)**	250 Mg 500 Mg	15200 22300
63876 [7631-86-9]	**Zeolite - Mesoporous Silica Nanopowder (3D-Cubic MCM-48 Type)**	250 Mg 500 Mg	15200 22300
79132	**Zeolite - Molecular Sieve SAPO-11**	5 Gms	10200
52201	**Zeolite - Ultrastable Y**	5 Gms	10800
74601 [7440-66-6]	**Zinc Nanopowder, 50nm, 99.9%**	5 Gms	5800
50841 [12063-19-3]	**Zinc Iron Oxide Nanopowder, 50nm, 98.5%** (Zinc Ferrite)	10 Gms	8900
91148 [1314-13-2]	**Zinc Oxide Nanopowder (Type I), 30nm, 99.9%**	10 Gms 25 Gms 100 Gms 500 Gms	2300 4600 6900 28400

Code	Name	Packing	Price	Code	Name	Packing	Price
60345 [1314-13-2]	**Zinc Oxide Nanopowder (Type II), 240nm, 99.999%**	5 Gms 25 Gms 100 Gms 500 Gms	2300 4600 6900 28400				
64168 [1314-13-2]	**Zinc Oxide Nanodispersion Type A-Nonionic (70nm), 50% in water**	10 Gms	10600				
80467 [1314-13-2]	**Zinc Oxide Nanodispersion Type B-Anionic (70nm), 50% in water**	10 Gms	10200				
74963 [1314-13-2]	**Zinc Oxide Nanodispersion Type C-Cationic (70nm), 50% in water**	10 Gms	10200				
99997 [12036-43-0]	**Zinc Titanate Nanopowder, 90 nm, 99.5%**	1 Gms	17500				
38501 [1314-23-4]	**Zirconium Dioxide Nanopowder, 45nm, 99.5%** (Zirconia, Zirconium (IV) Oxide)	25 Gms 100 Gms	3700 9400				

Part D

Dehydrated Culture Media Products

- Product Listing
- Media Supplements
- ChroMed™
- Media Kits

Harmonized Culture Media

In recent years, much of the focus has started shifing to newer methods of microbiological testing, led by the leading international pharmaceutical bodies across the world, namely European Pharmacopoeia, Japanese Pharmacopoeia and US Pharmacopoeia. These new methods include harmonization in microbiological testing and this will greatly simplify drug manufacturing norms.

SRL Harmonized Culture Media are manufactured, formulated, labelled and performance tested in compliance with the new harmonized standards.

SRL is ready to face the challenges of the future, be in Biotech, Pharmaceutical R & D, Diagnostics, or any application or research in Microbiology. We offer this new range of products to cater to the new demands that we will be fulfilling through this range.

Important highlights of our Harmonized Media:

- The usage of only the Finest Quality of Bacteriological ingredients
- Strict quality control of Raw Material testing
- Highly trained microbiologists ensure that our Harmonized media are comparable to the best international products
- Very sensitive and precise instrumentation and media manufacturing equipments make sure that each batch is consistent in quality
- Constant R & D and customer interaction helps us to design appropriate products

With our dedicated team of highly skilled microbiologists, an ISO 9001-2008 quality management system and our CE certification.





Conformité Européenne

All **SRL** culture media products are CE certified and now meet EU consumer safety, health and environmental requirements.

All the manufacturing, testing parameters and packaging are done according to the CE norms.

Code	Name	Packing	Price
38515	✔ **A1 Broth (A-1 Medium)** For the detection of faecal coliforms in foods, treated wastewater and seawater by Most Probable Number (MPN) Method.	100 Gms 500 Gms	1800 7620
26408	✔ **AATCC Bacteriostasis Agar (FDA Agar)** (FDA Agar, ATCC Medium 182) APHA, FDA recommended medium for detecting plate counts of microorganisms in food, water and waste water.	100 Gms 500 Gms	1800 5160
49025	✔ **AATCC Bacteriostasis Broth (FDA Broth)** For testing antibacterial activities of antiseptics and disinfectants.	100 Gms 500 Gms	1460 5040
59235	✔ **Acetamide Agar** For the differentiation of non-fermentative gram negative bacteria especially Pseudomonas aeruginosa.	100 Gms 500 Gms	1460 5270
20605	✔ **Acetamide Broth** For the confirmation of Pseudomonas aeruginosa in bottle water.	100 Gms 500 Gms	1350 4930
64773	✔ **Acetate Differential Agar** (Sodium Acetate Agar) For the differentiation of Shigella species from those of the Escherichia genus.	100 Gms 500 Gms	1570 5040
31257	✔ **Acetobacter Agar (Glucose)** Used for the cultivation and maintenance of glucose positive Acetobacter species.	100 Gms 500 Gms	1570 5380
85514	✔ **Acid Broth** For the isolation of acid tolerant bacteria from canned foods.	100 Gms 500 Gms	1460 4710
76557	✔ **Actidione Agar Base w/o Actidione** Used for the enumeration and detection of bacteria in specimens containing large numbers of yeasts and moulds.	100 Gms 500 Gms	1910 6050
15126	🌡 **Actinomyces Broth** Used for the maintenance and cultivation of Actinomyces species.	100 Gms 500 Gms	1460 4710
24729	✔ **Actinomycete Isolation Agar** Used for the isolation and cultivation of Actinomyces species from soil and water.	100 Gms 500 Gms	1570 5160
31047	✔ **Adams Agar** Used for examining sporulation in yeasts.	100 Gms 500 Gms	1350 4930

Code	Name	Packing	Price
59640	✔ **A. K. Agar No. 2 (Sporulating Agar, Arret & Kirshbaum Medium)** Used for production of spores of Bacillus subtilis ATCC 6633.	100 Gms 500 Gms	1680 5600
55418	✔ **Aeromonas Isolation Medium Base** For the selective differentiation and isolation of Aeromonas hydrophila and other Aeromonas species from clinical specimens and foods.	100 Gms 500 Gms	2020 6720
63013	🌡 **Aeromonas Selective Supplement (2.5mg/vl)** For the selective isolation and cultivation of Aeromonas species.	5 vl	1570
34609	✔ **Aero Pseudo Selective Agar (GSP Agar)** Used for detection of Aeromonas and Pseudomonas in food, as well as in equipments and water supplies of food industry.	100 Gms 500 Gms	1570 6160
77981 [9002-18-0]	**Agar powder BactoBio for microbiology** Specially manufactured to provide clear, high gel strength agar, suitable for use in preparation of microbiological media.	100 Gms 500 Gms 2.5 Kg	2360 8850 42120
35799	✔ **Aleksandrow Agar** Widely used media for isolation and identification of Potassium solubilizing bacteria from soil samples.	100 Gms 500 Gms	1570 5720
72439	✔ **Algae Culture Broth** Used for the isolation and cultivation of algae from soil, water and sewage.	100 Gms 500 Gms	1240 5040
36191	✔ **Algae Culture Agar** Used for isolation and cultivation of algae from soil, water and sewage.	100 Gms 500 Gms	1350 5160
45712	✔ **Alicyclobacillus Detection Agar (Alicyclobacillus Agar)** Used for detection of Alicyclobacillus in fruit juices.	100 Gms 500 Gms	2800 9750
80501	✔ **Alkaline Peptone Water** For the cultivation of Vibrio species from sea foods and clinical specimens.	100 Gms 500 Gms	1010 3480
84893	✔ **Alkaline Peptone Water (ISO Basis)** Used for the enrichment of Vibrio species.	100 Gms 500 Gms	1010 3480
83863	✔ **Alkaline Peptone Water (BIS Basis)** Used for the enrichment of Vibrio species.	100 Gms 500 Gms	1010 3480

Code	Name	Packing	Price
58784	✔ **Alkaline Saline Peptone Water (ISO Basis)** Used for the enrichment of Vibrio species from food and water samples.	100 Gms 500 Gms	1240 4260
44860	✔ **AC Agar (All Culture Agar)** For growth of aerobic, anaerobic and microaerophilic organisms. Recommended for sterility testing of solutions and materials not containing mercurial preservatives.	100 Gms 500 Gms	1240 5040
87091	✔ **AC Broth (All Culture Broth)** For cultivation and isolation of aerobic, anaerobic and microaerophilic organisms. Recommended for sterility testing of solutions and materials not containing mercurial preservatives.	100 Gms 500 Gms	1120 4600
52023	✔ **Alternative Thioglycollate Medium (IP, USP Basis)** (NIH Thioglycollate Broth) For the sterility testing of biological products which are turbid or viscous as mentioned in pharmacopeia.	100 Gms 500 Gms	680 2020
39422	✔ **Alternative Thioglycollate Medium, Sterile** Gamma irradiated sterile powder used for sterility testing of biological products which are turbid or viscous.	100 Gms 500 Gms	790 2920
98349	✔ **Amies Transport Medium w/ Charcoal** For the transport of swab specimens to prolong the survival of microorganisms, especially Neisseria gonorrhoeae.	100 Gms 500 Gms	1460 4370
24193	✔ **Amies Transport Medium w/o Charcoal** For the transport of swab specimens to prolong the survival of microorganisms.	100 Gms 500 Gms	1460 4370
90600	✔ **Ammonium Phosphate Agar** For the detection of microorganisms that can utilize ammonium phosphate as nitrogen source.	100 Gms 500 Gms	1570 6950
69084	✔ **Ampicillin Dextrin Agar Base** Recommended for the selective isolation, differentiation and quantitative enumeration of Aeromonas species from water, food and environmental samples by the membrane filter method.	100 Gms 500 Gms	1570 6280
50747	✔ **Ampicillin Dextrin Broth Base** Recommended for the selective isolation, differentiation and quantitative enumeration of Aeromonas species from water, food and environmental samples by the membrane filter method.	100 Gms 500 Gms	1350 5490
24891	🌡 **Ampicillin Dextrin Selective Supplement** Used for selective isolation and differentiation of Aeromonas species.	5 vl 5 x5 vl	1350 4480
78162	🌡 **Ampicillin Supplement (10mg/vl)** For the selective isolotaion of Aeromonas species.	5 vl	1800
12213	✔ **Anaerobic Agar** For the cultivation of anaerobic bacteria, especially Clostridium species and other anaerobic organisms fr	100 Gms 500 Gms	1460 5720
20911	✔ **Anaerobic Agar (Brewer)** For the isolation and sensitivity testing of facultative and obligatory anaerobes.	100 Gms 500 Gms	1680 5720
97394	✔ **Anaerobic Blood Agar Base** Recommended for isolation and cultivation of Group A and Group B Streptococci from throat cultures and other clinical samples.	100 Gms 500 Gms	1680 5830
59821	✔ **Anaerobic Agar w/o Dextrose** For cultivation of variety of anaerobic bacilli, especially Clostridium species.	100 Gms 500 Gms	1570 5380
19422	✔ **Anaerobic Egg Yolk Agar Base** For the detection of Clostridium species, especially from food samples.	100 Gms 500 Gms	1680 5940
85198	✔ **Anaerobic Tryptone Soya Agar** Used for screening of anaerobes in cosmetics (Talcum powder).	100 Gms 500 Gms	1570 5160
47702	✔ **Andrade Lactose Peptone Water** For the study of fermentation reactions of members of Enterobacteriaceae.	100 Gms 500 Gms	1010 3810
41449	✔ **Andrade Peptone Water** A basal medium to study fermentation reactions by adding various carbohydrates	100 Gms 500 Gms	1120 5040
90295	✔ **Andrade Peptone Water BioVeg** A basal medium to study fermentation reactions by adding various carbohydrates.	100 Gms 500 Gms	1350 5380

Code	Name	Packing	Price
70197	✔ **Antibiotic Assay Medium No. 5 (Streptomycin Assay Agar w/ Yeast Extract)** For antibiotic assay testing. For the streptomycin assay using Bacillus subtilis as the test organism.	100 Gms 500 Gms	1460 6160
24301	✔ **Antibiotic Assay Medium No. 8 (Base Agar w/Low pH)** For microbiological assay testing of Tetracycline,Vancomycin,Oxytetracycline,etc.	100 Gms 500 Gms	1460 6160
52812	✔ **Antibiotic Assay Medium No. 9 (Polymyxin Base Agar)** Used as base agar for the plate assay of Carbenicillin, Colistimethate and Polymixin B.	100 Gms 500 Gms	1350 5040
26751	✔ **Antibiotic Assay Medium No. 10 (Polymyxin Seed Agar)** Used as seed agar for plate assay of Carbenicillin, Colistimethate and Polymixin B.	100 Gms 500 Gms	1680 5040
87878	✔ **Antibiotic Assay Medium No. 1 (Seed Agar)** (Penassay Seed Agar, Seed Agar, Agar Medium A) For antibiotic assay testing, detection of antibiotics in milk and determination of the antimicrobial effectiveness of antibiotics.This medium is in accordance to specifications mentioned in the United States Pharmacopoeia -National Formulatory, 2022.	100 Gms 500 Gms	1350 5940
50741	✔ **Antibiotic Assay Medium No. 2 (Base Agar, Penassay Base Agar))** (Penassay Base Agar, Base Agar) For use as a base layer in antibiotic assay testing.This medium is in accordance to formulation mentioned in USP.	100 Gms 500 Gms	1460 6160
66729	✔ **Antibiotic Medium No. 3 (Assay Broth, Penassay Broth) (USP, IP Basis)** For microbiological assay of antibiotics.	100 Gms 500 Gms	1240 3920
10892	✔ **Antifungal Assay Agar** For assaying antifungal activity.	100 Gms 500 Gms	1570 5160
91104	■ **Ampicillin (AMP) AST Disc, 10mcg**	1 Vial (100 Discs)	280
14976	■ **Ampicillin (AMP) AST Disc, 25mcg**	1 Vial (100 Discs)	280
86985	■ **Azithromycin (AZM) AST Disc, 15mcg**	1 Vial (100 Discs)	280
73758	■ **Cefazolin (CZ) AST Disc, 30mcg**	1 Vial (100 Discs)	280
55518	■ **Cefixime (CFM) AST Disc, 5mcg**	1 Vial (100 Discs)	280
53629	■ **Cefoperazone (CPZ) AST Disc, 30mcg**	1 Vial (100 Discs)	280
49236	■ **Cefotaxime (CTX) AST Disc, 30mcg**	1 Vial (100 Discs)	280
39776	■ **Ceftazidime (CAZ) AST Disc, 30mcg**	1 Vial (100 Discs)	280
63746	■ **Chloramphenicol (CMP) AST Disc, 30mcg**	1 Vial (100 Discs)	280
45449	■ **Ciprofloxacin (CIP) AST Disc, 5mcg**	1 Vial (100 Discs)	280
30238	■ **Colistin (CL) AST Disc, 10mcg**	1 Vial (100 Discs)	280
65158	■ **Erythromycin (E) AST Disc, 15mcg**	1 Vial (100 Discs)	280
40955	■ **Gentamicin (GEN) AST Disc, 10mcg**	1 Vial (100 Discs)	280
24572	■ **Gentamicin (GEN) AST Disc, 120mcg**	1 Vial (100 Discs)	280
93264	■ **Nalidixic Acid (NA) AST Disc, 30mcg**	1 Vial (100 Discs)	280
77759	■ **Netilmicin (NET) AST Disc, 30mcg**	1 Vial (100 Discs)	280
70150	■ **Polymyxin B (PB) AST Disc, 50Units**	1 Vial (100 Discs)	280
47075	■ **Tetracycline (TE) AST Disc, 30mcg**	1 Vial (100 Discs)	280
32890	■ **Vancomycin (VA) AST Disc, 30mcg**	1 Vial (100 Discs)	280
12797	✔ **APRY Agar Base** For the cultivation of acid resistant yeasts.	100 Gms 500 Gms	1350 5040
26126	✔ **APT Agar (APHA)** Used for the cultivation of heterofermentative lactobacilli requiring high thiamine content.	100 Gms 500 Gms	1460 5270
31990	✔ **APT Broth (APHA)** Used for the cultivation of heterofermentative lactobacilli requiring high thiamine content.	100 Gms 500 Gms	1460 5040
22599	✔ **Aqua Test E.Coli Medium** Aqua Test Ecoli Medium, is used for detection of faecal indicator bacteria in drinking water, especially in low-income countries and outside controlled laboratory settings.	100 Gms	9300
41524	✔ **Arabinose Agar Base** Used for selective isolation of Enterococcus faecium from faeces, sewage and water samples.	100 Gms 500 Gms	1570 6050
85908	✔ **Arginine Dihydrolase Broth** Used for the detection of arginine dihydrolase producing microorganisms.	100 Gms 500 Gms	1350 5380

Code	Name	Packing	Price
15577	✔ **Aseptic Commissioning Medium (Aseptic Packing Line Medium)** It is an ideal medium used for validation of aseptic packing and filling lines.	100 Gms 500 Gms	1460 5160
64778	✔ **Ascospore Agar** Ascospore Agar is a solid medium recommended for use in qualitative procedures for the cultivation of ascomycetous yeasts.	500 Gms	5160
88277	🌡 **Ashby's Glucose Agar** For the cultivation of Azotobacter species that can utilize glucose as carbon and atmospheric nitrogen as nitrogen source.	100 Gms 500 Gms	1570 5600
26569	✔ **Ashby's Mannitol Agar** For the cultivation of Azotobacter species that can utilize mannitol as carbon and atmospheric nitrogen as nitrogen source.	100 Gms 500 Gms	1570 5600
75830	✔ **ASLA Agar** Used for isolation and cultivation of Propionibacterium species from foods.	100 Gms 500 Gms	1460 5270
54470	🌡 **ASLA Enrichment Supplement** Recommended growth supplement for cultivation of Propionibacteria	5 vl	1070
57624	✔ **Aspergillus Differentiation Medium Base** Recommended for detection of aflatoxin producing Aspergillus species from food samples.	100 Gms 500 Gms	1800 6720
53903	✔ **Asparagine Nitrate Medium** Recommended media for isolation and cultivation of denitrifying bacteria.	100 Gms	3920
69318	✔ **Asparagine Proline Broth** For the cultivation of Pseudomonas aeruginosa by membrane filter technique.	100 Gms 500 Gms	1460 5270
95670	✔ **Asparagine Proline Broth BioVeg** For cultivation of Pseudomonas aeruginosa by membrane filter technique.	100 Gms 500 Gms	1570 5830
92301	✔ **Ayers and Johnson Agar (Stock Culture Agar)** Used for maintaining stock cultures of bacteria, particularly streptococci.	100 Gms 500 Gms	1570 5600
94692	✔ **Azide Blood Agar Base** Used for isolation of Streptococcus sps.i and Staphylococcus sps., with added blood may be used for determination of hemolytic reactions.	100 Gms 500 Gms	1570 5380
47479	✔ **Bacillus Differentiation Agar** Used for the differentiation between Bacillus cereus and Bacillus subtilis based on mannitol fermentation.	100 Gms 500 Gms	2020 6280
51747	✔ **Bacillus Medium** Used for the cultivation of Bacillus licheniformis.	100 Gms 500 Gms	1120 4260
74664	✔ **Bacillus Cereus Agar** Used for the isolation and enumeration of Bacillus cereus.	100 Gms 500 Gms	1350 5380
98947	✔ **Bacteroides Bile Esculin (BBE) Agar Base** Used for the selective isolation, identification and cultivation of Bacteroides fragilis group.	100 Gms 500 Gms	3810 15240
83141	🌡 **Bacteroides Selective Supplement (50mg/vl)** Used for the selective isolation of Bacteroides species.	5 vl	1460
50222	✔ **Baird-Parker Agar Base (USP, IP, ISO Basis)** Used as a selective and diagnostic medium for the isolation and enumeration of coagulase-positive staphylococci in foods.When used with added supplements;egg yolk emulsion and potassium tellurite,this medium meets the formulation and performance criteria mentioned in ISO 6888 parts 1 and 3 and ISO 11133 :2014.	100 Gms 500 Gms	1240 5040
24562	✔ **Baird-Parker Agar Base (BIS Basis)** Recommended medium for isolation and enumeration of coagulase positive staphylococci from food and other specimens.This medium is in accordance to specifications given by BIS under IS : 5887 (Part II)1976.	100 Gms 500 Gms	1240 5040
21748	✔ **Baird-Parker Agar Base w/Sulpha** Recommended medium for isolation and enumeration of coagulase-positive staphylococci from food and other specimens.	100 Gms 500 Gms	1800 6280

Code	Name	Packing	Price
26267	✔ **B.C. Motility Medium** Used for enumeration of motility of Bacillus cereus strains.	100 Gms 500 Gms	1350 4480
97676	✔ **B.C.P. - D.C.L.S. Agar** (Bromocresol Purple-Deoxycholic Acid-Citrate-Lactose-Sucrose Agar) Used for the selective isolation of Salmonella and Shigella species.	100 Gms 500 Gms	1570 5720
69574	✔ **Beef Extract Agar (ex. Buffalo)** For the cultivation and maintenance of a wide variety of microorganisms.	100 Gms 500 Gms	1240 4930
39359	✔ **Beef Extract Broth (ex. Buffalo)** For the cultivation of a wide variety of microorganisms.	100 Gms 500 Gms	1120 4260
60856 [68990-09-0]	✔ **Beef Extract ex. Buffalo BactoBio for bacteriology** A dehydrated extract of bovine tissue, used generally to compliment the nutritive properties of peptone. It is source of nitrogen, vitamins, amino acids and carbon.	500 Gms 5 Kg	2360 20610
57542	✔ **Beijerinckia Medium** Used for the isolation of Beijerinckia sps.	100 Gms 500 Gms	1570 6050
29141	✔ **BIGGY Agar (Nickerson Medium)** (Nickerson Agar) For the detection and isolation of Candida species.	100 Gms 500 Gms	1350 5600
49319	✔ **Bile Broth Base** For the culture of blood clots from patients with suspected enteric fever.	100 Gms 500 Gms	1570 6160
66013	✔ **Bile Esculin Agar w/o azide (ISO Basis)** For the isolation and presumptive identification of group D Streptococci.	100 Gms 500 Gms	4150 16800
66572	✔ **Bile Esculin Agar Base** For the isolation and presumptive identification of group D Streptococci.	100 Gms 500 Gms	3360 13440
34963	✔ **Bile Esculin Azide Agar** For the isolation and presumptive identification of group D Streptococci.	100 Gms 500 Gms	2920 8960
79630	✔ **B.D.G.Broth, Hajna** For the selective enrichment and cultivation of enteric bacilli from treated drinking water.	100 Gms 500 Gms	1460 4710
64104	✔ **Bile Salts Agar** For the isolation and cultivation of bile tolerant enteric bacilli.	100 Gms 500 Gms	2360 9190
50362	✔ **Bile Salts BactoBio for bacteriology** Used in microbiological culture media as a selectively inhibitory agent.	100 Gms 500 Gms	3300 15000
18276	✔ **Bile Salts No. 3 BactoBio for bacteriology** Bile salts No. 3 is a refined bile salt for use as a selective inhibitory agent in bacteriological culture media. It consist of a specially modified fraction of bile acid salts (comprising mainly of sodium cholate and sodium deoxycholate).	100 Gms 500 Gms	9640 44690
69577	✔ **Bile Salts Brilliant Green Starch Agar** Used for isolation and identification of Aeromonas hydrophila from food and environmental samples.	100 Gms 500 Gms	1350 4480
56175	✔ **Bile Salts Mixture BactoBio for bacteriology** Used as an inhibitory agent for gram positive organisms in microbiological media.	100 Gms 500 Gms	9190 26550
18037	✔ **Bile Salts Technical for bacteriology** It is used in bacteriological culture media as an inhibitory agent for gram positive bacteria, making it selective for isolation and identification of pathogens.	100 Gms 500 Gms	1680 6390
79665	✔ **Bile Tolerant Gram negative bacteria Test Kit (Pre-Weighed Media)** Used for isolation, identification and cultivation of Bile tolerant gram negative bacteria.	1 Kit	1460
92875	✔ **Biochemical Test Kit - I (Pre-Weighed Media)** Used for biochemical identification of gram negative enteric bacteria based on tests like carbohydrate utilization, H2S production, phenylalanine deaminase and urease production.	1 Kit	3140
65896	✔ **Biochemical Test Kit - II (Pre-Weighed Media)** Used for biochemical identification of microorganisms based on tests - Nitrate reduction,Malonate utilisation,Gelatin liquefaction and Starch hydrolysis.	1 Kit	2920

Part D – Dehydrated Culture Media

Code	Name	Packing	Price
55785	✔ **Bio Peptone (Mixture of Casein & Meat Peptone) BactoBio for bacteriology** Biopeptone is a mixture of enzymic digest of casein and meat peptone. It is recommended for the cultivation of fastidious microorganisms in Columbia Broth Base, Dextrose Agar and Casman Agar Base.	500 Gms 5 Kg 10 Kg	3480 27220 51970
55369	✔ **Bismuth Sulphite Agar (ISO Basis)** (Bismuth Sulfite Agar (I)) For the selective isolation and identification of Salmonella typhi and other enteric bacilli.	100 Gms 500 Gms	1120 4370
48794	✔ **Blood Agar Base** (Infusion Agar) For the preparation of blood agar by the addition of sheep blood.	100 Gms 500 Gms	1120 4370
90935	✔ **Blood Agar Base (w/o sheep blood)** For the isolation, cultivation and detection of haemolytic activity of streptococci and other fastidious microorganisms.	100 Gms 500 Gms 5 Kg	1120 4370 40660
67644	✔ **Blood Agar Base No.2 (w/ 1.2% Agar)** For the isolation, cultivation and detection of haemolytic activity of streptococci, pneumococci and other particularly fastidious microorganisms.	100 Gms 500 Gms	1240 4370
74999	✔ **Blood Agar Base w/ Low pH (w/o sheep blood)** For the isolation and growth of a wide variety of microorganisms. The slightly acid pH of this medium enhances distinct hemolytic reactions.	100 Gms 500 Gms	1240 4480
81083	✔ **Blood Free Campylobacter Selectivity Agar Base** Used for the selective isolation and differentiation of Campylobacter species from food and animal feeding stuffs. This medium is recommended by ISO 10272-1&2:2017.	100 Gms 500 Gms	1800 6390
71476	✔ **BG-11 Broth (Blue-Green Medium)** For culture and maintenance of cyanobacteria.	100 Gms 500 Gms	1800 5940
85308	✔ **Bordet - Gengou Agar Base (w/o Rabbit Blood)** For the detection and isolation of Bordetella species from clinical specimens.	100 Gms 500 Gms	1570 5940
15097	🌡 **B.P. Sulpha Supplement (50mg/vl)** Used for suppressing the growth of Proteus species on Baird Parker Agar Base.	5 vl	1460
91661	✔ **BPL Agar (Brilliant Green Phenol Red Agar)** For isolation and identification of Salmonellae with the exception of S.typhi in faeces, urine, milk,meat and other materials.	100 Gms 500 Gms	1570 5380
46530	✔ **Brain Heart CC Agar** (BHI CC Agar) Recommended for isolation and cultivation of fastidious pathogenic fungi and saprophytic fungi from specimens heavily contaminated with bacteria.	100 Gms 500 Gms	2690 11200
64944	✔ **Brain Heart Infusion Powder (BHI Powder) BactoBio for bacteriology** Used for cultivation of highly fastidious microorganisms like Streptococci, Pneumococci, Meningococci, Haemophilus, Neisseria sp. and pathogenic fungi. A highly nutritious mixture of calf brain infusion and beef heart infusion.	500 Gms 5 Kg	3590 33270
94124	✔ **Brain Heart Infusion Agar (BHI Agar)** For the cultivation of wide variety of fastidious bacteria, yeasts and molds from clinical and non-clinical specimens.	100 Gms 500 Gms 5 Kg	800 2400 21000
87864	✔ **Brain Heart Infusion Broth (BHI Broth)** For the cultivation of variety of fastidious and non-fastidious aerobic and anaerobic microorganisms.	100 Gms 500 Gms 5 Kg	750 2200 19000
82856	✔ **Brain Heart Infusion Agar, Modified (BHI Agar, Modified)** For the cultivation of wide variety of fastidious bacteria, yeast and molds from clinical and non-clinical specimens.	100 Gms 500 Gms 5 Kg	900 2700 23000
84560	✔ **Brain Heart Infusion w/ 0.1% Agar (BHI w/ 0.1% Agar)** For the cultivation of fastidious pathogenic microorganisms.	100 Gms 500 Gms	900 2700
26671	✔ **Brain Heart Infusion w/ PABA and Agar (BHI Agar w/PABA)** Used for the examination of blood samples from patients who have recieved sulfonamide therapy,Improved growth of anaerobes is seen on addition of agar.	100 Gms 500 Gms	1000 4000

Code	Name	Packing	Price
25357	✔ **Brain Heart Infusion Broth w/6.5% NaCl (BHI Broth w/6.5% NaCl)** Nutritive medium used for cultivation of salt tolerant micro-organisms.	100 Gms 500 Gms	1240 4480
33355	✔ **Brewer Thioglycollate Medium** For sterility testing of biological products and other materials.	100 Gms 500 Gms	900 3360
65774	✔ **Brilliant Green Agar** For the selective isolation of Salmonella other than S.typhi from faeces, food and dairy products.	100 Gms 500 Gms	1120 4150
22374	✔ **Brilliant Green Agar w/ 1.2% Agar** Recommended as an enrichment medium for isolation of Salmonellae, from faeces, urine and pathological samples.	100 Gms 500 Gms	1240 5160
14980	✔ **Brilliant Green Agar, Modified (Brilliant Green Agar Base w/ Phosphates) (ISO Basis)** For the selective isolation of Salmonellae from faeces, water, food and dairy products.	100 Gms 500 Gms	1120 4150
18545	✔ **Lactose Bile Brilliant Green Broth 2% (ISO, BIS Basis)** Recommended medium for detection and confirmation of coliforms from water, wastewater, foods, milk and dairy products.This medium conforms to specifications given by ISO 4831:2006 , ISO 4832 : 2006,BIS committee under the specifications IS:5401-1:2012.	100 Gms 500 Gms	1350 4260
98383	✔ **Brilliant Green Bile Agar** For the detection and enumeration of coliform bacteria in water, sewage and foods.	100 Gms 500 Gms	1240 5040
39284	✔ **Brilliant Green Sulfa Agar (B.G. Sulpha Agar)** For the selective isolation of Salmonella species, other than Salmonella typhi from food, dairy products, eggs, egg products and feed.	100 Gms 500 Gms	1240 5380
32445	✔ **Bolton Broth Base** A medium for the selective pre-enrichment of Campylobacter species in food samples.	100 Gms 500 Gms	1350 5830
55582	🌡 **Bolton Selective Supplement** Used for the selective enrichment of Campylobacter species from food samples.	5 vl 5 x5 vl	2130 9190
66006	🌡 **Bordetella Selective Supplement (20mg/vl)** Used for the selective isolation of Bordetella species.	5 vl	1680
31732	✔ **Bromocresol Purple Agar** A non-selective medium, used for the differentiation of Enterobacteriaceae. Recommended medium for bacteriological examination of containers and recipients used in water bottling.	500 Gms	4500
49977	✔ **Bromocresol Purple Azide Broth** For the confirmation tests for presence of faecal streptococci in water and waste water.	100 Gms 500 Gms	1350 4930
19202	✔ **Bromocresol Purple Lactose Broth (BCP Lactose Broth, Yeast Fermentation Broth Base)** Used as presumptive medium for the detection of coliform bacteria in water.	100 Gms 500 Gms	1240 4480
53213	✔ **Bromocresol Purple Lactose Broth w/o Sodium Chloride**	100 Gms 500 Gms	1240 4260
26384	🌡 **Bromocresol Purple Supplement (15mg/vl)** Used for isolation and identification of faecal streptococci.	5 vl	1120
46172	🌡 **Bromothymol Blue Supplement (20mg/vl)** Used for differentiating microflora from urine sample.	5 vl	790
68035	✔ **Bryant and Burkey Agar** For detection and enumeration spores of lactate fermenting Clostridium species in milk and dairy products.	100 Gms 500 Gms	1460 5380
61882	✔ **Bryant and Burkey Medium** For the detection and enumeration of lactate fermenting Clostridial species in milk and dairy products.	100 Gms 500 Gms	1350 5040
19814	✔ **B.T.B. Lactose Agar** (Bromothymol Blue Lactose Agar) For the isolation and cultivation of pathogenic Staphylococci.	100 Gms 500 Gms	1460 5270
28540	🌡 **B12 Culture Agar (USP Basis)** For the cultivation and maintenance of Lactobacillus leichmannii ATCC 7830 to be used as the test organism in the vitamin B12 assay.	100 Gms	2130

Code	Name	Packing	Price
79328	**B12 Inoculum Broth (USP Basis)** For the preparation of inoculum cultures of Lactobacillus leichmannii ATCC 7830, which is used as the test organism in the vitamin B12 assay.	100 Gms 500 Gms	2020 8520
21700	**Buffered Glycerol Saline Base** With addition of glycerol is used in the collection and transportation of faecal specimens.	100 Gms	1120
72993	**Buffered Peptone Water** A pre-enrichment medium for the isolation of Salmonella species, especially injured microorganisms from food sources.	100 Gms 500 Gms	900 3480
88680	**Buffered Peptone Water (ISO, BIS Basis)** Recommended medium for pre-enrichment of injured Salmonella species from foods.This medium is in accordance to specifications given by ISO 6579-1:2017, ISO 21528-1:2004 and ISO 22964-4:2017 and BIS ,IS :5887 -1999.	100 Gms 500 Gms	1010 3480
12995	**Buffered Peptone Water w/ NaCl** An isotonic agent recommended for use as a diluent in performing microbial examination of non-sterile products.	100 Gms 500 Gms	1240 3480
22197	**Buffered Peptone Water BioVeg** A pre-enrichment medium, manufactured from animal free bases for the isolation of Salmonella, especially injured microorganisms from various food sources.	100 Gms 500 Gms	1010 3480
93439	**Buffered Peptone Water w/ Casein** A pre-enrichment medium, for the isolation of Salmonella, especially from cocoa (chocolates and confectionary) etc. prior to selective enrichment and isolation.	100 Gms 500 Gms	1120 4150
85268	**Buffered Tryptone Glucose Yeast Extract Broth** Used for the isolation of Clostridium perfringens from food samples.	100 Gms 500 Gms	1350 4600
70054	**Buffered Yeast Agar** For the cultivation of yeast and molds and for controlling bottle washing operations in the soft drinks and related industries.	100 Gms 500 Gms	1240 4600
89637	**Buffered Sodium Chloride-Peptone Solution pH 7.0,Harmonized** Used as diluent in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1010 3480
73461	**Burk`s Medium** Used for isolation and cultivation of nitrogen fixing bacteria such as Azotobacter species from soil.	100 Gms 500 Gms	1570 4930
32639	**Bushnell - Haas Agar** For examining fuels for microbial contamination and for studying hydrocarbon utilization by microorganisms.	100 Gms 500 Gms	1240 5040
14509	**Bushnell - Haas Broth** For examining fuels for microbial contamination and for studying hydrocarbon utilization by microorganisms.	100 Gms 500 Gms	1010 4260
16574	**CAE Agar Base (Citrate Azide Enterococcus)** Used for identification of Enterococci in meat, dairy and other food products.	100 Gms 500 Gms	1350 5040
35512	**CAL Agar (Cellobiose Arginine Lysine Agar)** Used for the isolation and characterization of Yersinia enterocolitica from faecal specimens.	100 Gms	8740
95845	**Calcium Caseinate Agar** For the detection and enumeration of proteolytic microorganism in foodstuff and other materials.	100 Gms 500 Gms	1800 5830
25867	**Campylobacter Agar Base** Used with addition of blood and antibiotics in isolation and cultivation of Campylobacter species.	100 Gms 500 Gms	960 3590
16092	**Campylobacter Supplement-I (Blaser-Wang)** An antibiotic supplement used for the selective isolation of Campylobacter species.	5 vl	1570
32048	**Campylobacter Supplement-II (Skirrow)** An antibiotic supplement used for the selective isolation of Campylobacter species.	5 vl	1570
36198	**Candida BCG Agar Base** For use with added Neomycin for isolation and differentiation of Candida species from primary specimens.	100 Gms 500 Gms	1570 5720
95160	**Candida Isolation Agar (Candida Agar)** For the isolation and differentiation of Candida albicans.	100 Gms 500 Gms	1570 5720

Code	Name	Packing	Price
45568	✔ **Candida Albicans Test Kit (Pre-Weighed Media)** Used for isolation, identification and cultivation of Candida albicans.	1 Kit	1460
92558	✔ **Carbohydrate Fermentation Test Kit - I (Pre-Weighed Media)** For the differentiation of microorganisms on the basis of various carbohydrate fermentation reactions.	1 Kit	4480
40071	✔ **Carbohydrate Fermentation Test Kit - II (Pre-Weighed Media)** For the differentiation of microorganisms on the basis of various carbohydrate fermentation reactions.	1 Kit	8630
54951	✔ **Carbon Utilisation Agar** Used for characterization of Streptomyces on the basis of carbon utilization studies.	500 Gms	2800
24715	✔ **Cary-Blair Transport Medium Base (w/o Charcoal)** For maintenance as a holding medium or transport medium of clinical specimens during collection or shipment.	100 Gms 500 Gms	1350 5270
68806 [65072-00-6]	✔ **Acicase (Casein Acid Hydrolysate, Casamino Acids) BactoBio for bacteriology** Acid digest of casein, a rich source of nitrogen recommended for use in microbiological culture media that require a completely hydrolyzed protein as a nitrogen source.	500 Gms 2 Kg	7730 25880
95305 [65072-00-6]	✔ **Acicase (Casein Acid Hydrolysate,Casamino Acids) Technical for bacteriology** Acid digest of casien; is a rich source of nitrogen.	500 Gms 2 Kg	5940 21170
68211	✔ **Casein Magnesium Broth (NZM Broth)** Used in the cultivation of recombinant strains of Escherichia coli.	100 Gms 500 Gms	1680 6500
61796	✔ **Casein Peptone BactoBio for bacteriology** An enzymatic digest of casein, rich in amino acids and peptides. Used as a nitrogen source in microbiological media.	500 Gms 5 Kg	2920 26550
72689	✔ **Casitone (Casitose) BactoBio for bacteriology** A hydrolyzed casein used as an ingredient in microbiological culture media requiring free amino acids, especially for toxin and vaccine production.	500 Gms	2920
49946	🌡 **Casitose Agar w/ 1.5% Agar** Used as a general purpose medium for cultivation of microorganisms.	500 Gms 5 Kg	4930 44580
79033	✔ **Casitose Soya Agar, Modified (Casein Soya Agar, Modified)** For the growth of wide variety of organisms.	100 Gms 500 Gms	1570 5040
13927	✔ **Casitose Soya Blood Agar Base (Casein Soya Blood Agar Base)** (Trypticase Soya Blood Agar Base) Recommended for the cultivation of fastidious microorganisms and detection of their hemolytic activity.	100 Gms 500 Gms	1570 5380
91254	✔ **Casman Agar Base** Used with addition of blood, for isolation and cultivation of fastidious microrganisms from clinical specimes.	100 Gms 500 Gms	1350 5040
18871	🌡 **CBI Supplement** Used for selective isolation of Clostridium botulinum.	5 vl	1800
90465	✔ **Cetrimide Agar Base, Harmonized** Used for the selective isolation and identification of Pseudomonas aeruginosa in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1120 4260
45470	✔ **Cetrimide Agar Base** For the selective isolation, cultivation and identification of Pseudomonas aeruginosa.	100 Gms 500 Gms	1120 4260
93927	✔ **Cetrimide Agar BioVeg** Made from Veg.based ingredients ; it is used for the selective cultivation and identification of Pseudomonas aeruginosa.	100 Gms 500 Gms	1120 4260
31199	✔ **Cetrimide Broth (USP Basis)** Used for the cultivation of Pseudomonas aeruginosa.	100 Gms 500 Gms	1120 3920
94532	🌡 **Cetrinix Supplement** Used for the selective isolation of Pseudomonas species.	5 vl 5 x5 vl	2240 5380
70984	🌡 **CFC Supplement** Used for selective isolation of Pseudomonas species.	5 vl 5 x5 vl	2020 5940
11507	✔ **Chapman Stone Agar** For the isolation of staphylococci from a variety of specimens.	100 Gms 500 Gms	1350 4600

Part D – Dehydrated Culture Media

Code	Name	Packing	Price
58323	✔ **Charcoal Agar** For the cultivation and maintenance of fastidious microorganisms, especially Bordetella pertussis, for production of vaccines.	100 Gms 500 Gms	1460 5040
53641	✔ **Charcoal Agar w/ Niacin** Used for cultivation and isolation of Bordetella pertussis and Haemophilus influenzae.	100 Gms 500 Gms	1680 5830
52057	✔ **China Blue Lactose Agar** Used for differentiating lactose fermenters and non fermenters and also in enumeration of bacteria in milk.	100 Gms 500 Gms	1350 5380
52665	🌡 **Chloramphenicol Yeast Glucose Agar BioVeg** Used for the isolation and enumeration of molds in milk and milk products.	100 Gms 500 Gms	1010 4260
27880	🌡 **ChroMed Aureus Agar** A chromogenic medium used for the identification of Staphylococcus aureus.	100 Gms 500 Gms	2690 11650
78355	🌡 **ChroMed Bacillus Cereus Agar** A chromogenic medium used for the isolation of Bacillus species from a mixed culture in foods, clinical and non-clinical samples.	100 Gms 500 Gms	1680 6500
98478	🌡 **ChroMed Candida Agar Base** Chromogenic medium for the isolation and identification of major clinical significant Candida species.	100 Gms 500 Gms	5830 23520
42436	🌡 **ChroMed Chromogenic Coliform Agar (CCA) (ISO Basis)** For detection of Escherichia coli and coliforms in water samples as per ISO standards 9308-1:2014.	50 Gms 100 Gms 500 Gms	3700 6840 32150
19440	🌡 **ChroMed Chromogenic Coliform Agar w/ SLS (CCA w/ SLS)** For detection of Escherichia coli and coliforms in water samples.	50 Gms 100 Gms 500 Gms	3030 5380 20610
96599	🌡 **ChroMed ECC Agar** Differential medium recommended for presumptive identification of Escherichia coli and other coliforms in food and environmental samples.	100 Gms 500 Gms	5040 19040
75955	🌡 **ChroMed E.coli Agar** Recommended chromogenic medium, used for the detection and enumeration of Escherichia coli and coliforms from food samples.	100 Gms 500 Gms	3480 15240
46706	🌡 **ChroMed Klebsiella Selective Agar Base** Used for the selective isolation and detection of Klebsiella species from water and other sources. This medium can also be used in membrane filtration procedure.	100 Gms 500 Gms	2580 10200
45615	🌡 **ChroMed Listeria Agar (ISO Basis)** A selective and differential chromogenic medium recommended for rapid and direct identification of listeria species.	100 Gms 500 Gms	2130 9080
38284	🌡 **ChroMed Listeria Agar Base, Modified** For the selective identification and differentiation of Listeria monocytogenes.	100 Gms 500 Gms	2130 8740
19589	🌡 **ChroMed Listeria Ottaviani -Agosti Agar Base (ISO basis)** Recommended for the selective and differential isolation of Listeria monocytogenes from food , animal feeds and environmental samples in food manufacturing and handling facilities.The formulation and performance criteria are in accordance to specification given by ISO 11290- parts 1 & 2:2017 and ISO 11133:2014.	100 Gms 500 Gms	2470 11760
65621	🌡 **ChroMed MeReSA Agar** A chromogenic media recommended for the isolation and selective identification of Methicillin Resistant Staphylococcus aureus (MRSA) from clinical isolates.	100 Gms 500 Gms	3030 9410
75213	🌡 **Chromed M-Lauryl Sulfate Agar (MLGA)** (Membrane Lactose Glucuronide Agar) For the differentiation and enumeration of Escherichia coli and other coliforms by membrane filter technique.	100 Gms 500 Gms	6160 24640
90219	🌡 **ChroMed OGYE Agar Base (ISO Basis)** Selective and differential medium for isolation of yeasts and moulds from milk and milk products.	100 Gms 500 Gms	4480 19600
92671	🌡 **ChroMed Salmonella Agar** A chromogenic medium used for the differentiation of Salmonella species.	100 Gms 500 Gms	12210 57010

Code	Name	Packing	Price
52496	**ChroMed TBX Agar (E.coli Agar Modified)** A selective chromogenic medium used for the detection and enumeration of Escherichia coli in food samples,animal feeds and water.	100 Gms 500 Gms	2130 7730
42979	**ChroMed UTI Agar** (Chromogenic Orientation Agar) A chromogenic medium used for the presumptive identification and differentiation of microorganisms causing urinary tract infections.	100 Gms 500 Gms	3590 16920
42308	**ChroMed VRE Agar Base** Recommended medium for detection of Vancomycin Resistant Enterococci (VRE) from clinical specimens.	100 Gms 500 Gms	2580 11760
96238	**Chromogenic Listeria Lipase-C Supplement (ISO basis)** (LP Enrichment Supplement) Selective supplement for the isolation, enumeration and detection of Listeria monocytogenes.	5 vl	9190
80559	**CHO Medium Base** (Carbohydrate Medium Base) Used as a basal medium to which carbohydrates are added for fermentation studies of anaerobic bacteria.	100 Gms 500 Gms	1240 4480
28850	**Chocolate Agar Base** For the isolation and identification of a variety of fastidious microorganisms.	100 Gms 500 Gms	1460 5040
45553	**Cholera Medium Base** Used for selective isolation of Vibrio species from specimens grossly contaminated with Enterobacteriaceae.	100 Gms 500 Gms	1460 5490
93008	**Cholera Medium Base BioVeg** Veg. based (non-animal based) medium, used for selective isolation of Vibrio species from specimens grossly contaminated with Enterobacteriaceae.	100 Gms 500 Gms	1570 5720
14602	**Chloramphenicol Selective Supplement (50mg/vl)** Used for the selective isolation and cultivation of yeasts and moulds.	5 vl 5 x5 vl	1460 4150
60822	**Chloramphenicol Yeast Glucose Agar (ISO Basis)** For selective enumeration of fungi, yeasts and molds in milk and milk products.This medium is in accordance to specifications given by ISO 7218 :2013 (Amd.),ISO/DIS 7904 and ISO 5403 : 1999.	100 Gms 500 Gms	1010 4370
93306	**Christensen Citrate Agar** For the differentiation of enteric pathogens based on their ability to utilize citrate as a carbon source.	100 Gms 500 Gms	1240 5270
55769	**CLED Agar** (Brolacin Agar) For the isolation, enumeration and presumptive identification of microorganisms from urine.	100 Gms 500 Gms	1240 4370
88366	**CLED Agar w/ Andrade Indicator** For differentiating various microorganisms depending on their colony characteristics.	100 Gms 500 Gms	1240 4370
40980	**Clostridium Botulinum Isolation Agar Base** For the isolation, cultivation and differentiation of Clostridium botulinum from food samples.	100 Gms 500 Gms	1350 5270
26836	**Clostridium Difficile Agar Base** Recommended for selective isolation of Clostridium difficile from food and certain pathological specimens	100 Gms 500 Gms	1680 5830
93834	**Clostridium Difficile Supplement** Used for the selective isolation of Clostridium difficile.	5 vl 5 x5 vl	1570 4600
12655	**Clostridium Test Kit (Pre-Weighed Media)** Used for isolation, identification and cultivation of Clostridium species.	1 Kit	1460
64078	**CNA Supplement** Recommended for the selective isolation of anaerobic gram-positive cocci.	5 vl	2130
19156	**Coagulase Mannitol Agar Base** For the cultivation, differentiation of Staphylococcus aureus from other Staphylococcus species, based on coagulase production and mannitol fermentation.	100 Gms 500 Gms	1350 5600
52655	**Coliform Medium (CM, Coliform Broth)** For the isolation & cultivation of coliform organisms from raw milk, cream, yoghurt.	100 Gms 500 Gms	1240 4260
21493	**Columbia Agar Base, Harmonized** Used for the isolation of fastidious and non-fastidious microorganisms in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1120 4370

Part D – Dehydrated Culture Media

Code	Name	Packing	Price
45481	✔ **Columbia Agar Base** For the isolation and cultivation of fastidious and non-fastidious microorganisms.	100 Gms 500 Gms	1120 4370
97789	✔ **Columbia C.N.A. Agar Base** For selective isolation and cultivation of pathogenic Gram-positive cocci from clinical and non-clinical specimens.	100 Gms 500 Gms	3030 11990
23588	✔ **Conn`s Agar** Used for the cultivation of fungi.	100 Gms 500 Gms	2690 10760
16340	✔ **Corn Meal Agar w/ Dextrose** Cultivation of phytopathological and other fungi.	100 Gms 500 Gms	1240 5160
26998	✔ **Corn Meal Agar** For the chlamydospore production by Candida albicans and maintenance medium for fungal stock cultures.	100 Gms 500 Gms	1240 5160
60277	✔ **Corn Meal Peptone Yeast Agar**	100 Gms 500 Gms	1460 5040
73398	✔ **Cooke Rose Bengal Agar Base** Used for the selective cultivation and isolation of fungi.	100 Gms 500 Gms	2020 6950
45455	✔ **CPC Agar Base** For cultivation and differentiation of Vibrio species from food samples.	100 Gms	7510
65165	🌡 **CPC Supplement** For isolation and enumeration of Vibrio species from foods	5 vl	2800
48743	✔ **Craigs Medium**	100 Gms 500 Gms	1570 4600
70952	✔ **CRAMP Agar Base** For the cultivation of Yersinia species with plasmids.	100 Gms	6280
41775	✔ **Cystine Tellurite Agar Base** Recommended for selective isolation and differentiation of Corynebacterium diphtheriae types.	100 Gms 500 Gms	1570 5270
27793	✔ **Cystine Tryptone Agar** (Cystine Tryptic Agar) Used for the maintenance, subculturing and detection of motility of various bacteria.	100 Gms 500 Gms	1240 4370
77853	✔ **Crystal Violet Agar** (Crystal Violet Lactose Agar) For the differentiation of pathogenic and non- pathogenic staphylococci.	100 Gms 500 Gms	1460 5160
57485	✔ **Crystal Violet Lactose Broth** For the detection of coliforms in water filtration control works.	100 Gms 500 Gms	1680 4480
91022	✔ **Crystal Violet Tetrazolium Agar Base** Used for the detection of gram-negative psychrotrophic bacteria causing food spoilage.	100 Gms 500 Gms	1350 5380
29212	✔ **Czapek - Dox Agar** For the cultivation of fungi and bacteria capable of utilizing sodium nitrate as the sole source of nitrogen.	100 Gms 500 Gms	1120 3920
93812	✔ **Czapek - Dox Agar Modified** For the cultivation and maintenance of numerous fungal species, also for chlamydospore production by Candida albicans.	100 Gms 500 Gms	1460 5040
74367	✔ **Czapek - Dox Broth** For the cultivation of fungi and bacteria capable of utilizing sodium nitrate as the sole source of nitrogen.	100 Gms 500 Gms	1010 3480
95417	✔ **Czapek Malt Agar** For the isolation and cultivation of saprophytic fungi.	100 Gms 500 Gms	1120 4370
61195	🌡 **Dairy samples Testing Kit (Pre-Weighed Media)** Used for identification and enumeration of microorganisms from dairy products.	1 Kit	2580
96640	✔ **Decarboxylase Test Medium Base (Falkow)** (Decarboxylase Basal Medium, Decarboxylase Media, Falkow) For the differentiation of bacteria based on their ability to decarboxylate the amnio acid.	100 Gms 500 Gms	1240 4370
47129	✔ **Deoxycholate Agar** Differential medium for direct counts of coliforms from dairy products and for the isolation of enteric microorganisms from clinical and non-clinical specimens.	100 Gms 500 Gms	900 4040

Code	Name	Packing	Price
20023	✔ **Deoxycholate Citrate Agar (IP Basis)** Used for the selective isolation and cultivation of enteric pathogens, especially Salmonella and Shigella species.The composition and specification of this medium is in accordance with Indian Pharmacopoeia .	100 Gms 500 Gms	1680 8520
38695	✔ **Deoxycholate Lactose Agar** For the isolation and differentiation of Gram-negative enteric bacilli and for enumeration of coliforms from water, waste water, milk and dairy products.	100 Gms 500 Gms	1460 5160
18752	✔ **Dermatophyte Test Agar Base (D.T.M. Agar Base)** For the isolation and cultivation of dermatophytic fungi.	100 Gms 500 Gms	1460 5380
68650	✔ **DEV Nutrient Agar** To determine total germ count in water acc. To the "Deutschen Einheitsverfahren" regulations of drinking water.	100 Gms 500 Gms	1350 4600
22438	✔ **DEV Tryptophan Broth** Used for sub-cultivation for coliform differentiation and for indole testing in the bacteriological examination of water as per German Standard methods.	100 Gms 500 Gms	1010 3480
15923	✔ **Dextrose Agar** For the cultivation of a wide variety of microorganisms and used as a base for blood agar and for general laboratory purposes.	100 Gms 500 Gms	1120 4480
64098	✔ **Dextrose Broth** For the isolation and enrichment of fastidious microorganisms and for the study of glucose (dextrose) fermentation.	100 Gms 500 Gms	1010 3250
83607	✔ **Dextrose Mannitol Agar (Gillies Agar No.1)** Recommended medium for primary isolation of Salmonella and Shigella species, by detection of urease production, dextrose and mannitol fermentation.	100 Gms 500 Gms	1570 5600
98009	✔ **Dextrose Proteose Peptone Agar Base** Used in combination with blood and tellurite for the isolation of Corynebacterium diphtheriae.	100 Gms 500 Gms	1350 4600
11581	✔ **Dextrose Salt Agar** For the enumeration of yeasts and molds in butter and dairy products.	100 Gms 500 Gms	1350 5160
87461	✔ **Dextrose Salt Broth** For the enumeration of yeasts and molds in butter and dairy products.	100 Gms 500 Gms	1240 4260
44486	✔ **Dextrose Starch Agar** For the cultivation and maintenance of Neisseria gonorrhoeae and other fastidious microorganisms.	100 Gms 500 Gms	1460 5040
50223	✔ **Dextrose Tryptone Agar** (Dextrose Casein Peptone Agar) For the isolation and cultivation of mesophilic and thermophilic aerobic microorganisms related to food spoilage.	100 Gms 500 Gms	1350 4930
78209	✔ **Dey-Engley Neutralizing Broth** Recommended for the testing and neutralization of antiseptics and disinfectants.	100 Gms 500 Gms	1460 5270
73318	✔ **Diamalt Agar**	100 Gms 500 Gms	1680 6160
29412	✔ **Dichloran Glycerol Medium Base** Used for isolation and enumeration of viable osmophilic yeasts and xerophilic molds in food ,animal feed products and clinical samples.	100 Gms 500 Gms	1120 5380
63110	✔ **Dichloran Medium Base w/ Rose Bengal** Recommended for selective isolation and enumeration of yeasts and moulds associated with food spoilage.	100 Gms 500 Gms	1680 7060
66473	✔ **Differential Reinforced Clostridial Agar** For the cultivation and enumeration of sulfite - reducing clostridia from water samples.	100 Gms 500 Gms	1460 5040
27548	✔ **Differential Reinforced Clostridial Broth Base** For the cultivation of Clostridia from water samples.	100 Gms 500 Gms	1010 3360
52978	✔ **Differential Reinforced Clostridial Broth Base (ISO Basis)** For the cultivation of Clostridia from water samples.This medium is in accordance to specifications given by ISO 6461-1:1986.	100 Gms 500 Gms	1120 3920
88167	✔ **Differential Reinforced Clostridial Broth Base BioVeg** For the cultivation of Clostridia from water samples.	100 Gms 500 Gms	1120 4150
63899	✔ **Dihydrolase Broth Base** Used for studying dihydrolase reaction of Vibrio parahaemolyticus.	100 Gms 500 Gms	1570 4930

Part D – Dehydrated Culture Media

Code	Name	Packing	Price
14912	✔ **Diluting Fluid A** Used for sterility testing of pharmaceuticals.	100 Gms 500 Gms	1010 3590
23121	✔ **Diluting Fluid D** Used for sterility testing of pharmaceutical products in accordance with USP.	100 Gms 500 Gms	1120 4040
53277	✔ **Diluting Fluid K** Used as diluting or rinsing fluid, when performing. sterility testing of pharmaceuticals.	100 Gms 500 Gms	1240 3590
82965	✔ **Disinfectant Test Broth, AOAC** For the determination of phenol coefficient of disinfectants.	100 Gms 500 Gms	1240 4480
82359	✔ **DNase Test Agar Base** Used for differentiation of microorganisms based on their production of deoxyribonuclease.	100 Gms 500 Gms	4370 18260
78317	✔ **DNase Test Agar w/Methyl Green** Used for differentiation of microorganisms based on deoxyribonuclease activity, especially for identification of pathogenic Staphylococci from clinical samples.	100 Gms 500 Gms	5040 19040
88941	✔ **Drakes Medium 10** Recommended for isolation of Pseudomonas species from water samples.	100 Gms 500 Gms	1800 6050
53292	✔ **DRBC Agar (Dichloran Rose Bengal Chloramphenicol Agar) Base** Used for the selective isolation and enumeration of yeasts and moulds of significance in food spoilage.	100 Gms 500 Gms	1600 5700
84127	✔ **Drigalski Lactose Agar**	100 Gms 500 Gms	2020 5940
56107	✔ **Drigalski Litmus Lactose Agar**	100 Gms 500 Gms	1570 5270
76861	🌡 **Doyle`s Antibiotic Supplement** Used for enrichment of Campylobacter species.	5 vl 5 x5 vl	1350 4710
26609	✔ **Doyle`s Enrichment Broth Base** Used for enrichment of Campylobacter species	100 Gms 500 Gms	1460 5490
80907	✔ **Doyle`s Enrichment Broth Base BioVeg** Used for enrichment of Campylobacter species	100 Gms 500 Gms	1680 5830
49763	✔ **E. Coli Test Kit (Pre-Weighed Media)** Used for isolation, identification and cultivation of Escherichia coli.	1 Kit	1680
25549	✔ **EC Medium (EC Broth)** For the differentiation and enumeration of coliforms in water, wastewater, shellfish and foods.	100 Gms 500 Gms	1120 4260
66389	✔ **EC Medium w/MUG** Used for the fluorogenic detection of Escherichia coli in water, food and milk samples.	100 Gms 500 Gms	1680 6610
54723	✔ **Elliker Broth (Lactobacilli Broth)** Used for cultivation of streptococci and lactobacilli, especially in dairy procedures.	100 Gms 500 Gms	1240 4480
11983	✔ **Eijkman Lactose Medium (Eijkman Lactose Broth)** For the cultivation and differentiation of Escherichia coli from other coliform organisms based on their ability to ferment lactose and produce gas.	100 Gms 500 Gms	1350 4260
95514	✔ **EMB Agar** (Eosin Methylene Blue Agar) For the isolation, cultivation and differentiation of gram negative enteric bacteria based on lactose fermentation.	100 Gms 500 Gms	1010 4150
46675	✔ **EMB Agar, Levine** (Eosin Methylene Blue Agar, Levine) For the isolation, cultivation and differentiation of Gram-negative enteric bacteria based on lactose fermentation.	100 Gms 500 Gms	1010 4150
51207	✔ **EMB Broth** (Eosin Methylene Blue Broth) For the differentiation of gram negative enteric bacteria from clinical and non clinical specimens.	100 Gms 500 Gms	1010 3810
65512	✔ **Emerson Agar** Used for isolation and cultivation of Actinomycetaceae, Streptomycetaceae, fungi and moulds.	100 Gms 500 Gms	1570 5160
85046	✔ **Endo Agar** For the selective isolation, cultivation and differentiation of coliforms and other enteric microorganisms based on their ability to ferment lactose.	100 Gms 500 Gms	1120 4040

Code	Name	Packing	Price
81496	✔ **Endo DEV Agar** Used for selection and differentiation of E.coli acc. To the "Deutschen Einheitsverfahren" regulations for drinking water.	100 Gms 500 Gms	1350 4600
95904	✔ **Enriched Thioglycollate Broth** (Thioglycollate Medium, Enriched, Thioglycollate Medium, Base) For the isolation, cultivation and identification of a wide variety of obligate anaerobic bacteria.	100 Gms 500 Gms	1010 4370
33945	✔ **Enteric Fermentation Media** For the cultivation and differentiation of a variety of bacteria based on their ability to ferment different carbohydrates.	100 Gms 500 Gms	1240 3480
71611	✔ **Enterobacteria Enrichment Broth Mossel (EE Broth Mossel), Harmonized** Used for the selective enrichment of Enterobacteriaceae in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1240 4820
58444	✔ **Enterococcus Confirmatory Agar** For the identification of enterococci from water supplies and other sources by the confirmatory test.	100 Gms 500 Gms	1240 4710
64819	✔ **Enterococcus Confirmatory Broth** For the identification of enterococci from water by the confirmatory test.	100 Gms 500 Gms	1120 4370
47082	🌡 **Enterococcus Selective Supplement** Used for the selective isolation and cultivation of Enterococci.	5 vl 5 x5 vl	2020 6500
60312	✔ **Esculin Iron Agar** For cultivation and identification of Enterococci based on their ability to hydrolyze esculin.	100 Gms 500 Gms	4930 21060
27131	🌡 **Esculin Supplement (0.5gm/vl)** Used for detection of group D streptococci on the basis of esculin hydrolysis.	5 vl	2240
69404	✔ **Ethyl Violet Azide Broth (E.V.A. Broth)** For isolation, cultivation and enumeration of Enterococci from water and other specimens.	100 Gms 500 Gms	1350 4150
48357	✔ **Eugonic Agar (BIS Basis)** Used cultivation and maintenance of a variety of fastidious microorganisms.	100 Gms 500 Gms	1460 5380
66070	✔ **Eugonic Broth** Used cultivation and maintenance of a variety of fastidious microorganisms.	100 Gms 500 Gms	1460 4480
78507	🌡 **FGTC Antibiotic Supplement (2.5mg/vl)** Used for the selective isolation and cultivation of Enterococci.	5 vl	1460
93466	✔ **Fluid Sabouraud Medium** For the cultivation of yeasts, molds and aciduric microorganisms.	100 Gms 500 Gms	790 2690
93023	✔ **Fluid Thioglycollate Medium, Clear (USP Basis)** A clear medium used for sterility testing of variety of samples in pharmaceutical testing.	100 Gms 500 Gms	790 2240
42917	✔ **Fluid Thioglycollate Medium (USP, IP Basis)** Used for the cultivation of anaerobic, microaerophilic and aerobic organisms. For use in sterility testing of a variety of samples in compliance with USP and IP.	100 Gms 500 Gms	790 2240
46990	✔ **Fluid Thioglycollate Medium, Harmonized** Used for sterility testing of variety of samples in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	790 2240
76037	✔ **FNA Medium (Fluorescein Denitrification Agar)** Used for differentiation of Pseudomonas sps. from other bacilli by their ability to reduce nitrates or nitrites to nitrogen gas (denitrification) and detection of fluorescein pigment.	100 Gms 500 Gms	1350 4930
46142	✔ **Forget Fredette Agar**	100 Gms 500 Gms	1460 5380
33696	✔ **Formate Ricinoleate Broth**	100 Gms 500 Gms	1460 5720
37593	✔ **Fraser Broth Base** Recommended selective and enrichment medium for the isolation of Listeria species.	100 Gms 500 Gms	1240 5040
97539	**Fraser Supplement (250mg/vl)**	5 vl	1460
52071	🌡 **Fraser Selective Supplement** Used for selective isolation and cultivation of Listeria monocytogenes from food and environmental specimens.	5 vl	1800
83831	✔ **Fuchsin Lactose Broth**	100 Gms 500 Gms	1240 4600

Part D – Dehydrated Culture Media

Code	Name	Packing	Price
37346	✔ **Fungal Agar (Mycological Agar)** Used for selective isolation and maintenance of pathogenic fungi.	100 Gms 500 Gms	1460 5600
93411	✔ **Fungal Agar w/ low pH (Mycological Agar with Low pH)**	100 Gms 500 Gms	1460 5600
68653	✔ **Fungal Agar with Low pH (Mycological Agar w/ low pH)** For the selective isolation, cultivation and maintenance of pathogenic fungi.	100 Gms 500 Gms	1460 5720
62729	✔ **Fungal Broth (Mycological Broth)** For cultivation and maintenance of fungi.	500 Gms	5160
90195	✔ **Fungi Cultivation Kit (Pre-Weighed Media)** Used for cultivation and isolation of fungi,yeast and moulds.	1 Kit	1680
94016	✔ **Gassner Lactose Agar**	100 Gms 500 Gms	1350 5270
72464	✔ **GC Agar Base** Used for isolation and cultivation of Gonococci.	100 Gms 500 Gms	1120 4040
10474	✔ **Gelatin Agar DEV** Used in water and waste water testing for the determination of total microbial count and for the detection of gelatin liquefying microorganisms according to the German methods for the examination of water and the German drinking water regulations (1990).	100 Gms 500 Gms	1240 4480
23942	✔ **Gelatin Iron Agar** For the detection of gelatin liquefaction and H2S production.	100 Gms 500 Gms	1240 5160
75750 [91079-43-5]	✔ **Gelatin Peptone BactoBio for bacteriology** Gelatin peptone,prepared by enzymatic digest of Gelatin; is used as an ingredient in media for fermentation studies and to support growth of non-fastidious microorganisms. It is characterised by low Cystine,Tryptophan and Carbohydrate content .	500 Gms 5 Kg	2100 20300
43172 [9000-70-8]	✔ **Gelatin Powder Special ex. Porcine (Type B), 120g Bloom BactoBio for bacteriology** Gelatin is a collagenous protein which finds its application for the solidification of cultural media. It is also used for detection and differentiation of certain proteolytic bacteria. A bacteriological grade, readily soluble in water.	500 Gms	1800
24426 [9000-70-8]	✔ **Gelatin Powder Special Clear ex. Bovine (Type B), 200-250g Bloom BactoBio for bacteriology** Gelatin for bacteriology is a special bacteriological grade product, readily soluble in water to give a clear solution.	500 Gms	2690
78334	🌡 **Gentamicin Selective Supplement (20mg/vl)** An antibiotic supplement used for the isolation and detection of Clostridium perfringens.	5 vl	1240
50084	🌡 **Genta-Oxy Selective Supplement** Used for the selective isolation and enumeration of yeasts and moulds.	5 vl 5 x5 vl	2020 4480
39238	✔ **General Purpose Kit - I (Pre-Weighed Media)** Nutrient agar cultivates a wide variety of microorganism and MacConkey agar is used in isolation of gram negative enteric bacteria.	1 Kit	1350
84905	✔ **General Purpose Kit - II (Pre-Weighed Media)** Nutrient agar cultivates a wide variety of microorganisms,MacConkey agar aids in isolation of gram negative enteric bacteria and Blood Agar Base is used in detecting hemolytic of fastidious microorganisms.	1 Kit	1680
15628	✔ **Giollotti Cantoni Broth Base** Used for the cultivation and enrichment of Staphylococcus aureus from foods.	100 Gms 500 Gms	1570 4710
68382	✔ **Gluconate Test Medium**	100 Gms 500 Gms	1120 4260
35070	✔ **Glucose Salt Teepol Broth (Twin Pack)** Liquid medium,used for the detection of Vibrio parahemolyticus in foodstuffs.	100 Gms 500 Gms	1240 5160
58174	✔ **Glucose Yeast Extract Agar** For the cultivation of lactobacilli, especially in pharmaceutical preparations.	100 Gms 500 Gms	1460 5270
17991	✔ **Glucose Yeast Peptone Agar** For the isolation of yeasts from soil specimens.	100 Gms 500 Gms	1460 5270
79127	✔ **GN Broth, Hajna** (Gram-Negative Broth, Hajna) For the selective cultivation of gram negative microorganisms.	100 Gms 500 Gms	1010 3480
30963	🌡 **GTC Supplement (1.25mg/vl)** Used for the selective isolation and cultivation of Enterococci.	5 vl	1120

Code	Name	Packing	Price
11266	✔ **Glutamate Starch Phenol Red Agar Base** Used for detection of Pseudomonas and Aeromonas in foodstuffs, wastewater and equipment in food industry.	100 Gms 500 Gms	1460 4710
61655	✔ **Gum Listeria Medium** Recommended for the isolation of Listeria monocytogenes from clinical and non-clinical samples.	100 Gms 500 Gms	1350 5160
72275	✔ **H Broth** For the preparation of H agglutination antigen used in the differentiation and identification of Salmonella species types and subtypes.	100 Gms 500 Gms	1460 4710
85012	✔ **Hartley`s Digest Broth**	100 Gms 500 Gms	1120 4260
41861	🌡 **Haemoglobin Powder** Haemoglobin powder is used as an enrichment supplement . 2.0% w/v aq.solution of this supplement is autoclavable. Used with G.C.Agar base ,in preparation of Chocolate Agar.	50 Gms	1680
36277	🌡 **Haemoglobin Supplement**	50 Gms	1800
47704	✔ **Heart Infusion Agar (HI Agar)**	100 Gms 500 Gms	1350 5380
75134	✔ **Heart Infusion Broth (HI Broth)** Used for isolation and cultivation of a wide variety of fastidious microorganisms. Also used as a base for preparation of Blood Agar in determining haemolytic reactions.	100 Gms 500 Gms	1240 4930
28551	✔ **Heart Infusion Powder (HI Powder) BactoBio for bacteriology** It is rich in amino acids, peptides and other nutrients required by highly fastidious microorganisms like Streptococci, Pneumococci, Meningococi and Brucella species.	500 Gms 2.5 Kg	3030 13110
68100	✔ **Hektoen Enteric Agar** Used for the isolation and cultivation of gram negative enteric microorganisms from a variety of clinical and non clinical specimens.	100 Gms 500 Gms	1460 6500
20852	✔ **Hemmes Medium Base**	100 Gms 500 Gms	2800 8520
52740	✔ **Hofers Alkaline Medium**	100 Gms	1800
67658	✔ **Horie Arabinose Ethyl Violet Broth**	100 Gms 500 Gms	1570 6050
46946	✔ **Hoyle Medium Base** For the isolation and differentiation of Corynebacterium diphtheriae.	100 Gms 500 Gms	1240 4710
79035	✔ **High Plate Count Agar** Recommended medium for obtaining higher colony counts by spread plate/pour plate or membrane filteration technique.	100 Gms 500 Gms	1460 5720
64068	✔ **High Salt Nutrient Agar (ISO Basis)** For the isolation and cultivation of salt tolerant Vibrio species.	100 Gms 500 Gms	1240 4600
37353	✔ **High Salt Nutrient Agar (BIS Basis)**	100 Gms 500 Gms	1240 4600
39147	✔ **High Salt Peptone Yeast Extract Agar (ISO Basis)**	100 Gms 500 Gms	1240 4480
77707	✔ **HS Medium (Hydrosulphite of Sodium Medium)** Used for the cultivation of aerobic and anaerobic bacteria. It is also used for sterility testing.	100 Gms 500 Gms	1240 3920
26221	✔ **Hugh Leifson Medium (BIS Basis)** Used to distinguish between anaerobic and aerobic breakdown of carbohydrate (Dextrose).	100 Gms 500 Gms	1350 5160
79544	✔ **HYA Agar**	100 Gms 500 Gms	2020 7510
21128	✔ **Hydrolysate BioVeg BactoBio for bacteriology** Prepared by enzymatic hydrolysis of vegetable proteins under controlled conditions.Used as a nutritive source for growth of bacteria and fungi.	500 Gms 5 Kg	3030 23300
36873	✔ **IMViC Test Kit (Pre-Weighed Media)** Used for biochemical identification of gram negative enteric bacteria based on Indole production, Methyl Red and Voges Proskauer reactions and Citrate Utilization test.	1 Kit	1460
68788	✔ **IMViC Test Kit (Pre-Weighed Media) w/ Reagents** Used for biochemical identification of gram negative enteric bacteria based on Indole production, Methyl Red and Voges Proskauer reactions and Citrate Utilization test.	1 Kit	5940
12139	✔ **Indole Nitrate Medium** (Tryptone Nitrate Medium) For the identification of microorganisms by means of nitrate reduction and indole tests.	100 Gms 500 Gms	1120 3810

Code	Name	Packing	Price
54336	✔ **ISP Medium No. 1 (Tryptone Yeast Extract Broth)** Used for characterizing Streptomyces species according the International Streptomyces project.	100 Gms 500 Gms	1570 6610
76920	✔ **ISP Medium No.4** (Inorganic Salt Starch Agar) Used for characterizing Streptomyces species according the International Streptomyces project.	100 Gms 500 Gms	1800 7060
91875	✔ **ITC Broth**	100 Gms 500 Gms	1570 5830
42882	✔ **IUT Medium Base**	100 Gms 500 Gms	1680 7510
70916	✔ **Jensen`s Medium (Jensen`s Agar)** Used for the cultivation of nitrogen fixing bacteria.	100 Gms 500 Gms	1240 4600
87806	✔ **Jensens Broth** Recommended for the detection and cultivation of nitrogen fixing bacteria.	100 Gms 500 Gms	1570 5600
55686	✔ **Jensen Seedling Agar** Used for growing plants for collection of root exudates and to study the nodulating ability of Rhizobium isolates.	100 Gms 500 Gms	1350 4930
44485	✔ **Kanamycin Esculin Azide Agar** Used for the isolation of enterococci in foodstuffs.	100 Gms 500 Gms	2020 7510
50136	✔ **Kanamycin Esculin Azide Broth Base** Recommended for selective isolation and identification of group D Streptococci in foodstuffs	100 Gms	4480
60649	🌡 **Kanamycin Sulphate Selective Supplement (10mg/vl)** Used for the selective isolation and cultivation of Enterococci.	5 vl	1570
96343	✔ **Kaper`s Medium**	100 Gms 500 Gms	1460 5830
33613	✔ **Kenknight and Munaier`s Medium** Used for the isolation of Actinomyces species from soil samples.	100 Gms 500 Gms	1460 5720
75772	✔ **Ketogluconate Broth**	100 Gms 500 Gms	1460 4260
64342	✔ **KF Streptococcus Agar Base (APHA)** For the selective cultivation and enumeration of faecal streptococci.	100 Gms 500 Gms	1240 5040
66947	✔ **KG Agar Base** Recommended medium for promoting fast and free sporulation to differentiate between Bacillus cereus and Bacillus thuringiensis .	100 Gms 500 Gms	1800 6280
41924	✔ **Kimmig Agar Base** (Fungi Agar Base Kimmig)	100 Gms 500 Gms	1460 5720
42820	🌡 **Kimmig Selective Supplement** Used for selective isolation of fungi.	5 vl	2690
11411	✔ **King`s Medium A Base** For the non selective isolation, cultivation and pigment production of Pseudomonas.	100 Gms 500 Gms	1120 4600
21694	✔ **King's Medium B Base** For the non-selective isolation, cultivation and pigment production of Pseudomonas species.	100 Gms 500 Gms	1120 4600
29760	🌡 **Klebsiella Selective Supplement (25mg/vl)** Used for the selective isolation and detection of Klebsiella species.	5 vl	1570
33943	✔ **Kligler Iron Agar** For the differentiation and identification of gram-negative bacilli based upon the fermentation of dextrose and lactose and production of H2S.	100 Gms 500 Gms	1120 4150
61860	✔ **Kligler Iron Agar (ISO Basis)** For the differentiation and identification of gram-negative bacilli based upon the fermentation of dextrose and lactose and production of H2S. For identification of Pseudomonas species.	100 Gms 500 Gms	1120 4600
25844	✔ **Kligler Iron Agar BioVeg** For the differentiation and identification of gram-negative bacilli based upon the fermentation of dextrose and lactose and production of H2S.	100 Gms 500 Gms	1120 4150
44000	✔ **Kohn Two Tube Medium No. 1 Base** Used for identification of Enterobacteriaceae based on dextrose fermentation, mannitol fermentation and urease production.	100 Gms 500 Gms	3030 7960
58578	✔ **Kohn Two Tube Medium No. 2 Base** Used for identification of Enterobacteriaceae based on sucrose fermentation, salicin fermentation, motility, H2S production and indole production.	100 Gms 500 Gms	7510 27560

Code	Name	Packing	Price
65373	✔ **Koser Citrate Medium** For the differentiation of Escherichia coli and Enterobacter aerogenes based on the basis of citrate utilization.	100 Gms 500 Gms	1120 4040
77760	✔ **Kundrat Agar**	100 Gms 500 Gms	1570 5600
70663	🌡 **KV Supplement**	5 vl	2130
52911	✔ **Lachica`s Medium Base**	100 Gms	4930
63026 [91079-44-6]	✔ **Lactalbumin Hydrolysate BactoBio for bacteriology** A hydrolysate obtained by enzymatic digestion of protein, Lactalbumin. It is used primarily in the preparation of tissue culture media, fermentation media, media for the growth of Lactobacilli, Campylobacter sps, etc..	500 Gms	5160
24806	🌡 **L-Lactic Acid 10% Solution (10ml/vl)** Used for adjusting pH in acidic range.	5 vl 5 x5 vl	1350 4150
37985	✔ **Lactic Agar** Recommended for enumeration and identification of lactic Streptococci and Lactobacilli by pour plate method.	100 Gms	2020
68888	✔ **Lactic Bacteria Differential Agar**	100 Gms 500 Gms	1460 5380
84282	🌡 **Lactic Supplement**	5 vl	1680
39785	🌡 **Lactobacilli Supplement (50mg/vl)** An antibiotic supplement recommended for the selective isolation of Lactobacilli from wine.	5 vl 5 x5 vl	2020 5940
79562	🌡 **Lactobacillus MRS Agar (ISO Basis)** Used for the isolation and enumeration of lactic acid bacteria from meat and meat products.This medium is in accordance to specifications given by ISO 13721 : 1995.	100 Gms 500 Gms	1240 4370
50871	✔ **Lactobacillus MRS Agar BioVeg** Used for the isolation and enumeration of lactic acid bacteria from meat and meat products.	100 Gms 500 Gms	1240 5160
46812	✔ **Lactobacillus MRS Broth w/o Dextrose** Used for the selective isolation and enumeration of lactobacilli from various sources including feces and dairy.	100 Gms 500 Gms	1120 4150
47192	✔ **Lactobacilli Agar, AOAC** For the cultivation and maintenance of stock cultures used in microbiological assays of vitamins and amino acids.	100 Gms 500 Gms	2360 6840
18659	🌡 **Lactobacillus Selection Agar** Used for isolation, cultivation and enumeration of lactobacilli.	100 Gms 500 Gms	1460 5600
66749	✔ **Lactobacillus Selection Oxgall Agar Base (LBS Oxgall Agar)** Used for the cultivation and enumeration of lactobacilli.	100 Gms 500 Gms	1570 5380
49190	🌡 **Lactobacillus MRS Broth (MRS Broth)** For the cultivation of lactobacilli.	100 Gms 500 Gms	1120 4150
27750	✔ **Lactobacillus MRS Broth BioVeg** Used for the cultivation of lactobacilli.	100 Gms 500 Gms	1120 4370
97354	✔ **Lactose Broth** For the detection of lactose fermenting gram-negative coliforms, as pre-enrichment broth for Salmonella species.	100 Gms 500 Gms	790 3030
77121	✔ **Lactose Gelatin Medium (ISO Basis)** For the cultivation of Clostridium perfringens.	100 Gms 500 Gms	1240 4600
24171	✔ **Lactose Lecithin Agar** Recommended for isolation and differentiation of histotoxic Clostridia from clinical specimens.	100 Gms 500 Gms	1680 5720
88539 [10039-26-6, 64044-51-5]	**Lactose Monohydrate sterile (gamma irradiated)** Lactose monohydrate, Sterile (Gamma Irradiated) is used for Media Fill Runs in dry injectable in pharmaceutical industry.	500 Gms 5 Kg	1240 10760
62367	✔ **Lactose Peptone Broth** Used for the detection of coliform organisms in water.	100 Gms 500 Gms	1240 4150
71135	✔ **Lactose Sulphite Broth Base (ISO Basis)** Recommended for detection and enumeration of Clostridium perfringens in pharmaceutical products.	100 Gms 500 Gms	1240 4260

Code	Name	Packing	Price
70637	✔ **Lactose TTC Agar (w/ Tergitol 7) (ISO Basis)** (Lactose Triphenyltetrazolium Chloride Agar) Selective differential medium for the detection and enumeration of Escherichia coli and coliform bacteria in water using the membrane filtration method.	100 Gms 500 Gms	1240 4370
92343	✔ **Lauryl Sulfate Broth w/MUG**	100 Gms 500 Gms	1460 5380
44956	✔ **Lauryl Sulfate Broth (Lauryl Tryptose Broth) (ISO, BIS Basis)** For the selective enrichment and presumptive detection of coliform bacteria from food , animal feed,water and other materials.This culture medium is in accordance with the specifications given by ISO 4831, ISO 7251 and APHA.	100 Gms 500 Gms	1010 3480
94906	✔ **Lauryl Tryptose Mannitol Broth w/Tryptophan** For the detection of Escherichia coli in water samples.	100 Gms 500 Gms	1460 4930
10517	✔ **Lead Acetate Agar** For the cultivation and differentiation of gram negative coliform bacteria based on H2S production.	100 Gms 500 Gms	1350 4930
82632	✔ **Leptospira Medium Base** Used for the cultivation of Leptospira species.	100 Gms 500 Gms	1240 4710
59326	✔ **Letheen Agar** To determine the phenol coefficient of quaternary ammonium compounds using Escherichia coli or Staphylococcus aureus ATCC 6538.	500 Gms	6280
42264	✔ **Letheen Broth (USP Basis)** Recommended medium for determination of bactericidal activity of quaternary ammonium compounds.This medium is in accordance to specifications given by United States Pharmacopoeia-National Formulatory (USP-NF), 2022.	100 Gms 500 Gms	1500 5500
27680	✔ **Letheen Broth, Modifed (Twin Pack)** Recommended medium for screening cosmetic products for microbial contamination.	100 Gms 500 Gms	1800 6000
28685	✔ **Listeria Enrichment Medium Base (UVM Medium)** Recommended medium for use in qualitative procedures for the isolation of Listeria monocytogenes from food samples.	100 Gms 500 Gms	1800 8180
96310	✔ **Listeria Identification Agar Base (PALCAM)** Used for isolation of Listeria sp. from clinical and non-clinical samples.	100 Gms 500 Gms	2020 8850
96155	✔ **Listeria Motility Medium (ISO Basis)** Listeria Motility Medium is a semisolid medium recommended for testing motility of Listeria monocytogenes	100 Gms 500 Gms	1350 5160
15650	🌡 **Listeria Moxalactum Supplement** Recommended for selective isolation and cultivation of Listeria monocytogenes from food and environmental specimens.	5 vl	2020
76366	✔ **Listeria Oxford Medium Base (Oxford Agar)** Used for isolation of Listeria species from pathological samples.	100 Gms 500 Gms	1350 6050
90345	🌡 **Listeria Selective Supplement (PALCAM)** Used for the isolation of Listeria monocytogenes and other species of Listeria from food products.	5 vl	1460
51864	🌡 **Listeria UVM Supplement I**	5 vl	1350
72547	🌡 **Listeria UVM Supplement II**	5 vl	1350
59797	✔ **Litmus Lactose Agar** For the differentiation of lactose fermenting and non-lactose fermenting bacteria.	100 Gms 500 Gms	1350 5040
96658	✔ **Litmus Milk (Litmus SM Broth)** For the maintenance of lactic acid bacteria and the differentiation of several bacteria, based on their action on milk.	100 Gms 500 Gms	1350 4710
50037	✔ **Liquoid Broth**	100 Gms 500 Gms	1460 4930
54004	✔ **Littman Oxgall Agar Base** With added streptomycin, used for isolation of fungi , especially dermatophytes	100 Gms 500 Gms	1460 5160
87566	✔ **Liver Extract Powder BactoBio for bacteriology** Nutritive characteristic are good in Liver extract so used for cultivation of anaerobic bacteria such as Brucellae and Clostridia, also used for cultivation of pathogenic bacteria.	500 Gms 5 Kg	2240 18370

Code	Name	Packing	Price
81322	✔ **Liver Hydrolysate BactoBio for bacteriology** (Peptone from Liver, Liver Peptone) Used as source of nutrients in microbiology culture media.Water soluble; being derived from liver contains relatively high levels of iron.	500 Gms	3480
68861	✔ **Liver Infusion Agar (LI Agar)**	100 Gms 500 Gms	1240 4480
96453	✔ **Liver Infusion Broth (LI Broth)** Used for cultivation of Brucella and other pathogenic anaerobic bacteria.	100 Gms 500 Gms	1120 3810
98858	✔ **Liver Meat Glucose Cysteine Broth**	100 Gms 500 Gms	1460 5600
79612	✔ **LM Agar (Meat Liver Agar)** Used for the cultivation of anaerobic microorganisms, specially for the sulphite reducing Clostridia and Clostridium perfringens in water and milk.	500 Gms	5040
69741	✔ **Lowenstein-Jensen Medium Base** (LJ Medium) For the cultivation and differentiation of Mycobacterium species.	100 Gms 500 Gms	1010 3360
27888	✔ **LPM Agar Base** Used for isolation and cultivation of Listeria monocytogenes from food and dairy products.	100 Gms	3920
46502	✔ **Luria Bertani Agar, Lennox** (Luria Agar, LB Agar, Lennox) Used for maintaining and cultivating recombinant strains of Escherichia coli.	100 Gms 500 Gms	1240 4710
47436	✔ **Luria Bertani Agar, Miller** Used for maintaining and propagating Escherichia coli in molecular microbiology procedures.	100 Gms 500 Gms	1120 4480
14593	✔ **Luria Bertani Broth, Lennox (Luria Broth)** (Luria Broth, LB Broth, Lennox) Used for maintaining and propagating Escherichia coli in molecular biology procedures.	100 Gms 500 Gms 2 Kg	1120 4040 12880
29817	✔ **Luria Bertani Broth, Miller** (LB Broth, Miller) Used for maintaining and propagating Escherichia coli in molecular microbiology procedures.	100 Gms 500 Gms 2 Kg	1120 4040 12880
40826	✔ **Luria Bertani Agar, Miller BioVeg** (LB Agar, Miller BioVeg) Used for maintaining and propagating Escherichia coli in molecular microbiology procedures.	100 Gms 500 Gms	1120 4480
17320	✔ **Luria Bertani Broth, Miller BioVeg** (LB Broth, Miller BioVeg) Used for maintaining and propagating Escherichia coli in molecular microbiology procedures.	100 Gms 500 Gms 2 Kg	1120 4040 12880
97970	✔ **Luria Bertani Agar w/o Sodium Chloride** Used for maintaining and propagating Escherichia coli in molecular microbiology procedures.	100 Gms 500 Gms	1570 5380
22006	✔ **Luria Bertani Broth w/o Sodium Chloride (LB Growth Medium w/o Sodium Chloride)** Used for maintaining and propagating Escherichia coli in molecular microbiology procedures.	100 Gms 500 Gms	1120 4480
90769	✔ **Lysine Arginine Iron Agar (LAI Agar)** For the cultivation and differentiation of bacteria based on their ability to decarboxylate lysine, arginine and produce H2S.Recommended by APHA for isolation and identification of Yersinia from milk and milk products.	100 Gms 500 Gms	1120 4040
16624	✔ **Lysine Decarboxylase Broth w/o Peptone (ISO Basis)** For distinguishing Salmonella serotype arizonae from the Bethesda Ballerup group of Enterobacteriaceae.This medium is in accordance to specifications given by ISO 6579-1:2017 and ISO 10273:2003.	100 Gms 500 Gms	1010 3810
71125	✔ **Lysine Iron Agar** For the differentiation of enteric organisms, especially Salmonella, based upon lysine decarboxylation/ deamination and H2S production.	100 Gms 500 Gms	1120 3810
82502	✔ **M2 Agar** A semisynthetic culture medium, used as a general purpose plate count agar.	100 Gms 500 Gms	1350 5040
24415	✔ **M17 Agar Base** For improved growth of lactic streptococci and their bacteriophages.	100 Gms 500 Gms	1350 4930
13297	✔ **M17 Broth Base** For improved growth of lactic streptococci and their bacteriophages.	100 Gms 500 Gms	1120 4480

Code	Name	Packing	Price
76875	✔ **MacConkey Agar** For the selective isolation, cultivation and differentiation of coliforms and enteric pathogens based on the ability to ferment lactose.	100 Gms 500 Gms	1010 3920
96871	✔ **MacConkey Agar BioVeg** For the selective isolation, cultivation and differentiation of coliforms and enteric pathogens based on ability to ferment lactose.	100 Gms 500 Gms	1010 3920
37775	✔ **MacConkey Agar, Harmonized** Used for the selective isolation and differentiation of Escherichia coli and coliform bacteria in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1010 3920
80302	✔ **MacConkey Agar w/ MUG**	100 Gms 500 Gms	1800 7280
71117	✔ **MacConkey Agar (BIS Basis)**	100 Gms 500 Gms	1010 3920
83870	✔ **MacConkey Agar Base** Used for studying carbohydrate fermentation reactions of coliforms by adding desired carbohydrates.	100 Gms 500 Gms	1010 3920
26823	✔ **MacConkey Agar Medium** For the isolation and differentiation of enteric bacteria on the basis of lactose fermentation and also for isolation of faecal streptococci.	100 Gms 500 Gms	1120 4040
95199	✔ **MacConkey Agar, Modified**	100 Gms	3590
98384	✔ **MacConkey Agar w/ Bromothymol Blue**	100 Gms 500 Gms	1240 4370
58370	✔ **MacConkey Agar w/ 0.15% Bile salts, C.V. and NaCl (USP Basis)** Selective medium for gram-negative organisms by suppressing the growth of a number of gram positive bacteria.	100 Gms 500 Gms	1010 4260
33759	✔ **MacConkey Agar w/o C.V., w/ 1.2% Agar** For the selective isolation and differentiation of enteric bacteria on the basis of lactose fermentation.	100 Gms 500 Gms	1240 4600
50372	✔ **MacConkey Agar w/o C.V., NaCl w/ 0.5% Sodium Taurocholate** Differential medium for detection and isolation of enteric organisms by restricting the swarming of Proteus species.	100 Gms 500 Gms	1460 4600
79447	✔ **MacConkey Agar w/o C.V. w/ 0.15% Bile Salts** For the selective isolation of members of Enterobacteriaceae on the basis of lactose fermentation.	100 Gms 500 Gms	1010 3480
64300	✔ **MacConkey Agar w/o C.V. w/ 0.5% Bile Salts** For isolation and detection of coliforms and enteric pathogens.	100 Gms 500 Gms	1010 4040
64024	✔ **MacConkey Agar w/o C.V., NaCl, w/ 0.0075% NR, 1.2% Agar BioVeg**	100 Gms 500 Gms	1120 3920
64030	✔ **MacConkey Agar w/o C.V., NaCl w/ 0.5% Bile Salts** For cultivation and differentiation of enteric bacteria on which most Proteus species will not swarm thus preventing overgrowth on the plate and also restricting growth of pathogenic gram positive organisms.	100 Gms 500 Gms	1010 4040
72074	✔ **MacConkey Agar No.3** For the selective isolation, cultivation and differentiation between coliforms and non -lactose fermentors.	100 Gms 500 Gms	1460 4710
74726	✔ **MacConkey Broth (MacConkey Broth w/ Neutral Red)** For the selective isolation and cultivation of coliforms in milk and water.	100 Gms 500 Gms	1010 3360
44448	✔ **MacConkey Broth, Harmonized** Used for the identification of coliforms in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1010 3360
54056	✔ **MacConkey Broth w/ Neutral Red (BIS Basis)**	100 Gms 500 Gms	1010 3360
20309	✔ **MacConkey Broth Double Strength w/ Neutral Red** For the selective isolation and cultivation of coliforms in milk and water from large samples (water and waste water samples).	100 Gms 500 Gms	900 3140
92933	✔ **M-MacConkey Broth** For detecting and enumerating lactose fermenting enteric bacteria from milk and water by the membrane filter method.	100 Gms 500 Gms	1010 3920
42944	✔ **MacConkey Broth Purple** Selective medium for isolation and cultivation of coliforms in milk and water.	100 Gms 500 Gms	1010 3480
38989	✔ **MacConkey Broth Purple w/BCP (ISO Basis)** For presumptive identification of coliforms from water, milk and food samples.This medium is in accordance to the specifications given by ISO 9328-2:1990.	100 Gms 500 Gms	900 3140

Code	Name	Packing	Price
90084	✔ **MacConkey Broth Purple w/BCP (IP Basis)** For the selective identification of E. coli from pharmaceutical products .This medium is in accordance with Indian Pharmacopoeia 2014.	100 Gms 500 Gms	1120 3590
51545	✔ **MacConkey Broth Purple w/BCP (BIS Basis)** For presumptive identification of coliforms from water, milk, food etc.The composition and specification are in accordance to BIS ,IS : 5401-1969.	100 Gms 500 Gms	900 3140
57608	✔ **MacConkey Broth Purple Double Strength w/ BCP** For the presumptive identification of coliforms from large water and wastewater samples.	100 Gms 500 Gms	1120 3810
98664	✔ **MacConkey Sorbitol Agar (Sorbitol Agar)** For the isolation and identification of pathogenic Escherichia coli.	100 Gms 500 Gms	1570 5160
57340	✔ **MacConkey Sorbitol Agar Base (ISO Basis)** For the selective isolation and detection of Escherichia coli 0157:H7 from food and animal feeds stuff. The formulation and composition are in accordance to ISO 16654:2001 ,(Amd 2:2023)	100 Gms 500 Gms	1460 5040
99130	✔ **Malachite Green Broth** Used for cultivation of Pseudomonas aeruginosa.	100 Gms 500 Gms	1240 3920
51517	✔ **M-Brilliant Green Broth** For the selective isolation and differentiation of Salmonella from polluted water by the membrane filter method.	100 Gms 500 Gms	1460 4710
10054	✔ **McBride Listeria Agar Base** Moderately selective medium for isolation of Listeria monocytogenes.	100 Gms 500 Gms	1460 5720
70317	🌡 **McBride Listeria Supplement** An antibiotic supplement used for the selective isolation of Listeria species.	5 vl	1680
50642	✔ **McClung Toabe Agar Base**	100 Gms 500 Gms	1460 5830
87015	✔ **M-EC Test Agar**	100 Gms 500 Gms	1350 4710
66805	✔ **M-Endo Agar LES** For the enumeration of coliform in water by the membrane filter method.	100 Gms 500 Gms	1240 4480
57906	✔ **M-Endo Broth**	100 Gms 500 Gms	1240 4600
33919	✔ **M-Enrichment Broth**	100 Gms 500 Gms	1240 4600
52295	✔ **M-FC Agar Base** For the cultivation and enumeration of faecal coliforms by membrane filter technique at elevated temperatures.	100 Gms 500 Gms	1460 5720
76439	✔ **M-FC Broth Base** For the cultivation and enumeration of faecal coliforms by membrane filter technique.	100 Gms 500 Gms	1460 4710
51072	✔ **M-HPC Agar Base** (Heterotrophic Plate Count Agar, Base) For the cultivation and enumeration of heterotrophic organisms from potable water sources, swimming pools and other water specimens by the membrane filter method.	100 Gms 500 Gms	1240 5040
14750	✔ **M-Lauryl Sulfate Broth (ISO Basis)**	100 Gms 500 Gms	1460 4710
21051	✔ **M-Staphylococcus Broth** Used for detection and isolation of Staphylococci by memebrane filter technique.	100 Gms 500 Gms	1240 4150
51416	✔ **MLCB Agar Mannitol Lysine Crystal Violet Brilliant Green Agar** (Mannitol Lysine Agar)	100 Gms 500 Gms	1460 5600
99505	✔ **Malonate Broth** Used for differentiating Enterobacter and Escherichia coli based on malonate utilization.	100 Gms	4930
76507	✔ **Malt Agar** For the cultivation of yeasts and molds from dairy products, foods and other materials.	100 Gms 500 Gms	1680 6500
84980	✔ **Malt Extract Agar** For the isolation, cultivation and enumeration of yeasts and molds.	100 Gms 500 Gms	1460 5600
42936	✔ **Malt Extract Agar Base, Modified as per Thom and Church** Used for isolation, detection and enumeration of yeast and molds.	100 Gms 500 Gms	1570 5940

Part D – Dehydrated Culture Media

Code	Name	Packing	Price
83358 [8002-48-0]	✔ **Malt Extract BactoBio for bacteriology** Malt extract contain high concentration of carbohydrates particularly maltose. It also provides an acidic environment, carbon, protein and nutrients for the isolation and cultivation of yeasts and molds, restricting the growth of bacteria and also for the cultivation of yeasts and molds.	500 Gms 5 Kg	1460 13440
67973	✔ **Malt Extract Broth Base** For the isolation, cultivation and enumeration of yeasts and molds.	100 Gms 500 Gms	1570 6160
40828 [69-65-8]	**D-Mannitol Sterile (Gamma Irradiated)** For use in microbiological applications using D-Mannitol (gamma irradiated) and in Media Fill Runs in dry injectables.	500 Gms 5 Kg	4260 39090
48330	✔ **Mannitol Agar w/ Prilion**	100 Gms 500 Gms	5720 22070
75485	✔ **Mannitol Lysine Agar**	100 Gms 500 Gms	1350 5270
62331	✔ **Mannitol Selenite Broth Base** (Selenite Mannitol Broth) A modification of Selenite F Broth used for selective enrichment of Salmonellae.	100 Gms 500 Gms	1460 5720
92493	✔ **Mannitol Motility Test Medium** Used for the detection of motility and mannitol fermentation in microorganisms.	100 Gms 500 Gms	1350 4370
48939	✔ **Mannitol Nitrate Motility Medium** Used for the rapid identification of Enterobacteriacae based on motility, mannitol fermentation and nitrate reduction.	100 Gms 500 Gms	1240 4260
60062	✔ **Mannitol Salt Agar (USP Basis)** Used for the selective isolation, cultivation and enumeration of staphylococci from clinical and non-clinical specimens.This medium is recommended in U.S.P. for use in the performance of microbial limit test.	100 Gms 500 Gms	1010 3480
23626	✔ **Mannitol Salt Agar, Harmonized** Used for the isolation and cultivation of staphylococci in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1010 3480
55066	✔ **Mannitol Salt Broth** Used for the selective isolation and cultivation of staphylococci.	100 Gms 500 Gms	1120 3810
25511	✔ **Marine Agar 2216 (Zobell)** Used for cultivation of heterotrophic marine bacteria.	100 Gms 500 Gms	1460 5720
54174	✔ **Marine Broth 2216 (Zobell)** Used for cultivation of heterotrophic marine bacteria.	100 Gms 500 Gms	1350 5380
34832	✔ **Maximum Recovery Diluent** Physiologically isotonic and protective medium for maximum recovery of microorganisms from a variety of sources.	100 Gms 500 Gms	1240 4040
78784 [100085-61-8]	✔ **Meat Extract BactoBio for bacteriology** Meat extract,manufactured from lean meat used to complement nutritive properties of peptone.	500 Gms 5 Kg	3360 28680
70117	✔ **Meat Extract Powder BioVeg BactoBio for bacteriology** Meat extract, (animal free) manufactured from vegetable/plant source, is used to complement nutritive properties of peptone. It is equivalent to animal based Meat extract. Its distribution of various peptides meets the requirements for replacing Meat extract (animal based) to yield optimum growth of microorganisms.	500 Gms 5 Kg	3030 23300
17585	✔ **Meat Extract w/ Peptone** (Pepted Meat Broth)	100 Gms 500 Gms	1350 4930
11980	✔ **Meat Peptone BactoBio for bacteriology** A highly nutritious, enzymatic digest of lean meat. Used as a ingredient in various general purpose media for cultivation and isolation of wide variety of microorganisms.	500 Gms	3480
91980	✔ **Meat Peptone Type P BactoBio for bacteriology** It is obtained by proteolytic digestion of meat with pepsin to retain nutritive properties of peptone. Suitable for growth of fastidious microorganisms.	500 Gms	4600
11828	✔ **Meat Peptone Type T BactoBio for bacteriology** It is obtained by proteolytic digestion of meat with Trypsin. Suitable for mass scale cultivation of microorganisms for production of antibiotics, vitamins, enzymes etc.	500 Gms	4710
28238	✔ **Micro Inoculum Broth** (Micro Vitamin Test Inoculum Broth) For the preparation of inocula in microbiological assay and cultivation of lactobacilli.	100 Gms 500 Gms	1350 6160

Code	Name	Packing	Price
63261	✔ **Microbial Content Test Agar** For use in the microbial content test of water soluble cosmetic products. Also used for determining the efficiency of sanitization of containers, equipment and environmental surfaces	100 Gms 500 Gms	1800 7280
61200	✔ **Milk Agar (Modified Plate Count Agar)** For the enumeration of microorganisms in milk, icecreams, whey, etc.	100 Gms 500 Gms	1570 5830
18789	✔ **Milk Agar w/Cetrimide (Twin Pack) (Modified Pseudomonas Selective Agar w/Cetrimide)** For detection and enumeration of Pseudomonas aeruginosa in water samples.	100 Gms 500 Gms	1570 6390
93957	✔ **Milk Medium w/Reducing Agent (Modified Litmus Medium w/Reducing Agent)**	100 Gms 500 Gms	1240 5270
17595	✔ **Milk Salt Agar Base (Modified Salt Agar Base for Staphylococci)**	100 Gms 500 Gms	1570 5600
20767	✔ **Minimal Agar, Davis** For the isolation, cultivation and characterization of nutritional mutants of Escherichia coli.	100 Gms 500 Gms	1570 6050
39821	✔ **Minimal Broth, Davis** For the isolation, cultivation and characterization of nutritional mutants of Escherichia coli.	100 Gms 500 Gms	1350 4370
24911	✔ **Minimal Broth, Davis w/o Dextrose** Used for isolation and characterization of nutritional mutants of Escherichia coli and Bacillus subtilis.	100 Gms 500 Gms	1350 4370
88736	✔ **MIO Medium (Motility Indole Ornithine Medium)** For differentiation of Enterobacteriaceae based on their motility, indole production and ornithine decarboxylase activity.	100 Gms 500 Gms	1240 4710
46584	✔ **Mitis Salivaris Agar Base** Selective medium for the isolation of streptococci and enterococci from clinical specimens.	500 Gms	5380
16309	✔ **MIU Medium Base (Motility Indole Urea Medium)** Used for detecetion of motility, urease and indole production.	100 Gms 500 Gms	1460 5270
79667	✔ **M-TEC Agar** Recommended medium for isolation, differentiation and enumeration of Escherichia coli from water by membrane filtration.	100 Gms 500 Gms	1460 5040
62811	✔ **M-Nutrient Broth** For enumerating bacteria by membrane filter technique.	100 Gms 500 Gms	1460 4480
85769	✔ **Modified Rappaport Vassiliadis Medium** Used as a selective enrichment medium for the isolation of Salmonella species from food and environmental specimens.	100 Gms 500 Gms	1570 5600
82453	✔ **Modified Skim Milk Agar (Modified SM Agar)** For cultivation and enumeration of microorganisms encountered in dairy industries.	100 Gms 500 Gms	2130 6280
89414	✔ **Modified Tergitol-7 Agar Base (ISO Basis)** (Tergitol-7 Agar Base Modified)	100 Gms 500 Gms	1570 5830
30066	✔ **Moeller Decarboxylase Broth Base (Decarboxylase Broth Base, Moeller)** With the addition of L-Lysine/ L-Arginine/ L-Ornithine used for the differentiation of bacteria on the basis of their ability to decarboxylate the amino acid.	100 Gms 500 Gms	2580 9970
18124	✔ **Moeller Decarboxylase Broth w/ Lysine Hydrochloride** Used for differentiation of bacteria on the basis of their ability to decarboxylate L-lysine hydrochloride.	100 Gms 500 Gms	3480 16470
27510	✔ **Moeller Decarboxylase Broth w/ Ornithine Hydrochloride**	100 Gms 500 Gms	3590 17030
69483	✔ **Moeller Decarboxylase Broth w/ Arginine Hydrochloride** Used for differentiation of bacteria on the basis of their ability to decarboxylate L-Arginine hydrochloride.	100 Gms 500 Gms	3480 16470
34452	✔ **Monsur Medium Base**	100 Gms 500 Gms	1240 4930
68957	✔ **M9CA Medium** For the cultivation of recombinant strains of Escherichia coli.	100 Gms 500 Gms	3250 12880
23184	✔ **M9 Minimal Salts 5X** For the cultivation of recombinant strains of Escherichia coli.	100 Gms 500 Gms	2240 7510

Part D – Dehydrated Culture Media

Code	Name	Packing	Price
42853	✔ **Motility GI Medium** Used for detection of motility of microorganisms.	100 Gms 500 Gms	2920 11200
90174	✔ **Motility Medium S Base** (Motility Nitrate Medium) For the determination of bacterial motility by means of TTC reduction.	100 Gms 500 Gms	1120 4260
21740	✔ **Motility Sulphide Medium** For the determination of bacterial motility and the ability of bacteria to produce H2S from L- Cystine.	100 Gms 500 Gms	1240 4710
64421	✔ **Motility Test Medium** For the determination of bacterial motility.	100 Gms 500 Gms	1240 4480
17624	✔ **MPH Agar (Milk Protein Hydrolysate Agar)** For the enumeration of bacteria in water and dairy products.	100 Gms 500 Gms	1680 6720
14833	✔ **MR-VP Medium (Glucose Phosphate Broth, Buffered Glucose Broth)** (MRVP Broth) For the differentiation of bacteria based on acid production (Methyl Red test) and acetoin production (Voges-Proskauer reaction).	100 Gms 500 Gms	900 3360
14443	✔ **Multipurpose Test Kit (Pre-Weighed Media)** Used for isolation, idetification and cultivation of a wide variety of microorganism.	1 Kit	4710
24756	✔ **Mueller Hinton Agar (MH Agar)** For antimicrobial susceptibility testing of a variety of microorganisms. For the cultivation and maintenance of Moraxella osloensis and Neisseria meningitidis.	100 Gms 500 Gms 2 Kg	1120 4260 13890
26351	✔ **Mueller Kauffman Tetrathionate Broth Base** An improved enrichment medium for the isolation of Salmonellae and the suppression of Proteus species.	100 Gms 500 Gms	1120 3360
76135	✔ **Mueller Kauffman Tetrathionate Novobiocin Broth Base (MKTTn)** Used for selective enrichment for the isolation of Salmonella, when used with MKTT Novobiocin supplement.	100 Gms 500 Gms	1350 4710
29465	✔ **Mucate Broth**	100 Gms 500 Gms	7840 36180
38800	✔ **Mucate Control Broth** Used for identification of enteropathogenic Escherichia coli and Salmonella species from milk and milk products.	500 Gms	4260
49550	✔ **Mueller Hinton Broth (MH Broth)** For the cultivation of a wide variety of organisms and for antimicrobial susceptibility testing.	100 Gms 500 Gms 2 Kg	1120 3920 11650
24698	🌡 **MUG EC Broth** For detection of Escherichia coli in water and food samples by fluorogenic method.	100 Gms 500 Gms	2920 8070
26834	🌡 **MUG Supplement (50mg/vl)** A fluoregenic substrate, used for the detection of Escherichia coli .	5 vl	2580
26706	✔ **Mycological Peptone BactoBio for bacteriology** Used as a nutritive source for the isolation, cultivation and indentification of saprophytic, dermatophytic fungi-yeast and moulds.	500 Gms 5 Kg	3480 24870
66739	✔ **MYP Agar Base (Phenol Red Egg Yolk Polymyxin Agar Base) (BIS Basis)** For the isolation and identification of pathogenic staphylococci and bacilli.	100 Gms 500 Gms	1350 5040
74997	✔ **MYP Agar Base (Phenol Red Egg Yolk Polymyxin Agar Base) (ISO Basis)** For the enumeration of presumptive Bacillus cereus from food and animal feed as well as from environmental samples.The composition and specification of this medium is as given by ISO 7932 : 2004, ISO 21871 :2006.	100 Gms 500 Gms	1350 5040
67124	🌡 **Nalidixic Selective Supplement (15mg/vl)** Used for the selective isolation of Pseudomonas aeruginosa.	5 vl	1460
77316	🌡 **Neomycin Supplement (30mg/vl)** Used for selective isolation of various microorganisms.	5 vl	1460
44697	✔ **Neopeptone BactoBio for bacteriology** (Peptone Special) Neopeptone, an enzymatic digest of protein is well suited for the growth of many delicate and fastidious bacteria like Group A Streptococci. It also facilitates the development of mucoid colonies of the hemolytic Streptococci.	500 Gms	2920
79989	✔ **Neutral Red Chalk Lactose Agar**	100 Gms 500 Gms	1240 5380

Code	Name	Packing	Price
65466	✔ **NIH Agar**	100 Gms	1120
		500 Gms	4710
23527	✔ **Nitrate Agar**	100 Gms	1680
		500 Gms	7060
	For the differentiation of aerobic and facultative Gram-negative organisms based on their ability to reduce nitrate.		
48541	✔ **Nitrate Broth (ISO, BIS Basis)**	100 Gms	1240
		500 Gms	4710
	For the differentiation of aerobic and facultative gram negative organisms based on their ability to reduce nitrate.The medium is in accordance to ISO 7932 :1993.		
44388	✔ **Nitrofurantoin Broth Base**	100 Gms	1350
		500 Gms	5270
	Used for the enrichment and isolation of Pseudomonas species.		
34350	🌡 **Non-Spore Anaerobic Supplement**	5 vl	1800
	Used for selective isolation of non-sporing anaerobic bacteria.		
31601	🌡 **Novobiocin MKTT Supplement(40mg/vl)**	5 vl	2240
	Recommended for selective enrichment and isolation of Salmonellae.		
18352	🌡 **Novobiocin Supplement(5mg/vl)**	5 vl	1570
		5 x5 vl	4260
	Used to enhance inhibition of Gram-positive organisms.		
63971	✔ **Nutrient Agar**	100 Gms	1010
		500 Gms	3920
	For the cultivation and maintenance of a wide variety of microorganisms.		
18953	✔ **Nutrient Agar (ISO Basis)**	100 Gms	1120
		500 Gms	4150
	For the cultivation and maintenance of a wide variety of microorganisms.		
73478	✔ **Nutrient Agar BioVeg**	100 Gms	1010
		500 Gms	3920
	General purpose (animal origin free) medium for the cultivation of micro organisms.		
18876	✔ **Nutrient Agar 1.5% (ISO Basis)**	100 Gms	1350
		500 Gms	4710
	Used for the cultivation of a wide variety of micro organisms. This medium is in accordance to formulation and specifications given by ISO 21528-1:2017 and ISO 21528-2:2017.		
84955	✔ **Nutrient Agar 1.5%**	100 Gms	1120
		500 Gms	4480
	Used for the cultivation of a wide variety of microorganisms.		
48275	✔ **Nutrient Agar for Oxidase (ISO Basis)**	100 Gms	1350
		500 Gms	4710
	For confirmation of presence of oxidase in microorganisms in water.		
75435	✔ **Nutrient Agar w/ Manganese**	100 Gms	1240
		500 Gms	4600
73153	✔ **Nutrient Agar w/MUG**	100 Gms	2580
		500 Gms	9970
	Used for the detection of Escherichia coli by fluorogenic procedure especially from water and food samples.		
50427	✔ **Nutrient Agar w/ 1% Peptone**	100 Gms	1120
		500 Gms	3920
	A general purpose medium used for the cultivation of a wide variety of microorganisms, by addition of blood or other biological fluids.		
70955	✔ **Nutrient Agar No. 2 (BIS Basis)**	100 Gms	1120
		500 Gms	4370
	Used for the cultivation and maintenance of a wide variety of microorganisms.This medium is in accordance to composition and specification given by BIS ,IS :5887 (PART -I, II, V :1986.		
17137	✔ **Nutrient Agar No. 2**	100 Gms	1120
		500 Gms	4370
	For the cultivation and maintenance of a wide variety of bacteria.		
58561	✔ **Nutrient Agar, pH 6.8**	100 Gms	1240
		500 Gms	4370
	For the cultivation and enumeration of organisms in water, sewage, faeces and other materials.		
96434	✔ **Nutrient Agar, pH 7.0 (ISO Basis)**	100 Gms	1240
		500 Gms	4600
	For the cultivation of a wide variety of microorganisms,including Salmonella species.This medium is in accordance to specifications given by ISO 6579, ISO 10273 and ISO 21528.		
38009	✔ **Nutrient Agar, pH 7.0 (BIS Basis)**	100 Gms	1240
		500 Gms	4600
67948	✔ **Nutrient Agar pH 6 w/ 0.8 % NaCl**	100 Gms	1240
		500 Gms	4480
55427	✔ **Nutrient Broth**	100 Gms	900
		500 Gms	3140
	For the cultivation and maintenance of a wide variety of microorganisms.		
39905	✔ **Nutrient Broth BioVeg**	100 Gms	900
		500 Gms	3140
	An animal origin free media, for the cultivation and maintenance of a wide variety of microorganisms .		

Part D – Dehydrated Culture Media

Code	Name	Packing	Price
97327	✔ **Nutrient Broth No. 2** For the cultivation and maintenance of less fastidious bacteria.	100 Gms 500 Gms	900 3360
63679	✔ **Nutrient Broth w/ 1% Peptone (BIS Basis)** Used as a general purpose culture medium.	100 Gms 500 Gms	900 3360
93063	✔ **Nutrient Broth w/ 1% Peptone (IP, BP Basis)** For microbial limit test and sterility testing.	100 Gms 500 Gms	900 3360
66363	✔ **Nutrient Broth, pH 6.9 w/o NaCl** A general purpose medium for cultivation of a wide variety of microorganisms.	100 Gms 500 Gms	900 3360
94007	✔ **Nutrient Broth w/ 1% Yeast Extract** For the cultivation and maintenance of a wide variety of microorganisms.	100 Gms 500 Gms	1240 4480
22583	✔ **Nutrient Gelatin** For detection of gelatin liquefaction by proteolytic microorganisms	100 Gms 500 Gms	1120 3810
50540	✔ **N Z Broth**	100 Gms 500 Gms	6720 11430
73438	✔ **Oatmeal Agar** Used for the cultivation of fungi,particularly for macrospore formation.	100 Gms 500 Gms	1350 6160
25916 [134134-86-4]	✔ **Oatmeal Powder BactoBio for bacteriology** Oatmeal powder is obtained from oat grains and manufactured under controlled conditions. It is highly nutritious source for growth of fungi in media.	100 Gms 500 Gms	1240 5830
83951	✔ **OF Basal Medium** For use with added carbohydrate for differentiation of gram negative microorganisms based on oxidation-fermentation patterns.	100 Gms 500 Gms	1120 4260
13566	✔ **OGYE Agar Base (Oxytetracycline Glucose Yeast Extract Agar Base)** With added Oxytetracycline used for selective isolation of yeast and molds in foods.	100 Gms 500 Gms	1680 4710
14427	✔ **Orange Serum Agar** Used for the isolation and enumeration of spoilage organisms of citrus products and for examining fruit beverages.	100 Gms 500 Gms	1460 5270
57925	✔ **Orchid Agar**	100 Gms 500 Gms	1460 5830
34235	✔ **Ornithine Decarboxylase Broth** Recommended medium for the differentiation of gram-negative enteric bacilli based on ornithine decarboxylation.	100 Gms	3920
14792	✔ **L-Ornithine Decarboxylase Broth (ISO Basis)**	100 Gms	4040
42167	🌡 **Ottaviani-Agosti Selective Supplement (ISO basis)** Selective supplement for the isolation, enumeration and detection of Listeria species.	5 vl	2690
57253	🌡 **Oxacillin Resistant Selective Supplement** Recommended for screening of oxacillin resistant microorganisms.	5 vl	1680
99455	✔ **Ox Bile (Ox Gall) Powder BactoBio for bacteriology** Dehydrated ox bile, used is preparation of microbiological culture media as a selectively inhibitory agent.It inhibits bacteria not growing in enteric environment and favours growth of enteric microorganisms.	100 Gms 500 Gms	4200 19000
42116	🌡 **Oxford Listeria Supplement** Used for selective isolation of Listeria species.	5 vl	2130
51521	🌡 **Oxytetra Selective Supplement (50mg/vl)** Used for the selective isolation of yeasts and moulds.	5 vl 5 x5 vl	1240 3810
35973	✔ **Pagano Levin Base** (Pagano Levin Candida Test Medium)	100 Gms 500 Gms	1120 4040
39992	✔ **Pantothenate Culture Agar**	100 Gms 500 Gms	2020 8850
50769	✔ **PE 2 Medium**	100 Gms 500 Gms	1460 5040
89858	🌡 **PEMBA Supplement (10KU/vl)** An antibiotic supplement recommended for the selective cultivation of Bacillus cereus.	5 vl	1240
95292 [73049-73-7]	✔ **Peptone BactoBio for bacteriology** Used in various general and diagnostic media, also for large-scale production of enzymes, vaccines, antibiotics, etc. Peptic digest of animal tissue also used in microbiological culture media.	500 Gms 5 Kg 25 Kg	1500 14000 60000

Code	Name	Packing	Price
51960	**Peptone BioVeg BactoBio for bacteriology** Peptone derived from vegetable/plant source; is alternative to traditional animal based peptones.It contains polypeptides and amino acids which provides readily available source of nitrogen and minerals for the growth of bacteria and fungi.	500 Gms 5 Kg 10 Kg	2690 21060 40660
29002	**Peptone BioVeg Sterile for bacteriology (Gamma Irradiated)** Peptone,derived from vegetable source of protein, which is gamma irradiated (Sterile). Recommended for use in manufacturing of vaccines, enzymes, antibiotics,etc.	100 Gms 500 Gms	1010 3810
23136	**Peptone Granulated BactoBio for bacteriology** A peptic digest of animal tissue used in various general and diagnostic media also for large scale production of enzyme, vaccines, antibiotics, etc.	500 Gms 5 Kg	3360 31590
95672	**Peptone Special BactoBio for bacteriology** Peptic digest of animal tissue used in microbiological culture media.	500 Gms 5 Kg 10 Kg	3030 27220 50400
65059	**Peptone Iron Agar** For the cultivation and differentiation of microorganisms based on their ability to produce hydrogen sulfide.	100 Gms 500 Gms	2800 7510
21025	**Peptone Sorbitol Bile Broth** For the selective enrichment of Yersinia enterocolitica	100 Gms 500 Gms	1570 5270
55010	**Peptone Water** For the cultivation of non-fastidious microorganisms, for carbohydrate fermentation tests and also for performing the indole tests.	100 Gms 500 Gms	1010 3360
11522	**Peptone Water BioVeg** General purpose growth medium and as the base of carbohydrate fermentation media.	100 Gms 500 Gms	1120 3920
56562	**Peptone Water w/ Phenol Red** For cultivation of non-fastidious microorganisms and as a base for carbohydrate fermentation tests.	100 Gms 500 Gms	1010 3920
38357	**0.1% Peptone Salt Solution Dehydrated** Recommended as a diluent ,in certain test methods of food,dairy and water samples.	100 Gms 500 Gms	1350 3480
70376	**Peptonized Milk (Peptonized SM Powder) BactoBio for microbiology** Peptonized milk is an enzymatic digest of fat free milk solids.It contains peptones, polypeptides and amino acids derived from milk protein as well as other nutritive components of skim milk. Used for isolation and cultivation of Lactobacilli and Streptococci, especially in dairy products.	100 Gms 500 Gms	1240 3920
17511	**Peptonized Milk Agar (Peptonized SM Agar)**	100 Gms 500 Gms	1570 6160
66916	**Perfringens Agar Base (O.P.S.P.)** For presumptive identification of Clostridium perfringens in foods.	100 Gms 500 Gms	1570 5830
11976	**Perfringens Supplement-I (50mg/vl)** Used for selective isolation of Clostridium perfringens	5 vl 5 x5 vl	1240 3140
62598	**Perfringens Supplement-II** Used for selective isolation of Clostridium perfringens.	5 vl 5 x5 vl	1570 4150
92096	**Perfringens S.F.P. Supplement (S.F.P. Supplement)** Used for the selective isolation of clostridium perfringens.	5 vl 4 x5 vl	1460 4710
55688	**Perfringens T.S.C. Supplement (T.S.C Supplement) (200mg/vl)** Used for selective isolation of Clostridium perfringens.	5 vl 4 x5 vl	1570 4370
23504	**Pikovskaya`s Agar** Used for the detection of phosphate solubilizing soil microorganism.	100 Gms 500 Gms	1350 5270
52133	**Pikovskaya`s Broth** Used for detection of phosphate solubilizing microorganisms.	100 Gms 500 Gms	1240 5160
59048	**Phenol Red Agar Base** Used as a basal medium for the determination of fermentation reactions of microrganisms after addition of various carbohydrates.	100 Gms 500 Gms	1350 4930
14517	**Phenol Red Broth Base** For the study of carbohydrate fermentations.	100 Gms 500 Gms	1460 4150

Code	Name	Packing	Price
79054	✔ **Phenol Red Dextrose Broth**	100 Gms 500 Gms	1240 4710
35850	✔ **Phenol Red Lactose Broth** Used for determination of fermentation of lactose in the differentiation of microorganisms.	100 Gms 500 Gms	1350 4710
15572	✔ **Phenol Red Tartrate Agar** (Jordan Tartrate Agar)	100 Gms 500 Gms	4040 18600
12168	✔ **Phenolphthalein Phosphate Agar** For identification of phosphatase positive Staphylococcus aureus.	100 Gms 500 Gms	1460 5720
92335	✔ **Phenylalanine Agar** (Phenylalanine Deaminase Medium) For the differentiation of enteric gram-negative bacilli on the basis of their ability to produce phenyl pyruvic acid from phenylalanine.	100 Gms 500 Gms	1350 5380
63532	✔ **Phenylalanine Malonate Broth (Shaw and Clarke)** Used as a test medium to detect the utilization of malonate and phenylalanine deamination.	100 Gms 500 Gms	5720 13000
34115	✔ **P.E.A. (Phenylethyl Alcohol) Agar** For the selective isolation of gram positive organisms like Staphylococci and Streptococci.	100 Gms 500 Gms	3030 13330
21418	✔ **Phosphate Buffer pH 7.2 (APHA)** Used for preparation of dilutions and blanks in the examination of specimens like water, dairy products, foods, etc.	100 Gms 500 Gms	1010 2800
69049	✔ **Phosphate Buffer, pH 8.0** Recommended for preparation of dilutions and blanks as per methods specified in pharmacopoeia.	100 Gms 500 Gms	1120 4480
10453	✔ **Plate Count Agar** For the enumeration of viable bacteria in milk and dairy products, heterotrophic bacteria in water by pour plate technique.	100 Gms 500 Gms	1240 5040
52192	✔ **Plate Count Agar BioVeg** For the enumeration of viable bacteria in milk and dairy products, heterotrophic bacteria in water by pour plate technique.	100 Gms 500 Gms	1240 5040
36879	✔ **Plate Count Agar w/MUG**	100 Gms 500 Gms	2580 10310
21699	✔ **Plate Count Agar (ISO Basis)** (Standard Methods Agar) Recommended medium for detecting plate counts of microorganisms in food, water and waste water .The formula conforms to ISO 4833-1&2:2013, ISO 17410:2001, (APHA), (AOAC) and FDA Bacteriological Analytical Manual.	100 Gms 500 Gms	1350 5600
88393	✔ **Plate Count Agar (BIS Basis)** Used for determination of plate counts of microorganisms in milk and dairy products by pour plate technique or spread plate technique.This medium is in accordance to specifications mentioned in BIS, IS 5402:2012.	100 Gms 500 Gms	1350 5270
78609	✔ **Plate Count Agar, Special** For the enumeration of viable bacteria in raw milk and other dairy products as recommended by Netherlands Dairy Association.	100 Gms 500 Gms	2130 8290
45218	✔ **Plate Count Agar w/ BCP**	100 Gms 500 Gms	1460 5270
25496	✔ **PM Indicator Agar (Penicillin in Milk Indicator Agar)** (Penicillin in Milk Indicator Agar)	100 Gms 500 Gms	1570 5380
35980	🌡 **Polymyxin B Selective Supplement (50KU/vl)** Used for the selective isolation of various microorganisms.	5 vl 5 x5 vl	1240 4040
81798	✔ **Polysorbate 80 Agar (Twin Pack)** A nutritional medium containing neutralizing agent, used commonly for cultivating microorganisms.	100 Gms 500 Gms	1240 4600
14644	🌡 **Potassium Lactate 50% (10ml/vl)** Used for isolation and enumeration of yeasts.	5 vl	1800
12370	🌡 **Potassium Sorbate 10% (10 ml/vl)** Recommended for the detection and isolation of acid resistant yeasts.	5 vl	1570
91488	🌡 **Potassium Tellurite 1% Supplement (1ml/vl)** Used for Selective Isolation for Cornybacteria and Staphylococci.	5 vl 5 x5 vl	1120 3480
99482	🌡 **Potassium Tellurite 3.5% Supplement (1ml/vl)** Used for selective isolation of cornybacteria and staphylococci	5 vl	1120

Code	Name	Packing	Price
15344	✔ **Potato Carrot Agar** Used to enhance sporulation of a variety of fungi. Also recommended for sporangial production of Phytophthora species.	100 Gms 500 Gms	1460 6050
17235	✔ **Potato Carrot Broth** Recommended broth medium for the reproduction of Pyronema domesticum.	100 Gms 500 Gms	1350 5830
71788	✔ **Potato Dextrose Agar (PDA)** For the isolation, cultivation and enumeration of yeasts and molds from dairy and food products.Used to induce sporulation in many fungi.	100 Gms 500 Gms 2 Kg 5 Kg	1010 3810 13890 33270
94111	✔ **Potato Dextrose Agar (PDA), Harmonized** Used for the cultivation of yeasts and molds in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1120 3920
91338	✔ **Potato Dextrose Agar (PDA) BioVeg** For the isolation, cultivation and enumeration of yeasts and molds from dairy and food products. Used to induce sporulation in many fungi.	100 Gms 500 Gms	1010 3810
63506	✔ **Potato Dextrose Agar w/ Chloramphenicol (0.1g/lit) (PDA w/ 0.01% Chloramphenicol)** Recommended as a general purpose medium for cultivation and isolation of yeasts and molds from foods, dairy and cosmetics.	100 Gms 500 Gms	1350 5160
72704	✔ **Potato Dextrose Agar w/2% Glucose and 60% Sucrose (PDA w/2% Glucose and 60% Sucrose)** For the isolation of Zygosaccharomyces rouxii from chocolate syrup.	100 Gms 500 Gms	1460 5720
35208	✔ **Potato Dextrose Broth (PDB)** For the isolation, cultivation and enumeration of yeasts and molds from dairy and food products.	100 Gms 500 Gms	1010 3480
74474	✔ **Potato Dextrose Broth (PDB) BioVeg** For the cultivation and enumeration of yeasts and molds.	100 Gms 500 Gms	1240 3920
69190	✔ **Potato Dextrose Rose Bengal Agar** Used for promoting ascospore production.	100 Gms 500 Gms	1120 4930
60040	✔ **Potato Infusion Powder (PI Powder) BactoBio for bacteriology** Potato infusion, is a high quality infusion powder for use as an ingredient in microbiological culture media, for enhanced sporulation and growth of fungal microorganisms.	100 Gms 500 Gms	1120 3250
43921	✔ **Potato Malt Agar** For the cultivation of fungi and other aciduric microorganisms.	100 Gms 500 Gms	1350 4710
74426	✔ **Presence - Absence Broth** (PA Broth) For the detection of coliform bacteria in water from treatment plants or distribution systems using the presence-absence coliform test.	100 Gms 500 Gms	1240 4260
41406 [91079-38-8]	✔ **Proteose Peptone BactoBio for bacteriology** Enzymatic digest of protein, high in proteoses. Used in microbiological culture media and in bacterial toxin production and also in microbiological media for bacteria requiring varying nutritional properties.	500 Gms 2 Kg	3030 9190
93788	✔ **Protose BactoBio for bacteriology** A specially manufactured biological ingredient, very high in proteoses. It`s composition of various amino acids ,along with proteoses and peptones ,make it suitable ingredient for fermentation and vaccine manufacturing process.	100 Gms 500 Gms	1240 3590
90912	✔ **Pseudomonas Agar Base (ISO Basis)** Used for the selective isolation of Pseudomonas species.The composition and specification of this medium is as given by ISO 16266-1:2006	100 Gms 500 Gms	1120 4710
49698	✔ **Pseudomonas Agar for Fluorescein** For the cultivation and observation of fluorescein production in Pseudomonas species.	100 Gms 500 Gms	1240 4370
80358	✔ **Pseudomonas Agar for Pyocyanin** (Pseudomonas Agar-P) For the isolation, cultivation and differentiation of Pseudomonas species on the basis of pyocyanin production.	100 Gms 500 Gms	1120 4040
40976	✔ **Pseudomonas Asparagine Broth (APHA)** Used in the presumptive examination of Pseudomonas aeruginosa from natural & recreational waters.	100 Gms 500 Gms	4150 13330

Code	Name	Packing	Price
55523	✔ **Pseudomonas Isolation Agar Base** For the isolation and cultivation of Pseudomonas species.	100 Gms 500 Gms	1240 5040
36124	✔ **Pseudomonas Isolation Agar BioVeg** For the isolation and cultivation of Pseudomonas species.	100 Gms 500 Gms	1240 5040
47758	✔ **Pseudomonas Test Kit (Pre-Weighed Media)** Used for isolation, identification and cultivation of Pseudomonas species.	1 Kit	4150
69827	✔ **Purple Agar Base**	100 Gms 500 Gms	1570 6390
98775	✔ **Purple Broth Base** Used in the preparation of carbohydrate media especially for identification of pure cultures of members of Enterobacteriaceae, based on fermentation studies	100 Gms 500 Gms	1240 5040
10666	✔ **Pyridoxine Assay Medium** For the microbiological assay of pyridoxine using Neurospora sitophila as the test organism.	100 Gms	4370
29051	✔ **PYR Agar** Used for the isolation and identification of Streptococcus pyogenes.	100 Gms 500 Gms	1680 6160
50846	✔ **R2A Agar** Used for enumerating heterotrophic bacteria from drinking water.	100 Gms 500 Gms	1570 5600
34629	✔ **R2A Agar BioVeg** Used for enumeration of heterotrophic bacteria from drinking water	100 Gms 500 Gms	1460 5830
16437	✔ **R3A Agar** For the cultivation and maintenance of heterotrophic bacteria from potable waters.	100 Gms 500 Gms	1800 5830
39167	✔ **Rappaport Vassiliadis R 10 Broth** Used for the selective enrichment of Salmonella from meat and dairy products, faeces and sewage polluted water.	100 Gms 500 Gms	1240 4710
59912	✔ **Rappaport Vassiliadis Salmonella Enrichment Broth, Harmonized** Used for selective enrichment of Salmonella species in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1120 4260
22140	✔ **Rappaport Vassiliadis Soyabean Meal Broth (RVSM Broth)** Used for enrichment and isolation of Salmonella species.	100 Gms 500 Gms	1460 4930
18530	✔ **Razi`s Medium**	100 Gms 500 Gms	1460 5720
78777	✔ **Reinforced Clostridial Agar** Used for the cultivation and enumeration of Clostridium species and other anaerobes from clinical specimens and foods.	100 Gms 500 Gms	1120 4260
11999	✔ **Reinforced Medium for Clostridia, Harmonized** Used for the cultivation and enumeration of Clostridia in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1010 3480
71467	✔ **Rhizobium Medium** Used for the cultivation of Rhizobium species.	100 Gms 500 Gms	1680 5160
63380	✔ **Ringer Salt Solution Powder** Used as an isotonic diluent during food, milk and dairy products microbiological examination.	100 Gms 500 Gms	2020 8400
57672	🌡 **Rosolic Acid Supplement (0.1gm/vl)** Used for the selective isolation of coliforms.	5 vl	1570
76126	✔ **SA Agar Base** Used for isolation and differentiation of Aeromonas hydrophila from foods based on starch hydrolysis (as per APHA).	100 Gms 500 Gms	2240 6840
99752	✔ **S. Aureus Test Kit (Pre-Weighed Media)** Used for isolation, identification and cultivation of Staphylococcus aureus.	1 Kit	4710
67506	✔ **SABHI Agar Base** Used for isolation and cultivation of pathogenic fungi,especially dermatophytes.	100 Gms 500 Gms	1680 5380
44970	✔ **Sabouraud Chloramphenicol Agar** For the cultivation of yeasts and molds.	100 Gms 500 Gms	1350 5040
55941	✔ **Sabouraud Chloramphenicol Agar BioVeg** For the cultivation of yeasts and molds.	100 Gms 500 Gms	1350 5040
19427	✔ **Sabouraud Dextrose Agar (USP, IP Basis)** (Sabouraud Glucose Agar) For the cultivation of dermatophytes, yeasts and filamentous fungi.	100 Gms 500 Gms	1010 3920

Code	Name	Packing	Price
28055	✔ **Sabouraud Dextrose Agar** Used for the cultivation of dermatophytes, yeasts & filamentous fungi.	100 Gms 500 Gms	1010 3810
95430	✔ **Sabouraud Dextrose Agar Cycloheximide and Chloramphenicol (Mycobiotic Agar)** For the selective isolation and cultivation of pathogenic fungi.	50 Gms 100 Gms	5040 9190
24613	✔ **Sabouraud Dextrose Agar, Harmonized** For the cultivation and isolation of yeasts, molds and aciduric microorganisms in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1010 3920
88999	✔ **Sabouraud Dextrose Agar BioVeg** Used for the cultivation of dermatophytes, yeasts and filamentous fungi.	100 Gms 500 Gms	1010 3920
24835	✔ **Sabouraud Dextrose Broth, Harmonized** For the cultivation and isolation of yeasts, molds and aciduric microorganisms in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	790 2360
64402	✔ **Sabouraud Dextrose Maltose Broth** Used for the cultivation of moulds, yeasts and aciduric organisms as well as testing of antimycotic substances.	100 Gms 500 Gms	1460 5600
14506	✔ **Sabouraud Glucose Broth** (Sabouraud Dextrose Broth) For the cultivation of dermatophytes, yeasts and filamentous fungi.	100 Gms 500 Gms	680 1800
57791	✔ **Sabouraud Maltose Agar** For the cultivation and maintenance of yeasts, molds and a variety of fungi.	100 Gms 500 Gms	1460 5380
33387	✔ **Sabouraud Maltose Broth** Used for the cultivation of a wide variety of fungi.	100 Gms 500 Gms	1240 4930
36370	✔ **Saline Agar** Used for detection of alpha -toxin in Clostridium perfringens.	100 Gms 500 Gms	1460 5040
39930	✔ **SS Agar (Salmonella Shigella Agar)** For the selective isolation and differentiation of pathogenic enteric bacilli especially Salmonella and Shigella species.	100 Gms 500 Gms	1120 4260
36856	✔ **SS Agar, Modified (Salmonella-Shigella Agar, Modified)** (Salmonella Shigella Agar) A modified medium, which provides better growth of Salmonella and Shigella species.	100 Gms 500 Gms	1350 5040
55521	✔ **Salmonella Test Kit (Pre-Weighed Media)** Used for isolation, identification and cultivation of Salmonella species.	1 Kit	2690
48589	✔ **Salt Polymyxin Broth Base (ISO Basis)** (Polymyxin Salt Broth)	100 Gms 500 Gms	1240 4600
22887	✔ **Salt Polymyxin Broth Base (BIS Basis)**	100 Gms 500 Gms	1240 4600
58249	✔ **Schaedler Agar w/o Thioglycollate** (Schaedler Anaerobic Agar w/o Thioglycollate)	100 Gms 500 Gms	2020 5720
16164	✔ **Schaedler Agar** (Schaedler Anaerobic Agar) Used for the enumeration of various aerobic and anaerobic bacterial species present in the gastrointestinal tract.	100 Gms 500 Gms	1350 5160
94766	✔ **Schuberts Arginine Broth** Used as an enrichment medium for testing swimming pool water for isolation of chlorine damaged Pseudomonas aeruginosa.	100 Gms 500 Gms	2800 9500
79432	✔ **SDS Agar Base (ISO basis)** Recommended for isolation and enumeration of Vibrio vulnificus from seafood samples in accordance with APHA, ISO /TS 21872-2 and ISO 11133 : 2014.	100 Gms 500 Gms	1460 5720
87213	✔ **Selenite Broth w/ Dulcitol (Twin Pack)** Used for enrichment of Salmonella species.	500 Gms	8520
16709	✔ **Selenite Cystine Broth Base w/o Selenite**	100 Gms 500 Gms	1010 3920
38926	✔ **Selenite F Broth Base w/ Selenite (Twin Pack)** With addition of selenite, used for isolation and cultivation of Salmonellae species from faeces, and other pathological specimens	100 Gms 500 Gms	1350 5270
83179	✔ **Selenite F Broth Base w/o Selenite** With addition of selenite, used for isolation and cultivation of Salmonella species from faeces, dairy products and other specimens.	100 Gms 500 Gms	1010 3920
17458	✔ **Semi-Solid Nutrient Agar (ISO Basis)**	100 Gms 500 Gms	1350 5270

Code	Name	Packing	Price
75917	✔ **Semi-Solid Nutrient Agar (BIS Basis)**	100 Gms 500 Gms	1350 5270
37963	✔ **Semi-solid Rappaport Vassiliadis (MSRV) Medium Base** A semi-solid medium used for the detection of motile Salmonella species from food and environmental samples	100 Gms 500 Gms	1460 5160
84809	✔ **Sensitest Agar** For the performance of antibiotic sensitivity test.	100 Gms 500 Gms	2580 10080
11580	✔ **S.F.P Agar Base** (Shahidi Ferguson Perfringens Agar Base) Used for the isolation and enumeration of Clostridium perfringens from foods.	100 Gms 500 Gms	1350 6050
94074	✔ **Shigella Broth Base** Used for the isolation and cultivation of Shigella species from food.	100 Gms 500 Gms	1350 4710
45290	🌡 **Shigella Selective Supplement (0.55mg/vl)** Used for the selective isolation and cultivation of Shigella species.	5 vl	1460
91694	✔ **Slanetz and Bartley Medium** For the detection and enumeration of enterococci by membrane filtration method.	100 Gms 500 Gms	1350 5160
49810	✔ **SIM Medium (Sulfite Indole Motility Medium)** For the differentiation of members of the Enterobacteriaceae, based on H2S production, indole production and motility.	100 Gms 500 Gms	1120 4150
57309	✔ **SIM Medium BioVeg** Used for the differentiation of members of the Enterobacteriaceae, based onH2S production, indole production and motility.	100 Gms 500 Gms	1120 4150
59437	✔ **Simmon`s Citrate Agar** For the differentiation of gram-negative bacteria on the basis of citrate utilization.	100 Gms 500 Gms	800 2700
28582	✔ **Skim Milk Powder (SM Powder) BactoBio for bacteriology** Specially made from Fat free milk,high in protein.Suitable as ingredient in microbiology media.	500 Gms 1 Kg	1050 1800
69372	✔ **Skim Milk Plate Count Agar (SM Plate Count Agar)** For the enumeration of bacteria in milk and dairy products.	100 Gms 500 Gms	2920 6840
41023	✔ **Sorbitol Iron Agar**	100 Gms 500 Gms	1680 6050
77892	✔ **Soyabean Casein Digest Agar (SCDA) (Tryptone Soya Agar, TSA), Harmonized** For isolation and cultivation of fastidious and non-fastidious microorganisms in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms 2 Kg	1010 3920 14120
33291	✔ **Soyabean Casein Digest Agar (SCDA) (Tryptone Soya Agar, TSA)** Used for the isolation and cultivation of fastidious & non fastidious microorganisms.	100 Gms 500 Gms 2 Kg	1010 3920 14120
10935	✔ **Soyabean Casein Digest Medium (SCDM) (Tryptone Soya Broth, TSB), Harmonized** Used for cultivation of a wide variety of fastidious and non-fastidious microorganisms in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms 2 Kg	790 2130 6950
47404	✔ **Soyabean Casein Digest Medium (SCDM) (Tryptone Soya Broth, TSB), Sterile (Gamma Irradiated)** Gamma irradiated sterile powder used for the growth of wide variety of bacteria and fungi.	500 Gms	2690
51993	✔ **Soyabean Casein Digest Medium (SCDM) (Tryptone Soya Broth, TSB) BioVeg** (Soyabean Medium BioVeg) For the cultivation of a wide variety of microorganisms.	100 Gms 500 Gms	790 2360
44520	✔ **Soyabean Casein Digest Medium (SCDM) (Tryptone Soya Broth, TSB), BioVeg Sterile (Gamma Irradiated)** A nutritious general purpose medium, prepared by replacing animal based peptones with vegetable based peptones. It is gamma irradiated sterile powder, used in evaluation of sterility in manufacturing process.	500 Gms	3250
28501	✔ **Soyabean Casein Digest Medium (SCDM) (Tryptone Soya Broth, TSB)** A nutritious general purpose medium used for the growth of a wide variety of bacteria and fungi.	100 Gms 500 Gms 2 Kg	790 2130 6950
89320 [8002-43-5]	✔ **Soya Lecithin (30%) BactoBio for bacteriology** A natural phospholipid made from non-GMO soya beans.	100 Gms 500 Gms	2240 4600

Code	Name	Packing	Price
20626 [91079-46-8]	**Soyatone BactoBio for bacteriology** (Soya peptone) Papaic digest of soyabean meal phytone, a peptone of plant origin. Contains high concentration of carbohydrates and nitrogen constituents, hence making it suitable for growth of fastidious microorganisms.	500 Gms 5 Kg	1680 14120
24205	**Special YM Medium**	100 Gms 500 Gms	1350 4930
39778	**Spirit Blue Agar** For the detection, enumeration and study of lipolytic microorganisms.	100 Gms 500 Gms	1570 6280
54880	**Sporulation Broth** Recommended medium for promoting sporulation in Bacillus subtilis.	100 Gms 500 Gms	1240 4480
31363	**SPS Agar** For the detection and enumeration of Clostridium perfringens in foods.	100 Gms 500 Gms	1460 5040
14041	**Standard Methods Agar w /Tween 80 and Lecithin** For the determination of sterility of surfaces	100 Gms 500 Gms	3140 12320
92654	**Standard Nutrient Agar No.1** Used for the cultivation of fastidious bacteria and also for the isolation and enumeration of bacteria.	100 Gms 500 Gms	1350 5270
72899	**Standard Nutrient Agar No.2** Standard Nutrient Agar No. 2 is recommended for the cultivation and enrichment of less fastidious bacteria. It can be used in detection of inhibitors in bacteriological examination of meat.	100 Gms 500 Gms	1570 6500
68818	**Standard Nutrient Broth (H.S. Vaccine Medium)** Recommended nutritive medium for mass cultivation of bacteria for vaccine production.	100 Gms 500 Gms	1120 3140
18247	**Staph-Strepto Selective Supplement** Used for the selective isolation of Staphylococcus and Streptococcus.	5 vl 5 x5 vl	2020 4710
39726	**Staphylococcus Medium 110** (Stone Gelatin Agar, Gelatin Mannitol Salt Agar) For the isolation and enumeration of Staphylococci based on mannitol fermentation, pigment formation and gelatinase activity.	100 Gms 500 Gms	1240 4150
37547	**Starch Agar** For detection of starch hydrolyzing microorganisms.	100 Gms 500 Gms	1350 5270
78486 [9005-84-9]	**Starch Soluble Special (Cold Water Soluble) BactoBio for microbiology** A special grade of Starch, suitable for media preparations in microbiology since it is cold water soluble due to its special preparation method. The solution characteristics being clear, make it especially suitable for such use in labs.	100 Gms 500 Gms	2130 4600
76043	**Sterility Test Kit (Pre-Weighed Media)** Used for sterility testing of products as per pharmacoepial standards.	1 Kit	1350
11475	**Strepto Supplement** Used for the selection of Streptococcus species.	5 vl 5 x5 vl	1570 4930
35140	**Streptococcus Selective Supplement** Used for the selective isolation and cultivation of Streptococcus species.	5 vl 5 x5 vl	2130 5040
31025	**Streptococcus Thermophilus Isolation Agar** For determining the ratio of Streptococcus thermophilus and Lactobacillus bulgaricus in yoghurt.	100 Gms 500 Gms	1570 5380
77304	**Streptomyces Agar** For the cultivation and maintenance of Streptomyces.	100 Gms 500 Gms	1350 5160
43048	**Stuart Transport Medium** For the preservation of fastidious organisms and Neisseria species during their transportation from clinic to laboratory.	100 Gms 500 Gms	1680 6840
14840	**Sucrose Salicin Agar (Gillies Agar No.2)**	100 Gms	6950
69237	**Sugar Free Agar (Count Agar Sugar-free)** For the enumeration of psychotropic and mesophilic gram- negative rods in butter and other dairy products as per formula described by the International Dairy Federation.	100 Gms 500 Gms	1460 5040
49778	**Sulphite Agar** (Sulfite Agar)	100 Gms 500 Gms	1800 5600
84524	**Iron Sulphite Agar** (Sulfite Iron Agar) Used for the detection of Clostridia in the examination of meat and meat products.	100 Gms 500 Gms	1240 5040

Part D – Dehydrated Culture Media

Code	Name	Packing	Price
34050	**Sulpha Supplement** Used for the selective isolation of Salmonella species.	5 vl	2130
68041	**Syncase Broth**	100 Gms 500 Gms	1460 5720
23635	**T1N1 Agar (Tryptone Salt Agar)** (Tryptone Salt Agar)	100 Gms 500 Gms	1460 5380
65235	**T1N1 Broth (Tryptone Salt Broth)**	100 Gms 500 Gms	900 3480
15600	**T.A.T. Broth Base** (Tryptone-Azolectin-Tween Broth Base) Used for the cultivation of microorganisms from highly viscous or gelatinous materials.	100 Gms 500 Gms	1120 4930
86626	**TCBS Agar** (Thiosulphate Citrate Bile Salts Sucrose Agar) For the selective isolation of Vibrio cholerae and Vibrio parahaemolyticus from a variety of clinical specimens.	100 Gms 500 Gms	1240 4930
50165	**TCBS Agar (BIS Basis)** (Thiosulfate Citrate Bile Salts Sucrose Agar) Used for selective isolation of Enteropathogenic Vibrio species.	100 Gms 500 Gms	1240 4930
68331	**TCBS Agar BioVeg** For the selective isolation of Vibrio cholerae and Vibrio parahaemolyticus from clinical specimens.	100 Gms 500 Gms	1350 5380
23746	**Tellurite Blood Agar Base** Used for the selective isolation and cultivation of Corynebacterium species.	100 Gms 500 Gms	1460 5720
15297	**Tellurite Cefixime Supplement** Used for the selective isolation and cultivation of Escherichia coli 0157:H7.	5 vl	1680
30431	**Tellurite Glycine Agar Base**	100 Gms 500 Gms	1240 5830
97630	**Tergitol- 7 Agar Base (BIS Basis)** (T7 Agar) Used for selective enumeration and identification of coilform organisms	100 Gms 500 Gms	1350 6050
67662	**Tergitol-7 Agar H Base** (T7 Agar H) For the isolation and differentiation of enteric bacteria from urine specimens.	100 Gms 500 Gms	1680 6720
99371	**Tergitol-7 Broth** (T7 Broth) A selective and differential medium for detection and enumeration of coliforms	100 Gms 500 Gms	1240 4930
30086	**Terrific Broth (Tartoff-Hobbs Broth)** For the cultivation and maintenance of recombinant strains of Escherichia coli.	100 Gms 500 Gms	1120 3920
40873	**Tetracycline Selective Supplement** Used for the selective isolation and cultivation of yeasts and moulds.	5 vl	1800
13112	**Tetrathionate Broth Base w/o Iodine & BG** (TT Broth Base) For the enrichment of Salmonella species.	100 Gms 500 Gms	900 3030
81902	**Tetrathionate Broth Base (ISO Basis)** (TT Broth Base) Recommended for enrichment of Salmonellae from samples suspected to be contaminated with Salmonellae in accordance with IP.	100 Gms 500 Gms	1010 3030
13727	**TGB Agar (Tryptone Glucose Beef Extract Agar)** Recommended medium for enumeration of bacteria in water, milk and other dairy products.	100 Gms 500 Gms	1240 5160
44521	**Thayer-Martin Medium Base** (Thayer-Martin Agar, Modified) For the isolation of Neisseria species from specimens containing mixed flora of bacteria and fungi.	100 Gms 500 Gms	1120 4370
86908	**Thiobacillus Agar**	100 Gms 500 Gms	2240 6950
65489	**Thioglycollate Agar** Used for the growth and cultivation of anaerobic microorganisms.	100 Gms 500 Gms	1350 5160
38136	**Thioglycollate Broth w/ Liver Extract (B.Q. Vaccine Medium)** Recommended medium for cultivation of anaerobes on large scale, especially for vaccine manufacturing.	100 Gms 500 Gms	1120 2920
91690	**Thioglycollate Medium w/o Indicator (Diagnostic Thioglycollate Medium)** For the cultivation of anaerobic, microaerophilic and aerobic microorganisms. For use in sterility testing of a variety of specimens.	100 Gms 500 Gms	1120 3360
63270	**Thiol Broth**	100 Gms 500 Gms	1570 4480
67819	**Ticarcillin Supplement (1mg/vl)** Used for enrichment of Yersinia enterocolitica.	5 vl	1570
29536	**TN Agar**	100 Gms 500 Gms	1460 5380

Code	Name	Packing	Price
11483	✔ **Todd Hewitt Broth** For the cultivation of Group A Streptococci used in serological typing and for the cultivation of variety of pathogenic microorganisms.	100 Gms 500 Gms	1120 4370
95197	✔ **Toluidine Blue DNA Agar** (Dnase Test Agar with Toluidine Blue)	100 Gms 500 Gms	8850 42120
74398	✔ **Tomato Juice Agar** Used for the cultivation of lactobacilli.	100 Gms 500 Gms	2130 5940
57639	✔ **TPEY Agar Base** Used for the selective isolation and enumeration of Staphylococci from food samples.	100 Gms 500 Gms	1460 5720
33544	✔ **Transport Charcoal Medium** For transportation of clinical specimens.	100 Gms 500 Gms	1460 5940
87704	🌡 **Tributyrin Supplement (10ml/vl)** Used for the detection of lipolytic microorganisms.	5 vl 5 x5 vl	2130 8180
30842	✔ **Tributyrin Agar Base w/o Tributyrin** For the cultivation and enumeration of lipolytic microorganisms.	100 Gms 500 Gms	1460 5830
12081	✔ **Triple Sugar Iron Agar (ISO Basis)** (TSI Agar) Used for the differentiation of members of Enterobacteriaceae on the basis of fermenting of glucose, lactose, sucrose and H2S production.This medium is in accordance to specifications given by ISO 6579 and ISO 19250.	100 Gms 500 Gms	800 2500
55720	✔ **Triple Sugar Iron Agar (BIS Basis)** Used for the differentiation of Enterobacteriaceae on the basis of fermenting of glucose, lactose, sucrose and H2S production. This medium is in accordance to specifications given by BIS ; IS: 5887 (Part V)-1976.	100 Gms 500 Gms	800 2500
69309	✔ **Triple Sugar Iron Agar BioVeg** Used for the differentiation of members of Enterobacteriaceae on the basis of fermenting of glucose, lactose, sucrose and H2S production.	100 Gms 500 Gms	900 3000
80175	✔ **Tryptone Agar BioVeg** Used for the detection and enumeration of Escherichia coli in foods using a modified direct plating method	100 Gms 500 Gms	1680 5830
32007 [91079-40-2]	✔ **Tryptone BactoBio for bacteriology** (Casitose) Pancreatic Digest Of Casein used for cultivation of fastidious and non-fastidious bacteria and fungi.	500 Gms 5 Kg	2580 22180
10396	✔ **Tryptone Bile Agar (ISO Basis)** Used for rapid detection and enumeration of Escherichia coli in foods.	100 Gms 500 Gms	1350 5380
62600	✔ **Tryptone Broth (Tryptone Water)** For the cultivation of production of Indole by microorganisms.	100 Gms 500 Gms	900 3480
11947	✔ **Tryptone Broth BioVeg (Tryptone Water BioVeg)** For the cultivation of production of indole by microorganisms.	100 Gms 500 Gms	1010 3480
62087	✔ **Tryptone Glucose Yeast Extract Broth** Recommended for the enumeration of microorganisms from food samples by MPN method.	100 Gms 500 Gms	1350 4370
45227	✔ **Tryptone Phosphate Broth** For the cultivation of enteropathogenic Escherichia coli.	100 Gms 500 Gms	900 3360
73825	✔ **Tryptone Salt Broth** Used as diluent to prepare stock solutions of dairy and other food products, cosmetic products for microbiological analysis.	100 Gms 500 Gms	900 3480
53226	✔ **Tryptone Soya Agar (TSA) w/Lecithin and Polysorbate 80 (SCDA Lecithin & Tween 80)** (Soyabean Casein Digest Agar w/Lecithin & Tween 80, SCDA Lecithin & Tween 80) Tryptone Soya Agar w/ Lecithin and Polysorbate 80 is used for validation of cleanliness on surfaces of containers, equipment surfaces and water miscible cosmetics.	100 Gms 500 Gms	1350 5600
18039	✔ **Tryptone Soya Broth (TSB) w/ 6.5% Sodium Chloride (Soyabean Casein Digest Medium (SCDM) w/ NaCl)** Used for differentiation between Salt tolerant and salt -intolerant species.	100 Gms 500 Gms	1570 3250
21091	✔ **Tryptone Sucrose Tetrazolium Agar Base (ISO Basis)** (TSTA) Recommended medium with added TTC for isolation of Vibrio species.	100 Gms 500 Gms	1350 5270

Code	Name	Packing	Price
20407	✔ **Tryptone Soya Agar (TSA) (Soyabean Casein Digest Agar, SCDA)** For the cultivation and maintenance of a wide variety of microorganisms.	100 Gms 500 Gms 2 Kg	1010 3920 14120
32068	✔ **Tryptone Soya Agar (TSA) (Soyabean Casein Digest Agar, SCDA) BioVeg** Vegetable based (Non-animal) medium,used for cultivation of wide variety of microorganisms from various samples with or without addition of blood.	100 Gms 500 Gms	1240 4370
21794	✔ **Tryptone Soya Agar (TSA) w/MUG (Soyabean Casein Digest Agar w/MUG, SCDA MUG)** Used for the cultivation of fastidious and non fastidious microorganisms by fluorogenic procedure.	100 Gms 500 Gms	1680 7280
24392	✔ **Tryptone Soya Broth (TSB) (Soyabean Casein Digest Medium, SCDM)** For the cultivation of a variety of microorganisms.	100 Gms 500 Gms 2 Kg	790 2130 6950
32347	✔ **Tryptone Soya Broth (TSB) /SPS (Soyabean Casein Digest Medium w/SPS, SCDM w/SPS)** For the preparation of blood culture and cultivation of fastidious and non-fastidious microorganisms such as H. influenzae, N. meningitidis, S. pneumonia,etc..	100 Gms 500 Gms	1120 4150
16684	✔ **Tryptone Soya Broth (TSB) (Soyabean Casein Digest Medium, SCDM) BioVeg** For the cultivation of a variety of microorganisms.	100 Gms 500 Gms	790 2360
45511	✔ **Tryptone Soya Yeast Extract Broth** A highly nutritious non-selective medium for the cultivation of wide variety of microorganisms especially Listeria monocytogenes.	100 Gms 500 Gms	1240 4150
83009	✔ **Tryptone Tellurite Agar Base**	100 Gms 500 Gms	1350 5720
18528 [91079-40-2]	✔ **Tryptone Type I (Casitose) BactoBio for bacteriology** (Casein Enzymatic Hydrolysate (Type I)) Tryptone Type I is an enzymatic digest of milk protein, rich in amino nitrogen source; used in microbiological and fermentation media.	500 Gms	3480
20483	✔ **Tryptone Type II (Casitose) BactoBio for bacteriology** (Casein Enzymatic Hydrolysate (Type II)) Tryptone Type II is an enzymatic digest of lactose free milk protein, rich in amino nitrogen source used especially in Antibiotic assay medium.	500 Gms	3480
69165	✔ **Tryptone Water Broth w/BCP** Recommended for the cultivation of Salmonella species from food samples.	100 Gms 500 Gms	1240 3920
79223	✔ **Tryptone Yeast Extract Agar w/BCP (ISO, BIS Basis)** Recommended medium for isolation and enumeration of Enterobacteriaceae and Bacillus cereus.The specification and composition of this medium is as given by ISO 7402 :1993.	100 Gms 500 Gms	1460 5940
15965	✔ **Tryptophan Broth (Tryptophan Medium) (ISO Basis)** Used for differentiation of bacteria on the basis of indole production.The formulation of this medium is in accordance to specifications mentioned in ISO/DIS 16654.	100 Gms 500 Gms	1010 3360
40626	✔ **Tryptophan Broth (Tryptophan Medium) BioVeg** Used for the differentiation of bacteria on the basis of indole production.	100 Gms 500 Gms	1120 3810
47253	✔ **Tryptose Blood Agar Base** For the cultivation and maintenance of a wide variety of fastidious microorganisms.	100 Gms 500 Gms	1570 6050
12514 [84843-69-6]	✔ **Tryptose BactoBio for bacteriology** Mixed enzymatic hydrolysate of protein, especially used in media to grow fastidious microorganisms.	500 Gms	2690
82302	✔ **Tryptose Cycloserine Dextrose Agar Base** Used for isolation of mesophilic spore forming anaerobes in food spoilage.	100 Gms 500 Gms	1570 6050
48200	✔ **Tryptose Phosphate Broth** For the cultivation of variety of fastidious bacteria.	100 Gms 500 Gms	1010 3810
45989	✔ **Tryptose Phosphate Broth BioVeg** For the cultivation of a variety of fastidious bacteria.	100 Gms 500 Gms	1120 3920

Code	Name	Packing	Price
57498	✔ **Tryptose Sulfite Cycloserine Agar (TSC Agar) (ISO basis)** Used for the detection and enumeration of Clostridium perfringens from food samples.It is in accordance to guidelines mentioned in ISO 7937 ,ISO 11133 and ISO 14189.	100 Gms 500 Gms	1180 5600
55052	🌡 **TTC Solution 1% (10ml/vl)** Used for detecting microorganisms on the basis of TTC reduction.	5 vl 5 x5 vl	1460 4260
18339	✔ **Urea Agar Base (BIS, ISO Basis)** Used for the detection of urease production by Proteus species and for identification of other members of Enterobacteriaceae.	100 Gms 500 Gms	1100 3300
88708	✔ **Urea Agar Base** Used for the detection of urease production by Proteus species and for identification of other members of Enterobacteriaceae	100 Gms 500 Gms	900 3000
42465	✔ **Urea Agar Base (Christensen) BioVeg** Used for the detection of urease production by Proteus species and for identification of other members of Enterobacteriaceae.	100 Gms 500 Gms	1100 3300
97438	✔ **Urea Broth Base** For identifying bacteria on the basis of urea utilization, for the differentiation of Proteus species from Salmonella and Shigella species.	100 Gms 500 Gms	1120 4370
58028	✔ **Urea Indole Medium** Used for the differentiation of Enterobacteria on the basis of urease and indole production and the transdeamination of tryptophan.	100 Gms 500 Gms	1350 5160
80192	🌡 **40% Urea Solution (5ml/vl)** Used for the detection of urease activity.	5 vl 5 x5 vl	790 2690
26117	🌡 **Vancomycin Supplement (3mg/vl)** Used for the selective isolation of Enterococci.	5 vl	2130
34123	🌡 **VCN Supplement** Used for selective isolation of Neisseria species.	5 vl 5 x5 vl	3030 10310
32249	🌡 **VCNT Supplement** Used for selective isolation of Neisseria species.	5 vl 5 x5 vl	3030 10310
78604	✔ **V Infusion Agar (Veal Infusion Agar)**	100 Gms 500 Gms	1680 5040
86771	✔ **V Infusion Broth (Veal Infusion Broth)** Recommended for cultivation of fastidious pathogenic bacteria.	100 Gms 500 Gms	1350 4600
67389	✔ **Veg Extract Agar BioVeg**	500 Gms	5270
55605	✔ **Veg Extract Broth BioVeg**	500 Gms	5270
53818	✔ **Violet Red Bile Agar (BIS Basis)** Recommended medium for selective isolation, detection and enumeration of coliforms in water, milk and other dairy products.This medium is in compliance with BIS specification IS:5401-1969.	100 Gms 500 Gms	1120 4260
57081	✔ **Violet Red Bile Agar (1.2%) (ISO Basis)** Recommended medium for selective isolation, detection and enumeration of coliforms in water, milk and other dairy products.This medium is in accordance to specifications given by ISO 5541-1, ISO 4832.	100 Gms 500 Gms	1350 4370
33690	✔ **Violet Red Bile Agar w/ MUG** For the detection of coliforms and Escherichia coli in food and dairy products by fluorogenic procedure.	100 Gms 500 Gms	2020 8520
88966	✔ **Violet Red Bile Broth** For detection and enumeration of coliforms from water and foods.	100 Gms 500 Gms	1120 3920
43459	✔ **Violet Red Bile Dextrose Agar BioVeg** Used for detection and enumeration of Enterobacteriaceae from foods.	100 Gms 500 Gms	1460 4370
55373	✔ **Violet Red Bile Glucose Agar (ISO Basis)** (Violet Red Bile Glucose Agar w/o Lactose) Used for detection and enumeration of Enterobacteriaceae from foods.The formulation and performance criteria of this medium is in accordance to ISO 21528 :2017 and ISO 11133 :2014.	100 Gms 500 Gms	1120 4150
89912	✔ **Violet Red Bile Glucose Agar, Harmonized** Used for cultivation and isolation of Enterobacteriaceae in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1120 4040

Code	Name	Packing	Price
34326	**Vitamin K1 Supplement (10mg/vl)** A vitamin growth supplement used for the isolation of anaerobic organisms.	5 vl	1460
38907	✔ **Vogel and Johnson Agar Base (USP, IP Basis)** (V.J. Agar, Base) For the detection of coagulase positive Staphylococcus aureus.	100 Gms 500 Gms	1010 3810
56055	✔ **Water Agar** Used for cultivation and observation of sporulation of some fungi and also for the enumeration of fungi according to MPN method.	100 Gms 500 Gms	1460 4710
23892	✔ **Water (Faecal coliforms) Testing Kit (Pre-Weighed Media)** For the enumeration and identification of faecal coliforms from water and waste water samples.	1 Kit	2240
19990	✔ **Water Plate Count Agar (ISO basis)** Used for the enumeration of culturable microorganisms by colony count technique from water samples .It is a non-selective medium;which conforms to ISO 6222:1999.	100 Gms 500 Gms	1240 5940
55014	✔ **Wilkins-Chalgren Agar** (Anaerobic Agar) For cultivation and maintenance of anaerobic bacteria. In susceptibility testing to determine minimum inhibitory concentrations (MICs) of antibiotics for anaerobic bacteria.	100 Gms 500 Gms	1460 5270
91753	✔ **Willis and Hobbs Medium Base** Used for isolation and identification of Clostridium from food samples.	100 Gms 500 Gms	1460 5720
31560	✔ **WL Differential Agar** Used in isolating bacteria encountered in industrial fermentation processes.	100 Gms 500 Gms	1350 4710
45086	✔ **WL Nutrient Agar** (Wallerstein Laboratory Medium) For the detection of bacteria & yeasts in industrial fermentation processes, particularly from beer processing	100 Gms 500 Gms	1460 4930
78176	✔ **WL Nutrient Broth** (Wallerstein Laboratory Nutrient Broth) For cultivation of bacteria and yeasts in industrial fermentation processes, particularly from beer processing.	100 Gms 500 Gms	1460 4260
47617	✔ **Wort Agar** For the cultivation and enumeration of yeasts. Its low pH inhibits bacterial growth.	100 Gms 500 Gms	1240 5270
26434	✔ **Wurtz Medium** Used for isolation and differentiation of lactose fermenting bacteria.	100 Gms 500 Gms	1350 5040
53556	✔ **Xylose Lysine Deoxycholate Agar (XLD Agar) (USP, IP Basis)** Used for the isolation and differentiation of enteric pathogens, especially Shigella and Providencia.	100 Gms 500 Gms	1120 4480
70791	✔ **Xylose Lysine Deoxycholate Agar (XLD Agar)** Used for the isolation and identification of Salmonellae and Shigellae.	100 Gms 500 Gms	1120 4480
89660	✔ **XLT 4 Agar Base** Used for isolation and identification of Salmonellae from various samples.	100 Gms 500 Gms	1570 6500
73111	✔ **Xylose Lysine Deoxycholate Agar (XLD Agar), Harmonized** For cultivation and isolation of Salmonella species in pharmaceutical testing according to harmonized methods.	100 Gms 500 Gms	1350 5160
67146 [84604-16-0]	✔ **Yeast Autolysate (Yeast Peptone) BactoBio for microbiology** Water soluble portion of autolysed yeast especially rich in B-complex vitamins, used in microbiological and fermentation media.	100 Gms 500 Gms	1010 3360
22176	**Yeast Autolysate Supplement** Used for enrichment of Neisseria species.	5 vl 5 x5 vl	1460 4600
51089	✔ **Yeast and Mould Agar** Used for isolation of yeasts, molds & other aciduric microorganisms.	100 Gms 500 Gms	1460 5600
34266 [8013-01-02]	✔ **Yeast Extract BactoBio for bacteriology** Water soluble portion of autolysed yeast especially rich in B-complex vitamins used in microbiological and fermentation media.	100 Gms 500 Gms 5 Kg 10 Kg	1120 3030 28680 54880
89463 [8013-01-02]	✔ **Yeast Extract for bacteriology** Low pH, slightly turbid, extract of autolysed yeast cells used generally for technical purposes.	100 Gms 500 Gms 5 Kg	900 1460 13560

Code	Name	Packing	Price
52184	✔ **Yeast Extract for cell culture** Water soluble portion of autolysed yeast, suitable for cell culture applications. Many media containing Yeast extract have been recommended for use in fermentations, cell culture media, molecular genetic media, and culture media for cultivation and resuscitation of a large variety of fastidious microorganisms and antibiotics assays.	100 Gms 500 Gms	1460 3810
36504	✔ **Yeast Extract Agar** For the enumeration of microorganisms in potable and fresh water samples.	100 Gms 500 Gms	2020 6610
62332	✔ **Yeast Lactose Agar** For cultivation of soil microorganisms like Rhizobium species.	100 Gms 500 Gms	2690 10640
49723	✔ **Yeast Malt Agar** (YM Agar) For the cultivation of fungi including yeasts and other aciduric organisms.	100 Gms 500 Gms	1350 5380
63973	✔ **Yeast Malt Extract Broth** (YM Broth) For the cultivation of fungi, yeasts and other aciduric microorganisms.	100 Gms 500 Gms	1120 4260
40965	✔ **Yeast Mannitol Agar** For the isolation and enumeration of soil microorganisms like Rhizobium species.	100 Gms 500 Gms	1460 6500
56342	✔ **Yeast Mannitol Broth** For the cultivation of soil microorganisms like Rhizobium species.	100 Gms 500 Gms	1120 4480
54250	✔ **YPD Agar (Yeast Peptone Dextrose Agar, YEPD Agar)** For maintaining and propagating yeasts in molecular biology studies.	100 Gms 500 Gms	1240 4930
51830	✔ **Yersinia Isolation Agar** For the isolation of Yersinia species from foods.	100 Gms 500 Gms	1570 5600
27810	✔ **Yersinia Selective Agar Base** For the isolation of Yersinia enterocolitica from clinical specimens and food samples.	100 Gms 500 Gms	1570 5600
15061	🌡 **Yersinia Selective (C.I.N.) Supplement** Used for isolation of Yersinia enterocolitica.	5 vl 5 x5 vl	1800 5270

Code	Name	Packing	Price
83385	✔ **Yeast Glucose Chloramphenicol Agar (YGC Agar, Chloramphenicol Agar)** Used for the isolation and enumeration of molds in milk and dairy products.	100 Gms 500 Gms	1010 4370
66995	✔ **YPD Broth (Yeast Peptone Dextrose Broth, YEPD Broth)** For maintenance and propagation of yeasts in molecular microbiology procedures.	100 Gms 500 Gms	1120 4480
98271	✔ **YT Broth (2X YT Broth)** Used for the cultivaton of recombinant strains of Escherichia coli.	100 Gms 500 Gms	1240 4370
19510	✔ **YXT Agar Base** Used for detection of yeasts and molds, with or without added Tetracycline.	100 Gms 500 Gms	1350 5160

Microcentrifuge Tube and Screw Cap

Our gamma irradiated microcentrifuge tubes are highly recommended in the areas of Life Science, Industry or Chemistry for use in research and critical reagent packaging requirements.



Our tubes are designed with conical bottoms allowing for complete sample retrieval. Screw caps with O-rings guard against leakage by forming a positive seal against the rim of the tube. External ribs on the free standing microcentrifuge tubes allow for single-handed opening and closing when placed in racks with serrated holes. Screw Cap Microcentrifuge Tubes work with many automated systems. The O-ring is manufactured from durable ethylene polypropylene dimonomere compound.

Product Features:

- Capacity: 0.5 ml & 1.5ml sizes
- Plastic (Polypropylene) clear colorless tube
- Self Standing Tubes with Screw Cap & O-Ring
- Sterile
- Operating range: –196°C to +121°C
- Autoclavable at 121°C
- Withstands maximum RCF of 20000× G
- No Graduations
- DNase & RNase Free

Prdt. Code	Product Name	Packing Unit
31635	Microcentrifuge Tube and Screw Cap (0.5ml), Self Standing, Conical, O-Ring, Clear, Non-Hinged, Sterile, Molecular Grade	500 Nos.
60838	Microcentrifuge Tube and Screw Cap (1.5ml), Self Standing, Conical, O-Ring, Clear, Non-Hinged, Sterile, Molecular Grade	500 Nos.

* Bulk packs available on request

Nucleic Acid Purification Aids:

Spin Column Set

A purification aide utilizing silica membrane in a column form which is used for purifying/extracting nucleic acids (DNA or RNA) of high quality. These columns are compatible with in-house buffers as well as store-brought ones.



Applications:

- Purification of Gel / PCR products
- Extraction of Genomic DNA
- Extraction of Plasmid DNA

Prdt. Code	Product Name	Packing Unit
39738	Mini Silica Spin Column Set A (2ml SC w/CT each)	50 Nos, 100 Nos
89887	Mini Silica Spin Column Set B (2ml SC w/CT each)	50 Nos, 100 Nos
69516	Mini Silica Spin Column Set C (2ml SC w/CT each)	50 Nos, 100 Nos
37351	RNA Mini Spin Column Set (2ml SC w/CT each)	50 Nos, 100 Nos

Collection Tubes

Collection tube of 2ml capacity to be used in purification of DNA / RNA samples.
They are used along with mini silica spin column sets A/B/C & RNA spin column sets.

Prdt. Code	Product Name	Packing Unit
41425	Collection Tubes (CT) for DNA, RNA purification (2ml)	50 Nos, 500 Nos

* Bulk packs available on request

Part E

BioLit®

(DNA & Protein Tools and Kits)

Seraloses™

Media for Gel Filtration and Affinity Chromatography



Gel Filtration (Molecular Sieving or Size Exclusion Chromatography) and Affinity Chromatography are rapidly becoming indispensable to biochemists and molecular biologists. Gel Filtration is used for the separation of macromolecules from smaller molecules, for fractionation of mixtures of various biopolymers on the basis of molecular size, for characterizing macromolecules and determining molecular weights and molar masses, for pathological screening of body fluids and for structure determination studies. Affinity Chromatography, besides being put to many of the above uses, is mainly employed for the isolation and purification of biologically active substances like enzymes, antibodies, antigens, lectins, viruses and toxins. Affinity Chromatography has the advantage of being rapid and highly selective.

SRL Seraloses™	Particle size in wet state (microns)	Fractionation Range (as molecular weight)		Approximate Exclusion Limits (as molecular wight)		
		Proteins	Polysaccarides	Proteins	Polysaccarides	DNA (basis pairs-bp)
Seralose™ 2B Seralose™ Cl-2B	60-250	7.5×10^4 – 45×10^6	1×10^5 – 20×106	45×10^6	20×10^6	1353 bp
Seralose™ 4B Seralose™ Cl-4B	40-190	5×10^4 – 20×10^6	3×10^4 – 5×10^6	20×10^6	5×10^6	872 bp
Seralose™ 6B Seralose™ Cl-6B	40-190	1×10^4 – 5×10^6	1×10^4 – 1×10^4 –	5×10^6	10×10^6	194 bp

Detailed literature available on request

Seralite Products

Ion Exchange Resins

Our important range of Seralite products is very critical for ion exchange applications. These are acidic/basic exchange resins consisting of different polymers and copolymers. The particle size of these resins is between 20 – 50 mesh and are excellent for cross-linking. Typical application of these resins are for water conditioning, softening and deionionisation, removal of recovery of high molecular weight anions, fractional separation of organic acids, amino acids, metals, peptides, inorganic and rare earth metals. They can also be used for th treatment of sugar solutions, amino acids and enzymes. SRL Ion exchange resins are highly recommended products for critical separation work.

Product Name	Particle Size	Ion Exchange Capacity	pH Range	Ionoc Form
Seralite SRA-400	20 – 50 mesh BSS	3.0-3.5meq/g	0-14	Cl-
Seralite SRC-120	20 – 50 mesh BSS	min 4.5meq/g	0-14	H+
Seralite SRA-400	20 – 50 mesh BSS	min 10meq/g	5-14	H+

Our Seralite series is equivalent to the Amberlite series of the Rohm and Haas Company

code	product name	application/description	unit	price

BioLit® PCR Products (Taq Polymerase, Buffers & Kits)

BioLit® TAQ Polymerases

Taq DNA Polymerases (offered in different concentrations) are enzymes used in Routine PCR. They are thermostable recombinant DNA polymerases derived from thermophilic bacterium *Thermus aquaticus*. They can amplify DNA target up to 5 kb (simple template) and have a 5' to 3' polymerase activity but lack 3' to 5' exonuclease activity that result in a 3'-dA overhang PCR product.

BioLit® Routine PCR

Applications - PCR amplification of DNA fragments as long as 5 kb, DNA labeling, DNA sequencing, PCR for cloning.

61524	**TAQ DNA Polymerase (1U/µl)**	Includes Taq DNA Polymerase (1U/µl) 10X PCR Buffer ($MgCl_2$ free), 25mM $MgCl_2$	500 units 1000 units 5000 units	2510.00 4760.00 23120.00
32976	**TAQ DNA Polymerase (3U/µl)**	Includes Taq DNA Polymerase (3U/µl) 10X PCR Buffer ($MgCl_2$ free), 25mM $MgCl_2$	500 units 1000 units 5000 units	2510.00 4760.00 23120.00
82487	**TAQ DNA Polymerase (5U/µl)**	Includes Taq DNA Polymerase (5U/µl) 10X PCR Buffer ($MgCl_2$ free), 25mM $MgCl_2$	500 units 1000 units 5000 units	2510.00 4760.00 23120.00
12816	**TAQ DNA Polymerase (5U/µl) (w/6X Buffer)**	Includes Taq DNA Polymerase (5U/µl) 10X PCR Buffer (Mg^{2+} Plus) 6X Loading Buffer	500 units 1000 units 5000 units	3130.00 5380.00 24370.00
51613	**TAQ DNA Polymerase (5U/µl) (w/6X Buffer, w/dNTP)**	Includes Taq DNA Polymerase (5U/µl) 10X PCR Buffer (Mg^{2+} Plus) 6X Loading Buffer dNTP mix (each 2.5 mM)	500 units 1000 units 5000 units	3380.00 5640.00 25610.00

A Thermostable recombinant DNA polymerase, useful in amplification of DNA fragments upto 6 kb from genomic DNA with high purity and specificity. Sources of template - plant genomic DNA , mouse cDNA, bacterial colony DNA , bacterial plasmid DNA , viral DNA and zygocity PCR templates. No loading buffers or tracking dyes required. Samples may be added directly to an agarose gel after PCR and visualized.

57694	**Red TAQ DNA Polymerase (1U/µl)**	Includes Red Taq DNA Polymerase (1U/µl) 10X Red Taq Buffer, 25 mM $MgCl_2$	500 units 1000 units	4380.00 7630.00
12192	**Roseate Taq DNA (5U/mcl)**		500 units 1000 units	3340.00 5480.00

Roseate Taq DNA Polymerase is a thermostable recombinant DNA polymerase derived from thermophilic bacterium Thermus aquaticus. Its molecular weight is 94 kDa. Taq DNA Polymerase can amplify DNA targets up to 5 kb (simple template). The elongation velocity is 0.9~1.2kb/min (70~75°C). It has 5' to 3' polymerase activity but lacks of 3' to 5' exonuclease activity that result ion a 3'-dA overhangs PCR product.

The buffer given is 10X Rosette PCR Buffer with Mg2+. This buffer has a roseate coloured dye enabling the researcher to load the sample directly after PCR into an agarose gel for electrophoresis.

Other different polymerases described below are useful when there is a need for higher specificity for primer incorporation, high fidelity with proof reading capability, amplifying long and complex templates & need for efficient and fast elongation, and DNA Polymerase which is suitable for Hot Start PCR. The polymerases here each have one of these features and can be chosen to suit the needs of the researcher.

BioLit® Hot Start PCR

This product is a DNA Polymerase coated with antibody, which is suitable for Hot Start PCR. There is no activity of this enzyme at room temperature, and the antibody can be inactivated rapidly at 70°C or above. This ensures that the enzyme is active only during PCR and prolonged heating or denaturation steps can be avoided.

81646	**Hot Start Taq Polymerase (5U/µl)**	Includes Includes HotStart Taq Polymerase (5U/µl) 10X HotStart buffer with Mg^{2+}, 2.5mM dNTP Mix	250 units 1000 units	8510.00 32490.00

code	product name	application/description	unit	price

BioLit® High Specificity PCR

Can amplify DNA target up to 5 kb. The elongation velocity is 0.9~1.2kb/min. It has 5' to 3' polymerase activity but lacks 3' to 5' exonuclease activity that results in a 3'-dA overhang PCR product.

Applications - PCR amplification of DNA fragments as long as 5 kb, DNA labeling, DNA sequencing, PCR for cloning.

46033	**HSP TAQ DNA Polymerase (5U/µl)**	Includes HSP Taq DNA Polymerase (5U/µl) 10X HSP PCR Buffer (with Mg^{2+}) 6×Loading Buffer	500 units 1000 units	5260.00 9750.00

BioLit® High Fidelity PCR

Pfu DNA polymerase, derived from the hyperthermophilic archae Pyrococcus furiosus, has superior thermostability and proofreading properties. It can amplify DNA target up to 2 kb (simple template), possesses 3' to 5' exonuclease proofreading activity, generates PCR fragments will have fewer errors, results in blunt-ended PCR products, which are ideal for cloning into blunt-ended vectors
Applications - High-fidelity PCR and primer-extension reactions, High fidelity PCR for cloning into blunt-ended vectors, Site-directed mutagenesis.

22774	**Pfu DNA Polymerase (5U/µl)**	Includes Pfu DNA Polymerase (5U/µl) 10X Pfu Buffer (Mg^{2+}) 6X Loading Buffer	250 units 1000 units	7130.00 25000.00
49069	**Pfu DNA Polymerase (5U/µl) (w/dNTP)**	Includes Pfu DNA Polymerase (5U/µl) 10X Pfu Buffer (Mg^{2+}) 6X Loading Buffer dNTP mix (each 2.5mM)	250 units 1000 units	7390.00 26240.00

Taq Plus DNA Polymerase is a mixture of Taq and Pfu DNA polymerase that allows amplification of higher fidelity and longer templates Ideal for amplifying complex template, such as GC-rich template.
Applications - Amplification of long template up to 20kb, Amplification of complex template, High fidelity PCR

85737	**TAQ Plus DNA Polymerase (2.5U/µl)**	Includes Taq plus DNA Polymerase (2.5U/µl) 10X PCR Buffer (Mg^{2+} Plus) 6X Loading Buffer	250 units 1000 units	4500.00 16250.00
93561	**TAQ Plus DNA Polymerase (2.5U/µl) (w/dNTP)**	Includes Taq plus DNA Polymerase (2.5U/µl) 10X PCR Buffer (Mg^{2+}) 6X Loading Buffer dNTP mix (each 2.5mM)	250 units 1000 units	4760.00 16750.00

BioLit® Long PCR

Long Taq DNA Polymerase, a combination of two thermostable DNA polymerases, Taq and Pfu, is a special formulation designed for amplifying large fragments. This specially formulated Long Taq was shown to amplifiy long templates from λ phage genome of up to 40 kb. It is also a better choice for amplifying complex template, such as GC-rich template. Long Taq is suitable as a direct replacement for ordinary Taq Polymerase in most applications. Using Long Taq in your PCR results in PCR products with 3´-dA overhangs which can be used in TA clone.

15367	**Long TAQ DNA Polymerase (5U/µl)**	Includes Long Taq DNA Polymerase (5U/µl) 10X Long Taq Buffers I & II (Mg^{2+}) 6X Loading Buffer PCR Enhancer	250 units 1000 units	3630.00 13200.00
94126	**Long TAQ DNA Polymerase (5U/µl) (w/dNTP)**	Includes Long Taq DNA Polymerase (5U/µl) 10X Long Taq Buffers I & II (Mg^{2+}) 6X Loading Buffer PCR Enhancer dNTP mix (each 2.5mM)	250 units 1000 units	5500.00 16500.00

code	product name	application/description	unit	price

BioLit® Fast PCR

FSP Taq DNA polymerase is the latest generation Taq-based DNA polymerase. It possesses high amplification efficiency as Taq polymerase and with fast elongation ability, it can be used in a variety of PCR. The FSP PCR Buffer, designed for FSP Taq DNA polymerase, can be used in fast amplification reactions.
The elongation rate of FSP Taq DNA polymerase is 2-fold higher than the one of regular Taq DNA polymerase, which shortens the amplification time by half.

58235	**FSP TAQ DNA Polymerase (5U/µl)**	Includes FSP Taq DNA Polymerase (5U/µl) 10X FSP Taq Buffer (with Mg^{2+}) 6X Loading Buffer	250 units 1000 units	4010.00 12500.00
70073	**FSP TAQ DNA Polymerase (5U/µl) (w/dNTP)**	Includes FSP Taq DNA Polymerase (5U/µl) 10X FSP Taq Buffer (with Mg^{2+}) 6X Loading Buffer dNTP mix (each 2.5 mM)	250 units 1000 units	4260.00 13010.00

BioLit® PCR Mixes

Routine PCR Master Mixes

Ready-to-use solutions containing the specific DNA Polymerase, dNTP, Mg^{2+} and Reaction Buffer at optimal concentrations for efficient amplification of DNA templates by PCR. This pre-mixed formulation saves time and reduces contamination due to the fewer pipetting steps. Includes Specific DNA Polymerase Mix, Water Nuclease Free.

50900	**Taq Mix (2X) (PCR Master Mix (2X))**	1ml (40 rxns) 5×1ml (200 rxns) 5×2.5ml (500 rxns)	2750.00 8250.00 16500.00
24184	**Taq Mix (2X) (PCR Master Mix (2X))(w/o Tracking Dye)**	1ml (40 rxns) 5×1ml (200 rxns) 5×2.5ml (500 rxns)	2530.00 7810.00 15680.00
46503	**Roseate PCR Master Mix (2X) (Taq Mix (2X))**	1ml (40 rxns) 5×1ml (200 rxns) 5×2.5ml (500 rxns)	2970.00 9020.00 18700.00
94531	**Pfu Mix (2X)**	1ml (40 rxns) 5×1ml (200 rxns)	8070.00 31790.00
84729	**FSP Taq Mix Direct for blood (2X)**	1ml (40 rxns) 5×1ml (200 rxns)	8130.00 33730.00

CRISPR/Cas9 Products

90326	**Cas9 Nickase (50pmol)**	Suitable for in vitro gRNA-guided specific single-strand breaks in double-stranded DNA, or double-strand DNA breaks, etc.	50 mcl	11500.00
72316	**Cas9 Nuclease (5mg/ml)**	CRISPR-Cas9 - Clustered Regularly Interspaced Short Palindromic Repeats- associated Protein 9 is an RNA-guided endonuclease that catalyzes site- specific cleavage of double stranded DNA.	100 mcg	47200.00

code	product name	application/description	unit	price

Realtime PCR Master Mix (qPCR)

2X Master Mix for Real-Time PCR, contains all components, except the primer & template. This reagent is applicable for intercalation assay with FluroGreen dye which is similar to SYBR Green I. Reagent can be used in various detection devices. Also available with passive reference dye of low or high concentrations.
Includes FluroGreen qPCR Mix(2X)w/wt passive refrence dye, Nuclease Free Water .

60159	**BioLit FluroGreen qPCR Master Mix (2X)**	0.5ml (20 amps) 1ml (40 amps) 5×1ml (200 amps)	2750.00 4950.00 19800.00
78797	**BioLit FluroGreen qPCR Master Mix (Low CAR) (2X)**	0.5ml (20 amps) 1ml (40 amps) 5×1ml (200 amps)	3300.00 5500.00 20900.00
38151	**BioLit FluroGreen qPCR Master Mix (High CAR) (2X)**	0.5ml (20 amps) 1ml (40 amps) 5×1ml (200 amps)	3300.00 5500.00 20900.00
45304	**One Step RT-qPCR Master Mix (2X)**	500 mcl	7500.00

BioLit® PCR Accessory Reagents

Ready-to-use reagents for PCR for the user's convenience.

61318	**10X PCR Buffer with KCl (Mg^{2+}free)**		4×1.25ml	1610.00
41212	**10X PCR Buffer with 15mM Mg^{2+}**		4×1.25ml	1610.00
43386	**10X PCR Buffer with Mg^{2+} set**		6×1.25ml	1980.00
90508	**10X Roseate PCR Buffer with Mg^{2+}**	Ready-to-use reagent having two tracking dyes along with the buffer & MgCl2. After PCR, samples can be loaded directly onto the agarose gel for electrophoresis.	4×1.25ml	1560.00
86687	**25mM $MgCl_2$**		4×1.25ml	1820.00
66886	**DEPC-Treated Water**	Suitable in any application where RNA is involved and all molecular biology applications.	100 ml 500 ml	610.00 1890.00
96370	**Water, nuclease-free**	Suitable for use in all molecular biology applications.	5×1ml 100ml 500 ml	320.00 610.00 1890.00
20870	**PCR Enhancer**	Useful to optimize PCR of problematic and/or GC-rich templates. Offers higher primer specificity,and improved themostability of DNA polymerases.	4 x0.1 ml 4×1ml	880.00 3080.00
52531	**PCR Sample Preparation Solution**	Simple, all-inclusive sample preparation solutions for varied samples suitable for direct PCR techniques.	50 preps 200 preps	1540.00 5500.00
48044	**Glycogen (20mg/ml) Aq. Solution Ultrapure**		2×0.25 ml 1 ml	4950.00 9350.00
79326	**Shrimp Alkaline Phosphatase (rSAP) (1U/ul)**		500 units 1000 units	8660.00 16230.00

Shrimp Alkaline Phosphatase (r-SAP) is a heat labile alkaline phosphatase purified from a recombinant source. rSAP is identical to the native enzyme and contains no affinity tags or other modifications. rSAP nonspecifically catalyzes the dephosphorylation of 5´ and 3´ ends of DNA and RNA phosphomonoesters.

code	product name	application/description	unit	price

BioLit® NTP and dNTP Sets, Mixes & Solutions

NTP Set & Solutions (100mM each)

The set and individual solutions consist of 100mM aqueous solutions of ATP, TTP, CTP and GTP each supplied in separate vials. These ready-to-use solutions offer maximum flexibility for various molecular biology applications *viz. in vitro* transcription, RNA amplification.
Applications - for use in PCR, long PCR, RT-PCR, cDNA synthesis, primer extension, DNA sequencing, DNA labeling.

65240	**NTP Set Ultrapure (100mM each)**	4x0.25ml	6050.00
89889	**ATP Solution Ultrapure (100mM)**	0.25ml	2010.00
71636	**CTP Solution Ultrapure (100mM)**	0.25ml	2010.00
88705	**GTP Solution Ultrapure (100mM)**	0.25ml	2010.00
93633	**UTP Solution Ultrapure (100mM)**	0.25ml	2010.00

dNTP Sets (100mM, 10mM)

The sets consists of 10mM/100mM aqueous solutions of dATP, dTTP, dCTP and dGTP, pH 7.0, each supplied in a separate vial. The dNTP Set offers maximum flexibility in preparation of reaction mixes for different molecular biology applications.
Applications - for use in PCR, long PCR, RT-PCR, cDNA synthesis, primer extension, DNA sequencing, DNA labeling.

94903	**dNTP Set Ultrapure (100mM each)**	4×0.25ml	8630.00
		4×1ml	29990.00
69478	**dNTP Set Ultrapure(10mM each)**	4×0.1ml	1890.00
		4×0.5ml	5010.00
		5x(4x0.2ml)	10000.00

dNTP Mixes

dNTP Mix is an aqueous solution at pH 7.0 containing dATP, dCTP, dGTP and dTTP, each at a final concentration of 2.5 mM/10mM/25mM. They are ready-to-use solutions designed to save time, provide higher reproducibility in PCR and other applications, reduce the number of pipetting steps and the risk of reaction set-up errors.
Applications - for use in PCR, long PCR, RT-PCR, cDNA synthesis, primer extension, DNA sequencing, DNA labeling.

70596	**dNTP Mix Ultrapure (2.5mM each)**	1ml	2630.00
		5×1ml	9380.00
18821	**dNTP Mix Ultrapure (10mM each)**	1ml	3510.00
		5×1ml	10630.00
14464	**dNTP Mix Ultrapure (25mM each)**	1ml	4130.00
		5×1ml	16250.00

dNTP Solutions

Each individual dNTP is supplied as a 100mM solution in sterile, deionized water at pH 7.0 and is prepared from monosodium salts of the free nucleotides.
Applications - for use in PCR, long PCR, RT-PCR, cDNA synthesis, primer extension, DNA sequencing, DNA labeling.

72009	**dATP (100mM Solution)**	0.2 ml	2750.00
		0.5 ml	5010.00
19809	**dCTP (100mM Solution)**	0.2 ml	2750.00
		0.5 ml	5010.00
47435	**dGTP (100mM Solution)**	0.2 ml	2750.00
		0.5 ml	5010.00
91474	**dTTP (100mM Solution)**	0.2 ml	2750.00
		0.5 ml	5010.00
71411	**dUTP (100mM Solution)**	0.25ml	12380.00

BioLit® RNA Accessory Reagents

75044	**Oligo (dT) 15 primer (50μM)**	A synthetic single-stranded 15-mer oligonucleotide with 5'- and 3'-hydroxyl ends. Supplied as an aqueous solution. Used for first strand cDNA synthesis synthesis in presence of reverse transcriptase	50 mcl	2260.00
53704	**Random Hexamer Primer 6-mer (50μM)**	A mixture of single-stranded random hexanucleotides with 5'- and 3'-hydroxyl ends. Used during cDNA synthesis to produce cDNAs randomly scattered through mRNA resulting in short strands of cDNA.	50 mcl	2260.00

code	product name	application/description	unit	price
96710	**Ribonuclease A Solution ex. Bovine Pancreas (RNAse A Soln) (10mg/ml) (DNase, Protease Free)**	RNase A is an endoribonuclease extracted from bovine pancreas that specifically degrades single-stranded RNA.	3 x1 ml 10 ml	4010.00 12500.00
38080	**Ribonuclease H (RNase H) (5U/ul)**	RNase H is a recombinant endoribonuclease that hydrolyses specifically the phosphodiester bonds of RNA hybridized to DNA and produces 5´ phosphate-terminated oligoribonucleotides and single-stranded DNA.	50 units 250 units	8340.00 25000.00
92030	**RNasin (RNase inhibitor) (40U/µl)**	Ribonuclease inhibitor purified from human placenta. It inhibits ribonuclease (RNase) activity of common eukaryotic enzymes such as RNaseA, RNase B, RNase C. Inhibits RNA degradation in the following procedures: in vitro transcription, cDNA synthesis, in vitro translation, isolation of mammalian cell fractions that contain mRNA-protein complex etc.	1000 units	4500.00
71969	**RNA Stabilizing Solution**	An aqueous tissue storage reagent that stabilizes and protects cellular RNA in intact, unfrozen tissue samples. It eliminates the need to immediately process tissue samples or to freeze samples in liquid nitrogen. It does not jeopardize the quality or quantity of RNA in the sample.	5 x1 ml 100 ml 500 ml	2780.00 3700.00 8670.00
20310	**RNA Loading Buffer (5X)**	It is used to load RNA sample on a formaldehyde agarose gel. It consists of two tracking dyes, Bromophenol blue and Xylene cyanol.	1 ml 5 x1ml	1320.00 6050.00
98343	**T7 RNA Polymerase (50U/ul) (w/10X IVT Buffer)**	T7 RNA polymerase is a monomeric bacteriophage encoded DNA directed RNA polymerase which catalyzes the formation of RNA in the 5'-3' direction.	1000 units 5000 units	8800.00 42900.00
76516	**M-MLV (Recombinant) (200U/µl)**	A Reverse Transcriptase enzyme that generates first PCR strand cDNA for use in PCR, real-time. Includes 10X Reaction Buffer	5000 units 10000 units	5720.00 9900.00
94837	**RT-PCR Kit (First Strand cDNA Synthesis Kit)**	Reverse Transcriptase -PCR Kit is optimized to synthesize first-strand cDNA from purified poly (A)+ or total RNA. Includes - M-MLV (200U/µl), RNasin (40U/µl), Oligo(dT)15 (50µM), Random hexamer primer (50µM), 5X First-strand Buffer, RNase-free Water, dNTP mix (10mM each)	20 preps 40 preps 100 preps	7810.00 13130.00 31230.00
41066	**TRIzol-C Reagent**	Ready-to-use reagent for the isolation of high yields of total RNA from cells and tissues. As it is differently buffered from TRIzol-T reagent, the isolated RNA may interact with these metal ions if involved in downstream processes.	20 ml 100 ml	2200.00 7370.00
99359	**TRIzol-T Reagent**	Ready-to-use reagent for the isolation of high yields of total RNA from cells and tissues. Differs in the buffer of TRIzol-T reagent, and is innocuous to all downstream processes.	20 ml 100 ml	1870.00 7150.00

BioLit® Tracking dyes & Buffers

6x Gel Loading Dyes are pre-mixed loading buffers with tracking dye/s to be added to the nucleic acid samples for agarose and non-denaturing poylacrylamide gel electrophoresis.

78591	**6X Blue Loading Dye**	Contains Bromophenol blue as tracking dye.	5×1 ml	3050.00
83184	**6X DNA Loading Dye**	Contains Bromophenol blue and Xylene Cyanol FF as tracking dyes.	5×1 ml	3050.00
56885	**6X DNA Loading Dye & SDS solution**	Contains Bromophenol blue and Xylene Cyanol FF as tracking dyes. The buffer also contains SDS which results in sharper bands as some restriction enzymes remain bound to DNA following cleavage	5×1 ml	3440.00
46200	**6X Orange Loading Dye**	Contains Orange G and Xylene Cyanol FF as tracking dyes	5×1 ml	3050.00
55895	**6X TriColor Loading Dye**	Contains Orange G, Bromophenol Blue and Xylene Cyanol FF as tracking dyes	5×1 ml	3050.00

code	product name	application/description	unit	price
73530	**20X DoubleQuick DNA Electrophoresis Buffer**	An exceptionally 'swift' agarose gel tank buffer that electrophoreses the DNA samples in less than 10 minutes. Three times more effective as compared to TAE or TBE buffer.	250 ml 500 ml	4580.00 8290.00

BioLit® Staining Dyes for Nucleic Acids

Nucleic acid & protein staining dyes.

50358	**BioLit FluroBronze Stain**	A non-carcinogenic nucleic acid stain that is eco-friendly and replaces the use of ethidium bromide. Used for viewing nucleic acids on agarose gel.	50 mcl 250 mcl 1 ml	990.00 3300.00 9350.00
	Advantage: The stain is not added to the gel matrix or the running buffers, but is mixed as a 6X loading buffer with the nucleic acid samples prior to loading onto the gel.			
21156	**SafeDye Green Stain (20000X in Water)**		100 mcl 500 mcl 1 ml	2750.00 8800.00 15400.00
25834	**BioLit SafeDye Stain**	SafeDye is a new, safe nucleic acid stain for the detection of double-stranded DNA, single-stranded DNA and RNA in Agarose gels. This dye replaces Ethidium Bromide (toxic, potential mutagen) commonly used in Agarose gel electrophoresis and assists in determining the size of PCR product, insert and restriction analysis.	50 mcl 250 mcl 1 ml	880.00 3080.00 8800.00
	SafeDye emits green fluorescence when bound to dsDNA, ssDNA and RNA. This stain has a fluorescence excitation maxima when bound to nucleic acid at approx. 290-320nm and emits at 490nm.			

BioLit® Reference Digests/Markers

87978	**ProxiB Lambda HindIII Digest Marker** Size range: 2.0 to 23kb	Lambda DNA is used routinely for measurement of restriction enzyme activities. Lambda-Hind III digestion is a useful demonstration experiment for restriction analysis and it is a substrate for checking DNA ligase. This is a ready-to-load marker and contains Bromophenol Blue as the tracking dye. (bands at 0.1 & 0.5kb although present, are very faint.)	50 mcg 100 mcg	2090.00 3960.00
90023	**ProxiB Lambda-EcoRI Digest Marker** Size range: 3.5 to 21kb	Lambda-EcoRI digest is a popular teaching tool to demonstrate the activity of restriction enzymes. This is a ready-to-load marker and contains Bromophenol Blue as the tracking dye.	50 mcg 100 mcg	2090.00 3960.00
53625	**ProxiB Lambda Double Digest Marker** Size range: 0.5 to 23kb	Lambda DNA double digest of Hind III/EcoRI is a popular teaching tool to demonstrate double-digests (i.e., the simultaneous use of two restriction enzymes). This is a ready-to-load marker and contains Bromophenol Blue as the tracking dye.	50 mcg 100 mcg	2530.00 4180.00
46658	**ProxiB Lambda BstEII Digest Marker** Size range: 0.1-8.4 kb		50 mcg	2860.00
93663	**Lambda PstI Digest Marker** Size range: 0.25-11.5 kb		50 mcg	2860.00

code	product name	application/description	unit	price

BioLit DNA Reagents

64438	**DH5α Competent Cells for Subcloning**		2 x0.1 ml 4 x0.5 ml	2900.00 20790.00
67291	**Lambda DNA**	Lambda DNA is a double stranded linear DNA of 48,502 base pairs in length. It is isolated from bacteriophage Lambda. The Lambda Phage is a bacterial virus (bacteriophage), that infects the bacterial species E. coli. The phage is isolated from the heat inducible lysogen *E. coli* strain. This NA is useful as size markers, in restriction enzyme assays and cloning. Concentration — 250ng/µl	100 mcg 250 mcg	4620.00 11000.00
12658	**Autoinduction Medium (AIM) Solution**	AIM-LB medium w/o Trace elements is used to grow IPTG-inducible expression in bacterial strains	500 ml	1760.00
69359	**DNA Methylation Reagent**	The modified DNA is ideal for PCR amplification for downstream analyses, including restriction endonuclease digestion, sequencing, and microarrays.	5 ml 10 ml	1900.00 3500.00
78116	**DNA Quick Extraction Reagent**	The isolated DNA can be used for downstream applications such as Restriction Digestion, Southern blot analysis, Polymerase chain reaction (PCR) and Molecular cloning.	20 ml 100 ml	2000.00 6900.00
79326	**Shrimp Alkaline Phosphatase (rSAP) (1U/ul)**		500 Units 1000 Units	9530.00 17860.00
56002	**Uracil DNA Glycosylase (5KU/ml)**	Monofunctional DNA glycosylase that catalyzes the hydrolysis of the N-glycosidic bond from deoxyuridine to release Uracil. Active on ss and dsDNA.	100 Units 500 Units	3500.00 9200.00

BioLit® Cloning Vectors

A cloning vector is a small piece of DNA into which a foreign DNA fragment can be inserted. The insertion of the fragment into the cloning vector is carried out by treating the vehicle and the foreign DNA with the same restriction vector, then ligating the fragments together. There are many types of cloning vectors. Genetically engineered plasmids are most commonly used for this purpose.

55912	**pBR322** pBR322 was the first artificial plasmid. It is a 4361 bp long plasmid and was originally developed for cloning. This low copy number plasmid contains tetracycline and ampicllin resistance-conferring genes. pBR322 remains a popular teaching aid for demonstration of plasmid isolation, restriction digestion and cloning.	50 mcg 100 mcg	4180.00 6600.00
93116	**pUC18** pUC18 vector is a circular double stranded DNA and is 2686 base pairs long. This plasmid possesses – i) ampicillin resistance gene ii) the origin of replication derived from pBR322 and iii) the lac Z gene derived from *E. coli*. Within the lac region is also found a polylinker sequence having unique restriction sites. It is a popular teaching aid for demonstrating plasmid DNA isolation, restriction digestion and cloning.	50 mcg 100 mcg	4180.00 6600.00
48402	**pUC19** pUC19 vector is also a circular double stranded DNA and is 2686 base pairs long. It is similar to pUC18, but possesses different restriction sites. This plasmid possesses – i) ampicillin resistance gene ii) the origin of replication derived from pBR322 and iii) the lac Z gene derived from *E. coli*. Within the lac region is also found a polylinker sequence having unique restriction sites. It is a popular teaching aid for demonstrating plasmid DNA isolation, restriction digestion and cloning.	50 µg 100 µg	4180.00 6600.00

code	product name	application/description	unit	price

BioLit® Ready-to-use DNA Molecular Weight Standards

BioLit DNA markers come in many useful ranges to cover most of the requirements of researchers and students. The markers are provided in a ready-to-use format for convenience and high reproducibility. We recommend a volume of 10 µL (~1 µg total DNA) per load to ensure that distinct, prominent DNA bands are observed even with a less sensitive UV transilluminator or in a less than optimal dark chamber. Users might be able to load a volume less than that recommended (without compromising results) if a high-sensitivity gel documentation system is available.

The markers have been carefully prepared under sterile conditions and stored in 10 mM Tris-HCl, pH 7.5, 1 mM EDTA and 1 mg/ml sodium azide and are unstained markers. Stocks will be stable for several months if properly stored frozen or even refrigerated.

10344	**Low-Range DNA Marker: 0.1 – 1kb**	This marker range has been designed to assist in determining the size of PCR product, insert and restriction analysis. These ready-to-load DNA markers contain Bromophenol Blue as the tracking dye. Size range: 103 to 1074bp	15 mcg 50 mcg 100 mcg	2260.00 6380.00 12010.00
12447	**Mid-Range DNA Marker: 0.5 – 5kb**	This marker range should facilitate the analysis of inserts cloned into commonly used plasmids and restriction digests. These ready-to-load DNA markers contain Bromophenol Blue as the tracking dye. Size range: 448 to 5024bp	15 mcg 50 mcg 100 mcg	2260.00 6380.00 12010.00
45587	**Hi-Range DNA Marker: 2.0 – 25kb**	This marker range is useful for experiments involving mega plasmids or genomic analysis. These ready-to-load DNA markers contain Bromophenol Blue as the tracking dye. Size range: 2.0 to 25kb	15 mcg 50 mcg 100 mcg	2260.00 3510.00 6380.00
48234	**ProxiO DNA Marker 2 : 0.1-1.2kb**	Size range: 100 to 1200bp These ready-to load DNA markers contain Bromophenol Blue and Xylene Cyanol as the tracking dye. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg	9970.00

BioLit® DNA Ladders

These ready-to-load ladders are designed with precise DNA fragments for accurate quantification by agarose gel electrophoresis. The concentration of each band is kept uniform to give sharp bands. A single or multiple bands in the ladders are made intense to orient the user to the other bands in the ladder. Recommended Loading volume 3-5µl. The products are unstained DNA Ladders.

44897	**ProxiO 10bp DNA Ladder**	Size range: 80 to 300 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg	11550.00
98628	**ProxiO 20bp DNA Ladder**	Size range: 60 to 300 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg	10340.00
68905	**ProxiO Low DNA Ladder**	Size range: 25 to 700 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg 5X 50 mcg	8690.00 42570.00

code	product name	application/description	unit	price
82208	**ProxiO 50bp DNA Ladder**	Size range: 50 to 500 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg 5X 50 mcg	8250.00 26950.00
88417	**ProxiO 50bp DNA Ladder Plus**	Size range: 50 to 1,000 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg 5X 50 mcg	7920.00 27390.00
92972	**ProxiO 100bp DNA Ladder**	Size range: 100 to 1,500 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg 5X 50 mcg	6600.00 25410.00
84628	**ProxiB 100bp DNA Ladder**	Size range: 100 to 1,500 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes.	50 mcg 5X 50 mcg	6600.00 25410.00
52075	**ProxiO 100bp DNA Ladder Plus**	Size range: 100 to 3,000 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg 5X 50 mcg	6600.00 25410.00
50255	**ProxiB 100bp DNA Ladder Plus**	Size range: 100 to 3,000 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes.	50 mcg 5X 50 mcg	6600.00 25410.00
37217	**ProxiO 200bp DNA Ladder**	Size range: 200 to 4,000 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg	10120.00
93121	**ProxiO 500bp DNA Ladder**	Size range: 500 to 5,000 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg	10560.00
48290	**ProxiO 1kb DNA Ladder**	Size range: 500 to 10,000 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg 5X 50 mcg	6490.00 27500.00
93406	**ProxiO 1kb DNA Ladder Plus**	Size range: 100 to 10,000 bp Contains Bromophenol Blue and Xylene Cyanol as the tracking dyes. Accompanied by a vial of 6X Orange Loading Dye.	50 mcg 5X 50 mcg	6490.00 27500.00
77305	**PCR Ladder**	This size marker covers the most common range encountered in PCR and shows an aesthetically pleasing ladder of five bands (50, 150, 300, 500, and 800 bp) that separate uniformly as shown. Contains Bromophenol Blue as the tracking dye.	10 mcg 50 mcg	2200.00 6930.00

BioLit® FluroBronze™ DNA Ladder Packs

The varying range of specially prepared DNA ladder packs (ready-to-load) contain in addition to the tracking dyes, a non-carcinogenic nucleic acid stain (BioLit FluroBronze Stain) for detection of the ladder bands. This nucleic acid stain replaces ethidium bromide in agarose gel preparation thus simplifying the preparation of agarose gel sans a potentially mutagenic dye.

97958	**FluroBronze™ 100bp Ladder**	Size range: 100 to 1,500 bp Contains Bromophenol Blue & Xylene Cyanol as the tracking dyes along with FluroBronze stain as the fluorescent dye. Also accompanied by a vial of 6X FluroBronze Stain.	50 mcg	9570.00

code	product name	application/description	unit	price
72601	**FluroBronze™ 100bp Ladder Plus**	Size range: 100 to 3,000 bp Contains Bromophenol Blue & Xylene Cyanol as the tracking dyes along with FluroBronze stain as the fluorescent dye. Also accompanied by a vial of 6X FluroBronze Stain.	50 mcg	9570.00
54896	**FluroBronze™ 500bp Ladder**	Size range: 500 to 5,000 bp Contains Bromophenol Blue & Xylene Cyanol as the tracking dyes along with FluroBronze stain as the fluorescent dye. Also accompanied by a vial of 6X FluroBronze Stain.	50 mcg	16610.00
23416	**FluroBronze™ 1kb Ladder Plus**	Size range: 100 to 10,000 bp Contains Bromophenol Blue & Xylene Cyanol as the tracking dyes along with FluroBronze stain as the fluorescent dye. Also accompanied by a vial of 6X FluroBronze Stain.	50 mcg	11880.00

BioLit® SafeDye™ DNA Ladder Packs

The varying range of DNA ladder packs (ready-to-load) contain DNA ladders accompanied by the BioLit SafeDye Stain, a non-carcinogenic nucleic acid stain mixed in a 6X loading buffer, for detection of the ladder bands. This nucleic acid stain replaces ethidium bromide in agarose gel preparation thus eliminating a potentially mutagenic dye.

64935	**SafeDye™ PCR Ladder**	Size range: 50 to 800 bp Contains Bromophenol Blue as the tracking dye. Also accompanied by a vial of 6X SafeDye Stain.	50 mcg	11990.00
19262	**SafeDye™ 50bp Ladder Plus**	Size range: 50 to 1,000 bp Contains Bromophenol Blue & Xylene Cyanol as the tracking dyes. Also accompanied by a vial of 6X SafeDye Stain.	50 mcg	12320.00
33902	**SafeDye™ 100bp Ladder**	Size range: 100 to 1,500 bp Contains Bromophenol Blue & Xylene Cyanol as the tracking dyes. Also accompanied by a vial of 6X SafeDye Stain.	50 mcg	12540.00
53261	**SafeDye™ 100bp Ladder Plus**	Size range: 100 to 3,000 bp Contains Bromophenol Blue & Xylene Cyanol as the tracking dyes. Also accompanied by a vial of 6X SafeDye Stain.	50 mcg	10560.00
27107	**SafeDye™ 500bp Ladder**	Size range: 500 to 5,000 bp Contains Bromophenol Blue & Xylene Cyanol as the tracking dyes. Also accompanied by a vial of 6X SafeDye Stain.	50 mcg	14960.00
29877	**SafeDye™ 1kb Ladder Plus**	Size range: 100 to 10,000 bp Contains Bromophenol Blue & Xylene Cyanol as the tracking dyes. Also accompanied by a vial of 6X SafeDye Stain.	50 mcg	11880.00

Surface Decontamination Spray

69341	**DNA Decontam Surface Decontamination Spray**	Eliminate unwanted DNA and DNase from glassware and plastic ware without affecting subsequent DNA samples and plastic ware without affecting subsequent DNA samples. Supplied in a spray bottle	250ml	3850.00
27248	**RNAse Decontam Surface Decontamination Spray**	A surface decontamination solution that destroys RNases on contact. Spray the Solution onto surfaces to be decontaminated and rinse off with RNase-free-water. Supplied in a spray bottle format.	250ml	4400.00

code	product name	application/description	unit	price

BioLit® Restriction Enzymes

BioLit® Restriction Enzymes are highly purified and ensure that the integrity of your DNA is not compromised, making them the enzymes of choice for even the most demanding applications. All the Restriction Enzymes are accompanied with their specific reaction buffer to enable 100% enzyme activity.

85844	**AluI (10U/µl)**	100 unit	6050.00
51914	**BamHI (10U/µl)**	2500 unit	3630.00
		10000 unit	12100.00
32338	**BglI (10U/µl)**	1500 unit	7070.00
24848	**BglII (10U/µl)**	250 unit	2420.00
		500 unit	4400.00
99402	**DraI (10U/µl)**	1000 unit	9350.00
91290	**EcoRI (10U/µl)**	4000 unit	3850.00
		10000 unit	9240.00
90279	**EcoRV (Eco321) (10U/µl)**	500 unit	3300.00
		1500 unit	8250.00
61338	**HaeIII (BsuRI, BshFI) (10U/µl)**	3000 unit	11000.00
		10000unit	35200.00
44261	**HincII (HindII) (10U/µl)**	250 units	3850.00
		1000 units	13750.00
77976	**HindIII (10U/µl)**	2500 units	4180.00
		5000 units	8250.00
36322	**HinfI (10U/µl)**	750 unit	3960.00
15100	**KspAI (HpaI) (10U/µl)**	200 unit	5390.00
		500 unit	12650.00
47167	**KpnI (10U/µl)**	1000 unit	2420.00
		5000 unit	11000.00
65301	**MboI (10U/µl)**	300 unit	13200.00
		1100 unit	44000.00
79971	**NcoI (10U/µl)**	125 unit	3190.00
		300 unit	6880.00
75760	**NdeI (10U/µl)**	500 unit	2420.00
		2500 unit	9900.00
97361	**NotI (10U/µl)**	100 unit	5280.00
		250 unit	12100.00
94651	**PstI (10U/µl)**	1500 unit	3300.00
		4500 unit	8800.00
72649	**PvuII (10U/µl)**	500 unit	2040.00
		2000 unit	6380.00
97505	**RsaI (10U/µl)**	250 unit	3630.00
		1000 unit	13200.00
27301	**SacI (10U/µl)**	750 unit	6600.00
		1250 unit	9900.00
10190	**SalI (10U/µl)**	500 unit	3960.00
		750 unit	5720.00

code	product name	application/description	unit	price
74952	**Sau3AI (Bsp143I) (10U/µl)**		300 unit	3960.00
			1500 unit	18700.00
15231	**SmaI (10U/µl)**		500 unit	6600.00
49365	**SpeI (BcuI) (10U/µl)**		200 unit	4070.00
			400 unit	7480.00
40653	**SspI (10U/µl)**		250 unit	4070.00
			500 unit	7700.00
96600	**XbaI (10U/µl)**		750 unit	3410.00
			1250 unit	5060.00
			4000 unit	14850.00
71245	**XhoI (SalI) (10U/µl)**		1000 unit	3300.00
			5000 unit	13200.00
46448	**Universal Restriction Enzyme Buffer (10X)**		5 x1 ml	1320.00

BioLit® QuickRxn Restriction Enzymes

These restriction enzymes reduce digestion time and speed up the downstream processes.

81536	**QuickRxn BamHI (20U/mcl)**	250 units	660.00
		1000 units	1870.00
15918	**QuickRxn EcoRI (20U/mcl)**	250 units	550.00
		1000 units	1320.00
15723	**QuickRxn HindIII (20U/mcl)**	250 units	550.00
		1000 units	1650.00
38616	**QuickRxn KpnI (20U/mcl)**	250 units	1050.00
		1000 units	2640.00
54781	**QuickRxn NcoI (20U/mcl)**	100 units	3300.00
		250 units	5500.00
57426	**QuickRxn NdeI (20U/mcl)**	250 units	1650.00
		1000 units	4620.00
60399	**QuickRxn NotI (20U/mcl)**	100 units	5610.00
		250 units	12980.00
64479	**QuickRxn PstI (20U/mcl)**	250 units	770.00
		1000 units	2420.00
37862	**QuickRxn SalI**	250 units	2420.00
49604	**QuickRxn SmaI (20U/mcl)**	250 units	3520.00
25887	**QuickRxn XhoI (20u/mcl)**	250 units	1540.00
		1000 units	3960.00

code	product name	application/description	unit	price

BioLit® Ligases

This DNA ligase, an ATP-dependent recombinant enzyme, isolated from Escherichia coli, is used for DNA cloning

85315	**T4 DNA Ligase VLC (20CEU/µl or 0.1 Weiss Unit/µl)**	T4 DNA ligase is a polypeptide that catalyzes the formation of phosphodiaster bonds between adjacent 3'-hydroxyl and 5'- phosphate termini in DNA. The enzyme is active on double-stranded DNA, nicked DNA and RNA substrates. One Weiss unit is equivalent to approximately 200 cohesive end ligation units (CEU) (Includes Reaction Buffer)	10000 Units (CEU) 40000 Units (CEU)	2750.00 5500.00
41686	**T4 DNA Ligase LC (200CEU/µl or 1 Weiss Unit/µl)**	T4 DNA ligase is a polypeptide that catalyzes the formation of phosphodiaster bonds between adjacent 3'-hydroxyl and 5'- phosphate termini in DNA. The enzyme is active on double-stranded DNA, nicked DNA and RNA substrates. One Weiss unit is equivalent to approximately 200 cohesive end ligation units (CEU) (Includes Reaction Buffer)	200 Units (Weiss) 1000 Units (Weiss) 5000 Units (Weiss)	2640.00 4620.00 13200.00
67428	**T4 DNA Ligase RC (1000CEU/µl or 5 Weiss Unit/µl)**	T4 DNA ligase is a polypeptide that catalyzes the formation of phosphodiaster bonds between The enzyme is active on double-stranded DNA, nicked DNA and RNA substrates. One Weiss unit is equivalent to approximately 200 cohesive end ligation units (CEU) (Includes Reaction Buffer)	250 Units (Weiss) 1000 Units (Weiss)	4950.00 15400.00
32581	**Tth DNA Ligase Swift (5U/µl) Thermostable**	An enzyme isolated from Escherichia Coli strain containing plasmid carrying the Thermus thermophilus DNA ligase, ligates DNA fragments within 15 minutes. Works on sticky ends only.	100 Units (Weiss) 200 Units (Weiss)	12100.00 24200.00

BioLit® Ready-to-use Unstained Protein Molecular Weight Standards

BioLit protein molecular weight standards cover the useful range of protein sizes that a researcher or student requires for routine experiments. The recommended loading volume is 20 µl, giving approximately 5 µg of each protein sufficient for detection by Coomassie blue staining. The products are unstained markers. The proteins have been dissolved and denatured in the standard SDS-PAGE sample-loading buffer. Prepared under sterile conditions, these markers have appropriate additives to ensure long-term stability when they are stored in a refrigerator or freezer and even with brief exposure to room temperature. Although boiling is not required prior to loading, doing so will not affect results.

58511	**Mid-Range Protein Marker: 14-95kDa**	This ready-to-use marker covers a slightly higher range of protein sizes of 14 - 95 kDa and contains Bromophenol Blue as the tracking dye. It has eight bands of varying molecular sizes that a researcher or student requires for routine experiments.	0.5 ml 1 ml 2.5 ml	2690.00 4840.00 10190.00
59893	**Low-Range Protein Marker: 3-40kDa**	This protein molecular weight standard includes proteins in the range of 2.5 - 40 kDa and contains Bromophenol Blue as the tracking dye. With the advent of proteomics and identification of many low molecular weight proteins and domains being involved in important functions, there is an emerging interest in polypeptides in this range. This is also a useful range for electrophoresis based peptide mapping, especially using Tricine SDS - PAGE.	0.5 ml 1 ml 2.5 ml	2690.00 4840.00 10190.00

code	product name	application/description	unit	price
99625	**High-Range Protein Marker: 14-220kDa**	This molecular range marker covers a larger range of molecular sizes - 14 to 220kDa and is helpful to researchers working on higher molecular weight proteins. It contains Bromophenol Blue as the tracking dye.	0.5 ml 1 ml	3300.00 6050.00

Ready-to-use Prestained Protein Molecular Weight Standards

68894	**PreStained Blue Protein Marker: 14 - 66kDa**	Blue Prestained Protein Molecular Weight Standard, contains a mix of proteins, coupled with a blue chromophore. This marker has 6 bands in range of 14 to 66kDa.	0.5 ml	2810.00
27839	**PreStained High Range 1 Tri-Coloured Protein Ladder: 10-180kDa**	This Protein Ladder is a three-color protein standard with 10 bands in the range of 10 to 180 kDa. Proteins are covalently coupled with a blue chromophore except for two reference bands (one green and one red band at 25 kDa and 75 kDa respectively) when separated on SDS-PAGE (Tris-glycine buffer).	0.25 ml 0.5 ml 1 ml	7200.00 14000.00 26000.00
67792	**PreStained High Range 2 Tri-Coloured Protein Ladder: 10-245kDa**	This Protein Ladder is a three-color protein standard with 12 bands in the range of 10 to 245 kDa. Proteins are covalently coupled with a blue chromophore except for two reference bands (one green and one red band at 25 kDa and 75 kDa respectively) when separated on SDS-PAGE (Tris-glycine buffer).	0.25 ml 0.5 ml 1 ml	8000.00 15000.00 28000.00
49772	**PreStained High Range 3 Tri-Coloured Protein Marker: 10-315kDa**	A ready-to-use Prestained protein marker giving 12 discrete fragments (in kDa): 315 (orange), 250, 180, 130, 95 (Blue) , 72 (orange), 52 (green), 43, 34, 26, 17 (blue), and 10 (green).	0.25 ml	35000.00

BioLit® Staining Dyes for Proteins

25437 Store at 2-8°C	**BioLit SpryBlue Stain for SDS PAGE**	A snappy stain for Polyacrylamide gels that accelerates the staining of proteins and eliminates the need for destains.	250 ml 500 ml	880.00 1430.00

Buffers for Proteins

85074	**Borate Buffer (20X) (Amine-free)**	Ideal buffer solution that is used in protein modification procedures requiring amine-free buffer at alkaline pH.	100 ml 500 ml	2420.00 8250.00
88370	**RIPA Protein Extraction Reagent** For RIPA assay (Radioimmunoprecipitation assay)	A mild, non-denaturing ready-to-use lysis buffer with optimised pH, designed for quick and effective preparation and isolation of nuclear and cytoplasmic proteins	50 ml 100 ml 500 ml	5800.00 9200.00 24200.00
88480	**SDS PAGE Sample Buffer (2X) (Reducing) (Laemmli Buffer 2X)**	It is the most commonly used sample buffer for SDS-PAGE of denatured proteins it is a ready-to-use 2X solution. It can be used for SDS-PAGE protein loading of conventional proteins.	5 ml 20 ml	430.00 1650.00
99573	**SDS PAGE Sample Buffer (5X) (Reducing) (Laemmli Buffer 5X)**	It is the most commonly used sample buffer for SDS-PAGE of denatured proteins it is a ready-to-use 5X solution. It can be used for SDS-PAGE protein loading of conventional proteins.	5 ml 20 ml	530.00 1980.00
98786	**SDS PAGE Sample Buffer (5X) (Reducing) w/ DTT w/ BPB**	It is the most commonly used sample buffer for SDS-PAGE of denatured proteins it is a ready-to-use 5X solution. It can be used for SDS-PAGE protein loading of conventional proteins.	5 ml	2750.00
35974	**SDS PAGE Sample Buffer (5X) (Non-Reducing) w/o DTT w/ BPB**	It is the most commonly used sample buffer for SDS-PAGE of denatured proteins it is a ready-to-use 5X solution. It can be used for SDS-PAGE protein loading of conventional proteins.	5 ml	2750.00

code	product name	application/description	unit	price

BioLit® SDS PAGE Accessory Reagents

Ready-to-use Separating gel mix in a various percentage ranges and Stacking gel mix to simplify the process of SDS PAG preparation and eliminate pipetting errors.

91574	**Ready-To-Use Stacking Gel mix - 4%**		250 ml	2530.00
86325	**Ready-To-Use Stacking Gel mix - 4% (Coloured)**		100 ml	1210.00
24723	**Ready-To-Use Separating Gel mix - 10%**		100 ml	2860.00
80306	**Ready-To-Use Separating Gel mix - 12%**		100 ml	2860.00
88394	**Ready-To-Use Separating Gel mix - 15%**		100 ml	2860.00

BioLit® Protease for Recombinant Protein Fusion Tag Cleavage

86851	**TEV Protease (1U/µl)**	A highly sequence-specific cysteine protease from Tobacco Etch Virus (TEV). Removes affinity tags from purified recombinant fusion proteins and its specificity gives controlled cleavage. Effective over a broad range of temperature	1000 Units 5000 Units	11880.00 45980.00

Antibodies

41262	**Express Coating Buffer (1X, pH 6.4) for Antigens & Antibodies**	A special coating buffer for rapid adsorptive immobilization of proteins, enzymes, antigens or antibodies on any plastic surface including ELISA plates. It coats the surface within 10-15mins.	100 ml 500 ml	2750.00 8800.00
91754	**Express Coating Buffer (1X, pH pH 8.5) for Antigens & Antibodies**	A special coating buffer for rapid adsorptive immobilization of proteins, enzymes, antigens or antibodies on any plastic surface including ELISA plates. It coats the surface within 10-15mins.	100 ml 500 ml	3850.00 9900.00
37971	**Goat anti-Human IgG (affinity purified) (2mg/ml)**	Antiserum is produced in Goat using highly purified human IgG as the immunogen. Whole antiserum is fractionated and then further purified by affinity chromatography to provide the IgG fraction of antiserum.	1 ml	5610.00
77924	**Goat anti-Human IgG-HRP (1mg/ml)**	Purified IgG conjugated to Horseradish Peroxidase (HRP). Pooled antisera from goats hyperimmunized with human IgG	1 ml	14850.00
57404	**Goat anti-Mouse IgG-HRP (1mg/ml)**	Polyclonal Goat antibodies conjugated with HRP. Pooled antisera from goats hyperimmunized with Mouse IgG	1 ml	14850.00
54537	**Goat anti-Rabbit IgG (whole serum) (3mg/ml)**	Antiserum is produced in goat using highly purified rabbit IgG as the immunogen. Whole antiserum is fractionated and then further purified by affinity chromatography to provide the IgG fraction of antiserum.	1 ml 5 ml	3740.00 14850.00
88495	**Goat anti-Rabbit IgG-HRP (1mg/ml)**	Polyclonal Goat antibodies conjugated with HRP. Pooled antisera from goats hyperimmunized with Rabbit IgG.	1 ml	14850.00
68977	**Rabbit anti-Bovine IgG (whole serum) (2mg/ml)**	Antiserum is produced in Rabbit using highly purified bovine IgG as the immunogen. Whole antiserum is fractionated and then further purified by affinity chromatography to provide the IgG fraction of antiserum.	1 ml	2970.00
21105	**Rabbit anti-Bovine IgG-HRP (whole serum) (1mg/ml)**	Anti-Bovine secondary antibodies conjugated with HRP. Pooled antisera from cows hyperimmunized with rabbit IgG	1 ml	14850.00

code	product name	application/description	unit	price
57919	**Rabbit anti-Chicken IgY (whole serum)**		1 ml 5 ml	4950.00 22330.00
26249	**Rabbit Anti-Human IgG (whole serum)**		1 ml 5 ml	4950.00 22330.00

BioLit® Research Kits

A column based quick and easy method for DNA Extraction from bacteria/ tissue/whole blood yielding pure DNA useful for down-stream processes.

51826	**BioLit Genomic DNA Extraction Mini Kit (Research)** Kit contents : ● Spin Columns ● Buffers ● Proteinase K	Isolation of high quality genomic DNA using the silica spin columns from various sources - plasma, serum, body fluids, cultured cells, sputum and other body fluids.	25 Extractions 50 Extractions	8700.00 16600.00
66865	**BioLit Plasmid DNA Extraction Miniprep Kit (Research)**		50 Preps	9100.00
56860	**BioLit Plant DNA Extraction Mini Kit (Research)** Kit contents : ● Spin Columns ● Buffers	Kit works to efficiently extract genomic DNA from various plant and tissue samples rapidly. Mini silica columns deliver high quality genomic DNA without use of expensive reagents and time consuming processes.	25 Extractions 50 Extractions	6100.00 11000.00
75718	**BioLit Tissue DNA Extraction Mini Kit (Research)** Kit contents : ● Spin Columns ● Buffers ● Proteinase K	Isolation of high quality tissue DNA using the silica spin columns from soft and hard tissues samples like liver, spleen, thymus, heart, kidney, brain, etc	25 Extractions 50 Extractions	11000.00 18800.00
94733	**BioLit Whole Blood DNA Extraction Mini Kit (Research)** Kit contents : ● Spin Columns ● Buffers	Isolation of high quality DNA using the silica spin columns from human whole blood.	25 Extractions 50 Extractions	10400.00 19400.00
33981	**BioLit Tissue RNA Extraction Mini Kit (Research)** Kit contents : ● Spin Columns ● Collection Tubes ● Buffers ● Carrier RNA -Poly A	Tissue RNA can be purified from soft and hard tissues samples like liver, spleen, thymus, heart, kidney, brain, etc	25 Extractions 50 Extractions	12500.00 23600.00
87883	**BioLit Viral RNA Extraction Mini Kit (Research)** Kit contents : ● Spin Columns ● Collection Tubes ● Buffers ● Carrier RNA -Poly A	Viral RNA can be purified from plasma, serum and other cell body fluids.	25 Extractions 50 Extractions 100 Extractions	7100.00 12100.00 22000.00
76108	**BioLit Whole Blood RNA Extraction Mini Kit (Research)** Kit contents : ● Spin Columns ● Collection Tubes ● Buffers ● ß-mercaptoethanol	A fast, easy method for the preparation of total cellular RNA from human whole blood. Multiple blood samples can be processed simultaneously in less than 1 hour.	25 Extractions 50 Extractions	13000.00 25200.00

code	product name	application/description	no. of experiments	unit	price
91210	**BioLit Animal Tissue Direct PCR Kit (Research)** Kit contents : ● 2X Taq mix ● Buffers ● Sterile Water	Amplifies DNA directly from a wide variety of animal tissues including mice, fish, birds and insects.		50 Preps	31400.00
74578	**BioLit Gel Extraction Kit (Research)** Kit contents : ● Silica Spin Columns● Buffers	Kit for quick and easy extraction of DNA fragment from Low melting agarose. Useful for various downstream processes.		50 Preps	12800.00
23512	**BioLit PCR and DNA Fragment Purification Kit (Research)**			50 Preps 100 Preps	13400.00 26000.00
62606	**BioLit Total RNA Extraction Kit from Plant (Research)** Kit contents : ● Spin Columns● Buffers	Quick Total RNA extraction kit using silica gel columns. Purifies RNA from Plants and Fungi.		50 Preps	27900.00
86792	**BioLit Total RNA extraction Kit From Various Sources (Research)** Kit contents : ● Spin Columns● Buffers	Isolates RNA from various sources - whole blood, animal and plant tissue, cultured cells and bacteria.		50 Preps	27900.00

BioLit® Teaching Kits

66602	**Agarose Gel Electrophoresis Kit (Teaching)** Kit contents : ● DNA Samples ● Loading Dye ● Agarose ● Ethidium Bromide ● Tank Buffer ● DNA Size Marker	The overall goal is to view different plasmid DNAs by their movement through a matrix viz agarose gel.	15	1 Kit	11000.00
64279	**Bacterial Transformation Kit - 1 (Teaching)** Kit contents : ● Host Cells ● Sterile LB Media ● LB agar mix ● Ampicillin ● $CaCl_2$ ● Plasmid DNA	Transformation is central to recombinant DNA technology. This highly self-contained Transformation Kit provides all the necessary materials including starter culture, sterile media, antibiotic & LB agar mix to perform competent cell preparation by $CaCl_2$ treatment and transformation within a day. Overnight incubation will reveal transformed colonies. Competent cells are provided as internal control.	15 30	1 Kit 1 Kit	9700.00 18200.00
28660	**Bacterial Transformation Kit – 2 (with blue white screening) (Teaching)** Kit contents : ● Host Cells ● Sterile LB Media ● LB agar mix ● Ampicillin ● X-gal ● IPTG ● $CaCl_2$ ● Plasmid DNA	Transformation is central to recombinant DNA technology. This highly self-contained Transformation Kit provides all the necessary materials including starter culture, sterile media, antibiotic & LB agar mix to perform competent cell preparation by $CaCl_2$ treatment and transformation within a day. Overnight incubation will reveal transformed colonies. Competent cells are provided as internal control. Blue white screening is also seen with the use of the X-Gal	5 15 30	1 Kit 1 Kit 1 Kit	7300.00 20600.00 36300.00
95043	**DNA Ligation Kit (Teaching)** Kit contents : ● Lambda DNA Digested with Appropriate Restriction Enzymes ● T4 DNA Ligase ● Ligase Buffer ● Agarose ● TAE buffer ● Gel Loading Dye ● Ethidium Bromide ● DNA Molecular Size Marker	This DNA ligation Kit contains reagents for ligation of two Lambda digests using T4 DNA ligase These experiments constitute essential learning experiences for a student interested in cloning and recombinant DNA strategies.	10	1 Kit	5500.00

code	product name	application/description	no. of experiments	unit	price
76941	**FITC Conjugation Kit (Teaching)**	Fluorescein isothiocyanate (FITC)is a fluorescent dye used in labeling antibodies and proteins. FITC is a derivative of fluroscein and reacts with the amines and sulphydryl groups on proteins. This kit provides components for conjugation of fluorescein isothiocyanate with antibodies or other proteins in a quick and easy manner. Kit contents : ● FITC isomer I ● Buffers	5	1 Kit	6700.00
19317	**BioLit Gene Amplification Kit - PCR (Teaching)**	PCR or Polymerase Chain Reaction has become a mainstay in modern recombinant DNA technology. BioLit's PCR kit has been specifically designed and standardized to eliminate the initial pitfalls so that even a beginner can be successful in the first try itself. The kit allows students to vary cycle conditions and find their effects on yield of the amplified product, an E.coli gene Kit contents : ● Template ● Reaction Mix: containing dNTP mix, Buffer, primers and $MgCl_2$ ● TAQ DNA Polymerase ● Enzyme Diluent Buffer ● Dye ● Agarose ● TAE buffer ● Ethidium Bromide ● Marker	10 25	1 Kit 1 Kit	14800.00 21000.00
83462	**BioLit Genomic DNA Extraction (from bacteria) Kit (Teaching)**	This genomic DNA extraction kit uses a rapid and simple technique. Kit contents : ● Bacterial Cell Pellets ● Lysis Buffer ● Saturated Phenol ● Chloroform ● DNA Precipitation Solution ● TE Buffer ● Ethidium Bromide ● Agarose ● Gel Loading Dye ● Marker	15	1 Kit	9100.00
61498	**BioLit Genomic DNA Extraction (from leaves) Kit (Teaching)**	This genomic DNA extraction kit uses a rapid and simple technique ideal for school and college students. Kit contents : ● Bacterial Cell Pellets ● Lysis Buffer ● Saturated Phenol ● Chloroform ● DNA Precipitation Solution ● TE Buffer ● Ethidium Bromide ● Agarose ● Gel Loading Dye ● Marker	15	1 Kit	12700.00
87377	**Onion Root Tip Mitosis Kit (Teaching)**	Demonstrates the different stages of mitosis present in the meristamatic cells located in the onion root tips. Mitosis is the process in which a eukaryotic cell nucleus splits in two, followed by division of the parent cell into two daughter cells.	10 50	1 Kit 1 Kit	3700.00 11000.00
27725	**Plasmid Isolation Kit from Bacterial cells (Teaching)**	This plasmid mini-prep kit facilitates the rapid (~ 2 h) isolation of plasmid DNA from bacterial cells already transformed with a plasmid. The detailed instructions which accompanythe kit will allow typical school and college students to use the kit with no additional assistance. Kit contents : ● Bacterial Cell Pellets ● Alkaline lysis Reagents ● DNA Precipitation Solution ● Wash Solution ● TE Buffer, Gel Loading Dye ● Agarose ● TAE Buffer ● Ethidium Bromide ● Plasmid Control	15 25	1 Kit 1 Kit	8300.00 11200.00

code	product name	application/description	no. of experiments	unit	price
18450	Restriction Digestion Kit (Teaching)	This kit contains materials to teach the use of restriction enzymes HindIII and EcoRI & demonstrate restriction mapping of lambda DNA mapping. These experiments constitute essential learning experiences for a student interested in cloning and genetic mapping.	10 50	1 Kit 1 Kit	4900.00 17000.00
		Kit contents : ● Lambda DNA ● Agarose ● Restriction Enzymes ● Reaction Buffer ● Ethidium Bromide ● Gel Loading Dye ● TAE buffer ● DNA Molecular Size Marker			
52046	Quick Restriction Digestion Kit (Teaching)	Quick Restriction Digestion Kit is designed for faster, more efficient, and easier to perform DNA digestions, while maintaining the high quality and specificity as expected from conventional restriction enzymes. QuickRxn enzymes significantly shorten digestion times, often requiring only 5-15 minutes to complete the test.	20	1 Kit	10400.00
		Kit contents : ● Lambda DNA ● Agarose ● QuickRxn Restriction Enzymes ● Reaction Buffer ● Ethidium Bromide ● Gel Loading Dye ● TAE Buffer ● DNA Molecular Size Marker			
83775	SDS PAGE Kit (Teaching)	This kit teaches the technique of separation of proteins based on their size under denaturing conditions using sodium dodecyl sulphate (SDS) and beta mercaptoethanol (BME). SDS cleaves non ionic interactions while BME cleaves disulphide bonds. SDS also imparts a strong negative charge to the polypeptides thereby their charge to mass ratio is rendered uniform. This allows the polypeptides to be separated on the basis of only their size, which is read by comparison with a standard marker	5 15	1 Kit 1 Kit	6100.00 10400.00
		Kit contents : ● Bis Acrylamide Solution ● Separating and Stacking Buffer Solutions ● APS ● TEMED ● Protein Marker ● Protein Samples ● Stain and Destain Solutions			
31105	Protein Silver Staining Kit (Teaching)	Silver staining is a rapid, easy and sensitive method used for staining proteins, DNA and RNA after polyacrylamide gel electrophoresis. Silver staining is a far more sensitive method of staining than Coomassie Brilliant blue and can detect protein & even polypeptides as low as 0.05 ng in concentration. Reagents sufficient for 15 reactions are enclosed.	5 15	1 Kit 1 Kit	6100.00 12100.00
		Kit contents : ● Fixing solution ● Silver nitrate ● Developing solution ● Stop solution			
46037	Southern Blotting Kit (Teaching)	E. M. Southern invented a technique to transfer separated DNA molecules to a membrane, which is capable of binding the DNA, and immobilizing it for further analysis. The technique named as Southern blotting, after its inventor, involves placing either a nitrocellulose or nylon membrane in contact with the agarose gel and facilitating the transfer process by capillary action. By drawing buffer from the gel to the membrane, the DNA fragments are forced to move with the buffer and deposit on the membrane. (Hydbridization not included)	5	1 Kit	10400.00
		Kit contents : ● NC Membranes ● SSPE Buffer ● Filter Paper Wicks and Strips ● DNA Samples ● Marker ● Agarose ● Ethidium Bromide			
52166	SpinCol Plasmid Miniprep kit from Bacterial cells (Teaching)	A new teaching kit for isolation of plasmid DNA from E.coli cultures based on silica - membrane technology on a mini spin columns. This is a fast and efficient method and yields high quality plasmid DNA.	20 50	1 Kit 1 Kit	6100.00 8700.00
		Kit contents : ● Spin column ● Lysis buffer ● precipitating buffer ● wash buffer ● elution buffer			

Part L
Laboratory Accessories

code	product name	unit	price
Buffer capsules			
87767	Buffer capsules pH 4	10 caps	260
95323	Buffer capsules pH 6	10 caps	260
18533	Buffer capsules pH 7	10 caps	260
74399	Buffer capsules pH 9.2	10 caps	260
Indicator Papers			
74183	Brilliant Yellow Indicator Papers	10 Books	380
11902	Bromothymol Blue Indicator Papers	10 Books	360
84795	Chlorine Test Indicator Papers	10 Books	410
54494	Congo Red Indicator Papers	10 Books	350
72874	Indicator Papers pH 1.0 - 14.0	10 Books	420
12794	Indicator Papers pH 1.0 - 10.0	4 Books	830
47603	Indicator Papers pH 2.0 - 4.5	10 Books	460
27671	Indicator Papers pH 2.0 - 10.5	10 Books	400
84794	Indicator Papers pH 3.5 - 6.0	10 Books	400
82811	Indicator Papers pH 3.8 - 5.3	10 Books	420
98044	Indicator Papers pH 3.8 - 5.4	10 Books	420
36120	Indicator Papers pH 5.0 - 7.5	10 Books	400
61169	Indicator Papers pH 6.5 - 9.0	10 Books	400
94744	Indicator Papers pH 8.0 - 10.5	10 Books	440
76345	Litmus Blue Indicator Papers	10 Books	330
67562	Litmus Red Indicator Papers	10 Books	330
48758	Phenolphthalein Indicator Papers	10 Books	380
27893	Starch Indicator Papers	10 Books	520
75672	Starch Iodide Indicator Papers	10 Books	440
45438	Turmeric Indicator Papers	10 Books	480
Barrier Pen			
93071	**PAP Pen Mini (Hydrophobic Barrier Pen, 1mm, Blue)**	1 nos	7560
59164	**PAP Pen Regular (Hydrophobic Barrier Pen, 3mm, Blue)**	1 nos	12600

Petri Plates

High quality disposable Petri Dishes, suitable for all microbiological applications manufactured as per ISO 24998. Dimensional accuracy of 90 mm ensures compatibility with most automated systems. Extremely high optical clarity for easy viewing on a microscope. Petri Dishes have elevated lids (vented) that allows for superior gas exchange and evaporation. Convenient & easy to stack.

90226	**Petri Plates (90mm)**	100 nos	1420
		500 nos	6990

code	product name	unit	price

Spin Columns

A purification aid set comprising of Silica membrane spin column and a collection tube. Silica membrane in a column form is used for purifying/ extracting nucleic acids (DNA or RNA) of high quality. These columns are compatible with in-house buffers as well as Store-bought ones.

39738	**Mini Silica Spin Column Set A (2ml SC w/ CT each)** Binding capacity — 45µg Supplied with 2ml Collection Tube (CT) Used in purifying nucleic acids of 65bp to 10kb size from PCR amplifications and agarose gel extractions.	50 nos 100 nos	5510 10250
89887	**Mini Silica Spin Column Set B (2ml SC w/ CT each)** Binding capacity — 50µg Supplied with 2ml Collection Tube (CT) Used in purifying nucleic acids of less than 10kb size from lysed cell cultures.	50 nos 100 nos	5510 10250
69516	**Mini Silica Spin Column Set C (2ml SC w/ CT each)** Binding capacity — 55µg Supplied with 2ml Collection Tube (CT) Used in purifying nucleic acids of less than 10kb size from lysed cell cultures.	50 nos 100 nos	5880 11000
37351	**RNA Mini Spin Column Set (2ml SC w/ CT each)** Binding capacity — 40µg Supplied with 2ml Collection Tube (CT) Used in purifying RNA from agarose gel extractions.	50 nos 100 nos	7410 12530

Tubes and Screw Caps

31635	**Microcentrifuge Tube and Screw Cap (0.5ml), Self Standing, Conical, O-Ring, Clear, Non-Hinged, Sterile, Molecular Grade** • Capacity - 0.5 ml • Self Standing Tube With Screw Cap & O-Ring • Sterile - Gamma Irradiated. • Operating range: -196 °C to +121°C. • Autoclavable at 121 °C • Withstands Maximum RCF of 20000X G • No Graduation	500 Nos	7720
60838	**Microcentrifuge Tube and Screw Cap (1.5ml), Self Standing, Conical, O-Ring, Clear, Non-Hinged ,Sterile, Molecular Grade** • Capacity - 1.5 ml • Self Standing Tubes With Screw Cap & O-Ring • Sterile - Gamma Irradiated. • Operating range: -196°C to +121°C • Autoclavable at 121°C • Withstands maximum of 20000X G • No Graduations	500 Nos	8280
41425	**Collection Tubes (CT) for DNA, RNA purification (2ml)**	50 nos 500 nos	470 4320
11876	**Centrifuge Tube and Screw Cap (15ml)**	50 nos 250 nos 500 nos	900 4200 7560
11048	**Centrifuge Tube and Screw Cap (50ml)**	50 nos 250 nos 500 nos	1050 5150 10080

code	product name	unit	price

Crucible Tall Form

These crucibles are made of high quality Porcelain for basic laboratory use to heat chemical compounds. Each product is well tested for its strength, thermal performance, and enhanced chemical resistance with low water absorbance. Widely suitable for burning solid matter, evaporation, concentration or crystallization of solutions.

45866	**Crucible Porcelain, 5 ml, Tall with Lid**	6 Nos	760
14903	**Crucible Porcelain, 10 ml, Tall with Lid**	6 Nos	790
61275	**Crucible Porcelain, 15 ml, Tall with Lid**	4 Nos	690
55162	**Crucible Porcelain, 30 ml, Tall with Lid**	4 Nos	790
94977	**Crucible Porcelain, 40 ml, Tall with Lid**	1 Nos	240



Capacity	Top OD	Height	Wall Thickness	Weight	Lid Dia
5 ml	24 mm	20 mm	1.7 mm	8 gm	26 mm
10 ml	31 mm	26 mm	1.8 mm	12 gm	34 mm
15 ml	35 mm	29 mm	1.8 mm	13 gm	38 mm
30 ml	45 mm	36 mm	1.8 mm	18.5 gm	47 mm
40 ml	50 mm	40 mm	1.8 mm	23 gm	53 mm

Crucible Wide Form

These crucibles are made of high quality Porcelain for basic laboratory use to heat chemical compounds. Each product is well tested for its strength, thermal performance, and enhanced chemical resistance with low water absorbance. Widely suitable for burning solid matter, evaporation, concentration or crystallization of solutions.

16024	**Crucible Porcelain, 2 ml, Wide with Lid**	6 Nos	780
97090	**Crucible Porcelain, 8 ml, Wide with Lid**	6 Nos	870
97897	**Crucible Porcelain, 17 ml, Wide with Lid**	4 Nos	720
71099	**Crucible Porcelain, 30 ml, Wide with Lid**	4 Nos	810
13057	**Crucible Porcelain, 50 ml, Wide with Lid**	1 Nos	270
		4 Nos	950



Capacity	Top OD	Height	Wall Thickness	Weight	Lid Dia
2 ml	18 mm	12 mm	1.3 mm	3 gm	26 mm
8 ml	32 mm	20 mm	1.6 mm	10 gm	34 mm
17 ml	41 mm	26 mm	1.6 mm	18 gm	47 mm
30 ml	50 mm	31 mm	1.8 mm	25 gm	53 mm
50 ml	61 mm	37 mm	1.8 mm	35 gm	69 mm

code	product name	unit	price

Mortar & Pestle

Mortar & Pestle are made from high quality Porcelain for basic laboratory application such as crushing, grinding, etc. They are designed with a spout which enables easy pouring.

44249	**Mortar & Pestle Set, Porcelain, 60 ml**	1 Nos	280
41278	**Mortar & Pestle Set, Porcelain, 100 ml**	1 Nos	320
13968	**Mortar & Pestle Set, Porcelain, 150 ml**	1 Nos	530
51686	**Mortar & Pestle Set, Porcelain, 275 ml**	1 Nos	630



Mortar						Pestle		
Capacity	Top OD	Top ID	Outer Height	Wall Thickness	Weight	Length	Head OD	Weight
60 ml	80 mm	60 mm	44 mm	7.5 mm	200 gm	98 mm	27 mm	50 gm
100 ml	90 mm	75 mm	48 mm	10 mm	280 gm	98 mm	27 mm	50 gm
150 ml	105 mm	85 mm	72 mm	9 mm	490 gm	114 mm	30 mm	90 gm
275 ml	125 mm	100 mm	70 mm	11.5 mm	650 gm	130 mm	36 mm	130-140 gm

Evaporating Dish

These evaporating dishes are made from high quality Porcelain with spout for easy pouring of liquids. These have wide applications in scientific laboratories such as evaporating, drying & crystallizing. They are highly suitable for critical laboratory use.

11443	**Dish Evaporating Porcelain, 20 ml**	6 Nos	840
42418	**Dish Evaporating Porcelain, 40 ml**	6 Nos	880
96704	**Dish Evaporating Porcelain, 100 ml**	6 Nos	1050



Capacity	Top OD	Height	Wall Thickness	Weight
20 ml	50 mm	22 mm	1.3 mm	19 gm
40 ml	60 mm	27 mm	1.3 mm	25 gm
100 ml	75 mm	35 mm	2 mm	51 gm

code	product name	unit	price

Buchner Funnel

Each Buchner funnel is autoclavable, made from high quality Porcelain with fixed perforated plate & a thick stem. It is commonly used in chemistry labs for filtration purposes to eliminate unwanted particles.

53092	**Funnel Buchner, Porcelain, 6 ml**	1 Nos	300
64702	**Funnel Buchner, Porcelain, 30 ml**	1 Nos	320
38635	**Funnel Buchner, Porcelain, 87 ml**	1 Nos	340
59425	**Funnel Buchner, Porcelain, 186 ml**	1 Nos	470



Capacity	Top OD	FP Dia	Height	Perforated Dia	Wall Thickness	Weight
6 ml	35 mm	20 mm	72 mm	23 mm	3 mm	32 gm
30 ml	50 mm	35-40 mm	89 mm	23 mm	3 mm	48 gm
87 ml	60 mm	50-55 mm	110 mm	38 mm	3 mm	110 gm
186 ml	84 mm	70 mm	143 mm	62 mm	3.5 mm	180 gm


Bio-Indenta®
Global searching station for
Unlisted Compounds

Bio-Indenta®

Our Indenting division, "Bio-Indenta" has completed its third successful year of operation. Headquartered in Mumbai, we have been able to master the art of supply chain and global sourcing for our valued scientific researchers. We have contributed to the scientific community by making available complex niche chemicals and biochemicals to our Institutional and Industrial customers. We have efficient logistics systems in place and offer you the transparency and trust you need to make most well informed decisions for your R&D & production needs.

Bio-Indenta exists to connect you to your Principal product provider. We offer you the option of

– importing the products yourself or

– allowing Bio-Indenta to import them and supply them to you.

Transactions can be conducted in Indian Rupees or in any foreign currency.

Through Bio-Indenta, we have established Partnerships and Distributorships with selected pioneer biotech companies for their complete range of biochemicals, diagnostics and fine chemicals to serve scientific and R&D community.

Our partners in this endeavor are,

1. **Carbosynth** (UK)
2. **Biosynth** (Switzerland)
3. **Bachem** (Switzerland)
4. **Elicityl SA** (France)
5. **ACS Materials**(USA)

And Multiple Chinese and European specialty chemical companies who can supply a range of specialty chemicals and intermediates with all required documentations. Now you can choose your most preferred supplier. Take command of your sourcing...

Bio-Indenta®

All your inquiries can be mailed to bio-indenta@srlchem.com for immediate quotes and delivery schedules. Our offices are standing by to assist you in your requirements.

Authorized Distributors for

Carbosynth
UK
(www.carbosynth.com)

Carbosynth

Carbosynth Limited was established to provide the most comprehensive product range in the fields of carbohydrates and nucleosides for the international Life Science R&D communities. Carbosynth aims to support its clients from discovery to launch by offering cost effective solutions to their carbohydrate and nucleoside requirements.

Carbosynth offers over 2500 carbohydrates and nucleosides. This range includes monosaccharides, enzyme substrates, D- and L-sugars, oligosaccharides, detergents and nucleosides. The Carbosynth catalogue offers quantities for R&D, but many are produced in bulk. For example methyl-a-D-glucopyranoside, n-octyl-b-D-glucopyranoside, 3,4,6-tri-acetyl-D-galactal, gulon's-1000's of kilos.

Carbosynth also offers a range of allied fine chemicals many of which are versatile building blocks or have specialist applications ranging from biochemical reagents to natural antioxidants. Included in their product range are coupling agents, such as EDAC, or activators in peptide and oligonucleotide synthesis.

Custom synthesis is one of Carbosynths' fastest growing business sectors and they are committed to providing in-time development and delivery of quality compounds. Carbosynth is able to draw upon a wide range of chemical reactions and processes especially in the fields of carbohydrates and nucleosides.

The Carbosynth product range broadly includes –

- **Carbohydrates and Derivatives of Carbohydrates.**
- **Unnatural Sugars and Carbohydrate Building Blocks.**
- **Nucleosides and Nucleotide Building Blocks.**
- **Oligosaccharides and a range of other Fine Chemicals and Enzyme Substrates.**

Carbosynth Prominent Products

Carbosynth Prod. Code	Product Name
EA06838	O-(2-Acetamido-2-deoxy-D-glucopyranosylidene)amino N-phenyl carbamate
OA08244	N-Acetyl-D-lactosamine
MA00746	N-Acetylneuraminic acid
FA02773	2-Amino-3-chloropyrazine
MA03786	1,6-Anhydro-b-D-galactopyranose
EB06680	5-Bromo-4-chloro-3-indolyl b-D-galactopyranoside
FC06230	6-Chloro-2-cyano-3-nitropyridine
MC02788	1-Chloro-2-deoxy-3,5-di-o-toluoyl-a-Dribofuranose
EC05990	2-Chloro-4-nitrophenyl a-L-fucopyranoside
ND04237	2'-Deoxy-2',2'-difluoro-D-cytidine
MD04977	2-Deoxy-D-ribose
FD06223	4,4'-Dibromo-2,2'-bipyridine
FD04039	R)-2,5-Dihydro-3,6-dimethoxy-2-isopropylpyrazine
FD06273	Diphenylcyclopropenone
FD05800	N-Ethyl-N'-(3-dimethylaminopropyl)carbodiimide HCl
MG05201	D-Galactose
MG05247	L-Glucose
MD05002	1,2:5,6-Di-O-isopropylidene-a-D-allofuranose
MD06795	1,2:5,6-Di-O-isopropylidene-a-D-glucofuranose
EN06363	2-Nitrophenyl b-D-galactopyranoside
DO05161	Octyl b-D-glucopyranoside
MP06689	1,2,3,4,6-Penta-O-acetyl-b-D-galactopyranose
FQ02772	R-(-)-3-Quinuclidinol
FQ06985	S-(+)-3-Quinuclidinol
FT04400	2,2,6,6-Tetramethylpiperidine 1-oxyl, free radical - solid melt
MT04086	2,3,4,6-Tetra-O-benzyl-D-galactopyranose
MT06691	2,3,4,6-Tetra-O-benzyl-D-glucopyranose
MT05853	3,4,6-Tri-O-acetyl-D-galactal
MT00681	3,4,6-Tri-O-acetyl-D-glucal

Please inquire with us for their latest product catalogues & quotations

Authorized Distributors for

Bachem
Switzerland
(www.bachem.com)

Bachem is an independent, technology-based, public biochemicals company providing full service to the pharma and biotech industry. Bachem is specialized in the process development and the manufacturing of peptides and complex organic molecules as active pharmaceutical ingredients (APIs), as well as innovative biochemicals for research purposes. With headquarters in Bubendorf, Switzerland, and affiliates in Europe and the US, Bachem works on a global scale and holds a leading position in the field of peptides.

Bachem are world leaders in Peptides and in large scale GMP Peptide production. The product range includes,

- **Peptides and Bio-chemicals — Bioactive peptides and proteins, Enzyme substrates and Inhibitors, Organic compounds.**
- **Building Blocks — Amino acids, Amino acid derivatives, Linkers, Resins and Reagents.**
- Immunochemicals — Antibodies, Antisera and Assay Kits.

Product Lines
Amyloid ß-Protein Fragments and Related Peptides
Adrenomedullins
Amylins and Fragments
Angiotensins and Related Peptides
Bradykinins, Analogs and Sequences
Caspase-1 Inhibitors
Corticotropin-Releasing Factor (CRF) and Analogs
Endothelins and Related Peptides
ƒ,…-Diaminobutyric Acid
Fibronectin Fragments and Analogs
Ghrelins
HIV Proteases
Hydroxyproline
Isoleucine, allo-Isoleucine
Luteinizing Hormone-Releasing Hormone (LHRH) and Analogs
Melanin-Concentrating Hormones (MCH) and Related Peptides
Neuropeptide Y (NPY), Analogs and Fragments
Orexins
Ornithine
p-Alkoxybenzyl Alcohol Resin (-R)-Linked Fmoc-Amino Acids
pTH, Sequences and Analogs
Somatostatin and Analogs
Threonine, allo-Threonine
Vasopressin and Analogs

Authorized Distributors for

Elicityl SA
France
(www.elicityl-oligotech.com)

Based near Grenoble, France, ELICITYL is a biotech company devoted to complex sugar engineering and production. The company mission is to supply oligosaccharides and polysaccharides for biological applications targeting plant health, pharmacy, nutraceutics and cosmetics.

A few highlighted products offered
ABO, Lewis & P Histo-Blood Group Antigens
Human milk oligosaccharides - HMOs
Neutral, Acidic, Minor HMO"s
Oligosaccharides from other Mammalian species (cow, goat, sheep)
Human & Bacterial Lectins
Functionalized Oligosaccharides
Functionalized Blood group & Lewis antigens
Functionalized Ganglio-Series
Functionalized Globo-Series antigens
Functionalized Fucosylated Oligosaccharides
Functionalized Core structure
Tumor-associated carbohydrate antigens - TACAs
Natural polysaccharides and Oligosaccharides extracted from Biomass
Glycobricks® Building Blocks (Biotin,...)
Neoglycoproteins
Aminoglycoside derivatives

Please inquire with us for their latest product catalogues & quotations

Authorized Distributors for

Biosynth
Switzerland
(www.biosynth.com)

Biosynth is leading swiss specialty reagents manufacturer operating out of its plants in Switzerland and China. It began operations in 1966 as one of the few biochemical manufacturers of the world. It manufactures specialty products for the life-science and biotech industry. It's main core business areas include the production of chromogenic & fluorogenic substrates, organic synthesis products, API's, molecular biology,microbiology, biochemistry, cell science and genetic research products. Bio-Indenta is an exclusive channel partner of Biosynth in India.

A few highlighted products offered
API's
Chromogenic & Fluorogenic Substrates
Luminogenic Substrates
Organic Synthesis Products
Biochemistry, Molecular Biology & Genetic Research Products
Microbiology, Cell Science & In-Vitro Diagnostic Reagents
Imagining Probes
Dyes For Forensic Research
Products for Biomedical Research
Veterinary Science Products
Plant Biotechnology Products
Nutrition Research Products
Cosmetic Research Products

Authorized Distributors for

ACS Materials
USA
(www.acsmaterial.com)

ACS Materials is a high-tech enterprise, based in Pasadena, California (USA), involved in advanced nanomaterials development and production. The Company provides high quality nanomaterial products, service or technical products for research and allied applications.

A few highlighted products offered
Graphene (single-layer graphene, graphene oxide, graphene nanoplates, carboxyl graphene, graphite oxide)
Molecular sieves (SBA-15, MCM-41, MCM-48, SBA-3, ZSM-5)
Carbon series (activated carbon, graphite nanopowder, order mesoporous carbon CMK-3, High-order 3-D mesporous carbon (Ia3d), porous carbon)
Metal nanomaterials (antioxidation dispersible copper nanoparticles, copper colloidals, nickel nanoparticles, silver nanoparticles, silver nanowires)
Battery Materials (Si/C Composite , Silicon Monoxide, Silicon Nanoparticles)

Cell Culture Reagents

Amino Acids

- L-Alanine
- L-Arginine
- L-Asparagine Monohydrate
- L-Cysteine
- L-Cystine
- L-Glutamic Acid
- L-Lysine (free base) Anhydrous
- L-Proline
- L-Serine
- L-Threonine
- L-Tryptophan
- L-Tyrosine
- L-Valine

Antibiotics

- Ampicillin Sodium Salt (AMP-Na)
- Blasticidin S Hydrochloride
- Carbenicillin Disodium Salt
- Chloramphenicol (CFP)
- ß-Cyclodextrin, base (BCD)
- Cycloheximide
- Erythromycin
- G-418 Sulphate (Geneticin)
- Gentamicin Sulphate (GM)
- Hygromycin B (HGR)
- Kanamycin Monosulphate (KM)
- Streptomycin Sulphate
- Tetracycline Hydrochloride (TC)

Buffers

- Glycine
- HEPES Buffer
- HEPES Sodium Salt
- Imidazole
- MES Hydrate Buffer
- PIPES Buffer
- Sodium Acetate Trihydrate

Carbohydrates & Derivatives

- Agar powder
- D-Fructose
- Gelrite Gellan Gum
- Lactose Monohydrate
- D-Maltose
- D-Mannose
- D-Ribose
- D-Sorbitol Powder
- Sucrose
- D-Trehalose Dihydrate

Dyes, Stains & Indicators

- Acridine Orange hemi (Zinc Chloride) Salt
- Alcian Blue
- Alizarin Red S
- DAPI Dihydrochloride (4,6-Diamidino-2 Phenylindole Dihydrochloride)
- New Methylene Blue N Zinc Chloride Double Salt
- Toluidine Blue O
- Trypan Blue
- Thiazolyl Blue Tetrazolium Bromide (MTT)

Enzyme Substrates

- 5-Bromo-4-Chloro-3-Indolyl-b-D Galactopyranoside (X-Gal)
- 5-Bromo-4-Chloro-3-Indolyl-B-D Glucuronide Cyclohexylammonium Salt
- 5-Bromo-4-Chloro-3-Indolyl-b-D Glucopyranoside (X-Glu, X-Glc)
- 5-Bromo-4-Chloro-3-Indolyl Phosphate Disodium Salt (BCIP)
- Nitro Blue Tetrazolium Chloride (Nitro BT) (NBT)

Inorganic Reagents

- Ammonium Chloride
- Ammonium Sulphate
- Calcium Chloride Dihydrate
- Ferrous Sulphate
- Magnesium Chloride Anhydrous
- Magnesium Chloride Hexahydrate
- Magnesium Sulphate Heptahydrate
- Potassium Acetate
- Potassium Chloride
- Potassium Dihydrogen Orthophosphate
- Potassium Phosphate Dibasic Anhydrous
- Potassium Sulphate
- Silver Nitrate
- Sodium Bicarbonate
- Sodium Chloride
- Sodium Phosphate Monobasic Anhydrous
- Zinc Chloride
- Zinc Sulphate Heptahydrate

Organic Reagents

- Adenine Sulphate
- Citric Acid Anhydrous
- Ferric Citrate Trihydrate
- Glutathione Oxidized (GSSG)
- Glutathione Reduced (GSH)
- Lactic Acid
- Polysorbate 20 (Tween 20)
- Putrescine Dihydryochloride
- Urea

Plant Growth Regulators

- Abscisic Acid
- DICAMBA (3,6-Dichloro-2 Methoxybenzoic Acid; 3,6-Dichloro-o Anisic Acid)
- 2,4-Dichlorophenoxyacetic Acid (2,4-D)
- EDTA Disodium Salt Dihydrate
- Gibberellic Acid (Ga3)
- Indole-3-Acetic Acid (IAA)
- Indole-3-Butyric Acid (IBA)
- Kinetin
- Naphthalene-1-Acetic Acid extrapure
- Thidiazuron (TDZ)
- Zeatin mixed Isomers

Proteins & Enzymes

- Albumin Bovine (pH 6.5-7.5) fraction V for tissue culture (Bovine Serum Albumin, BSA), 98%
- r-Albumin Human - Type LL ex. Yeast, 99% (Low Lipid)
- r-Albumin Human - Type Standard ACF ex. Plants, 97% (Animal Component Free)
- r-Albumin Human - Type EG ex. Yeast, 99% (Excipient use)
- r-Albumin Human - Type LE ex. Yeast, 98% (Low Endotoxin)
- r-CRP ex. E.Coli (Recombinant C-reactive protein)
- r-Insulin ex. Human
- r-Trypsin 1× Solution with EDTA in Phosphate Buffer for cell culture (AOF), 150U/ml
- r-Trypsin - Type Standard ACF ex. Porcine, 800U/mg

Vitamins

- L-Ascorbic Acid (Vitamin C)
- D-Biotin (Vitamin H)
- Folic Acid (Vitamin M)
- Inositol (Vitamin B8)
- Thiamine Hydrochloride (Vitamin B1 hydrochloride)
- DL-a-Tocopherol Acetate (Vitamin E Acetate)
- Vitamin B12 (Cyanocobalamin)



Pharmacopoeial Products

(IP/BP/USP/NF)

We take pride in introducing a select range of pharmacopoeial grade products according to IP/BP/USP/NF to meet the critical requirements of the pharmaceutical and biotech industries for drug manufacture and diagnostics use.

These products have stringent limits of purity and impurities, often more critical than specified in the pharmacopoeial standards. This will enable the niche pharmaceutical manufacturers to obtain very high quality raw material ingredients indigenously.

IP, BP, USP & NF Product List

A4P392	6-Aminocaproic Acid IP
A4P393	6-Aminocaproic Acid BP
A4P394	6-Aminocaproic Acid USP
G4P100	Glycine IP
G4P101	Glycine BP
G4P102	Glycine USP
P4P328	L-Proline BP
P4P329	L-Proline USP
P4P330	Propyl Gallate IP
P4P331	Propyl Gallate BP
P4P332	Propyl Gallate NF
S4P306	Salicylic Acid IP
S4P307	Salicylic Acid BP
S4P308	Salicylic Acid USP
S4P304	Sodium Caprylate BP
S4P305	Sodium Caprylate NF
S4P312	Sodium Lauryl Sulphate IP
S4P313	Sodium Lauryl Sulphate BP
S4P309	Sucrose IP
S4P310	Sucrose BP
S4P311	Sucrose USP
T4P289	Tris (Hydroxymethyl) Aminomethane USP (Tromethamine)(Tris Buffer)





Terms & Conditions of Sale

Prices: The Prices offered are net & subject to change without notice. Good shall be delivered subject to the availability of stocks & will be invoiced at prices ruling on the date of dispatch. The prices are unit prices exclusive of GST taxes & any levies.

Taxes: GST taxes & any other levies or charges in lieu of levies as applicable shall be charged extra as and where applicable at the prevailing rates. Customers seeking any exemption must submit the required form along with their orders. All Order must specify your GST registration number.

Delivery: Customers based in cities where SRL Depot's are located, are responsible to collect their ordered/invoiced goods from the respective SRL Depot. For all others customers in India, invoiced goods will be supplied C.P.T. destination by the cheapest mode of transport for invoice value of Rs. 20,000/- net or more. All invoice(s) below Rs. 20,000/- shall be executed on ex-godown basis, all expenses will be charged extra. In case of remote areas, please inform the nearest station of address where goods/ correspondence can be sent by courier. We reserve the right to make part supplies unless otherwise specified by the customers.

For invoices of temperature sensitive products below Rs. 5,000/- net, customers shall be charged additionally for dry ice, packing material, freight/courier, etc.

Special License Products: Methanol/Carbinol and all Pharmacopoeial labelled products of all grades required special permissions and are separately dispatched directly from our factory or respective warehouses which hold valid applicable licenses. Orders for Pharmacopoeial Grade Products must be accompanied with drug license no, date & validity. NMR Solvents require special End-Use certificates from customers and shall be supplied only against receipt of the same.

Shortage / Breakage: All products are packed with utmost care & forwarded at customer's risk. No claim for shortage / breakage, delay or damage will not be entertained if not informed to us in writing within seven days from the date of receipt of goods along with relevant evidences. Many of the products are very labile and likely to be damaged on prolonged storage. Many-a-times the products may get damaged during transit due to extreme condition of temperature, humidity, etc. Improper handling may also cause degradation of quality. Products should be stored in required conditions & protected from temperature, humidity, light, etc. Any complaint on quality of products should be intimated to us within 30 days from receipt of goods beyond which we shall not entertain complaints. No goods should be returned without our prior written consent. It is understood that buyers are aware of storage and handling condition of the products.

Amendment / Cancellation: Goods once sold will not be taken back. No amendment or cancellation of orders will be accepted once the goods have been packed for dispatch.

Warranty: In case of failure or rejection of any of our products for whatsoever reasons, our liability is expressly limited to the extent of its replacement or its cost and we shall not be responsible for any other damages or losses incurred directly or indirectly. All recommendation, information and description supplied with respect to any product are believed to be accurate and reliable, but do not constitute warranties.

Insurance: If expressly desired, goods can be insured at customer's cost at 1% premium on the invoice value.

Force Majeure Clause: Notwithstanding whatever is specified by the customers or mentioned herein, in the event due to reasons beyond our control such as war, riots, strikes, civil commotion, shortage of raw materials, natural calamities, Govt. restrictions, change in fiscal policies and budgetary changes by the Govt., etc. we or our authorised stockist/s will not be liable for any penalization for delay or non-supply of the products.

Products appearing in this catalogue are meant for laboratory use only.

All disputes are subject to Mumbai Jurisdiction.

Sisco Research Laboratories Pvt.Ltd.

The Electrophoresis Leaders of India

We are India's best company in the field of molecular biology since the last 40 years. Our expertise lies in the manufacture of reagents and buffers used in gel electrophoresis for nucleic acid, DNA/RNA and protein separation. Having the most extensive range of protein and nucleic acid electrophoresis applications, we are 100% sure that our products will be very important for you to conduct most critical molecular biology techniques such as electrophoresis of DNA and RNA, PFGE for large fragment analysis, PCR fragment analysis, cloning and recovery of DNA and RNA and various blotting techniques such as Northern & Southern blotting.



SRL is manufacturing few products in bulk scale to meet the following applications in molecular biology for critical applications:

- Gel Electrophoresis
- SDS Page Electrophoresis
- Pulse Field Gel Electrophoresis
- Immunoelectrophoresis
- Counter Current Electrophoresis
- Tissue Culture & Cell Culture Studies
- Forensic & Diagnostic Analysis
- Antigen Antibody Reaction Studies
- Molecular Weight Screening

SRL Special

Agarose Low-6 & Agarose Low-8 have been specially developed for molecular screening applications and improves resolution of small DNA fragments & PCR products

Product code	Agarose	Sulphate %	Gel strength (1.5% gel) gm/cm²	Gelling temperature (1.5% soln)	Melting temperature (1.5% gel)	Electro-endosmosis (EEO)
60645	High EEO for MB	<0.3	>800	36°C (±1.5)	85 - 90°C	0.25 - 0.30
10423	Medium EEO for MB	<0.2	>1500	36°C (±1.5)	85 - 90°C	0.16 - 0.19
91466	Low EEO Regular grade (100bp-25kb)	<0.15	>1200	36°C (±1.5)	85 - 90°C	0.09 - 0.13
36601	Low EEO for MB (100bp-25kb)	<0.15	>1200	36°C (±1.5)	85 - 90°C	0.09 - 0.13
23287	Low EEO Superior grade for MB (100bp-25kb)	<0.15	>1500	36°C (±1.5)	85 - 90°C	0.09 - 0.13
Product code	**Agarose (for Low Base Pairs)**	**Sulphate %**	**Gel strength (3% gel) gm/cm²**	**Gelling temperature (3% soln)**	**Melting temperature (3% gel)**	**Electro-endosmosis (EEO)**
32417	Low Melting for MB (10-1200bp)	<0.2	min 300	35°C (±1.5)	~ 65°C	max 0.1
27813	Metaphor (10-1200bp)	—	min 300	<35°C	<75°C	max 0.1
78711	Low-6 (10-1200bp)	<0.1	min 800	< 35°C	< 75°C	max 0.12
52734	Low-8 (10-1200bp)	<0.11	min 1500	<35.5°C	<80°C	max 0.12

Product code	Acrylamide	Solubility	Assay (Nitrogen)	pH (5% aq. soln)	A1% aq. 290nm
89314	1× Cryst. extrapure	30% aq. Clear & colourless	99.5%	min. 5.0	—
15657	3×Cryst. extrapure AR	60% aq. Clear & colourless	99.9%	min. 5.0	max. 0.5
61346	3× Cryst. for MB	30% aq. Clear & colourless	99.5%	min. 5.0	max. 0.5

Product code	Acrylamide/ Bis-acrylamide Premix	Solubility	pH (5% aq. soln)	Specific conductance (5% aq. soln)
24263	19:1	10% aq. Clear & colourless	5 – 8	min. 10μ mho/cm
10762	29:1	10% aq. Clear & colourless	5 – 8	min. 10μ mho/cm
74379	37.5:1	10% aq. Clear & colourless	5 – 8	min. 10μ mho/cm

Product code	N,N'-Methylene Bis Acrylamide	Solubility	Assay (CG)	pH (2.5% aq. soln)	A1% aq. 290nm
38516	3× Cryst. extrapure AR	2.5% aq.	99.5%	min. 5.0	max. 0.2
67320	3× Cryst. for MB	2.5% aq.	99.5%	min. 5.0	max. 0.1

SISCO RESEARCH LABORATORIES PVT. LTD.

Administrative Office:
608-B, Satellite Gazebo, Near Solitaire Corporate Park, Andheri Ghatkopar Link Road, Chakala, Andheri (East), Mumbai 400 099. India
Tel.: +91-22-4268 5800 • E-mail: info@srlchem.com, marketing@srlchem.com

Website: **www.srlchem.com**